

August 2015

virotype[®] ASFV PCR Kit Handbook



24 (catalog no. 281903)



96 (catalog no. 281905)

For detection of DNA from African swine
fever virus (ASFV)

Licensed in accordance with § 11 (2) of the German Animal Health Act
(Flu-B 670)

REF

281903, 281905



QIAGEN Leipzig GmbH, Deutscher Platz 5b,
04103 Leipzig, Germany

Contents

Kit Contents	3
Intended Use	3
Symbols	4
Storage	5
Safety Information	5
Quality Control	5
Introduction	6
Principle	7
DNA extraction	7
Equipment and Reagents to Be Supplied by User	9
Important Notes	10
General precautions.....	10
Negative Control	10
Positive Control.....	10
Extraction and Amplification Control	11
Protocol: Real-time PCR for detection of African swine fever virus	12
Data Analysis and Interpretation	15
Troubleshooting Guide.....	18
Ordering Information.....	19

Kit Contents

<i>virotype</i> ASFV PCR Kit	(24)	(96)
Catalog no.	281903	281905
Number of reactions	24	96
Master Mix (tube with orange cap), includes enzymes, primers and probes	1 x 500 µl	2 x 980 µl
Positive Control (tube with red cap)	1 x 25 µl	1 x 70 µl
Negative Control (tube with blue cap)	1 x 25 µl	1 x 70 µl
Handbook	1	1

Intended Use

The *virotype* ASFV PCR Kit is intended for the detection of DNA from African swine fever virus (ASFV) in serum, plasma, EDTA-blood, tissue, and swab samples from pigs and wild boar.

The kit is approved by the Friedrich-Loeffler-Institut and licensed in accordance with § 11 (2) of the German Animal Health Act (Flu-B 670) for use in Germany for veterinary diagnostic procedures.

For veterinary use only.

Symbols



<N>

Contains reagents for <N> tests



Legal manufacturer



Lot number



Use by date



Temperature limitations for storage



Handbook



Catalog number



Material number



Protect from light



For pig and wild boar samples

Storage

The components of the *virotype* ASFV PCR Kit should be stored at -30°C to -15°C and are stable until the expiration date stated on the label. Avoid repeated thawing and freezing ($>2\times$), as this may reduce assay sensitivity. Freeze the components in aliquots if they will only be used intermittently.

Safety Information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at www.qiagen.com/safety where you can find, view, and print the MSDS for each QIAGEN kit and kit component.

All sample residues and objects which have come into contact with samples must be decontaminated or disposed of as potentially infective material.

Quality Control

In accordance with QIAGEN's ISO-certified Quality Management System, each lot of *virotype* ASFV PCR Kit is tested against predetermined specifications to ensure consistent product quality.

Introduction

The *virotype* ASFV PCR Kit is a highly sensitive and specific solution for the detection of DNA from African swine fever virus (ASFV) in samples from pigs and wild boar.

African swine fever (ASF) is one of the most important infectious viral diseases of swine of all ages and causes a wide range of clinical signs characterized by a high rate of morbidity and mortality. The disease is notifiable to the World Organization for Animal Health (OIE).

The causative agent is a double-stranded DNA virus belonging to the family *Asfarviridae*, genus *Asfivirus*. ASF virus can be transmitted by vectors (soft ticks of the genus *Ornithodoros*) therefore classified as *Arbovirus* (arthropod-borne virus).

The high sensitivity of the *virotype* ASFV PCR Kit allows early detection of the pathogen in individual as well as in pooled samples.

Principle

Polymerase chain reaction (PCR) is based on the amplification of specific regions of the pathogen genome. In real-time PCR, the amplified product is detected using fluorescent dyes. These are usually linked to oligonucleotide probes that bind specifically to the amplified product. Monitoring the fluorescence intensities during the PCR run (i.e., in real time) allows detection of the accumulating product without the need to re-open the reaction tubes afterward.

The *virotype* ASFV PCR Kit contains all of the necessary reagents for the detection of ASFV DNA, including a positive and negative control.

An internal control excludes the possibility of false-negative results.

The kit uses two specific primer/probe combinations: one for ASFV DNA yielding FAM™ fluorescence and one for a housekeeping gene (β -actin DNA) present within the sample yielding HEX™ fluorescence. A Positive Control serves to verify the functionality of the reaction mix for the amplification of the ASFV DNA target.

DNA extraction

The *virotype* ASFV PCR Kit can be used for the detection of ASFV DNA from serum, plasma, EDTA-blood, tissue, and swab samples from swine. Due to the high sensitivity of the test individual or pooled samples can be tested. Pools of up to 20 individual serum,

plasma, EDTA-blood, or tissue samples can be used, provided that the sample quality is good. It is recommended to test dead wildlife samples on an individual basis.

Prior to real-time PCR, viral DNA must be extracted from the starting material. QIAGEN offers a range of products for DNA extraction from animal samples

- QIAamp® Viral RNA Mini Kit*
- QIAamp *cador*® Pathogen Mini Kit
- DNeasy® Blood & Tissue Kit
- QIAamp DNA Blood Mini Kit

If real-time PCR is not performed immediately after extraction, store the DNA at -30°C to -15°C .

DNA extraction using kits based on spin-column technology can be automated using the QIAcube®.

* Extraction kit allows simultaneous extraction of ASFV DNA and CSFV RNA.

Equipment and Reagents to Be Supplied by User

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, consult the appropriate material safety data sheets (MSDSs), available from the product supplier.

- Pipets
- Nuclease-free aerosol-resistant pipet tips with filters
- Sterile 1.5 ml Eppendorf® tubes
- Nuclease-free (RNase/DNase-free) consumables. Special care should be taken to avoid nuclease contamination of all reagents and consumables used to set up PCR for sensitive identification of viral nucleic acids
- Benchtop centrifuge with rotor for 1.5 ml tubes
- Cooling device or ice
- Rotor-Gene® Q or 96-well plate real-time cycler with appropriate fluorescent channels
- Rotor-Gene Q software version 1.7.94 or higher, or appropriate software for chosen 96-well plate cycler
- Strip Tubes and Caps, 0.1 ml, for use with Rotor-Gene Q (cat. no. 981103 or 981106) or 96-well optical microplate with optical sealing film or cover for chosen 96-well plate real-time cycler

Important Notes

General precautions

The user should always pay attention to the following:

- Use nuclease-free pipet tips with filters.
- Store and extract positive materials (specimens, positive controls, and amplicons) separately from all other reagents, and add them to the reaction mix in a spatially separated facility.
- Thaw all components on ice before starting an assay.
- When thawed, mix the components by inverting and centrifuge briefly.
- Do not use components of the test kit past the expiration date.
- Keep samples and controls on ice or in a cooling block during the setup of reactions.

Negative Control

At least one negative control reaction should be included in each PCR run. This enables assessment of contamination in the reaction.

Positive Control

When performing PCR on unknown samples, it is recommended to perform a positive control reaction in the PCR run, containing a sample that is known to include the targeted viral DNA. A positive

control serves to prove the functionality of the pathogen assay, for example, the correct setup of the reaction mix. Use 5 µl of the Positive Control provided with the *virotype* ASFV PCR Kit to test for successful amplification of the target.

Extraction and Amplification Control

For increased process safety and convenience, an extraction and amplification control assay is included in the form of a second primer/probe set that detects a housekeeping gene present within the sample. This allows both extraction and amplification to be monitored.

Protocol: Real-time PCR for detection of African swine fever virus

Important points before starting

- Please read “Important Notes” on page 10 before starting.
- Include at least one positive control (Positive Control) and one negative control (Negative Control) per PCR run.
- Before beginning the procedure, read through the protocol and ensure that you are familiar with the operation of the chosen real-time PCR cycler.
- Perform the protocol without interruption.

Things to do before starting

- Thaw all reagents on ice and protect from light.
- Maintain reagents on ice during PCR setup.
- Before use, spin the reagents briefly.

Procedure

1. Pipet 20 μ l of the Master Mix into each reaction tube. Then add 5 μ l of the sample DNA (Table 1).

Include positive and negative control reactions.

Positive Control: Use 5 μ l of the positive control (Positive Control) instead of sample DNA.

Negative Control: Use 5 μ l of the negative control (Negative Control) instead of sample DNA.

Table 1. Preparation of reaction mix

Component	Volume
Master Mix	20 µl
Sample	5 µl
Total volume	25 µl

2. Close the reaction tubes with the corresponding caps.
3. Set the filters for the reporter dyes in the software of your thermal cycler according to Table 2. Select the green and yellow channels on the Rotor-Gene Q.

Important: Set a fixed gain of +4 in the green and +1 in the yellow channels to ensure optimal fluorescence gains for the pathogen and the Internal Control assays when using the Rotor-Gene Q.

Table 2. Filter settings for the reporter

Pathogen/Internal Control	Reporter	Rotor-Gene Q
ASFV	FAM	green
Internal Control	HEX/JOE™*	yellow
Passive reference†	ROX™	

* Use the option appropriate for your thermal cycler.

† Internal reference for use with Applied Biosystems® ABI PRISM® Sequence Detection Systems.

4. Run the real-time PCR protocol according to Table 3 if running only the *virotype* ASFV PCR Kit.

Table 3. Real-time PCR Protocol for ASFV

Temperature	Time	Number of cycles
95°C	15 min	1
95°C	15 s	40
60°C*	60 s	

* Fluorescence data collection. Approximate run time (Rotor-Gene Q): 96min.

5. Run the real-time RT-PCR protocol according to Table 4 if running the *virotype* CSFV assay simultaneously.

Table 4. Real-time RT-PCR Protocol for parallel amplification of ASFV and CSFV[†]

Temperature	Time	Number of cycles
45°C	10 min	1
95°C	10 min	1
95°C	15 s	40
57°C‡	30 s	
72°C	35 s	

[†] Valid for *virotype* CSFV RT-PCR Kit only.

[‡] Fluorescence data collection. Approximate run time (Rotor-Gene Q): 118 min

Data Analysis and Interpretation

Interpretation of results

For the assay to be valid, the Positive Control must give a signal in both the FAM and HEX channels with a $C_T^* < 35$. The Negative Control must give no signal.

The following results are possible if working with unknown samples. The possible sample results are also summarized in Table 5 on page 17.

The sample is positive for ASFV, and the assay is valid, if the following criteria are met:

- The sample yields a signal in both the FAM and HEX[†] channels
- The Positive Control yields a signal in both the FAM and HEX channels
- The Negative Control does not yield a signal in the FAM and HEX channels

Note that very high concentrations of ASFV DNA in the sample may lead to a reduced HEX signal or no HEX signal due to competition with the internal control.

* Threshold cycle (C_T) — cycle at which the amplification plot crosses the threshold, i.e., there is the first clearly detectable increase in fluorescence.

† Green and yellow on the Rotor-Gene Q.

The sample is negative for ASFV, and the assay is valid, if the following criteria are met:

- The sample yields a signal in the HEX channel but not in the FAM channel
- The Positive Control yields a signal in both the FAM and HEX channels
- The Negative Control does not yield a signal in the FAM and HEX channels

A positive HEX signal means that extraction and amplification were successful as the housekeeping gene within the sample is amplified. However, if the C_T value of the internal control is >35 , pooled or individual samples could be partially inhibited. In such cases it is recommended that the respective individual samples are diluted (e.g., diluted 1:5) in nuclease free water and retested.

The sample results are inconclusive, and the assay is invalid, if the following occurs:

- The sample yields no signal in the FAM and HEX channel

If no signal is detected in both the FAM (pathogen) and the HEX (Internal Control) channel, the result is inconclusive. The absence of a signal for the housekeeping gene indicates PCR inhibition and/or other malfunctions.

To check for inhibition, we recommend 1:5 dilution of the sample DNA in nuclease free water, to repeat the DNA extraction, or repeat the whole test procedure starting with new sample material.

Check that there is a fluorescence signal in the FAM channel for the positive control reaction (Positive Control). Absence of a signal for the Positive Control indicates an error, which could be due to incorrect setup of the reaction mix or incorrect cycling conditions.

Table 5. Results interpretation table*

Sample result	Reporter	
	FAM (pathogen)	HEX (IC)
ASFV positive	X	X
ASFV positive (strong positive)	X	
ASFV negative		X
Inconclusive result		

* Interpretation of sample results can be determined provided positive and negative control reactions are performed. The positive control must yield a signal in both the FAM and HEX channels. The negative control must yield no signal in the FAM and HEX channels. For a complete explanation of possible sample results, please refer to “Data Analysis and Interpretation” on page 15.

Troubleshooting Guide

The scientists in QIAGEN Technical Services are always happy to answer any questions you may have about either the information and/or protocols in this handbook or sample and assay technologies (for contact information, visit www.qiagen.com).

Ordering Information

Product	Contents	Cat. no.
<i>virotype</i> ASFV PCR Kit (24)	For 24 reactions: Master Mix, Positive Control, Negative Control	281903
<i>virotype</i> ASFV PCR Kit (96)	For 96 reactions: Master Mix, Positive Control, Negative Control	281905
Related products		
<i>virotype</i> CSFV RT-PCR Kit (96)*	For 96 reactions: Master Mix, Positive Control, Negative Control	281805
<i>virotype</i> PEDV/TGEV RT-PCR Kit (96)*	For 96 reactions: Master Mix, Positive Control, Negative Control	283605
<i>virotype</i> BTV RT-PCR Kit (96)*	For 96 reactions: PCR Mix, Enzyme Mix, Positive Control, Negative Control	280435
<i>virotype</i> BTV pan/8 RT-PCR Kit (96)*	For 96 reactions: Master Mix, Positive Control, Negative Control	280445
<i>virotype</i> BVDV RT-PCR Kit (96)*	For 96 reactions: PCR Mix, Enzyme Mix, Positive Control, Negative Control	280375
<i>virotype</i> SBV RT-PCR Kit (96)*	For 96 reactions: Master Mix, Positive Control, Negative Control	281605
<i>virotype</i> PRRSV RT-PCR Kit (96)*	For 96 reactions: Master Mix, Positive Control, Negative Control	282305
<i>virotype</i> Influenza A RT-PCR Kit (96)*	For 96 reactions: Master Mix, Positive Control, Negative Control	282605

* Other kit sizes are available; see www.qiagen.com.

Product	Contents	Cat. no.
<i>bactotype</i> [®] Mycoplasma Mg/Ms PCR Kit (96)*	For 96 reactions: Master Mix, Positive Control, Negative Control	288105
<i>bactotype</i> MAP PCR Kit (96)*	For 96 reactions: Master Mix, Internal Control DNA, Positive Control, Negative Control	285905
QIAamp <i>cador</i> Pathogen Mini Kit (50)*	For 50 preps: 50 QIAamp Mini Spin Columns, Carrier RNA, Proteinase K, Collection Tubes (2 ml), RNase-free Buffers	54104
QIAamp Viral RNA Mini Kit (50)*	For 50 RNA preps: 50 QIAamp Mini Spin Columns, Carrier RNA, Collection Tubes (2 ml), RNase-free buffers	52904
DNeasy Blood & Tissue Kit (50)*	For 50 preps: 50 DNeasy Mini Spin Columns, Proteinase K, Buffers, Collection Tubes (2 ml)	69504
QIAamp DNA Blood Mini Kit (50)*	For 50 DNA preps: 50 QIAamp Mini Spin Columns, QIAGEN Protease, Reagents, Buffers, Collection Tubes (2 ml)	51104
Rotor-Gene Q 5plex Platform	Real-time PCR cycler with 5 channels (green, yellow, orange, red, crimson), laptop computer, software, accessories, 1-year warranty on parts and labor	9001570

* Other kit sizes are available; see www.qiagen.com.

QIAGEN offers a range of ELISA kits and real-time PCR and real-time RT-PCR kits for the detection of animal pathogens. Visit **www.qiagen.com/Animal-and-Veterinary-Testing** for more information about *bactotype*, *cador*, *cattletype*[®], *flocktype*[®], *pigtype*[®] and *virotype* products.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at **www.qiagen.com** or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN[®], QIAamp[®], QIAcube[®], *bactotype*[®], *cador*[®], *cattletype*[®], DNeasy[®], *flocktype*[®], *pigtype*[®], Rotor-Gene[®], Sample to Insight[®], *virotype*[®] (QIAGEN Group); Applied Biosystems[®], ABI PRISM[®] (Applied Biosystems); FAM[™], HEX[™], JOE[™], ROX[™] (Life Technologies Corporation); Eppendorf[®] (Eppendorf-Netheler-Hinz GmbH). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Limited License Agreement for *virotype* ASFV PCR Kit

Use of this product signifies the agreement of any purchaser or user of the product to the following terms:

1. The product may be used solely in accordance with the protocols provided with the product and this handbook and for use with components contained in the kit only. QIAGEN grants no license under any of its intellectual property to use or incorporate the enclosed components of this kit with any components not included within this kit except as described in the protocols provided with the product, this handbook, and additional protocols available at www.qiagen.com. Some of these additional protocols have been provided by QIAGEN users for QIAGEN users. These protocols have not been thoroughly tested or optimized by QIAGEN. QIAGEN neither guarantees them nor warrants that they do not infringe the rights of third-parties.
2. Other than expressly stated licenses, QIAGEN makes no warranty that this kit and/or its use(s) do not infringe the rights of third-parties.
3. This kit and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. QIAGEN specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the kit agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. QIAGEN may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the kit and/or its components.

For updated license terms, see www.qiagen.com.

The purchase of this product allows the purchaser to use it for amplification and identification of nucleic acid sequences for laboratory use only. No general patent or other license of any kind other than this specific right of use from purchase is granted hereby.

1096310 HB-1881-002 © 2015 QIAGEN, all rights reserved.

Notes

www.qiagen.com