

# *ipsogen*<sup>®</sup> BCR-ABL1 mbcR Kit Handbook



Version 1

**IVD**

Quantitative in vitro diagnostics

For use with Rotor-Gene<sup>®</sup> Q, ABI PRISM<sup>®</sup>, LightCycler<sup>®</sup>, and SmartCycler<sup>®</sup> instruments



**REF**

670023



QIAGEN GmbH, QIAGEN Strasse 1, 40724 Hilden, GERMANY

**R2**

**MAT**

1072506EN



## **QIAGEN Sample and Assay Technologies**

QIAGEN is the leading provider of innovative sample and assay technologies, enabling the isolation and detection of contents of any biological sample. Our advanced, high-quality products and services ensure success from sample to result.

### **QIAGEN sets standards in:**

- Purification of DNA, RNA, and proteins
- Nucleic acid and protein assays
- microRNA research and RNAi
- Automation of sample and assay technologies

Our mission is to enable you to achieve outstanding success and breakthroughs. For more information, visit [www.qiagen.com](http://www.qiagen.com).

# Contents

<b>Intended Use</b>	<b>5</b>
<b>Summary and Explanation</b>	<b>5</b>
<b>Principle of the Procedure</b>	<b>6</b>
<b>Materials Provided</b>	<b>8</b>
Kit contents	8
<b>Materials Required but Not Provided</b>	<b>9</b>
<b>Warnings and Precautions</b>	<b>10</b>
General precautions	10
<b>Reagent Storage and Handling</b>	<b>11</b>
<b>Procedure</b>	<b>12</b>
Sample RNA preparation	12
Protocols	
■ Recommended standardized EAC reverse transcription	12
■ qPCR on Rotor-Gene Q MDx 5plex HRM or Rotor-Gene Q 5plex HRM instruments with 72-tube rotor	15
■ qPCR on ABI PRISM 7000, 7700, and 7900HT SDS, and LightCycler 480 Instrument	19
■ qPCR on LightCycler 1.2 and 2.0 Instruments	24
■ qPCR on the SmartCycler instrument	28
<b>Interpretation of Results</b>	<b>31</b>
Data analysis principle	31
Results	32
Troubleshooting guide	34
<b>Quality Control</b>	<b>37</b>
<b>Limitations</b>	<b>38</b>
<b>Performance Characteristics</b>	<b>38</b>
Nonclinical studies	38
Clinical studies	41
<b>References</b>	<b>43</b>
<b>Symbols</b>	<b>44</b>
<b>Contact Information</b>	<b>45</b>



## Intended Use

The *ipsogen* BCR-ABL1 mbc Kit is intended for the quantification of BCR-ABL p190 transcripts in bone marrow or peripheral blood samples of Ph-positive acute lymphoblastic leukemia (ALL) patients previously diagnosed with a BCR-ABL mbc fusion gene (FG) event. The results obtained are intended to monitor efficacy of treatment in patients undergoing therapy and for minimal residual disease (MRD) follow-up to monitor disease relapse.

## Summary and Explanation


The Philadelphia (Ph) chromosome is the most frequent karyotypic aberration in adults with ALL. It occurs in 20–30% of adult patients with ALL overall, with the incidence rising to more than 50% in patients aged 50 years or older.

In this translocation, the 3' segment of ABL proto-oncogene on chromosome 9 is juxtaposed with the 5' segment of the BCR gene on chromosome 22. The BCR-ABL FG is the product of the Ph chromosome and is a constitutively active tyrosine kinase protein.

Breaks in the ABL gene typically occur in the first intron. Breaks in the BCR gene generally occur in one of the following 3 regions: a 5.8 kb region spanning exons 12–16, called the major breakpoint cluster region (Mbc), a 55 kb sequence of the first intron, called the minor breakpoint cluster region (mbc), and the micro breakpoint cluster region ( $\mu$ -bc).

Breakpoints occurring in mbc join exon 1 (e1) with the second exon of the ABL gene (a2) resulting in a smaller fusion transcript, e1a2, that encodes a 190 kDa (p190) chimeric protein (Figure 1). The p190 BCR-ABL protein is only observed in Ph+ ALL while the p210 BCR-ABL protein is common to 20–40% of patients with Ph+ ALL and nearly all patients with Ph+ chronic myelogenous leukemia (CML).

All forms of BCR-ABL fusion proteins display an increased and deregulated tyrosine kinase activity, and the p190 form has been shown to have more transforming potential than p210. Moreover, this chimeric protein seems to deregulate the normal cytokine-dependent signal transduction pathways, leading to the inhibition of apoptosis or growth factor independent growth.



**Figure 1. Schematic diagram of the BCR-ABL mbcR FG transcript covered by the qPCR primers and probe set: ENF402–ENP541–ENR561.** The number under the primers and probe refers to their nucleotide position in the normal gene transcript.

Therapy of Ph+ ALL patients has been optimized by the introduction of tyrosine kinase inhibitors, which significantly improved the survival of these patients (for a review, see reference 1). For these patients, monitoring of the MRD is required. The current methodology to measure the MRD level involves using real-time quantitative polymerase chain reaction (qPCR), whereby the BCR-ABL transcript numbers are related to transcript numbers of a control gene. The *ipsogen* BCR-ABL1 mbcR Kit is based on this technique.

## Principle of the Procedure


qPCR permits the accurate quantification of PCR products during the exponential phase of the PCR amplification process. Quantitative PCR data can be rapidly obtained, without post-PCR processing, by real-time detection of fluorescent signals during and/or subsequent to PCR cycling, thereby drastically reducing the risk of PCR product contamination. At present, 3 main types of qPCR techniques are available: qPCR analysis using SYBR® Green I Dye, qPCR analysis using hydrolysis probes, and qPCR analysis using hybridization probes.

This assay exploits the qPCR double-dye oligonucleotide hydrolysis principle. During PCR, forward and reverse primers hybridize to a specific sequence. A double-dye oligonucleotide is contained in the same mix. This probe, which consists of an oligonucleotide labeled with a 5' reporter dye and a downstream, 3'quencher dye, hybridizes to a target sequence within the PCR product. qPCR analysis with hydrolysis probes exploits the 5'→3' exonuclease activity of the *Thermus aquaticus* (*Taq*) DNA polymerase. When the probe is intact, the proximity of the reporter dye to the quencher dye results in suppression of the reporter fluorescence primarily by Förster-type energy transfer.

During PCR, if the target of interest is present, the probe specifically anneals between the forward and reverse primer sites. The 5'→3' exonuclease activity of the DNA polymerase cleaves the probe between the reporter and the quencher

only if the probe hybridizes to the target. The probe fragments are then displaced from the target, and polymerization of the strand continues. The 3' end of the probe is blocked to prevent extension of the probe during PCR (Figure 2). This process occurs in every cycle and does not interfere with the exponential accumulation of product.

The increase in fluorescence signal is detected only if the target sequence is complementary to the probe and hence amplified during PCR. Because of these requirements, nonspecific amplification is not detected. Thus, the increase in fluorescence is directly proportional to the target amplification during PCR.



**Figure 2. Reaction principle.** Total RNA is reverse-transcribed, and the generated cDNA is amplified by PCR using a pair of specific primers and a specific internal double-dye probe (FAM™–TAMRA™). The probe binds to the amplicon during each annealing step of the PCR. When the *Taq* DNA polymerase extends from the primer bound to the amplicon, it displaces the 5' end of the probe, which is then degraded by the 5'→3' exonuclease activity of the *Taq* DNA polymerase. Cleavage continues until the remaining probe melts off the amplicon. This process releases the fluorophore and quencher into solution, spatially separating them and leading to an increase in fluorescence from the FAM and a decrease in fluorescence from the TAMRA.

# Materials Provided

## Kit contents

<b>ipsogen BCR-ABL1 mbc Kit</b>		<b>(24)</b>
<b>Catalog no.</b>		<b>670023</b>
<b>Number of reactions</b>		<b>24</b>
ABL Control Gene Standard Dilution (10 <sup>3</sup> copies/5 µl)	C1-ABL	50 µl
ABL Control Gene Standard Dilution (10 <sup>4</sup> copies/5 µl)	C2-ABL	50 µl
ABL Control Gene Standard Dilution (10 <sup>5</sup> copies/5 µl)	C3-ABL	50 µl
BCR-ABL mbc Fusion Gene Standard Dilution (10 <sup>1</sup> copies/5 µl)	F1-BCR-ABL e1a2 mbc	50 µl
BCR-ABL mbc Fusion Gene Standard Dilution (10 <sup>2</sup> copies/5 µl)	F2-BCR-ABL e1a2 mbc	50 µl
BCR-ABL mbc Fusion Gene Standard Dilution (10 <sup>3</sup> copies/5 µl)	F3-BCR-ABL e1a2 mbc	50 µl
BCR-ABL mbc Fusion Gene Standard Dilution (10 <sup>5</sup> copies/5 µl)	F4-BCR-ABL e1a2 mbc	50 µl
BCR-ABL mbc Fusion Gene Standard Dilution (10 <sup>6</sup> copies/5 µl)	F5-BCR-ABL e1a2 mbc	50 µl
Primers and Probe Mix ABL*	PPC-ABL 25x	90 µl
Primers and Probe Mix BCR-ABL mbc Fusion Gene†	PPF-mbc 25x	110 µl
ipsogen <i>BCR-ABL1 mbc Kit Handbook</i> (English)		1

\* Mix of specific reverse and forward primers for the ABL control gene plus a specific FAM-TAMRA probe.

† Mix of specific reverse and forward primers for the BCR-ABL mbc fusion gene plus a specific FAM-TAMRA probe.

**Note:** Briefly centrifuge the standard dilutions and the primers and probe mixes before use.

## Materials Required but Not Provided

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, consult the appropriate safety data sheets (SDSs), available from the product supplier.

### Reagents

- Nuclease-free PCR grade water
- Reagents for reverse transcription: The validated reagent is Superscript® II (or Superscript) Reverse Transcriptase, includes 5x first-strand buffer, 100 mM DTT (Life Technologies, cat. no. 18064-022)
- RNase inhibitor: The validated reagent is RNaseOUT™ (Life Technologies, cat. no. 10777-019)
- Set of dNTPs, PCR grade
- Random hexamer
- MgCl<sub>2</sub>
- Buffer and Taq DNA polymerase: The validated reagents are TaqMan® Universal PCR Master Mix (Master Mix PCR 2x) (Life Technologies, cat. no. 4304437) and LightCycler TaqMan Master (Master Mix PCR 5x) (Roche, cat. no. 04535286001)

### Consumables

- Nuclease-free aerosol-resistant sterile PCR pipet tips with hydrophobic filters
- 0.5 ml or 0.2 ml RNase- and DNase-free PCR tubes
- Ice

### Equipment

- Microliter pipet\* dedicated for PCR (1–10 µl; 10–100 µl; 100–1000 µl)
- Benchtop centrifuge\* with rotor for 0.2 ml/0.5 ml reaction tubes (capable of attaining 10,000 rpm)
- Real-time PCR instrument:\* Rotor-Gene Q MDx 5plex HRM or other Rotor-Gene instrument; LightCycler 1.2, 2.0, or 480; ABI PRISM 7000, 7700, or 7900HT SDS; or SmartCycler instrument; and associated specific material
- Thermal cycler\* or water bath\* (reverse transcription step)

\* Ensure that instruments have been checked and calibrated according to the manufacturer's recommendations.

## Complementary reagents

- *ipsogen* BCR-ABL1 mbcR Controls Kit (cat. no. 670091), consisting of cell lines with negative, high, and low positive expression of the BCR-ABL mbcR fusion gene for the qualitative validation of the RNA extraction and the reverse transcription

## Warnings and Precautions

For in vitro diagnostic use

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at [www.qiagen.com/safety](http://www.qiagen.com/safety) where you can find, view, and print the SDS for each QIAGEN kit and kit component.

Discard sample and assay waste according to your local safety regulations.

## General precautions

qPCR tests require good laboratory practices, including equipment maintenance, that are dedicated to molecular biology and compliant with applicable regulations and relevant standards.

This kit is intended for in vitro diagnostic use. Reagents and instructions supplied in this kit have been validated for optimal performance. Further dilution of the reagents or alteration of incubation times and temperatures may result in erroneous or discordant data. PPC and PPF reagents may be altered if exposed to light. All reagents are formulated specifically for use with this test. For optimal performance of the test, no substitutions should be made.

Determining transcript levels using qPCR requires both the reverse transcription of the mRNA and the amplification of the generated cDNA by PCR. Therefore, the entire assay procedure must be performed under RNase-/DNase-free conditions.

Use extreme caution to prevent:

- RNase/DNase contamination, which might cause degradation of the template mRNA and the generated cDNA
- mRNA or PCR carryover contamination resulting in false positive signal

We therefore recommend the following.

- Use nuclease-free labware (e.g., pipets, pipet tips, reaction vials) and wear gloves when performing the assay.
- Use fresh aerosol-resistant pipet tips for all pipetting steps to avoid cross-contamination of the samples and reagents.

- Prepare pre-PCR master mix with dedicated material (pipets, tips, etc.) in a dedicated area where no DNA matrixes (cDNA, DNA, plasmid) are introduced. Add template in a separate zone (preferably in a separate room) with specific material (pipets, tips, etc.).
- Handle the standard dilutions (C1–3 and F1–5) in a separate room.

## Reagent Storage and Handling

The kits are shipped on dry ice and must be stored at  $-30^{\circ}\text{C}$  to  $-15^{\circ}\text{C}$  upon receipt.

- Minimize exposure to light of the primers and probe mixes (PPC and PPF tubes).
- Gently mix and centrifuge the tubes before opening.
- Store all kit components in original containers.

These storage conditions apply to both opened and unopened components. Components stored under conditions other than those stated on the labels may not perform properly and may adversely affect the assay results.

Expiration dates for each reagent are indicated on the individual component labels. Under correct storage conditions, the product will maintain performance until the expiration date printed on the label.

There are no obvious signs to indicate instability of this product. However, positive and negative controls should be run simultaneously with unknown specimens.

## Procedure

### Sample RNA preparation

RNA preparation from patient samples (blood or bone marrow) must have been performed with a validated procedure. The quality of the assay is largely dependent on the quality of input RNA. We therefore recommend qualifying the purified RNA by agarose\* gel electrophoresis or by using Agilent® Bioanalyzer® prior to analysis.

### Protocol: Recommended standardized EAC reverse transcription

#### Things to do before starting

- Prepare dNTPs, 10 mM each. Store at –20°C in aliquots.

#### Procedure

1. Thaw all necessary components and place them on ice.
2. Incubate 1  $\mu\text{g}$  of RNA (1–4  $\mu\text{l}$ ) for 10 minutes at 70°C and immediately cool on ice for 5 minutes.
3. Centrifuge briefly (approximately 10 seconds, 10,000 rpm, to collect the liquid in the bottom of the tube). Then keep on ice.
4. Prepare the following RT mix according to the number of samples being processed (Table 1).

\* When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles.

**Table 1. Preparation of RT mix**

<b>Component</b>	<b>Volume per sample (<math>\mu</math>l)</b>	<b>Final concentration</b>
First-Strand Buffer (supplied with Superscript II Reverse Transcriptase), 5x	4.0	1x
MgCl <sub>2</sub> (50 mM)	2.0	5 mM
dNTPs (10 mM each, to be prepared previously and stored at -20°C in aliquots)	2.0	1 mM
DTT (100 mM, supplied with Superscript II Reverse Transcriptase)	2.0	10 mM
RNase inhibitor (40 U/ $\mu$ l)	0.5	1 U/ $\mu$ l
Random hexamer (100 $\mu$ M)	5.0	25 $\mu$ M
Superscript II or Superscript Reverse Transcriptase (200 U/ $\mu$ l)	0.5	5 U/ $\mu$ l
Heated RNA sample (to be added in step 5)	1.0–4.0	50 ng/ $\mu$ l
Nuclease-free PCR grade water (to be added in step 5)	0.0–3.0	–
Final volume	20.0	–

- 5. Pipet 16  $\mu$ l of RT mix into each PCR tube. Then add 1–4  $\mu$ l (1  $\mu$ g) RNA (from step 3), and adjust the volume to 20  $\mu$ l with nuclease-free PCR grade water (see Table 2).**

**Table 2. Preparation of reverse transcription reaction**

<b>Component</b>	<b>Volume (<math>\mu</math>l)</b>
RT mix	16
Heated sample RNA (1 $\mu$ g)	1–4
Nuclease-free PCR grade water	0–3
Final volume	20

- 6. Mix well and centrifuge briefly (approximately 10 seconds, 10,000 rpm, to collect the liquid in the bottom of the tube).**
- 7. Incubate at 20°C for 10 minutes.**
- 8. Incubate at 42°C on a thermal cycler for 45 minutes, then immediately at 99°C for 3 minutes.**
- 9. Cool on ice (to stop the reaction) for 5 minutes.**
- 10. Briefly centrifuge (approximately 10 seconds, 10,000 rpm, to collect the liquid in the bottom of the tube). Then keep on ice.**
- 11. Dilute the final cDNA with 30  $\mu$ l of nuclease-free PCR grade water so that the final volume is 50  $\mu$ l.**
- 12. Carry out PCR according to the following protocols, according to your qPCR instrument.**

## Protocol: qPCR on Rotor-Gene Q MDx 5plex HRM or Rotor-Gene Q 5plex HRM instruments with 72-tube rotor


Using this instrument, we recommend performing all measurements in duplicate, as indicated in Table 3.

**Table 3. Number of reactions for Rotor-Gene Q instruments with 72-tube rotor**

<b>Samples</b>	<b>Reactions</b>
<b>With the ABL primers and probe mix (PPC-ABL)</b>	
n cDNA samples	n x 2 reactions
ABL standard	2 x 3 reactions (3 dilutions, each one tested in duplicate)
Water control	2 reactions
<b>With the BCR-ABL mbcr primers and probe mix (PPF-mbcr)</b>	
n cDNA samples	n x 2 reactions
mbcr standard	2 x 5 reactions (5 dilutions, each one tested in duplicate)
Water control	2 reactions

### Sample processing on Rotor-Gene Q instruments with 72-tube rotor

We recommend testing at least 8 cDNA samples in the same experiment to optimize the use of the standards and the primers and probe mixes.



**Figure 3. Suggested rotor setup for each experiment with the ipsogen BCR-ABL1 mbc Kit.** F1–5: BCR-ABL mbc standards; C1–3: ABL standards; S: cDNA sample; H<sub>2</sub>O: water control.

**Note:** Take care to always place a sample to be tested in position 1 of the rotor. Otherwise, during the calibration step, the instrument will not perform calibration, and incorrect fluorescence data will be acquired.

Fill all other positions with empty tubes.

### qPCR on Rotor-Gene Q instruments with 72-tube rotor

**Note:** Perform all steps on ice.

#### Procedure

1. Thaw all necessary components and place them on ice.
2. Prepare the following qPCR mix according to the number of samples being processed.

All concentrations are for the final volume of the reaction.

Table 4 describes the pipetting scheme for the preparation of one reagent mix, calculated to achieve a final reaction volume of 25  $\mu$ l. A pre-mix can be prepared, according to the number of reactions, using the same primer and probe mix (either PPC-ABL or PPF-mbcr). Extra volumes are included to compensate for pipetting error.

**Table 4. Preparation of qPCR mix**

<b>Component</b>	<b>1 reaction (<math>\mu</math>l)</b>	<b>ABL: 24+1 reactions (<math>\mu</math>l)</b>	<b>BCR-ABL mbcr: 28+1 reactions (<math>\mu</math>l)</b>	<b>Final concentration</b>
TaqMan Universal PCR Master Mix, 2x	12.5	312.5	362.5	1x
Primers and probe mix, 25x	1	25	29	1x
Nuclease- free PCR grade water	6.5	162.5	188.5	–
Sample (to be added at step 4)	5	5 each	5 each	–
Total volume	25	25 each	25 each	–

3. Dispense 20  $\mu$ l of the qPCR pre-mix per tube.
4. Add 5  $\mu$ l of the RT product (cDNA, 100 ng RNA equivalent) obtained in the reverse transcription (see “Protocol: Recommended standardized EAC reverse transcription”, page 12) in the corresponding tube (total volume 25  $\mu$ l).
5. Mix gently, by pipetting up and down.
6. Place the tubes in the thermal cycler according to the manufacturer recommendations.
7. Program the Rotor-Gene Q instrument with the thermal cycling program as indicated in Table 5.

**Table 5. Temperature profile**

<b>Mode of analysis</b>	Quantitation
<b>Hold</b>	Temperature: 50 deg Time: 2 mins
<b>Hold 2</b>	Temperature: 95 deg Time: 10 mins
<b>Cycling</b>	50 times 95 deg for 15 secs 60 deg for 1 min with acquisition of FAM fluorescence in channel Green: Single

- 8. For Rotor-Gene Q instruments, select "Slope Correct" for the analysis. We recommend setting the threshold at 0.03. Start the thermal cycling program, as indicated in Table 5.**

## Protocol: qPCR on ABI PRISM 7000, 7700, and 7900HT SDS, and LightCycler 480 Instrument


Using 96-well-plate qPCR equipment, we recommend performing all measurements in duplicate, as indicated in Table 6.

**Table 6. Number of reactions using 96-well-plate qPCR equipment**

<b>Samples</b>	<b>Reactions</b>
<b>With the ABL primers and probe mix (PPC-ABL)</b>	
n cDNA samples	n x 2 reactions
ABL standard	2 x 3 reactions (3 dilutions, each one tested in duplicate)
Water control	2 reactions
<b>With the BCR-ABL mbcr primers and probe mix (PPF-mbcr)</b>	
n cDNA samples	n x 2 reactions
mbcr standard	2 x 5 reactions (5 dilutions, each one tested in duplicate)
Water control	2 reactions

### Sample processing on ABI PRISM 7000, 7700, and 7900 SDS, and LightCycler 480 Instrument

We recommend testing at least 8 cDNA samples in the same experiment to optimize the use of the standards and the primers and probe mixes. The plate scheme in Figure 4 shows an example of such an experiment.



**Figure 4. Suggested plate setup for one experiment.** S: cDNA sample; F1–5: BCR-ABL mbc standards; C1–3: ABL standards; H<sub>2</sub>O: water control.

## qPCR on ABI PRISM 7000, 7700, and 7900 SDS, and LightCycler 480 Instrument

**Note:** Perform all steps on ice.

### Procedure

1. Thaw all necessary components and place them on ice.
2. Prepare the following qPCR mix according to the number of samples being processed. If using 96-well-plate qPCR equipment, we recommend performing all measurements in duplicate.

All concentrations are for the final volume of the reaction.

Table 7 describes the pipetting scheme for the preparation of one reagent mix, calculated to achieve a final reaction volume of 25  $\mu$ l. A pre-mix can be prepared, according to the number of reactions, using the same primer and probe mix (either PPC-ABL or PPF-mbcr). Extra volumes are included to compensate for pipetting error.

**Table 7. Preparation of qPCR mix**

<b>Component</b>	<b>1 reaction (<math>\mu</math>l)</b>	<b>ABL: 24+1 reactions (<math>\mu</math>l)</b>	<b>BCR-ABL mbr: 28+1 reactions (<math>\mu</math>l)</b>	<b>Final concentration</b>
TaqMan Universal PCR Master Mix, 2x	12.5	312.5	362.5	1x
Primers and probe mix, 25x	1	25	29	1x
Nuclease- free PCR grade water	6.5	162.5	188.5	–
Sample (to be added at step 4)	5	5 each	5 each	–
Total volume	25	25 each	25 each	–

3. Dispense 20  $\mu$ l of the qPCR pre-mix per well.
4. Add 5  $\mu$ l of the RT product (cDNA, 100 ng RNA equivalent) obtained in the reverse transcription (see “Protocol: Recommended standardized EAC reverse transcription”, page 12) in the corresponding well (total volume 25  $\mu$ l).
5. Mix gently, by pipetting up and down.
6. Close the plate and briefly centrifuge (300 x g, approximately 10 seconds).
7. Place the plate in the thermal cycler according to the manufacturer recommendations. Program the thermal cycler with the thermal cycling program as indicated in Table 8 for ABI PRISM 7000, 7700, and 7900HT SDS, or Table 9 for the LightCycler 480 Instrument.

**Table 8. Temperature profile for ABI PRISM 7000, 7700, and 7900HT SDS**

<b>Mode of analysis</b>	Standard Curve — Absolute Quantitation
<b>Hold</b>	Temperature: 50°C Time: 2 minutes
<b>Hold 2</b>	Temperature: 95°C Time: 10 minutes
<b>Cycling</b>	50 times 95°C for 15 seconds 60°C for 1 minute with acquisition of FAM fluorescence; quencher: TAMRA

**Table 9. Temperature profile for LightCycler 480 Instrument**

<b>Mode of analysis</b>	Absolute Quantification (“Abs Quant”)
<b>Detection formats</b>	Select “Simple Probe” in the Detection formats window
<b>Hold</b>	Temperature: 50°C Time: 2 minutes
<b>Hold 2</b>	Temperature: 95°C Time: 10 minutes
<b>Cycling</b>	50 times 95°C for 15 seconds 60°C for 1 minute with acquisition of FAM fluorescence corresponding to (483–533 nm) for LC version 01 and (465–510 nm) for LC version 02

**8. For the ABI PRISM 7000, 7700, and 7900HT SDS, follow step 8a. For the LightCycler 480 Instrument, follow step 8b.**

**8a. ABI PRISM 7000, 7700, and 7900HT SDS: We recommend a threshold set at 0.1 as described in the EAC protocol in the analysis step on the ABI PRISM SDS and a baseline set between cycles 3 and 15. Start the cycling program, as indicated in Table 8.**

**8b. LightCycler 480 Instrument: We recommend a Fit point analysis mode with background at 2.0 and threshold at 2.0. Start the thermal cycling program, as indicated in Table 9.**

## Protocol: qPCR on LightCycler 1.2 and 2.0 Instruments


Using capillary instruments, we recommend measuring samples in duplicate and controls only once, as indicated in Table 10.

**Table 10. Number of reactions for LightCycler 1.2 and 2.0 Instruments**

<b>Samples</b>	<b>Reactions</b>
<b>With the ABL primers and probe mix (PPC-ABL)</b>	
n cDNA samples	n x 2 reactions
ABL standard	1 x 3 reactions (3 standard dilutions, each one tested once)
Water control	1 reaction
<b>With the BCR-ABL mbcpr primers and probe mix (PPF-mbcpr)</b>	
n cDNA samples	n x 2 reactions
mbcpr standard	1 x 5 reactions (5 standard dilutions, each one tested once)
Water control	1 reaction

### Sample processing on LightCycler 1.2 and 2.0 Instruments

We recommend testing at least 5 cDNA samples in the same experiment to optimize the use of the standards and primers and probe mixes. The capillary scheme in Figure 5 shows an example of an experiment.



**Figure 5. Suggested rotor setup for each experiment with the *ipsogen* BCR-ABL1 mbc Kit.** F1–5: BCR-ABL mbc standards; C1–3: ABL standards; S: unknown DNA sample to be analyzed; H<sub>2</sub>O: water control.

## qPCR on LightCycler 1.2 and 2.0 Instruments

**Note:** Because of particular technological requirements, LightCycler experiments must be performed using specific reagents. We recommend to use the LightCycler TaqMan Master and to follow the manufacturer's instructions to prepare the Master Mix 5x.

**Note:** Perform all steps on ice.

### Procedure

1. Thaw all necessary components and place them on ice.
2. Prepare the following qPCR mix according to the number of samples being processed.

All concentrations are for the final volume of the reaction.

Table 11 describes the pipetting scheme for the preparation of one reagent mix, calculated to achieve a final reaction volume of 20  $\mu$ l. A pre-mix can be prepared, according to the number of reactions, using the same primers and probe mix (either PPC-ABL or PPF-mbcr). Extra volumes are included to compensate for pipetting error.

**Table 11. Preparation of qPCR mix**

<b>Component</b>	<b>1 reaction (<math>\mu</math>l)</b>	<b>ABL: 14+1 reactions (<math>\mu</math>l)</b>	<b>BCR-ABL mbr: 16+1 reactions (<math>\mu</math>l)</b>	<b>Final concentration</b>
Freshly prepared LightCycler TaqMan Master Mix, 5x	4.0	60	68.0	1x
Primers and probe mix, 25x	0.8	12	13.6	1x
Nuclease-free PCR grade water	10.2	153	173.4	–
Sample (to be added at step 4)	5.0	5 each	5.0 each	–
Total volume	20.0	20 each	20.0 each	–

3. Dispense 15  $\mu$ l of the qPCR pre-mix per capillary.
4. Add 5  $\mu$ l of the RT product (cDNA, 100 ng RNA equivalent) obtained in the reverse transcription (see "Protocol: Recommended standardized EAC reverse transcription", page 12) in the corresponding tube (total volume 20  $\mu$ l).
5. Mix gently, by pipetting up and down.
6. Place the capillaries in the adapters provided with the apparatus, and briefly centrifuge (700 x g, approximately 10 seconds).
7. Load the capillaries into the thermal cycler according to the manufacturer recommendations.
8. Program the LightCycler 1.2 or 2.0 Instrument with the thermal cycling program as indicated in Table 12.

**Table 12. Temperature profile**

<b>Mode of analysis</b>	Quantification
<b>Hold</b>	Temperature: 95°C Time: 10 minutes Ramp: 20
<b>Cycling</b>	50 times 95°C for 10 seconds; ramp: 20 60°C for 1 minute; ramp: 20; with acquisition of FAM fluorescence: Single
<b>Hold 2</b>	45°C for 1 minute; ramp: 20

9. For the LightCycler 1.2, follow step 9a. For the LightCycler 2.0, follow step 9b.
- 9a. LightCycler 1.2: The F1/F2 and “2<sup>nd</sup> derivative analysis” mode is recommended. Start the thermal cycling program, as indicated in Table 12.
- 9b. LightCycler 2.0: We recommend using Automated (F''max) analysis on LightCycler 2.0 Software version 4.0 to obtain reproducible results. Start the thermal cycling program, as indicated in Table 12.

## Protocol: qPCR on the SmartCycler instrument


Using this instrument, we recommend measuring samples in duplicate and controls only once, as indicated in Table 13.

**Table 13. Number of reactions for the SmartCycler instrument**

Samples	Reactions
<b>With the ABL primers and probe mix (PPC-ABL)</b>	
n cDNA samples	n x 2 reactions
ABL standard	1 x 3 reactions (3 standard dilutions, each one tested once)
Water control	1 reaction
<b>With the BCR-ABL mbc r primers and probe mix (PPF-mbc r)</b>	
n cDNA samples	n x 2 reactions
mbcr standard	1 x 5 reactions (5 standard dilutions, each one tested once)
Water control	1 reaction

### Sample processing on the SmartCycler instrument

We recommend testing at least 5 cDNA samples in the same experiment to optimize the use of the standards and primers and probe mixes. The two-block scheme in Figure 6 shows an example.



All the assays on this first block are performed with PPC-ABL.

All the assays on this second block are performed with PPF-mbc r.

**Figure 6. Suggested plate setup for one experiment.** S: cDNA sample; F1–5: BCR-ABL mbc r standards; C1–3: ABL standards; H<sub>2</sub>O: water control.

### qPCR on the SmartCycler instrument

**Note:** Perform all steps on ice.

## Procedure

1. **Thaw all necessary components and place them on ice.**
2. **Prepare the following qPCR mix according to the number of samples being processed.**

All concentrations are for the final volume of the reaction.

Table 14 describes the pipetting scheme for the preparation of one reagent mix, calculated to achieve a final reaction volume of 25  $\mu$ l. A pre-mix can be prepared, according to the number of reactions, using the same primer and probe mix (either PPC-ABL or PPF-mbcr). Extra volumes are included to compensate for pipetting error.

**Table 14. Preparation of qPCR mix**

<b>Component</b>	<b>1 reaction (<math>\mu</math>l)</b>	<b>ABL: 14+1 reactions (<math>\mu</math>l)</b>	<b>BCR-ABL mbcr: 16+1 reactions (<math>\mu</math>l)</b>	<b>Final concentration</b>
TaqMan Universal PCR Master Mix, 2x	12.5	187.5	212.5	1x
Primers and probe mix, 25x	1	15	17	1x
Nuclease- free PCR grade water	6.5	97.5	110.5	–
Sample (to be added at step 4)	5	5 each	5 each	–
Total volume	25	25 each	25 each	–

3. **Dispense 20  $\mu$ l of the qPCR pre-mix per well.**
4. **Add 5  $\mu$ l of the RT product (cDNA, 100 ng RNA equivalent) obtained in the reverse transcription (see “Protocol: Recommended standardized EAC reverse transcription”, page 12) in the corresponding tube (total volume 25  $\mu$ l).**

5. Mix gently, by pipetting up and down.
6. Load the samples into the thermal cycler according to the manufacturer recommendations.
7. Program the SmartCycler instrument with the thermal cycling program as indicated in Table 15.

**Table 15. Temperature profile**

<b>Hold</b>	Temperature: 50°C Time: 2 minutes
<b>Hold 2</b>	Temperature: 95°C Time: 10 minutes
<b>Cycling</b>	50 times 95°C for 15 seconds 60°C for 1 minute with acquisition: Single


8. We recommend a threshold set at 30. Start the thermal cycling program, as indicated in Table 15.

# Interpretation of Results

## Data analysis principle


Using TaqMan technology, the number of PCR cycles necessary to detect a signal above the threshold is called the threshold cycle ( $C_T$ ) and is directly proportional to the amount of target present at the beginning of the reaction.

Using standards with a known number of molecules, one can establish a standard curve and determine the precise amount of target present in the test sample. The *ipsogen* standard curves are plasmid-based; we use 3 plasmid standard dilutions for the CG, and 5 standard dilutions for the FG, in order to ensure accurate standard curves. Figures 7 and 8 show an example of TaqMan amplification curves obtained with the *ipsogen* BCR-ABL mbcR Kit.



**Figure 7. Detection of ABL standards (C1, C2, C3). 10<sup>3</sup>, 10<sup>4</sup>, and 10<sup>5</sup> copies/5  $\mu$ l.**

- m-bcr  $10^1$
- m-bcr  $10^2$
- m-bcr  $10^3$
- m-bcr  $10^5$
- m-bcr  $10^6$




**Figure 8. Detection of BCR-ABL mbc standards (F1–F5).  $10^1$ ,  $10^2$ ,  $10^3$ ,  $10^5$ ,  $10^6$  copies/5  $\mu$ l.**

## Results

### Standard curve and quality criteria

Raw data can be pasted into an Excel® file for analysis.

For each gene (ABL and BCR-ABL), raw  $C_T$  values obtained from plasmid standard dilutions are plotted according to the log copy number (3, 4, and 5 for C1, C2, and C3; 1, 2, 3, 5, and 6 for F1, F2, F3, F4, and F5). Figure 9 shows an example of the theoretical curve calculated on 5 standard dilutions.



**Figure 9. Theoretical curve calculated from the 5 standard dilutions.** A linear regression curve ( $y = ax + b$ ) is calculated for each gene (ABL and BCR-ABL), where  $a$  is the slope of the line and  $b$  is the y-intercept, which is the y-coordinate of the point where the line crosses the y axis. Its equation and coefficient of determination ( $R^2$ ) are printed on the graph.

As standards are tenfold dilutions, the theoretical slope of the curve is  $-3.3$ . A slope between  $-3.0$  and  $-3.9$  is acceptable as long as  $R^2$  is  $>0.95$  (2). However, a value for  $R^2 >0.98$  is desirable for precise results (3).

### Normalized copy number (NCN)

The ABL standard curve equation should be used to transform raw  $C_T$  values (obtained with PPC-ABL) for the unknown samples into ABL copy numbers ( $ABL_{CN}$ ).

The BCR-ABL standard curve equation should be used to transform raw  $C_T$  values (obtained with PPF-mbcr) for the unknown samples, into BCR-ABL copy numbers ( $BCR-ABL\ mbc_{CN}$ ).

The ratio of these CN values gives the normalized copy number (NCN):

$$NCN = \frac{BCR-ABL\ mbc_{CN}}{ABL_{CN}} \times 100$$

### MRD value

The minimal residual disease (MRD) value is the ratio between the CG normalized expression of the FG in follow-up ( $(FG_{CN}/CG_{CN})_{FUP}$ ) and diagnostic samples ( $(FG_{CN}/CG_{CN})_{DX}$ ).

$$MRD\ value\ (MRD_v) = \frac{(FG_{CN}/CG_{CN})_{FUP}}{(FG_{CN}/CG_{CN})_{DX}}$$

## Sensitivity

The sensitivity (SENS<sub>v</sub>) is calculated according to the relative expression of the FG at diagnosis (FG<sub>CN</sub>/CG<sub>CN</sub>)<sub>DX</sub> and CG expression (CG<sub>CN,FUP</sub>) in the follow-up sample.

$$\text{Sensitivity (SENS}_v\text{)} = \frac{\text{CG}_{\text{CN,DX}}}{\text{CG}_{\text{CN,FUP}} \times \text{FG}_{\text{CN,DX}}}$$

## Quality control on ABL values

Poor quality of the RNA or problems during the qPCR steps result in low ABL<sub>CN</sub>. We recommend discarding results from samples giving ABL<sub>CN</sub> < 1318 (lower value of the 95% CI from patient samples in the EAC study, reference 4).

## Reproducibility between replicates

The variation in C<sub>T</sub> values between replicates should be < 2, corresponding to a fourfold change in copy number values.

Variation in C<sub>T</sub> values between replicates is generally < 1.5 if the mean C<sub>T</sub> value of the replicates is < 36 (2).

**Note:** Each user should measure their own reproducibility in their laboratory.

## Water controls

Negative controls should give zero CN.

A positive water control results from a cross-contamination. See "Troubleshooting guide", below, to find a solution.

## Troubleshooting guide

This troubleshooting guide may be helpful in solving any problems that may arise. For more information, see also the Frequently Asked Questions page at our Technical Support Center: [www.qiagen.com/FAQ/FAQList.aspx](http://www.qiagen.com/FAQ/FAQList.aspx). The scientists in QIAGEN Technical Services are always happy to answer any questions you may have about either the information and protocol in this handbook or sample and assay technologies (for contact information, see "Contact Information", page 45).

## Comments and suggestions

---

### Negative result for the control gene (ABL) and BCR-ABL mbcr in all the samples — standard okay

- a) Poor RNA quality                      Always check the RNA quality and concentration before starting.  
Run a cell line RNA positive control (High positive control in *ipsogen* BCR-ABL1 mbcr Controls Kit, cat. no. 670091) in parallel.
- b) Failure of reverse transcription step                      Always check the RNA quality and concentration before starting.  
Run a cell line RNA positive control (*ipsogen* BCR-ABL1 mbcr Controls Kit, cat. no. 670091) in parallel.

### Negative result for the control gene (ABL) in the samples — standard okay

- a) Poor RNA quality                      Always check the RNA quality and concentration before starting.  
Run a cell line RNA positive control (*ipsogen* BCR-ABL1 mbcr Controls Kit, cat. no. 670091) in parallel.
- b) Failure of reverse transcription step                      Always check the RNA quality and concentration before starting.  
Run a cell line RNA positive control (*ipsogen* BCR-ABL1 mbcr Controls Kit, cat. no. 670091) in parallel.

### Standard signal negative

- a) Pipetting error                      Check pipetting scheme and the setup of the reaction.  
Repeat the PCR run.
- b) Inappropriate storage of kit components                      Store the *ipsogen* BCR-ABL1 mbcr Kit at –15 to –30°C and keep primers and probe mixes (PPC and PPF) protected from light. See “Reagent Storage and Handling”, page 11.  
Avoid repeated freezing and thawing.  
Aliquot reagents for storage.

## Comments and suggestions

---

### Negative controls are positive

Cross-contamination	Replace all critical reagents. Repeat the experiment with new aliquots of all reagents. Always handle samples, kit components, and consumables in accordance with commonly accepted practices to prevent carry-over contamination.
---------------------	--

### No signal, even in standard controls

a) Pipetting error or omitted reagents	Check pipetting scheme and the setup of the reaction. Repeat the PCR run.
b) Inhibitory effects of the sample material, caused by insufficient purification	Repeat the RNA preparation.
c) LightCycler: Incorrect detection channel chosen	Set Channel Setting to F1/F2 or 530 nm/640 nm.
d) LightCycler: No data acquisition programmed	Check the cycle programs. Select acquisition mode "single" at the end of each annealing segment of the PCR program.

### Absent or low signal in samples but standard controls okay

a) Poor RNA quality or low concentration	Always check the RNA quality and concentration before starting. Run a cell line RNA positive control ( <i>ipsogen</i> BCR-ABL1 mbcR Controls Kit, cat. no. 670091) in parallel.
b) Failure of reverse transcription step	Always check the RNA quality and concentration before starting. Run a cell line RNA positive control ( <i>ipsogen</i> BCR-ABL1 mbcR Controls Kit, cat. no. 670091) in parallel.

## Comments and suggestions

---

### Fluorescence intensity too low

- a) Inappropriate storage of kit components      Store the *ipsogen* BCR-ABL1 mbcR Kit at –15 to –30°C and keep primers and probe mixes (PPC and PPF) protected from light. See “Reagent Storage and Handling”, page 11.  
Avoid repeated freezing and thawing.  
Aliquot reagents for storage.
- b) Very low initial amount of target RNA      Increase the amount of sample RNA.  
**Note:** Depending of the chosen method of RNA preparation, inhibitory effects may occur.

### LightCycler: Fluorescence intensity varies

- a) Pipetting error      Variability caused by so-called “pipetting error” can be reduced by analyzing data in the F1/F2 or 530 nm/640 nm mode.
- b) Insufficient centrifugation of the capillaries      The prepared PCR mix may still be in the upper vessel of the capillary, or an air bubble could be trapped in the capillary tip.  
Always centrifuge capillaries loaded with the reaction mix as described in the specific operating manual of the apparatus.
- c) Outer surface of the capillary tip dirty      Always wear gloves when handling the capillaries.

### LightCycler: Error of the standard curve

- Pipetting error      Variability caused by so-called “pipetting error” can be reduced by analyzing data in the F1/F2 or 530 nm/640 nm mode.

## Quality Control

Quality control of the complete kit has been performed on a LightCycler 480 Instrument. This kit is manufactured according to ISO 13485:2003 standard. Certificates of analysis are available on request at [www.qiagen.com/support/](http://www.qiagen.com/support/).

## Limitations

The users must be trained and familiar with this technology prior to the use of this device.

Any diagnostic results generated must be interpreted in conjunction with other clinical or laboratory findings. It is the user's responsibility to validate system performance for any procedures used in their laboratory which are not covered by the QIAGEN performance studies.

Attention should be paid to expiration dates printed on the box and labels of all components. Do not use expired components.

**Note:** The kit has been designed according to the "Europe Against Cancer" (EAC) studies (4), and is compliant with the updated international recommendations (3, 5). It should be used following the instructions given in this manual, in combination with validated reagents and instruments (see "Materials Required but Not Provided", page 9). Any off-label use of this product and/or modification of the components will void QIAGEN's liability.

## Performance Characteristics


### Nonclinical studies

#### Materials and methods

Performance evaluation was performed on an ABI PRISM 7700 SDS, in combination with reagents listed in "Materials Required but Not Provided", page 9. Equivalence studies validated its use on the following instruments: ABI PRISM 7000 and 7900HT SDS, LightCycler 1.2 and 480 Instruments, Rotor-Gene 3000, and SmartCycler instrument (6).

Nonclinical studies were conducted to establish the analytical performance of the *ipsogen* BCR-ABL1 mbcR Kit. These nonclinical laboratory studies were performed on total RNA from TOM1 cell line diluted in a constant final amount of MV4-11 cell line total RNA.

To determine the repeatability of the assay, 5 different concentrations of TOM1 total RNA (5 ng, 500 pg, 50 pg, 5 pg, and 0.5 pg) diluted in MV4-11 total RNA, in a constant final total amount of 1000 ng, were analyzed in 5 replicates per run and in 4 different runs (Figure 10).



**Figure 10. Amplification plots of  $5 \times 10^{-3}$  (5 ng),  $5 \times 10^{-4}$  (0.5 ng), and  $5 \times 10^{-5}$  (0.05 ng) dilutions of TOM1 total RNA in MV4-11 negative total RNA.**

### Analytical data

Tables 16–19 show the inter-assay analyses with the mean threshold cycle ( $C_T$ ), standard deviation (SD), number of samples (n), coefficient of variation (CV), mean copy number (CN), and mean normalized copy number (NCN).

**Table 16. Inter-assay analysis — cell lines mbcr and ABL**

Cell line	Dilution	Mean $C_T$	SD	n	CV (%)
mbcr	$5 \times 10^{-3}$ (5 ng/1 $\mu$ g)	29.19	0.26	20	0.88
	$5 \times 10^{-4}$ (0.5 ng/1 $\mu$ g)	33.70	0.48	20	1.47
	$5 \times 10^{-5}$ (0.05 ng/1 $\mu$ g)	37.03	1.16	20	3.15
ABL	–	25.01	0.87	100	3.46

**Table 17. Inter-assay analysis — plasmids**

Gene	Plasmid	Mean C <sub>T</sub>	SD	n	CV (%)
mbcr	F1 (10 <sup>1</sup> copies)	35.19	0.90	11	2.57
	F2 (10 <sup>2</sup> copies)	31.87	0.64	12	1.99
	F3 (10 <sup>3</sup> copies)	28.41	0.71	12	2.50
	F4 (10 <sup>5</sup> copies)	21.48	0.59	12	2.76
	F5 (10 <sup>6</sup> copies)	18.37	0.71	12	3.89
ABL	C1 (10 <sup>3</sup> copies)	29.68	0.85	12	2.86
	C2 (10 <sup>4</sup> copies)	26.01	0.51	12	1.96
	C3 (10 <sup>5</sup> copies)	22.53	0.42	12	1.86

**Table 18. Inter-assay analysis — cell lines BCR-ABL mbcr and ABL (mean CN)**

Cell line	Dilution	Mean CN	SD	n	CV (%)
BCR-ABL mbcr	5 x 10 <sup>-3</sup> (5 ng/1 μg)	587.30	194.10	20	33.05
	5 x 10 <sup>-4</sup> (0.5 ng/1 μg)	57.84	20.38	20	35.23
	5 x 10 <sup>-5</sup> (0.05 ng/1 μg)	4.39	2.73	20	62.35
ABL	—	22,038.22	9459.17	100	42.92

**Table 19. Inter-assay analysis — cell line BCR-ABL mbcr (mean NCN)**

Cell line	Dilution	Mean NCN*	SD	n	CV (%)
BCR-ABL mbcr	5 x 10 <sup>-3</sup> (5 ng/1 μg)	267.46	93.22	20	34.85
	5 x 10 <sup>-4</sup> (0.5 ng/1 μg)	23.54	7.36	20	31.28
	5 x 10 <sup>-5</sup> (0.05 ng/1 μg)	2.60	2.80	20	107.66

\* For these study results only, the NCN is given as  $\frac{\text{BCR-ABL mbcr}_{\text{CN}}}{\text{ABL}_{\text{CN}}} \times 10,000$ .

## Clinical studies

Performance evaluation was performed on an ABI PRISM 7700 SDS, in combination with reagents listed in “Materials Required but Not Provided”, page 9. Equivalence studies validated its use on the following instruments: ABI PRISM 7000 and 7900HT SDS, LightCycler 1.2 and 480 Instruments, Rotor-Gene 3000, and SmartCycler instrument (6).

A group of 26 laboratories, in 10 European countries, organized in a Europe Against Cancer (EAC) concerted action, used plasmids provided by IPSOGEN to establish a standardized protocol for qPCR analysis of the major leukemia-associated fusion genes in the clinical setting. The BCR-ABL p190 transcript was one of the fusion genes (FG) included in the study. We present here a summary of this validation study; full results have been published in 2003 (4, 7).

### Inter-laboratory reproducibility for CG and FG plasmid standards

Eleven laboratories performed an inter-laboratory reproducibility experiment to assess variability in the measurement of CG and FG plasmid standard dilutions. Dilutions were performed in duplicate at each facility. Table 20 reports the mean, standard deviation, and CV (%) for each dilution.

**Table 20. Inter-laboratory reproducibility for CG and FG plasmid standards**

Gene	Dilution	Mean	C <sub>T</sub> SD	CV (%)
ABL control gene	C1	29.04	0.53	1.82
	C2	25.64	0.47	1.84
	C3	22.10	0.34	1.55
BCR-ABL mbcr fusion gene	F1	35.99	1.18	3.28
	F2	32.05	0.74	2.32
	F3	28.43	0.65	2.29
	F4	21.60	0.59	2.72
	F5	18.24	0.46	2.57

### Expression values of the BCR-ABL mbcr FG transcript

Tables 21 and 22 show the expression values of the BCR-ABL mbcr FG transcript and ABL CG, for TOM1 cell line, ALL patients at diagnosis, and normal patients.

**Table 21. Expression values of the BCR-ABL mbc<sub>r</sub> FG transcript and ABL CG — C<sub>T</sub> values**

	C <sub>T</sub> values (95% range)	
	BCR-ABL mbc <sub>r</sub>	ABL
<b>TOM1 cell line</b>	22.8	21.8
<b>ALL patient samples</b>		
BM (n = 17)	24.7 (21.3–27.1)	24.5 (21.7–27.1)
PB (n = 7)	23.3 (21.7–29.1)	22.5 (21.0–27.0)
<b>Negative patient samples</b>		
BM (n = 26)	–	25.35 (24.68–26.02)
PB (n = 74)	–	25.15 (24.83–25.48)

**Table 22. Expression values of the BCR-ABL mbc<sub>r</sub> FG transcript and ABL CG — CN and NCN values**

	CN values (95% range)		NCN values (95% range)
	BCR-ABL mbc <sub>r</sub>	ABL	CN BCR-ABL mbc <sub>r</sub> /CN ABL
<b>ALL patient samples</b>			
BM (n = 17)	9550 (1738–97,724)	11,912 (5012–70,795)	0.8 (0.35–1.38)
PB (n = 7)	91,201 (1905–208,930)	134,896 (4786–114,815)	0.68 (0.4–1.82)
<b>Negative patient samples</b>			
BM (n = 26)	–	19,201 (12,922–25,480)	–
PB (n = 74)	–	21,136 (17,834–24,437)	–

ABL C<sub>T</sub> values did not differ significantly between normal and leukemic samples, nor between samples types (PB or BM) or leukemia samples (ALL, AML, CML).

## False positive and false negative rates

False negative and false positive rates were computed using the following controls.

- Positive controls: TOM1 cells, a cell line well known for its positivity for BCR-ABL p190 fusion gene; patients samples already assessed for p190 positivity
- Negative controls: Negative RNA samples, no amplification controls (NAC) made of *E. coli* RNA instead of human RNA to check for PCR contamination, and no template controls (NTC), which contained water instead of human RNA

Amplification on RNA samples of the FG was done in triplicate and in duplicate for the CG.

A false-negative sample was defined as a positive RNA sample with less than 50% of positive wells (0/2, 0/3, or 1/3).

A false-positive sample was defined as a negative sample with at least 50% of positive wells (1/2, 2/3, or 3/3).

Table 23 shows the number and percentage of false negative and false positive samples.

**Table 23. False negative and false positive samples**

False negativity		False positivity	
10 <sup>-3</sup>	10 <sup>-4</sup>	FG negative control	NAC/NTC
0% (0/54)	4% (3/75)	4.8% (6/126)	5.8% (7/120)

## References

QIAGEN maintains a large, up-to-date online database of scientific publications utilizing QIAGEN products. Comprehensive search options allow you to find the articles you need, either by a simple keyword search or by specifying the application, research area, title, etc.

For a complete list of references, visit the QIAGEN Reference Database online at [www.qiagen.com/RefDB/search.asp](http://www.qiagen.com/RefDB/search.asp) or contact QIAGEN Technical Services or your local distributor.

### Cited references

1. Thomas, D.A. (2007) Philadelphia chromosome positive acute lymphocytic leukemia: a new era of challenges. *Hematology Am. Soc. Hematol. Educ. Program* **2007**, 435.

2. van der Velden, V.H., Hochhaus, A., Cazzaniga, G., Szczepanski, T., Gabert, J., and van Dongen, J.J. (2003) Detection of minimal residual disease in hematologic malignancies by real-time quantitative PCR: principles, approaches, and laboratory aspects. *Leukemia* **17**, 1013.
3. Branford, S. et al. (2006) Rationale for the recommendations for harmonizing current methodology for detecting BCR-ABL transcripts in patients with chronic myeloid leukaemia. *Leukemia* **20**, 1925.
4. Gabert, J. et al. (2003) Standardization and quality control studies of 'real-time' quantitative reverse transcriptase polymerase chain reaction of fusion gene transcripts for residual disease detection in leukemia — a Europe Against Cancer program. *Leukemia* **17**, 2318.
5. Hughes, T. et al. (2006) Monitoring CML patients responding to treatment with tyrosine kinase inhibitors: review and recommendations for harmonizing current methodology for detecting BCR-ABL transcripts and kinase domain mutations and for expressing results. *Blood* **108**, 28.
6. Silvy, M., Mancini, J., Thirion, X., Sigaux, F., and Gabert, J. (2005) Evaluation of real-time quantitative PCR machines for the monitoring of fusion gene transcripts using the Europe against cancer protocol. *Leukemia* **19**, 305.
7. Beillard, E. et al. (2003) Evaluation of candidate control genes for diagnosis and residual disease detection in leukemic patients using 'real-time' quantitative reverse-transcriptase polymerase chain reaction (RQ-PCR) - a Europe against cancer program. *Leukemia* **17**, 2474.

## Symbols

The following symbols may appear on the packaging and labeling:



<N>

Contains reagents sufficient for <N> reactions



Use by



In vitro diagnostic medical device



Catalog number



Lot number



Material number



Global Trade Item Number



Temperature limitation



Manufacturer



Consult instructions for use

## Contact Information

For technical assistance and more information, please see our Technical Support Center at [www.qiagen.com/Support](http://www.qiagen.com/Support), call 00800-22-44-6000, or contact one of the QIAGEN Technical Service Departments or local distributors (see back cover or visit [www.qiagen.com](http://www.qiagen.com)).

## Ordering Information

Product	Contents	Cat. no.
<i>ipsogen</i> BCR-ABL1 mbcR Kit (24)	For 24 reactions: ABL Control Gene Standards, BCR-ABL mbcR Fusion Gene Standards, Primer and Probe Mix ABL, Primer and Probe Mix BCR-ABL mbcR Fusion Gene	670023
<b>Rotor-Gene Q MDx — for IVD-validated real-time PCR analysis in clinical applications</b>		
Rotor-Gene Q MDx 5plex HRM Platform	Real-time PCR cycler and High Resolution Melt analyzer with 5 channels (green, yellow, orange, red, crimson) plus HRM channel, laptop computer, software, accessories, 1-year warranty on parts and labor, installation and training not included	9002032
Rotor-Gene Q MDx 5plex HRM System	Real-time PCR cycler and High Resolution Melt analyzer with 5 channels (green, yellow, orange, red, crimson) plus HRM channel, laptop computer, software, accessories, 1-year warranty on parts and labor, installation and training	9002033
<b><i>ipsogen</i> BCR-ABL1 mbcR Controls Kit — for qualitative validation of RNA extraction and reverse transcription of the BCR-ABL mbcR fusion gene</b>		
<i>ipsogen</i> BCR-ABL1 mbcR Controls Kit	Cell lines with negative, high, and low positive expression of the BCR-ABL mbcR fusion gene	670091

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

This product is intended for in vitro diagnostic use. *ipsogen* products may not be resold, modified for resale or used to manufacture commercial products without written approval of QIAGEN.

Information in this document is subject to change without notice. QIAGEN assumes no responsibility for any errors that may appear in this document. This document is believed to be complete and accurate at the time of publication. In no event shall QIAGEN be liable for incidental, special, multiple, or consequential damages in connection with, or arising from the use of this document.

*ipsogen* products are warranted to meet their stated specifications. QIAGEN's sole obligation and the customer's sole remedy are limited to replacement of products free of charge in the event products fail to perform as warranted.

Trademarks: QIAGEN<sup>®</sup>, *ipsogen*<sup>®</sup>, Rotor-Gene<sup>®</sup> (QIAGEN Group); ABI PRISM<sup>®</sup>, FAM<sup>™</sup>, RNaseOUT<sup>™</sup>, SuperScript<sup>®</sup>, SYBR<sup>®</sup>, TAMRA<sup>™</sup> (Life Technologies Corporation); Agilent<sup>®</sup>, Bioanalyzer<sup>®</sup> (Agilent Technologies, Inc.); Excel<sup>®</sup> (Microsoft Corporation); LightCycler<sup>®</sup>, TaqMan<sup>®</sup> (Roche Group); SmartCycler<sup>®</sup> (Cepheid).

### Limited License Agreement

Use of this product signifies the agreement of any purchaser or user of the *ipsogen* BCR-ABL1 mbcR Kit to the following terms:

1. The *ipsogen* BCR-ABL1 mbcR Kit may be used solely in accordance with the *ipsogen BCR-ABL1 mbcR Kit Handbook* and for use with components contained in the Kit only. QIAGEN grants no license under any of its intellectual property to use or incorporate the enclosed components of this Kit with any components not included within this Kit except as described in the *ipsogen BCR-ABL1 mbcR Kit Handbook* and additional protocols available at [www.qiagen.com](http://www.qiagen.com).
2. Other than expressly stated licenses, QIAGEN makes no warranty that this Kit and/or its use(s) do not infringe the rights of third-parties.
3. This Kit and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. QIAGEN specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the Kit agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. QIAGEN may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the Kit and/or its components.

For updated license terms, see [www.qiagen.com](http://www.qiagen.com).

HB-1357-002 © 2013–2015 QIAGEN, all rights reserved.

---

**www.qiagen.com**

**Australia** ■ techservice-au@qiagen.com

**Austria** ■ techservice-at@qiagen.com

**Belgium** ■ techservice-bnl@qiagen.com

**Brazil** ■ suportetecnico.brasil@qiagen.com

**Canada** ■ techservice-ca@qiagen.com

**China** ■ techservice-cn@qiagen.com

**Denmark** ■ techservice-nordic@qiagen.com

**Finland** ■ techservice-nordic@qiagen.com

**France** ■ techservice-fr@qiagen.com

**Germany** ■ techservice-de@qiagen.com

**Hong Kong** ■ techservice-hk@qiagen.com

**India** ■ techservice-india@qiagen.com

**Ireland** ■ techservice-uk@qiagen.com

**Italy** ■ techservice-it@qiagen.com

**Japan** ■ techservice-jp@qiagen.com

**Korea (South)** ■ techservice-kr@qiagen.com

**Luxembourg** ■ techservice-bnl@qiagen.com

**Mexico** ■ techservice-mx@qiagen.com

**The Netherlands** ■ techservice-bnl@qiagen.com

**Norway** ■ techservice-nordic@qiagen.com

**Singapore** ■ techservice-sg@qiagen.com

**Sweden** ■ techservice-nordic@qiagen.com

**Switzerland** ■ techservice-ch@qiagen.com

**UK** ■ techservice-uk@qiagen.com

**USA** ■ techservice-us@qiagen.com

