

# **virotype<sup>®</sup> BTV pan/8 RT-PCR Kit Handbook**



24 (catalog no. 280443)



96 (catalog no. 280445)



480 (catalog no. 280447)\*

For detection of RNA from bluetongue virus

Registered in accordance with § 17c of the German Law on Animal Diseases (Flu-B 539)

**REF**

280443, 280445, 280447\*



QIAGEN Leipzig GmbH, Deutscher Platz 5b,  
04103 Leipzig, Germany



\* Available only on request.

## **QIAGEN Sample and Assay Technologies**

QIAGEN is the leading provider of innovative sample and assay technologies, enabling the isolation and detection of contents of any biological sample. Our advanced, high-quality products and services ensure success from sample to result.

### **QIAGEN sets standards in:**

- Purification of DNA, RNA, and proteins
- Nucleic acid and protein assays
- microRNA research and RNAi
- Automation of sample and assay technologies

Our mission is to enable you to achieve outstanding success and breakthroughs. For more information, visit [www.qiagen.com](http://www.qiagen.com).

In addition, QIAGEN provides high-quality, easy-to-use, and sensitive molecular solutions to enable veterinary pathogen detection and animal pathogen research. The QIAGEN veterinary portfolio includes a broad range of pathogen-specific PCR-assays and an extensive and growing ELISA portfolio. For more information, visit [www.qiagen.com/Animal-and-Veterinary-Testing](http://www.qiagen.com/Animal-and-Veterinary-Testing).

## Contents

<b>Kit Contents</b>	<b>4</b>
<b>Intended Use</b>	<b>4</b>
<b>Symbols</b>	<b>5</b>
<b>Storage</b>	<b>5</b>
<b>Safety Information</b>	<b>6</b>
<b>Quality Control</b>	<b>6</b>
<b>Introduction</b>	<b>7</b>
Principle	7
RNA extraction	8
<b>Equipment and Reagents to Be Supplied by User</b>	<b>10</b>
<b>Important Notes</b>	<b>11</b>
General precautions	11
Negative Control	11
Positive Control	11
Extraction and amplification control	12
<b>Protocols:</b>	
■ <b>Real-time RT-PCR using the Rotor-Gene Q</b>	<b>13</b>
■ <b>Real-time RT-PCR using 96-well plate real-time cyler</b>	<b>16</b>
<b>Data Analysis and Interpretation</b>	<b>19</b>
<b>Troubleshooting Guide</b>	<b>22</b>
<b>Ordering Information</b>	<b>23</b>

## Kit Contents

<b><i>virotype</i> BTV pan/8 RT-PCR Kit</b>	<b>(24)</b>	<b>(96)</b>	<b>(480)</b>
<b>Catalog no.</b>	<b>280443</b>	<b>280445</b>	<b>280447*</b>
<b>Number of reactions</b>	<b>24</b>	<b>96</b>	<b>480</b>
Master Mix (tube with orange cap) includes enzymes, primers, and probes	1 x 500 µl	2 x 980 µl	6 x 1625 µl
Positive Control (tube with red cap)	1 x 25 µl	1 x 70 µl	2 x 50 µl
Negative Control (tube with blue cap)	1 x 25 µl	1 x 70 µl	2 x 50 µl
Handbook	1	1	1

\* Available only on request

## Intended Use

The *virotype* BTV pan/8 RT-PCR Kit is intended for the detection of bluetongue virus RNA in ruminant whole blood (preferred with anticoagulants, for example EDTA-blood) and tissue samples (spleen, lymph nodes) from cattle, sheep, and goats. The kit is approved by the Friedrich-Loeffler-Institut and registered in accordance with § 17c of the German Law on Animal Diseases (FluB 539) for use in Germany for veterinary diagnostic procedures. For veterinary use only.

## Symbols



Contains reagents for <N> tests



Legal manufacturer



Lot number



Use by date



Temperature limitations for storage



Handbook



Catalog number



Material number



Protect from light



For cattle, sheep, and goat samples

## Storage

The components of the *virotype* BTV pan/8 RT-PCR Kit should be stored at  $-15$  to  $-30^{\circ}\text{C}$  and are stable until the expiration date stated on the label. Avoid repeated thawing and freezing ( $>2x$ ), as this may reduce assay sensitivity. Freeze the components in aliquots if they will only be used intermittently.

## Safety Information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at [www.qiagen.com/safety](http://www.qiagen.com/safety) where you can find, view, and print the SDS for each QIAGEN kit and kit component.

All sample residues and objects which have come into contact with samples must be decontaminated or disposed of as potentially infective material.

## Quality Control

In accordance with QIAGEN's ISO-certified Quality Management System, each lot of *virotype* BTV pan/8 RT-PCR Kit is tested against predetermined specifications to ensure consistent product quality.

## Introduction

Bluetongue is an infectious, non-contagious disease of ruminants. The agent is the bluetongue virus (BTV), a double-stranded RNA virus of the genus *Orbivirus* of the family *Reoviridae* which includes at least 25 known serotypes. BTV is widely distributed around the world. Sheep, cattle, and goats are mainly affected by the disease. Clear clinical signs are usually seen only in sheep. In severe cases the tongue may show intense hyperemia and become cyanotic (Bluetongue). BTV serotype 8 is of epidemiological importance in Central Europe and cause of recent major Bluetongue Disease outbreaks. The virus is transmitted by certain midges of the genus *Culicoides*. Furthermore, the virus can be spread by contaminated needles and surgery equipment.

## Principle

Polymerase chain reaction (PCR) is based on the amplification of specific regions of the pathogen genome. In real-time RT-PCR, the amplified product is detected using fluorescent dyes. These are usually linked to oligonucleotide probes that bind specifically to the amplified product. Monitoring the fluorescence intensities during the PCR run (i.e., in real time) allows the detection of the accumulating product without the need to re-open the reaction tubes afterward.

The *virotype* BTV pan/8 RT-PCR Kit contains all of the necessary reagents for the detection of BTV RNA, including a positive and negative control. With this kit, both reverse transcription and PCR are performed in one reaction tube, reducing the risk of contamination.

The *virotype* BTV pan/8 RT-PCR Kit uses three specific primer/probe combinations: one for the RNA of all 25 known BTV serotypes yielding FAM™ fluorescence, one for BTV-8 RNA yielding Cy5™ fluorescence, and one for a housekeeping gene ( $\beta$ -actin mRNA), present within the sample, yielding HEX™ fluorescence.

The Positive Control contains BTV-8 RNA and allows the control of the denaturation step since the successful denaturation of the viral double-stranded RNA is a prerequisite for amplification.

## RNA extraction

*virotype* BTV pan/8 can be used for the detection of BTV RNA from ruminant whole blood (preferred with anticoagulants, for example EDTA-blood) and tissue samples (spleen, lymph nodes). Due to the high sensitivity of the test, pools of up to 10 individual blood samples may be analyzed. However, the optimum pool size depends on the regional prevalence for BTV.

Prior to real-time RT-PCR, viral RNA must be extracted from the starting material. QIAGEN offers a range of products for RNA extraction from animal samples.

- QIAamp® *cador*® Pathogen Mini Kit
- QIAamp Viral RNA Mini Kit
- RNeasy® Fibrous Tissue Mini Kit for tissue
- RNeasy Mini Kit

If real-time RT-PCR is not performed immediately after extraction, store the RNA at -20°C or at -70°C for longer storage.

RNA extraction using kits based on spin-column technology can be automated using the QIAcube®.

## Equipment and Reagents to Be Supplied by User

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, consult the appropriate safety data sheets (SDSs), available from the product supplier.

- Pipets
- Nuclease-free, aerosol-resistant pipet tips with filters
- Sterile 1.5 ml Eppendorf® tubes
- Nuclease-free (RNase/DNase-free) consumables. Special care should be taken to avoid nuclease contamination of all reagents and consumables used to set up PCR for sensitive identification of viral nucleic acids.
- Cooling device and ice or liquid nitrogen
- 96-well plate standard PCR cycler
- Benchtop centrifuge with rotor for 1.5 ml tubes
- Rotor-Gene® Q or 96-well plate real-time cycler with appropriate fluorescent channels
- Rotor-Gene Q software version 1.7.94 or higher, or appropriate software for chosen 96-well plate cycler
- Strip Tubes and Caps, 0.1 ml, for use with Rotor-Gene Q (cat. no. 981103 or 981106) or PCR tubes and Caps, 0,2 ml or 96-well optical microplate with optical sealing film or cover for chosen 96-well plate real-time cycler

## Important Notes

### General precautions

The user should always pay attention to the following:

- Use nuclease-free pipet tips with filters.
- Store and extract positive materials (specimens, positive controls, and amplicons) separately from all other reagents, and add them to the reaction mix in a spatially separated facility.
- Thaw all components on ice before starting an assay.
- When thawed, mix the components by inverting and centrifuge briefly.
- Do not use components of the test kit past the expiration date.
- Keep samples and controls on ice or in a cooling block during the setup of reactions.

### Negative Control

At least one negative control reaction should be included in each PCR run. This enables assessment of contamination in the reaction.

### Positive Control

When performing PCR on unknown samples, it is recommended to perform a positive control reaction in the PCR run, containing a sample that is known to include the targeted viral RNA. A positive control serves to prove the functionality of the pathogen assay, for example, the correct setup of the reaction mix. Use 5  $\mu$ l of the Positive Control provided with the *virotype* BTV pan/8 RT-PCR Kit to test for successful amplification of the target.

## **Extraction and amplification control**

For increased process safety and convenience, an extraction and amplification control assay is included in the form of an additional primer/probe set that detects a housekeeping gene present within the sample. This allows both extraction and amplification to be monitored.

# Protocol: Real-time RT-PCR using the Rotor-Gene Q

## Important points before starting

- Please read “Important Notes” on page 11 before starting.
- Include at least one positive control (Positive Control) and one negative control (Negative Control) per PCR run.
- Before beginning the procedure, read through the protocol and ensure that you are familiar with the operation of the chosen real-time PCR cycler.
- RNA is unstable. Perform the protocol without interruption.

## Things to do before starting

- Thaw all reagents on ice and protect from light.
- Maintain reagents on ice during PCR setup.
- Before use, spin the reagents briefly.

## Procedure

- 1. Pipet at least 7  $\mu$ l of RNA samples or Positive Control into individual 0.2 ml PCR reaction tubes. Cover the reaction tubes (e.g., with PCR sealing film).**  
Include positive control reactions.  
Positive control: Use at least 7  $\mu$ l of the positive control (Positive Control) instead of sample RNA.
- 2. Denature the samples for 5 min at 98°C in a 96-well plate standard cycler with a heated lid.**
- 3. Immediately cool down on ice water or liquid nitrogen for at least 20 s. Then store on ice or cooling device.**

- Pipet 5  $\mu$ l of RNA samples, Positive Control, and Negative Control into individual Strip Tubes and Caps, 0.1 ml, for use with Rotor-Gene Q. Use 5  $\mu$ l of the negative control (Negative Control) instead of sample RNA.
- Add 20  $\mu$ l of the Master Mix into each reaction tube. Thus the final volume is 25  $\mu$ l (Table 1).

**Table 1. Preparation of reaction mix**

<b>Component</b>	<b>Volume</b>
Master Mix	20 $\mu$ l
Sample	5 $\mu$ l
<b>Total volume</b>	25 $\mu$ l

- Close the reaction tubes with the corresponding caps.
- Set the filters for the reporter dyes in the software of your thermal cycler according to Table 2. Select the green, red, and, yellow channels on the Rotor-Gene Q.

**Table 2. Filter settings for reporter**

<b>Pathogen/internal control</b>	<b>Reporter</b>
BTV	FAM
BTV-8	Cy5
Internal Control	HEX/JOE™*

\* Use the option appropriate for your thermal cycler

8. Run the real-time RT-PCR protocol according to Table 3 if running only the *virotype* BTV/pan8 RT-PCR Kit.

**Table 3. Real-time RT-PCR protocol for BTV pan/8**

Temperature	Time	Number of cycles
50°C	10 min	1
95°C	10 min	1
95°C	15 s	40
60°C*	60 s	

\* Fluorescence data collection.

9. Run the real-time RT-PCR protocol according to Table 4 if running other *virotype* assays simultaneously (i.e., *virotype* PRRSV, *virotype* BVDV, *virotype* CSFV, *virotype* SBV and/or *virotype* Influenza A).

**Table 4. Real-time RT-PCR protocol for simultaneous assays**

Temperature	Time	Number of cycles
50°C	20 min	1
95°C	15 min	1
95°C	30 s	40
57°C†	45 s	
68°C	45 s	

\* Fluorescence data collection.

## **Protocol: Real-time RT-PCR using 96-well plate real-time cycler**

Please read “Important Notes”, page 11 and “Important points before starting” and “Things to do before starting”, page 13.

### **Procedure**

- 1. Pipet 5  $\mu$ l of RNA samples, Positive Control, and Negative Control into individual reaction tubes. Cover the reaction tubes (e.g., with PCR sealing film).**

Include positive and negative control reactions.

Positive control: Use 5  $\mu$ l of the positive control (Positive Control) instead of sample RNA.

Negative control: Use 5  $\mu$ l of the negative control (Negative Control) instead of sample RNA.

- 2. Denature the samples for 5 min at 98°C in a 96-well plate standard cycler with a heated lid.**
- 3. Immediately cool down on ice water or liquid nitrogen for at least 20 s. Then store on ice or in a cooling device.**
- 4. Pipet 20  $\mu$ l of the Master Mix into each reaction tube. Thus the final volume of a test is 25  $\mu$ l (Table 5).**

**Table 5. Preparation of reaction mix**

<b>Component</b>	<b>Volume</b>
Master Mix	20 $\mu$ l
Sample	5 $\mu$ l
<b>Total volume</b>	<b>25 <math>\mu</math>l</b>

5. Close the reaction tubes with the corresponding caps.
6. Set the filters for the reporter dyes in the software of your thermal cycler according to Table 6.

**Table 6. Filter settings for reporter**

<b>Pathogen/internal control</b>	<b>Reporter</b>
BTV	FAM
BTV-8	Cy5
Internal Control	HEX/JOE*
Passive reference <sup>†</sup>	ROX

\* Use the option appropriate for your thermal cycler

<sup>†</sup> Internal reference for use with the Applied Biosystems<sup>®</sup> ABI PRISM<sup>®</sup> Sequence Detection Systems.

7. Run the real-time RT-PCR protocol according to Table 7 if running only the *virotype* BTV/pan8 RT-PCR Kit.

**Table 7. Real-time RT-PCR protocol for BTV pan/8**

Temperature	Time	Number of cycles
50°C	10 min	1
95°C	10 min	1
95°C	15 s	40
60°C <sup>†</sup>	60 s	

<sup>†</sup> Fluorescence data collection.

8. Run the real-time RT-PCR protocol according to Table 8 if running other *virotype* assays simultaneously (i.e., *virotype* PRRSV, *virotype* BVDV, *virotype* CSFV, *virotype* SBV and/or *virotype* Influenza A).

**Table 8. Real-time RT-PCR protocol for simultaneous assays**

Temperature	Time	Number of cycles
50°C	20 min	1
95°C	15 min	1
95°C	30 s	40
57°C*	45 s	
68°C	45 s	

\* Fluorescence data collection.

## Data Analysis and Interpretation

### Interpretation of results

For the assay to be valid the Positive Control must give a signal in the FAM, Cy5, and HEX channels with a  $C_T^* < 35$ . If no FAM and no Cy5 signals of the Positive Control are measured the denaturation and cooling steps were insufficient and the testing should be repeated. The Negative Control must give no signal.

The following results are possible if working with unknown samples. The possible sample results are also summarised in Table 9 on page 21.

**The sample is positive for BTV and BTV-8, and the assay is valid, if the following criteria are met:**

- The sample yields a signal in the FAM, Cy5 and HEX<sup>†</sup> channels
- The Positive Control yields a signal in all channels
- The Negative Control yields no signal

Note that very high concentrations of BTV RNA in the sample may lead to a reduced HEX signal or no HEX signal due to competition with the internal control.

\* Threshold cycle ( $C_T$ ) – cycle at which the amplification plot crosses the threshold, i.e., there is the first clearly detectable increase in fluorescence.

<sup>†</sup> Green, red, and yellow on the Rotor-Gene Q.

**The sample is positive for BTV, negative for BTV-8, and the assay is valid, if the following criteria are met:**

- The sample yields a signal in the FAM and HEX channels but not in the Cy5 channel
- The Positive Control yields a signal in all channels
- The Negative Control yields no signal

Note that very high concentrations of BTV RNA in the sample may lead to a reduced HEX signal or no HEX signal due to competition with the internal control.

**The sample is negative for both BTV and BTV-8, and the assay is valid, if the following criteria are met:**

- The sample yields a signal only in the HEX channel
- The Positive Control yields a signal in all channels
- The Negative Control yields no signal

A positive HEX signal rules out the possibility of PCR inhibition and/or incorrect RNA extraction as the internal control is amplified.

**The sample results are inconclusive, and the assay is invalid, if the following occurs:**

- The sample yields no signal in any of the fluorescence channels

The PCR was inhibited or the sample extraction was incorrect. It is recommended to retest the respective individual samples in nuclease free water (e.g., diluted 1:5), to repeat the RNA

extraction, or repeat the whole test procedure starting with new sample material.

Check that there is a fluorescence signal in the all channels for the positive control reaction (Positive Control). Absence of a signal for the Positive Control indicates an error, which could be due to incorrect RNA denaturation, RNA extraction failure,, or incorrect cycling conditions.

Repeat RNA extraction or repeat the whole procedure starting with new sample material.

**Table 9. Results interpretation table\***

Fluorescence signal	Pathogen			
	BTV	BTV-8	Negative	Invalid
FAM	X	X		
Cy5		X		
HEX	(X)	(X)	X	

\* Interpretation of sample results can be determined provided positive and negative control reactions are performed. The positive control must yield a signal in the FAM, Cy5, and HEX channels. The negative control must yield no signal. For a complete explanation of possible sample results please refer to "Data Analysis and Interpretation" on page 19.

## Troubleshooting Guide

The scientists in QIAGEN Technical Services are always happy to answer any questions you may have about either the information and/or protocols in this handbook or sample and assay technologies (for contact information, see back cover or visit [www.qiagen.com](http://www.qiagen.com)).

## Ordering Information

Product	Contents	Cat. no.
<i>virotype</i> BTV pan/8 RT-PCR Kit (24)	For 24 reactions: Master Mix, Positive Control, Negative Control	280443
<i>virotype</i> BTV pan/8 RT-PCR Kit (96)	For 96 reactions: Master Mix, Positive Control, Negative Control	280445
<i>virotype</i> BTV pan/8 RT-PCR Kit (480)*	For 480 reactions: Master Mix, Positive Control, Negative Control	280447
<b>Related products</b>		
<i>bactotype</i> MAP PCR Kit (24) <sup>†</sup>	For 24 reactions: Master Mix, Internal Control DNA, Positive Control, Negative Control	285903
<i>virotype</i> ASFV PCR Kit (96)	For 96 reactions: Master Mix, Positive Control, Negative Control	281905
<i>virotype</i> BTV RT-PCR Kit (96) <sup>†</sup>	For 96 reactions: PCR Mix, Enzyme Mix, Positive Control, Negative Control	280435
<i>virotype</i> PRRSV RT-PCR Kit (96) <sup>†</sup>	For 96 reactions: Master Mix, Positive Control, Negative Control	282305

\* Available only on request.

<sup>†</sup> Other kit sizes are available; see [www.qiagen.com](http://www.qiagen.com).

<b>Product</b>	<b>Contents</b>	<b>Cat. no.</b>
<i>virotype</i> BVDV RT-PCR Kit (96)*	For 96 reactions: PCR Mix, Enzyme Mix, Positive Control, Negative Control	280375
<i>virotype</i> CSFV RT-PCR Kit (96)*	For 96 reactions: Master Mix, Positive Control, Negative Control	281805
<i>virotype</i> SBV RT-PCR Kit (96)*	For 96 reactions: Master Mix, Positive Control, Negative Control	281605
<i>virotype</i> Influenza A RT-PCR Kit (96)*	For 96 reactions: Master Mix, Positive Control, Negative Control	282605
<i>bactotype</i> Mycoplasma Mg/Ms PCR Kit (96)*	For 96 reactions: Master Mix, Positive Control, Negative Control	288105
QIAamp <i>cador</i> Pathogen Mini Kit (50)*	For 50 preps: 50 QIAamp Mini Spin Columns, Carrier RNA, Proteinase K, Collection Tubes (2 ml), RNase-free Buffers	54104
QIAamp Viral RNA Mini Kit (50)*	For 50 RNA preps: 50 QIAamp Mini Spin Columns, carrier RNA, Collection Tubes (2 ml), RNase-free buffers	52904

\* Other kit sizes are available; see [www.qiagen.com](http://www.qiagen.com).

<b>Product</b>	<b>Contents</b>	<b>Cat. no.</b>
RNeasy Fibrous Tissue Mini Kit (50)	For 50 preps: 50 RNeasy Mini Spin Columns, Collection Tubes (1.5 ml and 2 ml), Proteinase K, RNase-free DNase I, RNase-free Reagents and Buffers	74704
RNeasy Mini Kit (50)*	For 50 preps: 50 RNeasy Mini Spin Columns, Collection Tubes (1.5 ml and 2 ml), RNase-free Reagents and Buffers	74104
Rotor-Gene Q 5plex Platform	Real-time PCR cycler with 5 channels (green, yellow, orange, red, crimson), laptop computer, software, accessories, 1-year warranty on parts and labor	9001570

\* Other kit sizes are available; see [www.qiagen.com](http://www.qiagen.com).

QIAGEN offer a range of ELISA kits and real-time PCR and real-time RT-PCR kits for the detection of animal pathogens. Visit [www.qiagen.com/Animal-and-Veterinary-Testing](http://www.qiagen.com/Animal-and-Veterinary-Testing) for more information about the *bactotype*<sup>®</sup>, *cador*<sup>®</sup>, *cattletype*<sup>®</sup>, *flocktype*<sup>®</sup>, *pigtype*<sup>®</sup>, and *virotype* products.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, QIAamp®, QIAcube®, *bactotype*®, *cador*®, *cattletype*®, *flocktype*®, *pigtype*®, RNeasy®, Rotor-Gene®, *virotype*® (QIAGEN Group); Applied Biosystems®, ABI PRISM®, FAM™, HEX™, JOE™, ROX™, TAMRA™ (Applied Biosystems Corporation or its subsidiaries); Eppendorf® (Eppendorf-Netheler-Hinz GmbH); Cy™ (GE Healthcare). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

The purchase of this product allows the purchaser to use it for amplification and detection of nucleic acid sequences for providing veterinary in vitro diagnostics. No general patent or other license of any kind other than this specific right of use from purchase is granted hereby.

#### Limited License Agreement for *virotype* BTV pan/8 RT-PCR Kit

Use of this product signifies the agreement of any purchaser or user of the product to the following terms:

1. The product may be used solely in accordance with the protocols provided with the product and this handbook and for use with components contained in the kit only. QIAGEN grants no license under any of its intellectual property to use or incorporate the enclosed components of this kit with any components not included within this kit except as described in the protocols provided with the product, this handbook, and additional protocols available at [www.qiagen.com](http://www.qiagen.com). Some of these additional protocols have been provided by QIAGEN users for QIAGEN users. These protocols have not been thoroughly tested or optimized by QIAGEN. QIAGEN neither guarantees them nor warrants that they do not infringe the rights of third-parties.
2. Other than expressly stated licenses, QIAGEN makes no warranty that this kit and/or its use(s) do not infringe the rights of third-parties.
3. This kit and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. QIAGEN specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the kit agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. QIAGEN may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the kit and/or its components.

For updated license terms, see [www.qiagen.com](http://www.qiagen.com).

© 2013–2014 QIAGEN, all rights reserved.

---

**www.qiagen.com**

**Australia** ■ techservice-au@qiagen.com

**Austria** ■ techservice-at@qiagen.com

**Belgium** ■ techservice-bnl@qiagen.com

**Brazil** ■ suportetecnico.brasil@qiagen.com

**China** ■ techservice-cn@qiagen.com

**Denmark** ■ techservice-nordic@qiagen.com

**Finland** ■ techservice-nordic@qiagen.com

**France** ■ techservice-fr@qiagen.com

**Germany** ■ techservice-de@qiagen.com

**Hong Kong** ■ techservice-hk@qiagen.com

**India** ■ techservice-india@qiagen.com

**Ireland** ■ techservice-uk@qiagen.com

**Italy** ■ techservice-it@qiagen.com

**Japan** ■ techservice-jp@qiagen.com

**Korea (South)** ■ techservice-kr@qiagen.com

**Luxembourg** ■ techservice-bnl@qiagen.com

**Mexico** ■ techservice-mx@qiagen.com

**The Netherlands** ■ techservice-bnl@qiagen.com

**Norway** ■ techservice-nordic@qiagen.com

**Singapore** ■ techservice-sg@qiagen.com

**Sweden** ■ techservice-nordic@qiagen.com

**Switzerland** ■ techservice-ch@qiagen.com

**UK** ■ techservice-uk@qiagen.com

