

# Qproteome<sup>®</sup> Mitochondria Isolation Kit

Lysis Buffer of the Qproteome Mitochondria Isolation Kit (cat. no. 37612) should be stored at  $-30$  to  $-15^{\circ}\text{C}$  upon arrival. All other buffers and Protease Inhibitor Solution (100x) should be stored at  $2$ – $8^{\circ}\text{C}$ .

## Further information

- *Qproteome Mitochondria Isolation Handbook*: [www.qiagen.com/HB-0870](http://www.qiagen.com/HB-0870)
- Safety Data Sheets: [www.qiagen.com/safety](http://www.qiagen.com/safety)
- Technical assistance: [support.qiagen.com](http://support.qiagen.com)

## Notes before starting

- This protocol is for standard preparation of mitochondria from eukaryotic cell lysates. For high-purity preparation and preparation from tissues, refer to the handbook.
  - All steps are performed at  $4^{\circ}\text{C}$ . Use pre-cooled buffers. Separated protein fractions should be snap-frozen in liquid nitrogen and stored at  $-80^{\circ}\text{C}$ .
  - For downstream applications (e.g., SDS-PAGE analysis), elution fractions should be pooled and concentrated (e.g., by acetone precipitation).
  - Immediately before use, supplement Lysis Buffer and Disruption Buffer with 1/100 volume of Protease Inhibitor Solution (100x) (i.e., if using 2 ml Disruption Buffer, add 20  $\mu\text{l}$  Protease Inhibitor Solution [100x]).
1. Transfer a cell suspension containing approximately  $5 \times 10^6$  to  $2 \times 10^7$  cells into a 15 ml conical tube and centrifuge at  $500 \times g$  for 10 min at  $4^{\circ}\text{C}$ . Carefully remove and discard the supernatant.
  2. Wash the cells using 1 ml 0.9% sodium chloride solution.

3. Resuspend the cell pellet in 1 ml ( $<1 \times 10^7$  cells) or 2 ml ( $\geq 1 \times 10^7$  cells) ice-cold Lysis Buffer by pipetting up and down using a 1 ml pipet tip. Incubate for 10 min at 4°C on an end-over-end shaker.
4. Centrifuge the lysate at 1000 x g for 10 min at 4°C.
5. Carefully remove the supernatant.  
This fraction primarily contains cytosolic proteins.
6. Resuspend the cell pellet in 1.5 ml ice-cold Disruption Buffer by pipetting up and down using a 1 ml pipet tip. Complete cell disruption by using a blunt-ended needle and a syringe (not provided). Draw the lysate slowly into the syringe and eject with one stroke. Repeat 10 times. Alternatively, disrupt cells using a Dounce or Potter homogenizer.
7. Centrifuge the lysate at 1000 x g for 10 min at 4°C and carefully transfer the supernatant to a clean 1.5 ml tube.  
The pellet contains nuclei, cell debris and unbroken cells. If desired, proteins can be re-extracted from the cell pellet by repeating steps 6 and 7 using 500 µl ice-cold Disruption Buffer. The supernatants from each extraction should be combined before the next step.
8. Centrifuge the supernatant(s) from step 7 at 6000 x g for 10 min at 4°C.
9. Carefully remove the supernatant.  
The pellet contains mitochondria. The supernatant constitutes the microsomal fraction.
10. Wash the mitochondrial pellet from step 9 with 1 ml Mitochondria Storage Buffer by carefully pipetting up and down using a 1 ml pipet tip. Centrifuge at 6000 x g for 20 min at 4°C.
11. Resuspend the mitochondrial pellet in Mitochondria Storage Buffer or a buffer of choice for further analysis.



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. Trademarks: QIAGEN®, Sample to Insight®, Qproteome® (QIAGEN Group). 1101320 02/2016 HB-1298-002 © 2016 QIAGEN, all rights reserved.