



Results

Average DNA yield for all 48 samples was $5.4 \mu\text{g} \pm 0.17 \mu\text{g}$, which corresponds to a deviation of approximately 3% (Figure 1). Average DNA purity was consistently high ($A_{260}/A_{280}=1.87$; S.D.=0.01). Agarose gel electrophoresis showed consistent high-quality DNA (Figure 2). In addition, clean, strong bands were observed for the single-copy PCR amplification (Figure 3).

Reproducible Yields




Figure 1. DNA yields from 48 blood samples from the same donor using the BioRobot M48 workstation.

Reproducible High-Quality DNA




Figure 2. High-quality genomic DNA purified from 48 blood samples using the BioRobot M48 workstation. 10 μl of each 200 μl eluate was analyzed by agarose gel electrophoresis.

Reproducible Performance in Sensitive PCR Analysis




Figure 3. Amplification of a single-copy gene (MECL-1) using template DNA purified from whole blood. Samples 1–6 and 43–48 are from the first and last sample prep plate positions on the BioRobot M48 workstation, respectively. All blood samples were obtained from the same donor. **M:** 100 bp DNA ladder (100 ng); **+**: positive control; **-**: negative control. 5 μl (10%) of each PCR was loaded on the agarose gel.



Conclusions

These results clearly demonstrate high reproducibility in both yield and quality of DNA purified using the BioRobot M48 workstation in combination with MagAttract magnetic particle technology.

Purification of DNA from whole blood using the BioRobot M48 workstation and MagAttract technology results in:

- **Highly reproducible yields** — average $5.4 \mu\text{g} \pm 0.17 \mu\text{g}$ DNA from 200 μl blood
- **High-quality DNA** — consistently strong, smear-free bands from all 48 samples
- **High-performance PCR** — strong, clean bands even from single-copy genes



Ordering Information

| Product | Description | Order No. |
|--|---|-----------|
| BioRobot M48 workstation | Robotic workstation for automation of magnetic-particle purification technology | 9000708 |
| MagAttract DNA Blood Mini M48 Kit (192)* | MagAttract Suspension and reagents for purification of genomic DNA from 200 µl whole blood samples using the BioRobot M48 workstation | 951336 |

The BioRobot M48 is intended as a microtiter diluting and dispensing device. No claim or representation is intended for its use in identifying any specific organism or for a specific clinical use (diagnostic, prognostic, therapeutic, or blood banking). It is the user's responsibility to validate the performance of the BioRobot M48 for any particular use, since its performance characteristics have not been validated for any specific organism. The BioRobot M48 may be used in clinical diagnostic laboratory systems after the laboratory has validated their complete system as required by CLIA '88 regulations in the U.S. or equivalents in other countries.

* MagAttract Kits are intended as general-purpose devices that may be used in clinical diagnostic laboratory systems after the laboratory has validated their complete system as required by CLIA '88 regulations in the U.S. or equivalents in other countries.

Contact QIAGEN today to discover more about standardized purification of DNA from blood.

QIAGEN Robotic Systems are not available in all countries; please inquire.

The PCR process is covered by U.S. Patents 4,683,195 and 4,683,202 and foreign equivalents owned by Hoffmann-La Roche AG.

Trademarks: BioRobot®, MagAttract®, QIAGEN®, (QIAGEN Group).

© 2003 QIAGEN, all rights reserved.

www.qiagen.com

QIAGEN:

Australia
Tel. 03-9489-3666
Fax 03-9489-3888

Canada
Tel. 800-572-9613
Fax 800-713-5951

France
Tel. 01-60-920-930
Fax 01-60-920-925

Germany
Tel. 02103-29-12400
Fax 02103-29-22000

Italy
Tel. 02-33430411
Fax 02-33430426

Japan
Tel. 03-5547-0811
Fax 03-5547-0818

Switzerland
Tel. 061-319-30-31
Fax 061-319-30-33

UK and Ireland
Tel. 01293-422-999
Fax 01293-422-922

USA
Tel. 800-426-8157
Fax 800-718-2056

Distributors:

Argentina Tecnolab S.A. (011) 4555 0010 **Austria** VWR International GmbH (01) 576 00 0 **Belgium/Luxemburg** Westburg b.v. 0800-1-9815 **Brazil** Uniscience do Brasil 011 3622 2320 **China** Gene Company Limited (852)2896-6283 **Cyprus** Scientronics Ltd (02) 357 22 765416 **Czech Republic** BIOCONSUIT spol. s.r.o. (420) 2 417 29 792 **Denmark** VWR International A/S 43 86 87 88 **Egypt** Clinilab 52 57 212 **Finland** VWR International Oy (09) 804 551 **Greece** BioAnalytica S.A. (210)640 03 18 **India** Genetix (011)542 1714 or (011)515 9346 **Israel** Westburg (Israel) Ltd. 08 6650813/4 or 1-800 20 22 20 **Korea** LRS Laboratories, Inc. (02) 924-86 97 **Malaysia** RESEARCH BIOLABS SDN. BHD. (603)-8070 3101 **Mexico** Quimica Valaner S.A. de C.V. (55) 55 25 57 25 **The Netherlands** Westburg b.v. (033)4950094 **New Zealand** Biolab Ltd (09) 980 6700 or 0800 933 966 **Norway** VWR International AS 22 90 00 00 **Poland** Syngen Biotech Sp.z.o.o. (071) 351 41 06 or 0601 70 60 07 **Portugal** IZASA PORTUGAL, LDA (21) 424 7312 **Singapore** Research Biolabs Pte Ltd 273 1066 **Slovak Republic** BIOCONSUIT Slovakia spol. s.r.o. (02) 5022 1336 **South Africa** Southern Cross Biotechnology (Pty) Ltd (021) 671 5166 **Spain** IZASA, S.A. (93) 902.20.30.90 **Sweden** VWR International AB (08) 621 34 00 **Taiwan** TAIGEN Bioscience Corporation (02) 2880 2913 **Thailand** Theera Trading Co. Ltd. (02) 412-5672 **In other countries contact:** QIAGEN, Germany

