

August 2011

Multi-Analyte ELISArray Handbook

For the simultaneous detection of up to 12
cytokines or chemokines in multiple samples



Sample & Assay Technologies

QIAGEN Sample and Assay Technologies

QIAGEN is the leading provider of innovative sample and assay technologies, enabling the isolation and detection of contents of any biological sample. Our advanced, high-quality products and services ensure success from sample to result.

QIAGEN sets standards in:

- Purification of DNA, RNA, and proteins
- Nucleic acid and protein assays
- microRNA research and RNAi
- Automation of sample and assay technologies

Our mission is to enable you to achieve outstanding success and breakthroughs. For more information, visit www.qiagen.com.

Product Use Limitations

Multi-Analyte ELISArray Kits are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

CONTENTS

I.	Background and Introduction	4	
II.	Materials Provided	6	
III.	Additional Materials Required	6	
IV.	Precautions	6	
V.	Complementary Products	6	
VI.	Protocol	7	
	A.	Sample Preparation and Handling	7
	B.	Reagent Preparation	8
	C.	Assay Procedure	9
	D.	Data Analysis	11
VII.	Troubleshooting and Frequently Asked Questions	11	

I. Background and Introduction

The Multi-Analyte ELISArray Kits are designed to simultaneously profile the level of multiple cytokines and/or chemokines using the conventional and simple sandwich-based enzyme-linked immunosorbent assay (ELISA) technique. The 96-well ELISA microplate has been coated with a panel of twelve target-specific capture antibodies, one in each eight-well strip allowing you to obtain qualitative relative profiling results from up to six samples. Figure 1 displays the typical layout of a catalogued ELISArray Kit. Each kit also includes the corresponding detection antibodies, Antigen Standards, and a complete set of detection reagents for a colorimetric ELISA.

The ELISArray Kits provide a rapid, simple and cost-effective solution for the assessment of protein expression changes for multiple protein targets in your experimental samples. We have screened commercially available antibodies to identify the capture and detection antibodies with the best sensitivity, best linearity and lowest background. All assays are optimized under uniform conditions to allow simultaneous detection with the same development or incubation time without compromising performance. Unlike other array-based technologies, the ELISArray Kits enable you to analyze multiple protein targets without the need for special equipment.

The Multi-Analyte ELISArray Kit uses a standard ELISA technique. Figure 2 displays an overview of the assay protocol. Incubation allows the capture antibodies to bind their specific protein of interest. After washing away unbound protein, biotinylated detection antibodies added to the wells also bind the captured analyte. After washing again to remove unbound material, an avidin-horseradish peroxidase conjugate is added. The wells are again washed and a colorimetric substrate solution is added, which produces a blue color in direct proportion to the amount of protein analyte present in the initial sample. The color development is stopped by adding stop solution, and the absorbance at 450 nm is read and compared across your samples.

Benefits of the Multi-Analyte ELISArray Kit:

- **Multi-Protein Flexibility:**
Detect up to twelve cytokines or chemokines at once using the same incubation and development times.
- **High Performance:**
Profile with a high level of sensitivity and linearity with screened and verified capture and detection antibodies.
- **Ease of Use:**
No special equipment is required, only the ELISA plate reader already in your lab.

		Cytokines											
		1	2	3	4	5	6	7	8	9	10	11	12
Negative Control	A												
Sample 1	B												
Sample 2	C												
Sample 3	D												
Sample 4	E												
Sample 5	F												
Sample 6	G												
Positive Control	H												

Figure 1: Layout of Cataloged Multi-Analyte ELISArray Kits

All wells of an eight-well strip are coated with the same capture antibody for the same cytokine or chemokine. A set of twelve strips in one ELISArray microplate therefore represents twelve cytokines or chemokines. See the information for the specific kit that you ordered for the list of proteins represented. The corresponding lettered wells of each row will characterize the same biological sample. Each ELISArray microplate therefore characterizes up to six biological samples plus positive and negative controls.

1. Prepare all reagents.
Set up experimental samples, positive controls, and negative control.
2. Add 50 µl assay buffer into each well of ELISArray plate.
Transfer 50 µl of samples and control samples into the appropriate wells of the ELISArray plate.
Incubate 2 hours.
3. Wash three times.
4. Add 100 µl Detection Antibody Solution. Incubate 1 hour.
5. Wash three times.
6. Add 100 µl Avidin – HRP. Incubate 30 minutes.
7. Wash four times.
8. Add 100 µl Development Solution. Incubate 15 minutes in the dark.
9. Add 100 µl Stop Solution. Read OD 450 nm within 30 minutes.

**Figure 2: Overview of the Multi-Analyte ELISArray Kit procedure.
Total Time ~ 4.5 hours**

II. Materials Provided:

Component / Description	Quantity
BOX 1: Shipped on dry ice or blue ice packs. Store at -20 °C.	
Avidin-HRP Conjugate	1.5 ml tube
10% BSA	15 ml bottle
Donkey Serum	15 ml bottle
BOX 2: Shipped at ambient temperature. Store at 4 °C.	
96-well pre-coated Capture Antibody microplate	One plate carrying 12 x 8-well strips
Detection Antibody Dilution Tube Strip	One strip of 12 tubes
Sample Dilution Buffer Stock	60 ml bottle
Assay Buffer Stock	60 ml bottle
Wash Buffer (10x Concentrate)	125 ml bottle
Development Solution	25 ml bottle
Stop Solution	60 ml bottle
BOX 3: Shipped on dry ice or blue ice packs. Store at -20°C.	
Antigen Standards (1 µg/mL)	12 x 1.5 ml tubes
Detection Antibodies	12 x 1.5 ml tubes

Storage Conditions:

Box 1 is shipped on dry ice or blue ice packs and should be stored in **non-frost-free freezer** at -20°C upon receipt.

Box 2 is shipped at ambient temperature and should be stored at 4°C upon receipt.

Box 3 is shipped on dry ice or blue ice packs and should be stored in **non-frost-free freezer** at -20°C upon receipt.

Do not use kit beyond the expiration date printed on the label.

III. Additional Materials Required:

1. Standard ELISA Microplate Reader
Capable of measuring 450 nm absorbance with a 570 nm correction wavelength
2. Calibrated Multi-Channel Pipettor
3. Wash Bottle
4. Microcentrifuge Tubes
5. Laboratory Timer
6. Culture Tubes

IV. Precautions

1. The Development Solution is toxic if inhaled or swallowed. Avoid contact with skin. Keep container tightly closed when not in use.
2. Stop Solution is an acidic solution. Wear eye, hand, face, and clothing protection when using this material.

V. Complementary Products:

Single-Analyte ELISArray Kits

VI. Protocol:

Please read through this entire protocol before beginning your experiment.

A. Sample Preparation and Handling:

High quality input material is **ESSENTIAL** for obtaining good results.

The most important prerequisite for any ELISA analysis experiment is consistent, high-quality experimental samples. Therefore, the sample preparation and handling procedures are critical to the success of the experiment. Specimens should be clear, because residual traces of particulate matter, heme, lipids or other contaminants will interfere with the performance of the ELISA.

1. Recommended Preparation Methods:

You will need roughly 650 μl of each sample to complete the experiment.

High quality protein samples for the ELISA experiment must be prepared using one of the following methods, each specific for your biological sample:

a. Cell Culture Supernatants:

Remove any particulate material by centrifugation for 10 minutes at 1000 x g and assay immediately, or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze / thaw cycles.

b. Serum:

Use a serum separator tube and allow samples to clot for 30 minutes before centrifugation for 15 minutes at 1000 x g . Remove serum and assay immediately or store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze / thaw cycles.

c. Plasma:

Collect plasma using citrate, EDTA, or heparin as an anticoagulant. Within 30 minutes of collection, centrifuge for 10 minutes at 1000 x g . Assay immediately or store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze / thaw cycles.

B. Reagent Preparation:

Use only polypropylene tubes for the antigen and detection antibody steps.

Briefly centrifuge any 1.5 ml tube in the kit before opening it in order to collect its contents at the bottom of the tube.

1. ELISArray Kit Reagents:

- a. Thaw the 10% BSA and Donkey Serum at room temperature. Store on ice.
- b. Bring the Wash Buffer Concentrate, Assay Buffer Stock, Sample Dilution Buffer Stock and ELISArray plate to room temperature. Do not remove the ELISArray plate from its pouch at this time.
- c. Once finished with each reagent, return it to its proper storage conditions.
- d. Keep all other kit components at their recommended storage temperature until recommended by the protocol.

2. Wash Buffer

Visually inspect the Wash Buffer Concentrate to ensure that all components are in solution. If any precipitates are visible, briefly shake the bottle to suspend any precipitate. Dilute 50 ml of Wash Buffer Concentrate into de-ionized or distilled water (dH₂O) to a final volume of 500 ml. Transfer to wash bottle. Keep at room temperature.

3. Assay Buffer

Dilute 0.6 ml of 10% BSA into a final volume of 30 ml with Assay Buffer Stock. Keep at room temperature.

4. Sample Dilution Buffer:

For cell culture supernatant samples, prepare Sample Dilution Buffer 1:

Dilute 2 ml of 10% BSA into a final volume of 20 ml with Sample Dilution Buffer Stock. Store on ice.

For serum or plasma samples, prepare Sample Dilution Buffer 2:

Dilute 6 ml of Donkey Serum into a final volume of 20 ml with Sample Dilution Buffer Stock. Store on ice.

5. Sample Preparation

If you wish to characterize dilutions of your samples to insure that the results are in the linear dynamic range of the assay, perform dilutions of your samples using the appropriate Sample Dilution Buffer.

C. Assay Procedure:

NOTE: *Very accurate and precise pipetting is critical to the reliability and consistency of any ELISA experiment. Make sure that all of your single- and multi-channel pipettors are calibrated before starting this protocol. Maintain incubation times and temperatures as consistent as possible across arrays and experiments for best results.*

1. Generate an Antigen Standard Cocktail:
 - a. Just before use, thaw the Antigen Standards on ice for 20 minutes. Vortex gently to mix well.
 - b. In a single tube, prepare the Concentrated Antigen Standard Cocktail containing all 12 of the Antigen Standards by pipeting 10 μ l of each Antigen Standard into the same 880 μ l volume of the appropriate Sample Dilution Buffer to yield 1 ml of a Concentrated Antigen Standard Cocktail.

NOTE: *For Mix-N-Match ELISArrays see your Product Specification Sheet for Generating Antigen Standard Cocktail*

- c. Dilute 200 μ l of the concentrated Standard Antigen Cocktail into 800 μ l of the appropriate Sample Dilution Buffer for the final Antigen Standard Cocktail.
2. Remove the ELISArray plate from its pouch. Using a multi-channel pipettor, add 50 μ l of Assay Buffer into each well of the ELISArray plate.

NOTE: *For best results, add all samples and reagents carefully by touching the pipette tips to the side of the wells just above the bottom and allowing the dispensed volume to run down the well wall. Do not blow out the last drop of volume dispensed to avoid introducing bubbles into the wells.*

3. Add 50 μ l of the appropriate Sample Dilution Buffer into each well of Row A in the ELISArray plate to setup the negative control.

NOTE: *Change your pipette tips with **every** addition to avoid cross-contamination of the wells.*

4. Add 50 μ l of each experimental sample or their dilutions to each well of their respective rows (B through G) in the ELISArray plate.
5. Add 50 μ l of the final Antigen Standard Cocktail into each well of Row H in the ELISArray plate to set up the positive control.

6. Cover the plate and gently shake or tap the plate for 10 seconds to mix. Incubate for 2 hours at room temperature on your bench-top.

7. Allow the Development Solution and Stop Solution to warm to room temperature in lab bench drawer protected from light.

NOTE: *Avoid prolonged exposure of the Development Solution to light or contact with water, air, or extreme temperature.*

8. Dilute the Detection Antibodies just before the next step:
 - a. Thaw the Detection Antibodies on ice for 30 minutes.
 - b. Add 855 μ l of Assay Buffer to each tube of Detection Antibody.
 - c. Mix well but gently.
 - d. Transfer each detection antibody to its own empty tube of the Detection Antibody Dilution Tube Strip.

- e. Be sure to add the Detection Antibodies to the tubes in the same left-to-right order as the capture antibodies in columns 1 through 12. See the Product Specification Sheet for the correct orientation.
9. Washing the ELISArray Plate:
 - a. Decant or aspirate the contents of the wells.
 - b. Wash wells by filling each well with 1x Wash Buffer (350 μ l per well).
 - c. Gently shake or tap the plate for 10 seconds to mix and follow by decanting or aspirating the solution.
 - d. Blot the plate upside down on absorbent paper to remove any residual buffer.
 - e. Repeat twice for a total of three washes.
 10. Using a 12-channel pipettor, transfer 100 μ l of the dilute Detection Antibodies from the Dilution Tube Strip to the appropriate rows of the ELISArray plate. Cover the plate and gently shake or tap the ELISArray plate for 10 seconds to mix. Incubate for 1 hour at room temperature on your bench-top.

NOTE: *Be sure to correctly orient your Detection Antibody Dilution Tube Strip above row A of the ELISArray plate to transfer the correct Detection Antibody to its correct wells.*
 11. Prepare Avidin-HRP:

Just before use, thaw the Avidin-HRP Conjugate on ice for 20 minutes. Vortex gently to mix well. Add 11 μ l of the Avidin-HRP Conjugate to 11 ml of Assay Buffer.

NOTE: *Avoid prolonged exposure of the Avidin-HRP Solution to light.*
 12. Wash ELISA wells as in step 9.
 13. Add 100 μ l of dilute Avidin-HRP into all wells. Cover the plate and gently shake or tap the ELISArray for 10 seconds to mix. Incubate for 30 minutes at room temperature in the dark.
 14. Wash ELISA wells as in step 9, except for a total of 4 washes.
 15. Just before use, transfer 12 ml of the Development Solution from its original bottle into a clean multi-channel pipettor reservoir.

NOTE: *Use a different clean multi-channel reservoir for each solution to be dispensed.*
 16. Add 100 μ l of the Development Solution to each well. Incubate the plate for 15 minutes at room temperature in the dark.
 17. Add 100 μ l of Stop Solution to each well in the same order as the Development Solution was added. The color in the wells should change from blue to yellow.
 18. Read absorbance at 450 nm within 30 minutes of stopping the reaction. If wavelength correction is available, subtract readings at 570 nm from the reading at 450 nm.

NOTE: *Subtracting the A450 reading by the A570 reading corrects for any minor optical imperfections in the ELISA plate. Uncorrected A450 readings may yield artificially high signals.*

D. Data Analysis:

1. Processing the Raw Data

- a. Typical absorbance values should range from 0.00 to 2.50.

Absorbance values greater than 2.50 are not within the linear range of the assay.

Also, do not interpret absorbance values less than two times the negative control absorbance values for each antigen.

- b. For each antigen, subtract the observed absorbance by the absorbance of the negative control to obtain the corrected absorbance values.

2. Profiling Experiment

- a. For each antigen separately, compare the corrected absorbance values between samples to determine which cytokine or chemokine changes its protein level most or least dramatically. Be sure to account for any dilution of the samples.
- b. Follow up the profiling experiment with a more detailed analysis using the corresponding Single Analyte ELISArray Kit.

VII. Troubleshooting and FAQs

For technical assistance and more information, please see our Technical Support Center at www.qiagen.com/Support or call one of the QIAGEN Technical Service Departments or local distributors (see back cover or visit www.qiagen.com).

Ordering Information

Product	Contents	Cat. no.
Multi-Analyte ELISArray Kits	Antigen Standards, Detection Antibodies, Avidin-HRP Conjugate, 10% BSA, Donkey Serum, 96-Well Precoated Capture Antibody Microplate, Detection Antibody Dilution Tube Strip, Sample Dilution Buffer Stock, Assay Buffer Stock, 10x Wash Buffer, Development Solution, Stop Solution	Varies

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor

Notes

Notes

Trademarks: QIAGEN® (QIAGEN Group).

LIMITED PRODUCT WARRANTY

This warranty limits our liability to replace this product in the event the product fails to perform due to any manufacturing defect. QIAGEN makes no other warranties of any kind, expressed or implied, including without limitation, warranties of merchantability or fitness for a particular purpose. QIAGEN shall not be liable for any direct, indirect, consequential or incidental damages arising out of the use, the results of use or the inability to use this product.

NOTICE TO PURCHASER

Purchase of an ELISArray Kit does not grant rights to use the kit components for reproduction of any primer pair mix, to modify kit components for resale or to use the ELISArray Kit to manufacture commercial products without written approval of QIAGEN. No other license, expressed, implied or by estoppels, is granted. U.S. patents may cover the use of certain antibodies included in the ELISArray Kit. Presently, it is not clear under U.S. laws whether commercial users must obtain licenses from the owners of the rights to these U.S. patents before using the ELISArray Kit.

Limited License Agreement

Use of these products signifies the agreement of any purchaser or user of the Multi-Analyte ELISArray Kits to the following terms:

1. The Multi-Analyte ELISArray Kits may be used solely in accordance with the *Multi-Analyte ELISArray Handbook* and for use with components contained in the Kits only. QIAGEN grants no license under any of its intellectual property to use or incorporate the enclosed components of these Kits with any components not included within these Kits except as described in the *Multi-Analyte ELISArray Handbook* and additional protocols available at www.qiagen.com.
2. Other than expressly stated licenses, QIAGEN makes no warranty that these Kits and/or their use(s) do not infringe the rights of third-parties.
3. These Kits and their components are licensed for one-time use and may not be reused, refurbished, or resold.
4. QIAGEN specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of these Kits agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. QIAGEN may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to these Kits and/or their components.

For updated license terms, see www.qiagen.com.

© 2011 QIAGEN, all rights reserved.

www.qiagen.com

Australia ■ Orders 1-800-243-800 ■ Fax 03-9840-9888 ■ Technical 1-800-243-066

Austria ■ Orders 0800-28-10-10 ■ Fax 0800-28-10-19 ■ Technical 0800-28-10-11

Belgium ■ Orders 0800-79612 ■ Fax 0800-79611 ■ Technical 0800-79556

Brazil ■ Orders 0800-557779 ■ Fax 55-11-5079-4001 ■ Technical 0800-557779

Canada ■ Orders 800-572-9613 ■ Fax 800-713-5951 ■ Technical 800-DNA-PREP (800-362-7737)

China ■ Orders 86-21-3865-3865 ■ Fax 86-21-3865-3965 ■ Technical 800-988-0325

Denmark ■ Orders 80-885945 ■ Fax 80-885944 ■ Technical 80-885942

Finland ■ Orders 0800-914416 ■ Fax 0800-914415 ■ Technical 0800-914413

France ■ Orders 01-60-920-926 ■ Fax 01-60-920-925 ■ Technical 01-60-920-930 ■ Offers 01-60-920-928

Germany ■ Orders 02103-29-12000 ■ Fax 02103-29-22000 ■ Technical 02103-29-12400

Hong Kong ■ Orders 800 933 965 ■ Fax 800 930 439 ■ Technical 800 930 425

Ireland ■ Orders 1800 555 049 ■ Fax 1800 555 048 ■ Technical 1800 555 061

Italy ■ Orders 800-789-544 ■ Fax 02-334304-826 ■ Technical 800-787980

Japan ■ Telephone 03-6890-7300 ■ Fax 03-5547-0818 ■ Technical 03-6890-7300

Korea (South) ■ Orders 080-000-7146 ■ Fax 02-2626-5703 ■ Technical 080-000-7145

Luxembourg ■ Orders 8002-2076 ■ Fax 8002-2073 ■ Technical 8002-2067

Mexico ■ Orders 01-800-7742-639 ■ Fax 01-800-1122-330 ■ Technical 01-800-7742-436

The Netherlands ■ Orders 0800-0229592 ■ Fax 0800-0229593 ■ Technical 0800-0229602

Norway ■ Orders 800-18859 ■ Fax 800-18817 ■ Technical 800-18712

Singapore ■ Orders 1800-742-4362 ■ Fax 65-6854-8184 ■ Technical 1800-742-4368

Spain ■ Orders 91-630-7050 ■ Fax 91-630-5145 ■ Technical 91-630-7050

Sweden ■ Orders 020-790282 ■ Fax 020-790582 ■ Technical 020-798328

Switzerland ■ Orders 055-254-22-11 ■ Fax 055-254-22-13 ■ Technical 055-254-22-12

UK ■ Orders 01293-422-911 ■ Fax 01293-422-922 ■ Technical 01293-422-999

USA ■ Orders 800-426-8157 ■ Fax 800-718-2056 ■ Technical 800-DNA-PREP (800-362-7737)

