

## Quantitative RT-PCR analysis of human monocytes using total RNA purified using the BioRobot® M48 workstation

Ovstebo R., Haug K.B.F., and Kierulf P.


Department of Clinical Chemistry, The R&D Group, Ullevaal University Hospital, Oslo, Norway.

This study describes the suitability of total RNA purified using the BioRobot® M48 workstation for sensitive RT-PCR applications. Numbers of gene transcripts correlated well between samples containing different numbers of cell equivalents, indicating highly pure RNA and consistent purification efficiency. Ratios of target transcripts relative to a housekeeping gene correlated well, independent of cell numbers.

### Introduction

Examining the differential transcription of genes in populations of cells helps decipher the complex mechanisms involved in cellular interactions, migration, and differentiation. To accurately measure gene expression in cell populations, it is important to use a reliable and standardized purification method that gives representative yields of RNA from a range of cell counts. Purified RNA must also be of high quality, and free from contaminating genomic DNA.

#### Linear Increase in 3 Different mRNA Transcripts With Cell Input



**Figure 1** Copy numbers of 3 different mRNA transcripts from LPS-stimulated human monocytes. Copy numbers from duplicate PCR correlate well with numbers of cells within the original lysate. **A**  $\beta$ -actin. **B** Tissue Factor. **C** MCP-1.




## Materials and methods

Human monocytes were extracted from whole blood by density gradient centrifugation, quantified by flow cytometry, and activated by incubation with LPS (1). Cells were lysed according to the MagAttract<sup>®</sup> RNA Cell Mini M48 Kit protocol and the lysate diluted to yield aliquots equivalent to  $1 \times 10^5$ ,  $2.5 \times 10^5$ , and  $5 \times 10^5$  human monocytes. Purification of total RNA from lysates (400  $\mu$ l) was fully automated using the BioRobot M48 workstation with the MagAttract RNA Cell Mini M48 Kit. DNase digestion was performed on the workstation and total RNA was eluted in 200  $\mu$ l.

cDNA was transcribed in 20  $\mu$ l reactions, combining 8  $\mu$ l RNA and 12  $\mu$ l RT mastermix. DNA was amplified from 2  $\mu$ l cDNA in a 20  $\mu$ l PCR using primers specific for MCP-1, Tissue Factor, and  $\beta$ -actin (2). PCR was performed in duplicate and control PCR was performed from sham RT reactions to ensure the absence of genomic DNA. Transcript numbers were read from transcript-specific standard curves (2).

### Consistent Analysis of Relative Gene Expression with Different Cell Inputs



**Figure 2** Ratios of mRNA copy numbers relative to  $\beta$ -actin mRNA in increasing numbers of LPS-stimulated human monocytes. Ratios from duplicate PCR and also from increasing numbers of cell equivalents, correlate well. **A** Tissue Factor:  $\beta$ -actin. **B** MCP-1:  $\beta$ -actin.

## Results

RT-PCR analyses showed accurate quantification of mRNA transcript numbers that accurately correlate to numbers of cell equivalents in the original sample (Figure 1). Transcript copy numbers were calculated from a calibration curve generated using defined standards of known copy-number (2). Ratios of MCP-1 and Tissue Factor transcripts relative to the housekeeping gene  $\beta$ -actin correlate well, independent of cell input (Figure 2).

## Conclusion

The BioRobot M48 workstation in combination with the MagAttract RNA Cell Mini M48 Kit provides standardized, fully automated purification of total RNA from cells.

- RNA is purified from a wide range of cell inputs with consistent efficiency and purified RNA performs reliably in sensitive RT-PCR analysis.
- Highly pure RNA enables reliable RT-PCR, which allows sensitive and accurate gene expression analyses.

## References

1. Lund, P.K., Joo, G.B., Westvik, A.B., Ovstebo, R., Kierulf, P. (2000) Isolation of monocytes from whole blood by density gradient centrifugation and counter-current elutriation followed by cryopreservation: six years' experience [In Process]. *Scand J Clin Lab Invest.* **60**: 357.
2. Ovstebo, R., Haug, K.B., Lande, K., Kierulf, P. (2003) PCR-based calibration curves for

## Ordering Information

Product	Contents	Cat. no.
MagAttract RNA Cell Mini M48 Kit (192)	MagAttract Suspension E, Buffers, RNase-free DNase,	958336
BioRobot M48 workstation	Robotic workstation for automation of magnetic particle technology	9000708
<b>Related Products</b>		
MagAttract RNA Tissue Mini M48 Kit (192)	MagAttract Suspension E, Buffers, RNase-free DNase,	959336
QuantiTect Probe RT-PCR Kit (200)	For 200 x 50 µl reactions: 3 x 1.7 ml QuantiTect Probe RT-PCR Master Mix (providing a final concentration of 4 mM MgCl <sub>2</sub> ), 100 µl QuantiTect RT Mix, 2 x 2 ml RNase-free water	204443
Omniscript RT Kit (50)	For 50 reverse-transcription reactions: 200 units Omniscript Reverse Transcriptase, 10x Buffer RT, dNTP Mix (containing 5 mM each dNTP), RNase-free water	205111

Visit us at [www.qiagen.com](http://www.qiagen.com) or contact **QIAGEN** today to discover how automated RNA preps enable sensitive and accurate gene expression analysis!

### Trademarks and disclaimers

QIAGEN®, BioRobot®, MagAttract®, Omniscript®, QuantiTect™ (QIAGEN Group).

QIAGEN Robotic Systems are not available in all countries; please inquire.

The BioRobot EZ1 and EZ1 Kits are intended as general-purpose devices. No claim or representation is intended for their use in identifying any specific organism or for a specific clinical use (diagnostic, prognostic, therapeutic, or blood banking). It is the user's responsibility to validate the performance of BioRobot workstations and QIAGEN Kits for any particular use, since their performance characteristics have not been validated for any specific organism. BioRobot workstations and QIAGEN Kits may be used in clinical diagnostic laboratory systems after the laboratory has validated their complete system as required by CLIA '88 regulations in the U.S. or equivalents in other countries.

The PCR process is covered by U.S. Patents 4,683,195 and 4,683,202 and foreign equivalents owned by Hoffmann-La Roche AG

Purchase of QIAGEN products for PCR containing Taq DNA Polymerase, HotStarTaq DNA Polymerase, or ProofStart DNA Polymerase is accompanied by a limited license to use them in the Polymerase Chain Reaction (PCR) process for research and development activities in conjunction with a thermal cycler whose use in the automated performance of the PCR process is covered by the up-front license fee, either by payment to Applied Biosystems or as purchased, i.e. an authorized thermal cycler. The PCR process is covered by U.S. Patents 4,683,195 and 4,683,202 and foreign equivalents owned by Hoffmann-La Roche AG. \*

The 5' nuclease process is covered by patents owned by Roche Molecular Systems, Inc. and F. Hoffmann-La Roche Ltd.

1025399 08/2003 © 2003 QIAGEN, all rights reserved.

**Australia** = QIAGEN Pty Ltd = Orders 03-9489-3666 = Fax 03-9489-3888 = Technical 1-800-243-066

**Canada** = QIAGEN Inc. = Orders 800-572-9613 = Fax 800-713-5951 = Technical 800-DNA-PREP (800-362-7737)

**France** = QIAGEN S.A. = Orders 01-60-920-920 = Fax 01-60-920-925 = Technical 01-60-920-930

**Germany** = QIAGEN GmbH = Orders 02103-29-12000 = Fax 02103-29-22000 = Technical 02103-29-12400

**Italy** = QIAGEN S.p.A. = Orders 02-33430411 = Fax 02-33430426 = Technical 02-33430414

**Japan** = QIAGEN K.K. = Telephone 03-5547-0811 = Fax 03-5547-0818 = Technical 03-5547-0811

**Switzerland** = QIAGEN AG = Orders 061-319-30-30 = Fax 061-319-30-33 = Technical 061-319-30-31

**United Kingdom** = QIAGEN Ltd. = Orders 01293-422-911 = Fax 01293-422-922 = Technical 01293-422-999

**USA** = QIAGEN Inc. = Orders 800-426-8157 = Fax 800-718-2056 = Technical 800-DNA-PREP (800-362-7737)

