

AdnaTest BreastCancerDetect and AdnaTest ER/PR-Detect

AdnaTest BreastCancerDetect Box 2 (Box 2 of cat. no. 396012), containing the AdnaTest PrimerMixes and AdnaTest Positive Controls, and AdnaTest ER/PR-Detect (cat. no. 396062), containing the AdnaTest PrimerMixes and AdnaTest Positive Controls, must be stored separately at -30 to -15°C . To prevent possible contamination and repeated temperature changes, aliquot the primer mix. The components must not be used beyond the expiration date.

Further information

- *AdnaTest BreastCancerSelect, BreastCancerDetect and ER/PR-Detect Handbook:* www.qiagen.com/HB-2453
- Safety Data Sheets: www.qiagen.com/safety
- Technical assistance: toll-free 00800-22-44-6000 or support.qiagen.com

Protocol

1. Equilibrate AdnaTest Lysis/Binding Buffer and RNA Purification Buffers A and B to room temperature and place Tris-HCL Buffer on ice.
2. Wash 20 μl Oligo(dT)₂₅ Beads per sample twice with 20 μl AdnaTest Lysis/Binding Buffer per sample.
3. Add 20 μl washed Oligo(dT)₂₅ Beads to each sample.
4. Incubate for 10 min at room temperature under tilting and rotation at approximately 5 rpm.
5. Place the reaction tube in the AdnaMag-S and remove supernatant.
6. Wash beads with 2 x 100 μl RNA Purification Buffer A.

Important: To avoid any loss of beads, rinse lid and tube wall thoroughly.

7. Resuspend beads in 100 μl RNA Purification Buffer B and transfer into a new 1.5 ml tube (provided).
8. Wash beads with 1 x 100 μl RNA Purification Buffer B.
9. Wash beads with 1 x 100 μl ice cold Tris-HCL Buffer.
10. Resuspend beads in 29.5 μl RNase-free water.
11. Incubate for 5 min at 50°C and place on ice for at least 2 min.



12. Continue with reverse transcription; see Table 1 and Table 2.

Table 1. Reverse transcription reaction setup

| Component | Volume |
|--|-------------------------------|
| RT master mix | |
| 10x Buffer RT | 4.0 μ l |
| dNTP Mix (5 mM each dNTP) | 4.0 μ l |
| RNase inhibitor, 40 U/ μ l (Promega) | 0.5 μ l |
| Sensiscript Reverse Transcriptase | 2.0 μ l |
| Template RNA* | |
| mRNA/bead complex or RNase-free water | 29.5 μ l |
| Total volume | 40.0 μl |

* As RT control add 29.5 μ l of RNase-free water instead of mRNA/bead complex. The volume of the mRNA/bead complex may vary slightly. Always use the total volume of this in the reverse transcription reaction.

Table 2. Reverse transcription program

| Step | Time | Temperature |
|-----------------------|----------|-------------|
| Reverse transcription | 60 min | 37°C |
| Denaturation | 5 min | 93°C |
| Cooling | ∞ | 4°C |

13. Continue with multiplex PCR for AdnaTest BreastDetect (Table 3 and Table 4) and/or AdnaTest ER/PR-Detect (Table 5 and Table 6) or store cDNA at -20°C for a maximum of 4 weeks.

Table 3. Preparation of the multiplex PCR (AdnaTest BreastDetect)

| Component | Volume |
|---|-----------------------------|
| Multiplex PCR master mix | |
| HotStarTaq Master Mix | 25 μ l |
| RNase-free water | 13 μ l |
| PrimerMix BreastDetect | 4 μ l |
| cDNA or RT control or Negative control (RNase-free water) or Positive Control (C+) each: | 8 μ l |
| Total volume | 50 μl |

Table 4. PCR cycling program (AdnaTest BreastDetect)

| | Temperature | Time |
|-----------------------------------|-------------|--------|
| Initial activation step | 95°C | 15 min |
| 3-step cycling (35 cycles) | | |
| Denaturation: | 94°C | 30 s |
| Annealing: | 60°C | 30 s |
| Extension: | 72°C | 60 s |
| Number of cycles: | 35 | |
| Final extension: | 72°C | 10 min |
| Cool-down: | 4°C | ∞ |

Table 5. Preparation of the multiplex PCR (AdnaTest ER/PR-Detect)

| Component | Volume |
|--|--------------|
| Multiplex PCR master mix | |
| HotStarTaq Master Mix | 25 µl |
| RNase-free water | 13 µl |
| AdnaTest PrimerMix ER/PR-Detect | 4 µl |
| cDNA or RT control or Negative control (RNase-free water) or AdnaTest Positive Control (C+) each: | 8 µl |
| Total volume | 50 µl |

Table 6. PCR cycling program

| | Temperature | Time |
|-----------------------------------|-------------|--------|
| Initial activation step | 95°C | 15 min |
| 3-step cycling (37 cycles) | | |
| Denaturation: | 94°C | 30 s |
| Annealing: | 60°C | 30 s |
| Extension: | 72°C | 30 s |
| Number of cycles: | 37 | |
| Final extension: | 72°C | 5 min |
| Cool-down: | 4°C | ∞ |

14. For fragment analysis, use an Agilent® 2100 Analyzer or alternative analysis system.

For evaluation of the results, please refer to the *AdnaTest BreastCancerSelect*, *BreastCancerDetect* and *ER/PR-Detect Handbook*.

IMPORTANT: If the protocol is not followed exactly, this may result in false-negative or false-positive results.



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN®, Sample to Insight® (QIAGEN Group); Agilent® (Agilent Technologies, Inc.). 1106857 03/2017 HB-2455-001
© 2017 QIAGEN, all rights reserved.