

Quick-Start Protocol

RNeasy[®] Midi Kit

March 2016

The RNeasy Midi Kit (cat. no. 75142) can be stored for at least 9 months at room temperature (15–25°C) if not otherwise stated on label.

Further information

- *RNeasy Midi/Maxi Handbook*: www.qiagen.com/HB-0436
- Safety Data Sheets: www.qiagen.com/safety
- Technical assistance: support.qiagen.com

Notes before starting

- If purifying RNA from cell lines rich in RNases, or from tissue, add either 10 µl β-mercaptoethanol (β-ME), or 20 µl 2 M dithiothreitol (DTT),* to 1 ml Buffer RLT before use. Buffer RLT containing β-ME or DTT can be stored at room temperature for up to 1 month.
- Add 4 volumes of ethanol (96–100%) to Buffer RPE for a working solution.
- Remove RNAlater[®] stabilized tissue from the reagent using forceps.
- For RNeasy Protect Midi Kit (cat. no. 75154), please start with the *Quick-Start Protocol RNAlater RNA Stabilization Reagent, RNAlater TissueProtect Tubes, and RNeasy Protect Kits*.

* This option not included in handbook; handbook to be updated.

1. **Cells:** Harvest a maximum of 1×10^8 cells. Add the appropriate volume of Buffer RLT and homogenize the lysate (see Table 1).

Tissues: Disrupt and homogenize the tissue (≤ 250 mg) in the appropriate volume of Buffer RLT (see Table 1). Centrifuge the lysate for 10 min at 3000–5000 x g. Remove the supernatant by pipetting and use it in step 2.

2. Add 1 volume of 70% ethanol to the lysate, and mix well by shaking vigorously. Do not centrifuge. Proceed immediately to step 3.

3. Transfer the sample (maximum 4 ml) to an RNeasy Midi column placed in a 15 ml centrifuge tube (supplied). Close the tube gently, and centrifuge for 5 min at 3000–5000 x g. Discard the flow-through.
4. Add 4 ml Buffer RW1 to the RNeasy spin column. Close the lid gently, and centrifuge for 5 min at 3000–5000 x g. Discard the flow-through.

Optional: For on-column DNase digestion, follow steps in Appendix F of the *RNeasy Midi/Maxi Handbook*.

5. Add 2.5 ml Buffer RPE to the RNeasy spin column. Close the lid gently, and centrifuge for 2 min at 3000–5000 x g. Discard the flow-through.
6. Add 2.5 ml Buffer RPE to the RNeasy spin column. Close the lid gently, and centrifuge for 5 min at 3000–5000 x g to dry the RNeasy silica membrane.
7. To elute, transfer the RNeasy column to a new 15 ml collection tube (supplied). Add the appropriate volume of RNase-free water (see Table 1) directly to the spin column membrane. Close the tube gently. Let it stand for 1 min, then centrifuge for 3 min at 3000–5000 x g.
8. Repeat step 7 as described with a second volume of RNase-free water.

Table 1. Volumes for sample homogenization and elution of RNA

Sample	Amount	Buffer RLT (ml)	Disruption and homogenization	Expected RNA yield and elution volume
Animal cells	5×10^6 – 3×10^7	2	Add Buffer RLT, use QIAshredder, TissueRuptor [®] ; or needle and syringe	≤150 µg in 150 µl
	3×10^7 – 5×10^7	2 or 4*		150 µg–1 mg in 250 µl
	5×10^7 – 1×10^8	4		
Animal tissues	20–75 mg	2	TissueLyser II; TissueRuptor; or mortar and pestle followed by QIAshredder or needle and syringe	≤150 µg in 150 µl
	75–130 mg	2 or 4†		150 µg–1 mg in 250 µl
	130–250 mg	4		

* If expected RNA yield is higher than 150 µg per 10^7 cells.

† If using difficult-to-lyse tissue samples >75 mg.



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. Trademarks: QIAGEN[®], Sample to Insight[®], RNeasy[®], TissueRuptor[®] (QIAGEN Group). *RNAiater[™] is a trademark of AMBION, Inc., Austin, Texas and is covered by various U.S. and foreign patents. 1101296 03/2016 HB-0569-002 © 2016 QIAGEN, all rights reserved.