

High-throughput QIAxcel[®] based subtyping: Enhanced epidemiologic surveillance of *Campylobacter jejuni*



Eduardo Taboada¹, Steven Mutschall¹, and Mirjana Kozulic²

¹ Campylobacter Genomics Laboratory, Division of Science and Technology, Laboratory for Foodborne Zoonoses, Public Health Agency of Canada, Lethbridge, Canada

² QIAGEN Instruments AG, Garstligweg 8, CH-8634 Hombrechtikon, Switzerland




Public Health Agency of Canada
Agence de la santé publique du Canada

C. jejuni — a significant public health concern

- *C. jejuni* is the leading cause of foodborne bacterial enteritis in Canada and worldwide.
- Approximately 90% of human foodborne infections are caused by *C. jejuni*.
- Campylobacteriosis causes symptoms such as fever, abdominal pain, nausea, vomiting, and diarrhea.
- In Canada, there are probably ~300,000 cases a year, costing in total ~\$300 million.

Challenges in *C. jejuni* subtyping


- Pulsed-field gel electrophoresis (PFGE) is a gold standard method for many bacterial foodborne pathogens, but not optimal for *C. jejuni* due to fingerprint "instability" (Barton et al., 2006).
- *flaA* typing through RFLP or sequencing of *flaA* alleles is rapid and inexpensive. However, high rates of recombination in *C. jejuni* can lead to erroneous genotypic clusters.
- Multi locus sequence typing (MLST) is a gold standard for *C. jejuni* subtyping and has several advantages, such as unambiguous and highly portable sequence data. However, MLST is a low-throughput method, lacks resolution in some epidemiological contexts, and is costly to implement.



Campylobacter jejuni.

Comparative genomic fingerprinting (CGF)


- Pan-genome = core genes + accessory genes of a species
- Accessory genes are responsible for many 'major phenotypic traits' of importance
- The accessory genome houses most of the genetic variability. Its diversity can be utilized for strain subtyping.
- The 'true' gold standard is the whole-genome phylogeny and a good subtyping method should compare favorably to the whole-genome phylogeny.



The pan-genome concept.



C. jejuni diversity and the CGF scheme

- The method is based on presence/absence of accessory genes.
- Each marker has intermediate carriage in the *C. jejuni* population = 'random' binary marker.
- Multi-marker profile: 40 markers → 2⁴⁰ possible fingerprints.
- SNP-free PCR primer design → multiplex PCR → 8 x multiplex PCRs = a comparative genomic fingerprint.



Genomic diversity in *C. jejuni*, containing 16 hypervariable regions with accessory genes.


CGF workflow



CGF workflow timeline. Complete CGF fingerprint analysis of 12 isolates in less than 5 hours, with minimal hands-on time.

QIAxcel Advanced improves CGF data analysis and workflow

Improved accuracy of CGF band sizing analysis: fewer falsely detected or misclassified bands.



As opposed to the standard gel electrophoresis, the QIAxcel platform enables automated streamlined analysis of CGF mPCRs, thereby facilitating high-throughput analysis and improving the discriminatory power.

Advantages:

- Automated analysis
- Automated data reporting
- Decreased time for sample-answer (~48 isolates per work day)
- Reduced costs (~\$7 USD per isolate)

Conclusion

- There are a number of challenges with existing methods of subtyping *Campylobacter jejuni*.
- CGF exploits significant differences in accessory genome content in *C. jejuni* and thereby has a greater discriminatory power than MLST.
- Adoption of the QIAxcel Advanced System has improved CGF data analysis and facilitated method assessment/deployment, resulting in:
 - Improved sensitivity and accuracy
 - A simplified, automated workflow
 - Automated data reporting
 - Increased throughput.
- A common method for subtyping has made it feasible to envision national networks for campylobacter surveillance, involving medical research and environmental settings, as well as the food production and agricultural industries.
- Comprehensive human, farm-to-fork, and source-to-tap surveillance (under a national database), will lead to:
 - Improved understanding of campylobacter epidemiology
 - The development and assessment of intervention strategies to mitigate risk.

The applications presented here are for molecular biology applications. They are not intended for the diagnosis, prevention or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN[®], QIAxcel[®] (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law. © 2014 QIAGEN, all rights reserved.