

# QuantiFast<sup>®</sup> Probe PCR Kit

The QuantiFast Probe PCR Kit (cat. nos. 204254, 204256 and 204257) should be stored immediately upon receipt at  $-30^{\circ}\text{C}$  to  $-15^{\circ}\text{C}$  and protected from light. 2x QuantiFast Probe PCR Master Mix can also be stored protected from light at  $2-8^{\circ}\text{C}$  for up to 2 months, depending on the expiration date.

## Further information

- *QuantiFast Probe PCR Handbook*: [www.qiagen.com/HB-0458](http://www.qiagen.com/HB-0458)
- Safety Data Sheets: [www.qiagen.com/safety](http://www.qiagen.com/safety)
- Technical assistance: [support.qiagen.com](http://support.qiagen.com)

## Notes before starting

- This protocol is optimized for quantification of gDNA and cDNA targets using dual-labeled probes (e.g., TaqMan<sup>®</sup> probes) on any cycler from Applied Biosystems except the Applied Biosystems<sup>®</sup> 7500 and ViiA<sup>™</sup> 7. PCR is carried out in the presence of ROX passive reference dye, which is included in 2x QuantiFast Probe PCR Master Mix.
- For the highest efficiency in real-time PCR using sequence-specific probes, targets should ideally be 70–200 bp in length.
- Always start with the cycling conditions specified in this protocol, even if using previously established primer–probe systems.
- The PCR must start with an initial incubation step of 3 min at  $95^{\circ}\text{C}$  to activate HotStarTaq<sup>®</sup> Plus DNA Polymerase.

1. Thaw 2x QuantiFast Probe PCR Master Mix, template gDNA or cDNA, primer and probe solutions and RNase-free water. Mix the individual solutions.
2. Prepare a reaction mix according to Table 1. Due to the hot start, it is not necessary to keep samples on ice during reaction setup or while programming the real-time cycler.

**Note:** We strongly recommend starting with the  $Mg^{2+}$  concentration as provided in 2x QuantiFast Probe PCR Master Mix.

3. Mix the reaction mix thoroughly, and dispense appropriate volumes into PCR tubes, PCR capillaries or the wells of a PCR plate.
4. Add template gDNA or cDNA to the individual PCR tubes, capillaries or wells containing the reaction mix.

**Note:** For two-step RT-PCR, the volume of the cDNA added (from the undiluted reverse-transcription reaction) should not exceed 10% of the final PCR volume.

**Table 1. Reaction setup**

Component	Volume/reaction		
	96-well block	384-well block	Final concentration
<b>Reaction mix</b> 2x QuantiFast Probe PCR Master Mix	12.5 µl	5 µl	1x
Primer A	Variable	Variable	0.4 µM
Primer B	Variable	Variable	0.4 µM
Probe	Variable	Variable	0.2 µM
RNase-free water	Variable	Variable	–
<b>Template gDNA or cDNA</b> (added at step 4)	Variable	Variable	≤200 ng/reaction
<b>Total reaction volume</b>	25 µl	10 µl	

5. Program the real-time cycler according to Table 2.

**Note:** Data acquisition should be performed during the combined annealing/extension step.

6. Place the PCR tubes, capillaries or plates in the real-time cycler, and start the cycling program.

**Table 2. Cycling conditions**

Step	Time	Temperature	Ramp rate
<b>PCR initial heat activation</b>	3 min	95°C	Maximal/fast mode
<b>2-step cycling:</b>			
Denaturation	3 s*	95°C	Maximal/fast mode
Combined annealing/extension	30 s*	60°C	Maximal/fast mode
Number of cycles	35–40†		

\* If your cycler does not accept this short time, choose the shortest time possible (e.g., 5 s denaturation for the ABI PRISM® 7700 or 31 s annealing/extension for the ABI PRISM 7000 or Applied Biosystems 7300).

† The number of cycles depends on the amount of template DNA.



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN®, Sample to Insight®, HotStarTaq®, QuantiFas® (QIAGEN Group); ABI PRISM®, Applied Biosystems®, ViiA™ (Applied Biosystems Corporation or its subsidiaries); TaqMan® (Roche Group). 1101262 03/2016 HB-0733-003 © 2016 QIAGEN, all rights reserved.