

QIAGEN Supplementary Protocol


ABI PRISM[®] 7000 software setup for the QuantiFast[®] Probe RT-PCR Kit

This protocol shows the necessary parameters that need to be entered into the ABI PRISM 7000 software (version 1.2 or higher) when using the QuantiFast Probe RT-PCR Kit.

IMPORTANT: Please read the *QuantiFast Probe RT-PCR Handbook*, paying careful attention to the “Safety Information” section, before beginning this procedure.

Procedure

1. Launch the software and start a thermal cycler protocol.
2. In the “Thermal Profile” tab, set up the parameters for RT-PCR as shown below.



3. Load your PCR plate and start the program.



QIAGEN handbooks can be requested from QIAGEN Technical Service or your local QIAGEN distributor. Selected handbooks can be downloaded from www.qiagen.com/literature.

Material safety data sheets (MSDS) for any QIAGEN product can be downloaded from www.qiagen.com/Support/MSDS.aspx.

The QuantiFast Probe RT-PCR Kit is intended for research use. No claim or representation is intended to provide information for the diagnosis, prevention, or treatment of a disease.

Trademarks: QIAGEN®, QuantiFast® (QIAGEN Group); ABI PRISM® (Applied Biosystems or its subsidiaries). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

The QuantiFast Probe RT-PCR Kit is an Authorized 5' Nuclease Core Kit without Licensed Probe. Its purchase price includes a limited, non-transferable immunity from suit under certain patents owned by Roche Molecular Systems, Inc. or F. Hoffmann-La Roche Ltd, for using only this amount of the product in the practice of the 5' nuclease process solely for the purchaser's own internal research when used in conjunction with Licensed Probe. No right under any other patent claims (such as apparatus or system claims) and no right to use this product for any other purpose is hereby granted expressly, by implication or by estoppel. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

Purchase of the QuantiFast Probe RT-PCR Kit is accompanied by a limited, non-transferable license under RT and Reverse Transcription-PCR patents owned by Roche Molecular Systems, Inc. and F. Hoffmann-La Roche Ltd to use it for the purchaser's own internal research. No real-time patent rights of any kind, no right under any other patent claims (such as apparatus or system claims), and no right to use this product for any other purpose is hereby granted expressly, by implication or by estoppel.

PCR94 Jan-09 © 2009 QIAGEN, all rights reserved.

www.qiagen.com

Canada = 800-572-9613

China = 0086-21-3865-3865

Denmark = 80-885945

Australia = 1-800-243-800

Finland = 0800-914416

Austria = 0800/281010

France = 01-60-920-930

Belgium = 0800-79612

Germany = 02103-29-12000

Brazil = 0800-557779

Hong Kong = 800 933 965

Ireland = 1800 555 049

Italy = 800-787980

Japan = 03-6890-7300

Korea (South) = 1544 7145

Luxembourg = 8002 2076

Mexico = 01-800-7742-639

The Netherlands = 0800 0229592

Norway = 800-18859

Singapore = 65-67775366

Spain = 91-630-7050

Sweden = 020-790282

Switzerland = 055-254-22-11

UK = 01293-422-911

USA = 800-426-8157

