

---

June 2020

# QIAseq<sup>®</sup> Multimodal Panel Handbook

Consolidated targeted next-generation  
sequencing of DNA and RNA

# Contents

|   |    |
|---|----|
| Kit Contents .....  | 4  |
| Storage .....   | 10 |
| Intended Use.....   | 10 |
| Safety Information.....   | 10 |
| Quality Control.....  | 11 |
| Introduction.....   | 12 |
| Principle and procedure .....   | 13 |
| Important Notes.....  | 20 |
| Total nucleic acid isolation .....  | 20 |
| Simultaneous purification of DNA and RNA into separate eluates .....        | 20 |
| Specific recommendations for FFPE samples (total nucleic acid or DNA) ..... | 21 |
| Nucleic acid input amount and sequencing depth .....                        | 21 |
| Sample multiplexing recommendations for Illumina sequencing platforms ..... | 23 |
| Equipment and Reagents to Be Supplied by User .....                         | 28 |
| Protocol: Nucleic Acid Fragmentation, Standard Samples.....                 | 29 |
| Protocol: Nucleic Acid Fragmentation, FFPE Samples.....                     | 32 |
| Protocol: RNA Polyadenylation.....  | 35 |
| Protocol: DNA Ligation.....   | 37 |
| Protocol: Reverse Transcription .....                                       | 40 |
| Protocol: Combined Targeted DNA+RNA Enrichment in a Single Tube .....       | 43 |
| Protocol: qPCR Determination of Universal PCR Cycles.....                   | 51 |
| Protocol: Universal PCR .....   | 54 |

---

|  |    |
|--|----|
| Recommendations: Library QC and Quantification .....                           | 62 |
| NGS library QC .....   | 62 |
| Preferred library quantification method .....                                  | 62 |
| Protocol: Sequencing Setup on Illumina MiSeq and NextSeq .....                 | 63 |
| Sequencing preparations for MiSeq.....   | 64 |
| Sequencing preparations for NextSeq .....                                      | 66 |
| Protocol: Data Analysis Using the QIAseq Multimodal Data Analysis Portal ..... | 68 |
| Troubleshooting Guide .....  | 71 |
| Appendix A: FFPE Sample Quality and Quantity .....                             | 72 |
| Appendix B: Separated Targeted DNA and RNA Enrichment in Separate Tubes .....  | 74 |
| Ordering Information .....   | 81 |
| Document Revision History.....   | 85 |

# Kit Contents

| <b>QIAseq Multimodal Panel</b>      | <b>(12)</b>   | <b>HC (12)</b> | <b>(96)</b>   | <b>HC (96)</b> | <b>Custom (96)</b> |
|-------------------------------------|---------------|----------------|---------------|----------------|--------------------|
| <b>Catalog no.</b>                  | <b>333932</b> | <b>333942</b>  | <b>333935</b> | <b>333945</b>  | <b>333955</b>      |
| <b>Number of samples</b>            | <b>12</b>     | <b>12</b>      | <b>96</b>     | <b>96</b>      | <b>96</b>          |
| Multimodal DHS Panel (DNA)          | 120 µl        | 120 µl         | 960 µl        | 960 µl         | 960 µl             |
| Multimodal VHS Panel (RNA)          | 96 µl         | 96 µl          | 768 µl        | 768 µl         | 768 µl             |
| Fragmentation Buffer, 10x           | 40 µl         | 40 µl          | 192 µl        | 192 µl         | 192 µl             |
| Fragmentation Enzyme Mix            | 90 µl         | 90 µl          | 384 µl        | 384 µl         | 384 µl             |
| FERA Solution                       | 15 µl         | 15 µl          | 60 µl         | 60 µl          | 60 µl              |
| Side Reaction Reducer               | 48 µl         | 48 µl          | 192 µl        | 192 µl         | 192 µl             |
| FG Solution                         | 170 µl        | 170 µl         | 170 µl        | 170 µl         | 170 µl             |
| Ligation Buffer, 5x                 | 160 µl        | 160 µl         | 1250 µl       | 1250 µl        | 1250 µl            |
| DNA Ligase                          | 75 µl         | 75 µl          | 600 µl        | 600 µl         | 600 µl             |
| Ligation Solution                   | 125 µl        | 125 µl         | 970 µl        | 970 µl         | 970 µl             |
| DNA Ligation Adapter                | 34 µl         | 34 µl          | 270 µl        | 270 µl         | 270 µl             |
| ATP Solution                        | 36 µl         | 36 µl          | 290 µl        | 290 µl         | 290 µl             |
| PAP Enzyme                          | 12 µl         | 12 µl          | 96 µl         | 96 µl          | 96 µl              |
| PAP Dilution Buffer, 10x            | 24 µl         | 24 µl          | 192 µl        | 192 µl         | 192 µl             |
| T4 Polynucleotide Kinase            | 12 µl         | 12 µl          | 125 µl        | 125 µl         | 125 µl             |
| EZ Reverse Transcriptase            | 15 µl         | 15 µl          | 150 µl        | 150 µl         | 150 µl             |
| RNase Inhibitor                     | 12 µl         | 12 µl          | 96 µl         | 96 µl          | 96 µl              |
| Multimodal RT Primer                | 12 µl         | 12 µl          | 96 µl         | 96 µl          | 96 µl              |
| Multimodal RT Buffer, 5x            | 60 µl         | 60 µl          | 480 µl        | 480 µl         | 480 µl             |
| Multimodal RT Enhancer              | 12 µl         | 12 µl          | 48 µl         | 48 µl          | 48 µl              |
| TEPCR Buffer, 5x                    | 60 µl x 2     | 60 µl x 2      | 500 µl x 2    | 500 µl x 2     | 500 µl x 2         |
| HotStarTaq® DNA Polymerase (6 U/µl) | 60 µl         | 60 µl          | 480 µl        | 480 µl         | 480 µl             |
| Nuclease-free Water                 | 1 tube        | 1 tube         | 10 ml         | 10 ml          | 10 ml              |
| QIAseq Beads                        | 10 ml         | 10 ml          | 38.4 ml x 2   | 38.4 ml x 2    | 38.4 ml x 2        |
| QIAseq Bead Binding Buffer          | 10.2 ml       | 10.2 ml        | 10.2 ml x 2   | 10.2 ml x 2    | 10.2 ml x 2        |

## Indexing for combined targeted DNA+RNA enrichment in a single tube

**Note:** Follow “Protocol: Combined Targeted DNA+RNA Enrichment in a Single Tube”, page 43.

|  |               |
|--|---------------|
| <b>QIAseq Multimodal Index I<br/>(12 DNA + 12 RNA sample indexes for 12 samples on Illumina® platform)</b>   | <b>(12)</b>   |
| <b>Catalog no.</b>   | <b>333962</b> |
| <b>Number of samples</b>   | <b>12</b>     |
| Multimodal N7 Plate (12):  | 1             |
| Each plate allows N7 indexing of 12 samples: 12 for DNA and 12 for RNA. Each well in the plate is single use. <b>In each indicated well of the cuttable plate, there are dried N7 index primers for both DNA and RNA, in the same well.</b> The plates can be cut in columns to enable indexing of the desired number of samples.  |               |
| Multimodal S5 Plate (12):  | 1             |
| Each plate allows S5 indexing of 12 samples: 12 for DNA and 12 for RNA. Each well in the plate is single use. SQDIB001 to SQDIB012 are mixed with universal DNA primer for DNA library amplification. SQDIB0049 to SQDIB060 are mixed with universal RNA primer for RNA library amplification. The plates can be cut in columns to enable indexing of the desired number of samples. |               |
| UPCR Buffer, 5x  | 60 µl x 3     |
| DNA qPCR AMP Set   | 12 µl         |
| RNA qPCR AMP Set   | 12 µl         |
| QIAseq A Read 1 Primer I (100 µM)  | 24 µl         |
| Multimodal Read 2 Primer (100 µM)  | 24 µl         |
| Optical Thin-wall 8-cap Strips   | 24 strips     |

|  |               |
|--|---------------|
| <b>QIAseq Multimodal Index I Set A</b><br><b>(Two identical sets of 48 DNA + 48 RNA sample indexes for Illumina platforms)</b>   | <b>(96)</b>   |
| <b>Catalog no.</b>   | <b>333965</b> |
| <b>Number of samples</b>   | <b>96</b>     |
| Multimodal N7 Plate Set A (48):  | 2             |
| Each plate allows N7 indexing of 48 samples: 48 for DNA and 48 for RNA. Each well in the plate is single use. <b>In each indicated well of the cuttable plate, there are dried N7 index primers for both DNA and RNA, in the same well.</b> The plates can be cut in columns to enable indexing of the desired number of samples.  |               |
| Multimodal S5 Plate Set A (48):  | 2             |
| Each plate allows S5 indexing of 48 samples: 48 for DNA and 48 for RNA. Each well in the plate is single use. SQDIB001 to SQDIB048 are mixed with universal DNA primer for DNA library amplification. SQDIB0049 to SQDIB096 are mixed with universal RNA primer for RNA library amplification. The plates can be cut in columns to enable indexing of the desired number of samples. |               |
| UPCR Buffer, 5x  | 500 µl x 3    |
| DNA qPCR AMP Set   | 96 µl         |
| RNA qPCR AMP Set   | 96 µl         |
| QIAseq A Read 1 Primer I (100 µM)  | 4 x 24 µl     |
| Multimodal Read 2 Primer (100 µM)  | 4 x 24 µl     |
| Optical Thin-wall 8-cap Strips   | 48 strips     |

|   |               |
|---|---------------|
| <b>QIAseq Multimodal Index I Set B</b><br><b>(Two identical sets of 48 DNA + 48 RNA sample indexes for Illumina platforms)</b>  | <b>(96)</b>   |
| <b>Catalog no.</b>  | <b>333975</b> |
| <b>Number of samples</b>  | <b>96</b>     |
| Multimodal N7 Plate Set B (48):   | 2             |
| <p>Each plate allows N7 indexing of 48 samples: 48 for DNA and 48 for RNA. Each well in the plate is single use. <b>In each indicated well of the cuttable plate, there are dried N7 index primers for both DNA and RNA, in the same well.</b> The plates can be cut in columns to enable indexing of the desired number of samples.</p>  |               |
| Multimodal S5 Plate Set B (48):   | 2             |
| <p>Each plate allows S5 indexing of 48 samples: 48 for DNA and 48 for RNA. Each well in the plate is single use. SQDIB097 to SQDIB144 are mixed with universal DNA primer for DNA library amplification. SQDIB0145 to SQDIB192 are mixed with universal RNA primer for RNA library amplification. The plates can be cut in columns to enable indexing of the desired number of samples.</p> |               |
| UPCR Buffer, 5x   | 500 µl x 3    |
| DNA qPCR AMP Set  | 96 µl         |
| RNA qPCR AMP Set  | 96 µl         |
| QIAseq A Read 1 Primer I (100 µM)   | 4 x 24 µl     |
| Multimodal Read 2 Primer (100 µM)   | 4 x 24 µl     |
| Optical Thin-wall 8-cap Strips  | 48 strips     |

## Indexing for separated targeted DNA and RNA enrichment in separate tubes

**Note:** Follow “Appendix B: Separated Targeted DNA and RNA Enrichment in Separate Tubes”, page 74.

| QIAseq Multimodal Index I Set A SW<br>(Two identical sets of 48 DNA + 48 RNA sample indexes for Illumina platforms)<br>Catalog no.<br>Number of samples   | (96)<br>333985<br>96 |
|---|----------------------|
| Multimodal N7 Plate Set A SW (48):<br><br>Each plate allows N7 indexing of 48 samples: 48 for DNA and 48 for RNA. Each well in the plate is single use. <b>In each indicated well of the cuttable plate, there are dried N7 index primers for DNA and RNA libraries, in separate wells.</b> The plates can be cut in columns to enable indexing of the desired number of samples.   | 2                    |
| Multimodal S5 Plate Set A (48):<br><br>Each plate allows S5 indexing of 48 samples: 48 for DNA and 48 for RNA. Each well in the plate is single use. SQDIB001 to SQDIB048 are mixed with universal DNA primer for DNA library amplification. SQDIB0049 to SQDIB096 are mixed with universal RNA primer for RNA library amplification. The plates can be cut in columns to enable indexing of the desired number of samples. | 2                    |
| UPCR Buffer, 5x   | 500 µl x 3           |
| DNA qPCR AMP Set  | 96 µl                |
| RNA qPCR AMP Set  | 96 µl                |
| QIAseq A Read 1 Primer I (100 µM)   | 4 x 24 µl            |
| Multimodal Read 2 Primer (100 µM)   | 4 x 24 µl            |
| Optical Thin-wall 8-cap Strips  | 48 strips            |

|  |               |
|--|---------------|
| <b>QIAseq Multimodal Index I Set B SW</b><br><b>(Two identical sets of 48 DNA + 48 RNA sample indexes for Illumina platforms)</b>  | <b>(96)</b>   |
| <b>Catalog no.</b>   | <b>333995</b> |
| <b>Number of samples</b>   | <b>96</b>     |
| Multimodal N7 Plate Set B SW (48):   | 2             |
| Each plate allows N7 indexing of 48 samples: 48 for DNA and 48 for RNA. Each well in the plate is single use. <b>In each indicated well of the cuttable plate, there are dried N7 index primers for DNA and RNA libraries, in separate wells.</b> The plates can be cut in columns to enable indexing of the desired number of samples.  |               |
| Multimodal S5 Plate Set B (48):  | 2             |
| Each plate allows S5 indexing of 48 samples: 48 for DNA and 48 for RNA. Each well in the plate is single use. SQDIB097 to SQDIB144 are mixed with universal DNA primer for DNA library amplification. SQDIB0145 to SQDIB192 are mixed with universal RNA primer for RNA library amplification. The plates can be cut in columns to enable indexing of the desired number of samples. |               |
| UPCR Buffer, 5x  | 500 µl x 3    |
| DNA qPCR AMP Set   | 96 µl         |
| RNA qPCR AMP Set   | 96 µl         |
| QIAseq A Read 1 Primer I (100 µM)  | 4 x 24 µl     |
| Multimodal Read 2 Primer (100 µM)  | 4 x 24 µl     |
| Optical Thin-wall 8-cap Strips   | 48 strips     |

---

## Storage

QIAseq Multimodal Panels (except Ligation Solution, QIAseq Beads, and QIAseq Bead Binding Buffer) are shipped on dry ice and should be stored immediately upon receipt at  $-30$  to  $-15^{\circ}\text{C}$  in a constant-temperature freezer. The Ligation Solution, QIAseq Beads, and QIAseq Bead Binding Buffer are shipped on cold packs. Upon receipt, QIAseq Beads and Bead Binding Buffer should be stored at  $2$ – $8^{\circ}\text{C}$ , while the Ligation Solution should be stored at  $-30$  to  $-15^{\circ}\text{C}$  in a constant-temperature freezer.

QIAseq Multimodal Index kits are shipped on dry ice and should be stored at  $-30$  to  $-15^{\circ}\text{C}$  upon arrival.

## Intended Use

QIAseq Multimodal Panels and QIAseq Multimodal Index kits are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

All due care and attention should be exercised in the handling of the products. We recommend all users of QIAGEN® products to adhere to the NIH guidelines that have been developed for recombinant DNA experiments, or to other applicable guidelines.

## Safety Information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at [www.qiagen.com/safety](http://www.qiagen.com/safety) where you can find, view, and print the SDS for each QIAGEN kit and kit component.

---

# Quality Control

In accordance with QIAGEN's ISO-certified Quality Management System, each lot of QIAseq Multimodal Panels is tested against predetermined specifications to ensure consistent product quality.

---

# Introduction

Recent advancements in NGS have enabled the analysis of single nucleotide variants (SNVs), InDels, and copy number variants (CNVs) from DNA, and the analysis of fusions and gene expression levels from RNA. Existing solutions, however, only allow users to perform such DNA and RNA analyses using 2 separate workflows: one for DNA and one for RNA. Additionally, these solutions require separate inputs of DNA and RNA, making the sequencing of low-yielding samples very difficult.

To overcome the limitations of existing solutions, the QIAseq Multimodal Panels have been developed. QIAseq Multimodal Panels enable Sample to Insight<sup>®</sup>, simultaneous targeted next-generation sequencing (NGS) of DNA and RNA using total nucleic acids in a single-tube workflow. Resulting DNA and RNA libraries can be sequenced together for cost effectiveness. This highly optimized solution facilitates ultrasensitive DNA variant detection as well as fusions and gene expression detection from RNA using integrated unique molecular indices (UMIs) from cells, tissue, and biofluids. The starting material for QIAseq Multimodal can be total nucleic acid or separately isolated DNA and RNA.

The QIAseq Multimodal Panels use a targeted approach to sequencing by enriching specific genomic or transcriptomic regions. This enhances DNA and RNA NGS by enabling users to sequence specific regions of interest, which in turn effectively increases sequencing depth and sample throughput while minimizing cost. More importantly, QIAseq Multimodal Panels enable simultaneous enrichment for both DNA and RNA regions of interest using total nucleic acids as input, thereby saving precious biological material, maximizing library prep efficiency, and reducing handling errors. Using a robust chemistry with integrated UMIs, QIAseq Multimodal Panels enable sensitive detection of DNA and RNA analytes of interest. Furthermore, data analysis tools have been developed to perform all steps necessary to generate both a DNA sequence variant report as well as an RNA fusion and gene expression report from NGS data. Collectively, QIAseq Multimodal Panels are a Sample to Insight solution for consolidated targeted DNA and RNA analysis using NGS (Figure 1).



**Figure 1. Overview of the Sample to Insight NGS workflow with QIAseq Multimodal Panels.** The complete Sample to Insight procedure begins with total nucleic acid or AllPrep® (separate DNA and RNA) extractions. Next is library construction and target enrichment with QIAseq Multimodal Panels. Following NGS, data analysis is performed using the QIAseq Multimodal Panel Analysis Software pipeline in Genomics Workbench. Ultimately, detected variants can be interpreted with QIAGEN Clinical Insight Interpret (QCI™-I) for QIAseq.

## Principle and procedure


QIAseq Multimodal Panels enable the simultaneous enrichment and library prep of DNA+RNA, with up to 28,000 primers per panel (20,000 DNA + 8,000 RNA). For DNA, the recommended input range is 10–40 ng for fresh samples or 40–250 ng for FFPE samples. For RNA, the recommended amount input range is 10–250 ng for fresh samples or ≤250 ng for FFPE samples (up to 500 ng for “severely” fragmented FFPE sample, with “severely” being defined as samples that have less than 40% of fragments >200 nt by smear analysis on the Bioanalyzer®). When working with total nucleic acid samples, input amounts should be based on DNA, because RNA is usually in vast excess to DNA. Lower input amounts are possible; however, this will lead to fewer sequenced UMIs and reduced variant detection sensitivity. The following reactions occur in a streamlined, single-tube workflow (Figure 2).

### Nucleic acid fragmentation

RNA molecules are heat fragmented and DNA molecules are enzymatically fragmented, end repaired, and A-tailed within a single controlled multienzyme reaction.

## RNA polyadenylation

Specific to RNA, synthetic polyadenylation is performed to create a binding site for subsequent reverse transcription.



**Figure 2. QIaseq Multimodal Panels workflow, using Combined Targeted DNA+RNA Enrichment. Alternatively, using Separated Targeted DNA and RNA Enrichment can be performed using "Appendix B: Separated Targeted DNA and RNA Enrichment in Separate Tubes".**

---

## DNA ligation

Specific to DNA, UMI-containing adapters are ligated at the 3' ends of the molecules. The UMI is a 12-base fully random sequence, which statistically provides  $4^{12}$  possible sequences per adapter and ensures that each molecule receives a UMI sequence. In addition, this adapter contains a binding site for subsequent target enrichment.

## RNA reverse transcription and template switching

Specific to RNA, reverse transcription and template switching are performed. For reverse transcription, the anchored oligo-dT primer contains a 10-base fully random UMI sequence, and the template switching oligonucleotide also contains a 10-base fully random UMI sequence. This allows each RNA molecule to be tagged with a unique UMI, regardless of which strand it was derived from. Lastly, the reverse transcription and template switching oligonucleotides each contain the same binding site for subsequent target enrichment.

## Target enrichment

Two protocols are provided for target enrichment:

- "Protocol: Combined Targeted DNA+RNA Enrichment in a Single Tube", page 43.
- "Appendix B: Separated Targeted DNA and RNA Enrichment in Separate Tubes", page 74.

**IMPORTANT:** Table 1 outlines important points to consider when choosing between the 2 protocols.

**Table 1. Target enrichment options**

| Choose ...  | For ...   | Notes   |
|---|---|---|
| <p>“Protocol: Combined Targeted DNA+RNA Enrichment in a Single Tube”</p>  | <p><b>Maximal detection sensitivity of DNA and RNA variants,</b> because there is no splitting of samples</p> | <p>Points to note with the combined protocol are that RNA primers will amplify DNA, and DNA primers have the possibility of amplifying RNA. As a result, when considering sequencing read budgets, the panel size needs to be accounted for as “DNA+RNA” primers. For example, if you are working with a DNA panel of 10,000 primers and an RNA panel of 1,000 primers, the read budget ultimately needed for both the DNA and the RNA library must be based on 11,000 primers.</p> |
| <p>“Appendix B: Separated Targeted DNA and RNA Enrichment in Separate Tubes”</p>  | <p><b>Maximal panel specificity</b></p>   | <p>With this protocol, when considering sequencing read budgets, the DNA panel and the RNA panel can be considered separately. For example, if you are working with a DNA panel of 10,000 primers and an RNA panel of 1,000 primers, the read budget ultimately needed for the DNA library must be based on 10,000 primers and the read budget ultimately needed for the RNA library must be based on 1,000 primers.</p>  |
| <p><b>Note: Samples are split before targeted enrichment, so the input DNA recommendations would be doubled, compared to combined targeted enrichment workflow.</b></p> | <p>Flexibility for Custom Multimodal Primer Panel design</p>  | <p>Separated targeted DNA and RNA enrichment prevents the need to account for potential interactions (i.e. dimerization potential) between DNA and RNA primer pool</p>  |
|   | <p>Custom Multimodal Primer Panels with specific primer numbers</p>   | <p>When the number of DNA+RNA primers (assuming the number DNA primers is greater than the number of RNA primers) is <math>\geq 14,000</math> and at the same time the number of RNA primers is <math>\geq 2500</math>, “Appendix B: Separated Targeted DNA and RNA Enrichment in Separate Tubes” workflow is highly recommended</p>  |
|   | <p><b>QIAseq Multimodal Panel UHS-5000Z</b></p>   | <p>QIAseq Multimodal Pan Cancer Panel (UHS-5000Z) is designed exclusively for the “Separated Targeted DNA and RNA Enrichment in Separate Tubes” workflow</p>  |

---

For both DNA and RNA, target enrichment is performed post-UMI assignment to ensure that molecules containing UMIs are sufficiently enriched in the sequenced library. For enrichment, ligated DNA molecules and reverse-transcribed/template-switched cDNA molecules are subject to several cycles of targeted PCR using a single primer extension (SPE) approach. This reaction includes highly optimized chemistry to amplify traditionally difficult regions using pools of DNA and RNA region-specific primers. Universal primers complementary to the DNA adapter-binding sequence and reverse-transcription/template-switching oligonucleotides ensure specificity for DNA and RNA molecules and assign the N7 index, which is one of the 2 sample unique dual indexes (UDIs). The N7 indexes are listed on [QIAGEN.com/QIAseqMultimodalPanels](https://www.qiagen.com/QIAseqMultimodalPanels).

### Library amplification

A Universal PCR is ultimately carried out separately on DNA and RNA libraries to both optimally amplify each library as well as add the second UDI. Collectively, DNA and RNA libraries for a given sample have their own unique dual indexes. The S5 indexes are listed on [QIAGEN.com/QIAseqMultimodalPanels](https://www.qiagen.com/QIAseqMultimodalPanels).

**Important:** The sample recommendations for the UDIs should be maintained between the N7 and S5 indexes.

Next-generation sequencing  
QIAseq Multimodal Panels are compatible with Illumina NGS platforms including MiniSeq®, MiSeq®, NextSeq® 500/550, HiSeq® 2500, HiSeq 3000/4000, and NovaSeq™ 6000. The QIAseq Multimodal Panels cannot be used on Illumina's iSeq 100 platform. When using Illumina NGS systems, QIAseq Multimodal libraries require a custom sequencing primer for Read 1 (QIAseq A Read 1 Primer I), custom sequencing primer for Read 2 (Multimodal Read 2 Primer), and 149 bp paired-end reads.

## Data analysis

Data from QIAseq Multimodal Panels can be analyzed using the QIAGEN CLC Genomics Workbench, which allows you to optimize analysis parameters for your specific panels. The parameters can then be locked for routine use. All detected variants can be further interpreted using QCI for QIAseq.


Alternatively, the QIAseq Multimodal Data Analysis pipeline is available at [ngsdataanalysis2.qiagen.com/MultiModal](https://ngsdataanalysis2.qiagen.com/MultiModal). The pipeline automatically performs all steps necessary to generate a DNA sequence variant report, as well as an RNA fusion and differential expression analysis report from your raw NGS data. An explanation of the principles of UMI-directed variant detection and the features of the primary sequence analysis output can be found and downloaded from [doi.org/10.1093/bioinformatics/bty790](https://doi.org/10.1093/bioinformatics/bty790).<sup>\*</sup> All detected variants can be further interpreted using QCI for QIAseq.

## Principle of variant detection with UMIs

The principle of variant detection with UMIs is described in Figure 3. Due to intrinsic noise and sequence-dependent bias, indexed molecules may be amplified unevenly across the target regions. Target region coverage can be better achieved, however, by counting the number of UMIs rather than counting the number of total reads for each region. Sequence reads having different UMIs represent different original molecules, while sequence reads having the same UMIs are the result of PCR duplication from one original molecule. Errors from PCR amplification and from the sequencing process may also be present in final reads that lead to false positive variants in sequencing results. These artifactual variants can be

<sup>\*</sup> Xu, C., Nezami Ranjbar, M.R., Wu, Z., DiCarlo, J., Wang, Y. (2017) Detecting very low allele fraction variants using targeted DNA sequencing and a novel molecular barcode-aware variant caller. *BMC Genomics*. **18**, 5.

greatly reduced by calling variants across all reads within a unique UMI instead of picking up variants at the original read level.



**Figure 3. Principle of variant detection with UMIs.** Each original molecule is tagged by a UMI. True variants are those mutations present in the majority of reads within a UMI, while false positives are mutations present in only one or a few reads within a UMI. Description of the variant calling algorithm can be found and downloaded from [doi.org/10.1186/s12864-016-3425-4](https://doi.org/10.1186/s12864-016-3425-4). \*

\* Xu, C., Nezami Ranjbar, M.R., Wu, Z., DiCarlo, J., Wang, Y. (2017) Detecting very low allele fraction variants using targeted DNA sequencing and a novel molecular barcode-aware variant caller. *BMC Genomics*. **18**, 5.

---

## Important Notes

The most important prerequisite for sequence analysis is consistent, high-quality DNA and RNA from every experimental sample. Therefore, sample handling and isolation procedures are critical to the success of the experiment. Residual traces of proteins, salts, or other contaminants may either degrade the nucleic acids or decrease the efficiency of – if not block completely – the enzymatic activity necessary for optimal targeted enrichment.

### Total nucleic acid isolation

Supplementary protocols for the simultaneous isolation of total nucleic acid (DNA+RNA) from cells and tissue, blood, or FFPE samples are available at [www.qiagen.com/TotalNucleicAcid](http://www.qiagen.com/TotalNucleicAcid).

**Note:** For quantification of DNA and RNA from total nucleic acid samples, we recommend the high-sensitivity Quant-iT™ dsDNA Assay Kit (Thermo Fisher Scientific, cat. no. Q33120) for DNA and the Quant-iT RNA Assay Kit (Thermo Fisher Scientific, cat. no. Q33140) for RNA. When working with total nucleic acid samples, input amounts should be based on DNA because RNA is usually in vast excess to DNA.

### Simultaneous purification of DNA and RNA into separate eluates

The QIAGEN kits listed in Table 2 are recommended for the preparation of DNA and RNA samples from cells, tissues, and FFPE tissues. For whole blood, we recommend the PAXgene Blood DNA Kit (cat. no. 761133) and the PAXgene Blood RNA Kit (cat. no. 762174).

**Note:** If samples must be harvested from biological samples for which kits are not available, please contact QIAGEN Technical Services ([support.qiagen.com](mailto:support.qiagen.com)) for suggestions.

**Table 2. Recommended AllPrep kits for simultaneous purification of DNA and RNA into separate eluates**

| Kit                      | Starting material | Cat. no. |
|--------------------------|-------------------|----------|
| AllPrep DNA/RNA Mini Kit | Cells and tissue  | 80204    |
| AllPrep DNA/RNA FFPE Kit | FFPE samples      | 80234    |

## Specific recommendations for FFPE samples (total nucleic acid or DNA)

If FFPE samples are used for QIAseq Multimodal Panels, the QIAseq DNA QuantiMIZE kits (cat. no. 333404 or 333414) are strongly recommended for determining the quality of each FFPE sample. Appendix A provides detailed information for FFPE DNA quality assessment and input amount.

## Nucleic acid input amount and sequencing depth

The number of UMIs captured from the original DNA sample correlates with the DNA input amount and sequencing read depth. Adequate sequencing of captured UMIs requires relatively deep sequencing coverage. Table 3 provides guidance on variant detection with fresh DNA amounts at different depths of coverage. Additionally, the number of UMIs sequenced directly impacts the variant detection sensitivity. Therefore, low-frequency mutation detection usually requires more DNA input and sequencing at deeper coverage (i.e., more reads/UMI) to generate a sufficient amount of UMIs.

**Table 3. Suggested fresh DNA input amount and read depth coverage for variant detection\***

| Variant frequency | Input (ng)      | Read pairs/UMI | Mean read depth (coverage) |
|-------------------|-----------------|----------------|----------------------------|
| 5%                | 10 <sup>†</sup> | 4              | 7200x                      |
| 5%                | 20 <sup>†</sup> | 2              | 3640x                      |
| 1%                | 40 <sup>†</sup> | 4              | 25,600x                    |

\* Variant detection is based on 90% sensitivity on the entire region of the QIAseq Multimodal DHS (DNA) Panel.

<sup>†</sup> If performing Appendix B: Separated Targeted DNA and RNA Enrichment in Separate Tubes, double the input amounts listed.

As RNA expression levels for each transcript varies broadly between different samples, there is no definitive calculation for a required number of reads. Based on previous testing with RNA libraries, allocation of 5000 reads per primer is a reasonable starting point. Table 4a (Combined Targeted DNA+RNA Enrichment) and Table 4b (Separated Targeted

DNA and RNA Enrichment) provide recommendations for the number of reads that should be allocated for the prepared QIAseq Multimodal DNA and RNA libraries. As a note, the QIAseq Multimodal Pan Cancer Panel (UHS-5000Z) is designed exclusively for the Separated Targeted DNA and RNA Enrichment workflow.

**Table 4a. Read allocation for cataloged QIAseq Multimodal DNA and RNA libraries (Combined Targeted DNA+RNA Enrichment)**

|                 |       | Primer Number | UHS-005Z                                 | UHS-009Z | UHS-006Z |
|-----------------|-------|---------------|--|----------|----------|
|                 |       | DNA primer #  | 4149                                     | 6244     | 11243    |
|                 |       | RNA primer #  | 487                                      | 1116     | 665      |
| Per DNA library | Input | Coverage (X)  | Reads allocated to each DNA library (M)* |          |          |
| 5% VAF          | 10 ng | 7200          | 33                                       | 53       | 86       |
| 5% VAF          | 20 ng | 3640          | 17                                       | 27       | 43       |
| 1% VAF          | 40 ng | 25600         | 119                                      | 188      | 305      |
|                 |       |               | Reads allocated to each RNA library (M)* |          |          |
| Per RNA library |       | 5000          | 23                                       | 37       | 60       |

\* Reads allocated to each library (M) = Coverage X Primer number/10<sup>6</sup>

**Table 4b. Read allocation for cataloged QIAseq Multimodal DNA and RNA libraries (Separated Targeted DNA and RNA Enrichment)**

|                 |                    | Primer Number | UHS-005Z   | UHS-009Z | UHS-006Z | UHS-5000Z* |
|-----------------|--------------------|---------------|--|----------|----------|------------|
|                 |                    | DNA primer #  | 4149   | 6244     | 11243    | 19995      |
|                 |                    | RNA primer #  | 487  | 1116     | 665      | 2571       |
| Per DNA library | Input              | Coverage (X)  | Reads allocated to each DNA library (M) <sup>†</sup> |          |          |            |
| 5% VAF          | 20 ng <sup>‡</sup> | 7200          | 30   | 45       | 81       | 144        |
| 5% VAF          | 40 ng <sup>‡</sup> | 3640          | 15   | 23       | 41       | 73         |
| 1% VAF          | 80 ng <sup>‡</sup> | 25600         | 106  | 160      | 288      | 512        |
|                 |                    |               | Reads allocated to each RNA library (M)*             |          |          |            |
| Per RNA library |                    | 5000          | 2  | 6        | 3        | 13         |

\* QIAseq Multimodal Pan Cancer Panel (UHS-5000Z) is designed exclusively for the “Separated Targeted DNA and RNA Enrichment” workflow.

<sup>†</sup> Reads allocated to each library (M) = Coverage X Primer number/10<sup>6</sup>

<sup>‡</sup> Since samples are split before targeted enrichment, the recommend DNA input is doubled, compared to combined targeted enrichment workflow.

## Sample multiplexing recommendations for Illumina sequencing platforms

Sample multiplexing level is determined by the size of the panel, required depth of coverage, and sequencing platform total output. For the Illumina platforms, sample indexes are available to multiplex up to 96 samples, one targeted DNA, and one targeted RNA library per sample, per run. General guidelines are provided for the number of samples that can be multiplexed in different sequencing platforms, based on panel size and read depth calculated for QIAseq Multimodal DNA and RNA libraries prepared using the Combined Targeted DNA+RNA Enrichment workflow (Table 5a through Table 5d) or the Separated Targeted DNA and RNA Enrichment workflow (Table 6a through Table 6d). Fine-tuning the read depth is possible after the first run. See read allocation and sample multiplexing template for the custom panel on [www.qiagen.com/PROM-16466](http://www.qiagen.com/PROM-16466).

**Table 5a. Number of multiplexed samples (one DNA and one RNA library for each sample) based on a DNA+RNA panel size (Combined targeted DNA+RNA Enrichment workflow) with 500x mean coverage for the DNA library and 5000 reads/primer for the RNA library\***

|                      |                       |                                     | UHS-005Z            | UHS-009Z    | UHS-006Z    |              |
|----------------------|-----------------------|-------------------------------------|---------------------|-------------|-------------|--------------|
|                      |                       |                                     | <b>DNA primer #</b> | <b>4149</b> | <b>6244</b> | <b>11243</b> |
|                      |                       |                                     | <b>RNA primer #</b> | <b>487</b>  | <b>1116</b> | <b>665</b>   |
| <b>Instrument</b>    | <b>Version</b>        | <b>Capacity (paired-ends reads)</b> |                     |             |             |              |
| MiniSeq              | Mid output            | 16M                                 | N/A                 | N/A         | N/A         |              |
| MiniSeq              | High output           | 50M                                 | 1                   | 1           | N/A         |              |
| MiSeq                | v2                    | 30M                                 | 1                   | N/A         | N/A         |              |
| NextSeq 500          | Mid output            | 260M                                | 10                  | 6           | 3           |              |
| NextSeq 500          | High output           | 800M                                | 31                  | 19          | 12          |              |
| HiSeq 2500 rapid run | Dual Flowcell v2      | 1200M                               | 47                  | 29          | 18          |              |
| HiSeq 3000           | 8 lanes per flow cell | 5B                                  | 196                 | 123         | 76          |              |
| HiSeq 4000           | 8 lanes per flow cell | 10B                                 | 392                 | 247         | 152         |              |
| NovaSeq 6000         | SP (per flow cell)    | 1.6B                                | 62                  | 39          | 24          |              |

N/A: not applicable, no samples can be run.

\* Based on 2 x 149 bp paired-end reads.

**Table 5b. Number of multiplexed samples (one DNA and one RNA library for each sample) based on a DNA+RNA panel size (Combined targeted DNA+RNA Enrichment workflow) with 3640x mean coverage for the DNA library and 5000 reads/primer for the RNA library\***

|                      |                       |                              | UHS-005Z     | UHS-009Z | UHS-006Z |       |
|----------------------|-----------------------|------------------------------|--------------|----------|----------|-------|
|                      |                       |                              | DNA primer # | 4149     | 6244     | 11243 |
|                      |                       |                              | RNA primer # | 487      | 1116     | 665   |
| Instrument           | Version               | Capacity (paired-ends reads) |              |          |          |       |
| MiniSeq              | Mid output            | 16M                          | N/A          | N/A      | N/A      |       |
| MiniSeq              | High output           | 50M                          | 1            | N/A      | N/A      |       |
| MiSeq                | v2                    | 30M                          | N/A          | N/A      | N/A      |       |
| NextSeq 500          | Mid output            | 260M                         | 6            | 4        | 2        |       |
| NextSeq 500          | High output           | 800M                         | 19           | 12       | 7        |       |
| HiSeq 2500 rapid run | Dual Flowcell v2      | 1200M                        | 29           | 18       | 11       |       |
| HiSeq 3000           | 8 lanes per flow cell | 5B                           | 124          | 78       | 48       |       |
| HiSeq 4000           | 8 lanes per flow cell | 10B                          | 249          | 157      | 97       |       |
| NovaSeq 6000         | SP (per flow cell)    | 1.6B                         | 39           | 25       | 15       |       |

N/A: not applicable, no samples can be run.

\* Based on 2 x 149 bp paired-end reads.

**Table 5c. Number of multiplexed samples (one DNA and one RNA library for each sample) based on a DNA+RNA panel size (Combined targeted DNA+RNA Enrichment workflow) with 7200x mean coverage for the DNA library and 5000 reads/primer for the RNA library\***

|                      |                       |                              | UHS-005Z     | UHS-009Z | UHS-006Z |       |
|----------------------|-----------------------|------------------------------|--------------|----------|----------|-------|
|                      |                       |                              | DNA primer # | 4149     | 6244     | 11243 |
|                      |                       |                              | RNA primer # | 487      | 1116     | 665   |
| Instrument           | Version               | Capacity (paired-ends reads) |              |          |          |       |
| MiniSeq              | Mid output            | 16M                          | N/A          | N/A      | N/A      |       |
| MiniSeq              | High output           | 50M                          | N/A          | N/A      | N/A      |       |
| MiSeq                | v2                    | 30M                          | N/A          | N/A      | N/A      |       |
| NextSeq 500          | Mid output            | 260M                         | 4            | 2        | 1        |       |
| NextSeq 500          | High output           | 800M                         | 14           | 8        | 5        |       |
| HiSeq 2500 rapid run | Dual Flowcell v2      | 1200M                        | 21           | 13       | 8        |       |
| HiSeq 3000           | 8 lanes per flow cell | 5B                           | 88           | 55       | 34       |       |
| HiSeq 4000           | 8 lanes per flow cell | 10B                          | 176          | 111      | 68       |       |
| NovaSeq 6000         | SP (per flow cell)    | 1.6B                         | 28           | 17       | 11       |       |

N/A: not applicable, no samples can be run.

\* Based on 2 x 149 bp paired-end reads.

**Table 5d. Number of multiplexed samples (one DNA and one RNA library for each sample) based on a DNA+RNA panel size (Combined targeted DNA+RNA Enrichment workflow) with 25,600x mean coverage for the DNA library and 5000 reads/primer for the RNA library\***

|                      |                       |                              | UHS-005Z     | UHS-009Z | UHS-006Z |       |
|----------------------|-----------------------|------------------------------|--------------|----------|----------|-------|
|                      |                       |                              | DNA primer # | 4149     | 6244     | 11243 |
|                      |                       |                              | RNA primer # | 487      | 1116     | 665   |
| Instrument           | Version               | Capacity (paired-ends reads) |              |          |          |       |
| MiniSeq              | Mid output            | 16M                          | N/A          | N/A      | N/A      | N/A   |
| MiniSeq              | High output           | 50M                          | N/A          | N/A      | N/A      | N/A   |
| MiSeq                | v2                    | 30M                          | N/A          | N/A      | N/A      | N/A   |
| NextSeq 500          | Mid output            | 260M                         | 1            | 1        | N/A      | N/A   |
| NextSeq 500          | High output           | 800M                         | 5            | 3        | 2        |       |
| HiSeq 2500 rapid run | Dual Flowcell v2      | 1200M                        | 8            | 5        | 3        |       |
| HiSeq 3000           | 8 lanes per flow cell | 5B                           | 35           | 22       | 13       |       |
| HiSeq 4000           | 8 lanes per flow cell | 10B                          | 70           | 44       | 27       |       |
| NovaSeq 6000         | SP (per flow cell)    | 1.6B                         | 11           | 7        | 4        |       |

N/A: not applicable, no samples can be run.

\* Based on 2 x 149 bp paired-end reads.

**Table 6a. Number of multiplexed samples (one DNA and one RNA library for each sample) based on DNA or RNA panel size (Separated targeted DNA and RNA Enrichment workflow) with 500x mean coverage for DNA library and 5000 reads/primer for the RNA library\***

|                      |                       |                              | UHS-005Z     | UHS-009Z | UHS-006Z | UHS-5000Z |       |
|----------------------|-----------------------|------------------------------|--------------|----------|----------|-----------|-------|
|                      |                       |                              | DNA primer # | 4149     | 6244     | 11243     | 19995 |
|                      |                       |                              | RNA primer # | 487      | 1116     | 665       | 2571  |
| Instrument           | Version               | Capacity (paired-ends reads) |              |          |          |           |       |
| MiniSeq              | Mid output            | 16M                          | 3            | 1        | 1        | N/A       |       |
| MiniSeq              | High output           | 50M                          | 11           | 5        | 5        | 2         |       |
| MiSeq                | v2                    | 30M                          | 6            | 3        | 3        | 1         |       |
| NextSeq 500          | Mid output            | 260M                         | 57           | 29       | 29       | 11        |       |
| NextSeq 500          | High output           | 800M                         | 177          | 91       | 89       | 35        |       |
| HiSeq 2500 rapid run | Dual Flowcell v2      | 1200M                        | 266          | 137      | 134      | 52        |       |
| HiSeq 3000           | 8 lanes per flow cell | 5B                           | 1108         | 574      | 558      | 218       |       |
| HiSeq 4000           | 8 lanes per flow cell | 10B                          | 2217         | 1149     | 1117     | 437       |       |
| NovaSeq 6000         | SP (per flow cell)    | 1.6B                         | 354          | 183      | 178      | 70        |       |

N/A: not applicable, no samples can be run.

\* Based on 2 x 149 bp paired-end reads.

**Table 6b. Number of multiplexed samples (one DNA and one RNA library for each sample) based on DNA or RNA panel size (Separated targeted DNA and RNA Enrichment workflow) with 3640x mean coverage for DNA library and 5000 reads/primer for the RNA library\***

|                         |                       |                              | UHS-005Z     | UHS-009Z | UHS-006Z | UHS-5000Z |       |
|-------------------------|-----------------------|------------------------------|--------------|----------|----------|-----------|-------|
|                         |                       |                              | DNA primer # | 4149     | 6244     | 11243     | 19995 |
|                         |                       |                              | RNA primer # | 487      | 1116     | 665       | 2571  |
| Instrument              | Version               | Capacity (paired-ends reads) |              |          |          |           |       |
| MiniSeq                 | Mid output            | 16M                          | N/A          | N/A      | N/A      | N/A       | N/A   |
| MiniSeq                 | High output           | 50M                          | 2            | 1        | 1        | N/A       | N/A   |
| MiSeq                   | v2                    | 30M                          | 1            | 1        | N/A      | N/A       | N/A   |
| NextSeq 500             | Mid output            | 260M                         | 14           | 9        | 5        | 3         |       |
| NextSeq 500             | High output           | 800M                         | 45           | 28       | 18       | 9         |       |
| HiSeq 2500<br>rapid run | Dual Flowcell v2      | 1200M                        | 68           | 42       | 27       | 14        |       |
| HiSeq 3000              | 8 lanes per flow cell | 5B                           | 285          | 176      | 112      | 58        |       |
| HiSeq 4000              | 8 lanes per flow cell | 10B                          | 570          | 353      | 225      | 116       |       |
| NovaSeq 6000            | SP (per flow cell)    | 1.6B                         | 91           | 56       | 36       | 18        |       |

N/A: not applicable, no samples can be run.

\* Based on 2 x 149 bp paired-end reads.

**Table 6c. Number of multiplexed samples (one DNA and one RNA library for each sample) based on DNA or RNA panel size (Separated targeted DNA and RNA Enrichment workflow) with 7200x mean coverage for DNA library and 5000 reads/primer for the RNA library\***

|                         |                       |                              | UHS-005Z     | UHS-009Z | UHS-006Z | UHS-5000Z |       |
|-------------------------|-----------------------|------------------------------|--------------|----------|----------|-----------|-------|
|                         |                       |                              | DNA primer # | 4149     | 6244     | 11243     | 19995 |
|                         |                       |                              | RNA primer # | 487      | 1116     | 665       | 2571  |
| Instrument              | Version               | Capacity (paired-ends reads) |              |          |          |           |       |
| MiniSeq                 | Mid output            | 16M                          | N/A          | N/A      | N/A      | N/A       | N/A   |
| MiniSeq                 | High output           | 50M                          | 1            | N/A      | N/A      | N/A       | N/A   |
| MiSeq                   | v2                    | 30M                          | N/A          | N/A      | N/A      | N/A       | N/A   |
| NextSeq 500             | Mid output            | 260M                         | 8            | 5        | 3        | 1         |       |
| NextSeq 500             | High output           | 800M                         | 24           | 15       | 9        | 5         |       |
| HiSeq 2500<br>rapid run | Dual Flowcell v2      | 1200M                        | 37           | 23       | 14       | 7         |       |
| HiSeq 3000              | 8 lanes per flow cell | 5B                           | 154          | 98       | 59       | 31        |       |
| HiSeq 4000              | 8 lanes per flow cell | 10B                          | 309          | 197      | 118      | 63        |       |
| NovaSeq 6000            | SP (per flow cell)    | 1.6B                         | 49           | 31       | 18       | 10        |       |

N/A: not applicable, no samples can be run.

\* Based on 2 x 149 bp paired-end reads.

**Table 6d. Number of multiplexed samples (one DNA and one RNA library for each sample) based on DNA or RNA panel size (Separated targeted DNA and RNA Enrichment workflow) with 25,600x mean coverage for DNA library and 5000 reads/primer for the RNA library\***

|                         |                       |                              | UHS-005Z     | UHS-009Z | UHS-006Z | UHS-5000Z |
|-------------------------|-----------------------|------------------------------|--------------|----------|----------|-----------|
|                         |                       |                              | DNA primer # | 6244     | 11243    | 19995     |
|                         |                       |                              | RNA primer # | 1116     | 665      | 2571      |
| Instrument              | Version               | Capacity (paired-ends reads) |              |          |          |           |
| MiniSeq                 | Mid output            | 16M                          | N/A          | N/A      | N/A      | N/A       |
| MiniSeq                 | High output           | 50M                          | N/A          | N/A      | N/A      | N/A       |
| MiSeq                   | v2                    | 30M                          | N/A          | N/A      | N/A      | N/A       |
| NextSeq 500             | Mid output            | 260M                         | 2            | 1        | N/A      | N/A       |
| NextSeq 500             | High output           | 800M                         | 7            | 4        | 2        | 1         |
| HiSeq 2500<br>rapid run | Dual Flowcell v2      | 1200M                        | 11           | 7        | 4        | 2         |
| HiSeq 3000              | 8 lanes per flow cell | 5B                           | 46           | 30       | 17       | 9         |
| HiSeq 4000              | 8 lanes per flow cell | 10B                          | 92           | 60       | 34       | 19        |
| NovaSeq 6000            | SP (per flow cell)    | 1.6B                         | 14           | 9        | 5        | 3         |

N/A: not applicable, no samples can be run.

\* Based on 2 x 149 bp paired-end reads.

# Equipment and Reagents to Be Supplied by User

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, consult the appropriate safety data sheets (SDSs), available from the product supplier.

In addition to the QIAseq Multimodal Panels and the QIAseq Multimodal Index kit, the following are required:

- EvaGreen®, 20x in water (Biotium, cat. no. 31000-T or 31000)
- 80% ethanol (made fresh daily) \*
- Nuclease-free pipette tips and tubes
- 1.5 ml LoBind® tubes (Eppendorf, cat. no. 022431021)
- PCR tubes (0.2 ml individual tubes [VWR, cat. no. 20170-012] or tube strips [VWR, cat. no. 93001-118]) or plates
- Ice
- Microcentrifuge
- Thermal cycler
- Magnet for bead cleanups:
  - **Tubes:** MagneSphere® Technology Magnetic Separation Stand (Promega, cat. no. Z5342)
  - **Plates:** DynaMag™-96 Side Magnet (Thermo Fisher Scientific Inc., cat. no. 12331D)
- 2100 Bioanalyzer (Agilent, cat. no. G2939BA)
- Agilent® High Sensitivity DNA Kit (Agilent, cat. no. 5067-4626)

\* Do not use denatured alcohol, which contains other substances, such as methanol or methylethylketone.

# Protocol: Nucleic Acid Fragmentation, Standard Samples

## Important points before starting

- This protocol describes fragmentation of nucleic acids from “standard samples” (i.e., cells or tissue). For fragmentation of FFPE samples, please refer to “Protocol: Nucleic Acid Fragmentation, FFPE Samples”, page 32.
- This protocol is designed to work with either total nucleic acid eluates (those containing DNA+RNA) or separate DNA and RNA eluates.
- When performing “Protocol: Combined Targeted DNA+RNA Enrichment in a Single Tube”, the recommended amount of DNA is 10–40 ng.
- When performing “Appendix B: Separated Targeted DNA and RNA Enrichment in Separate Tubes”, the recommended amount of DNA is 20-80 ng.
- The recommended amount of RNA is 10 ng to 250 ng total RNA. When working with total nucleic acid samples, we recommend basing the input on the amount of quantified DNA.
- Set up reactions on ice.
- Do not vortex any reagents or reactions.

## Procedure

1. Thaw nucleic acid samples on ice. Gently mix, briefly centrifuge to collect residual liquid from the sides of the tubes, and then return to ice.
2. Prepare the reagents required for the fragmentation.
  - 2a. Thaw Fragmentation Buffer, 10x, and FERA Solution at room temperature (15–25°C).
  - 2b. Mix by flicking the tube, and then centrifuge briefly.

**Note:** Side Reaction Reducer and Fragmentation Enzyme Mix should be removed from the freezer just before use and placed on ice. After use, immediately return the enzymes to the freezer.

- On ice, prepare the fragmentation mix according to Table 7. Briefly centrifuge, mix by pipetting up and down 10–12 times, and then briefly centrifuge again.

**Note:** If setting up more than one reaction, prepare a volume of master mix 10% greater than what is required for the total number of reactions.

**Table 7. Reaction mix for nucleic acid fragmentation**

| Component   | Volume/reaction                                    |
|---|--|
| DNA (see input recommendation in the “Important points before starting” section)* | Variable A   |
| DNA (see input recommendation in the “Important points before starting” section)* | Variable B   |
| Fragmentation Buffer, 10x   | 2 $\mu$ l  |
| FERA Solution   | 0.6 $\mu$ l  |
| Side Reaction Reducer   | 1.6 $\mu$ l  |
| Fragmentation Enzyme Mix  | 4 $\mu$ l  |
| Nuclease-free Water   | 11.8 $\mu$ l – variable A (DNA) – variable B (RNA) |
| <b>Total</b>  | <b>20 <math>\mu</math>l</b>                        |

\* Instead of adding DNA and RNA separately, total nucleic acid containing both DNA and RNA can be added. If adding total nucleic acid, base the input amount on DNA.

- Program the thermal cycler according to Table 8. Use the instrument’s heated lid.

**Table 8. Incubation conditions for nucleic acid fragmentation**

| Step | Incubation temperature | Incubation time |
|------|------------------------|-----------------|
| 1    | 4°C                    | 1 min           |
| 2    | 32°C                   | 24 min          |
| 3    | 72°C                   | 30 min          |
| 4    | 4°C                    | Hold            |

- Prior to adding the tubes/plate to a thermal cycler, start the program. When the thermal cycler reaches 4°C, pause the program.

Important: The thermal cycler must be prechilled and paused at 4°C.

- 
6. Transfer the tubes/plate prepared in step 3 to the prechilled thermal cycler and resume the program.
  7. Upon completion, allow the thermal cycler to return to 4°C.
  8. Place the samples on ice and immediately proceed to “Protocol: RNA Polyadenylation”, page 35.

# Protocol: Nucleic Acid Fragmentation, FFPE Samples

## Important points before starting

- This protocol describes fragmentation of nucleic acids from FFPE samples. For fragmentation of “standard samples” (i.e., cells or tissue), please refer to “Protocol: Nucleic Acid Fragmentation, Standard Samples”, page 29.
- This protocol is designed to work with either total nucleic acid eluates (those containing DNA+RNA) or separate DNA and RNA eluates.
- The recommended amount of FFPE DNA is up to 250 ng DNA if QIAseq QuantiMIZE kits have been used (See “Appendix A: FFPE Sample Quality and Quantity”, page 72). If an alternative method was used to determine the concentration of FFPE DNA, then up to 100 ng DNA can be used. The recommended amount of FFPE RNA is 250 ng total RNA (up to 500 ng for “severely” fragmented FFPE sample, with “severely” being defined as samples that have less than 40% of fragments >200 nt by smear analysis on the Bioanalyzer). When working with total nucleic acid samples, we recommend basing the input on the amount of quantified DNA.
- Set up reactions on ice.
- Do not vortex any reagents or reactions.

## Procedure

1. Thaw nucleic acid samples on ice. Gently mix, briefly centrifuge to collect residual liquid from the sides of the tubes, and then return to ice.
2. Prepare the reagents required for the fragmentation.
  - 2a. Thaw Fragmentation Buffer, 10x, and FERA Solution at room temperature.
  - 2b. Mix by flicking the tube, and centrifuge briefly.

**Note:** Side Reaction Reducer and Fragmentation Enzyme Mix should be removed from the freezer just before use and placed on ice. After use, immediately return the enzymes to the freezer.

3. On ice, prepare the fragmentation mix according to Table 9. Briefly centrifuge, mix by pipetting up and down 10–12 times, and then briefly centrifuge again.

**Note:** If setting up more than one reaction, prepare a volume of master mix 10% greater than what is required for the total number of reactions.

**Table 9. Reaction mix for nucleic acid fragmentation**

| Component   | Volume/reaction                               |
|---|---|
| DNA (see input recommendation in the “Important points before starting” section)* | Variable A                                    |
| RNA (see input recommendation in the “Important points before starting” section)* | Variable B                                    |
| Fragmentation Buffer, 10x   | 2 µl  |
| FERA Solution   | 0.6 µl  |
| Side Reaction Reducer   | 1.6 µl  |
| Nuclease-free Water   | 11.8 µl – variable A (DNA) – variable B (RNA) |
| <b>Total</b>  | <b>16 µl</b>                                  |

\* Instead of adding DNA and RNA separately, total nucleic acid containing both DNA and RNA can be added. If adding total nucleic acid, base the input amount on DNA.

4. Incubate for 15 min at 37°C and then place on ice.

5. Add 4 µl of Fragmentation Enzyme Mix to each reaction. Briefly centrifuge, mix by pipetting up and down 10–12 times (do not vortex), and then briefly centrifuge again.

**Important:** Keep the reaction tubes/plate on ice during the entire reaction setup.

6. Program the thermal cycler according to Table 10. Use the instrument’s heated lid.

**Table 10. Incubation conditions for nucleic acid fragmentation**

| Step | Incubation temperature | Incubation time |
|------|------------------------|-----------------|
| 1    | 4°C                    | 1 min           |
| 2    | 32°C                   | 14 min          |
| 3    | 72°C                   | 30 min          |
| 4    | 4°C                    | Hold            |

---

7. Prior to adding the tubes/plate to a thermal cycler, start the program. When the thermal cycler reaches 4°C, pause the program.

**Important:** The thermal cycler must be prechilled and paused at 4°C.

8. Transfer the tubes/plate prepared in step 5 to the prechilled thermal cycler and resume the program.

9. Upon completion, allow the thermal cycler to return to 4°C.

10. Place the samples on ice, and immediately proceed to “Protocol: RNA Polyadenylation”, page 35.

---

# Protocol: RNA Polyadenylation

## Important points before starting

- The product from “Protocol: Nucleic Acid Fragmentation, Standard Samples”, page 29, or “Protocol: Nucleic Acid Fragmentation, FFPE Samples”, page 32, is the starting material for this protocol.
- Set up reactions on ice.
- Do not vortex any reagents or reactions.

## Procedure

1. Prepare the reagents required for the polyadenylation.
  - 1a. Thaw PAP Dilution Buffer, 10x, and ATP Solution at room temperature.
  - 1b. Mix by flicking the tube, and then centrifuge briefly.  
**Note:** T4 Polynucleotide Kinase and PAP Enzyme should be removed from the freezer just before use and placed on ice. After use, immediately return the enzymes to the freezer.
2. Prepare 1x PAP Dilution Buffer by diluting 2  $\mu\text{l}$  of the 10x PAP Dilution Buffer with 18  $\mu\text{l}$  Nuclease-free Water.
3. Use the 1x PAP Dilution Buffer to dilute an aliquot of the PAP Enzyme from 5 U/ $\mu\text{l}$  to 2 U/ $\mu\text{l}$ . Briefly centrifuge, mix by pipetting up and down 10–12 times, and then briefly centrifuge again.
4. Prepare the RNA polyadenylation mix according to Table 11. Briefly centrifuge, mix by pipetting up and down 10–12 times, and then briefly centrifuge again.  
**Note:** If setting up more than one reaction, prepare a volume of master mix 10% greater than what is required for the total number of reactions.

**Table 11. Reaction mix for RNA polyadenylation**

| Component                                | Volume/reaction             |
|--|-----------------------------|
| Fragmentation reaction (already in tube) | 20 $\mu$ l                  |
| ATP Solution                             | 1.25 $\mu$ l                |
| T4 Polynucleotide Kinase                 | 1 $\mu$ l                   |
| Diluted PAP Enzyme (2 U/ $\mu$ l)*       | 1 $\mu$ l                   |
| Nuclease-free Water                      | 1.75 $\mu$ l                |
| <b>Total</b>                             | <b>25 <math>\mu</math>l</b> |

\* Ensure PAP Enzyme has been diluted from its stock 5U/ $\mu$ l concentration to 2U/ $\mu$ l using 1x PAP Dilution Buffer.

5. Incubate the reactions in a thermal cycler according to Table 12. Use the instrument's heated lid.

**Table 12. Incubation conditions for RNA polyadenylation**

| Step | Incubation temperature | Incubation time |
|------|------------------------|-----------------|
| 1    | 4°C                    | 1 min           |
| 2    | 30°C                   | 10 min          |
| 3    | 4°C                    | Hold            |

6. Upon completion, place the reactions on ice and proceed to "Protocol: DNA Ligation", page 37.

---

# Protocol: DNA Ligation

## Important points before starting

- The product from “Protocol: RNA Polyadenylation”, page 35, is the starting material for this protocol.
- Set up reactions on ice.
- Do not vortex any reagents or reactions.
- QIAseq Beads are used for all reaction cleanups.
- **Important:** Prepare fresh 80% ethanol daily.
- Ensure that the QIAseq Beads are thoroughly mixed at all times. This necessitates working quickly and resuspending the beads immediately before use. If a delay in the protocol occurs, simply vortex the beads.

## Procedure

1. Prepare the reagents required for the DNA ligation.
  - 1a. Thaw DNA Ligation Adapter; Ligation Buffer, 5x; and Ligation Solution at room temperature.
  - 1b. Mix by flicking the tube, and then centrifuge briefly.

**Note:** DNA Ligase should be removed from the freezer just before use and placed on ice. After use, immediately return the enzyme to the freezer.
2. Prepare the DNA ligation mix according to Table 13. Briefly centrifuge, mix by pipetting up and down 10–12 times, and then briefly centrifuge again.

**Important:** Pipet slowly to mix. The reaction mix is very viscous. Do not vortex.

**Note:** If setting up more than one reaction, prepare a volume of master mix 10% greater than what is required for the total number of reactions.

**Table 13. Reaction mix for DNA ligation**

| Component                                      | Volume/reaction             |
|--|-----------------------------|
| RNA polyadenylation reaction (already in tube) | 25 $\mu$ l                  |
| Ligation Buffer, 5x                            | 10 $\mu$ l                  |
| DNA Ligation Adapter                           | 2.8 $\mu$ l                 |
| DNA Ligase                                     | 5 $\mu$ l                   |
| Ligation Solution*                             | 7.2 $\mu$ l                 |
| <b>Total</b>                                   | <b>50 <math>\mu</math>l</b> |

\* Ligation Solution is very viscous. It should be added into each reaction individually and not premixed with other components for a master mix. Do not coat the outside of the pipette tip with Ligation Solution, because in doing so, excess volume may be added.

3. Incubate the reactions in a thermal cycler according to Table 14.

**Important:** Do not use the heated lid.

**Table 14. Incubation conditions for DNA ligation**

| Step | Incubation temperature | Incubation time |
|------|------------------------|-----------------|
| 1    | 4°C                    | 1 min           |
| 2    | 20°C                   | 15 min          |
| 3    | 4°C                    | Hold            |

4. Add 50  $\mu$ l of Nuclease-free Water to bring each sample to 100  $\mu$ l.

5. Add 130  $\mu$ l QIAseq Beads, and then mix by vortexing.

6. Incubate for 5 min at room temperature.

7. Place the tubes/plate on a magnetic rack for 2 min (for tubes) or 10 min (for plates).

8. After the solution has cleared, with the beads still on the magnetic stand, carefully remove and discard the supernatant.

**Important:** Do not discard the beads, because they contain the DNA of interest.

9. Add 80  $\mu$ l of Nuclease-free Water to resuspend the beads and then 128  $\mu$ l of QIAseq NGS Bead Binding Buffer. Mix by vortexing and incubate for 5 min at room temperature.

10. Place the tubes/plate on a magnetic rack for 2 min (for tubes) or 10 min (for plates). Once the solution has cleared, with the beads still on the magnetic stand, carefully remove and discard the supernatant.
11. With the beads still on the magnetic stand, add 200  $\mu$ l of 80% ethanol. Rotate the tube 2–3 times or move the plate side-to-side between the 2 column positions of the magnet to wash the beads. Carefully remove and discard the wash.
12. Repeat the ethanol wash.  
**Important:** Completely remove all traces of the ethanol after the second wash. Remove the ethanol with a 200  $\mu$ l pipette tip first, spin down briefly and then use a 10  $\mu$ l pipette tip to remove any residual ethanol.
13. With the beads still on the magnetic stand, air-dry at room temperature for 5–10 min.  
**Note:** Visually inspect the pellet to confirm that it is completely dry. If the pellet is not completely dry, additional drying time may be required.
14. Remove the beads from the magnetic stand, and then elute the DNA from the beads by adding 19  $\mu$ l Nuclease-free Water. Mix well by pipetting.
15. Return the tube/plate to the magnetic rack until the solution has cleared.
16. Transfer 16.62  $\mu$ l of the supernatant to clean tubes/plate.  
Proceed to “Protocol: Reverse Transcription”, page 40. Alternatively, the samples can be stored at  $-30$  to  $-15^{\circ}\text{C}$  in a constant-temperature freezer.

---

# Protocol: Reverse Transcription

## Important points before starting

- The 16.62  $\mu$ l product from “Protocol: DNA Ligation”, page 37, is the starting material for this protocol.
- Set up reactions on ice.
- Do not vortex any reagents or reactions.
- QIAseq Beads are used for all reaction cleanups.
- **Important:** Prepare fresh 80% ethanol daily.
- Ensure that the QIAseq Beads are thoroughly mixed at all times. This necessitates working quickly and resuspending the beads immediately before use. If a delay in the protocol occurs, simply vortex the beads.

## Procedure

1. Prepare the reagents required for the reverse transcription.
  - 1a. Thaw the Multimodal RT Primer; Multimodal RT Buffer, 5x; and Multimodal RT Enhancer at room temperature.
  - 1b. Mix by flicking the tube, and then centrifuge briefly.

**Note:** The RNase Inhibitor and EZ Reverse Transcriptase should be removed from the freezer just before use and placed on ice. After use, immediately return the enzymes to the freezer.
2. Prepare the reverse transcription mix according to Table 15. Briefly centrifuge, mix by pipetting up and down 10–12 times, and then briefly centrifuge again.

**Note:** If setting up more than one reaction, prepare a volume of master mix 10% greater than what is required for the total number of reactions.

**Table 15. Reaction mix for reverse transcription**

| Component                                       | Volume/reaction             |
|---|-----------------------------|
| Sample (from "Protocol: DNA Ligation", page 37) | 16.62 $\mu$ l               |
| Multimodal RT Primer                            | 1 $\mu$ l                   |
| Multimodal RT Buffer, 5x                        | 5 $\mu$ l                   |
| Multimodal RT Enhancer                          | 0.5 $\mu$ l                 |
| RNase Inhibitor                                 | 0.63 $\mu$ l                |
| EZ Reverse Transcriptase                        | 1.25 $\mu$ l                |
| <b>Total</b>                                    | <b>25 <math>\mu</math>l</b> |

3. Incubate the reactions in a thermal cycler according to Table 16. Use the instrument's heated lid.

**Table 16. Incubation conditions for reverse transcription**

| Step | Incubation temperature | Incubation time |
|------|------------------------|-----------------|
| 1    | 4°C                    | 1 min           |
| 2    | 25°C                   | 10 min          |
| 3    | 42°C                   | 45 min          |
| 4    | 70°C                   | 15 min          |
| 5    | 4°C                    | Hold            |

4. Add 75  $\mu$ l of Nuclease-free Water to bring each sample to 100  $\mu$ l.
5. Add 130  $\mu$ l QIAseq Beads and mix by vortexing or by pipetting up and down several times.
6. Incubate for 5 min at room temperature.
7. Place the tubes/plate on a magnetic rack for 2 min (for tubes) or 10 min (for plates).
8. After the solution has cleared, with the beads still on the magnetic stand, carefully remove and discard the supernatant.

**Important:** Do not discard the beads, because they contain the DNA of interest.

9. With the beads still on the magnetic stand, add 200  $\mu$ l of 80% ethanol. Rotate the tube 2–3 times or move the plate side-to-side between the 2 column positions of the magnet to wash the beads. Carefully remove and discard the wash.
10. Repeat the ethanol wash.  
**Important:** Completely remove all traces of the ethanol after this second wash. Remove the ethanol with a 200  $\mu$ l pipette tip first, spin down briefly, and then use a 10  $\mu$ l pipette tip to remove any residual ethanol.
11. With the beads still on the magnetic stand, air-dry at room temperature for 5–10 min.  
**Note:** Visually inspect the pellet to confirm that it is completely dry. If the pellet is not completely dry, additional drying time may be required. Ethanol carryover to the next universal PCR step will affect PCR efficiency.
12. Remove the beads from the magnetic stand, and then elute the DNA from the beads by adding 15  $\mu$ l Nuclease-free Water.  
**Important: If performing separated DNA and RNA target enrichment (Appendix B, page 74), elute by adding 22.4  $\mu$ l Nuclease-free Water.**
13. Return the tube/plate to the magnetic rack until solution the solution has cleared.
14. Transfer 12.4  $\mu$ l of the supernatant to clean tubes/plate.  
**Important: If performing separated DNA and RNA target enrichment (Appendix B), transfer 10.2  $\mu$ l of the eluate to each of 2 tubes and proceed to the Appendix B protocol.**
15. Proceed to “Protocol: Combined Targeted DNA+RNA Enrichment in a Single Tube”, page 43. Alternatively, the samples can be stored at  $-30$  to  $-15^{\circ}\text{C}$  in a constant-temperature freezer.

# Protocol: Combined Targeted DNA+RNA Enrichment in a Single Tube

## Important points before starting

- The 12.4 µl product from “Protocol: Reverse Transcription”, page 40, is the starting material for this protocol.
- Set up reactions on ice.
- Do not vortex any reagents or reactions.
- The final library dual sample index is determined by the combination of the QIAseq Multimodal N7 Plate and the QIAseq Multimodal S5 Plate. QIAseq Beads are used for all reaction cleanups.

**IMPORTANT:** To use this protocol, one of the following is required:

- **QIAseq Multimodal Index I (12)** (cat. no. 333962)
- **QIAseq Multimodal Index I Set A (96)** (cat. no. 333965)
- **QIAseq Multimodal Index I Set B (96)** (cat. no. 333975)

These plates come in either 12- or 48-reaction formats allowing the indexing of 12 DNA and RNA samples or 48 DNA and RNA samples (using one plate of Set A or Set B). In each indicated well of the cuttable plate, there are dried N7 index primers for both DNA and RNA. The plates can be cut in columns to enable indexing of the desired number of samples.

Two plates of Set A or Set B are included for each index kit for making a total of 96 DNA and 96 RNA libraries. By combining Set A and Set B, up to 96 DNA and 96 RNA libraries can be multiplexed. Each well in the plate is single use.

- **Important:** The required combinations of indexes are described in the sequencing sample setup sheets:
  - Sample Sheet Multimodal UDI Set A: [www.qiagen.com/PROM-15281](http://www.qiagen.com/PROM-15281)
  - Sample Sheet Multimodal UDI Set B: [www.qiagen.com/PROM-15282](http://www.qiagen.com/PROM-15282)
  - Sample Sheet Multimodal UDI Set A and Set B: [www.qiagen.com/PROM-15283](http://www.qiagen.com/PROM-15283)
- **Important:** Prepare fresh 80% ethanol daily.
- Ensure that the QIAseq Beads are thoroughly mixed at all times. This necessitates working quickly and resuspending the beads immediately before use. If a delay in the protocol occurs, simply vortex the beads.

## Procedure

1. Prepare the reagents required for target enrichment.
  - 1a. Thaw TEPCR Buffer, 5x; Multimodal DHS Panel (DNA); and Multimodal VHS Panel (RNA); and bring QIAseq Multimodal N7 Plate to room temperature.
  - 1b. Mix by flicking the tube, and then centrifuge briefly.

**Note:** QIAseq Multimodal N7 Plate only needs to be centrifuged, not mixed.

**Note:** HotStarTaq DNA Polymerase should be removed from the freezer just before use and placed on ice. After use, immediately return the enzyme to the freezer.
2. Prepare the target enrichment mix according to Table 17. Briefly centrifuge, mix by pipetting up and down 10–12 times, and then briefly centrifuge again.

**Note:** If setting up more than one reaction, prepare a volume of master mix 10% greater than what is required for the total number of reactions.

**Table 17. Reaction mix for target enrichment**

| <b>Component</b>   | <b>Volume/reaction</b>      |
|--|-----------------------------|
| Sample (from "Protocol: Reverse Transcription", page 40) | 12.4 $\mu$ l                |
| TEPCR Buffer, 5x   | 8 $\mu$ l                   |
| Multimodal DHS Panel (DNA)                               | 10 $\mu$ l                  |
| Multimodal VHS Panel (RNA)                               | 8 $\mu$ l                   |
| HotStarTaq DNA Polymerase (6 U/ $\mu$ l)                 | 1.6 $\mu$ l                 |
| <b>Total</b>   | <b>40 <math>\mu</math>l</b> |

3. Add the 40  $\mu$ l target enrichment reaction mix into a well of a QIAseq Multimodal N7 Plate (Table 18, Table 19, and Table 20), which are cuttable plates that contain predispensed, dried index primer pairs for both DNA and RNA samples in the same well.

**Note:** The plates can be cut in columns to enable indexing of the desired number of samples.

**Table 18. Layout of QIAseq Multimodal N7 Plates, 12 reactions**

|          | 1                            | 2                             | 3     | 4     | 5     | 6     | 7     | 8     | 9     | 10    | 11    | 12    |
|----------|------------------------------|-------------------------------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| <b>A</b> | S1<br>DNAp-M001<br>RNAp-M049 | S9<br>DNAp-M009<br>RNAp-M067  | Empty |
| <b>B</b> | S2<br>DNAp-M002<br>RNAp-M050 | S10<br>DNAp-M010<br>RNAp-M059 | Empty |
| <b>C</b> | S3<br>DNAp-M003<br>RNAp-M051 | S11<br>DNAp-M011<br>RNAp-M059 | Empty |
| <b>D</b> | S4<br>DNAp-M004<br>RNAp-M052 | S12<br>DNAp-M012<br>RNAp-M060 | Empty |
| <b>E</b> | S5<br>DNAp-M005<br>RNAp-M053 | Empty                         | Empty | Empty | Empty | Empty | Empty | Empty | Empty | Empty | Empty | Empty |
| <b>F</b> | S6<br>DNAp-M006<br>RNAp-M054 | Empty                         | Empty | Empty | Empty | Empty | Empty | Empty | Empty | Empty | Empty | Empty |
| <b>G</b> | S7<br>DNAp-M007<br>RNAp-M055 | Empty                         | Empty | Empty | Empty | Empty | Empty | Empty | Empty | Empty | Empty | Empty |
| <b>H</b> | S8<br>DNAp-M008<br>RNAp-M056 | Empty                         | Empty | Empty | Empty | Empty | Empty | Empty | Empty | Empty | Empty | Empty |

**Table 19. Layout of QIAseq Multimodal N7 Plates, 48 reactions Set A**

|   | 1                            | 2                             | 3                             | 4                             | 5                             | 6                             | 7     | 8     | 9     | 10    | 11    | 12    |
|---|------------------------------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|-------|-------|-------|-------|-------|-------|
| A | S1<br>DNAp-M001<br>RNAp-M049 | S9<br>DNAp-M009<br>RNAp-M057  | S17<br>DNAp-M095<br>RNAp-M095 | S25<br>DNAp-M026<br>RNAp-M073 | S33<br>DNAp-M033<br>RNAp-M081 | S41<br>DNAp-M041<br>RNAp-M086 | Empty | Empty | Empty | Empty | Empty | Empty |
|   | S2<br>DNAp-M002<br>RNAp-M050 | S10<br>DNAp-M010<br>RNAp-M058 | S18<br>DNAp-M016<br>RNAp-M096 | S26<br>DNAp-M026<br>RNAp-M074 | S34<br>DNAp-M034<br>RNAp-M082 | S42<br>DNAp-M042<br>RNAp-M090 | Empty | Empty | Empty | Empty | Empty | Empty |
| B | S3<br>DNAp-M003<br>RNAp-M051 | S11<br>DNAp-M009<br>RNAp-M059 | S19<br>DNAp-M019<br>RNAp-M097 | S27<br>DNAp-M027<br>RNAp-M075 | S35<br>DNAp-M035<br>RNAp-M083 | S43<br>DNAp-M043<br>RNAp-M091 | Empty | Empty | Empty | Empty | Empty | Empty |
|   | S4<br>DNAp-M004<br>RNAp-M052 | S12<br>DNAp-M012<br>RNAp-M060 | S20<br>DNAp-M020<br>RNAp-M098 | S28<br>DNAp-M028<br>RNAp-M076 | S36<br>DNAp-M036<br>RNAp-M084 | S44<br>DNAp-M044<br>RNAp-M092 | Empty | Empty | Empty | Empty | Empty | Empty |
| C | S5<br>DNAp-M005<br>RNAp-M053 | S13<br>DNAp-M005<br>RNAp-M061 | S21<br>DNAp-M021<br>RNAp-M099 | S29<br>DNAp-M029<br>RNAp-M077 | S37<br>DNAp-M037<br>RNAp-M085 | S45<br>DNAp-M045<br>RNAp-M093 | Empty | Empty | Empty | Empty | Empty | Empty |
|   | S6<br>DNAp-M006<br>RNAp-M054 | S14<br>DNAp-M014<br>RNAp-M062 | S22<br>DNAp-M024<br>RNAp-M100 | S30<br>DNAp-M030<br>RNAp-M078 | S38<br>DNAp-M038<br>RNAp-M088 | S46<br>DNAp-M046<br>RNAp-M094 | Empty | Empty | Empty | Empty | Empty | Empty |
| D | S7<br>DNAp-M007<br>RNAp-M055 | S15<br>DNAp-M015<br>RNAp-M063 | S23<br>DNAp-M023<br>RNAp-M101 | S31<br>DNAp-M031<br>RNAp-M079 | S39<br>DNAp-M039<br>RNAp-M087 | S47<br>DNAp-M047<br>RNAp-M095 | Empty | Empty | Empty | Empty | Empty | Empty |
|   | S8<br>DNAp-M008<br>RNAp-M056 | S16<br>DNAp-M016<br>RNAp-M064 | S24<br>DNAp-M024<br>RNAp-M102 | S32<br>DNAp-M032<br>RNAp-M080 | S40<br>DNAp-M040<br>RNAp-M089 | S48<br>DNAp-M048<br>RNAp-M096 | Empty | Empty | Empty | Empty | Empty | Empty |

**Table 20. Layout of QIAseq Multimodal N7 Plates, 48 reactions Set B**

|          | 1                                    | 2                      | 3                      | 4                      | 5                                    | 6                      | 7     | 8     | 9     | 10    | 11    | 12    |
|----------|--------------------------------------|------------------------|------------------------|------------------------|--------------------------------------|------------------------|-------|-------|-------|-------|-------|-------|
| <b>A</b> | <b>S49</b><br>DNAp-M097<br>RNAp-M146 | DNAp-M105<br>RNAp-M153 | DNAp-M113<br>RNAp-M161 | DNAp-M121<br>RNAp-M169 | <b>S81</b><br>DNAp-M129<br>RNAp-M177 | DNAp-M137<br>RNAp-M185 | Empty | Empty | Empty | Empty | Empty | Empty |
|          | <b>S50</b><br>DNAp-M098<br>RNAp-M146 | DNAp-M106<br>RNAp-M154 | DNAp-M114<br>RNAp-M162 | DNAp-M122<br>RNAp-M170 | <b>S82</b><br>DNAp-M130<br>RNAp-M178 | DNAp-M138<br>RNAp-M186 | Empty | Empty | Empty | Empty | Empty | Empty |
| <b>C</b> | <b>S51</b><br>DNAp-M096<br>RNAp-M147 | DNAp-M107<br>RNAp-M155 | DNAp-M115<br>RNAp-M163 | DNAp-M123<br>RNAp-M171 | <b>S83</b><br>DNAp-M131<br>RNAp-M179 | DNAp-M139<br>RNAp-M187 | Empty | Empty | Empty | Empty | Empty | Empty |
|          | <b>S52</b><br>DNAp-M100<br>RNAp-M148 | DNAp-M108<br>RNAp-M156 | DNAp-M116<br>RNAp-M164 | DNAp-M124<br>RNAp-M172 | <b>S84</b><br>DNAp-M132<br>RNAp-M180 | DNAp-M140<br>RNAp-M188 | Empty | Empty | Empty | Empty | Empty | Empty |
| <b>E</b> | <b>S53</b><br>DNAp-M101<br>RNAp-M149 | DNAp-M109<br>RNAp-M157 | DNAp-M117<br>RNAp-M165 | DNAp-M125<br>RNAp-M173 | <b>S85</b><br>DNAp-M133<br>RNAp-M181 | DNAp-M141<br>RNAp-M189 | Empty | Empty | Empty | Empty | Empty | Empty |
|          | <b>S54</b><br>DNAp-M102<br>RNAp-M150 | DNAp-M110<br>RNAp-M158 | DNAp-M118<br>RNAp-M166 | DNAp-M126<br>RNAp-M174 | <b>S86</b><br>DNAp-M134<br>RNAp-M182 | DNAp-M142<br>RNAp-M190 | Empty | Empty | Empty | Empty | Empty | Empty |
| <b>G</b> | <b>S55</b><br>DNAp-M103<br>RNAp-M151 | DNAp-M111<br>RNAp-M159 | DNAp-M119<br>RNAp-M167 | DNAp-M127<br>RNAp-M175 | <b>S87</b><br>DNAp-M135<br>RNAp-M183 | DNAp-M143<br>RNAp-M191 | Empty | Empty | Empty | Empty | Empty | Empty |
|          | <b>S56</b><br>DNAp-M104<br>RNAp-M152 | DNAp-M112<br>RNAp-M160 | DNAp-M120<br>RNAp-M168 | DNAp-M128<br>RNAp-M176 | <b>S88</b><br>DNAp-M138<br>RNAp-M184 | DNAp-M144<br>RNAp-M192 | Empty | Empty | Empty | Empty | Empty | Empty |

4. Briefly centrifuge, mix by pipetting up and down 8 times, and then briefly centrifuge again.

**Note:** If only a column is used, cut that column from the cuttable plate and proceed to the next step.

5. Program a thermal cycler using the cycling conditions in Table 21 (DNA+RNA primers <1500) or Table 22 (DNA+RNA primers ≥1500).

**Table 21. Cycling conditions for target enrichment if DNA+RNA primers <1500**

| Step                 | Time   | Temperature |
|----------------------|--------|-------------|
| Initial denaturation | 13 min | 95°C        |
|                      | 2 min  | 98°C        |
| 8 cycles             | 15 sec | 98°C        |
|                      | 10 min | 68°C        |
| Hold                 | 5 min  | 72°C        |
|                      | ∞      | 4°C         |

**Table 22. Cycling conditions for target enrichment if number of primers ≥1500/tube**

| Step                 | Time (1500–12,000 primers/tube) | Time (>12,000 primers/tube) | Temperature |
|----------------------|---------------------------------|-----------------------------|-------------|
| Initial denaturation | 13 min                          | 13 min                      | 95°C        |
|                      | 2 min                           | 2 min                       | 98°C        |
| 6 cycles             | 15 s                            | 15 s                        | 98°C        |
|                      | 15 min                          | 30 min                      | 65°C        |
| 1 cycle              | 5 min                           | 5 min                       | 72°C        |
| Hold                 | 5 min                           | 5 min                       | 4°C         |
| Hold                 | ∞                               | ∞                           | 4°C         |

6. Place the target enrichment reaction in the thermal cycler and start the run.

7. Once the run has finished, add 60 µl of Nuclease-free Water to bring each sample to 100 µl.

8. Add 100 µl QIAseq Beads and mix by vortexing or by pipetting up and down several times.

9. Incubate for 5 min at room temperature.

10. Place the tubes/plate on a magnetic rack for 2 min (for tubes) or 10 min (for plates).
11. After the solution has cleared, with the beads still on the magnetic stand, carefully remove and discard the supernatant.

**Important:** Do not discard the beads, because they contain the DNA of interest.

**Important: When DNA+RNA primer  $\geq 12000$ , add 75  $\mu$ l of Nuclease-free Water to resuspend beads, and then add 75  $\mu$ l of QIAseq Bead Binding Buffer.** Mix by vortexing or pipetting up and down. Repeat step 9 to 11.
12. With the beads still on the magnetic stand, add 200  $\mu$ l of 80% ethanol. Rotate the tube 2–3 times or move the plate side-to-side between the 2 column positions of the magnet to wash the beads. Carefully remove and discard the wash.
13. Repeat the ethanol wash.

**Important:** Completely remove all traces of the ethanol after this second wash. Remove the ethanol with a 200  $\mu$ l pipette tip first, spin down briefly, and then use a 10  $\mu$ l pipette tip to remove any residual ethanol.
14. With the beads still on the magnetic stand, air-dry at room temperature for 5–10 min or longer.

**Note:** Visually inspect the pellet to confirm that it is completely dry. Ethanol carryover to the next universal PCR step will affect PCR efficiency.
15. Remove the beads from the magnetic stand, and then elute the DNA from the beads by adding 25  $\mu$ l Nuclease-free Water. Mix well by pipetting.
16. Return the tube/plate to the magnetic rack until the solution has cleared.
17. Transfer 24  $\mu$ l of the supernatant to clean tubes/plate. This will be used in the next 2 protocols.
18. Proceed to “Protocol: qPCR Determination of Universal PCR Cycles”, page 51. Alternatively, the samples can be stored at  $-30$  to  $-15^{\circ}\text{C}$  in a constant-temperature freezer.

# Protocol: qPCR Determination of Universal PCR Cycles

## Important points before starting

- Two microliters of the product from “Protocol: Combined Targeted DNA+RNA Enrichment in a Single Tube”, page 43, or “Appendix B: Separated Targeted DNA and RNA Enrichment in Separate Tubes”, page 74, is the starting material for each of the reaction mixes.
- **Important:** EvaGreen, 20x in water, is required for this procedure and must be purchased from Biotium (cat. no. 31000-T, 31000)
- Set up reactions on ice.
- Do not vortex any reagents or reactions.

## Procedure

1. Prepare the reagents required for the qPCR.
  - 1a. Thaw UPCR Buffer, 5x; DNA qPCR AMP Set; and RNA qPCR AMP. Set at room temperature.
  - 1b. Mix by flicking the tube, and then centrifuge briefly.

**Note:** HotStarTaq DNA Polymerase should be removed from the freezer just before use and placed on ice. After use, immediately return the enzyme to the freezer.
2. Prepare the qPCR reactions according to Table 23 for DNA library or Table 24 for RNA library. Briefly centrifuge, mix by pipetting up and down 10–12 times, and then briefly centrifuge again.

**Note:** If setting up more than one reaction, prepare a volume of master mix 10% greater than what is required for the total number of reactions.

**Table 23. Reaction mix for qPCR of DNA library**

| Component  | Volume/reaction             |
|--|-----------------------------|
| Sample (from "Protocol: Combined Targeted DNA+RNA Enrichment in a Single Tube", page 43)             |                             |
| or   | 2 $\mu$ l                   |
| DNA sample (from "Appendix B: Separated Targeted DNA and RNA Enrichment in Separate Tubes", page 74) |                             |
| UPCR Buffer, 5x  | 2 $\mu$ l                   |
| Nuclease-free Water  | 4.1 $\mu$ l                 |
| DNA qPCR AMP Set   | 1 $\mu$ l                   |
| HotStarTaq DNA Polymerase (6 U/ $\mu$ l)   | 0.4 $\mu$ l                 |
| EvaGreen, 20x in water*  | 0.5 $\mu$ l                 |
| <b>Total</b>   | <b>10 <math>\mu</math>l</b> |

\* Must be purchased from Biotium (cat. no. 31000-T, 31000).

**Table 24. Reaction mix for qPCR of RNA library**

| Component  | Volume/reaction             |
|--|-----------------------------|
| Sample (from "Protocol: Combined Targeted DNA+RNA Enrichment in a Single Tube", page 43)             |                             |
| or   | 2 $\mu$ l                   |
| RNA sample (from "Appendix B: Separated Targeted DNA and RNA Enrichment in Separate Tubes", page 74) |                             |
| UPCR Buffer, 5x  | 2 $\mu$ l                   |
| Nuclease-free Water  | 4.1 $\mu$ l                 |
| RNA qPCR AMP Set   | 1 $\mu$ l                   |
| HotStarTaq DNA Polymerase (6 U/ $\mu$ l)   | 0.4 $\mu$ l                 |
| EvaGreen, 20x in water*  | 0.5 $\mu$ l                 |
| <b>Total</b>   | <b>10 <math>\mu</math>l</b> |

\* Must be purchased from Biotium (cat. no. 31000-T, 31000).

3. Program a qPCR instrument using the cycling conditions in Table 25.

**Note:** No melting curve is required.

**Table 25. Reaction mix for qPCR of RNA library**

| Step                  | Time             | Temperature |
|-----------------------|------------------|-------------|
| <b>Hold</b>           | 13 min           | 95°C        |
|                       | 2 min            | 98°C        |
| <b>2-step cycling</b> |                  |             |
| Denaturation          | 15 s             | 98°C        |
| Annealing/Extension*  | 2 min            | 62°C        |
| <b>Cycle number</b>   | <b>30 cycles</b> |             |
| <b>Hold</b>           | ∞                | 4°C         |

\* Perform fluorescence data collection.

4. Following the reaction, determine the  $C_T$  values. Based on the  $C_T$  values, the number of universal PCR cycles is defined as  $C_T^{(qPCR)}+3$ , for both the DNA and RNA libraries. For example, if the DNA qPCR is  $C_T=19$ , then perform 22 cycles for DNA universal PCR. If the RNA qPCR is  $C_T=15$ , then perform 18 cycles for RNA universal PCR.

Alternative method:

When the run has finished, observe the amplification plot in “Log View” and define the baseline using “auto baseline”. Using the “Log View” of the amplification plot, determine the cycle in which the amplification curve reaches its Plateau Phase, and use 2 cycles fewer. For example, if the plateau phase is reached when the  $C_T$  is 18, then 16 is the required number of universal PCR amplification cycles.

5. Once the amplification cycles for universal PCR have been determined, proceed to “Protocol: Universal PCR”, page 54.

# Protocol: Universal PCR

## Important points before starting

- Nine microliters of the product from “Protocol: Combined Targeted DNA+RNA Enrichment in a Single Tube”, page 43, or “Appendix B: Separated Targeted DNA and RNA Enrichment in Separate Tubes”, page 74, is the starting material for each of the reaction mixes.
- The number of cycles required for amplification are determined in “Protocol: qPCR Determination of Universal PCR Cycles”, page 51.
- Set up reactions on ice.
- Do not vortex any reagents or reactions.
- QIAseq Multimodal S5 Plates are used and come in either 12- or 48-reaction formats allowing the indexing of 12 DNA and RNA samples or 48 DNA and RNA samples, respectively. Each well in a QIAseq Multimodal S5 Plates is single use. SQDIB001 to SQDIB048 and SQDIB097 to SQDIB144 are mixed with universal DNA primer for DNA library amplification. SQDIB0049 to SQDIB096 and SQDIB0145 to SQDIB192 are mixed with universal RNA primer for RNA library amplification. The S5 primers are expected to be used in pairs with N7 primers, with SQDIB001 being paired with DNAP-M001, SQDIB002 being paired with DNAP-M002, etc.; and SQDIB049 being paired with RNAP-M049 and SQDIB050 being paired with RNAP-M050, etc. The plates can be cut in columns to enable indexing of the desired number of samples.
- The final library dual sample index is determined by the combination of the QIAseq Multimodal N7 Plate and the QIAseq Multimodal S5 Plate.
- **Important:** The required combinations of indexes are described in the sequencing sample setup sheets:
  - Sample Sheet Multimodal UDI Set A: [www.qiagen.com/PROM-15281](http://www.qiagen.com/PROM-15281)
  - Sample Sheet Multimodal UDI Set B: [www.qiagen.com/PROM-15282](http://www.qiagen.com/PROM-15282)
  - Sample Sheet Multimodal UDI Set A and Set B: [www.qiagen.com/PROM-15283](http://www.qiagen.com/PROM-15283)

- QIAseq Beads are used for all reaction cleanups.
- **Important:** Prepare fresh 80% ethanol daily.
- Ensure that the QIAseq Beads are thoroughly mixed at all times. This necessitates working quickly and resuspending the beads immediately before use. If a delay in the protocol occurs, simply vortex the beads.

## Procedure

1. Prepare the reagents required for universal PCR.
  - 1a. Thaw UPCR Buffer, 5x, and bring QIAseq Multimodal S5 Plate to room temperature.
  - 1b. Mix by flicking the tube, and then centrifuge briefly.

**Note:** HotStarTaq DNA Polymerase should be removed from the freezer just before use and placed on ice. After use, immediately return the enzyme to the freezer.

2. Prepare the universal PCR according to Table 26 for DNA library or Table 27 for RNA library. Briefly centrifuge, mix by pipetting up and down 10–12 times, and then briefly centrifuge again.

**Note:** If setting up more than one reaction, prepare a volume of master mix 10% greater than what is required for the total number of reactions.

**Table 26. Reaction mix for universal PCR of DNA library**

| Component  | Volume/reaction             |
|--|-----------------------------|
| Sample (from "Protocol: Combined Targeted DNA+RNA Enrichment in a Single Tube", page 43)             |                             |
| or   | 9 $\mu$ l                   |
| DNA sample (from "Appendix B: Separated Targeted DNA and RNA Enrichment in Separate Tubes", page 74) |                             |
| UPCR Buffer, 5x  | 5 $\mu$ l                   |
| Nuclease-free Water  | 10 $\mu$ l                  |
| HotStarTaq DNA Polymerase (6 U/ $\mu$ l)   | 1 $\mu$ l                   |
| <b>Total</b>   | <b>25 <math>\mu</math>l</b> |

**Table 27. Reaction mix for universal PCR of RNA library**

| Component  | Volume/reaction             |
|--|-----------------------------|
| Sample (from "Protocol: Combined Targeted DNA+RNA Enrichment in a Single Tube", page 43)<br>or<br>RNA sample (from "Appendix B: Separated Targeted DNA and RNA Enrichment in Separate Tubes", page 74) | 9 $\mu$ l                   |
| UPCR Buffer, 5x  | 5 $\mu$ l                   |
| Nuclease-free Water  | 10 $\mu$ l                  |
| HotStarTaq DNA Polymerase (6U/ $\mu$ l)  | 1 $\mu$ l                   |
| <b>Total</b>   | <b>25 <math>\mu</math>l</b> |

3. To the QIAseq Multimodal S5 Plate (Table 28, Table 29, or Table 30), add the 25  $\mu$ l reaction mix for universal PCR of DNA library to the DNA wells (columns 1–6), and add the 25  $\mu$ l reaction mix for universal PCR of RNA library to the RNA wells (columns 7–12) of a QIAseq Multimodal S5 Plate.

**Note:** The QIAseq Multimodal S5 plates are cuttable plates that contain predispensed, dried index primers for both DNA and RNA samples in separate wells. The DNA and RNA primers are expected to be used in pairs: for example, sample 1 should use DNA primer SQDIB001 and RNA primer SQDIB049, sample 2 should use DNA primer SQDIB002 and RNA primer SQDIB050, etc.

**Note:** If only a column is used, cut that column from the cuttable plate and proceed to next step.

Table 28. Layout of QIAseq Multimodal S5 Plates, 12 reactions

|   | 1                       | 2                        | 3     | 4     | 5     | 6     | 7                       | 8                        | 9     | 10    | 11    | 12    |
|---|-------------------------|--------------------------|-------|-------|-------|-------|-------------------------|--------------------------|-------|-------|-------|-------|
| A | \$1 SOD1B001<br>uDNA Pr | \$9 SOD1B009<br>uDNA Pr  | Empty | Empty | Empty | Empty | \$1 SOD1B049<br>uRNA Pr | \$9 SOD1B057<br>uRNA Pr  | Empty | Empty | Empty | Empty |
| B | \$2 SOD1B002<br>uDNA Pr | \$10 SOD1B010<br>uDNA Pr | Empty | Empty | Empty | Empty | \$2 SOD1B050<br>uRNA Pr | \$10 SOD1B058<br>uRNA Pr | Empty | Empty | Empty | Empty |
| C | \$3 SOD1B003<br>uDNA Pr | \$11 SOD1B011<br>uDNA Pr | Empty | Empty | Empty | Empty | \$3 SOD1B051<br>uRNA Pr | \$11 SOD1B059<br>uRNA Pr | Empty | Empty | Empty | Empty |
| D | \$4 SOD1B004<br>uDNA Pr | \$12 SOD1B012<br>uDNA Pr | Empty | Empty | Empty | Empty | \$4 SOD1B052<br>uRNA Pr | \$12 SOD1B060<br>uRNA Pr | Empty | Empty | Empty | Empty |
| E | \$5 SOD1B005<br>uDNA Pr | Empty                    | Empty | Empty | Empty | Empty | \$5 SOD1B053<br>uRNA Pr | Empty                    | Empty | Empty | Empty | Empty |
| F | \$6 SOD1B006<br>uDNA Pr | Empty                    | Empty | Empty | Empty | Empty | \$6 SOD1B054<br>uRNA Pr | Empty                    | Empty | Empty | Empty | Empty |
| G | \$7 SOD1B007<br>uDNA Pr | Empty                    | Empty | Empty | Empty | Empty | \$7 SOD1B055<br>uRNA Pr | Empty                    | Empty | Empty | Empty | Empty |
| H | \$8 SOD1B008<br>uDNA Pr | Empty                    | Empty | Empty | Empty | Empty | \$8 SOD1B056<br>uRNA Pr | Empty                    | Empty | Empty | Empty | Empty |

Table 29. Layout of QIAseq Multimodal S5 Plates, 48 reactions Set A

|   | 1                         | 2                          | 3                          | 4                          | 5                          | 6                          | 7                          | 8                          | 9                          | 10                         | 11                         | 12                         |
|---|---------------------------|----------------------------|----------------------------|----------------------------|----------------------------|----------------------------|----------------------------|----------------------------|----------------------------|----------------------------|----------------------------|----------------------------|
| A | S1<br>SQDIB001<br>uDNA Pt | S5<br>SQDIB008<br>uDNA Pt  | S17<br>SQDIB017<br>uDNA Pt | S25<br>SQDIB025<br>uDNA Pt | S33<br>SQDIB033<br>uDNA Pt | S41<br>SQDIB041<br>uDNA Pt | S49<br>SQDIB049<br>uDNA Pt | S57<br>SQDIB057<br>uDNA Pt | S65<br>SQDIB065<br>uDNA Pt | S73<br>SQDIB073<br>uDNA Pt | S81<br>SQDIB081<br>uDNA Pt | S89<br>SQDIB089<br>uDNA Pt |
| B | S2<br>SQDIB002<br>uDNA Pt | S6<br>SQDIB009<br>uDNA Pt  | S18<br>SQDIB018<br>uDNA Pt | S26<br>SQDIB026<br>uDNA Pt | S34<br>SQDIB034<br>uDNA Pt | S42<br>SQDIB042<br>uDNA Pt | S50<br>SQDIB050<br>uDNA Pt | S58<br>SQDIB058<br>uDNA Pt | S66<br>SQDIB066<br>uDNA Pt | S74<br>SQDIB074<br>uDNA Pt | S82<br>SQDIB082<br>uDNA Pt | S90<br>SQDIB090<br>uDNA Pt |
| C | S3<br>SQDIB003<br>uDNA Pt | S7<br>SQDIB010<br>uDNA Pt  | S19<br>SQDIB019<br>uDNA Pt | S27<br>SQDIB027<br>uDNA Pt | S35<br>SQDIB035<br>uDNA Pt | S43<br>SQDIB043<br>uDNA Pt | S51<br>SQDIB051<br>uDNA Pt | S59<br>SQDIB059<br>uDNA Pt | S67<br>SQDIB067<br>uDNA Pt | S75<br>SQDIB075<br>uDNA Pt | S83<br>SQDIB083<br>uDNA Pt | S91<br>SQDIB091<br>uDNA Pt |
| D | S4<br>SQDIB004<br>uDNA Pt | S8<br>SQDIB011<br>uDNA Pt  | S20<br>SQDIB020<br>uDNA Pt | S28<br>SQDIB028<br>uDNA Pt | S36<br>SQDIB036<br>uDNA Pt | S44<br>SQDIB044<br>uDNA Pt | S52<br>SQDIB052<br>uDNA Pt | S60<br>SQDIB060<br>uDNA Pt | S68<br>SQDIB068<br>uDNA Pt | S76<br>SQDIB076<br>uDNA Pt | S84<br>SQDIB084<br>uDNA Pt | S92<br>SQDIB092<br>uDNA Pt |
| E | S5<br>SQDIB005<br>uDNA Pt | S9<br>SQDIB012<br>uDNA Pt  | S21<br>SQDIB021<br>uDNA Pt | S29<br>SQDIB029<br>uDNA Pt | S37<br>SQDIB037<br>uDNA Pt | S45<br>SQDIB045<br>uDNA Pt | S53<br>SQDIB053<br>uDNA Pt | S61<br>SQDIB061<br>uDNA Pt | S69<br>SQDIB069<br>uDNA Pt | S77<br>SQDIB077<br>uDNA Pt | S85<br>SQDIB085<br>uDNA Pt | S93<br>SQDIB093<br>uDNA Pt |
| F | S6<br>SQDIB006<br>uDNA Pt | S10<br>SQDIB013<br>uDNA Pt | S22<br>SQDIB022<br>uDNA Pt | S30<br>SQDIB030<br>uDNA Pt | S38<br>SQDIB038<br>uDNA Pt | S46<br>SQDIB046<br>uDNA Pt | S54<br>SQDIB054<br>uDNA Pt | S62<br>SQDIB062<br>uDNA Pt | S70<br>SQDIB070<br>uDNA Pt | S78<br>SQDIB078<br>uDNA Pt | S86<br>SQDIB086<br>uDNA Pt | S94<br>SQDIB094<br>uDNA Pt |
| G | S7<br>SQDIB007<br>uDNA Pt | S11<br>SQDIB014<br>uDNA Pt | S23<br>SQDIB023<br>uDNA Pt | S31<br>SQDIB031<br>uDNA Pt | S39<br>SQDIB039<br>uDNA Pt | S47<br>SQDIB047<br>uDNA Pt | S55<br>SQDIB055<br>uDNA Pt | S63<br>SQDIB063<br>uDNA Pt | S71<br>SQDIB071<br>uDNA Pt | S79<br>SQDIB079<br>uDNA Pt | S87<br>SQDIB087<br>uDNA Pt | S95<br>SQDIB095<br>uDNA Pt |
| H | S8<br>SQDIB008<br>uDNA Pt | S12<br>SQDIB015<br>uDNA Pt | S24<br>SQDIB024<br>uDNA Pt | S32<br>SQDIB032<br>uDNA Pt | S40<br>SQDIB040<br>uDNA Pt | S48<br>SQDIB048<br>uDNA Pt | S56<br>SQDIB056<br>uDNA Pt | S64<br>SQDIB064<br>uDNA Pt | S72<br>SQDIB072<br>uDNA Pt | S80<br>SQDIB080<br>uDNA Pt | S88<br>SQDIB088<br>uDNA Pt | S96<br>SQDIB096<br>uDNA Pt |

Table 30. Layout of QIAseq Multimodal S5 Plates, 48 reactions Set B

|   | 1                   | 2                   | 3                   | 4                   | 5                   | 6                   | 7                   | 8                   | 9                   | 10                  | 11                  | 12                  |
|---|---------------------|---------------------|---------------------|---------------------|---------------------|---------------------|---------------------|---------------------|---------------------|---------------------|---------------------|---------------------|
| A | SQDIB097<br>uDNA Pr | SQDIB105<br>uDNA Pr | SQDIB113<br>uDNA Pr | SQDIB121<br>uDNA Pr | SQDIB129<br>uDNA Pr | SQDIB137<br>uDNA Pr | SQDIB145<br>uRNA Pr | SQDIB153<br>uRNA Pr | SQDIB161<br>uRNA Pr | SQDIB169<br>uRNA Pr | SQDIB177<br>uRNA Pr | SQDIB185<br>uRNA Pr |
|   | S89                 | S85                 | S85                 | S73                 | S81                 | S88                 | S46                 | S57                 | S85                 | S73                 | S81                 | S89                 |
| B | SQDIB098<br>uDNA Pr | SQDIB106<br>uDNA Pr | SQDIB114<br>uDNA Pr | SQDIB122<br>uDNA Pr | SQDIB130<br>uDNA Pr | SQDIB138<br>uDNA Pr | SQDIB146<br>uRNA Pr | SQDIB154<br>uRNA Pr | SQDIB162<br>uRNA Pr | SQDIB170<br>uRNA Pr | SQDIB178<br>uRNA Pr | SQDIB186<br>uRNA Pr |
|   | S90                 | S86                 | S86                 | S74                 | S82                 | S90                 | S50                 | S38                 | S66                 | S74                 | S82                 | S90                 |
| C | SQDIB099<br>uDNA Pr | SQDIB107<br>uDNA Pr | SQDIB115<br>uDNA Pr | SQDIB123<br>uDNA Pr | SQDIB131<br>uDNA Pr | SQDIB139<br>uDNA Pr | SQDIB147<br>uRNA Pr | SQDIB155<br>uRNA Pr | SQDIB163<br>uRNA Pr | SQDIB171<br>uRNA Pr | SQDIB179<br>uRNA Pr | SQDIB187<br>uRNA Pr |
|   | S91                 | S87                 | S87                 | S75                 | S83                 | S81                 | S85                 | S87                 | S87                 | S76                 | S83                 | S81                 |
| D | SQDIB100<br>uDNA Pr | SQDIB108<br>uDNA Pr | SQDIB116<br>uDNA Pr | SQDIB124<br>uDNA Pr | SQDIB132<br>uDNA Pr | SQDIB140<br>uDNA Pr | SQDIB148<br>uRNA Pr | SQDIB156<br>uRNA Pr | SQDIB164<br>uRNA Pr | SQDIB172<br>uRNA Pr | SQDIB180<br>uRNA Pr | SQDIB188<br>uRNA Pr |
|   | S92                 | S88                 | S88                 | S76                 | S84                 | S92                 | S52                 | S90                 | S89                 | S76                 | S84                 | S92                 |
| E | SQDIB101<br>uDNA Pr | SQDIB109<br>uDNA Pr | SQDIB117<br>uDNA Pr | SQDIB125<br>uDNA Pr | SQDIB133<br>uDNA Pr | SQDIB141<br>uDNA Pr | SQDIB149<br>uRNA Pr | SQDIB157<br>uRNA Pr | SQDIB165<br>uRNA Pr | SQDIB173<br>uRNA Pr | SQDIB181<br>uRNA Pr | SQDIB189<br>uRNA Pr |
|   | S93                 | S89                 | S89                 | S77                 | S85                 | S93                 | S53                 | S61                 | S69                 | S77                 | S85                 | S93                 |
| F | SQDIB102<br>uDNA Pr | SQDIB110<br>uDNA Pr | SQDIB118<br>uDNA Pr | SQDIB126<br>uDNA Pr | SQDIB134<br>uDNA Pr | SQDIB142<br>uDNA Pr | SQDIB150<br>uRNA Pr | SQDIB158<br>uRNA Pr | SQDIB166<br>uRNA Pr | SQDIB174<br>uRNA Pr | SQDIB182<br>uRNA Pr | SQDIB190<br>uRNA Pr |
|   | S94                 | S82                 | S82                 | S78                 | S86                 | S94                 | S54                 | S63                 | S70                 | S78                 | S86                 | S94                 |
| G | SQDIB103<br>uDNA Pr | SQDIB111<br>uDNA Pr | SQDIB119<br>uDNA Pr | SQDIB127<br>uDNA Pr | SQDIB135<br>uDNA Pr | SQDIB143<br>uDNA Pr | SQDIB151<br>uRNA Pr | SQDIB159<br>uRNA Pr | SQDIB167<br>uRNA Pr | SQDIB175<br>uRNA Pr | SQDIB183<br>uRNA Pr | SQDIB191<br>uRNA Pr |
|   | S95                 | S83                 | S83                 | S79                 | S87                 | S95                 | S55                 | S63                 | S71                 | S79                 | S87                 | S95                 |
| H | SQDIB104<br>uDNA Pr | SQDIB112<br>uDNA Pr | SQDIB120<br>uDNA Pr | SQDIB128<br>uDNA Pr | SQDIB136<br>uDNA Pr | SQDIB144<br>uDNA Pr | SQDIB152<br>uRNA Pr | SQDIB160<br>uRNA Pr | SQDIB168<br>uRNA Pr | SQDIB176<br>uRNA Pr | SQDIB184<br>uRNA Pr | SQDIB192<br>uRNA Pr |
|   | S96                 | S84                 | S84                 | S80                 | S88                 | S96                 | S56                 | S64                 | S72                 | S80                 | S88                 | S96                 |

4. Program a thermal cycler using the cycling conditions in Table 31.

**Table 31. Cycling conditions for universal PCR**

| Step                  | Time  | Temperature |
|-----------------------|---|-------------|
| <b>Hold</b>           | 13 min  | 95°C        |
|                       | 2 min   | 98°C        |
| <b>2-step cycling</b> |   |             |
| Denaturation          | 15 s  | 98°C        |
| Annealing/Extension   | 2 min   | 62°C        |
| <b>Cycle number</b>   | <b>Based on "Protocol: qPCR Determination of Universal PCR Cycles", page 51</b> |             |
| <b>Hold</b>           | ∞   | 4°C         |

5. After the reaction is complete, add 75 µl of Nuclease-free Water to bring each sample to 100 µl.
6. Add 100 µl QIAseq Beads, and then mix by vortexing or pipetting up and down several times.
7. Incubate for 5 min at room temperature.
8. Place the tubes/plate on a magnetic rack for 2 min (for tubes) or 10 min (for plates) to separate the beads from the supernatant.
9. After the solution has cleared, with the beads still on the magnetic stand, carefully remove and discard the supernatant.

**Important:** Do not discard the beads, because they contain the DNA of interest.

10. With the beads still on the magnetic stand, add 200 µl of 80% ethanol. Rotate the tube 2–3 times or move the plate side-to-side between the 2 column positions of the magnet to wash the beads. Carefully remove and discard the wash.
11. Repeat the ethanol wash.


**Important:** Completely remove all traces of the ethanol after this second wash. Remove the ethanol with a 200 µl pipette tip first, spin down briefly, and then use a 10 µl pipette tip to remove any residual ethanol.

- 
12. With the beads still on the magnetic stand, air-dry at room temperature for 5–10 min.  
**Note:** Visually inspect the pellet to confirm that it is completely dry. If the pellet is not completely dry, additional drying time may be required.
  13. Remove the beads from the magnetic stand, and then elute the DNA from the beads by adding 20  $\mu$ l Nuclease-free Water. Mix well by pipetting.
  14. Return the tube/plate to the magnetic rack until the solution has cleared.
  15. Transfer 18  $\mu$ l of the supernatant to clean tubes/plate.
  16. Proceed to “Recommendations: Library QC and Quantification”, page 62. Alternatively, the samples can be stored at  $-30$  to  $-15^{\circ}\text{C}$  in a constant-temperature freezer.

# Recommendations: Library QC and Quantification

## NGS library QC

QC can be performed with the Agilent Bioanalyzer or TapeStation. Check for the correct size distribution of library fragments (~400–500 bp median size) and for the absence of adapters or adapter–dimers (~<200 bp) (Figure 4).



**Figure 4.** QIAseq Multimodal Targeted DNA (left) and Targeted RNA (right) libraries.

## Preferred library quantification method

The library yield measurements from the Bioanalyzer or TapeStation® rely on fluorescence dyes that intercalate into DNA or RNA. These dyes cannot discriminate between molecules with or without adapter sequences, yet only complete QIAseq Multimodal libraries with full adapter sequences will be sequenced. As a result, QIAGEN's QIAseq Library Quant Array Kit (cat. no. 333304) or QIAseq Library Quant Assay Kit (cat. no. 333314), which contains laboratory-verified forward and reverse primers, together with a DNA standard, is highly recommended for accurate quantification of the prepared QIAseq Multimodal library.

---

# Protocol: Sequencing Setup on Illumina MiSeq and NextSeq


## Important points before starting

- **Important:** Recommendations for library dilution concentrations and library loading concentrations are based on QIAseq Library Quant System.
- **Important:** QIAseq A Read 1 Primer I (Custom Read 1 Sequencing Primer) and Multimodal Read 2 Primer (Custom Read 2 Sequencing Primer) **must** be used when performing sequencing on Illumina platform.
- **Important:** QIAseq A Read 1 Primer I (Custom Read 1 Sequencing Primer) goes into the following specific reagent cartridge positions:
  - MiSeq position #18
  - NextSeq position #7
  - **Important:** Multimodal Read 2 Primer (the Custom Read 2 Sequencing Primer) goes into the following specific reagent cartridge positions:
    - MiSeq position #20
    - NextSeq position # 8
- **Important:** Paired-end sequencing **should** be used for QIAseq Multimodal libraries on Illumina platforms.
  - Read 1: 149 bp
  - Read 2: 149 bp
  - Custom Index 1: 10 bp
  - Custom Index 2: 10 bp
- For complete instructions on how to denature sequencing libraries, prepare custom index primers, and set up a sequencing run, please refer to the system-specific Illumina documents.
- Instrument-specific imagery is included to aid in sequencing preparations.

## Sequencing preparations for MiSeq

1. Download the appropriate template from the “Resource” tab of the QIAseq Multimodal Panel:
  - Sample Sheet Multimodal UDI Set A (also used for Index 1–12):  
**[www.qiagen.com/PROM-15281](http://www.qiagen.com/PROM-15281)**
  - Sample Sheet Multimodal UDI Set B: **[www.qiagen.com/PROM-15282](http://www.qiagen.com/PROM-15282)**
  - Sample Sheet Multimodal UDI Set A and Set B: **[www.qiagen.com/PROM-15283](http://www.qiagen.com/PROM-15283)**
2. On the template:
  - 2a. Modify **Investigator Name**, **Date**, **Sample\_ID**, and **Sample Name**.  
**Important:** We recommend adding **-DNA** in the Sample name of a DNA library and **-RNA** for an RNA library, to allow automatic parsing of the DNA and RNA libraries during data analysis. If the libraries are not labeled, they must be manually parsed into either the DNA or RNA box.
  - 2b. Delete any unused index pairs and save the sample sheet for uploading.
  - 2c. Read 1 is **149 bp**, Read 2 is **149 bp**, and each Index Read is **10 bp**.
3. **Sample dilution and pooling:** Dilute libraries to 2 or 4 nM for MiSeq. Then combine libraries with different sample indexes in equimolar amounts if similar sequencing depth is needed for each library.  
**Note:** If combining libraries with the same number of primers, pool equal volumes of the individual libraries at 4 nM together. If the libraries have different primer numbers, then combine the libraries at volume ratios according to their number of primers. For example, Library A has 5000 primers at 4 nM, and Library B has 600 primers at 4 nM; combining 50 µl Library A with 6 µl Library B will result in similar coverage depth for both libraries A and B in the same sequencing run.
4. **Library preparation and loading:** Prepare and load library on a MiSeq according to the *MiSeq System Denature and Dilute Libraries Guide*. The final library concentration is 10–12 pM on MiSeq.  
**Note:** Recommendations for library loading concentrations are based on the QIAseq Library Quant System.

5. **Custom sequencing primer for Read 1 preparation and loading:** Use 597  $\mu$ l HT1 (hybridization buffer) from the Illumina Sequencing Kit to dilute 3  $\mu$ l of QIAseq A Read 1 Primer I (provided) to obtain a final concentration of 0.5  $\mu$ M. Load 600  $\mu$ l of the diluted QIAseq A Read 1 Primer I to position 18 of the MiSeq reagent cartridge (Figure 5). For more details, please refer to Illumina’s *MiSeq System: Custom Primers Guide*.
6. **Custom sequencing primer for Read 2 preparation and loading:** Use 597  $\mu$ l HT1 (hybridization buffer) from the Illumina Sequencing Kit to dilute 3  $\mu$ l of Multimodal Read 2 Primer (provided) to obtain a final concentration of 0.5  $\mu$ M. Load 600  $\mu$ l of the diluted QIAseq Read 2 Primer to position 20 of the MiSeq reagent cartridge (Figure 5). For more details, please refer to Illumina’s *MiSeq System: Custom Primers Guide*.




**Figure 5. MiSeq reagent cartridge.** A: position 18 for Read 1 Custom Primer; C: position 20 for Read 2 Custom Primer.

7. Upon completion of the sequencing run, proceed to “Protocol: Data Analysis Using the QIAseq Multimodal Data Analysis Portal”, page 68.

## Sequencing preparations for NextSeq

- 1. Sample dilution and pooling:** Dilute libraries to 0.5, 1, 2, or 4 nM for NextSeq. Then combine libraries with different sample indexes in equimolar amounts if similar sequencing depth is needed for each library.  
**Note:** Recommendations for library dilution concentrations are based on the QIAseq Library Quant System.  
**Note:** If combining libraries with the same number of primers, pool equal volumes of the individual libraries at 4 nM together. If the libraries have different primer numbers, then combine the libraries at volume ratios according to their number of primers. For example, Library A has 5000 primers at 4 nM, and Library B has 600 primers at 4 nM; combining 50  $\mu$ l Library A with 6  $\mu$ l Library B will result in similar coverage depth for both libraries A and B in the same sequencing run.
- 2. Library preparation and loading:** Prepare and load library onto a NextSeq according to the *NextSeq System Denature and Dilute Libraries Guide*. The final library concentration is 1.2–1.5 pM on NextSeq.  
**Note:** Recommendations for library loading concentrations are based on the QIAseq Library Quant System.
- 3. Custom sequencing primer for Read 1 preparation and loading:** Use 1994  $\mu$ l HT1 (hybridization buffer) from the Illumina Sequencing Kit to dilute 6  $\mu$ l of QIAseq A Read 1 Primer I (provided) to obtain a final concentration of 0.3  $\mu$ M. Load 2 ml of the diluted QIAseq A Read 1 Primer I to position 7 of the NextSeq reagent cartridge (Figure 6).  
**Note:** All other steps refer to run setup workflow as described in the NextSeq 500 System Guide (part # 15046563) or NextSeq 550 System Guide (part # 15069765-02).

- Custom sequencing primer for Read 2 preparation and loading:** Use 1994  $\mu$ l HT1 (hybridization buffer) from the Illumina Sequencing Kit to dilute 6  $\mu$ l of Multimodal Read 2 Primer to obtain a final concentration of 0.3  $\mu$ M. Load 2 ml of the diluted Multimodal Read 2 Primer to position 8 of the NextSeq reagent cartridge (Figure 6).



| Position # | Custom Primer        |
|------------|----------------------|
| 7          | Custom Read 1 primer |
| 8          | Custom Read 2 primer |


**Figure 6. NextSeq reagent cartridge.**

- When working with the QIAseq Multimodal custom UDIs, use Local Run Manager (LRM) V2 on the instrument to upload sample sheet (see page 66 for downloading the appropriate template and modifying the template) and proceed with sequencing: Read 1 is 149 bp, Read 2 is 149 bp, and each Index Read is 10 bp.
- Upon completion, proceed to “Protocol: Data Analysis Using the QIAseq Multimodal Data Analysis Portal”, page 68.

# Protocol: Data Analysis Using the QIAseq Multimodal Data Analysis Portal

After sequencing, the results can be analyzed using QIAGEN's QIAseq Multimodal Data Analysis Portal. Our data analysis pipeline will perform mapping to the reference genome, UMI counting, read trimming (removing primer sequences), and variant identification. All detected variants can be further interpreted using QCIH.

1. Go to the QIAseq Multimodal Data Analysis Portal, [ngsdataanalysis2.qiagen.com/MultiModal/](https://ngsdataanalysis2.qiagen.com/MultiModal/)
2. Log in to the portal.
3. In the **Read Files** tab, select **BaseSpace** to upload files from BaseSpace, or select **Uploaded** > **Upload New Files** to upload files from your local hard drive.



**Note:** All files that have been uploaded to the portal are listed under the **Read Files** tab. Here it is possible to delete files that are no longer needed and share files with collaborators.

4. Select the boxes next to the files that will be analyzed, and then click **Select For Analysis**.

1 - 4 of 4 files    50 per page    [Upload New Files](#)    [Delete](#)    [Share](#)    [Refresh](#)    [Select For Analysis](#)

| File Name                          | File Size            | Uploaded At          | Status               | <input type="checkbox"/>            |
|------------------------------------|----------------------|----------------------|----------------------|-------------------------------------|
| <input type="text"/>               | <input type="text"/> | <input type="text"/> | <input type="text"/> | <input type="text"/>                |
| rna-265-1_S11_L001_R2_001.fastq.gz | 23.63 MB             | 2019/10/18 05:34:59  | Ready                | <input checked="" type="checkbox"/> |
| rna-265-2_S12_L001_R2_001.fastq.gz | 27.75 MB             | 2019/10/18 05:34:59  | Ready                | <input checked="" type="checkbox"/> |
| rna-265-1_S11_L001_R1_001.fastq.gz | 51.98 MB             | 2019/10/18 05:34:59  | Ready                | <input type="checkbox"/>            |
| rna-265-2_S12_L001_R1_001.fastq.gz | 60.49 MB             | 2019/10/18 05:34:59  | Ready                | <input type="checkbox"/>            |

5. Under the **Analysis Jobs** tab, configure the analysis per the drop-down menus and check boxes as described below:

[Read Files](#)    [Analysis Jobs](#)

[View All Jobs](#)

Read Files \*    2 read files selected    JSON Request   

Job Title \*   

Job Type   

Catalog # \*         Custom Catalog #

File Lanes        ⓘ

DNA smCounter Version        ⓘ

Copy Number Reference Job IDs        ⓘ

RNAseq Analysis Modes     Fusion Calling     SNP/Indel Variant Calling

[Match DNA RNA Files](#)    [ANALYZE](#)

- **Read Files:** Verify that the correct read files have been selected.
  - **Job Title:** Enter a title for the analysis job.
  - **Job Type:** Select **DNA & RNA**, **DNA Only**, or **RNA only**.
  - **Catalog #:** If using a catalog panel, select the number from the dropdown menu. If using a custom panel, enter the custom catalog number manually.
  - **File Lanes:** Choose **1-lane** if MiSeq/HiSeq/NovaSeq was used or if NextSeq was used with concatenated files. Choose **4-lane** if NextSeq was used without ultimately having the files concatenated.
  - **DNA smCounter Version:** Select the appropriate version, based on the experimental needs. For guidance, select **Information** .
  - **Copy Number Reference Job IDs:** For copy number analysis, normal sample(s) need to be analyzed with the portal before case samples are set up. Enter the job ID corresponding to your control samples for copy number analysis.
  - **RNAscan Analysis Modes:** Select **Fusion Calling** and/or **SNP/Indel Variant Calling**. By default, gene expression analysis is provided with each selection.
6. Click **Match DNA RNA Files** to manually drag the selected read files into either the **DNA Files** box or the **RNA Files** box.
- Note:** We recommend that you include **-DNA** at the end of the sample name for a DNA library and **-RNA** for an RNA library, to allow automatic parsing of the DNA and RNA libraries during data analysis.
7. Click **ANALYZE**. The analysis job status changes from “Queued” to “In progress”, and then to “Done successfully”.
8. Once the analysis is completed, output files can be downloaded by clicking **Download**.
- Note:** Ultimately, detected variants can be interpreted with QCI-I.

# Troubleshooting Guide

This troubleshooting guide may be helpful in solving any problems that may arise. For more information, see also the Frequently Asked Questions page at our Technical Support Center: [www.qiagen.com/FAQ/FAQList.aspx](http://www.qiagen.com/FAQ/FAQList.aspx). The scientists in QIAGEN Technical Services are always happy to answer any questions you may have about either the information and protocols in this handbook or sample and assay technologies (for contact information, visit [www.qiagen.com](http://www.qiagen.com)).

---

## Comments and suggestions

---

### Low library yield

- |   |  |
|---|--|
| a) Suboptimal reaction conditions due to low sample quality | Make sure to use high-quality samples to ensure optimal activity of the library enzymes.   |
| b) Inefficient targeted enrichment or universal PCR         | QIaseq beads need to be completely dried before elution. Ethanol carryover to targeted enrichment and universal PCR will affect PCR reaction efficiency. |

### Unexpected signal peaks

- |   |  |
|---|--|
| a) Short peaks <200 bp                      | These are primer–dimers from targeted enrichment or universal PCR (<200 bp). The presence of primer–dimers indicates either not enough DNA/RNA input or inefficient PCR reactions or handling issues with bead purifications.  |
| b) Larger DNA fragments after universal PCR | After the universal PCR, library fragments are larger than the intended peak and can be a PCR artifact due to overamplification of the DNA library. Overamplification of the library won't affect the sequencing performance. Decreasing the number of universal PCR cycle numbers can reduce overamplification. |

### Sequencing issues

- |  |   |
|--|---|
| a) Too low or too high cluster density | Accurate library quantification is the key for optimal cluster density on any sequencing instrument. PCR-based quantification method is recommended. Other methods may lead to incorrect quantification of the library, especially when there is overamplification.   |
| b) Very low clusters passing filter    | Make sure the library is accurately quantified and that the correct amount is loaded onto the sequencing instrument. In addition, the QIaseq A Read 1 Primer I (100 µM) Custom Read 1 Sequencing Primer and Custom Multimodal Read 2 Primer (100uM) <b>must</b> be used when sequencing on any Illumina platform. |

### Variant detection issues

- |                             |  |
|-----------------------------|--|
| Known variants not detected | Variant detection sensitivity is directly related to the input DNA and read depth. Check Table 3 (page 21), Table 4 (page 22), and Table 5 (page 23) to determine if the required input DNA, UMI numbers, and read depth are met for the specific variant detection application. |
|-----------------------------|--|

## Appendix A: FFPE Sample Quality and Quantity

High-sensitivity Quant-iT dsDNA Assay Kit (Thermo Fisher Scientific, cat. no. Q33120) is recommended for quantification of DNA from FFPE total nucleic acid samples. This kit or other methods like Nanodrop® or QIAxperi® can be used for quantification of FFPE DNA (in separate elutes with RNA).

Genomic DNA present in FFPE archives is usually damaged and fragmented to an uncertain extent. The QIAseq DNA QuantiMIZE System is a qPCR-based approach that determines the quantity and quality of the DNA amenable to PCR-based targeted enrichment prior to NGS. The system provides a cost-effective approach to qualify and quantify the DNA isolated from biological samples – mainly for FFPE samples. Please refer to the corresponding handbook for determining FFPE DNA quantity and quality with the QIAseq DNA QuantiMIZE System.

FFPE DNA input can be determined by the following: If FFPE DNA is defined as high quality (quality control [QC] score  $\leq 0.04$ ) by QuantiMIZE, then up to 100 ng of DNA can be used. If the DNA is determined as low quality (QC score  $>0.04$ ) then up to 250 ng of DNA can be used. The QC score of QuantiMIZE reflects the amount of amplifiable DNA present in the sample, therefore correlating with the number of UMIs that can be sequenced in the library (Figure 7).




Figure 7. Correlation between QIAseq QuantiMIZE QC Score and the number of UMIs.

---

Compared to the same amount of fresh DNA, only 10–50% of UMIs can be captured from FFPE DNA, depending on the quality. This is due to a lower amplifiable DNA amount present in the FFPE samples. Therefore, a higher input amount is recommended for FFPE DNA samples to ensure that enough UMIs can be sequenced for variant detection.

However, if the quality of the FFPE DNA is not assessed by QIAseq QuantiMIZE kits, up to 100 ng can be used. If the FFPE DNA quality is high, an input of more than 100 ng will potentially overload the QIAseq Multimodal Panels system.

---

# Appendix B: Separated Targeted DNA and RNA Enrichment in Separate Tubes

## Important points before starting

- Two 10.2 µl aliquots of the product from “Protocol: Reverse Transcription”, page 40, is the starting material for this protocol.
- Set up reactions on ice.
- Do not vortex any reagents or reactions.
- The final library dual sample index is determined by the combination of the QIAseq Multimodal N7 Plate SW and the QIAseq Multimodal S5 Plate. QIAseq Beads are used for all reaction cleanups.

**Important:** The required combinations of indexes are described in the sequencing sample setup sheets:

- Sample Sheet Multimodal UDI Set A: [www.qiagen.com/PROM-15281](http://www.qiagen.com/PROM-15281)
- Sample Sheet Multimodal UDI Set B: [www.qiagen.com/PROM-15282](http://www.qiagen.com/PROM-15282)
- Sample Sheet Multimodal UDI Set A and Set B: [www.qiagen.com/PROM-15283](http://www.qiagen.com/PROM-15283)
- **Important:** Prepare fresh 80% ethanol daily.
- Ensure that the QIAseq Beads are thoroughly mixed at all times. This necessitates working quickly and resuspending the beads immediately before use. If a delay in the protocol occurs, simply vortex the beads.

**IMPORTANT:** To use this protocol, one of the following is required:

- **QIAseq Multimodal Index I Set A SW (96)** (cat. no. 333985)
- **QIAseq Multimodal Index I Set B SW (96)** (cat. no. 333995)

These plates come in 48-reaction formats allowing the indexing of 48 DNA and 48 RNA libraries using one plate of Set A or Set B. In each indicated well of the cuttable plate, there are dried N7 index primers for either DNA or RNA. The plates can be cut in columns to enable indexing of the desired number of samples.

Two plates of Set A or Set B are included for each index kit for making a total of 96 DNA and 96 RNA libraries. By combining Set A and Set B, up to 96 DNA and 96 RNA libraries can be multiplexed. Each well in the plate is single use.

## Procedure

1. Prepare the reagents required for target enrichment.
  - 1a. Thaw TEPCR Buffer, 5x; Multimodal DHS Panel (DNA); and Multimodal VHS Panel (RNA); and bring QIAseq Multimodal N7 SW Plate to room temperature.
  - 1b. Mix by flicking the tube, and then centrifuge briefly.

**Note:** HotStarTaq DNA Polymerase should be removed from the freezer just before use and placed on ice. After use, immediately return the enzyme to the freezer.

2. Prepare the target enrichment mix according to Table 32. Briefly centrifuge, mix by pipetting up and down 10–12 times, and then briefly centrifuge again.

**Note:** If setting up more than one reaction, prepare a volume of master mix 10% greater than what is required for the total number of reactions.

**Table 32. Reaction mix for target enrichment**

| Component  | Separate DNA                | Separate RNA                |
|--|-----------------------------|-----------------------------|
| Sample (from "Protocol: Reverse Transcription", page 40) | 10.2 $\mu$ l                | 10.2 $\mu$ l                |
| TEPCR Buffer, 5x   | 4 $\mu$ l                   | 4 $\mu$ l                   |
| Multimodal DHS Panel (DNA)                               | 5 $\mu$ l                   | 0 $\mu$ l                   |
| Multimodal VHS Panel (RNA)                               | 0 $\mu$ l                   | 4 $\mu$ l                   |
| HotStarTaq Polymerase (6 U/ $\mu$ l)                     | 0.8 $\mu$ l                 | 0.8 $\mu$ l                 |
| Nuclease-free Water                                      | 0 $\mu$ l                   | 1 $\mu$ l                   |
| <b>Total</b>   | <b>20 <math>\mu</math>l</b> | <b>20 <math>\mu</math>l</b> |

3. Add the 20  $\mu$ l target enrichment reaction mix into a well of a QIAseq Multimodal N7 Plate SW (Table 33 and Table 34), which are cuttable plates that contain predispensed, dried N7 index primers for either DNA (columns 1–6) or RNA (columns 7–12).

**Note:** The plates can be cut in columns to enable indexing of the desired number of samples.

**Important:** Put the unused plate in the foil back and keep in  $-20^{\circ}\text{C}$  for long-term storage.

**Important:** Index primers for DNA and RNA should be used in pairs. For example: use well A1 for DNA and well A7 for RNA for sample 1 (S1); use well B1 for DNA and well B7 for RNA for sample 2 (S2), and so on.

**Table 33. Layout of QIAseq Multimodal N7 Plates SW, 48 reactions Set A**

|   | 1                | 2                | 3                | 4                | 5                | 6               | 7                | 8                | 9                | 10               | 11               | 12               |
|---|------------------|------------------|------------------|------------------|------------------|-----------------|------------------|------------------|------------------|------------------|------------------|------------------|
| A | DNAG-M001<br>S8  | DNAG-M005<br>S17 | DNAG-M017<br>S25 | DNAG-M025<br>S33 | DNAG-M033<br>S41 | DNAG-M041<br>S1 | RNAc-M049<br>S9  | RNAc-M057<br>S17 | RNAc-M065<br>S25 | RNAc-M073<br>S33 | RNAc-M081<br>S41 | RNAc-M089<br>S49 |
| B | DNAG-M002<br>S9  | DNAG-M006<br>S18 | DNAG-M018<br>S26 | DNAG-M026<br>S34 | DNAG-M034<br>S42 | DNAG-M042<br>S2 | RNAc-M050<br>S10 | RNAc-M058<br>S18 | RNAc-M066<br>S26 | RNAc-M074<br>S34 | RNAc-M082<br>S42 | RNAc-M090<br>S50 |
| C | DNAG-M003<br>S11 | DNAG-M007<br>S19 | DNAG-M019<br>S27 | DNAG-M027<br>S35 | DNAG-M035<br>S43 | DNAG-M043<br>S3 | RNAc-M051<br>S11 | RNAc-M059<br>S19 | RNAc-M067<br>S27 | RNAc-M075<br>S35 | RNAc-M083<br>S43 | RNAc-M091<br>S51 |
| D | DNAG-M004<br>S12 | DNAG-M008<br>S20 | DNAG-M020<br>S28 | DNAG-M028<br>S36 | DNAG-M036<br>S44 | DNAG-M044<br>S4 | RNAc-M052<br>S12 | RNAc-M060<br>S20 | RNAc-M068<br>S28 | RNAc-M076<br>S36 | RNAc-M084<br>S44 | RNAc-M092<br>S52 |
| E | DNAG-M005<br>S13 | DNAG-M009<br>S21 | DNAG-M021<br>S29 | DNAG-M029<br>S37 | DNAG-M037<br>S45 | DNAG-M045<br>S5 | RNAc-M053<br>S13 | RNAc-M061<br>S21 | RNAc-M069<br>S29 | RNAc-M077<br>S37 | RNAc-M085<br>S45 | RNAc-M093<br>S53 |
| F | DNAG-M006<br>S14 | DNAG-M010<br>S22 | DNAG-M022<br>S30 | DNAG-M030<br>S38 | DNAG-M038<br>S46 | DNAG-M046<br>S6 | RNAc-M054<br>S14 | RNAc-M062<br>S22 | RNAc-M070<br>S30 | RNAc-M078<br>S38 | RNAc-M086<br>S46 | RNAc-M094<br>S54 |
| G | DNAG-M007<br>S15 | DNAG-M011<br>S23 | DNAG-M023<br>S31 | DNAG-M031<br>S39 | DNAG-M039<br>S47 | DNAG-M047<br>S7 | RNAc-M055<br>S15 | RNAc-M063<br>S23 | RNAc-M071<br>S31 | RNAc-M079<br>S39 | RNAc-M087<br>S47 | RNAc-M095<br>S55 |
| H | DNAG-M008<br>S16 | DNAG-M012<br>S24 | DNAG-M024<br>S32 | DNAG-M032<br>S40 | DNAG-M040<br>S48 | DNAG-M048<br>S8 | RNAc-M056<br>S16 | RNAc-M064<br>S24 | RNAc-M072<br>S32 | RNAc-M080<br>S40 | RNAc-M088<br>S48 | RNAc-M096<br>S56 |



- Briefly centrifuge, mix by pipetting up and down 8 times, and then briefly centrifuge again.

**Note:** If only a column is used, cut that column from the cuttable plate and proceed to the next step.

- Program a thermal cycler using the cycling conditions based on the number of DNA or RNA primers (not combined DNA+RNA).

**Table 35. Cycling conditions for target enrichment if primers <1500**

| Step                 | Time   | Temperature |
|----------------------|--------|-------------|
| Initial denaturation | 13 min | 95°C        |
|                      | 2 min  | 98°C        |
| 8 cycles             | 15 sec | 98°C        |
|                      | 10 min | 68°C        |
| Hold                 | 5 min  | 72°C        |
|                      | ∞      | 4°C         |

**Table 36. Cycling conditions for target enrichment if number of primers ≥1500/tube**

| Step                 | Time (1500–12,000 primers/tube) | Time (>12,000 primers/tube) | Temperature |
|----------------------|---------------------------------|-----------------------------|-------------|
| Initial denaturation | 13 min                          | 13 min                      | 95°C        |
|                      | 2 min                           | 2 min                       | 98°C        |
| 6 cycles             | 15 s                            | 15 s                        | 98°C        |
|                      | 15 min                          | 30 min                      | 65°C        |
| 1 cycle              | 5 min                           | 5 min                       | 72°C        |
| Hold                 | 5 min                           | 5 min                       | 4°C         |
| Hold                 | ∞                               | ∞                           | 4°C         |

6. Place the target enrichment reaction in the thermal cycler and start the run.
7. Once the run has finished, add 80  $\mu$ l of Nuclease-free Water to bring each sample to 100  $\mu$ l.
8. Add 100  $\mu$ l QIAseq Beads and mix by vortexing or by pipetting up and down several times.
9. Incubate for 5 min at room temperature.
10. Place the tubes/plate on a magnetic rack for 2 min (for tubes) or 10 min (for plates).
11. After the solution has cleared, with the beads still on the magnetic stand, carefully remove and discard the supernatant.

**Important:** Do not discard the beads, because they contain the DNA of interest.

12. With the beads still on the magnetic stand, add 200  $\mu$ l of 80% ethanol. Rotate the tube 2–3 times or move the plate side-to-side between the 2 column positions of the magnet to wash the beads. Carefully remove and discard the wash.
13. Repeat the ethanol wash.

**Important:** Completely remove all traces of the ethanol after this second wash. Remove the ethanol with a 200  $\mu$ l pipette first, spin down briefly, and then use a 10  $\mu$ l pipette to remove any residual ethanol.

14. With the beads still on the magnetic stand, air-dry at room temperature for 5–10 min.

**Note:** Visually inspect the pellet to confirm that it is completely dry. If the pellet is not completely dry, additional drying time may be required. Ethanol carryover to the next universal PCR step will affect PCR efficiency.

15. Remove the beads from the magnetic stand, and then elute the DNA from the beads by adding 14  $\mu$ l Nuclease-free Water. Mix well by pipetting.
16. Return the tube/plate to the magnetic rack until the solution has cleared.
17. Transfer 12  $\mu$ l of the supernatant to clean tubes/plate.
18. Proceed to “Protocol: qPCR Determination of Universal PCR Cycles”, page 51.  
Alternatively, the samples can be stored at  $-30$  to  $-15^{\circ}\text{C}$  in a constant-temperature freezer.

# Ordering Information

| Product                              | Contents   | Cat. no. |
|--------------------------------------|--|----------|
| QIAseq Multimodal Panel (12)*        | Kit containing <b>all</b> reagents (except indexes) for multimodal (DNA and RNA) sequencing; fixed small panel for a total of 12 samples: 12 DNA and 12 RNA libraries                                | 333932   |
| QIAseq Multimodal Panel (96)*        | Kit containing <b>all</b> reagents (except indexes) for multimodal (DNA and RNA) sequencing; fixed small panel for a total of 96 samples: 96 DNA and 96 RNA libraries                                | 333935   |
| QIAseq Multimodal HC Panel (12)*     | Kit containing <b>all</b> reagents (except indexes) for multimodal (DNA and RNA) sequencing; fixed high content (HC) panel for a total of 12 samples: 12 DNA and 12 RNA libraries                    | 333942   |
| QIAseq Multimodal HC Panel (96)*     | Kit containing <b>all</b> reagents (except indexes) for multimodal (DNA and RNA) sequencing; fixed high content (HC) panel for a total of 96 samples: 96 DNA and 96 RNA libraries                    | 333945   |
| QIAseq Multimodal Custom Panel (96)* | Kit containing <b>all</b> reagents (except indexes) for multimodal (DNA and RNA) sequencing; custom panel for a total of 96 samples: 96 DNA and 96 RNA libraries                                     | 333955   |
| QIAseq Multimodal Index I (12)*      | Box containing oligos, enough to process a total of 12 samples, for indexing up to a total of 12 samples (12 for DNA and 12 for RNA libraries) for multimodal panel sequencing on Illumina platforms | 333962   |

\* Visit [www.qiagen.com/GeneGlobe](http://www.qiagen.com/GeneGlobe) to search for and order these products.

| <b>Product</b>                           | <b>Contents</b>   | <b>Cat. no.</b> |
|--|---|-----------------|
| QIAseq Multimodal Index I Set A (96)*    | Box containing oligos, enough to process a total of 96 samples, for indexing up to a total of 96 libraries (48 for DNA and 48 for RNA ) in one run for multimodal panel sequencing on Illumina platforms; Set A (of A and B) required for multiplexing 192 libraries (96 for DNA and 96 for RNA ) in one run  | 333965          |
| QIAseq Multimodal Index I Set B (96)*    | Box containing oligos, enough to process a total of 96 samples, for indexing up to a total of 96 libraries (48 for DNA and 48 for RNA ) in one run for multimodal panel sequencing on Illumina platforms; Set B (of A and B) required for multiplexing 192 libraries (96 for DNA and 96 for RNA ) in one run  | 333975          |
| QIAseq Multimodal Index I Set A SW (96)* | Box containing oligos, enough to process a total of 96 samples, for indexing up to a total of 96 libraries (48 for DNA and 48 for RNA ) in one run for multimodal panel sequencing on Illumina platforms using the “separated targeted enrichment” workflow; Set A (of A and B) required for multiplexing 192 libraries (96 for DNA and 96 for RNA ) in one run | 333985          |
| QIAseq Multimodal Index I Set B SW (96)* | Box containing oligos, enough to process a total of 96 samples, for indexing up to a total of 96 libraries (48 for DNA and 48 for RNA ) in one run for multimodal panel sequencing on Illumina platforms using the “separated targeted enrichment” workflow; Set B (of A and B) required for multiplexing 192 libraries (96 for DNA and 96 for RNA ) in one run | 333995          |

\* Visit [www.qiagen.com/GeneGlobe](http://www.qiagen.com/GeneGlobe) to search for and order these products.

| <b>Product</b>                           | <b>Contents</b>   | <b>Cat. no.</b> |
|--|---|-----------------|
| <b>Related products</b>                  |   |                 |
| QIAseq Library Quant Array Kit           | Reagents for quantification of libraries prepared for Illumina or Ion Torrent® platforms; array format  | 333304          |
| QIAseq Library Quant Assay Kit           | Reagents for quantification of libraries prepared for Illumina or Ion Torrent platforms; assay format   | 333314          |
| QIAseq DNA QuantitiMIZE Array Kit        | qPCR arrays for optimizing amount of input DNA and PCR cycling conditions for targeted enrichment of FFPE DNA   | 333404          |
| QIAseq DNA QuantitiMIZE Assay Kit        | qPCR assays for optimizing amount of input DNA and PCR cycling conditions for targeted enrichment of FFPE DNA   | 333414          |
| QIAamp® DNA Mini Kit (50)                | For 50 DNA preps: 50 QIAamp Mini Spin Columns, QIAGEN Proteinase K, collection tubes (2 ml), reagents and buffers   | 51304           |
| QIAamp Circulating Nucleic Acid Kit (50) | For 50 DNA preps: QIAamp Mini Columns, Tube Extenders (20 ml), QIAGEN Proteinase K, carrier RNA, buffers, VacConnectors, and collection tubes (1.5 ml and 2 ml) | 55114           |
| AllPrep DNA/RNA Mini Kit (50)            | For 50 minipreps: AllPrep DNA Spin Columns, RNeasy Mini Spin Columns, collection tubes, RNase-free water, and buffers   | 80204           |
| AllPrep DNA/RNA FFPE Kit (50)            | 50 RNeasy MinElute Spin Columns, 50 QIAamp MinElute Spin Columns, collection tubes, RNase-free reagents, and buffers  | 80234           |
| PAXgene Blood DNA Kit (25)               | Processing tubes and buffers for 25 preparations  | 761133          |

| Product                     | Contents   | Cat. no. |
|-----------------------------|--|----------|
| PAXgene Blood RNA Kit (50)  | 50 PAXgene Spin Columns, 50 PAXgene Shredder Spin Columns, processing tubes, RNase-Free DNase I, RNase-free reagents and buffers | Inquire  |
| GeneRead™ DNA FFPE Kit (50) | QIAamp MinElute® columns, proteinase K, UNG, collection tubes (2 ml), buffers, deparaffinization solution, RNase A               | 180134   |

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

# Document Revision History

| Date    | Changes  |
|---------|--|
| 10/2019 | Initial release  |
| 6/2020  | <p>Removed Index sequences from HB. Index sequences can be found on <a href="https://qiagen.com/QIAseqMultimodalPanels">qiagen.com/QIAseqMultimodalPanels</a>.</p> <p>Removed references to the Ingenuity Variant Analysis (IVA) tool.</p> <p>Updated volumes of clean up steps after DNA ligation</p> <p>Updated the following sections: Kit Contents (Optical Thin-wall 8-cap Strips), Storage, Principle and procedure, Nucleic acid input amount and Sequencing depth section, Protocol: Combined Targeted DNA+RNA Enrichment in a Single Tube Protocol: qPCR Determination of Universal PCR Cycles, Protocol: Universal PCR, and Protocol: Sequencing Setup on Illumina MiSeq and NextSeq.</p> <p>Layout and formatting changes</p> |

---

**This page is intentionally left blank.**

### Limited License Agreement for QIAseq Multimodal Panels

Use of this product signifies the agreement of any purchaser or user of the product to the following terms:

1. The product may be used solely in accordance with the protocols provided with the product and this handbook and for use with components contained in the kit only. QIAGEN grants no license under any of its intellectual property to use or incorporate the enclosed components of this kit with any components not included within this kit except as described in the protocols provided with the product, this handbook, and additional protocols available at [www.qiagen.com](http://www.qiagen.com). Some of these additional protocols have been provided by QIAGEN users for QIAGEN users. These protocols have not been thoroughly tested or optimized by QIAGEN. QIAGEN neither guarantees them nor warrants that they do not infringe the rights of third-parties.
2. Other than expressly stated licenses, QIAGEN makes no warranty that this kit and/or its use(s) do not infringe the rights of third-parties.
3. This kit and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. QIAGEN specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the kit agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. QIAGEN may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the kit and/or its components.

For updated license terms, see [www.qiagen.com](http://www.qiagen.com).

Trademarks: QIAGEN®, Sample to Insight®, QIAamp®, QIAseq®, QIAxpert®, AllPrep®, GeneRead™, HotStarTaq®, Ingenuity®, MinElute®, QCI™ (QIAGEN Group); Agilent®, Bioanalyzer®, TapeStation® (Agilent Technologies); Chrome®, Chromium® (Google LLC); DynaMag™, Ion Torrent®, NanoDrop®, Quanti-iT™ (Thermo Fisher Scientific Inc.); EvaGreen® (Biotium, Inc.); Illumina®, HiSeq®, MiniSeq®, MiSeq®, NextSeq®, NovaSeq™ (Illumina, Inc.); LoBind® (Eppendorf AG); MagneSphere® (Promega Corporation). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

06/2020 HB-2719-002 © 2020 QIAGEN, all rights reserved.

---

Ordering [www.qiagen.com/shop](http://www.qiagen.com/shop) | Technical Support [support.qiagen.com](http://support.qiagen.com) | Website [www.qiagen.com](http://www.qiagen.com)