

Protein expression analysis in Formalin-Fixed, Paraffin-Embedded (FFPE) tissue



Marcia Armstrong, Karl-Friedrich Becker*, Christina Schott*, Anja Erpenbach†, and Peter Porschewski†
 QIAGEN Inc, Germantown, MD, USA; * Technische Universität München, Munich, Germany;
 † QIAGEN GmbH, Hilden, Germany



Introduction

- As the standard method for histological tissue preparation, millions of Formalin-Fixed Paraffin-Embedded (FFPE) samples are classified and archived each year.
- FFPE archives constitute a huge library of clinically appraised specimens, whose analysis on the molecular level can be directly related to progression of a disease or therapy.
- Previously, due to the crosslinks formed during the fixation process, their analysis was limited to histochemical techniques.
- Here we present a method that opens this vast, highly valuable resource up to proteomic and protein expression analysis by enabling efficient extraction of full-length proteins suitable for western blotting and protein arrays.

Methods

- Unstained FFPE samples are deparaffinized and rehydrated by successive washes in xylene, 100%, 96%, and 70% ethanol and water. Proteins are extracted by incubation in an extraction buffer at 100°C for 20 min followed by a 2 h incubation at 80°C. After incubation, proteins are recovered by centrifugation.

Qproteome FFPE Tissue Procedure





Table 1. Protein Yields from Different Tissues as Starting Material


Tissue	Total size (mm ²)	No. of sections	Protein yield
Heart (rat)*	~150	3	40 µg
Liver (rat) *	~150	3	80 µg
Colon (rat) *	~100–150	3	20 µg
Brain (rat) *	~150	3	60 µg
Breast cancer (human)†	~100	1	25–80 µg

* Proteins were extracted from FFPE tissue sections (10 µm) directly cut from a FFPE tissue sample block.
 † Proteins were extracted from two areas with different morphological structures in the same FFPE tissue section (10 µm) mounted on a microscope slide.


Efficient and reproducible isolation of full-length proteins from FFPE tissues



Detection and identification of biomarkers



Proteins suitable for lysate arrays or western blot analysis



Summary

- Proteins obtained using the FFPE kit are full-length and suitable for array or western blot analysis
- Investigation of FFPE samples enables detection and identification of biomarkers on the protein level
- Extraction of proteins from FFPE tissue is efficient and reproducible, and yields are comparable with those obtained from frozen tissue
- Any kind of FFPE material can be processed, including unstained slide-mounted sections, freshly-cut sections, and tissue punches from FFPE blocks.