




**Clinical Sample Concentrators – QIAGEN Sample
Purification for In Vitro Diagnostic Use
11:00 AM, EST
December 14, 2010**

**Dr. Anke Homann-Wischinski, Ph.D.
Associate Director, QIAGEN R&D**





Clinical Sample Concentrator Kits IVD–class I

- Background on Clinical Sample Concentrators
 - Classification
 - Regulation
- QIAGEN solutions
 - Portfolio of Clinical Sample Concentrator Kits
 - Features and advantages
 - Performance characteristics
- Summary



Agenda



Clinical Sample Concentrator Kits IVD–class I

- **Background on Clinical Sample Concentrators**
 - **Classification**
 - **Regulation**

- **QIAGEN solutions**
 - **Portfolio of Clinical Sample Concentrator Kits**
 - **Features and advantages**
 - **Performance characteristics**

- **Summary**



Clinical Sample Concentrator Background I – Classification

The US Food and Drug Administration (FDA) requires that **sample preparation products** used for clinical diagnostic applications from human samples be classified as **Clinical sample concentrators** under the US Code of Federal Regulations pertaining to medical devices (21 CFR 862.2310).

FDA requires registration and listing of sample preparation products which are part of systems for human diagnostics, if such products are:

- Used by sponsor of a submission as front end
- Used by customers for clinical diagnostic use



Clinical sample concentrators (CSCs) are Class 1 Medical Devices that are exempt from FDA 510k pre-market notification submission requirements. CSCs require manufacture under an FDA compliant Quality System.



Clinical Sample Concentrator

Background II – regulation

TITLE 21 – FOOD AND DRUGS
CHAPTER I – FOOD AND DRUG ADMINISTRATION
DEPARTMENT OF HEALTH AND HUMAN SERVICES
SUBCHAPTER H – MEDICAL DEVICES
PART 862 – CLINICAL CHEMISTRY AND CLINICAL
TOXICOLOGY DEVICES
Subpart C – Clinical Laboratory Instruments

§ 862.2310 Clinical sample concentrator.

(a) Identification. A clinical sample concentrator is a device intended to concentrate (by dialysis, evaporation, etc.) serum, urine, cerebrospinal fluid, and other body fluids before the fluids are analyzed.

(b) Classification. Class I (general controls). The device is exempt from the premarket notification procedures in subpart E of part 807 of this chapter subject to the limitations in 862.9.





Clinical Sample Concentrator Kits IVD-class I

- Background on Clinical Sample Concentrators
 - Classification
 - Regulation

- **QIAGEN solutions**
 - **Portfolio of Clinical Sample Concentrator Kits**
 - **Features and advantages**
 - **Performance characteristics**

- Summary



QIAGEN's solutions

DSP Kits

Diagnostics Sample Preparation





QIAGEN's Clinical Sample Concentrator Offerings I

Available in November 2010

QIAamp[®] DSP Kits – manual

- | | |
|--|----------------|
| 1. QIAamp DSP Virus Kit (50) | Cat. no. 60704 |
| 2. QIAamp DSP Circulating NA Kit (50) | Cat. no. 61504 |
| 3. QIAamp DSP 96 DNA Blood Kit (12) | Cat. no. 61162 |
| 4. QIAamp DSP DNA FFPE Tissue Kit (50) | Cat. no. 60404 |



QIAamp DSP Kits – manual and automatable on QIAcube[®]

- | | |
|---------------------------------------|----------------|
| 5. QIAamp DSP DNA Blood Mini Kit (50) | Cat. no. 61104 |
| 6. QIAamp DSP DNA Mini Kit (50) | Cat. no. 61304 |
| 7. QIAamp DSP Virus Spin Kit (50) | Cat. no. 61704 |
| 8. QIAamp DSP Viral RNA Mini Kit (50) | Cat. no. 61904 |





QIAGEN's Clinical Sample Concentrator Offerings II

Available in November 2010

EZ1[®] DSP Kits for the EZ1 Advanced / Advanced XL

- | | |
|--------------------------------|----------------|
| 9. EZ1 DSP Virus Kit (48) | Cat. no. 62724 |
| 10. EZ1 DSP DNA Blood Kit (48) | Cat. no. 62124 |



QIASymphony[®] DSP Kits for the QIASymphony SP

- | | |
|--|-----------------|
| 11. QIASymphony DSP Virus/Pathogen Kit | |
| Mini Kit (192) | Cat. no. 937036 |
| Midi Kit (96) | Cat. no. 937055 |
| 12. QIASymphony DSP AXpH DNA Kit (192) | Cat. no. 937156 |





Features of Clinical Sample Concentrator Kits

Production standards

DSP Kits are ...

- Designed for in vitro diagnostic (IVD) use
- Manufactured in compliance with Quality System Regulations (QSR)
 - Produced under cGMP
- Expiry dated at kit and component level
- Tested by a functional kit end-test before release (QC)
- Registered and listed under 21 CFR 862.2310 (Clinical sample Concentrator)



Features of Clinical Sample Concentrator Kits

Regulatory compliance

DSP Kits are...

- Developed and manufactured in compliance with
 - ISO 13485 medical device: Quality management systems requirement for regulatory purposes
 - FDA Quality System Regulations 21 CFR 820

- In compliance with international IVD standards, e.g.
 - EN 13640 Stability testing of in vitro diagnostic reagents
 - ISO 14971 Application of risk management to medical devices



What does CSC mean to you?

Key benefits of DSP Kits

- Easy integration into existing diagnostic workflows
- GMP compliant for use in lab-developed assay and IVD applications
- Batch-to-batch consistency
- Robustness for FDA 510(K) clearance

As the *gold standard* in purification, QIAGEN is positioned to provide the best in class sample prep for in vitro diagnostic use

Performance evaluation

- ⇒ Limit of detection
- ⇒ Linear range
- ⇒ Accuracy
- ⇒ Precision
- ⇒ Robustness
- ⇒ Cross contamination
- ⇒ Stability testing
- ⇒ Eluate stability
- ⇒





Performance characteristics

Clinical Sample Concentrator

QIAamp DSP DNA Blood Kit system






Performance characteristics


QIAamp DSP DNA Blood Mini Kit (Cat. no. 61104) automated on QIAcube

Average absolute yields of DNA from blood samples collected using various collection tubes and anticoagulants





Mean absolute yield of DNA of 6 donor bloods per collection tube (QIAcube standard protocol with elution volume 200 μ l)

QIAamp DSP DNA Mini Kit (Cat. no. 61304) automated on QIAcube



DNA integrity –
genomic gel
(QIAcube manual lysis protocol)



Inhibition assay –
different PCR input volumes
(QIAcube standard protocol with 100 µl elution)




Performance characteristics

Clinical Sample Concentrator

QIAamp DSP Circulating NA Kit system



QIAamp DSP Circulating NA Kit (Cat. no. 61504)



Recovery of spike-in control in EDTA plasma

- 200 bp and 1000 bp fragments added to plasma before extraction
- Recovery measured by duplex real-time PCR
- Reproducible DNA yield among different donors



Performance characteristics

Clinical Sample Concentrator


EZ1 DSP DNA Blood Kit system





Performance characteristics

EZ1 Adv XL, Blood 200 μ / 200 μ														
1	MALE	W												
2	W	MALE												
3	W	W	W	W	W	W	W	W	W	W	W	W	W	W
4	MALE	W												
5	W	MALE												
6	W	W	W	W	W	W	W	W	W	W	W	W	W	W
7	MALE	W												
8	W	MALE												
9	W	MALE												



Typical amplification plot of male samples
 Mean quantification: 35,700 copies/PCR

Positive controls (10 copies/ μ l)

Negative samples
 – no SYBR sequence detected

No detectable cross-contamination



Performance characteristics

Clinical Sample Concentrator

QIASymphony DSP Virus/Pathogen Kit system






Performance characteristics

The QIA Symphony DSP Virus/Pathogen Mini and Midi Kits can be used for fully automated purification of nucleic acids from a broad range of DNA and RNA viruses as well as bacterial DNA from Gram-negative and Gram-positive bacteria from


- Plasma
- Serum
- CSF
- Respiratory samples (BAL, dried swabs, transport media, aspirates, sputum)
- Urogenital samples (urine, transport media)





Performance characteristics

QIASymphony DSP Virus/Pathogen Kits



Panel member	n	IU/ml	CV (%)	log IU/ml	SD (log IU/ml)
1	6	1835700	30.04	6.24	0.15
2	6	199931	26.99	5.28	0.13
3	5	13785	21.02	4.13	0.09
4	5	1363	17.49	3.13	0.09
5	6	642	24.82	2.79	0.12
6	6	294	31.12	2.44	0.16
7	6	123	23.25	2.08	0.11

Linearity

Linear range of yields using the Virus Cellfree 1000 protocol. The linear range of the Virus Cellfree protocol was determined. Standard deviations and coefficients of variations (CVs) were determined for different dilution series in the linear range of the appropriate downstream assays.



Clinical Sample Concentrator Kits IVD–class I

- Background on Clinical Sample Concentrators
 - Classification
 - Regulation
- QIAGEN solutions
 - Portfolio of Clinical Sample Concentrator Kits
 - Features and advantages
 - Performance characteristics
- Summary



Summary

Clinical Sample Concentrator Kits IVD–class I

QIAGEN now offers you a broad range of Diagnostic Sample Preparation Kits listed under CFR 862.2310 for manual and automated isolation of DNA and RNA from:

- Blood
- Plasma/serum
- CSF
- FFPE tissue
- Respiratory samples
- Urogenital samples

...to simplify your diagnostic workflow and help you remain in compliance with your regulatory requirements

Questions and Answers

