

DNeasy® PowerSoil® Pro Kits

Overview

The next-generation kits for the isolation of microbial DNA from all soil types

Procedure

Lyse Tough
Microbes

Remove
Inhibitors

Improve NGS
Results

Ordering
Information



Procedure

Overview

Procedure


Lyse Tough Microbes

Remove Inhibitors

Improve NGS Results

Ordering Information

DNeasy PowerSoil Pro Kit



Starting material:
up to 250 mg of soil

Elution Volume: 50–100 μ l

Automated on the QIAcube[®]

Shorter workflow compared to original DNeasy PowerSoil Kit:

- 4 fewer steps
- No incubation times
- Total time saved = 12 min



Download Quick-Start Protocol

Overview

Procedure


Lyse Tough Microbes

Remove Inhibitors

Improve NGS Results

Ordering Information


New bead tubes increase lysis efficiency for higher DNA yields



Various soil sample types (250 mg) were prepared using commercially available sample preparation solutions and compared to the new DNeasy PowerSoil Pro Kit. Yields were measured by fluorometric quantification (Qubit®).

Up to 8-fold higher DNA yields than alternative methods

Streamlined Inhibitor Removal Technology[®] (IRT) removes inhibitors for pure DNA



Various soil samples (250 mg) were prepared using commercially available sample preparation solutions and compared to the new DNeasy PowerSoil Pro Kit. DNA purity was measured via spectrophotometry and is presented as the ratios **A** A_{260}/A_{280} and **B** A_{260}/A_{230} .

Faster, one-step IRT results in highly pure DNA

Overview

Procedure

Lyse Tough Microbes

Remove Inhibitors

Improve NGS Results

Ordering Information

Better lysis of tough microbes results in increased alpha diversity in 16S rRNA gene sequencing

Overview


Procedure

Lyse Tough Microbes

Remove Inhibitors

Improve NGS Results

Ordering Information



DNA prepared from soil samples was isolated with different methods and enriched. Analysis of the 16S rRNA genes was done using the QIAseq[®] 1-Step Amplicon Kit and data analysis was done using the Microbial Genomics Pro suite (CLC workbench). Alpha diversity was determined as total number of bacterial operational taxonomic units (OTUs).

More bacteria (OTUs) identified than alternative methods

- Overview
- Procedure
- Lyse Tough Microbes
- Remove Inhibitors
- Improve NGS Results

Ordering Information

Product	Description	Cat. no.
DNeasy PowerSoil Pro (50)	For 50 DNA minipreps: Buffers, PowerBead Tubes, Spin Filters, 2 ml Collection Tubes	47014
DNeasy PowerSoil Pro Kit (250)	For 250 DNA minipreps: Buffers, PowerBead Tubes, Spin Filters, 2 ml Collection Tubes	47016

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, Sample to Insight®, QIAcube®, QIAseq®, DNeasy®, Inhibitor Removal Technology®, PowerSoil®, PowerLyzer® (QIAGEN Group); FastDNA® (MP Biomedicals); E.Z.N.A.® (Omega Bio-tek, Inc.); PureLink®, Qubit® (Thermo Fisher Scientific or its subsidiaries); ZymoBIOMICS® (Zymo Research Corp.).

© 2018 QIAGEN, all rights reserved.

Ordering www.qiagen.com/contact | Technical Support support.qiagen.com | Website www.qiagen.com

Visit Webpage

Ordering Information