

June 2020

Labware List — QIAAsymphony® DSP HPV Media Kit

Sample and eluate tubes/racks that can be used with the QIAAsymphony DSP HPV Media Kit and the QIAAsymphony SP (software version 5.0; labware package SOW-516-8)

This document is the QIAAsymphony DSP HPV Media Kit, version 1, *Labware List* R1 for labware package SOW-516-8.

General information

The QIASymphony DSP HPV Media Kit is intended for in vitro diagnostic use.

Important: Before using this Labware List, confirm that it is consistent with the labware package installed on your QIASymphony system.

Note: The QIASymphony DSP HPV Media Kit is not available in the United States. Please visit www.qiagen.com for product availability.

Legend

- Recommended: These tubes are for use in combination with this protocol.
- User-defined: User is responsible for defining and/or validating as the user deems appropriate.
- Not recommended: These tubes are not for use in combination with this protocol.

- x/■ Minimum sample volume (µl) required per sample per protocol (denoted by x)/clot detection possible.
- x/□ Minimum sample volume (µl) required per sample per protocol (denoted by x)/clot detection not possible.

Note: Be aware that other tubes are available that can be used in combination with other QIASymphony Kits, but are not recommended for any protocols that can be used with this kit.

“Sample” drawer, tube carrier

Supplier	Material	Cat. no.	Name in touchscreen	Insert	Protocol	
					PC2500_HC2	SP800_HC2
BD*	14 ml Falcon® polystyrene round-bottom tube 17 x 100 mm	352051	BD#352051 FalconPP 17x100	No insert required	x/□	
QIAGEN	2 ml Micro tube, PP, NON-SKIRTED	997055	QIA#997055 T2.0 Conical blue	03/3B		x/□

* BD was the previous supplier of this tube and Corning Inc. is now the new supplier.

“Eluate” drawer

Supplier	Material	Example cat. no.	Category	Name in touchscreen	Adapter on Elution slots	Protocol	
						PC2500_HC2	SP800_HC2
QIAGEN	Hybridization microplate	6000-1203	Micro Plate	QIA#6000-1203 *MTP96 RB	96-Well Round Bottom QS (cooling slot 1) Elution slots 2–4 are not used in this procedure		

* Indicates labware that can be cooled using a cooling adapter with or without bar code.

For ordering information, visit www.qiagen.com/goto/QIAsymphony.

Revision history

Document Revision History	
R1 06/2020	Release of new labware package

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN® kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, Sample to Insight®, QIAsymphony® (QIAGEN Group); BD™ (Becton, Dickinson and Company); Corning®, Falcon® (Corning, Inc.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.
06/2020 HB-1554-105-001 © 2020 QIAGEN, all rights reserved.

Ordering www.qiagen.com/shop | Technical Support support.qiagen.com | Website www.qiagen.com