

January 2011

Signal Finder Reporter Array Plate Format Handbook

For cell-based multi-pathway activity assays



Sample & Assay Technologies

QIAGEN Sample and Assay Technologies

QIAGEN is the leading provider of innovative sample and assay technologies, enabling the isolation and detection of contents of any biological sample. Our advanced, high-quality products and services ensure success from sample to result.

QIAGEN sets standards in:

- Purification of DNA, RNA, and proteins
- Nucleic acid and protein assays
- microRNA research and RNAi
- Automation of sample and assay technologies

Our mission is to enable you to achieve outstanding success and breakthroughs. For more information, visit www.qiagen.com.

Product Use Limitations

Signal Finder 10-Pathway Reporter Arrays are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

CONTENTS

I.	Introduction	3
II.	Product Contents and Descriptions	5
III.	Additional Materials Required	8
IV.	Protocol	9
	A. Before you begin	9
	B. Protocol	11
	Appendix: Signal Finder Multi-Pathway Array Product Descriptions	14
	Ordering Information	16

I. Introduction

The Signal Finder Multi-Pathway Reporter Arrays enable you to pinpoint the pathways regulated by the gene products or chemical compounds being studied in your laboratory. The Signal Finder Arrays consist of 10 or 45 dual-luciferase reporter assays, and are designed for use in one of four research areas. The targeted research areas are cancer, immunology, development, and toxicology. In this era of post-genomics life science research, many labs are investigating how diverse signal transduction pathways function on their own, and in combination, within the cell. The Signal Finder Arrays equip life science researchers to carry out such studies with speed and confidence.

These arrays are cell culture-ready 96-well plates. For the 10-pathway arrays, each of the twelve columns of the 96-well plate contains a pathway-focused reporter or control dried down in all eight wells. For the 45-pathway array, each pathway reporter assay is dried down in two wells, with the remaining wells being used for positive and negative controls. The reporter assays are reverse transfected into your cells.

Each pathway-focused dual-luciferase reporter encodes for the mammalian codon-optimized, non-secreted form of the firefly luciferase gene, carrying a protein-destabilizing sequence. Cells rapidly degrade the destabilized form of the firefly luciferase protein and hence the background luciferase activity (noise level) is greatly reduced. Due to low background activity, the magnitude of the response that can be measured (signal-to-noise ratio) as well as the speed of measuring changes in transcription are enhanced. The Signal dual-luciferase reporter assays provide outstanding reproducibility, sensitivity, specificity, and signal-to-noise ratio. The Signal reporters are useful assays for carrying out quantitative pathway regulation studies.

Benefits of Signal Finder Multi-Pathway Reporter Arrays

- **Multi-Pathway Analysis:** Profile the changes in the activities of ten or forty-five signaling pathways relevant to a specific biological process
- **High Performance:** Dual-luciferase assay provides high sensitivity, specificity, and reproducibility
- **Flexibility and Convenience:** Utilize a straightforward reverse transfection procedure with your favorite cell lines to rapidly generate valuable mechanism of action data




Figure 1: Overview of Signal Finder 10-Pathway Reporter Array Protocol.

II. Product Contents and Descriptions

A. Product Contents

Signal Finder 10-Pathway Reporter Array Contents:

Table 1: Signal Finder Reporter Array (plate format) Specifications

Component	Specification	Total DNA in each well
Each of the 10 Reporter Assays	A mixture of an inducible transcription factor responsive firefly luciferase reporter and constitutively expressing <i>Renilla</i> construct (20:1).	200 ng
Negative control	A mixture of non-inducible firefly luciferase reporter and constitutively expressing <i>Renilla</i> construct (20:1).	200 ng
Positive control	A mixture of a constitutively expressing GFP construct, constitutively expressing firefly luciferase construct, and constitutively expressing <i>Renilla</i> luciferase construct (20:1:1).	200 ng

Signal Finder 45-Pathway Reporter Array Contents:

Table 2: Signal Finder Reporter Array (plate format) Specifications

Component	Specification	Total DNA in each well
Each of the 45 Reporter Assays	A mixture of an inducible transcription factor responsive firefly luciferase reporter and constitutively expressing <i>Renilla</i> construct (20:1).	200 ng
Negative control	A mixture of non-inducible firefly luciferase reporter and constitutively expressing <i>Renilla</i> construct (20:1).	200 ng
Positive control	A mixture of a constitutively expressing GFP construct, constitutively expressing firefly luciferase construct, and constitutively expressing <i>Renilla</i> luciferase construct (20:1:1).	200 ng

NOTE: All constructs are **transfection-grade** and are ready for transient transfection. These constructs are specifically designed to inhibit transformation and are **NOT MEANT** for introduction and amplification in bacteria.

Each kit also includes a **white self-adhesive sealing tape** for each plate included in the kit. This tape should be affixed to the bottom of each plate immediately prior to reading the plate in a plate-reading luminometer, in order to maximize the signal-to-noise ratio of each reading.

B. Description of Individual Signal Reporter Assays:

Each Signal Reporter Assay Kit includes the following components:

1. **Reporter:** Each reporter is a mixture of an inducible transcription factor responsive construct and constitutively expressing *Renilla* luciferase construct (20:1). The inducible transcription factor-responsive construct encodes the firefly luciferase reporter gene under the control of a basal promoter element (TATA box) joined to tandem repeats of a specific Transcriptional Response Element (TRE; Figure 2A). This construct monitors both increases and decreases in the activity of a key transcription factor, which is a downstream target of a specific signaling pathway. The constitutively expressing *Renilla* construct encodes the *Renilla* luciferase reporter gene under the control of a CMV immediate early enhancer/promoter (Figure 2B) and acts as an internal control for normalizing transfection efficiencies and monitoring cell viability. It is also useful to confirm transfection and to verify active luciferase in the transfected culture.
2. **Negative control:** The negative control is a mixture of non-inducible reporter construct and constitutively expressing *Renilla* luciferase construct (20:1). The non-inducible reporter construct encodes firefly luciferase under the control of a basal promoter element (TATA box), without any additional transcriptional response elements (Figure 2C). The negative control is critical to identifying specific effects and determining background reporter activity.
3. **Positive control:** The positive control is a constitutively expressing GFP construct (Figure 2D), pre-mixed with a constitutively expressing firefly luciferase construct (Figure 2E), and a constitutively expressing *Renilla* luciferase construct (Figure 2B) (20:1:1). The positive control is necessary for visual confirmation of transfection. It is also useful for transfection optimization studies. The expression of the GFP from the positive control construct can be monitored by fluorescence microscopy using an excitation filter of 470 ± 20 nm (470 / 40 nm) and an emission filter of 515 nm (long pass).




Figure 2: Schematic representation of constructs involved in the Signal Reporter Assay. (A) The inducible transcription factor-responsive construct expressing firefly luciferase, (B) The constitutively expressing *Renilla* luciferase construct, (C) The non-inducible firefly luciferase reporter construct, (D) The constitutively expressing GFP construct, and (E) The constitutively expressing firefly luciferase construct.

IMPORTANT NOTE: There are a few reports in the literature of the CMV regulatory element being activated by certain stimuli (see below). We recommend that you confirm that the stimulus used in each Signal reporter assay does not induce the CMV regulatory element, in order to confirm that the CMV-*Renilla* construct is the appropriate normalization construct for the experiment. This can be done empirically by testing the impact of a stimulus on the Signal positive control reporters, which are each under the control of the CMV enhancer/promoter cassette. If stimulus is one of the very few reported activators of the CMV regulatory element, we advise contacting technical support.

- W. Bruening, B. Giasson, W. Mushynski, and H. D. Durham. 1998. *Nucleic Acids Research* 26(2):486-489. **Activation of stress-activated MAP protein kinases up-regulates expression of transgenes driven by the cytomegalovirus immediate/early promoter**

- Madhu S. Malo, Moushumi Mozumder, Alexander Chen, Golam Mostafa, Xiao Bo Zhang, Richard A. Hodin. 2006. Analytical Biochemistry 350:307-309. **pFRL7: An ideal vector for eukaryotic promoter analysis**

III. Additional Materials Required:

- Mammalian cell line cultured in the appropriate growth medium
- Cell culture medium and standard cell culture supplies
- Multi-channel pipettor and pipettor reservoirs
- Transfection reagent [Recommended reagent: Attractene Transfection Reagent (cat. no. 301005), however, other transfection reagents work equally well]
- Polystyrene test tubes (BD FALCON, Cat # 352099)
- Opti-MEM® I Reduced Serum Medium (Invitrogen, Cat. No. 31985-062)
- Fetal bovine serum (FBS)
- Non-essential amino acids (NEAA) (Invitrogen, Cat. No. 11140-050)
- Penicillin/Streptomycin
- Hemacytometer
- Dual-Luciferase® Assay System
 - Dual-Luciferase® Reporter Assay System (Promega, Cat. No. E1910)
This system requires cell lysis, and is well-suited for the rapid quantitation of both luciferase reporters when using luminometers with reagent auto-injectors.
 - Dual-Glo® Luciferase Assay System (Promega, Cat. No. E2920)
This system is used to assay for both luciferase reporters on intact cells in growth medium. This system can be used with any luminometer, including those without reagent auto-injectors.
- 96-well white opaque flat bottom microtiter plate
- Luminometer


IV. Protocol:

A. Before you begin:

1. Cell line selection: The Signal Reporter Assay may be used with various mammalian cell lines. Cell lines show a great deal of variation in the levels of signaling proteins. The transcriptional activator activities in the cell line used will determine the sensitivity of the assay. A cell line should be selected based on the functionality of the signal transduction pathway under investigation, as well as for the “transfectability” of the cell line (see below).
2. Transfection reagent selection: We recommend the use of Attractene Transfection Reagent (cat. no. 301005) as a transfection reagent. The Signal Reporter Assay, however, also performs equally well with other transfection reagents. When using alternative transfection reagents, please refer to the manufacturer’s instructions on the use of those reagents.
3. Optimization of transfection conditions: The sensitivity of the Signal Reporter Assay depends on the transfection efficiency. The transfection efficiency, in turn, primarily depends upon cell line used. Therefore, it is very important to optimize the transfection conditions for each cell type under study. Variables to consider, when optimizing the transfection conditions include cell density, cell viability, amount of DNA, ratio of DNA to transfection reagent, transfection complex formation time, and transfection incubation time (see the detailed protocols for our recommendations). The positive control construct included with each Signal Reporter Assay can be used for determining the optimal transfection conditions.
4. Optimization of assay condition: The response rate in the Signal Reporter Assay depends on the assay conditions (conditions of the experimental treatment). To obtain maximum response given by any stimulus, perform dosing and time-course studies. The optimal amount of stimulus and the time of treatment must be obtained empirically for each experiment (see different protocols for our recommendations).
5. Important recommendations for best results:
 - A. Perform all transfections in **triplicate** to minimize variability among treatment groups.
 - B. Include positive and negative controls in each experiment to obtain reliable results.
 - C. Use low-passage cells that are actively growing and are greater than 90% viable, for maximal transfection efficiencies.
 - D. Do not add antibiotics to media during transfection, as this may cause cell death.
 - E. Take care to always seed the same number of cells in each well, in order to maximize the reproducibility of your experiment.

- F. Serum induces various signaling pathways, leading to cross-talk and high background. Therefore, use reduced amounts of serum (0.5%) in the assay medium during the experimental treatment to minimize these serum effects.
6. **Transfection Protocols:** In order to use the Signal Finder 10 or 45-Pathway Arrays in the plate format, a reverse transfection method must be employed. This approach involves seeding the cell line of interest onto the transfection complexes in a one day procedure. This is in contrast to traditional transfection methods, in which cells are seeded on the first day of the experiment and transfection complexes are added to the cells the following day. Attractene Transfection Reagent has been specifically developed as a reverse transfection reagent. Optimized reverse transfection protocols using the Attractene Transfection Reagent are described throughout the Signal Finder Reporter Arrays User Manual. Utilizing reverse transfection procedures results in both a time savings as well as improved reproducibility, when compared to traditional forward transfection methods. Conditions for using transfection reagents from other vendors in reverse transfection protocols may also be developed. This will require initial process optimization studies. Below is a general protocol overview for reverse transfection of the Signal Finder 10- or 45-Pathway Reporter Arrays.

Reverse Transfection Protocol Overview (1 DAY PROCEDURE)



- Add 50 μL of Opti-MEM to each well of Signal Finder array plate to resuspend reporter constructs
- Dilute Attractene into Opti-MEM
- Add 50 μL of diluted Attractene to 50 μL of resuspended reporter constructs, mix well and incubate at room temperature (15–25°C) for 20 minutes

- Trypsinize (if necessary), count, and suspend cells to appropriate density
- Immediately seed 50 μ L of suspended cells to each well
- Replace growth medium after 16-24 hours of transfection

B. PROTOCOL

The following protocol is designed to reverse transfect an adherent cell line, HEK293, using Attractene Transfection Reagent (cat. no. 301005). *If you are using a transfection reagent other than Attractene follow their manufacturer's protocol for optimizing transfection.* **This is just a general guideline; the optimal conditions/amounts should be optimized according to the cell type and the study requirements.** Read the protocol completely before starting the experiment.

IMPORTANT: (1) Do not add antibiotics to media during transfection as this causes cell death.

1. Add 50 μ l of Opti-MEM® to each well of the Cignal Finder Array plate. Resuspend the reporter assay constructs by gently tapping the side of the plate, while slightly rocking the plate back and forth, then left to right, five times each and incubate it for 5 minutes at room temperature.
2. QIAGEN uses 0.6 μ l of Attractene in 50 μ l of Opti-MEM® per well for each individual transfection. In order to prepare sufficient Attractene for an entire 96-well plate, QIAGEN recommends diluting 64.8 μ l of Attractene into 5400 μ l of Opti-MEM® (sufficient for 108 transfections). Mix gently by inverting tube slowly and set the tube at room temperature for 5 minutes.
3. After the 5 minute incubation, add 50 μ l of diluted Attractene into each well containing 50 μ l of the diluted nucleic acids (1:1 ratio).
4. Mix by gently tapping the sides of the plate for at least 30 seconds and incubate for 20 minutes at room temperature to allow complex formation to occur.
5. Meanwhile, wash cells in a culture dish once with Dulbecco's PBS **without calcium and magnesium**, and treat with 1-3 ml trypsin-EDTA for 2-5 minutes at 37°C in a humidified atmosphere containing 5% CO₂. Suspend the cells in 7-9 ml of Opti-MEM® containing 5% of fetal bovine serum, then centrifuge the cells down, remove the supernatant, and resuspend the cells to 8×10^5 cells/ml in Opti-MEM® containing 10% of fetal bovine serum and 1% NEAA**. To ensure reproducible transfection results, it is important to accurately measure the cell density with a hemacytometer or an automated cytometry device.

6. After the 20 minute incubation for complex formation is completed, mix the cell suspension by several inversions of the tube containing the cells or by gentle pipeting of the cell suspension.

7. Add 50 μ l of prepared cell suspension (8×10^4 cells in Opti-MEM® containing 10% of fetal bovine serum) to each well containing constructs- Attractene complexes. This gives a final volume in each well of 150 μ l. Mix gently by rocking the plate back and forth, then left to right. Do not move the plate in a circular motion, as this may cause the cells to preferentially sediment around the edges of each well.

8. Incubate cells at 37°C in a 5% CO₂ incubator for 16-24 hours.

9. After 16-24 hours of transfection, change the medium to complete growth medium (DMEM with 10% FBS, 0.1mM NEAA, 1mM Sodium pyruvate, 100 U/ml penicillin and 100 μ g/ml streptomycin).

10. Carry out the **luciferase assay** using either the Dual-Luciferase Reporter Assay System or Dual-Glo Luciferase Assay System from Promega. Follow the manufacturer's protocol for developing the assay. Please see specific recommendations in the **Important Notes** section below, for some general recommendations on when to carry out the luciferase assays for different types of studies. Each Signal Finder Array plate comes along with a white self-adhesive sticker, which should be attached to the bottom of the plate before reading the luciferase activity. Using the sticker to cover the optical bottom of the 96-well plate helps to maximize the signal-to-noise ratio of each reading.

Important Notes: Listed below are general recommendations for different experimental designs.

1. To determine the **effect of siRNA/shRNA** on different cell signaling pathways, we recommend doing transient co-transfection of siRNA/shRNA and reporter constructs. For this one can add 2 pmol of siRNA or 200 ng of shRNA plasmid to the resuspended reporter construct in step 1 of the protocol. The luciferase assay can be developed 48-72 hours after the co-transfection. Please remember to include negative control siRNA/shRNA to assist in the interpretation of your results.
2. To determine the **effect of cDNA overexpression** on different cell signaling pathways, we recommend doing the transient co-transfection of experimental vector and reporter constructs. For this one can add 100-200 ng of experimental vector to the resuspended reporter construct in step 1 of the protocol. The luciferase assay can be developed 36-48 hours after the co-transfection. Please remember to include negative control vector(empty vector) to assist in the interpretation of your results.
3. To determine the **effect of recombinant protein or small peptide** on different cell signaling pathways, we recommend changing the cell medium to assay medium (Opti-MEM® containing 0.5% of fetal bovine serum, 1% NEAA, 100 U/ml Penicillin and 100 μ g/ml Streptomycin) instead of growth medium in step 9 and treating the

transfected cells with 3 or 4 different concentrations of recombinant protein or small peptide 6 to 24 hours prior to assay development.

4. To determine the **effect of small chemicals** on different cell signaling pathways, we recommend changing the cell medium to assay medium (Opti-MEM® containing 0.5% of fetal bovine serum, 1% NEAA, 100 U/ml Penicillin and 100 µg/ml Streptomycin) instead of growth medium in step 9 and treating the transfected cells with 3 or 4 different concentrations of small chemicals 6 to 24 hours prior to assay development.

*Cells that have been passed 1:3 or 1:4 the day before are generally more easily transfected than cells that have reached a confluent state at the time of use.

**In most cases, cells grow well in Opti-MEM® serum-reduced growth medium with 3-5% FBS due to extra growth factors and nutrients supplied in Opti-MEM®. Cell should reach ~50-90% confluence once attached to the wells, otherwise increase the cell numbers.

For any other troubleshooting or technical questions about the Signal Reporter Assay, please call one of our Technical Support representatives at 1-888-503-3187 or 301-682-9200 or email at support@SABiosciences.com.

Appendix:

Signal Finder 10- and 45- Pathway Reporter Arrays

Signal Finder Signal Transduction 45-Pathway Reporter Array

(Plate Format: CCA-901L)

See www.sabiosciences.com for Signal 45-Pathway Reporter Arrays

Signal Finder Cancer 10-Pathway Reporter Array

(Tube Format: CCA-001L; Plate Format: CCA-101L)

Pathway	Transcriptional Regulatory Element (TRE)	Transcription Factor
Wnt	TCF/LEF response element	TCF/LEF
Notch	RBP-J κ binding element	RBP-J κ
p53/DNA Damage	p53 response element	p53
TGF β	SMAD response element	SMAD2/SMAD3/SMAD4
Cell Cycle/pRb-E2F	E2F binding element	E2F/DP1
NF κ B	NF κ B binding element	NF κ B
Myc/Max	E-box binding element	Myc/Max
Hypoxia	HIF response element	Hypoxia-inducible factor-1 (HIF-1)
MAPK/ERK	Serum response element (SRE)	Elk-1/SRF
MAPK/JNK	AP-1 binding element	AP-1

Signal Finder Immune Signaling 10-Pathway Reporter Array

(Tube Format: CCA-008L; Plate Format: CCA-108L)

Pathway	Transcriptional Regulatory Element (TRE)	Transcription Factor
NF κ B	NF κ B binding element	NF κ B
Type I Interferon	Interferon stimulated response element (ISRE)	STAT1/STAT2
Interferon Gamma	Interferon gamma activation sequence (GAS)	STAT1/STAT1
IL-6	STAT3 binding element	STAT3
Interferon Regulation	IRF-1 binding element	IRF-1
TGF β	SMAD response element	SMAD2/SMAD3/SMAD4
cAMP/PKA	cAMP regulatory element (CRE)	CREB
PKC/Ca ⁺⁺	NFAT response element	NFAT
C/EBP	C/EBP binding element	C/EBP
Glucocorticoid Receptor	Glucocorticoid response element (GRE)	Glucocorticoid Receptor (GR)

Signal Finder Development 10-Pathway Reporter Array

(Tube Format: CCA-003L; Plate Format: CCA-103L)

Pathway	Transcriptional Regulatory Element (TRE)	Transcription Factor
Notch	RBP-J κ binding element	RBP-J κ
Wnt	TCF/LEF response element	TCF/LEF
Myc/Max	E-box binding element	Myc/Max
NF κ B	NF κ B binding element	NF κ B
TGF β	SMAD response element	SMAD2/SMAD3/SMAD4
Cell Cycle/pRb-E2F	E2F binding element	E2F/DP1

C/EBP	C/EBP binding element	C/EBP
cAMP/PKA	cAMP regulatory element (CRE)	CREB
MAPK/ERK	Serum response element (SRE)	Elk-1/SRF
MAPK/JNK	AP-1 binding element	AP-1

Signal Finder Stem Cell & Differentiation 10-Pathway Reporter Array

(Tube Format: CCA-006L; Plate Format: CCA-106L)

Pathway	Transcriptional Regulatory Element (TRE)	Transcription Factor
Oct4	Oct4 binding element	Oct4
Nanog	Nanog binding element	Nanog
KLF4	KLF4 binding element	KLF4
Sox2	Sox2 binding element	Sox2
Myc/Max	E-box binding element	Myc/Max
Hedgehog	Gli binding element	Gli
Notch	RBP-J κ binding element	RBP-J κ
Wnt	TCF/LEF response element	TCF/LEF
Pax6	Pax6 binding element	Pax6
MEF2	MEF2 binding element	MEF2

Signal Finder Nuclear Receptors 10-Pathway Reporter Array

(Tube Format: CCA-005L; Plate Format: CCA-105L)

Pathway	Transcriptional Regulatory Element (TRE)	Transcription Factor
Estrogen	Estrogen Response Element	Estrogen Receptor (ER)
Androgen	Androgen Response Element	Androgen Receptor (AR)
PPAR	PPARbinding element	PPAR
Retinoic Acid	Retinoic Acid Response Element	RAR
Vitamin D	Vitamin D Response Element	Vitamin D Receptor (VDR)
Glucocorticoid	Glucocorticoid Response Element	Glucocorticoid Receptor (GR)
Progesterone	Progesterone Response Element	Progesterone Receptor (PR)
Retinoid X	RXR binding element	RXR
Liver X	LXR binding element	LXR α
Hepatocyte Nuclear Factor 4	HNF4 binding element	HNF4

Signal Finder Stress & Toxicity 10-Pathway Reporter Array

(Tube Format: CCA-007L; Plate Format: CCA-107L)

Pathway	Transcriptional Regulatory Element (TRE)	Transcription Factor
Antioxidant Response	Antioxidant Response Element (ARE)	Nrf2/Nrf1
DNA Damage	p53 response element	p53
NF κ B	NF κ B binding element	NF κ B
Hypoxia	HIF response element	HIF-1 α
ER Stress	ER Stress Response Element (ERSE)	CBF/NF-Y/YY1
Heavy Metal Stress	MTF1 binding element	MTF1
Heat Shock	Heat Shock Response Element (HSE)	HSF-1
Glucocorticoid	Glucocorticoid response element (GRE)	Glucocorticoid Receptor (GR)
MAPK/JNK	AP-1 binding element	AP-1
Xenobiotic	Xenobiotic Response Element	AhR

Ordering Information

Product	Contents	Cat. no.
Signal Finder 10-or 45-Pathway Reporter Arrays (plate format)	10 or 45 dual-luciferase reporter assays in plate format	Varies

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN® (QIAGEN Group).

Limited License Agreement

Use of this product signifies the agreement of any purchaser or user of the Cignal Finder 10-Pathway Reporter Arrays (plate format) to the following terms:

1. The Cignal Finder 10-Pathway Reporter Arrays (plate format) may be used solely in accordance with the *Cignal Finder Reporter Array Plate Format Handbook* and for use with components contained in the Kit only. QIAGEN grants no license under any of its intellectual property to use or incorporate the enclosed components of this Kit with any components not included within this Kit except as described in the *Cignal Finder Reporter Array Plate Format Handbook* and additional protocols available at www.qiagen.com.
2. Other than expressly stated licenses, QIAGEN makes no warranty that this Kit and/or its use(s) do not infringe the rights of third-parties.
3. This Kit and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. QIAGEN specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the Kit agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. QIAGEN may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the Kit and/or its components.

For updated license terms, see www.qiagen.com.

Firefly and/or Renilla Luciferase and Monster Green Limited Use Label License

READ THIS FIRST BEFORE OPENING PRODUCT

For research use only. The terms of the limited license conveyed with the purchase of this product are as follows: Researchers may use this product in their own research and they may transfer derivatives to others for such research use provided that at the time of transfer a copy of this label license is given to the recipients and the recipients agree to be bound by the conditions of this label license. Researchers shall have no right to modify or otherwise create variations of the nucleotide sequence of the luciferase gene or Monster Green® gene except that Researchers may: (1) clone heterologous DNA sequences at either or both ends of said luciferase or Monster Green® gene so as to create fused gene sequences provided that the coding sequence of the resulting luciferase or Monster Green gene has no more than four deoxynucleotides missing at the affected terminus when compared to the intact luciferase or Monster Green® gene sequence, and (2) insert and remove nucleic acid sequences in furtherance of splicing research predicated on the inactivation or reconstitution of the luminescent activity of the encoded luciferase. In addition, Researchers must do one of the following: (1) use luminescent assay reagents purchased from Promega Corporation for all determinations of luminescence activity resulting from the research use of this product and its derivatives; or, (2) contact Promega Corporation to obtain a license for the use of the product and its derivatives. No other use or transfer of this product or its derivatives is authorized without the express written consent of Promega Corporation including, without limitation, Commercial Use. Commercial Use means any and all uses of this product and derivatives by a party for monetary or other consideration and may include, but is not limited to use in: (1) product manufacture; and (2) to provide a service, information or data; and/or resale of the product or its derivatives, whether or not such product or derivatives are resold for use in research. With respect to such Commercial Use, or any diagnostic, therapeutic or prophylactic uses, please contact Promega Corporation for supply and licensing information. If the purchaser is not willing to accept the conditions of this limited use statement, SABiosciences is willing to accept the return of the unopened product and provide the purchaser with a full refund. However, in the event the product is opened, then the purchaser agrees to be bound by the conditions of this limited use statement. The above license relates to Promega Corporation patents and/or patent applications on improvements to the luciferase and Monster Green® gene.

United States Patent No. 5,292,658 licensed from Millipore Corporation.

Dual-Glo, Dual-Luciferase and Monster Green are trademarks of Promega Corporation.

Opti-MEM is a registered trademark of Life Technologies.

© 2011 QIAGEN, all rights reserved.

www.qiagen.com

Australia ■ Orders 1-800-243-800 ■ Fax 03-9840-9888 ■ Technical 1-800-243-066

Austria ■ Orders 0800-28-10-10 ■ Fax 0800-28-10-19 ■ Technical 0800-28-10-11

Belgium ■ Orders 0800-79612 ■ Fax 0800-79611 ■ Technical 0800-79556

Brazil ■ Orders 0800-557779 ■ Fax 55-11-5079-4001 ■ Technical 0800-557779

Canada ■ Orders 800-572-9613 ■ Fax 800-713-5951 ■ Technical 800-DNA-PREP (800-362-7737)

China ■ Orders 86-21-3865-3865 ■ Fax 86-21-3865-3965 ■ Technical 800-988-0325

Denmark ■ Orders 80-885945 ■ Fax 80-885944 ■ Technical 80-885942

Finland ■ Orders 0800-914416 ■ Fax 0800-914415 ■ Technical 0800-914413

France ■ Orders 01-60-920-926 ■ Fax 01-60-920-925 ■ Technical 01-60-920-930 ■ Offers 01-60-920-928

Germany ■ Orders 02103-29-12000 ■ Fax 02103-29-22000 ■ Technical 02103-29-12400

Hong Kong ■ Orders 800 933 965 ■ Fax 800 930 439 ■ Technical 800 930 425

Ireland ■ Orders 1800 555 049 ■ Fax 1800 555 048 ■ Technical 1800 555 061

Italy ■ Orders 800-789-544 ■ Fax 02-334304-826 ■ Technical 800-787980

Japan ■ Telephone 03-6890-7300 ■ Fax 03-5547-0818 ■ Technical 03-6890-7300

Korea (South) ■ Orders 080-000-7146 ■ Fax 02-2626-5703 ■ Technical 080-000-7145

Luxembourg ■ Orders 8002-2076 ■ Fax 8002-2073 ■ Technical 8002-2067

Mexico ■ Orders 01-800-7742-639 ■ Fax 01-800-1122-330 ■ Technical 01-800-7742-436

The Netherlands ■ Orders 0800-0229592 ■ Fax 0800-0229593 ■ Technical 0800-0229602

Norway ■ Orders 800-18859 ■ Fax 800-18817 ■ Technical 800-18712

Singapore ■ Orders 1800-742-4362 ■ Fax 65-6854-8184 ■ Technical 1800-742-4368

Spain ■ Orders 91-630-7050 ■ Fax 91-630-5145 ■ Technical 91-630-7050

Sweden ■ Orders 020-790282 ■ Fax 020-790582 ■ Technical 020-798328

Switzerland ■ Orders 055-254-22-11 ■ Fax 055-254-22-13 ■ Technical 055-254-22-12

UK ■ Orders 01293-422-911 ■ Fax 01293-422-922 ■ Technical 01293-422-999

USA ■ Orders 800-426-8157 ■ Fax 800-718-2056 ■ Technical 800-DNA-PREP (800-362-7737)

