

September 2006

---

## Important Note

Dear AllPrep user,

### **SDS-PAGE of protein purified with the AllPrep RNA/Protein Kit**

We would like to inform you that the RNA-stabilizing agent in Buffer APL (lysis buffer) causes precipitation of sodium dodecyl sulfate (SDS). Therefore, an SDS-containing buffer should not be used to equilibrate the Protein Cleanup spin column in step 5 of the protocol in the handbook (page 13).

To avoid possible SDS precipitation in applications such as SDS-PAGE, protein purified using the AllPrep RNA/Protein Kit should be cleaned up by acetone precipitation, as described in Appendix D of the handbook (page 25), to remove any traces of the RNA-stabilizing agent.

If you have any questions, please call one of the QIAGEN Technical Service Departments or distributors listed on the back cover of the handbook.

### **Parallel purification of DNA, RNA, and denatured protein**

For experiments requiring DNA, RNA, and denatured protein from a single cell or tissue sample, we recommend using the AllPrep DNA/RNA Mini Kit (cat. no. 80204) in combination with a supplementary protocol ( [www.qiagen.com/literature/protocols/pdf/Ry22.pdf](http://www.qiagen.com/literature/protocols/pdf/Ry22.pdf) ). The purified protein is suitable for SDS-PAGE and western blotting.

Best regards,

QIAGEN

