

## Duplex assay for reliable detection of norovirus RNA and internal control

Lillian Roth, Denis Flügge, Sandy Leifholz, Mario Scherer, and Holger Engel  
QIAGEN GmbH, Hilden, Germany

Reliable detection of norovirus RNA by real-time RT-PCR in the presence of an internal control can be achieved without any optimization steps using the QuantiTect® Virus +ROX Vial Kit and the Rotor-Gene® Q cyclers. Purification of norovirus RNA for real-time RT-PCR detection can be successfully automated on the QIASymphony® SP using the QIASymphony Virus/Bacteria Mini Kit. This combination of QIAGEN instrumentation and chemistries streamlines the entire workflow for analyzing norovirus RNA, from purification to detection.

### Introduction

Duplex, real-time one-step RT-PCR provides a reliable method for detecting norovirus RNA. The norovirus RNA sample is spiked with internal control (IC) RNA, and both RNA targets are then reverse transcribed and amplified in the same tube. Detection of IC, but not of norovirus RNA, would indicate that real-time RT-PCR was successful and that norovirus RNA is either absent or below the limit of detection.

Success in duplex, real-time one-step RT-PCR often requires careful optimization of reaction and cycling conditions. Here we describe a method for duplex analysis that enables precise detection of norovirus RNA over a wide linear range without the need for PCR optimization. In addition, we demonstrate that this method works successfully in combination with automated purification of viral nucleic acids.

### Materials and methods

#### Detection of in vitro transcribed norovirus RNA in a duplex assay

Serial tenfold dilutions of in vitro transcribed norovirus RNA were spiked with in vitro transcribed IC. Duplex, real-time one-step RT-PCR was then performed using the QuantiTect Virus +ROX Vial Kit in combination with 2 primer–probe sets: one specific for norovirus RNA from genogroup II (1), and the other specific for the IC. Reactions were run in triplicate on the Rotor-Gene Q real-time PCR cycler according to the appropriate protocol in the *QuantiTect Virus Handbook*, without any optimization of PCR parameters.

#### Detection of norovirus RNA purified from medium using an automated platform


RNA was purified in triplicate from 5 viral transport medium samples (100 µl sample diluted with 100 µl phosphate-buffered saline) from a commercially available norovirus test panel using the QIASymphony SP in combination with the QIASymphony Virus/Bacteria Mini Kit. Samples 1–4 contained norovirus from genogroup II, genotype 4 (GII.4), originally isolated from stool; sample 5 was the negative control and contained no norovirus. All samples were spiked with in vitro transcribed IC, and duplex, real-time one-step RT-PCR was carried out as described above.



## Results

### Detection of in vitro transcribed norovirus RNA in a duplex assay


The QuantiTect Virus +ROX Vial Kit and Rotor-Gene Q enabled reliable detection of the in vitro transcribed norovirus RNA, and required no optimization of reaction and cycling conditions (Figure 1). Norovirus RNA was detected over a 6-log range with high linearity ( $R^2 = 0.995$ ) and precision (the  $C_T$  values for each set of triplicates showed an overall relative standard deviation of  $1.30 \pm 0.59\%$ ). The IC was also reliably detected, as demonstrated by the comparable  $C_T$  values for all samples (relative standard deviation of  $0.59 \pm 0.49\%$ ).



**Figure 1. Precise detection of in vitro transcribed norovirus RNA over a wide linear range.** Duplex, real-time one-step RT-PCR was carried out using the QuantiTect Virus +ROX Vial Kit and the Rotor-Gene Q according to the supplied protocol and without any PCR optimization. The template was serial tenfold dilutions of in vitro transcribed norovirus RNA (genogroup II; GGII) containing a fixed amount of IC. Each dilution was analyzed in triplicate. **A** Amplification plots for norovirus RNA (one replicate is shown for each template dilution). **B** Amplification plots for the IC (one replicate is shown for each norovirus RNA dilution analyzed). **C** Plot of log template amount versus mean  $C_T$  value, demonstrating the high linearity and precision in detection of norovirus RNA.

## Detection of norovirus RNA purified from medium using an automated platform

Norovirus RNA purified using the QIAAsymphony Virus/Bacteria Mini Kit on the QIAAsymphony SP was reliably detected using the QuantiTect Virus +ROX Vial Kit and the Rotor-Gene Q (Figure 2). For samples 1–4, norovirus RNA was detected with high precision, with the  $C_T$  values for each set of triplicates showing an overall relative standard deviation of  $1.16 \pm 0.86\%$ . As expected, norovirus RNA was not detected in sample 5, the negative control. The IC was detected with comparable  $C_T$  values in samples 1–5 (relative standard deviation of  $1.12 \pm 0.51\%$ ).



**Figure 2. Precise detection of different norovirus GGII RNA isolates.** Duplex, real-time one-step RT-PCR was carried out using the QuantiTect Virus +ROX Vial Kit and the Rotor-Gene Q. The template was viral RNA purified in triplicate from 5 viral transport medium samples on the QIAAsymphony SP using the QIAAsymphony Virus/Bacteria Mini Kit. After purification, each sample was spiked with a fixed amount of IC. The bar chart shows the  $C_T$  values obtained: samples 1–4 contain norovirus RNA; sample 5 is the negative control and contains no norovirus RNA.

## Conclusions

- The QuantiTect Virus +ROX Vial Kit and Rotor-Gene Q allowed precise detection of norovirus RNA over a wide linear range.
- Reliable detection of norovirus RNA and internal control in a duplex assay was achieved without any optimization of reaction and cycling conditions.
- The QIAAsymphony Virus/Bacteria Mini Kit and QIAAsymphony SP provided reliable purification of norovirus RNA from viral transport medium for real-time RT-PCR analysis.
- The combination of QIAAsymphony, Rotor-Gene, and QuantiTect technologies streamlined the entire workflow for purifying and detecting norovirus RNA.

## References

- (1). Höhne, M. and Schreier, E. (2004) Detection and characterization of norovirus outbreaks in Germany: Application of a one-tube RT-PCR using a fluorogenic real-time detection system. *J. Med. Virol.* **72**, 312.

## Ordering Information

Product	Contents	Cat. no.
QuantiTect Virus +ROX Vial Kit (50)*	For 50 x 50 µl reactions: 5x Master Mix (without ROX dye), ROX Dye Solution, 100x RT Mix, RNase-Free Water, and Nucleic Acid Dilution Buffer	211031
Rotor-Gene Q	Real-time PCR cycler with 2, 5, or 6 channels, and with or without high-resolution melting (HRM®) analyzer	Varies†
QIASymphony SP	QIASymphony sample prep module, 1-year warranty on parts and labor	9001297
QIASymphony Virus/Bacteria Mini Kit (192)	For 192 preps of 200 µl each: Includes 2 reagent cartridges and enzyme racks and accessories	931036

\* For use with Rotor-Gene cyclers, the Applied Biosystems® 7500, and cyclers from Bio-Rad, Cepheid, Eppendorf, Roche, and Stratagene. Larger kit sizes available; please inquire. Kits for use with other Applied Biosystems instruments also available; please inquire.

† For ordering information, visit [www.qiagen.com/Products/Rotor-GeneQ.aspx](http://www.qiagen.com/Products/Rotor-GeneQ.aspx).

The QuantiTect Virus +ROX Vial Kit is intended for research use. No claim or representation is intended to provide information for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

**Visit [www.qiagen.com/goto/VirusResearch](http://www.qiagen.com/goto/VirusResearch) to discover comprehensive solutions for virus research!**

Trademarks: QIAGEN®, QIASymphony®, QuantiTect® (QIAGEN Group); Applied Biosystems® (Applied Biosystems Corporation or its subsidiaries); HRM®, Rotor-Gene® (Corbett Research Pty Ltd). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

This product is an Authorized 5' Nuclease Core Kit without Licensed Probe. Its purchase price includes a limited, non-transferable immunity from suit under certain patents owned by Roche Molecular Systems, Inc. or F. Hoffmann-La Roche Ltd, for using only this amount of the product in the practice of the 5' nuclease process solely for the purchaser's own internal research when used in conjunction with Licensed Probe. No right under any other patent claims (such as apparatus or system claims) and no right to use this product for any other purpose is hereby granted expressly, by implication or by estoppel. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

Purchase of this product is accompanied by a limited, non-transferable license under RT and Reverse Transcription-PCR patents owned by Roche Molecular Systems, Inc. and F. Hoffmann-La Roche Ltd to use it for the purchaser's own internal research. No real-time patent rights of any kind, no right under any other patent claims (such as apparatus or system claims), and no right to use this product for any other purpose is hereby granted expressly, by implication or by estoppel.

Certain specific embodiments of the process of multiplex PCR may be covered by patents of third parties in certain countries and may require a license.

1060222 11/2009 © 2009 QIAGEN, all rights reserved.

[www.qiagen.com](http://www.qiagen.com)

Australia ■ 1-800-243-800

Austria ■ 0800-281011

Belgium ■ 0800-79612

Brazil ■ 0800-557779

Canada ■ 800-572-9613

China ■ 86-21-3865-3865

Denmark ■ 80-885945

Finland ■ 0800-914416

France ■ 01-60-920-930

Germany ■ 02103-29-12000

Hong Kong ■ 800 933 965

Ireland ■ 1800 555 049

Italy ■ 800-787980

Japan ■ 03-6890-7300

Korea (South) ■ 080-000-7145

Luxembourg ■ 8002 2076

Mexico ■ 01-800-7742-436

The Netherlands ■ 0800-0229592

Norway ■ 800-18859

Singapore ■ 1800-742-4368

Spain ■ 91-630-7050

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911

USA ■ 800-426-8157

