


QIAseq FastSelect –rRNA HMR, –Globin and Epidemiology Kits

Fast, single-kit solutions for efficient removal of rRNA, globin mRNA and/or bacterial rRNA



- Streamline RNA removal using the fastest kit on the market
- No limits on RNA quantity or quality – compatible with intact and degraded RNA (e.g., FFPE RNA), with input from 1 ng – 1 µg
- Works with any RNA library prep kit
- Nearly complete removal of rRNA and/globin mRNA

One-step QIAseq FastSelect –rRNA HMR and –Globin protocols completed in just 14 minutes



- Single pipetting step – simply add QIAseq FastSelect reagent to RNA sample

- QIAseq FastSelect works with or without fragmentation, which depends on library prep kit and sample type

- 14-minute incubation easily intergrates into existing workflow with no additional accessories or procedures

Single reagent for one-step rRNA and/or globin mRNA from human, mouse, rat and other mammalian species

Overview

One-step Protocol

Fast, 14-min Protocol

Species Compatibility

Robust Performance with FFPE Samples

Reproducibility

Works with any RNA-seq Protocol


Flexibility in Single or Dual-seq Applications

Effective Removal of Bacterial rRNA


Ordering Information

A 14-min, one-step protocol means greater efficiency and speed

A QIAseq FastSelect: one inline reaction step



B QIAseq FastSelect: 90% less time



Data provided by QIAGEN R&D

Save even more time and effort with a 30% faster protocol compared to previous kit version

A single-kit solution for various mammalian samples

QIAseq FastSelect –rRNA HMR Kit is a single solution that removes mitochondrial and cytoplasmic rRNA from human, mouse, rat and other mammalian samples

Species	Reads mapped to rRNA (%)		rRNA removed (%)
	No treatment	After QIAseq FastSelect –rRNA HMR	
Human	91.8	1.6	98.3
Mouse	92.9	2.2	97.6
Rat	93.0	2.6	97.2
Cow	70.0	1.2	98.3
Horse	74.1	1.8	97.6
Hamster	77.5	3.5	95.5
Sheep	87.3	4.0	95.4
Dog	87.8	9.1	89.6
Monkey	90.4	10.6	88.3
Chicken	88.3	16.3	81.5
Rabbit	88.8	16.5	81.4
Pig	87.6	17.3	80.3

Overview

One-step Protocol

Fast, 14-min Protocol

Species Compatibility

Robust Performance with FFPE Samples

Reproducibility

Works with any RNA-seq Protocol

Flexibility in Single or Dual-seq Applications

Effective Removal of Bacterial rRNA

Ordering Information

Reduce the number of wasted reads from FFPE samples

Overview

One-step Protocol

Fast, 14-min Protocol

Species Compatibility

Robust Performance with FFPE Samples

Reproducibility

Works with any RNA-seq Protocol

Flexibility in Single or Dual-seq Applications

Effective Removal of Bacterial rRNA


Ordering Information

Superior rRNA removal compared to other suppliers

(A)

Sample	QIAseq FastSelect rRNA remaining (%)	Ribo-Zero Gold rRNA remaining (%)
S34022 N	0.69	7.46
S34022 C	0.72	7.94
S34019 N	0.89	8.65
S34019 C	0.84	9.05
S33119 N	1.30	6.23
S33119 C	0.50	5.97

(B)



QIAseq FastSelect results in highly efficient removal of rRNA from fragmented samples. Total RNA was isolated from 5 µm normal and cancer lung FFPE sections using the miRNeasy FFPE Kit. The QIAseq Stranded Total RNA Lib Kit was used for RNA-seq library prep. Average % rRNA remaining is provided in **A** and plotted in **B**. Ribo-Zero Gold required substantially more amplification cycles than the QIAseq FastSelect libraries, suggesting some sample may be lost with Ribo-Zero Gold (Data provided by QIAGEN R&D).

Sequence what matters – efficient use of sequencing resources and capacity

Robust removal of rRNA and globin mRNA

Overview

One-step Protocol

Fast, 14-min Protocol

Species Compatibility

Robust Performance with FFPE Samples

Reproducibility

Works with any RNA-seq Protocol


Flexibility in Single or Dual-seq Applications

Effective Removal of Bacterial rRNA

Ordering Information


Increased exon mapping

Exon/intron/intergenic reads



Strong correlation: Replicates

Gene expression analysis



Substantial removal of rRNA and globin

Mapping metrics

Percentage of reads mapped to rRNA		Percentage of TPM mapped to Globin	
No treatment	QIAseq FastSelect -rRNA HMR +/-Globin	No treatment	QIAseq FastSelect -rRNA HMR +/-Globin
88.50	1.16	82.46	0.16

Experimental overview. As starting sample, 100 ng human whole blood total RNA was used. The QIAseq FastSelect -rRNA/Globin Kit was used for RNA removal. RNA-seq library prep was performed using the QIAseq Stranded Total RNA Lib Kit, with sequencing on the NextSeq 550 and read mapping using QIAGEN CLC Genomics Workbench.

Efficient rRNA and globin removal resulting in an increased percentage of reads mapped to exons, with gene expression values from FastSelect-treated samples highly correlative (Log₂ RPKM >0.3)

Overview

One-step Protocol

Fast, 14-min Protocol

Species Compatibility

Robust Performance with FFPE Samples

Reproducibility

Works with any RNA-seq Protocol

Flexibility in Single or Dual-seq Applications

Effective Removal of Bacterial rRNA


Ordering Information

Compatible with QIAGEN, Illumina, KAPA and NEB RNA-seq library prep workflows

A

Sample	rRNA remaining (%)	Protein coding (%)
100 ng UHRR	1.01	89.4
1 µg UHRR	1.09	89.7

B



QIAseq FastSelect integrated into Illumina TruSeq workflow (performance). **A** Average % rRNA remaining and % protein-coding reads shown. **B** rRNA remaining is shown when no rRNA removal is performed (no treatment) and when rRNA removal is performed using QIAseq FastSelect technology (blue = 100 ng and red = 1 µg). QIAseq FastSelect results in highly efficient removal of rRNA. Average % rRNA remaining is provided in **A** and plotted in **B**.

Consistent performance and high downstream reproducibility

Separate reagents for bacteria and human/mouse/rat for flexibility in single or dual-seq applications

QIAseq FastSelect **-rRNA HMR Kit** + QIAseq FastSelect **-5S/16S/23S Kit**
= QIAseq FastSelect Epidemiology Kits



Use with mixed samples for co-expression studies to reduce costs, expand mRNA insights and uncover more biology

Overview

One-step Protocol

Fast, 14-min Protocol

Species Compatibility

Robust Performance with FFPE Samples

Reproducibility

Works with any RNA-seq Protocol

Flexibility in Single or Dual-seq Applications

Effective Removal of Bacterial rRNA

Ordering Information

The QIAseq FastSelect Epidemiology Kit enables fast and effective removal of pan-bacterial 5S/16S/23S rRNA

Overview

One-step Protocol

Fast, 14-min Protocol

Species Compatibility

Robust Performance with FFPE Samples

Reproducibility


Works with any RNA-seq Protocol

Flexibility in Single or Dual-seq Applications

Effective Removal of Bacterial rRNA

Ordering Information

Sample	Reads mapped to total bacterial rRNA (%)		rRNA removed (%)
	No treatment	Treatment	
Supplier T – method ME	95.0	89.9	89.9
Supplier T – method RM	95.0	87.6	87.6
QIAseq FastSelect –5S/16S/23S	93.3	8.5	8.5



Total RNA was isolated from Gut Microbiome Whole cell Mix (ATCC) using the RNeasy PowerMicrobiome Kit (QIAGEN). Stranded transcriptome libraries were then prepared from 1 µg aliquots of the RNA using the QIAseq Stranded Total RNA Lib Kit. For rRNA removal, three methods were used: Supplier T – methods ME and RM, which used a hybrid capture workflow and QIAGEN’s QIAseq FastSelect –5S/16S/23S Kit. Following a bead cleanup, the remaining library prep steps were completed (starting with 1st strand synthesis). Sequencing was performed on a NextSeq 550, and data analysis was performed using QIAGEN CLC Genomics Workbench. (Data provided by QIAGEN R&D).

Removal of over 90% of bacterial rRNA, unlike hybrid capture workflows, which remove less than 10% of the same rRNAs

- Overview
- One-step Protocol
- Fast, 14-min Protocol
- Species Compatibility
- Robust Performance with FFPE Samples
- Reproducibility
- Works with any RNA-seq Protocol
- Flexibility in Single or Dual-seq Applications
- Effective Removal of Bacterial rRNA

Ordering Information

Ordering Information

Product	Species available	Targets	Cat. no.
<u>QIAseq FastSelect –rRNA HMR Kit (96)*</u>	Human, mouse, rat and other mammalian species	rRNA	334387
<u>QIAseq FastSelect –rRNA/Globin Kit (96)*</u>	Human, mouse, rat and other mammalian species	rRNA and globin mRNA	335377
<u>QIAseq FastSelect –Globin Kit (96)*</u>	Human, mouse, rat and other mammalian species	Globin mRNA	334377
<u>QIAseq FastSelect Epidemiology Kit (96)*</u>	Human, mouse, rat and other mammalian species, and bacterial samples	rRNA	333275

* Other kit sizes and/or formats available; see www.qiagen.com for details.

Learn more and see key publications at www.qiagen.com/QIAseqFastSelect/HMR and www.qiagen.com/QIAseqFastSelect/Epidemiology

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, Sample to Insight®, QIAseq® (QIAGEN Group); Illumina®, NextSeq®, Ribo-Zero®, TruSeq® (Illumina, Inc.); KAPA® (Roche Group); NEB®, NEBNext® (New England Biolabs, Inc.).

© 2021 QIAGEN, all rights reserved.

Ordering www.qiagen.com/contact | Technical Support support.qiagen.com | Website www.qiagen.com

Visit Webpage