

# miScript<sup>®</sup> SYBR<sup>®</sup> Green PCR Kit used with miScript miRNA PCR Arrays

The miScript SYBR Green PCR Kit (cat. nos. 218073, 218075, 218076) should be stored at –30 to –15°C upon arrival.

## Further information

- *miScript miRNA PCR Array Handbook*: [www.qiagen.com/HB-0903](http://www.qiagen.com/HB-0903)
- Safety Data Sheets: [www.qiagen.com/safety](http://www.qiagen.com/safety)
- Technical assistance: [support.qiagen.com](http://support.qiagen.com)

## Notes before starting

- This protocol is for real-time PCR using cDNA prepared with the miScript II RT Kit (using miScript HiSpec Buffer) and preamplified cDNA. The reverse-transcription reaction should be diluted as shown in Table 1.
- **IMPORTANT: Only cDNA prepared with miScript HiSpec Buffer should be used for real-time PCR with miScript miRNA PCR Arrays.**
- Prior to preparing real-time PCR reaction mixes, thaw all reagents (2x QuantiTect<sup>®</sup> SYBR Green PCR Master Mix, 10x miScript Universal Primer and RNase-free water) and the cDNA templates prepared using the miScript II RT Kit (and optionally the miScript PreAMP PCR Kit). Gently mix the contents in each tube and briefly centrifuge.
- The format of the miScript miRNA PCR Array is indicated by the last letter of the catalog number. The formats and the real-time cyclers that are supported are detailed in the *miScript miRNA PCR Array Handbook*.

**Table 1. cDNA dilution prior to PCR**

Application	Array	Reaction dilution
Preamplified cDNA	Any miScript Array	Refer to <i>miScript PreAMP Handbook</i>
Pathway profiling	Pathway-Focused miScript miRNA PCR Arrays	Add 200 µl RNase-free water to each 20 µl reverse-transcription reaction
Pathway profiling (high content)	miScript miRNA HC PCR Arrays	Add 90 µl RNase-free water to each 20 µl reverse-transcription reaction
Whole miRNome profiling	miRNome miScript miRNA PCR Arrays	Dilution depends on the number of plates/Rotor-Discs: 1 x 384-well plate or 4 x 96-well plates/Rotor-Discs: add 90 µl RNase-free water to the 20 µl reverse-transcription reaction 2 x 384-well plate or 8 x 96-well plates/Rotor-Discs: add 200 µl RNase-free water to the 20 µl reverse-transcription reaction 3 x 384-well plate or 12 x 96-well plates/Rotor-Discs: add 310 µl RNase-free water to the 20 µl reverse-transcription reaction 4 x 384-well plate or 16 x 96-well plates/Rotor-Discs: add 420 µl RNase-free water to the 20 µl reverse-transcription reaction

1. Prepare a reaction mix according to Table 2 (for Pathway-Focused miScript miRNA PCR Arrays) or Table 3 (for miRNome miScript miRNA PCR Arrays or miScript miRNA HC PCR Arrays), depending on your chosen miScript miRNA PCR Array.

**Table 2. Reaction mix for Pathway-Focused miScript miRNA PCR Arrays**

Array format: Component	384-well (4 x 96)* Formats E, G	96-well Formats A, C, D, F	Rotor-Disc® 100 Format R
2x QuantiTect SYBR Green PCR Master Mix	550 µl	1375 µl	1100 µl
10x miScript Universal Primer	110 µl	275 µl	220 µl
RNase-free water	340 µl	1000 µl	780 µl
Template cDNA†	100 µl	100 µl	100 µl
<b>Total volume</b>	<b>1100 µl</b>	<b>2750 µl</b>	<b>2200 µl</b>

\* Volumes shown are sufficient for one cDNA template. In total, 4 cDNA templates can be analyzed on one 384-well (4 x 96) Pathway-Focused miScript miRNA PCR Array.

† Provides 0.5–1 ng cDNA per well.

**Table 3. Reaction mix for miRNome miScript miRNA PCR Arrays\* and miScript miRNA HC PCR Arrays†**

<b>Array format: Component</b>	<b>384-well Formats E, G</b>	<b>96-well Formats A, C, D, F</b>	<b>Rotor-Disc 100 Format R</b>
2x QuantiTect SYBR Green PCR Master Mix	2050 µl	1375 µl	1100 µl
10x miScript Universal Primer	410 µl	275 µl	220 µl
RNase-free water	1540 µl	1075 µl	855 µl
Template cDNA‡	100 µl	25 µl	25 µl
<b>Total volume</b>	<b>4100 µl</b>	<b>2750 µl</b>	<b>2200 µl</b>

\* Volumes are for a single plate or Rotor-Disc associated with a whole miRNome set. Scale up volumes according to the number of plates/Rotor-Discs to be run.

† For miScript miRNA HC PCR Arrays, use the volumes shown for 384-well, Formats E, G.

‡ Provides 0.5–1 ng cDNA per well.

2. Carefully remove the miScript miRNA PCR Array from its sealed bag.

**Optional for 96-well and 384-well array formats:** If the reaction mix is in a tube, transfer to a loading reservoir, such as the RT<sup>2</sup> PCR Array Loading Reservoir (cat. no. 338162).

3. Add reaction mix to each well of the miScript miRNA PCR Array as follows:

For 384-well miScript miRNA PCR Array: add 10 µl per well.

For 96-well miScript miRNA PCR Array: add 25 µl per well.

For Rotor-Disc miScript miRNA PCR Array: add 20 µl per well.

4. Carefully, tightly seal the miScript miRNA PCR Array with Optical Thin-Wall 8-Cap Strips (Formats A and D), Optical Adhesive Film (Formats C, E, F and G) or Rotor-Disc Heat-Sealing Film (Format R).

5. Centrifuge for 1 min at 1000 x g at room temperature (15–25°C).

**Note:** This step is not necessary for Format R miScript miRNA PCR Arrays in Rotor-Discs.

6. Program the real-time cycler according to Table 4.

**Note:** Perform dissociation (melting) curve analysis of the PCR product(s) to verify their specificity and identity. Dissociation curve analysis is an analysis step built into the software of real-time cyclers. Follow the instructions provided by the supplier.

**Table 4. Cycling conditions**

Step	Time	Temperature	Additional comments
<b>Initial activation step</b>	15 min	95°C	HotStarTaq® DNA Polymerase is activated by this heating step.
<b>3-step cycling:***</b>			
Denaturation	15 s	94°C	
Annealing	30 s	55°C	
Extension <sup>§</sup>	30 s	70°C	Perform fluorescence data collection.
Cycle number	40 cycles <sup>†</sup>		Cycle number depends on the amount of template cDNA and abundance of the target.

\* For Bio-Rad® models CFX96™ and CFX384™: adjust the ramp rate to 1°C/s.

† For Eppendorf® Mastercycler® ep realplex models 2, 2S, 4 and 4S: for the Silver Thermoblock, adjust the ramp rate to 26%; for the Aluminum Thermoblock, adjust the ramp rate to 35%.

‡ If using a Roche® LightCycler® 480, adjust the ramp rate to 1°C/s.

§ Due to software requirements, the fluorescence detection step must be at least 30 s with the ABI PRISM® 7000 or 34 s with the Applied Biosystems® 7300 and 7500.

† If using a Roche LightCycler 480, use 45 cycles.

7. Place the plates in the real-time cycler and start the cycling program.



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN®, Sample to Insight®, QuantiTect®, HotStarTaq®, miScript®, Rotor-Disc® (QIAGEN Group); SYBR® (Molecular Probes, Inc.); ABI PRISM®, Applied Biosystems® (Applied Biosystems or its subsidiaries); Bio-Rad®, CFX96™, CFX384™ (Bio-Rad Laboratories, Inc.); Eppendorf®, Mastercycler® (Eppendorf AG); Roche®, LightCycler® (Roche Group). 1102264 04/2016 HB-0845-003 © 2016 QIAGEN, all rights reserved.