

# Labware List — QIASymphony® DSP HPV Media Kit

## Sample and eluate tubes/racks that can be used with the QIASymphony DSP HPV Media Kit and the QIASymphony SP (software version 4.0; labware package SOW-516-6)

### Version management

This document is the QIASymphony DSP HPV Media Kit Labware List for labware package SOW-516-6, Version 1, Revision 1.

**Important:** Before using this Labware List, confirm that it is consistent with the labware package installed on your QIASymphony system.

**Note:** The QIASymphony DSP HPV Media Kit is not available in the United States. Please visit [www.qiagen.com](http://www.qiagen.com) for product availability.

### Legend

- Required: These tubes are required for use in combination with this protocol.
- User defined: User is responsible for defining and/or validating as the user deems appropriate.
- Not recommended: These tubes are not for use in combination with this protocol.
- x/■ Minimum sample volume ( $\mu$ l) required per sample per protocol (denoted by x)/clot detection possible.
- x/□ Minimum sample volume ( $\mu$ l) required per sample per protocol (denoted by x)/clot detection not possible.

**Note:** Be aware that other tubes are available that can be used in combination with other QIASymphony Kits, but are not recommended for any protocols that can be used with this kit.

### "Sample" drawer, tube carrier

Supplier	Material	Cat. no.	Name in touchscreen	Insert	Protocol	
					PC2500_HC2	SP800_HC2
BD*	14 ml Falcon polystyrene round-bottom tube 17 x 100 mm	352051	BD#352051 FalconPP 17x100	No insert required	x/□	
QIAGEN	2 ml Micro tube, PP, NON-SKIRTED	997055	QIA#997055 T2.0 Conical blue	03/3B		x/□

\* BD was the previous supplier of this tube and Corning Inc. is now the new supplier.

January 2016



Sample & Assay Technologies

# Labware List — QIASymphony DSP HPV Media Kit

## "Eluate" drawer

						Protocol	
Supplier	Material	Example cat. no.	Category	Name in touchscreen	Adapter on Elution slots	PC2500_HC2	SP800_HC2
QIAGEN	Hybridization microplate	6000-1203	Micro Plate	QIA#6000-1203 *MTP96 RB	96-Well Round Bottom QS (cooling slot 1) Elution slots 2-4 are not used in this procedure		

\* Indicates labware that can be cooled using a cooling adapter with or without bar code.

For ordering information, visit [www.qiagen.com/goto/QIASymphony](http://www.qiagen.com/goto/QIASymphony).

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, QIASymphony® (QIAGEN Group); BD® (Becton, Dickinson and Company); Corning®, Falcon® (Corning, Inc.).  
Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Jan-16 HB-1554-L03-001 © 2013–2016 QIAGEN, all rights reserved.

