

## Quick-Start Protocol

# QIAseq<sup>®</sup> miRNA Library Kit

## Part 3b: Library amplification using tube indices

### Further information

- When using Illumina<sup>®</sup> NGS systems, refer to the *QIAseq miRNA Library Kit Handbook: Illumina NGS Systems*: [www.qiagen.com/HB-2157](http://www.qiagen.com/HB-2157)
- When using Thermo Fisher Scientific<sup>®</sup> NGS systems, refer to the *QIAseq miRNA Library Kit Handbook: Thermo Fisher Scientific NGS Systems*: [www.qiagen.com/HB-2573](http://www.qiagen.com/HB-2573)
- Safety Data Sheets: [www.qiagen.com/safety](http://www.qiagen.com/safety)
- Technical assistance: [support.qiagen.com](http://support.qiagen.com)

### Notes before starting

- Prepare the reagents according to the *QIAseq miRNA Library Kit* handbooks
  - Do not vortex reactions or reagents unless instructed
1. On ice, prepare the library amplification reaction using Table 1 for Illumina tube indices (cat. nos. 331592 and 331595) or Table 2 for Thermo Fisher Scientific tube indices (cat no. 331582).

**Table 1. Setup of library amplification reactions using Illumina tube indices (331592 and 331595)**

Component	Volume/rxn
Product from cDNA Cleanup protocol	15 µl
QIAseq miRNA NGS Library Buffer	16 µl
HotStarTaq® DNA Polymerase	3 µl
QIAseq miRNA NGS ILM Library Forward Primer	2 µl
QIAseq miRNA NGS ILM IPD1 through IPD48 (index primer)*	2 µl
Nuclease-free water	42 µl
<b>Total volume</b>	<b>80 µl</b>

\* Up to 48 different QIAseq miRNA NGS ILM IPD primers (index primers) are available for use.

**Table 2. Setup of library amplification reactions using Thermo Fisher Scientific tube indices (331582)**

Component	Volume/rxn
Product from cDNA Cleanup protocol	15 µl
QIAseq miRNA NGS Library Buffer	16 µl
HotStarTaq DNA Polymerase	3 µl
QMI TF IP1 through IP12 (index primer)*	2 µl
QMI TF Lib Rev Primer	2 µl
Nuclease-free water	42 µl
<b>Total volume</b>	<b>80 µl</b>

\* Up to 12 different QMI TF IP primers (index primers) are available for use.

2. Program the thermal cycler according to Table 3. The correct number of cycles depends on the original RNA input and is shown in Table 4.

**Table 3. Library amplification protocol**

Step	Time	Temperature
<b>Hold</b>	15 min	95°C
<b>3-step cycling</b>		
Denaturation	15 s	95°C
Annealing	30 s	60°C
Extension	15 s	72°C
Cycle number (see Table 4)		
<b>Hold</b>	2 min	72°C
<b>Hold</b>	∞*	4°C

\* At least 5 min.

**Table 4. Cycles of library amplification**

Original RNA input (total RNA)	Cycle number
500 ng	13
100 ng	16
10 ng	19
1 ng	22
Serum/plasma	22

3. Place the library amplification reaction in the thermal cycler and start the run.

**Important:** Upon completion of the protocol, hold at 4°C for at least 5 min.

4. Add 75 µl of QMN Beads into the tubes.

5. Briefly centrifuge the 80 µl library amplification reactions, and then transfer 75 µl of the library to the tubes containing the QMN Beads. Vortex for 3 sec, and then centrifuge briefly.

6. Incubate for 5 min at room temperature (15–25°C).

7. Place the tubes on a magnet stand for ~4 min or until the beads have fully migrated.

8. Transfer 145 µl of the supernatant to new tubes.

**Important:** Do not discard the supernatant at this step.

9. To the 145 µl supernatant, add 130 µl of QMN Beads. Vortex for 3 sec, and then briefly centrifuge.

10. Incubate for 5 min at room temperature.

11. Place the tubes on a magnet stand until the beads have fully migrated.

**Note:** Ensure beads have fully migrated before proceeding.

12. Discard the supernatant, keep the beads, and keep the tubes on the magnet stand.

13. Add 200 µl of 80% ethanol into the tubes, and then immediately remove and discard the ethanol wash.

---

14. Repeat the 80% ethanol wash: Add 200  $\mu$ l of 80% ethanol into the tubes again, and then immediately remove and discard this second ethanol wash.

**Note:** It is important to completely remove all traces of the ethanol after the second wash. To do this, briefly centrifuge the tubes, and then return them to the magnetic stand. Remove the ethanol with a 200  $\mu$ l pipette first, and then use a 10  $\mu$ l pipette to remove any residual ethanol.

15. With the bead pellets still on the magnetic stand, air dry them at room temperature for 10 min.

**Note:** Visually inspect to ensure that each pellet is completely dry.

16. With the bead pellets still on the magnetic stand, elute the DNA by adding 17  $\mu$ l of nuclease-free water into the tubes. Subsequently close the caps and remove the tubes from the magnetic stand.

17. Carefully pipet the mixture up and down until all the beads are thoroughly resuspended. Briefly centrifuge, and then incubate at room temperature for 2 min.

18. Place the tubes on the magnetic stand for ~2 min or until the beads have fully migrated.

**Note:** Ensure that the beads have fully migrated before proceeding.

19. Transfer 15  $\mu$ l of eluted DNA to new tubes. This is the miRNA sequencing library.

20. Store the completed miRNA sequencing library at  $-20^{\circ}\text{C}$  in a constant-temperature freezer, or proceed to the *miRNA Library Pre-Sequencing QC* protocol in the *QIAseq miRNA Library Kit Handbook* for either the Illumina NGS systems or the Thermo Fisher NGS systems.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN<sup>®</sup>, Sample to Insight<sup>®</sup>, QIAseq<sup>®</sup>, HotStarTaq<sup>®</sup> (QIAGEN Group); Illumina<sup>®</sup> (Illumina, Inc.); Thermo Fisher Scientific<sup>®</sup> (Thermo Fisher Scientific or its subsidiaries). 1114600 08/2018 HB-2577-001 © 2018 QIAGEN, all rights reserved.