

QIAamp® Viral RNA Mini Accessory Set

For use with multiple QIAamp Viral RNA Mini Kits on the QIAcube™

Set Contents

QIAamp Viral RNA Mini Accessory Set	
Catalog no.	1048147
Buffer AVL	2 x 31 ml
Buffer AVE	20 x 2 ml
Carrier RNA	310 µg
Product Sheet	1

Storage

All buffers can be stored at room temperature (15–25°C). Lyophilized carrier RNA is stable for up to 1 year when stored at room temperature (15–25°C). Carrier RNA can only be dissolved in Buffer AVE; dissolved carrier RNA should be immediately added to Buffer AL as described in the *QIAamp Viral RNA Mini Handbook*. This solution should be prepared fresh, and is stable at 2–8°C for up to 48 hours. Unused portions of carrier RNA dissolved in Buffer AVE should be frozen in aliquots at –20°C.

Product Use Limitations

The QIAamp Viral RNA Mini Accessory Set contains additional buffers and reagents for use with at least **11 x QIAamp Viral RNA Mini Kits (50)**, catalog number 52904, or **5 x QIAamp Viral RNA Mini Kits (250)**, catalog number 52906, on the QIAcube. Depending on run parameters, the set may be sufficient for use with more than 11 or 5 kits, respectively. For setting up the QIAcube for viral RNA purification using the QIAamp Viral RNA Mini Kit, follow the instructions in the Protocol Sheet, which can be downloaded from www.qiagen.com/MyQIAcube.

The QIAamp Viral RNA Mini Accessory Set is intended for molecular biology applications. This product is not intended for the diagnosis, prevention, or treatment of disease.



Safety Information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate material safety data sheets (MSDSs).

These are available online in convenient and compact PDF format at

www.qiagen.com/Support/MSDS.aspx where you can find, view, and print the MSDS for each QIAGEN® kit and kit component.

CAUTION: DO NOT add bleach or acidic solutions directly to waste containing Buffer AVL.

Buffer AVL contains guanidine salts, which can form highly reactive compounds when combined with bleach. If liquid containing Buffer AVL is spilt, clean with suitable laboratory detergent and water. If the spilt liquid contains potentially infectious agents, clean the affected area first with laboratory detergent and water, and then with 1% (v/v) sodium hypochlorite.

The following risk and safety phrases apply to components of the QIAamp Viral RNA Mini Accessory Set:

Buffer AVL

Contains guanidine thiocyanate: harmful. Risk and safety phrases:* R20/21/22-32, S13-26-36-46

24-hour emergency information

Emergency medical information in English, French, and German can be obtained 24 hours a day from:

Poison Information Center Mainz, Germany

Tel: +49-6131-19240

* R20/21/22: Harmful by inhalation, in contact with skin and if swallowed; R32: Contact with acids liberates very toxic gas; S13: Keep away from food, drink, and animal feedingstuffs; S26: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice; S36: Wear suitable protective clothing; S46: If swallowed, seek medical advice immediately and show container or label.

Trademarks: QIAGEN®, QIAamp®, QIAcube™ (QIAGEN Group).

1063054 04/2010 © 2007 QIAGEN, all rights reserved.

www.qiagen.com

Australia = 1800-243-800

Finland = 0800-914416

Korea (South) = 080-000-7146

Sweden = 020-790282

Austria = 0800-28-10-10

France = 01-60-920-920

Luxembourg = 8002-2076

Switzerland = 055-254-22-11

Belgium = 0800-79612

Germany = 02103-29-12000

Mexico = 01-800-7742-639

UK = 01293-422-911

Brazil = 55-11-5079-4000

Hong Kong = 800-933-965

The Netherlands = 0800-0229592

USA = 800-426-8157

Canada = 800-572-9613

Ireland = 1800-555-049

Norway = 800-18859

China = 800-988-0326

Italy = 02-33430411

Singapore = 1800-742-4362

Denmark = 80-885945

Japan = 03-5547-0811

Spain = +34-91-630-7050



Sample & Assay Technologies