

AllPrep[®] DNA/mRNA Nano Kit

Simultaneous DNA and mRNA purification from low-biomass samples



Overview

Procedure

Minimal Cell Requirement

HMW DNA

Long-Range PCR


Sequencing-Ready DNA

Sequencing-Ready RNA

Competitive Performance

Ordering Information

Genomic and transcriptomic profiling from the same sample



Magnetic-bead system maximizes mRNA capture and isolates HMW gDNA

Download Quick-Start Protocol

To learn more, go to our webpage

HMW: high molecular weight.

Overview

Procedure

Minimal Cell Requirement

HMW DNA

Long-Range PCR

Sequencing-Ready DNA

Sequencing-Ready RNA


Competitive Performance

Ordering Information

Requires minimal amounts of cells

Down to 5 cells

PIK3CA E545K detection



AF: allele frequency.

Pyrosequencing[®] for PIK3CA hotspot mutation E545K was done on gDNA isolated using the AllPrep DNA/mRNA Nano Kit.




Mutation detection is feasible down to **five MCF-7 cells** spiked in 5 ml of blood.



Download Handbook

- Overview
- Procedure
- Minimal Cell Requirement
- HMW DNA
- Long-Range PCR
- Sequencing-Ready DNA
- Sequencing-Ready RNA
- Competitive Performance
- Ordering Information

Extracts high-molecular-weight gDNA



≥60.000 bp

Genomic DNA was isolated from various numbers of MCF-7 cells and analyzed with an automated electrophoresis instrument.



Genomic DNA isolated with the AllPrep DNA/mRNA Nano Kit has high molecular weights of **60,000 bp or more**.

Overview

Procedure

Minimal Cell Requirement

HMW DNA

Long-Range PCR


Sequencing-Ready DNA

Sequencing-Ready RNA

Competitive Performance

Ordering Information

Highly suitable for long-range PCR



Genomic DNA was purified from ≤ 2000 cells and amplified using the UltraRun LongRange PCR Kit.

A 4 kb amplicon from P53 was amplified from as few as 20 cells.



The resulting PCR fragments are highly suitable for a variety of applications, such as cloning, genome mapping and sequencing.

Overview

Procedure

Minimal Cell Requirement

HMW DNA

Long-Range PCR


Sequencing-Ready DNA

Sequencing-Ready RNA

Competitive Performance

Ordering Information

Allows the detection of cell-line variants from low-biomass samples



VAF: variant allele frequency.



All variants described for the selected cell lines were **detected with the expected allele frequency**.

None of them were found in the healthy donor control.

mRNA and gDNA were isolated from different cell lines that were mixed and spiked into healthy background.

Isolated DNA was used for library preparation, and NGS was performed.

Overview

Procedure

Minimal Cell Requirement

HMW DNA

Long-Range PCR


Sequencing-Ready DNA

Sequencing-Ready RNA

Competitive Performance

Ordering Information

Isolates mRNA suitable for single-cell RNA sequencing



mRNA from single cells was purified using the AllPrep mRNA/DNA Nano procedure.


RNA sequencing was done with the QIAseq[®] UPX 3' Transcriptome Kit.



The UPX gene-expression profiles differentiate single cells from Breast CTC, Ovary CTC, and healthy donor controls.

- Overview
- Procedure
- Minimal Cell Requirement
- HMW DNA
- Long-Range PCR
- Sequencing-Ready DNA
- Sequencing-Ready RNA
- Competitive Performance
- Ordering Information

gDNA and mRNA yield significantly exceeded all other competitor procedures



mRNA was analyzed using the AdnaTest BreastCancerDetect Kit and yields were determined by densitometry.

Matched gDNA yield was determined by qPCR.

- Overview
- Procedure
- Minimal Cell Requirement
- HMW DNA
- Long-Range PCR
- Sequencing-Ready DNA
- Sequencing-Ready RNA
- Competitive Performance
- Ordering Information**

Ordering Information

Product	Contents	Cat. no.
AllPrep DNA/mRNA Nano Kit (12)	For 12 preps: Oligo-dT magnetic beads, DNA-isolation beads, collection tubes, RNase-free water and buffers	80272
AdnaMag-S	Magnetic rack for 8 tubes, 1.5 ml	399911

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, Sample to Insight®, QIAseq®, AllPrep®, Pyrosequencing® (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

© 2020 QIAGEN, all rights reserved. PROM-16719-001 10/2020

Ordering www.qiagen.com/shop | Technical Support support.qiagen.com | Website www.qiagen.com

[Visit Webpage](#)