

Quick-Start Protocol

Pf/Non-Pf Detection Assay

The Pf/Non-Pf Detection Assay kit (cat. no. 224113) is intended for molecular biology application for human epidemiological research using qPCR. Upon receipt, the kit should be stored at -30 to -15°C in a constant-temperature freezer and protected from light.

Further information

- *QIAprep& Plasmodium Kit Handbook*: www.qiagen.com/HB-3663
- Safety Data Sheets: www.qiagen.com/safety
- Technical assistance: support.qiagen.com
- *Pv/Pm/Po/Pk Detection Assay Kit Quick-Start Protocol*: www.qiagen.com/HB-3671

Notes before starting

- The assay is compatible with all qPCR cyclers with the mentioned dye detection capabilities.

Each tube contains primers and probes for detection of *Plasmodium falciparum* (*P. falciparum*) and *Plasmodium spp.*. Additionally, an assay for the detection of human RNase P is included as sampling and in-process PCR control.

- The *P. falciparum* signal can be detected in the FAM channel;
- Other *Plasmodium* species (*Plasmodium spp.*) can be detected in the Cy5 channel; and

- The signal for RNase P can be detected in the HEX channel.
- Each tube of the Pf/Non-Pf Detection Assay Kit contains a mixture of primers and probes at a 20x concentration; therefore, only 1 µL of the assay is needed per 20 µL reaction. Each tube contains sufficient reagents for 100 reactions.

Procedure

1. Ensure that the qPCR cyclers has the following detection channels: FAM, HEX, and Cy5.
2. Follow the reaction setup and cycling condition of the QIAprep& Plasmidium Kit when using this assay. Refer to www.qiagen.com/HB-3663
3. Select the relevant detection channels (FAM, HEX, and Cy5) in the PCR cyclers during run setup. See Table 1 for alternative channel names.

Table 1. Dyes, typical channel names and their alternative names

Dye	Alternative channel names in qPCR cyclers	
	QIAquant	RGQ
FAM	Blue channel	Green channel
HEX/JOE/VIC	Green channel	Yellow channel
Cy5	Red channel	Red channel

Result interpretation

Table 2. Result interpretation in the different channels from analyzed human samples

Signal in FAM* (<i>P. falciparum</i>)	Signal in HEX† (In-process control)	Signal in Cy5‡ (Non <i>falciparum</i>)	Status	Result
+	+	–	VALID	Positive signal for Pf
–	+	+	VALID	Positive signal for Non-Pf
+	–	–	VALID	Positive signal for Pf
–	–	+	VALID	Positive signal for Non-Pf
+	+	+	VALID	Sample with a mixed infection
+	–	+	VALID	Sample with a mixed infection
–	+	–	VALID	Negative signal for Plasmodium
–	–	–	INVALID	Repeat test

* A positive signal in the FAM channel indicates the presence of *Plasmodium falciparum* in the analyzed human sample

† A signal in the HEX channel indicates that the PCR reaction was successful and human sample material was present in the reaction.

‡ Non-*falciparum* species: *Plasmodium vivax*, *Plasmodium ovale*, *Plasmodium knowlesi*, and/or *Plasmodium malariae*.

Note: If no signal is observed, then the test is invalid and needs to be repeated.

Summary of the *in silico* investigation on potential cross-reactive species

Due to the close relationship among different *Plasmodium* species, the possibility of cross-reactivity was investigated. The following species might be detected by the assays based on an *in silico* analysis (Table 3 and Table 4):

Table 3. Detected species in Plasmodium falciparum assay*Plasmodium sp. gorilla**Plasmodium sp. chimpanzee**Plasmodium billcollinsi**Plasmodium lutzi***Table 4. Detected species in Non-falciparum assay**

<i>Plasmodium brasilianum</i>	<i>Plasmodium gonderi</i>	<i>Polychromophilus murinus</i>	<i>Plasmodium minuoviride</i>
<i>Plasmodium inui</i>	<i>Plasmodium sp. chimpanzee</i>	<i>Leucocytozoon caulleryi</i>	<i>Plasmodium azurophilum</i>
<i>Plasmodium cynomolgi</i>	<i>Plasmodium sp. gorilla</i>	<i>Polychromophilus melanipherus</i>	<i>Haemoproteus pyodactylii</i>
<i>Plasmodium fragile</i>	<i>Plasmodium hylobati</i>	<i>Plasmodium juxtannucleare</i>	<i>Plasmodium elongatum</i>
<i>Plasmodium simium</i>	<i>Plasmodium simiovale</i>	<i>Plasmodium leucocytica</i>	<i>Plasmodium billcollinsi</i>
<i>Plasmodium fieldi</i>	<i>Plasmodium coatneyi</i>	<i>Plasmodium circumflexum</i>	<i>Plasmodium lacertiliae</i>
<i>Plasmodium sp. pongo</i>	<i>Haemosporida</i>	<i>Plasmodium megalotrypa</i>	<i>Leucocytozoon sp.</i>
<i>Hepaticystis sp.</i>	<i>Polychromophilus sp.</i>	<i>Plasmodium koreafense</i>	<i>Plasmodium reichenowi</i>

Document Revision History

Date	Changes
02/2025	Initial release

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN®, Sample to Insight®, QIAcard®, QIAprep&®, Uni-Core® (QIAGEN Group); Whatman® (Whatman International Limited Corporation). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

02/2025 HB-3669-001 © 2025 QIAGEN, all rights reserved.