

EZ2[®] Connect Fx Recovery Procedure Instruction Manual

For use with EZ1&2[™] DNA Investigator[®] Kit



9003220



QIAGEN GmbH, QIAGEN Strasse 1, 40724 Hilden, GERMANY



1126342

Contents

1	Introduction	3
2	Equipment and Reagents to be Supplied By User (for Manual Recovery)	4
3	General Instructions	5
	3.1 System description	5
4	Index for Recovery Process	7
5	Manual Sample Recovery and Processing Procedure	8
	5.1 Manual Recovery from unprocessed sample lysates	8
	5.2 Manual Recovery 1	8
	5.3 Manual Recovery 2	9
	5.4 Manual Recovery 3	9
	5.5 Manual Recovery 4	9
	5.6 Manual Recovery from rinse step	10
	5.7 Recovery of Samples in the elution step	10
6	Semi-Automated Recovery Procedure	11
	6.1 Recovery from unprocessed sample lysates	11
	6.2 Recovery of samples in binding step	11
	6.3 Recovery of samples in wash 1–3	11
	6.4 Recovery of samples in the rinse step	11
	6.5 Recovery of samples in the elution step	12
	Document Revision History	13

1 Introduction

EZ2 Connect Fx recovery mechanism is available to ensure that sample recovery is possible after an unexpected protocol execution problem. This manual describes how to perform the recovery process (both manually and semi-automatically) by using the QIAamp® DNA Investigator Kit. The semi-automatic execution saves time and allows you to continue the EZ2 Connect Fx process in a few steps. The manual recovery process allows you to purify your samples in the usual QIAGEN® standard, even if the EZ2 Connect Fx instrument becomes inoperable.

2 Equipment and Reagents to be Supplied By User (for Manual Recovery)

When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles. For more information, consult the appropriate safety data sheets (SDSs), available from the product supplier.

- QIAamp DNA Investigator Kit (cat. no. 56504)
- Ethanol (96–100%)
- Thermomixer, heated orbital incubator, heating block, or water bath
- Microcentrifuge with rotor for 2 mL tubes

3 General Instructions

3.1 System description

Note: If the device was switched off unintentionally: Start the device. The recovery screen should appear. If it is not possible to start the device, please contact QIAGEN Technical Support.

Note: If the sample remains in the tip after the device switched off: Place a tube under the tip and remove the tip from the pipette head. The liquid will now run out of the tip. If you have any problems with this step, please contact QIAGEN Technical Support.

Before you proceed, please read these general instructions first.




Figure 1. Position key for sample location from recovery screen description.




Figure 2. Recovery screen (displayed after unintended protocol abortion).

Regardless of the respective recovery step, the following steps must be performed first:

1. The display message shows important information. Note the position of the sample and the step to be performed in the manual.
2. Open the hood.
3. Remove and keep the sample containing tubes/cartridges.
 - Important:** Label the sample tubes/cartridges and be careful not to mix up their order.
4. Proceed with the indicated manual recovery step using the table in Section 4.

4 Index for Recovery Process

Indicated Step	Manual Recovery Step	Semi-automatic recovery step	Process status
1	5.1	6.1	Sample untouched
2	5.2	6.2	Buffer MTL added to the sample, beads may or may not have been added
3	5.3	6.3	During wash 1
4	5.4	6.3	During wash 2
5	5.5	6.3	During wash 3
6	5.6	6.4	During rinse step
7	5.7	6.5	During elution step

5 Manual Sample Recovery and Processing Procedure

5.1 Manual Recovery from unprocessed sample lysates

1.
 - 1a. Recovery of Trace, Trace TD, Normalization, or Normalization TD protocols:**
Add 200 μL buffer AL to the sample lysate. Add 100 μL EtOH, close the lid, and mix by pulse-vortexing for 10 s.
 - 1b. Recovery of Large Volume protocols with 400 μL buffer MTL added to the sample tube:**
No addition of buffer AL required. Add 250 μL EtOH, close the lid, and mix by pulse-vortexing for 10 s.
 - 1c. Recovery of Large Volume protocols with no buffer MTL added:**
Add 500 μL buffer AL to the sample lysate. Add 250 μL EtOH, close the lid, and mix by pulse-vortexing for 10 s.
2. Carefully transfer the sample to a QIAamp MinElute[®] column (in a 2 mL collection tube) without wetting the rim. Large Volume recovery requires to load the sample in two steps, with a maximum of 700 μL applied to the column.
3. Centrifuge the samples at 6000 $\times g$ for 1 min.
4. Place the QIAamp MinElute column in a clean 2 mL collection tube and discard the collection tube containing the flow-through.
5. Add 500 μL Buffer AW1 from the QIAamp DNA Investigator Kit to the QIAamp MinElute column and continue with Step 2 of "Manual Recovery 1" section.

5.2 Manual Recovery 1

1. Carefully transfer 700 μL lysate from the recovered sample to a QIAamp MinElute column (in a 2 mL collection tube), without wetting the rim.
2. Close the lid and centrifuge at 6000 $\times g$ for 1 min.
3. Carefully discard the flow-through from the collection tube then place the QIAamp MinElute column back into the collection tube.
4. Carefully apply the remaining lysate to the QIAamp MinElute column without wetting the rim.
5. Close the lid and centrifuge at 6000 $\times g$ for 1 min.

Note: To save the lysates from Large Volume protocols, a third loading step might be required.

Note: Ensure the entire lysate, including magnetic particles, has been transferred to the QIAamp MinElute column.
6. Place the QIAamp MinElute column in a clean 2 mL collection tube and discard the collection tube containing the flow-through.
7. Add 500 μL Buffer AW1 from the QIAamp DNA Investigator Kit to the QIAamp MinElute column and continue with Step 2 of "Manual Recovery 2" section.

5.3 Manual Recovery 2

1. Carefully transfer 700 μL sample, including all magnetic particles to a QIAamp MinElute column (in a 2 mL collection tube), without wetting the rim.
2. Close the lid and centrifuge at 6000 $\times g$ for 1 min.
3. Carefully discard the flow-through from the collection tube and place the QIAamp MinElute column back into the collection tube.
4. Add 700 μL buffer AW2 from the QIAamp DNA Investigator Kit to the QIAamp MinElute column and continue with Step 2 of "Manual Recovery 3" section.

5.4 Manual Recovery 3

1. Carefully transfer 700 μL sample, including all magnetic particles to a QIAamp MinElute column (in a 2 mL collection tube), without wetting the rim.
2. Close the lid and centrifuge at 6000 $\times g$ for 1 min.
3. Carefully discard the flow-through from the collection tube and place the QIAamp MinElute column back into the collection tube.
4. Add 700 μL ethanol (96–100%) to the QIAamp MinElute column and continue with Step 2 of "Manual Recovery 4" section.

5.5 Manual Recovery 4

1. Carefully transfer 700 μL sample, including all magnetic particles to a QIAamp MinElute column (in a 2 mL collection tube), without wetting the rim.
2. Close the lid and centrifuge at 6000 $\times g$ for 1 min.
3. Carefully discard the flow-through from the collection tube and place the QIAamp MinElute column back into the collection tube.
4. Centrifuge at full speed (20,000 $\times g$) for 3 min. to dry the membrane completely.
5. Carefully open the lid of the QIAamp MinElute column and incubate at 56°C for 3 min.
6. Place the QIAamp MinElute column in a clean 1.5 mL microcentrifuge tube (not provided; without lid), and discard the collection tube containing the flow-through.
7. Add Elution Buffer ATE (i.e., 20–200 μL) to the QIAamp MinElute column.
8. Close the lid and incubate in a thermomixer or heated orbital incubator at 56°C, with shaking at 900 rpm for 5 min. Centrifuge at full speed (20,000 $\times g$) for 1 min.

5.6 Manual Recovery from rinse step

1. Split the 1200 μL rinse volume and beads into two 600 μL . Add 600 μL buffer AL and 300 μL EtOH to each. Close the lid, and mix by pulse-vortexing for 10 s.
2. Carefully transfer 750 μL sample to QIAamp MinElute columns (in a 2 mL collection tube), without wetting the rim.
3. Centrifuge at $6000 \times g$ for 1 min.
4. Repeat Steps 2 and 3 until the entire lysate of both tubes is loaded onto the QIAamp MinElute column.
5. Place the QIAamp MinElute column in a clean 2 mL collection tube and discard the collection tube containing the flow-through.
6. Add 500 μL Buffer AW1 from the QIAamp DNA Investigator Kit to the QIAamp MinElute column, and continue with Step 2 of "Manual Recovery 2" section.

5.7 Recovery of Samples in the elution step

1. Place samples in a shaker for 5 min at 900 rpm.
2. Separate magnetic beads by centrifugation, or by using a magnetic stand.
3. Transfer the eluate to a new tube.

6 Semi-Automated Recovery Procedure

6.1 Recovery from unprocessed sample lysates

1. Add new labware.
2. Restart the run.

6.2 Recovery of samples in binding step

Trace or Trace Tip Dance protocols

1. Recover the lysate and all beads, if already added. The expected volume is approximately 900 μL .
2. Use as samples in the Large Volume protocol (old version, with no MTL from well 10).

Large Volume protocols (new protocols using MTL from well 10)

1. Recover the lysate and all beads, if already added. The maximum expected volume is approximately 1600 μL . Split into two fractions of maximum 900 μL .
2. Use as samples in a Large Volume protocol (old version, with no MTL from well 10).

Normalization, or Normalization Tip Dance protocols:

1. No normalized recovery is possible. In case no sample material is available to re-run the normalization protocols, the recovery procedure for Trace and Trace Tip Dance protocols can be applied.

6.3 Recovery of samples in wash 1–3

1. Carefully discard the supernatant of the wash buffer in the Cartridge (beads should be covered with washing buffer).
2. Re-suspend the beads by pipetting up and down.
3. Transfer the remaining beads from the Cartridge into the elution tube.
4. To get a bead pellet, centrifuge the elution tube at $6000 \times g$ for 1 min.
5. Carefully discard wash buffer from the elution tube, until approximately 200 μL are left.
6. Resuspend all beads in a volume of 200 μL of the wash buffer and transfer buffer and beads to a new 2 mL sample tube.
7. Use sample in a Trace protocol.

6.4 Recovery of samples in the rinse step

1. Recover all beads and the entire water used for the rinse.
2. Split into two fractions of approximately 500 μL and use as samples in a Large Volume protocol. Pool eluates of split samples after the run.

6.5 Recovery of samples in the elution step

1. Place samples in a shaker for 5 min at 900 rpm.
2. Separate magnetic beads by centrifugation, or by using a magnetic stand.
3. Transfer the eluate to a new tube.

Document Revision History

Term	Description
06/2023	Updated the Semi-automatic recovery step column of the table in Section 4. Updated the subsection titles of Sections 5 and 6.
08/2022	Added Figure 1 and updated Figure 2. Updated steps 5 and 6 of Section 6.3 (Recovery of samples in wash 1–3 (Steps 3–5)).
11/2021	Initial draft of the <i>EZ2 Connect Fx Recovery Procedure Instruction Manual</i> .

Trademarks: QIAGEN®, Sample to Insight®, QIAamp®, EZ1 & 2™, EZ2®, Investigator®, MinElute® (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

HB-2981-003 06/2023 © 2023 QIAGEN, all rights reserved.

This page intentionally left blank

