EGFR Mutation Analysis – Performance Evaluation of Blood Collection Tubes for ccfDNA Stabilization


PreAnalytiX GmbH, Hombrechtikon, Switzerland


Standardized preanalytical workflows for circulating tumor DNA (ctDNA) analysis from blood are key to obtaining reliable results for cancer research. Some blood collection tubes can decrease assay sensitivity or accuracy by diluting ctDNA with genomic DNA released from white blood cells. Others chemically modify DNA, impairing quality and recovery. This research study compared how different stabilizing blood collection tubes containing either non-crosslinking or formaldehyde-releasing substances perform in a quantitative PCR (qPCR) assay workflow for detecting EGFR mutations.

Blood from healthy consented donors was collected into PAXgene[®] Blood ccfDNA Tubes⁺⁺ (PreAnalytiX[®], n = 20 + 10), Cell-Free DNA BCT[®] (Streck[®], n = 20) and LBgard[®] Blood Tubes (Biomatrica[®], n = 10). Sheared DNA with 500 copies of different EGFR mutations (T790M, L858R and Exon 19 deletion (DEL)) or 500 µl male EDTA plasma was spiked into subsets of samples. Blood was processed either directly after phlebotomy or stored for 7 days at 30°C. Manual extraction of ccfDNA was performed using spin column-based QIAGEN® extraction kits and protocols. Yield and stability of ccfDNA as well as EGFR mutations were analyzed by qPCR. For monitoring changes of 18S rDNA or male specific DYS14 fragments were analyzed in relation to sample replicates processed directly after blood draw. EGFR mutation occurrences were determined by building ΔC_{τ} of the wildtype assay from C_{τ} for the respective mutation. For each mutation, the assay supplier defined a ΔC_{τ} threshold which defines a valid assay result.

High analytical sensitivity in EGFR mutation detection

• Higher analytical sensitivity for EGFR mutation detection after storage with the PAXgene Blood ccfDNA Tube compared to alternative tubes - Stable control C_T values for PAXgene samples while control C_T values in many Streck tubes are lower





DNA containing a low copy number (500 copies) of different EGFR mutations (T790M, L858R and Exon 19 deletion (DEL)) was spiked into PAXgene Blood ccfDNA Tubes, Streck Cell-Free DNA BCT, and Biomatrica LBgard Blood Tubes directly after blood draw. Plasma was processed after blood storage for 7 days at 30°C. EGFR mutation variants were analyzed by qPCR. The result of a wildtype control PCR is part of the assay analysis. High levels of wildtype DNA can lead to negative results although mutant was detected. No signal in the qPCR reaction is indicated with a red symbol.

• The PAXgene Blood ccfDNA Tube outperforms alternative tubes in detection of low EGFR mutant copy numbers in ccfDNA generated from stabilized blood. • Recovery of spiked male ccfDNA in female blood confirmed the performance of the PAXgene Blood ccfDNA Tube.

*For Research Use Only. Not for use in diagnostic procedures.

[†]For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN or PreAnalytiX kit handbook or user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor. Trademarks: PAXgene[®], PreAnalytiX[®] (PreAnalytiX GmbH); QIAGEN[®], QIAsymphony[®] (QIAGEN Group); Streck[®], BCT[®] (Streck, Inc.); Biomatrica[®], LBgard[®] (Biomatrica, Inc.). © 2019 PreAnalytiX GmbH. PreAnalytiX, the PreAnalytiX logo and all other trademarks are property of PreAnalytiX GmbH, Hombrechtikon, CH.

Andrea Ullius, Daniel Groelz, Tomasz Krenz, Eric Provencher and Thorsten Voss

Background


Methods

Results

- Overall stable detection of spiked target molecules in PAXgene samples while samples

Conclusions

Disclaimer



- Storage of blood in Streck tubes resulted in an increased variance of total DNA release into plasma.
- The non-crosslinking stabilization reagent in the PAXgene Blood ccfDNA Tube enables detection of low abundant target molecules and ensures detection sensitivity in downstream reactions.