PartoSure® Test

Instructions For Use

Assess the Risk of Preterm Birth

EN English	3
BG Български	4
CS Český	5
DA Dansk	6
DE Deutsch	7
ES Español	8
FI Finish	9
FR Français	10
HR Hrvatski	11
ID Indonesian	12
IT Italiano	13
KK Kazakh	14
LT Lietuvių	15
LV Latviešu	16
MK Macedonian	17
NL Nederlands	18
NO Norsk	19
PL Polish	20
PT Portuguese (Basil berian)l	21
PT Português	22
RO Română	23
RU Русский	24
SK Slovenčina	25
SR Srpski	26
SV Svenska	27
TR Türkçe	28

QIAGEN
19300 Germantown Road
Germantown, MD 20874 USA
+1-800-426-8157
www.QIAGEN.com

QIAGEN GmbH
QIAGEN Strasse 1
40724 Hilden, GERMANY

TTDT-1-5-OUS
TTDT-1-20-OUS

1108422 Rev. 04

PATENT AND TRADEMARK INFORMATION

©2018 QIAGEN

The PartoSure Test and its use are covered by one or more of the following patents granted or licensed to QIAGEN and/or its subsidiaries; US Patent No. 7,709,272; 8,114,610; 9,891,233, 9,429,580; corresponding foreign patents and other patents pending.

K-Y is a registered trademark of Reckitt Benckiser.

Surgiliube is a registered trademark of Fougera.

Monistat is a registered trademark of Insight Pharmaceuticals.

Betadine is a registered trademark of Purdue Products.

PartoSure and associated logos are trademarks and/or registered trademarks of QIAGEN and/or its subsidiaries in the United States and/or other countries.

GLOSSARY OF SYMBOLS

CODE

EN In Vitro Diagnostic Medical Device
BG Медицинско изделие за in vitro диагностика
CS Průstředek zdravotnické techniky pro diagnostiku in vitro
DA Medicinsk udstyr til in vitro-diagnostik
DE Medizinisches In-vitro-Diagnostikum
ES Dispositivo médico para diagnóstico in vitro
FI Diagnostinen in vitro -lääkintälaite
FR Dispositif médical de diagnostic in vitro
HR Medicinski proizvod za in vitro dijagnostiku
ID Perangkat Medis Diagnostik In Vitro
IT Dispositivo medico-diagnostico in vitro
KK In Vitro диагностикалык медициналык куралдары
LT „In vitro“ diagnostinis medicinos prietaisas
LV In vitro diagnostikas medicīnas ierīce
MK Ин vitro дијагностичко медицинско средство
NL Medisch hulpmiddel voor in-vitrodiagnostiek
NO Medisinsk utstyr til in vitro-diagnostikk
PL Wyrób medyczny do diagnostyki in vitro
PT Dispositivo médico de diagnóstico in vitro
PT Dispositivo médico para diagnóstico in vitro
RO Dispozitiv medical pentru diagnostic in vitro
RU Медицинское устройство для диагностики in vitro
SK Zdravotnícke diagnostické zariadenie na použitie v podmienkach in vitro
SR Medicinsko sredstvo za in vitro dijagnostiku
SV In vitro-diagnostiskt medicinskt produkt
TR In Vitro Tanı Amaçlı Tıbbi Cihaz

Do Not Reuse
Да не се използва повторно
Nepoužívejte opakovaně
Må ikke genbruges
Nicht wiederverwenden
No reutilizar
Ei saa käyttää uudelleen
Ne pas réutiliser
Za jednu upotrebu
Jangan Digunakan Kembali
Non riutilizzare
Қайта пайдалану болмайды
Nenaudoti pakartotinai
Nelietot atkārtoti
He употребуйте повторно
Voor eenmalig gebruik
Ikke til flergangsbruk
Nie używać ponownie
Não reutilizar
Não reutilizar
A nu se refolosi
He использовать повторно
Nepoužívat opakovane
Nemojte da koristite više puta
Får ej återanvändas
Tkrar Kullanmayın

Use By Date
Годно до
Datum použitelnosti
Anvendes inden
Zur Verwendung bis
Fecha de caducidad
Viimeinen käyttöpäivämäärä
À utiliser avant le
Upotrebljivo do
Gunakan Hingga Tanggal
Utilizzare entro la data indicata
Жарамдылық мерзімі
Tinka iki datos
Dergigma terminş
Користете до датум
Uiterste gebruiksdatum
Siste forbruksdato
Termin ważności
Data de validade
Prazo de validade
Termen de valabilitate
Срок годности
Datum použitelnosti
Rok upotrebe
Utgångsdatum
Son Kullanma Tarihi

Caution Consult Accompanying Instructions for Use
Внимание! Вижте придружаващите инструкции за употреба
Uprozornění - přečtěte si příložené pokyny k použití
Advarsel: Læs den medfølgende brugsanvisning
Vorsicht! Beiliegende Gebrauchsanweisungen beachten
Precaución, consulte las instrucciones de uso adjuntas
Varoitus: Lue mukana tulevat käyttöohjeet
Attention, consultez le mode d'emploi joint
Opze: pročítajte prateče upute za upotrebu
Perhatian! Konsultasi Petunjuk yang Menyertai Penggunaan
Attenzione: consultare le istruzioni per l'uso fornite
Назар аударарңыз: бирге берилген пайдалану нускаларын караңыз
Désimos: žr. pridėdamas naudojimo instrukcijas
Uzmanību, iepazīstieties ar pievienoto lietošanas instrukciju
Внимание, разгледайте го придружното упатство за употреба
Let op: raadpleeg de meegeleverde gebruiksaanwijzing
Forsiktig, se bruksanvisningen
Uwaga: zapoznać się z dołączoną instrukcją obsługi
Cuidado, consulte as instruções anexas de uso
Atenção, consultar as Instruções de Utilização que acompanham o produto
Atenție! Consultați Instrucțiunile de utilizare însoțitoare
Внимание! Обратитесь к прилагаемой инструкции по применению
Uprozornenie: Prečítajte si príložený návod na použitie
Opze, pročítajte príložené uputstvo za upotrebu
Försiktig läs den medföljande bruksanvisning
Dikkat! Birlikte Verilen Kullanım Talimatlarına Başvurun

Consult Instructions for Use
Вижте инструкциите за употреба
Viz návod k použití
Læs brugsanvisningen
Gebrauchsanweisung beachten
Consulte las instrucciones de uso
Katsa käyttöohjeet
Consultez le mode d'emploi
Pročítajte upute za upotrebu
Lihat Petunjuk Penggunaan
Consultare le istruzioni per l'uso
Пайдалану нускаларын караңыз
Žr. naudojimo instrukcijas
Skatīt lietošanas norādījumus
Разгледайте го упатството за употреба
Raadpleeg de gebruiksaanwijzing
Se informasjonen i håndboken
Zapoznać się z instrukcją obsługi
Consulte as instruções de uso
Consultar as instruções de utilização
Consultați Instrucțiunile de utilizare
Обратитесь к инструкции по применению
Prečítajte si návod na použitie
Pročítajte uputstvo za upotrebu
Läs bruksanvisningen
Kullanım Talimatlarına Başvurun

Temperature Limitation
Температурни ограничения
Templotni omezení
Temperaturbegrensning
Zulässiger Temperaturbereich
Limitación de temperatura
Lämpötilarajoitus
Limite de temperatura
Ograničenja temperatura
Batasan Suhu
Limite di temperatura
Температуралық шектеуі
Temperatūros apribojimai
Temperatūros ierobežojums
Температурно ограничаване
Zakres temperatury
Temperaturbegrensninger
Limites de temperatura
Limites de temperatura
Limitări de temperatură
Ограничение по температуре
Templotno obmedzenia
Temperaturno ograničenje
Temperaturbegrensning
Sıcaklık Sınırlaması

CODE

Catalog Number
Каталожен номер
Katalogové číslo
Katalognummer
Katalognummer
Número de referencia
Luettelonumero
Numéro de catalogue
Kataloški broj
Nomor Katalog
Numero di catalogo
Katalogo nómír
Katalogo numeris
Kataloga numurs
Каталошки број
Catalogusnummer
Katalognummer
Numer katalogowy
Referência
Número de referência
Număr de catalog
Номер по каталогу
Katalogové číslo
Kataloški broj
Katalognummer
Katalog Numarasi

Manufacturer
Производител
Výrobce
Producent
Hersteller
Fabricante
Valmistaja
Fabriquant
Proizvodač
Pabrik Pembuat
Produttore
Өндіруші
Gamintojas
Ražotājs
Производител
Fabrikant
Produsent
Producent
Fabricante
Fabricante
Producător
Изготвитель
Výrobca
Proizvodač
Tillverkare
Üretici

Sufficient for "X" Tests
Достаточно про тесты "X"
Dostatečně pro testy "X"
Tilstrækkelig til "X" tests
Ausreichend für „X“ Tests
Suficiente para pruebas «X»
X testikertaa
En quantité suffisante pour « X » tests
Dovoljno za „X“ testova
Cukup untuk Uji "X"
Suficiente per "X" test
«X» сынақтары үшін жеткілікті
Galima atlikti „X“ tyrimus
Pietiekams "X" testiem
Довольно за „X“ тестови
Voldoende voor „X“ testen
Tilstrækkelig til «X» tester
Wystarcza do wykonania „X” testów
Suficiente para "X" testes
Suficiente para "X" testes
Suficient pentru „X” teste
Достаточно для выполнения «X» процедур анализа
Postaćujące pre „X” testow
Dovolljno za više testova (x)
Tillräckligt för "X" tester
"X" Testleri için yeterlidir

Batch Code
Код на партида
Číslo šarže
Batchcode
Chargennummer
Código de lote
Eräkoodi
Code de lot
Kód šarže
Kode Batch
Codice lotto
Бума коды
Partijos kodas
Sérijas kods
Број на серија
Batchcode
Partinummer
Kod partii
Código de lote
Código do lote
Cod lot
Код серии
Kód šarže
Šifra serije
Batchkod
Parti Kodu

Global Trade Item Number
Глобален номер на трговска единица
Globální číslo obchodní položky
Globalt handelsvarenummer
Global Trade Item Number
Número mundial de artículo comercial
GTIN-numero
Code article international
Global Trade Item Number (globalni broj trgovačke jedinice)
Nomor Item Perdagangan Global (GTIN)
Global Trade Item Number
Тауардың глобалдық сауда нөмірі
Visuotinis prekės numeris
Globālais tirdzniecības identifikācijas numurs (GTIN)
Меѓународен број за трговска единица
Global Trade Item Number
Globalt artikelnummer
Globalny numer jednostki handlowej
Número Global de Item Comercial (Global Trade Item Number)
Número global de item comercial
Număr de comercializare global articol
Глобалный торговый номер единицы товара
Identifikátor GTIN (Global Trade Item Number)
Globalni broj trgovinske jedinice
GS-artikelnnummer
Küresel Ticaret Parça Numarasi

PartoSure® Test

INSTRUCTIONS FOR *IN VITRO* DIAGNOSTIC USE

INTENDED USE

The PartoSure Test is a rapid, non-instrumented, qualitative immunochromatographic test for the *in vitro* detection of placental alpha microglobulin-1 (PAMG-1) in vaginal secretions of pregnant women. The device is designed as an aid to rapidly assess the risk of preterm delivery in ≤ 7 or ≤ 14 days from the time of cervicovaginal sample collection in pregnant women with signs and symptoms of early preterm labor, intact amniotic membranes, and minimal cervical dilatation (≤ 3 cm), sampled between 20 weeks, 0 days and 36 weeks, 6 days gestation.

SUMMARY AND EXPLANATION OF THE TEST

An accurate risk assessment of preterm birth is clinically important among pregnancies with threatened preterm labor. This is particularly true with respect to both the administration of corticosteroids, which have an optimal benefit within 7 days of administration,¹ as well as the transfer of patients to a tertiary care center capable of caring for the birth of a premature infant. Clinical evaluation alone, including the measurement of cervical length and dilatation, is not sufficiently predictive of imminent delivery.² Similarly, currently available biomarker tests, such as that based on the detection of fetal fibronectin, have extremely poor predictive values for imminent delivery.³ Thus, the accurate prediction of time to delivery in patients with threatened preterm labor remains an important clinical concern. The high negative and positive predictive values of the PartoSure Test may allow for the timely identification of patients that would benefit from available interventions such as antenatal corticosteroids. Additionally, a more accurate identification of these patients will likely reduce the costs associated with unnecessary admissions caused by the poor positive predictive values of existing methods. The test kit is a self-contained system that provides rapid, qualitative results without outside instrumentation or an invasive speculum examination.

PRINCIPLE OF THE TEST

The PartoSure Test is a lateral flow, immunochromatographic assay designed to identify the presence of human placental alpha microglobulin-1 (PAMG-1). The test employs monoclonal antibodies sufficiently sensitive to detect 1 ng/ml of PAMG-1. For the analysis, a sample of cervicovaginal discharge collected by vaginal swab is extracted into a solvent. The presence of PAMG-1 antigen is then detected by inserting a lateral-flow test strip into the vial. The sample flows from an absorbent pad to a nitrocellulose membrane, passing through a reactive area containing monoclonal anti-PAMG-1 antibodies conjugated to a gold particle. The antigen-antibody complex flows to the test region where it is immobilized by a second anti-PAMG-1 antibody. This event leads to the appearance of the test line. Unbound antigen-antibody complexes continue to flow along the test strip and are immobilized by a second antibody. This leads to the appearance of the internal control line.

REAGENTS AND COMPONENTS

The PartoSure Test kit includes the following components: the PartoSure Test strip in a foil pouch with desiccant, a sterile flocked vaginal swab, and a plastic vial with solvent solution (0.9% NaCl, 0.05% NaN₃, 0.01% Triton X100). No additional materials are required to perform the PartoSure Test.

PRECAUTIONS AND WARNINGS

- Use of any swab or solvent solution other than the one provided with the test kit is prohibited.
- Specimens should be collected prior to collection of culture specimens. Collection of vaginal specimens for microbiologic culture frequently requires aggressive collection techniques that may abrade the cervical or vaginal mucosa and may potentially interfere with sample preparation.
- Care must be taken not to contaminate the swab or cervicovaginal secretions with lubricants or antiseptics (e.g. K-Y® or Surgilube® lubricating jelly, Betadine® Cleanser). These substances may interfere with absorption of the specimen by the swab or with the antibody-antigen reaction of the PartoSure Test and lead to invalid test results.
- If it is suspected that the patient has applied a topical disinfectant (e.g. Monistat®, miconazole nitrate cream) to the vaginal area within 24 hours, delay specimen collection until 24 hours from application of the topical disinfectant have passed as these products can lead to false-negative test results.
- The PartoSure Test is not intended for use in women with moderate or gross vaginal bleeding. The presence of vaginal bleeding may contribute to difficulty in interpreting the PartoSure Test result. Testing a moderately to grossly bloody sample may lead to false-positive results. If upon visual examination you are concerned about the presence of moderate or gross vaginal blood, it is recommended that the sample be collected following the cessation of active vaginal bleeding.
- Specimens should not be obtained from patients with suspected or known placental abruption or placenta previa.
- Do not use the PartoSure Test after the expiration date, which is printed on the product packaging.
- PartoSure Test kit components are for single-use only.
- Specimens not tested within 24 hours of collection must be stored refrigerated at 2° to 8°C and tested within five days of collection.
- The PartoSure Test is for *in vitro* diagnostic use only and no component of the test should be taken internally.
- Do not use the kit if the swab or test strip package integrity is compromised or if the solvent vial has leaked.
- Do not bend or fold the test strip or the foil pouch with the test strip in it; doing so may damage the strip and lead to inaccurate results.
- The PartoSure Test kit includes a plastic vial with solvent solution (0.9% NaCl, 0.05% NaN₃, 0.01% Triton X100). Sodium azide may react with plumbing to form potentially explosive metal azides. Avoid contact with skin, eyes, and clothing. In case of contact with any of these reagents, wash area thoroughly with water. If disposing of this reagent, always flush the drain with large volumes of water to prevent azide build-up.
- Safety precautions should be observed when collecting, handling, and disposing of test samples.
- Used PartoSure Test kits are biohazardous; dispose of all components using necessary precautions.

STORAGE AND STABILITY

Store the PartoSure Test kit in a dry place at 15 to 25°C (59 to 77°F). The test should not be frozen. When stored in the foil pouch at the recommended temperature, the test is stable until the expiration date printed on the pouch. The PartoSure Test strip should be used within six (6) hours after removing the test strip from the foil pouch.

QUALITY CONTROL

The PartoSure Test strip contains an internal procedural control mechanism that ensures analytical functionality. The appearance of one or two lines in the results region of the test strip verifies the integrity of the test procedure and components.

TEST PROCEDURE

1. Take the solvent vial by its cap and ensure all liquid in the vial has dropped to the bottom. Open the solvent vial and place it in a vertical position.
2. To collect a sample from the vagina, use only the sterile flocked swab provided with the PartoSure Test kit. Remove the swab from its package following the instructions on the packaging. The tip of the swab should not touch anything prior to insertion into the vagina. Hold the swab by the middle of its shaft and, while the patient is lying on her back, carefully insert the tip of the swab into the vagina until the fingers contact the skin (no more than 5-7 cm deep). Withdraw the swab from the vagina after 30 seconds.
3. After the swab has been removed from the vagina, immediately place the tip into the provided solvent vial and rinse by rotating for 30 seconds.
4. Remove the swab from the vial and dispose of it.
5. Tear open the foil pouch at the tear notches and remove the PartoSure Test strip.
6. Insert the white end of the test strip (marked with arrows facing downward) into the vial with solvent.
7. Remove the test strip from the vial if two lines are clearly visible in the test region or after 5 minutes sharp. Read the results by placing the test strip on a clean, dry, and flat surface in a well-lit environment via either natural or fluorescent lighting. A positive result is indicated by two lines in the test region, while a negative result is indicated by a single line in the test region. Do not read or interpret the results after 10 minutes have passed since inserting the test strip into the vial.

Two lines: POSITIVE	One line: NEGATIVE	No lines: INVALID
Control Line ↓ Test Line	Control Line ↓	Test Line Only
Positive: Two Lines	Negative: One Control Line	Invalid: No Lines or Test Line Only
Imminent delivery within ≤ 7 or ≤ 14 days is highly likely	Imminent delivery within ≤ 7 or ≤ 14 days is highly unlikely	Results not valid; retest

The intensity of the lines may vary; the test result is valid even if the lines are faint or uneven. Do not interpret the test result based on the intensity of the lines.

LIMITATIONS OF THE TEST

- The PartoSure Test result should not be interpreted as absolute evidence for the presence or absence of a process that will result in delivery ≤ 7 or ≤ 14 days from specimen collection.
- The PartoSure Test result should always be used in conjunction with information available from the clinical evaluation of the patient and other diagnostic procedures such as cervical examination, assessment of uterine activity, and evaluation of other risk factors.
- These instructions for use must be followed exactly; failure to do so may lead to inaccurate results.
- PartoSure Test performance has been characterized from specimens taken from the vaginal cavity. Samples obtained from other locations should not be used. A speculum examination is not required.
- Results should be interpreted with caution when a specimen is obtained from a patient with unconfirmed gestational age.
- PartoSure Test results are qualitative and not quantitative. No quantitative interpretation should be made based on the strength of the test or control lines.
- The PartoSure Test should only be used in patients with signs and symptoms of preterm labor.
- Invalid test results may occur if lubricants or antiseptics (e.g. K-Y® or Surgilube® lubricating jelly, Betadine® Cleanser) have been used by the patient.

EXPECTED VALUES

PartoSure detects trace amounts of human PAMG-1. PAMG-1 concentrations greater than 1 ng/ml indicate an elevated risk of delivering within the next 7 days. The cutoff of 1 ng/ml was validated in an international multicenter clinical trial evaluating the association of PAMG-1 present in cervicovaginal secretions and the test-to-delivery time interval.

PERFORMANCE CHARACTERISTICS

Precision and reproducibility were determined using three lots of the PartoSure Test at three intended-use sites by three different intended users at each site. Five replicates of seven different PAMG-1 concentrations above and below the limit of detection of the PartoSure Test were used, including an absolute zero (0.0 ng/ml), a low-negative (0.2 ng/ml), a high-negative (0.5 ng/ml), the limit of detection (1.0 ng/ml), a low-positive (2.0 ng/ml), and two concentrations within the C₅-C₉ interval for the device (0.7 ng/ml and 0.9 ng/ml). The results of this study demonstrated 100% negative readings at ≤ 0.5 ng/ml and 100% positive readings at ≥ 1 ng/ml PAMG-1 levels by intended users.

The clinical performance of the PartoSure Test was determined in an international multi-site study. Patients who (i) reported signs, symptoms, or complaints suggestive of preterm labor, (ii) were between 20 and 36¹/₂ weeks of gestation, and (iii) had clinically intact membranes were invited to participate in the trial. In total, 203 patients were included in the final analysis. The following performance estimates with corresponding 95% confidence intervals were calculated:

Performance Measure	Delivery Within	
	≤ 7 Days	≤ 14 Days
Sensitivity (SN) (95% CI)*	80% (63.1% – 91.6%)	63% (47.4% – 76.1%)
Specificity (SP) (95% CI)*	95% (90.1% – 97.5%)	96% (90.9% – 98.2%)
Negative Predictive Value (NPV) (95% CI)*	96% (91.5% – 98.3%)	89% (83.4% – 93.5%)
Positive Predictive Value (PPV) (95% CI)*	76% (58.8% – 88.2%)	81% (64.8% – 92.0%)

*95% confidence intervals (CI) computed via the Clopper-Pearson procedure
Nikolova T, Bayev O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. *J Perinat Med.* 2015 Jul;43(4):395-402.

CROSS REACTIVITY

The PartoSure Test was evaluated using a panel of potentially cross-reactive protein substances likely to be found in vaginal specimens, including human chorionic gonadotropin, trophoblastic beta-2 glycoprotein, human placental lactogen, alpha-fetoprotein, IGFBP-3, and human serum albumin. Four (4) replicates of each sample containing the potentially cross-reactive substance were tested using a high-negative PAMG-1 sample (0.2 ng/ml) and a low-positive sample (2.0 ng/ml). Each potentially cross-reactive substance was tested at the highest concentration of substance that was considered clinically relevant. None of the potentially cross-reactive substances tested demonstrated cross-reactivity with the PartoSure Test.

INTERFERENCE STUDIES

Lubricants, Disinfectants, Antiseptics, Soaps & Creams
Four (4) replicates each of a high-negative sample (0.2 ng/ml) and a low-positive sample (2.0 ng/ml) were tested against lubricant jelly, miconazole nitrate disinfectant cream, antiseptic cleanser, body wash soap, and vaginal cream. While all high-negative replicates were correctly identified in the presence of miconazole nitrate, low-positive replicates were not. The presence of lubricating jelly led to several invalid results for both high-negative and low-positive replicates, though in cases where a negative or positive result was reported, the result was correct. The presence of antiseptic cleanser led to invalid test results in all replicates. Therefore, care must be taken not to contaminate the swab or cervicovaginal secretions with lubricants or antiseptics (e.g. K-Y® or Surgilube® lubricating jelly, Betadine® Cleanser). If it is suspected that the patient has applied a topical disinfectant (e.g. Monistat®, miconazole nitrate cream) to the vaginal area within 24 hours, delay specimen collection until 24 hours from application have passed as these products can lead to false-negative test results. The presence of soap or cream did not interfere with the PartoSure Test.

Pharmacological Agents

Four (4) replicates each of a high-negative sample (0.2 ng/ml) and a low-positive sample (2.0 ng/ml) were tested against 10 non-biologically based interfering substances including: 17-OH progesterone (50 µg/ml), ampicillin (152 µmol/L), cephalixin (337 µmol/L), erythromycin (81.6 µmol/L), gentamycin (21 µmol/L), dexamethasone (1.53 µg/ml), magnesium sulfate (50 µg/ml), oxytocin (58 µU/ml), terbutaline (1 mg/ml), and ritodrine (100 µg/ml). None of these non-biological substances interfered with the PartoSure Test.

Vaginal Bacterial Pathogens

Four (4) replicates each of a high-negative sample (0.2 ng/ml) and a low-positive sample (2.0 ng/ml) were tested against 3 vaginal bacterial pathogens, including Gardnerella vaginalis, Candida albicans, and Trichomonas vaginalis. None of these vaginal infection pathogens interfered with the PartoSure Test.

Maternal Bleeding

Four (4) replicates each of a high-negative sample (0.2 ng/ml) and a low-positive sample (2.0 ng/ml) were tested against ten (10) individual maternal blood samples at the three (3) lowest admixture levels determined to represent "trace," "moderate," or "gross" levels of maternal bleeding on the vaginal collection swab. Moderate or gross vaginal bleeding may therefore contribute to difficulty in interpreting the PartoSure Test result and may lead to false-positive results.

Semen & Urine

A high-negative (0.2 ng/ml) sample and a low-positive sample (2.0 ng/ml) were tested against 10 individual samples of semen and 10 individual samples of maternal urine. Neither semen nor maternal urine demonstrated interference with the PartoSure Test.

BIBLIOGRAPHY

1. Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins—Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. *Obstet Gynecol* 2012; 119:1308-1316.
3. Berghella V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Art. No.: CD006843. [Systematic review]

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1. Antenatal Corticosteroids Revised: Repeat Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. Obstet Gynecol 2012; 119:1308-1316.
3. Berghella V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. Cochrane Database of Systematic Reviews 2008, Issue 4. Art. No.: CD006843. [Systematic review]

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
*95% intervaly spolehlivosti (CI) byly vypočteny pomocí postupu Clopper-Pearson Nikolova T, Bayev O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. *J Perinat Med.* 2015 Jul;43(4):395-402.

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LITERATURA

- Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
- ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. *Obstet Gynecol* 2012; 119:1308-1316.
- Borghello V, Hayes E, Visintin J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Art. No.: CD006843. (Systematic review)

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Nikolova T, Bayev O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. *J Perinat Med.* Jul. 2015;43(4):395-402.

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LITTERATURLISTE

1. Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin nr. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. *Obstet Gynecol* 2012; 119:1308-1316.
3. Berghella V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, 4. udgave. Art. nr.: CD006843. (Systematic review)

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Nikolova T, Baye O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. J Perinat Med. 2015 Jul;43(4):395-402.

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LITERATUR

- Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
- ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. Obstet Gynecol 2012; 119:1308-1316.
- Borghello V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. Cochrane Database of Systematic Reviews 2008, Issue 4. Art. No.: CD006843. (Systematic review)

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*Intervalos de confianza (IC) del 95 % computados a través del procedimiento Clopper-Pearson
Nikolova T, Bayev O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. *J Perinat Med.* 2015 Jul;43(4):395-402.

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1. Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. *Obstet Gynecol* 2012; 119:1308-1316.
3. Berghello V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Art. No.: CD006843. (Systematic review)

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LÄHTEET

1. Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. Obstet Gynecol 2012; 119:1308-1316.
3. Berghella V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. Cochrane Database of Systematic Reviews 2008, Issue 4. Art. No.: CD006843. (Systematic review)

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*Intervalles de confiance (IC) 95 % calculés selon la procédure de Clopper-Pearson Nikolova T, Bayev O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. *J Perinat Med.* 2015 Jul;43(4):395-402.

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BIBLIOGRAPHIE

- Antenatal Corticosteroids Revisited: Repeat Courses. *NIH Consensus Statement*: NIH, 2000.
- ACOG Committee on Practice Bulletins-Obstetrics. *ACOG Practice Bulletin* no. 127: Management of Preterm Labor. *Clinical management guidelines for obstetrician-gynecologists*. *Obstet Gynecol* 2012; 119:1308-1316.
- Berghele V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Art. No.: CD006843. (Systematic review)

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BIBLIOGRAFIJA

- Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
- ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. Obstet Gynecol 2012; 119:1308-1316.
- Berghele V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. Cochrane Database of Systematic Reviews 2008, Issue 4. Art. No.: CD006843. (Systematic review)

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*95% interval kepercayaan (CI/Confidence Intervals) terkomputasi melalui prosedur Clopper-Pearson
Nikolova T, Baye O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. *J Perinat Med.* 2015 Jul;43(4):395-402.

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*Intervalli di confidenza (IC) del 95% calcolati con la procedura di Clopper-Pearson. Nikolova T, Bayev O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. *J Perinat Med.* 2015; *44*(4):395-402.

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BIBLIOGRAFIA

- Antenatal Corticosteroids Revisited: Repeat Courses. *NIH Consensus Statement: NIH*, 2000.
- ACOG Committee on Practice Bulletins-Obstetrics. *ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists.* *Obstet Gynecol* 2012; *119*:1308-1316.
- Borghello V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Art. No.: CD006843. (Systematic review)

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- Antenatal Corticosteroids Revisited: Repeat Courses. *NH Consensus Statement.* *NH.* 2000.
- ACOG Committee on Practice Bulletins-Obstetrics. *ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists.* *Obstet Gynecol.* 2012; 119:1308-1316.
- Berghella V, Hayes E, Vasiline J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Art. no.: CD006843. (Systematic review)

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- Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
- ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. Obstet Gynecol 2012; 119:1308-1316.
- Borghello V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. Cochrane Database of Systematic Reviews 2008, Issue 4. Art. No.: CD006843. [Systematic review]

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* betrouwbaarheidsinterval (CI) van 95% berekend met de Clopper-Pearson-methode Nikolova T, Bayev O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. J Perinat Med. 2015 Jul;43(4):395-402.

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1. Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins—Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. Obstet Gynecol 2012; 119:1308-1316.
3. Berghello V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. Cochrane Database of Systematic Reviews 2008, Issue 4. Art.Nr.: CD006843. (Systematic review)

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
^a95 % konfidensintervaller (CI) beregnet via Clopper-Pearson-prosedyren Nikolova T, Bayev O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. *J Perinat Med.* Juli 2015;43(4):395-402.

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BIBLIOGRAFI

1. Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. *Obstet Gynecol* 2012; 119:1308-1316.
3. Berghella V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Art. No.: CD006843. [Systematic review]

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BIBLIOGRAFIA

- Antenatal Corticosteroids Revisited: Repeat Courses. *NH Consensus Statement: NIH*, 2000.
- ACOG Committee on Practice Bulletins-Obstetrics. *ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists*. *Obstet Gynecol* 2012; 119:1308-1316.
- Berghello V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Art. No.: CD006843. (Systematic review)

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* Intervalos de confiança de 95% (CI) foram calculados através do procedimento de Clopper-Pearson
Nikolova T, Bayev O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. J Perinat Med. 2015 Jul;43(4):395-402.

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Nikolova T, Bayev O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. *J Perinat Med.* 2015 Jul;43(4):395-402.

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BIBLIOGRAFIA

- Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
- ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. *Obstet Gynecol* 2012; 119:1308-1316.
- Borghella V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Art. No.: CD006843. (Systematic review)

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*Interval de încredere (I^a) 95% calculat prin procedura Clapper-Pearson Nikolova T, Bayev O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. J Perinat Med. 2015 Jul;43(4):395-402.

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BIBLIOGRAFIE

1. Antenatal Corticosteroids Revised: Repeat Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. Obstet Gynecol 2012; 119:1308-1316.
3. Berghele V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. Cochrane Database of Systematic Reviews 2008, Issue 4. Art. No.: CD006843. [Systematic review]

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1. Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins—Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. *Obstet Gynecol* 2012; 119:1308–1316.
3. Berghella V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Art. No.: CD006843. [Systematic review]

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Nikolova T, Baye O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. *J Perinat Med.* 2015 Jul;43(4):395-402.

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1. Antenatal Corticosteroids Revised: Repeat Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. *Obstet Gynecol* 2012; 119:1308-1316.
3. Bergella V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, 4. vydanie. <http://dx.doi.org/10.1002/CD006843>. [Systematic review]

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

Nikolova P, Bayev O, Nikolova N, Di Renzo GC. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. *J Perinat Med.* 2015 Jul;43(4):395-402.

UNAKRSNA REAKTIVNOST

Procena PartoSure testa je izvršena pomoću panela potencijalno unakrsno reaktivnih proteina koji su česti u vaginalnim uzorcima, uključujući humani horionski gonadotropin, beta-2 glikoprotein trofoblasta, humani placentalni laktogen, alfafetoprotein, IGFBP-3 i humani serumski albumin. Četiri (4) replikata svakog uzorka sa potencijalno unakrsno reaktivnom supstancom su testirana pomoću uzorka PAMG-1 visoke negativne vrednosti (0,2 ng/ml) i uzorka niske pozitivne vrednosti (2,0 ng/ml). Svaka potencijalno unakrsno reaktivna supstanca je testirana pri najvećoj koncentraciji supstance koja je ocenjena kao klinički relevantna. Nijedna testirana potencijalno unakrsno reaktivna supstanca nije iskazala unakrsnu reaktivnost sa PartoSure testom.

STUDIJE SMETNJI

Lubrikanti, sredstva za dezinfekciju, antiseptici, sapuni i kreme Četiri (4) replikata uzorka visoke negativne vrednosti (0,2 ng/ml) i uzorka niske pozitivne vrednosti (2,0 ng/ml) je testirano na lubrikant, krem sa mikonazol-nitratom, antiseptik, sapun za telo i vaginalnu kremu. Dok je prisustvo svih replikata visoke negativne vrednosti utvrđeno pri korišćenju mikonazol-nitrata, to se nije desilo u slučaju replikata niske negativne vrednosti. Prisustvo lubrikanta proizvelo je nekoliko nevažnih rezultata za replikate sa visokom negativnom i niskom pozitivnom vrednošću, ali je rezultat bio tačan u slučaju prijavljivanja negativnog ili pozitivnog rezultata. Prisustvo antiseptika je proizvelo nevažne rezultate testa u svim replikatima. Dakle, vodite računa da ne kontaminirate štapić za uzorkovanje ili cervikovaginalni iscedak lubrikantima ili antisepticima (npr. lubrikant K-Y® ili Surgilube®, antiseptik Betadine®). Ako sumnjate da je pacijentkinja nanela lokalno sredstvo za dezinfekciju (npr. Monistat®, krem sa mikonazol-nitratom) u područje vagine tokom prethodna 24 sata, odložite prikupljanje uzorka 24 sata nakon primene jer ovi proizvodi mogu da proizvedu lažno negativne rezultate testa. Prisustvo sapuna ili kreme nije uticalo na PartoSure test.

Farmakološka sredstva

Četiri (4) replikata uzorka visoke negativne vrednosti (0,2 ng/ml) i uzorka niske pozitivne vrednosti (2,0 ng/ml) je testirano na 10 nebioloških interferentnih supstanci, uključujući: 17-OH progesteron (50 µg/ml), ampicilin (152 µmol/l), cefaleksin (337 µg/ml), eritromicin (81,6 µmol/l), gentamicin (21 µmol/l), deksametazon (1,53 µg/ml), magnezijum sulfat (50 µg/ml), oksitocin (58 µg/ml), terbutalin (1 mg/ml) i ritodrin (100 µg/ml). Nijedna od ovih nebioloških supstanci nije uticala na test.

Vaginalni bakterijski patogeni

Četiri (4) replikata uzorka visoke negativne vrednosti (0,2 ng/ml) i uzorka niske pozitivne vrednosti (2,0 ng/ml) je testirano na 3 vaginalna bakterijska patogena, uključujući Gardnerella vaginalis, Candida albicans i Trichomonas vaginalis. Nijedna od ovih vaginalnih infekcijskih patogena nije uticala na test.

Krvarenje trudnice

Četiri (4) replikata uzorka visoke negativne vrednosti (0,2 ng/ml) i uzorka niske pozitivne vrednosti (2,0 ng/ml) je testirano na deset (10) pojedinačnih uzorka krvi trudnice u tri (3) najniža nivoa mešavine koji su predstavljali „mestimično“, „umereno“ ili „obilno“ krvarenje trudnice na štapiću za uzorkovanje iz vagine. Umereno ili obilno vaginalno krvarenje može dovesti do poteškoća pri tumačenju rezultata PartoSure testa i posledično do lažno pozitivnih rezultata.

Semena tečnost i urin

Uzorak visoke negativne vrednosti (0,2 ng/ml) i uzorak niske pozitivne vrednosti (2,0 ng/ml) su testirani sa 10 pojedinačnih uzorka semene tečnosti i 10 pojedinačnih uzorka urina. Ni semena tečnost ni urin nisu uticali na rezultate PartoSure testa.

LITERATURA

- Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
- ACOG Committee on Practice Bulletins-Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. *Obstet Gynecol* 2012; 119:1308-1316.
- Borghello V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Art. No.: CD006843. (Systematic review)

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