



BioMag[®] Goat Anti-Human IgG Fc

BioMag[®] suspensions are a series of products containing superparamagnetic particles. BioMag particles are irregularly shaped, with an average diameter of 1.6 μm , and consist of an iron oxide core with a silane coating. The surface of the particle is coated with amine or carboxyl groups, facilitating the covalent attachment of proteins, glycoproteins, secondary antibodies, and other ligands, with retention of biological activity. The particles' irregular shape provides a large surface area, ensuring a high binding efficiency.

BioMag particles are superparamagnetic — they respond to magnetic fields but do not retain magnetic properties upon removal of the magnetic field. This inability to become magnetized permits magnetic extraction without magnetically induced aggregation. Rapid and efficient removal of BioMag particles from suspension is achieved by the application of an external magnetic field.

Specifications

Form	Magnetic particle suspension (5 mg/ml) in phosphate-buffered saline containing EDTA and 0.08% sodium azide.
BioMag concentration	5 mg/ml
Binding capacity	1 ml (5 mg) BioMag Goat Anti-Human IgG Fc binds >0.10 mg human IgG. The antibody will react with heavy chains on human IgG but not with any light chains.
Storage and stability	BioMag Goat Anti-Human IgG Fc can be stored at 2–8°C until the expiration date. Do not freeze. Do not dry. Centrifugation should only be used if it is the last step of a procedure, i.e., if resuspension of the BioMag particles is not required. Freezing, drying, and centrifugation result in extensive aggregation of the BioMag particles and a loss of binding activity.
Safety information	When working with chemicals always wear a suitable lab coat, disposable gloves, and protective goggles. For more information please consult the appropriate material safety data sheet, available online at www.qiagen.com/ts/msds.asp

Applications

BioMag Goat Anti-Human IgG Fc is a standard BioMag particle coated with polyclonal goat anti-human IgG Fc antibodies and is highly suited for use in cell sorting methods where a human IgG antibody is used as a primary antibody. BioMag Goat Anti-Human IgG Fc can be used to separate the cells of interest from a heterogeneous cell population using negative depletion. BioMag Goat Anti-Human IgG Fc can also be used as a secondary antibody in enzyme immunoassays and radioassays that utilize a human IgG primary monoclonal antibody.

Ordering Information

Product	Contents	Cat. No.
BioMag Goat Anti-Human IgG Fc (50 ml)	BioMag goat anti-human IgG (Fc-region specific) secondary antibody suspension (5 mg/ml)	310344

Accessories

Single Tube Magnet	Magnet for separating magnetic particles in a 1.5 ml or 2 ml tube	36910
12-Tube Magnet	Magnet for separating magnetic particles in 12 x 1.5 ml or 2 ml tubes	36912
15 ml/50 ml Tube Magnet	Magnet for separating magnetic particles in 5 x 15 ml and 3 x 50 ml tubes	36935
Flask Magnet	Magnet for separating magnetic particles in a cell culture flask	36937

Patented or patent-pending and/or registered or registration-pending trademark of the QIAGEN Group of Companies: QIAGEN®, BioMag®.
Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

© 2002 QIAGEN, all rights reserved.

www.qiagen.com

QIAGEN:

Australia Tel. 03-9489-3666 Fax 03-9489-3888	Canada Tel. 800-572-9613 Fax 800-713-5951	France Tel. 01-60-920-930 Fax 01-60-920-925
Germany Tel. 02103-29-12400 Fax 02103-29-22000	Italy Tel. 02-33430411 Fax 02-33430426	Japan Tel. 03-5547-0811 Fax 03-5547-0818
Switzerland Tel. 061-319-30-31 Fax 061-319-30-33	UK and Ireland Tel. 01293-422-999 Fax 01293-422-922	USA Tel. 800-426-8157 Fax 800-718-2056

Distributors:

Argentina Tecnolab S.A. (011) 4555 0010 Austria/Slovenia Merck Eurolab GmbH (01) 576 00 0 Belgium/Luxembourg Westburg b.v. 0800-1-9815 Brazil Uniscience do Brasil 011 3622 2320 China Gene Company Limited (852)2896-6283 Cyprus Scientronics Ltd (02) 765 416 Czech Republic BIO-CONSULT spol. s.r.o. (420) 2 417 29 792 Denmark Merck Eurolab A/S 43 86 87 88 Egypt Clinilab 52 57 212 Finland Merck Eurolab Oy (09)804 551 Greece BioAnalytica S.A. (01)-640 03 18 India Genetix (011)-542 1714 or (011)-515 9346 Israel Westburg (Israel) Ltd. 08 6650813/4 or 1-800 20 22 20 Korea IRS Laboratories, Inc. (02) 924-86 97 Malaysia RESEARCH BIOLABS SDN. BHD. (603)-8070 3101 Mexico Quimica Volante S.A. de C.V. (5) 525 57 25 The Netherlands Westburg b.v. (033)-4950094 New Zealand Biolab Scientific Ltd. (09) 980 6700 or 0800 933 966 Norway Merck Eurolab AS 22 90 00 00 Poland Syngene Biotech Sp.z o.o. (071) 351 41 06 or 0601 70 60 07 Portugal IZASA PORTUGAL, LDA (21) 424 7312 Singapore Research Biolabs Pte Ltd 2731066 Slovak Republic BIO-CONSULT Slovakia spol. s.r.o. (02) 5022 1336 South Africa Southern Cross Biotechnology (Pty) Ltd (021) 671 5166 Spain IZASA, S.A. (93) 902.20.30.90 Sweden Merck Eurolab AB (08) 621 34 00 Taiwan TAIGEN Bioscience Corporation (02) 2880 2913 Thailand Theera Trading Co. Ltd. (02) 412-5672 In other countries contact: QIAGEN, Germany



1019288 01/2002