

High-sensitivity screening of a large number of samples for *KRAS* and *PIK3CA* mutations using digital PCR



James Qin
QIAGEN Sciences LLC, 19300 Germantown Rd. Germantown, MD 20874, USA

Introduction and methods

Colorectal cancer (CRC) is the third most common cancer in both women and men, with ~1.9 million new cases and >900,000 deaths worldwide in 2022 (1)

KRAS and PIK3CA in CRC

Activating mutations in the proto-oncogene *KRAS* are among the most frequent in CRC. Of these, the *KRAS* G12 gene has 15 different point mutations (including G12S), and accounts for 6.5% of CRC cases (2). Somatic mutations in one of the PI3K catalytic subunit genes, *PIK3CA*, are also frequently found in numerous malignancies, including CRC.

There is a strong association between *KRAS* and *PIK3CA* mutations, with simultaneous mutations in *PIK3CA* exon 9 and *KRAS* exon 2 the most common (67.2%) (3). Compared with WT tumors, tumors with concomitant *KRAS* and *PIK3CA* mutations were significantly associated with aggressive clinically pathological features. For example, increased proportions of proximal colon cancer or prone to metastasize to distal organs at initial surgery (4, 5).

Methods


Here, we present a novel dPCR methodology and instrumentation to detect rare mutations, such as *PIK3CA* and *KRAS*, with consistent sensitivity, robustness and precision.

We selected two dPCR LNA Mutation Assays (QIAGEN® GeneGlobe®, cat. no. 250200; wet-lab validated) and germline cell-model systems (Horizon Discovery):


- *KRAS* G12S reference gDNA from SW48 cell line (GGT/AGT), with 50% mutation (cat. no. HD288)
- *PIK3CA* H1047R reference gDNA from PI3Ka(H1047/+) cell line, (CAT/CGT) with 50% mutation (cat. no. HD690)

These two model systems were used to demonstrate assay sensitivity and explore the minimum range of total input amount on the QIAcuity Digital PCR System (QIAcuity Four).


Results: Robust PIK3CA dPCR assay design and performance




Results: Highly sensitive PIK3CA mutation detection using dPCR




Results: High PIK3CA dPCR assay sensitivity




Results: Low sample input for PIK3CA mutation detection with dPCR



Results: Highly sensitive KRAS mutation detection using dPCR



Results: Low sample input for KRAS mutation detection with dPCR



Conclusions

- The QIAcuity Digital PCR System is highly consistent and sensitive, requiring no standards. The system can detect as low as 0.05–0.1% of mutation level reliably and consistently from 100–10 ng of total sample input volumes from two germline cell models:
 - *KRAS* (G12S/+) SW48 cell line (Horizon Discovery, cat. no. HD288)
 - *PIK3CA* (H1047/+) cell line (Horizon Discovery, cat. no. HD690).
- We have demonstrated the QIAcuity Four instrument and related kits are capable of mutation detection to the sub-0.1% level, making it highly suited for use in clinical cancer research with samples that are often scarce and/or poor quality (e.g., FFPE tumor samples and cell-free DNA)
- QIAcuity dPCR is easy to adapt from quantitative PCR assays, making the migration of laboratory-developed or commercial qPCR to dPCR assays robust and seamless.
- Mutational LNA probe assays from QIAGEN GeneGlobe database are wet-lab validated. This study demonstrates that assay performance is robust, sensitive and applicable for the detection of cancer targets in multiplexing.
- The sensitivity, robustness and low sample requirement of dPCR assays may be well suited for translational research topics, such as disease progression monitoring.

References
1. WHO International Agency for Research on Cancer. <https://www.iarc.who.int/cancer-type/colorectal-cancer/>. Accessed: March 17, 2025. **2.** Meng M, et al. *Biomed Pharmacother.* 2021;140:111717. **3.** Li H-T, et al. *Oncol Rep.* 2011;25(6):1691-7. **4.** Kosty C, et al. *PLoS One.* 2013;8(6):e65479. **5.** Luo Q, et al. *Transl Oncol.* 2020;13(12):10087.
The QIAcuity is intended for molecular biology applications. This product is not intended for the diagnosis, prevention or treatment of a disease. For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit instructions for use or user operator manual. QIAGEN instructions for use and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services (or your local distributor).
Trademarks: QIAGEN®, Sample to Insight®, QIAcuity®, GeneGlobe® (QIAGEN Group); FAM®, HEX® (Thermo Fisher Scientific Inc.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, may still be protected by law. QRC09921 © 2025 QIAGEN, all rights reserved.