AmniSure® ROM (Rupture of [fetal] Membranes) Test

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참고 문헌
1. Cousins LM et al. AmniSure Placental Alpha Microglobulin-1 Rapid Immunoassay versus Standard Diagnostic Methods for Detection of Rupture of Membranes. Am J Perinatol. 2005; (22):317-320. 2. Petrunin, D. Immunochemical identification of organ specific human placental alpha-globulin and its concentration in amniotic fluid. Akush Ginekolo (Mosk) 1977; (1):64-5. 3. Chen, F, Dudenhausen, J. Comparison of Two Rapid Strip Tests Based on IGFBP-1 and PAMG-1 for the Detection of Amniotic Fluid. Am J Perinatol 2008; (25):243-246. 4. Gaucherand P, Doret M., et al. Detection of Placental Alpha Microglobulin-1 versus Insulin-Like Growth Factor-Binding Protein-1 in Amniotic Fluid at Term: A Comparative Study. Am J Perinatol 2011; (6):489-94. 5. Tagore S, Kwek K. Comparative analysis of insulin-like growth factor binding protein-1 (IGFBP-1), placental alpha microglobulin-1 (PAMG-1) and nitrazine test to diagnose premature rupture of membranes in pregnancy. Am J Perinatol 2008; (25):243-6. 6. Albayrak M, Ozdemir I, Koc O, Ankarali H, & Oren O. Comparison of the diagnostic efficacy of the two rapid bedside immunoassays and combined clinical conventional diagnosis in prelabor rupture of membranes. Eur J Obstet Gynecol Reprod Biol. 2011 Oct; 158(2):179-82. 7. Di Renzo GC, Cabero Roura L, Facchinetti F & the EAPM-Study Group on "Preterm Birth". Guidelines for the management of spontaneous preterm labor: identification of spontaneous preterm labor, diagnosis of preterm premature rupture of membranes, and preventive tools for preterm birth. The Journal of Maternal-Fetal & Neonatal Medicine 2011; 24(5):659-667.

상표
K-Y® Jelly(소유주: Reckitt Benckiser, Slough, England)

Monistat® Yeast Infection Treatment(소유주: Insight Pharmaceuticals, White Plains, New York)

Replens® Feminine Moisturizer(소유주: Church and Dwight Co., Inc., Ewing, New Jersey)

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© 2015 QIAGEN
QIAGEN GmbH