AmniSure® ROM (ruptura fetalnih membrana) Test

UPUTSTVO ZA UPOTREBU U *IN VITRO* DIJAGNOSTICI

NAMENA

AmniSure ROM Test je brz, kvalitativni imunohromatografski test za *in vitro* otkrivanje amnijske tečnosti u vaginalnom isjetku trudnica. Test je namenjen zdravstvenim radnicima kako bi se pomoglo u otkrivanju rupture fetalnih membrana kod trudnica.¹

SAŽETAK I OBJAŠNJENJE TESTA

Preвременa ruptura fetalnih membrana (engl. PROM: Premature rupture of fetal membranes) se događa u oko 10% trudnoća i predstavlja jednu od najvažnijih terapijskih dilema iz svakodnevne ginekološke prakse.¹ Lečenje trudnica sa PPO ili prevremenom rupturom fetalnih membrana pre 37. nedelje trudnoće (engl. PPRO: pre-term PROM) je skupo, i predstavlja važnu perinatalnu dilemu kliničkog lekara koji mora pokušati da uravnoteži rizike produženja trudnoće i razvoja infekcije.¹

Rizici PROM-a u terminu povezani su sa ozbiljnim neonatalnim posledicama kao što su prevremeni porodaj, fetalni distres, prolaps vrpce, abruptio placentae i infekcija.¹ PPROM čini 20% do 40% slučajeva PROM-a, a povezan je sa 20% do 50% prevremenih porođaja, infektivnim morbiditetom kod majke i fetusa, plućnom hipoplazijom fetusa, prolapsom pupčane vrpce, razvojem fetalnih deformiteta i postnatalnim endometritisom.¹ Sve ove posledice značajno povećavaju morbiditet i smrtnost fetusa i majki.

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



LITERATURA

1. Cousins LM et al. AmniSure Placental Alpha Microglobulin-1 Rapid Immunoassay versus Standard Diagnostic Methods for Detection of Rupture of Membranes. Am J Perinatol. 2005; (22):317-320. 2. Petrunin, D. Immunochemical identification of organ specific human placental alpha-globulin and its concentration in amniotic fluid. Akush Ginekol (Mosk) 1977; (1):64-5.

3. Chen, F, Dudenhausen, J. Comparison of Two Rapid Strip Tests Based on IGFBP-1 and PAMG-1 for the Detection of Amniotic Fluid. Am J Perinatol 2008; (25):243-246. 4. Gaucherand P., Doret M., et al. Detection of Placental Alpha Microglobulin-1 versus Insulin-Like Growth Factor-Binding Protein-1 in Amniotic Fluid at Term: A Comparative Study. Am J Perinatol 2011; (6):489-94. 5. Tagore S, Kwek K. Comparative analysis of insulin-like growth factor binding protein-1 (IGFBP-1), placental alpha microglobulin-1 (PAMG-1) and nitrazine test to diagnose premature rupture of membranes in pregnancy. Am J Perinatol 2008; (25):243-6. 6. Albayrak M, Ozdemir I, Koc O, Ankarali H, & Oren O. Comparison of the diagnostic efficacy of the two rapid bedside immunoassays and combined clinical conventional diagnosis in prelabour rupture of membranes. Eur J Obstet Gynecol Reprod Biol. 2011 Oct; 158(2):179-82. 7. Di Renzo GC, Cabero Roura L, Facchinetti F & the EAPM-Study Group on "Preterm Birth". Guidelines for the management of spontaneous preterm labor: identification of spontaneous preterm labor, diagnosis of preterm premature rupture of membranes, and preventive tools for preterm birth. The Journal of Maternal-Fetal & Neonatal Medicine 2011; 24(5):659-667.

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 **QIAGEN**
19300 Germantown Road
Germantown, MD 20874, SAD

 **QIAGEN GmbH**
QIAGEN Strasse 1
40724 Hilden, NEMAČKA

 **FMRT-1-10-ML-RT**
FMRT-1-25-ML-RT

 **1135880 Rev. 07**
01/2025

ISTORIJA REVIZIJA DOKUMENTA

Revizija	Datum izdanja	Promene
01	mar. 2015.	Inicijalno izdanje
02	apr. 2016.	Ažurirano je vreme testiranja i očitavanja
03	jan. 2017.	Ažurirano je kako bi se prilagodilo francuskim prevodima
04	feb. 2017.	Dodati su holandski, norveški, švedski i danski prevodi
05	dec. 2019.	Dodati su portugalski (Brazil), finski, indonežanski, korejski (Južna Koreja), poljski prevodi
06	mar. 2020.	Uklonjen je mađarski prevod
07	jan. 2025.	Reorganizovana su uputstva da imaju jasnoću i bolju definiciju uputstava; uklonjeni su prevodi

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U.S. Patent 7,709,272. Izdati patenti i patenti na čekanju u drugim zemljama. © 2025 QIAGEN
Replens[®] Feminine Moisturizer (u vlasništvu kompanije Church and Dwight Co., Inc., Juimg, Nju Džerzi)