

Ensuring viral safety through multiplexed adventitious virus detection with QIAcuity dPCR

Kathrin Wolf, Anida Mesihovic Karamitsos, Daniela Waack, Nina Laureys, Jörg Krummheuer, Colin Donohoe
QIAGEN GmbH, Hilden, Germany

Introduction

Testing for adventitious viruses, including human viral contaminants, is critical to ensure the safety and integrity of all biotechnology products derived from cell lines. The risk of viral contamination must be minimized, as such events can lead to serious clinical consequences. Regulatory bodies such as the EMA and FDA adopted the guidance on viral safety evaluation of biotechnology products derived from cell lines of human or animal origin, known as ICH Q5A(R2). This document emphasizes the need for sensitive and validated methods to detect both known and unknown viral agents that may be introduced from the original source of the cell lines or during production (*International Council for Harmonisation of Technical Requirements for Pharmaceuticals for Human Use, 2023*). Adventitious virus testing is a key component of viral safety evaluation and helps manufacturers demonstrate process control and product safety throughout biopharmaceutical development. As part of a risk-based approach and in alignment with ICH Q5A(R2), Next-Generation Sequencing (NGS) is encouraged as a powerful tool for broad, untargeted and targeted screening, and may supplement or replace traditional *in vivo* and *in vitro* assays for adventitious virus detection. While NGS is well-suited for detecting unexpected viral sequences, its use for routine or selective quality control (QC) testing may vary depending on specific program needs, testing goals, available resources and timelines.

In this context, digital PCR (dPCR) offers an efficient alternative. It can be used to screen for known viral targets

in a focused panel format or to confirm NGS findings with high specificity. For QC labs and other viral safety testing contexts, dPCR provides several practical advantages:

- Rapid turnaround for decision-critical workflows
- Multiplexing capabilities that allow simultaneous detection of multiple viral targets with minimal sample input
- High sensitivity and precision for low-level viral detection
- Robust performance across diverse sample types and matrices
- Cost-efficiency compared to sequencing-based methods

Together, these attributes make dPCR an ideal tool for routine targeted adventitious virus testing, enabling reliable detection of human viral contaminants in biologics manufacturing environments. This application note demonstrates the feasibility of dPCR for multiplexed detection of a subset of relevant human viruses.

Performance evaluation

To demonstrate the feasibility of this approach, a panel targeting selected human viruses relevant to biologics manufacturing was tested using the [QIAcuity Digital PCR System](#). The panel was comprised of 13 viral targets, organized into two independent multiplex reactions: one for RNA viruses (including RNaseP as a human sampling control, where applicable) and one for DNA viruses. These represent a diverse set of virus families, including ►

Herpesviridae, *Retroviridae*, *Flaviviridae* and *Polyomaviridae*, which are commonly referenced in regulatory guidance and biosafety testing frameworks.

The ready-to-use assays are available through the [GeneGlobe Design & Analysis Hub](#), enabling streamlined setup and integration into digital PCR workflows using the QIAcuity system (Table 1).

This human virus panel was assembled to demonstrate the potential of digital PCR for multiplex detection of representative human viral contaminants in an adventitious virus testing context.

The testing strategy used combines classical fluorophore-based multiplexing with amplitude-based multiplexing. This strategy enables the simultaneous detection of multiple viral targets within the same optical channel by adjusting probe concentrations to generate distinct signal intensities. With the [QIAcuity High Multiplex Probe PCR Kit](#), it is possible to perform higher-order multiplexing (up to 12 targets) within a single reaction on the existing QIAcuity hardware. This approach efficiently uses limited sample input and minimizes errors while maintaining high specificity and clear signal separation. In validation experiments using synthetic DNA targets, each assay

demonstrated robust performance across both 5-plex (RNA panel) and 9-plex (DNA panel) configurations (Figure 1).

This multiplexing strategy illustrates the flexibility of the QIAcuity platform for high-multiplex applications. By combining classical with amplitude-based signal separation, the workflow achieves precise discrimination of multiple viral targets without compromising assay accuracy or throughput.

Each assay was tested in both singleplex and multiplex reactions. The results demonstrated comparable quantification across both formats, confirming that multiplexing does not compromise assay performance. This consistency highlights the robustness and reliability of digital PCR for high-multiplex applications in adventitious virus testing, where maintaining sensitivity and specificity across multiple targets is essential (Figure 1, Table 2).

Furthermore, non-template controls (NTCs) consistently showed no signal, confirming the absence of contamination and highlighting assay specificity. Together, these findings support the suitability of the QIAcuity dPCR platform for viral detection workflows in QC and biosafety settings.

Table 1. Human virus panel for adventitious agents detection

Target Name	Acronym	GeneGlobe ID	Dye	Assay Concentration	RNA/DNA panel
Human betaherpesvirus 6A	HHV-6A	DMA00837	FAM*	1x	DNA
Human betaherpesvirus 6B	HHV-6B	DMA00838	FAM*	1.5x	DNA
Hepatitis B	HBV	DMA00808	HEX*	0.5x	DNA
Human parvovirus B19	B19V	DMA00862	HEX*	1x	DNA
Betapolyomavirus macacae (SV40)	SV40	DMA00936	TAMRA	1x	DNA
Human betaherpesvirus 7	HHV-7	DMA00839	ROX*	1x	DNA
Human gammaherpesvirus 8	HHV-8	DMA00836	ROX*	1.5x	DNA
Human gammaherpesvirus 4	HHV-4 (EBV)	DMA00840	Cy5	1x	DNA
Human betaherpesvirus 5	HHV-5 (HCMV)	DMA00835	ATTO700	2x	DNA
Hepatitis C	HCV	DMA00809	FAM	1x	RNA
Hepatitis A	HAV	DMA00807	HEX	1x	RNA
Human T-cell leukemia virus type I	HTLV-1	DMA00867	TAMRA	1x	RNA
Human T-cell leukemia virus 2 (2)	HTLV-2	DMA00870	ROX	1x	RNA
RNAseP(g)	-	DMA00720	ATTO700	2x	RNA **

* Amplitude-based multiplexing

** Added to the RNA panel as a human sampling control target




Figure 1. The multiplexing capabilities of the QIAcuity High Multiplex Probe PCR Kit, combined with the QIAcuity Software Suite 3.1, provide significant benefits for virus detection. dPCR Microbial DNA Detection Assays were used in a 5-plex reaction for relevant adventitious RNA viruses (Multiplex A) and a 9-plex reaction for DNA viruses (Multiplex B). Each multiplex assay setup was used to detect an individual synthetic target DNA or a pool of all synthetic target DNAs on an 8.5k nanoplate. Specific amplification of each target was observed in both singleplex and multiplex reactions. The no template control (NTC) reactions were consistently negative. Amplitude-based multiplexing allowed for the detection of two distinct targets within the Green, Yellow and Red fluorescence channels by modulating probe concentrations. For these channels, multiple thresholds (red lines) can easily be set, resulting in clearly distinguishable signal clusters. Imaging settings were adjusted to increase fluorescence intensities. A custom crosstalk matrix was applied.

Table 2. Highly concordant quantification between singleplex and multiplex reactions. Absolute quantification (copies/ μ L) of human virus targets measured on the QIAcuity dPCR System shows comparable results between singleplex and multiplex reactions

Target	Singleplex (copies/ μ L)	RNA Multiplex (copies/ μ L)	Target	Singleplex (copies/ μ L)	DNA Multiplex (copies/ μ L)
HCV	819	840	HHV-6A	1134	1117
HAV	1294	1308	HHV-6B	950	1009
HTLV-1	921	938	HBV	558	541
HTLV-2	986	986	B19V	1050	932
RNAseP	1145	1148	SV40	1155	1060
			HHV-7	1183	1008
			HHV-8	1116	1116
			HHV-4 (EBV)	1218	1140
			HHV-5 (HCMV)	1102	1069

Conclusion

The tested panel effectively demonstrates the capability of dPCR for multiplexed detection of human viruses in the context of adventitious virus testing. The human virus panel described here is available for easy ordering and convenient integration into digital PCR workflows as two assay bundles (DNA and RNA) on the GeneGlobe platform, with additional human viral targets of interest on the bundle page.

As outlined in ICH Q5A(R2), the type and extent of testing depend on a risk assessment considering the specific risk factors associated with the cell substrate and the manufacturing process. In line with this principle, users are encouraged to adopt a mix-and-match approach, selecting viral assays that best reflect the risk profile of their product and production environment.

Additional common human viral targets, such as those targeting HSV-1 and HSV-2 (Assay IDs: DMA00724 and DMA00725), HIV (Assay ID: DMA00811), HPV 16 and 18 (Assay IDs: DMA00882 and DMA007883), BK (Assay ID: DMA00894) and JC (Assay ID: DMA00861) and others, are also available through the [dPCR Microbial DNA Detection Assays page](#) as part of a wet-lab tested assay portfolio that continues to expand to support evolving biosafety and quality control needs. Furthermore, for targets not currently included in the portfolio, users can

leverage the easy-to-use [Custom dPCR Microbial Assays design tool](#). This tool allows the creation of tailored assays for any pathogen of interest, including but not limited to human viruses, supporting diverse biosafety and microbial testing needs. It is important to note that multiplexing configurations involving additional targets require testing and optimization to ensure reliable performance.

Additionally, regulatory guidelines require these methods to be validated for their intended use, based on product-specific risk assessment. Performance characteristics, including sensitivity, specificity and robustness, should be demonstrated through appropriate validation studies.

In conclusion, dPCR is valuable for targeted screening, or confirmatory testing of adventitious agents. This initiative is part of our broader commitment to supporting biopharmaceutical development with high-performance molecular testing solutions, including assays, instruments and digital tools tailored to the needs of regulated environments.

The combined use of wet-lab tested assays, customizable assay design and scalable multiplexing capabilities positions the QIAcuity Digital PCR System as a robust and adaptable platform for advancing adventitious virus testing in modern viral safety workflows.

References

International Council for Harmonisation of Technical Requirements for Pharmaceuticals for Human Use. (November 1, 2023). *Viral Safety Evaluation of Biotechnology Products Derived from Cell Lines of Human or Animal Origin (ICH Q5A(R2))* [Guideline]. https://database.ich.org/sites/default/files/ICH_Q5A%28R2%29_Guideline_2023_1101.pdf

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit instructions for use or the user operator manual. QIAGEN instructions for use and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services (or your local distributor). All figures and underlying data are from the QIAGEN R&D Department.



Discover how the QIAcuity Digital PCR System enables rapid, multiplexed adventitious virus detection: <http://www.qiagen.com/qiacuitydpcr>



Trademarks: QIAGEN®, Sample to Insight®, QIAcuity® (QIAGEN Group).

Registered names, trademarks, etc. used in this document, even when not specifically marked as such, may still be protected by law.

© 2025 QIAGEN, all rights reserved. QPRO-11882 11/2025