

April 2019

Connectivity Package Quick Start Guide



Contents

Introduction	3
Installation Procedures	3
Connecting the Q-Base	4
Connecting the Q-Base	4
To define WLAN/Wi-Fi (WPAv2 only) setting	5
To define LAN (Static IP and DHCP) setting	5

Introduction

QIAGEN brings the QIAcube® Connect customer experience to the next level with remote instrument monitoring through the Connectivity Package with the included Q-Base (an IoT gateway) and the QIAcube Connect tablet Application. The gateway can be connected to one or more QIAcube Connect instruments and serves as an interface between the instrument and the QIAcube Connect tablet Application.

Installation Procedures

This section provides you with important information on how to install the Q-Base and tablet including the required instrument setup and network settings.

If you need detailed installing information, or if you need to use your laptop for configuration via cable instead of Bluetooth®/Wi-Fi adapter, refer to the respective user manual.

Note: Before using the Q-Base, it is essential that you read the supplier's user manual carefully and pay attention to the safety information. The safety information for Q-Base can be found at <https://www.eurotech.com/en/products/iot/multi-service-iot-edge-gateways/reliagate-10-12>.

Note: You might need the help of an IT professional to enter detailed network settings.

Connecting the Q-Base


The Q-Base must be connected to the same network as QIAcube Connect instrument and tablet. The Q-Base does not need to be in proximity with the QIAcube Connect instrument to be connected.

Connecting the Q-Base

1. Power on the Q-Base by plugging the power cord into an appropriate power socket.

Note: If you use a LAN cable to connect to your network, plug the cable into ethernet port ETH1 prior to powering on the Q-Base.

Note: Bluetooth® is automatically turned off after 10 minutes.

2. Switch on the tablet and activate Bluetooth in the settings.
3. Launch the QIAcube Connect App.
4. To proceed, read and accept the "Terms and Agreement".
5. Press the **Configuration** icon () to set up the network.

Note: For successful Bluetooth communication between Q-Base and tablet, the devices must be within a maximum distance of 10 meters.

6. After a few seconds, the Q-Base name for the first-time configuration will be listed under **Available Q-Base devices for configuration of network settings**. If the Q-Base cannot be found, refresh the Bluetooth search by pressing the **Configuration** icon again.
7. To pair Q-Base with the tablet via Bluetooth, select the Q-Base name from the list.
8. Enter the default login details. The **User Name** is Admin and the **Password** is Qiagen123.
9. Select **Login**. The tablet application prompts you to change the password.
10. Enter a new password of your choice according to the listed password rules. Press **Reset**.
11. Follow the procedures below to define the LAN or Wi-Fi settings of the Q-Base.

Note: Do not press **Finish Configuration** before you enter the IP address of Q-Base into the QIAcube Connect instrument.

To define WLAN/Wi-Fi (WPAv2 only) setting

- 11a. Press the **WLAN** tab in the QIAcube Connect App to configure the Wi-Fi connection.
- 11b. Scan to search for available Wi-Fi networks.
- 11c. Select the appropriate Wi-Fi network and enter the password. The server, SSID and security type fields are automatically filled.
- 11d. Connect and proceed to step 12.

Note: Q-Base supports connecting to Wi-Fi networks with visible SSIDs only.

To define LAN (Static IP and DHCP) setting

- 11e. Select the **LAN** tab from the **Configuration** menu, enter the local network settings and set up the LAN Network. **Static IP Manual** is selected by default.

Note: When using DHCP-configured LAN, the network settings for Q-Base will be obtained automatically from the DHCP server. No further input is required.

Note: You might need the help of an IT professional to enter detailed network settings.

12. On the QIAcube Connect instrument screen, select the **Configuration** menu and then select the **Q-Base** tab.
13. Check the **Communication enabled** box.
14. Enter the Q-Base IP address in the **Q-Base URL** field on the QIAcube Connect touchscreen. The IP address is shown in the QIAcube Connect App in the **Configuration** menu under the **Q-Base** tab. You can also scan the 2-D barcode.
Note: The test will fail if you press **Test** before pressing **Finish Configuration** on the tablet.
15. To continue, select **Apply**.
16. Press **Finish Configuration**.
Note: If you do not press **Finish Configuration**, a connection will be established after 10 minutes.
17. **Optional:** To test if the connection is successful, wait for 10 minutes and then press **Test**.
18. Reboot the QIAcube Connect instrument.

For further support with technical difficulties, please contact QIAGEN Technical Services at **www.qiagen.com**.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN user manual. QIAGEN kit handbooks and user manuals are available at **www.qiagen.com** or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN[®], Sample to Insight[®], QIAcube[®] (QIAGEN Group); Bluetooth[®] (Bluetooth SIG, Inc.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.
HB-2657-001 04/2019 1117327 © 2019 QIAGEN, all rights reserved.

