

Performance Characteristics

artus GBS QS-RGQ Kit, Version 1, **REF** 4576366

Version management

This document is the *artus* GBS QS-RGQ Kit Performance Characteristics, Version 1, R1.



Check availability of new electronic labeling revisions at www.qiagen.com/p/artus-GBS-QS-RGQ-Kit-CE before test execution.

Limit of detection

The limit of detection (LoD) was assessed for the *artus* GBS QS-RGQ Kit using two serotypes of group B Streptococcus (GBS): GBS Serotype Ia (ATCC[®] BAA-1177), and Serotype III (ATCC 12403). Each strain was serially diluted in Lim broth plus matrix (Lim cultures inoculated with vaginal/rectal swabs and found to be negative for GBS) using 2-fold dilutions and run in replicates of 20 on the QIASymphony RGQ. The LoD for each strain was determined by plotting the detection rate against the bacterial input titer and performing probit analysis (Figures 1 and 2).

The limit of detection of the *artus* GBS QS-RGQ Kit is 74.07 CFU/ml for GBS Serotype Ia and 24.29 CFU/ml for Serotype III. This means that there is a 95% probability that 74.07 CFU/ml of GBS Serotype Ia and 24.29 CFU/ml of GBS Serotype III will be detected.

May 2014






Figure 1. Probit analysis: GBS Serotype Ia. Limit of detection of the *artus* GBS QS-RGQ Kit.




Figure 2. Probit analysis: GBS Serotype III. Limit of detection of the *artus* GBS QS-RGQ Kit.

Analytical reactivity

The analytical reactivity of the *artus* GBS QS-RGQ Kit was assessed to determine whether the kit could detect an array of GBS serotypes and clinical isolates. A total of 21 GBS isolates representing various serotypes (Table 1) were diluted in Lim broth plus matrix to approximately 2–3x LoD and tested with the *artus* GBS QS-RGQ Kit. GBS target was detected in all strains tested.

Table 1. GBS strains tested in analytical reactivity studies

Serotype	Strain ID	Hemolysis	GBS target detected	Internal control
Ib	ATCC*: BAA 1174	Hemolytic	+	+
V	ATCC: BAA 23	Hemolytic	+	+
VI	CDC†: ABC 20028757	Hemolytic	+	+
VII	CDC: ABC20013332	Hemolytic	+	+
VIII	CDC: ABC20026068	Hemolytic	+	+
IV	ATCC: 49446	Hemolytic	+	+
Ia	ATCC: BAA1138	Hemolytic	+	+
V	ATCC: 700046	Hemolytic	+	+
III	ATCC: BAA22	Hemolytic	+	+
V	ATCC: 700048	Hemolytic	+	+
V	ATCC: 49447	Hemolytic	+	+
Ib	BEI‡: MNZ929	Hemolytic	+	+
II	BEI: MNZ933	Hemolytic	+	+
III	BEI: MNZ938	Hemolytic	+	+

* American Type Culture Collection.

† Centers for Disease Control and Prevention.

‡ BEI Resources.

Table continued on next page

Table 1. Continued

Serotype	Strain ID	Hemolysis	GBS target detected	Internal control
V	CDC*: ABC020033655	Hemolytic	+	+
II	ATCC†: BAA1175	Hemolytic	+	+
Ib	ATCC: 55193	Hemolytic	+	+
Ic	ATCC: 27591	Hemolytic	+	+
Unknown	York Hosp 130	Non-Hemolytic	+	+
Unknown	York Hosp 144	Non-Hemolytic	+	+
Unknown	Evanston 2	Non-Hemolytic	+	+

* Centers for Disease Control and Prevention.

† American Type Culture Collection.

Cross-reactivity and microbial interference

A panel of microorganisms that may be present in patient specimens was tested to determine whether these microorganisms interfered with the detection of GBS or were cross-reactive with the *artus* GBS QS-RGQ Kit. To assess microbial interference, GBS Serotype III was diluted in Lim broth plus matrix to a concentration of approximately 2–3x LoD. Potential interfering organisms were diluted in Lim broth plus matrix. Bacteria and fungi were diluted to approximate working concentrations of $\geq 1.0 \times 10^6$ CFU/ml and viruses were diluted to approximate concentrations of $\geq 1.0 \times 10^5$ TCID₅₀/ml or PFU/ml. Bacterial and viral nucleic acids were diluted to a final concentration of $\geq 1 \times 10^6$ genome copies/ml. Human genomic DNA was diluted to a final concentration of $\geq 1 \times 10^6$ genome copies/ml. Each potentially interfering organism was tested with the *artus* GBS QS-RGQ Kit in the presence of GBS Serotype III. To assess cross reactivity, each organism from the panel was tested in the absence of GBS, at the same concentrations described above.

None of the potential interfering organisms cross-reacted or interfered with the detection of GBS Serotype III by the *artus* GBS QS-RGQ Kit.

Table 2. Panel of organisms tested for interference and cross-reactivity

Organism tested	Source ID	Interference		Cross-reactivity	
		GBS target detected	Internal control	GBS target detected	Internal control
<i>Acinetobacter baumannii</i>	ATCC* 19606	+	+	-	+
Adenovirus 40	ZMC† 0810084CF	+	+	-	+
<i>Aeromonas hydrophila</i>	ATCC 7966	+	+	-	+
<i>Bacillus cereus</i>	ATCC 11778	+	+	-	+
<i>Bacteroides fragilis</i>	ZMC 0601533	+	+	-	+
<i>Campylobacter coli</i>	ATCC 33559	+	+	-	+
<i>Candida albicans</i>	ATCC 10231	+	+	-	+
<i>C. glabrata</i>	ZMC Z007	+	+	-	+
<i>C. guilliermondii</i>	ZMC Z008	+	+	-	+
<i>C. krusei</i>	ZMC Z009	+	+	-	+
<i>C. parapsilosis</i>	ZMC Z011	+	+	-	+

* American Type Culture Collection.

† Zeptomatrix Corporation.

Table continued on next page

Table 2. Continued

Organism tested	Source ID	Interference		Cross-reactivity	
		GBS target detected	Internal control	GBS target detected	Internal control
<i>Candida tropicalis</i>	ZMC* Z012	+	+	-	+
<i>Chlamydia trachomatis</i>	ATCC† VR885	+	+	-	+
<i>Citrobacter amalonaticus</i>	ATCC 25405	+	+	-	+
<i>C. freundii</i>	ATCC 8090	+	+	-	+
<i>C. koseri</i>	ATCC 27028	+	+	-	+
Coxsackie virus	ZMC 0810075CF	+	+	-	+
Cytomegalovirus	ZMC 0810003CF	+	+	-	+
Echovirus	ZMC 0810023CF	+	+	-	+
<i>Enterobacter aerogenes</i>	ATCC 13048	+	+	-	+
<i>E. cloacae</i>	ATCC 13047	+	+	-	+
<i>Enterococcus casseliflavus</i>	ATCC 25788	+	+	-	+
<i>E. faecalis</i>	ATCC 35550	+	+	-	+

* Zeptomatrix Corporation.

† American Type Culture Collection.

Table continued on next page

Table 2. Continued

Organism tested	Source ID	Interference		Cross-reactivity	
		GBS target detected	Internal control	GBS target detected	Internal control
<i>Enterococcus faecium</i>	ATCC* BAA-2319	+	+	-	+
<i>Escherichia coli</i>	ATCC 23511	+	+	-	+
Enterovirus Type 71	ZMC† 0810047CF	+	+	-	+
Herpes simplex virus 1	ATCC VR- 1493	+	+	-	+
Herpes simplex virus 2	ATCC VR- 734	+	+	-	+
HIV	ZMC 0801032CF	+	+	-	+
HPV 11	ATCC 45151D	+	+	-	+
HPV 16	ATCC 45113D	+	+	-	+
HPV 18	ATCC 45152	+	+	-	+
<i>Klebsiella pneumoniae</i>	ATCC 13883	+	+	-	+
<i>Lactobacillus acidophilus</i>	ATCC 4356	+	+	-	+
<i>Lactococcus lactis</i>	ATCC 11454	+	+	-	+

* American Type Culture Collection.

† Zeptomatrix Corporation.

Table continued on next page

Table 2. Continued

Organism tested	Source ID	Interference		Cross-reactivity	
		GBS target detected	Internal control	GBS target detected	Internal control
<i>Listeria monocytogenes</i>	ATCC* 19115	+	+	-	+
<i>Mobiluncus mulieris</i> BV 64-5	ATCC 35240D-5	+	+	-	+
<i>Moraxella catarrhalis</i>	ATCC 8176	+	+	-	+
<i>Mycoplasma hominis</i>	ATCC 23114D	+	+	-	+
<i>Neisseria gonorrhoeae</i>	ATCC 19424	+	+	-	+
<i>N. meningitides</i>	ATCC 13077	+	+	-	+
<i>Proteus mirabilis</i>	ATCC 25933	+	+	-	+
<i>P. penneri</i>	ATCC 35198	+	+	-	+
<i>Pseudomonas aeruginosa</i>	ATCC 35554	+	+	-	+
Rotavirus	ZMC 0810041CF	+	+	-	+
<i>Serratia marcescens</i>	ATCC 13880	+	+	-	+
<i>Streptococcus mutans</i> Z072	ZMC Z072	+	+	-	+

* American Type Culture Collection.

† Zeptomatrix Corporation.

Table continued on next page

Table 2. Continued

Organism tested	Source ID	Interference		Cross-reactivity	
		GBS target detected	Internal control	GBS target detected	Internal control
<i>Staphylococcus aureus</i>	ATCC* 29213	+	+	-	+
<i>S. epidermidis</i>	ATCC 51625	+	+	-	+
<i>Streptococcus dysgalactiae</i>	ATCC 43078	+	+	-	+
<i>S. mitis</i>	ZMC† Clinical isolate	+	+	-	+
<i>S. pneumoniae</i>	ATCC 33400	+	+	-	+
<i>S. pyogenes</i>	ATCC 49399	+	+	-	+
<i>Toxoplasma gondii</i>	ZMC 0810007CF	+	+	-	+
<i>Trichomonas vaginalis</i>	ATCC 30238	+	+	-	+
<i>Ureaplasma urealyticum</i>	ATCC 27618	+	+	-	+
<i>Vibrio parahaemolyticus</i>	ATCC 17802	+	+	-	+
<i>Yersinia enterocolitica</i>	ATCC 23715	+	+	-	+
Human gDNA	Promega G304A	+	+	-	+

* American Type Culture Collection.

† Zeptomatrix Corporation.

Precision

The precision of the *artus* GBS QS-RGQ Kit was assessed using a 7-member panel consisting of 2 serotypes of GBS: Serotype Ia (ATCC BAA-1177) and Serotype III (ATCC 12403). Panel members were formulated in Lim broth plus matrix with a single serotype present (Ia or III) at 3 concentrations; positive (approximately 2–3x LoD), low positive (1x LoD), and high negative (<1x LoD). A seventh panel member (negative) was prepared using Lim broth plus matrix only. The data obtained were used to determine the mean C_T , standard deviation (SD), and the coefficient of variation (%CV) for each target and the internal control.

For the within laboratory repeatability study, the seven-member panel was tested in replicates of 3, once a day for a total of 12 days. The testing was conducted by 2 alternating operators using one instrument system (QIASymphony RGQ) and one lot of the *artus* GBS QS-RGQ Kit.

Table 3. Summary of within laboratory repeatability for the *artus* GBS QS-RGQ Kit

Serotype	Panel member	Internal control			GBS		
		Mean C_T	SD	%CV	Mean C_T	SD	%CV
Ia	Positive	29.87	0.32	1.08	31.19	0.53	1.71
	LoD	29.77	0.30	1.00	32.85	1.04	3.15
	High negative	29.78	0.31	1.05	35.54	0.88	2.47
III	Positive	29.72	0.30	1.01	32.16	0.50	1.55
	LoD	29.82	0.27	0.90	33.35	0.82	2.45
	High negative	29.76	0.29	0.96	35.93	0.62	1.72
	Negative	30.09	0.25	0.84	n.a.*	n.a.	n.a.

* n.a.: not applicable.

Lot-to-lot reproducibility was assessed using 3 different lots of the *artus* GBS QS-RGQ Kit. A single run was performed for each of the 3 lots (for a total of 3 runs) using a single instrument system (QIASymphony RGQ), by a single operator. For each run, each panel member was tested in replicates of 6 (Table 4).

Table 4. Summary of lot-to-lot reproducibility for the *artus* GBS QS-RGQ Kit

Strain	Panel member	Lot	Internal control			GBS		
			Mean C _T	SD	%CV	Mean C _T	SD	%CV
Ia	Positive	1	29.50	0.42	1.41	29.88	0.27	0.90
		2	29.46	0.17	0.59	29.87	0.13	0.43
		3	30.00	0.50	1.65	30.10	0.32	1.07
		Overall	29.65	0.44	1.49	29.95	0.26	0.87
	LoD	1	29.44	0.20	0.69	31.30	0.46	1.46
		2	29.39	0.19	0.64	31.19	0.35	1.11
		3	30.15	0.38	1.28	31.73	0.39	1.24
		Overall	29.66	0.44	1.48	31.41	0.44	1.42
	High negative	1	29.65	0.13	0.44	n.a.*	n.a.	n.a.
		2	29.46	0.33	1.13	34.00	0.51	1.50
		3	29.90	0.35	1.16	36.16	n.a.	n.a.
		Overall	29.67	0.33	1.11	34.72	1.30	3.74
III	Positive	1	29.61	0.31	1.06	31.77	0.92	2.89
		2	29.49	0.35	1.17	32.03	0.45	1.41
		3	30.05	0.31	1.04	31.91	0.41	1.28
		Overall	29.72	0.39	1.33	31.90	0.61	1.90

* n.a.: not applicable.

Table continued next page

Table 4. Continued

Strain	Panel member	Lot	Internal control			GBS		
			Mean C _T	SD	%CV	Mean C _T	SD	%CV
III	LoD	1	29.49	0.12	0.42	33.28	0.49	1.47
		2	29.54	0.26	0.87	32.56	0.44	1.37
		3	29.96	0.36	1.20	33.62	0.69	2.04
		Overall	29.66	0.33	1.12	33.15	0.69	2.08
	High negative	1	29.46	0.25	0.83	35.72	0.12	0.34
		2	29.40	0.32	1.10	36.69	0.60	1.64
		3	29.99	0.39	1.29	36.05	n.a.	n.a.
		Overall	29.62	0.41	1.38	36.32	0.64	1.76
	Negative	1	29.79	0.27	0.89	36.01	n.a.	n.a.
		2	29.69	0.18	0.59	n.a.*	n.a.	n.a.
		3	30.20	0.22	0.73	n.a.	n.a.	n.a.
		Overall	29.90	0.31	1.04	36.01	n.a.	n.a.

* n.a.: not applicable.

To measure instrument reproducibility, the 7-member panel was run on 3 different QIASymphony RGQ instrument systems. Panel members were tested in replicates of 6 by a single operator with one lot of the *artus* GBS QS-RGQ Kit (see Table 5).

Table 5. Summary of instrument-to-instrument reproducibility for the *artus* GBS QS-RGQ Kit

Strain	Panel member	Instrument	Internal control			GBS		
			Mean C _T	SD	%CV	Mean C _T	SD	%CV
Ia	Positive	1	29.63	0.24	0.80	30.11	0.20	0.66
		2	29.89	0.20	0.68	29.96	0.19	0.64
		3	29.63	0.20	0.68	30.03	0.22	0.72
		Overall	29.72	0.24	0.80	30.04	0.20	0.67
	LoD	1	29.49	0.19	0.65	31.81	0.72	2.26
		2	29.87	0.17	0.58	31.66	0.37	1.18
		3	29.57	0.35	1.19	30.90	0.17	0.56
		Overall	29.64	0.29	0.99	31.46	0.61	1.93
	High negative	1	29.66	0.36	1.22	36.37	0.01	0.02
		2	30.00	0.30	1.00	35.88	n.a.*	n.a.
		3	29.55	0.29	0.98	34.68	1.74	5.02
		Overall	29.74	0.36	1.21	35.59	1.22	3.43
III	Positive	1	29.83	0.20	0.67	32.55	0.31	0.95
		2	29.93	0.37	1.23	32.77	0.55	1.67
		3	29.47	0.22	0.74	31.99	0.26	0.82
		Overall	29.74	0.33	1.10	32.44	0.50	1.54

* n.a.: not applicable.

Table continued next page

Table 5. Continued

Strain	Panel member	Instrument	Internal control			GBS		
			Mean C _T	SD	%CV	Mean C _T	SD	%CV
III	LoD	1	29.56	0.22	0.76	33.41	0.72	2.16
		2	30.10	0.32	1.05	34.17	1.24	3.63
		3	29.63	0.24	0.81	33.16	0.57	1.73
		Overall	29.76	0.35	1.17	33.58	0.95	2.82
	High negative	1	29.51	0.26	0.88	36.22	0.82	2.26
		2	30.10	0.24	0.81	36.12	1.44	3.99
		3	29.56	0.27	0.92	35.65	0.77	2.16
		Overall	29.72	0.37	1.23	36.04	1.12	3.10
	Negative	1	30.07	0.16	0.53	n.a.*	n.a.	n.a.
		2	30.95	0.48	1.54	n.a.	n.a.	n.a.
		3	30.15	0.40	1.31	n.a.	n.a.	n.a.
		Overall	30.39	0.54	1.77	n.a.	n.a.	n.a.

* n.a.: not applicable.

To measure site-to-site reproducibility, the 7-member panel was run by 2 users at each of 3 sites. Each of the 2 users performed 5 runs on alternating testing days. Panel members were tested in replicates of 3 that were randomized and blinded to the user. A single QIA Symphony RGQ instrument system and one lot of the *artus* GBS QS-RGQ Kit were used at each site to conduct the study (Table 6).

Table 6. Summary of site-to-site reproducibility for the *artus* GBS QS-RGQ Kit

Strain	Panel member	Site	Internal control			GBS		
			Mean C _T	SD	%CV	Mean C _T	SD	%CV
Ia	Positive	1	29.51	0.31	1.05	30.25	0.34	1.14
		2	29.15	0.35	1.19	30.21	0.78	2.59
		3	29.47	0.27	0.90	30.30	0.28	0.94
		Overall	29.38	0.34	1.17	30.25	0.51	1.70
	LoD	1	29.49	0.39	1.33	31.63	0.47	1.49
		2	29.21	0.41	1.39	31.40	0.93	2.96
		3	29.43	0.32	1.07	31.40	0.39	1.26
		Overall	29.38	0.39	1.32	31.48	0.65	2.05
	High negative	1	29.47	0.35	1.17	35.32	0.95	2.70
		2	29.24	0.67	2.28	34.66	0.77	2.23
		3	29.51	0.37	1.26	35.39	1.41	4.00
		Overall	29.40	0.49	1.68	35.13	1.07	3.05
III	Positive	1	29.37	0.33	1.14	31.45	0.33	1.06
		2	29.14	0.47	1.60	31.35	0.90	2.89
		3	29.55	0.34	1.14	31.92	0.45	1.40
		Overall	29.36	0.41	1.41%	31.57	0.65	2.07

Table continued next page

Table 6. Continued

Strain	Panel member	Site	Internal control			GBS		
			Mean C _T	SD	%CV	Mean C _T	SD	%CV
III	LoD	1	29.48	0.33	1.14	31.67	0.52	1.65
		2	29.19	0.41	1.42	31.79	0.74	2.32
		3	29.40	0.33	1.13	32.50	0.54	1.67
		Overall	29.36	0.38	1.29	31.99	0.70	2.20
	High negative	1	29.38	0.35	1.20	35.29	0.82	2.33
		2	29.06	0.29	0.99	35.59	1.19	3.34
		3	29.46	0.25	0.84	35.89	0.99	2.75
		Overall	29.30	0.34	1.17	35.57	1.03	2.88
	Negative	1	29.81	0.32	1.08	36.70	0.23	0.62
		2	29.44	0.33	1.12	n.a.*	n.a.	n.a.
		3	29.73	0.39	1.30	36.77	n.a.	n.a.
		Overall	29.66	0.38	1.28	36.72	0.17	0.45

* n.a.: not applicable.

Carryover

Less than 1% carryover (cross-contamination) between samples for the entire workflow was proven by the correct detection of >99% of high positive and negative samples in alternating positions. Contrived positive GBS and negative GBS samples were prepared in Lim broth plus matrix. High positive samples were formulated at a concentration of $\geq 2 \times 10^6$ CFU/ml. These samples were processed with the complete *artus* GBS QS-RGQ workflow.

Interfering substances

A panel of 34 substances that may be present in patient specimens was tested to determine whether these substances interfered with the performance of the *artus* GBS QS-RGQ Kit. GBS Serotype III was used for this study and diluted to approximately 2–3x LoD in Lim broth plus matrix. GBS was individually spiked with each potentially inhibitory substance and tested with the *artus* GBS QS-RGQ Kit. The samples and brands are shown in Table 7. None of the substances showed an inhibitory effect on the signals of the internal control and GBS.

Table 7. Substances tested for potential interference

Substance	Potential interfering substance	Concentration tested*
Aleve [®]	Naproxen sodium	0.5mg/ml
Anti-hemorrhoidal cream/gel	Phenylephedrine HCL	0.25%
Amniotic fluid	Amniotic fluid	10% v/v
Astringent	Witch Hazel	50%
Barium sulfate	Barium sulfate	12 mg/ml
Condoms	Nonoxynol-9	7% v/v
Ethanol	Ethanol	0.05% v/v
ex-lax [®]	Sennosides	0.6% v/v
Feces	Feces	100%
Gaviscon [®]	Aluminum hydroxide	42 µg/ml
Germicidal soap	Triclosan	0.015% v/v
Human hemoglobin	Hemoglobin	0.025mg/ml
Imodium [®]	Loperamide HCL	0.8 µg/ml
Isopropyl alcohol	Isopropyl alcohol	0.05% v/v
Leukocytes	Leukocytes	2.5% v/v
Methicillin	Methicillin	0.5mg/ml

* Represents physiologically relevant concentrations of substances

Table continued on next page

Table 7. Continued

Substance	Potential interfering substance	Concentration tested*
Mineral oil	Mineral oil	1.5% v/v
Moist towelettes	Benzalkonium chloride	0.30%
Monistat Derm	Miconazole nitrate cream	2%
Mucus	Mucin	1% v/v
Palmitic acid	Palmitic acid	0.025mg/ml
Pepto Bismol®	Bismuth subsalicylate	1% v/v
Personal deodorant spray	Personal deodorant spray	100%
Preparation H®	Hydrocortisone	1%
Prilosec®	Omeprazole magnesium	8 µg/ml
Sodium hypochlorite	Sodium hypochlorite	0.05% v/v
Stearic acid	Stearic acid	0.025mg/ml
Tagamet®	Cimetidine	80 µg /ml
Tums®	Calcium carbonate	200 µg/ml
Urine	Urine	10% v/v
Vancomycin	Vancomycin	0.5mg/ml
Vaseline®	Petroleum jelly	100%
Whole blood	Whole blood EDTA	5% v/v
Zinc oxide	Zinc oxide	20% w/w

* Represents physiologically relevant concentrations of substances

Diagnostic performance evaluation

The performance of the *artus* GBS QS-RGQ Kit was evaluated using clinical samples from 4 geographically diverse locations within the United States in 2014. A total of 369 Lim broth cultures were evaluated and results from the *artus* GBS QS-RGQ Kit were compared to results obtained from reference culture and bi-directional sequencing. Clinical samples were considered positive if the C_T values were ≤ 39.7 . With the clinical samples tested, the *artus* GBS QS-RGQ Kit showed an overall diagnostic sensitivity of 95% (with a positive predictive value of 87%) and specificity of 95% (with a negative predictive value of 98%) for GBS in comparison with enriched culture and bi-directional sequencing (Tables 8 and 9).

Table 8. GBS clinical agreement study results

		Culture + bi-directional sequencing		
		+	-	Total
artus GBS QS-RGQ Kit	+	90	13	103
	-	5	261	266
Total		95	274	369

Table 9. GBS clinical agreement study results

		95% CI
Sensitivity	95%	88–98%
Specificity	95%	92–97%
Positive predictive value	87%	80–92%
Negative predictive value	98%	96–99%
Prevalence	26%	22–30%
Agreement	95%	N/A

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, QIASymphony®, *artus*®, Rotor-Gene® (QIAGEN Group); Aleve® (Bayer HealthCare LLC); ATCC® (American Type Culture Collection Corporation); *ex-lax*® (Sandoz Pharmaceuticals Corporation); Gaviscon® (Aventis Holdings Inc); Imodium® (Johnson & Johnson Corporation); Preparation H® (American Home Products Corporation); Prilosec® (AstraZeneca AB Corporation); Tagamet®, Tums® (SmithKline Beecham Pharmaceuticals Co.); Vaseline® (Unilever Supply Chain, Inc.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

May-14 © 2014 QIAGEN, all rights reserved.

www.qiagen.com

Australia = 1-800-243-800
Austria = 0800-281011
Belgium = 0800-79612
Brazil = 0800-557779
Canada = 800-572-9613
China = 800-988-0325

Denmark = 80-885945
Finland = 0800-914416
France = 01-60-920-930
Germany = 02103-29-12000
Hong Kong = 800 933 965
India = 1-800-102-4114
Ireland = 1800 555 049

Italy = 800-787980
Japan = 03-6890-7300
Korea (South) = 080-000-7145
Luxembourg = 8002 2076
Mexico = 01-800-7742-436
The Netherlands = 0800 0229592
Norway = 800-18859

Singapore = 1800-742-4368
Spain = 91-630-7050
Sweden = 020-790282
Switzerland = 055-254-22-11
Taiwan = 0800-665-957
UK = 0808-234-3665
USA = 800-426-8157



Sample & Assay Technologies