

# miScript<sup>®</sup> Single Cell qPCR Kit

## Protocols: Cell Lysis, 3' Ligation, 5' Ligation, Reverse-Transcription

### Further information

- *miScript Single Cell qPCR System Handbook*: [www.qiagen.com/handbooks](http://www.qiagen.com/handbooks)
- Safety Data Sheets: [www.qiagen.com/safety](http://www.qiagen.com/safety)
- Technical assistance: toll-free 00800-22-44-6000, or [www.qiagen.com/contact](http://www.qiagen.com/contact)

### Notes before starting

- Add 2.5 µl Buffer FCPM to Buffer FCPL, mix well and store at 2–8°C. This only needs to be added once. Briefly shake or vortex the FCPL / FCPM mixture before each use.
- Add 750 µl RNase-free water to lyophilized gDNA Wipeout Buffer 2. Mix by gently inverting the vial, divide into single-use aliquots and store at –20°C.
- Prepare the reagents required for each protocol according to the *miScript Single Cell qPCR System Handbook*.
- Use 1–100 cells as starting material. The cells should be in 2 µl or less of neutral buffer like PBS and thawed on ice.
- The 3' ligation and 5' ligation reactions are very viscous. Pipet slowly and thoroughly (pipet up and down 12 times) to mix.
- Do not vortex reactions or reagents unless instructed.
- Ensure master mixes are centrifuged briefly before mixing, thoroughly mixed by pipetting up and down 12 times and briefly centrifuged once again.
- Ensure reactions are thoroughly mixed, prepared at recommended temperatures and incubated at recommended temperatures.



## Cell Lysis

1. Prepare the cell processing master mix on ice according to Table 1.

**Table 1. Preparation of cell processing master mix**

Component	Master mix 24 samples	Master mix 96 samples
Buffer FCPL/FCPM	124 $\mu$ l	496 $\mu$ l
gDNA Wipeout Buffer	8 $\mu$ l	32 $\mu$ l
<b>Total volume</b>	<b>132 <math>\mu</math>l</b>	<b>528 <math>\mu</math>l</b>

2. Add 3  $\mu$ l of cell processing master mix into each well of the PCR strip tubes or 96-well plates containing a cell in 2  $\mu$ l PBS. Do not mix the cell / cell processing master mix.

**Important:** If the cells are in less PBS, bring final volume to 5  $\mu$ l with additional mix.

3. Incubate for 5 min at room temperature, 5 min at 75°C and hold at 4°C.
4. Proceed to 3' Ligation. Alternatively, the lysed cells can be stored at -80°C.

## 3' Ligation

1. Thaw PCR strip tube(s) or 96-well plate containing cell lysate on ice. Centrifuge briefly.
2. Prepare the 3' ligation reaction master mix at room temperature according to Table 2.

**Table 2. Setup of 3' ligation reactions**

Component	Volume one sample	Master mix 24 samples	Master mix 96 samples
Cell lysate (already in the well of a PCR strip tube or 96-well plate)	5 $\mu$ l	–	–
miScript SC 3' Ligation Buffer	1 $\mu$ l	27 $\mu$ l	108 $\mu$ l
Nuclease-free water	0.5 $\mu$ l	13.5 $\mu$ l	54 $\mu$ l
miScript SC 3' RNA Ligase	0.5 $\mu$ l	13.5 $\mu$ l	54 $\mu$ l
miScript SC Ligation Activator (Step 4)	8 $\mu$ l	–	–
<b>Total volume</b>	<b>15 <math>\mu</math>l</b>	<b>54 <math>\mu</math>l</b>	<b>216 <math>\mu</math>l</b>

- At room temperature, aliquot 2  $\mu\text{l}$  of the 3' ligation master mix into each well of the PCR strip tubes or 96-well plates containing 5  $\mu\text{l}$  cell lysate.
- Slowly add miScript SC Ligation Activator to each well. Briefly centrifuge, mix by pipetting up and down 12 times and briefly centrifuge again.

**Important:** Pipet slowly to mix. miScript SC Ligation Activator is very viscous.

- Incubate for 1 h at 16°C and hold at 4°C.
- Proceed immediately to 5' Ligation.

### 5' Ligation

- Immediately prior to setting up the 5' ligation reactions, allow the completed 3' ligation reactions to equilibrate to room temperature for 2 min. Centrifuge briefly.
- Prepare the 5' ligation reaction master mix at room temperature according to Table 3.

**Table 3. Setup of 5' ligation reactions**

Component	Volume one sample	Master mix 24 samples	Master mix 96 samples
3' ligation reaction (already in the well of a PCR strip tube or 96-well plate)	15 $\mu\text{l}$	–	–
miScript SC 5' Ligation Buffer	4 $\mu\text{l}$	108 $\mu\text{l}$	432 $\mu\text{l}$
miScript SC 5' RNA Ligase	1 $\mu\text{l}$	27 $\mu\text{l}$	108 $\mu\text{l}$
<b>Total volume</b>	<b>20 <math>\mu\text{l}</math></b>	<b>135 <math>\mu\text{l}</math></b>	<b>540 <math>\mu\text{l}</math></b>

- At room temperature, aliquot 5  $\mu\text{l}$  of the 5' ligation master mix into each well of the PCR strip tubes or 96-well plates containing the completed 3' ligation reaction. Briefly centrifuge, mix by pipetting up and down 12 times and briefly centrifuge again.

**Important:** Pipet slowly to mix. The reactions are very viscous.

- Incubate for 5 min at 37°C, 15 min at 65°C and hold at 4°C.

**Important:** Hold at 4°C for at least 5 min.

5. Proceed immediately to *Reverse-Transcription*.

## Reverse-Transcription

1. Prepare the reverse-transcription reaction master mix on ice according to Table 4.

**Table 4. Setup of reverse-transcription reactions**

Component	Volume one sample	Master mix 24 samples	Master mix 96 samples
5' ligation reaction (already in the well of a PCR strip tube or 96-well plate)	20 $\mu$ l	–	–
5x miScript SC RT Buffer	8 $\mu$ l	216 $\mu$ l	864 $\mu$ l
10x miScript SC RT Nucleics	4 $\mu$ l	108 $\mu$ l	432 $\mu$ l
Nuclease-free water	7 $\mu$ l	189 $\mu$ l	756 $\mu$ l
miScript SC Reverse Transcriptase	1 $\mu$ l	27 $\mu$ l	108 $\mu$ l
<b>Total volume</b>	<b>40 <math>\mu</math>l</b>	<b>540 <math>\mu</math>l</b>	<b>2160 <math>\mu</math>l</b>

2. On ice, aliquot 20  $\mu$ l of the reverse-transcription master mix into each well of the PCR strip tubes or 96-well plates containing the completed 5' ligation reaction.

3. Incubate for 2 h at 37°C, 5 min at 95°C and hold at 4°C.

**Important:** Hold at 4°C for at least 5 min.

4. Proceed to the second Quick-Start protocol associated with the miScript Single Cell qPCR Kit (Protocols: "Beads + Bind" Mixture Preparation, cDNA Cleanup, Pre-amplification). Alternatively, the completed cDNA syntheses can be stored at –20°C in a constant-temperature freezer.