

EasyXpress® Disulfide *E. coli* Kit (96)

The EasyXpress Disulfide *E. coli* Kit (96) (cat. no. 32578), including buffers and reagents, should be stored immediately upon receipt at -80°C in a constant-temperature freezer.

For more information, please refer to the *EasyXpress Disulfide E. coli Handbook*, which can be found at www.qiagen.com/handbooks.

For technical assistance, please call toll-free 00800-22-44-6000, or find regional phone numbers at www.qiagen.com/contact.

Notes before starting

- This protocol is for the expression of multiple recombinant, disulfide-bonded proteins and protein complexes, especially Fab and scFv antibody fragments, using linear templates or plasmid DNA and the EasyXpress Disulfide *E. coli* Kit.
- For template design, refer to “Expression templates” in the *EasyXpress Disulfide E. coli Handbook*.
- A PCR cycler, stepper pipet, and centrifuge with a rotor for 96-well plates are required.
- The in vitro translation system is extremely sensitive to nuclease contamination. Wear gloves and use RNase- and DNase-free reaction tubes and pipet tips with filters.
- EasyXpress Disulfide *E. coli* Extract is provided in four individual tubes, each containing 1160 μl . Once thawed, use *E. coli* extracts within 4 h. If required, the EasyXpress Disulfide *E. coli* Kit (96) can be used for just 24 reactions. Use one tube of EasyXpress Disulfide *E. coli* Extract to prepare a volume of master mix for 24 reactions, including an excess of 10%.
- One 96-well plate and silicone mat are provided with the EasyXpress Disulfide *E. coli* Kit (96). For 24 reactions, place tape over the empty wells.
- Additional plates can be ordered from nerbe plus GmbH (cat. no. 040830140).
- Except for the transcription/translation incubation, all handling steps should be performed on ice. For protein synthesis reactions, add the components in the order given in the protocol and Table 1.

October 2011



- After incubation for 4 h, reactions are boosted with EasyXpress Disulfide *E. coli* Folding Buffer to enhance disulfide bond formation. Usually, a total incubation time of 6 h is sufficient for disulfide bond formation. In some cases, extending the incubation time for up to 20 h after the boost may help ensure completeness of disulfide bond formation.
- If using PCR products generated with EasyXpress Linear Template Kits as templates, thaw XE-Solution on ice. XE-Solution is provided with EasyXpress Linear Template Kits.

1. Dilute linear expression templates or plasmids for the number of reactions required in the 96-well plate provided. If necessary, the plate can be tightly sealed and stored at 4°C overnight.

Plasmid DNA template: For each reaction, prepare a plasmid DNA expression template solution containing 2 µg of plasmid in 24 µl of RNase-free water. If expressing Fab heavy and light chains from separate plasmids, mix 1 µg each of the heavy chain plasmid and the light chain plasmid in a total volume of 24 µl of RNase-free water.

PCR product template: For each reaction add 4 µl XE-Solution to ~2 µg DNA from the second PCR (performed using EasyXpress Linear Template Kits) in a final volume of 24 µl, and mix by pipetting up and down. If expressing Fab heavy and light chains from separate PCR products, add a total of 2.5 µg PCR product.

2. The EasyXpress Disulfide *E. coli* Kit (96) can be used to prepare a master mix for 24 samples. Use one tube of EasyXpress Disulfide *E. coli* Extract to prepare a volume of master mix for 24 reactions, including a 10% excess. Add the required volumes of EasyXpress Disulfide *E. coli* Reaction Buffer and RNase-free water as specified in Table 1.
3. Thaw and store EasyXpress Disulfide *E. coli* Extract and EasyXpress Disulfide *E. coli* Reaction Buffer on ice. Make sure not to invert the tube containing the EasyXpress Disulfide *E. coli* Extract. Thaw RNase-free water at room temperature (15–25°C). Gently mix the EasyXpress Disulfide *E. coli* Extract by pipetting before use. Make sure that the EasyXpress Disulfide *E. coli* Reaction Buffer is completely thawed and fully dissolved.
4. Calculate the amount of template, XE-Solution, and water needed to obtain a final volume of 24 µl for the template mix. Table 1 provides a pipetting scheme for reaction setup.

Table 1. Pipetting scheme for setup of EasyXpress Disulfide *E. coli* protein synthesis reactions

Reagent/component (for 24 reactions, including 10% excess)	Plasmid DNA template	PCR product/linear template	Positive control	No template control
Add plasmid, linear template, positive control, or RNase-free water for no template control to 96-well plate				
Template DNA	2 μ g/24 μ l per well	Approx. 2 μ g	10 μ l in one well	–
XE-Solution*	–	4 μ l per well [†]	–	–
Add RNase-free water to bring the volume to 24 μ l per well				
Preparation of master mix for 24 reactions				
EasyXpress Disulfide <i>E. coli</i> Extract		1160 μ l		
EasyXpress Disulfide <i>E. coli</i> Reaction Buffer		1188 μ l		
RNase-free water		712 μ l		
Total master mix volume		3060 μ l		
Volume of master mix per well		116 μ l		
Total reaction volume per well (before boost)		140 μ l		
Incubate for 4 h at 27°C				
EasyXpress Disulfide <i>E. coli</i> Folding Buffer		30 μ l per well		
Total		170 μ l per well		
Incubate for an additional 2–20 h at 27°C				

* XE-solution should be used only when PCR templates from EasyXpress Linear Template Kits are used. Mix XE-Solution with the second PCR product before adding to the reaction. It is important that XE-Solution is not added directly to the diluted *E. coli* extract.

[†] If expressing Fab heavy and light chain from separate PCR products, add a total amount of ~1.25 μ g PCR product per chain (approx. 8.5 μ l of a PCR product concentrated at 150 ng/ μ l).

EasyXpress Disulfide *E. coli* Positive Control DNA (yellow screw-cap) is supplied with the kit. The heavy chain of the synthesized anti-CD4 Fab that serves as a positive control contains a C-terminal 6xHis tag.

5. Add all components in the given order and mix the master mix by pipetting.
6. Using a stepper pipet, pipet the master mix into the 96-well plate containing the templates.
7. Incubate the reactions for 4 h at 27°C in a PCR cycler. Thaw EasyXpress Disulfide *E. coli* Folding Buffer at room temperature.
8. After incubation for 4 h, boost each expression reaction by adding 30 μ l of EasyXpress Disulfide *E. coli* Folding Buffer Mix and centrifuge briefly. Incubate for an additional 2–20 h at 27°C.
9. Proceed to “Sample analysis” or store samples at –80°C until ready for analysis.

Sample analysis

1. Centrifuge expression reactions for 1 h at 4000 x g at room temperature (15–25°C). If expression reactions have been stored at –80°C, thaw samples before centrifugation.
2. Transfer the supernatant (soluble protein fraction) to a fresh 96-well plate (not supplied). Store on ice, if necessary. If the insoluble protein fraction is to be analyzed, resuspend the pellet in 170 μ l PBS containing 0.5% Triton® X-100 by vigorous pipetting. For western blot analysis, use 5–10 μ l of supernatant or suspended pellet per gel lane. For ELISA assays, perform serial dilutions of the soluble fraction starting at 1:100.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN®, EasyXpress® (QIAGEN Group); Triton® (Union Carbide Corporation). 1070005 10/2011 © 2011 QIAGEN, all rights reserved.

