

## Type-it<sup>®</sup> CNV Probe PCR+ qC Kit

The Type-it CNV Probe PCR+ qC Kit (cat. nos. 206662 and 206664), including buffers and reagents, should be stored immediately upon receipt at  $-20^{\circ}\text{C}$  in a constant-temperature freezer and should be protected from exposure to light. The 2x Type-it Probe PCR Master Mix should also be protected from exposure to light and can be stored at  $2-8^{\circ}\text{C}$  for up to 1 month (depending on the expiration date), without showing any reduction in performance.

For more information, please refer to the *Type-it CNV Probe PCR+ qC Kit Handbook*, which can be found at [www.qiagen.com/handbooks](http://www.qiagen.com/handbooks).

For technical assistance, please call toll-free 00800-22-44-6000, or find regional phone numbers at [www.qiagen.com/contact](http://www.qiagen.com/contact).

### Notes before starting

- This protocol is optimized for relative quantification ( $\Delta\Delta\text{C}_T$ ) of DNA copy number in the human genome and is intended for use with the Type-it CNV Probe PCR+ qC Kit and TaqMan<sup>®</sup> probes on all real-time cyclers, including the Rotor-Gene<sup>®</sup> Q, and all instruments from Applied Biosystems, Bio-Rad, Cepheid, Eppendorf, Roche, and Agilent. Using this protocol, duplex or multiplex PCR is carried out in the presence of ROX<sup>™</sup> passive reference dye.
- For each gene of interest (GOI), we recommend using the provided Type-it CNV Reference Probe Assay in the same reaction well as a universal reference assay for reliable  $\Delta\Delta\text{C}_T$ -based quantification of the CNV in the human genome. The TaqMan probe of the reference assay is labeled with MAX<sup>™</sup> dye, which is detected on the HEX<sup>™</sup>/VIC<sup>®</sup> channel.
- We recommend preparing a 25x primer–probe mix containing specific primers and probe (recommended concentrations in the 25x primer–probe mix:  $10\ \mu\text{M}$  of each primer and  $5\ \mu\text{M}$  of the probe) for each GOI. See Tables 1 and 3 for details.
- It is recommended to use a FAM<sup>™</sup> labeled probe to detect the gene of interest by duplex PCR on any cycler. For multiplex analyses, we strongly recommend using dual-labeled probes with nonfluorescent quenchers (standard TaqMan or TaqMan MGB<sup>™</sup> Probes). For information on suitable combinations of reporter dyes for triplex or 4-plex PCR on various cyclers, please refer to *Type-it CNV*

February 2012



*Probe PCR+ qC Kit Handbook*, which can be found at [www.qiagen.com/handbooks](http://www.qiagen.com/handbooks).

- To reconstitute the Type-it CNV Reference Probe Assay (100) to a 25x working solution, briefly centrifuge the tube, add 110  $\mu$ l Buffer TE (pH 8.0; provided with the kit), and mix by vortexing the tube 4–6 times. If necessary, gently warm the tube to help the primers dissolve.
  - Always use 30 pg–30 ng template DNA (see Tables 1–3) and the optimized cycling conditions specified in Table 4.
  - The PCR must start with an initial incubation step of 5 min at 95°C to activate HotStarTaq® Plus DNA Polymerase.
1. Thaw the Type-it Probe PCR Master Mix, template DNA, Type-it CNV Reference Probe Assay 25x working solution, 25x primer–probe mix containing primers and probe for the gene of interest, and RNase-free water. Mix the individual solutions and place them on ice.
  2. Prepare a reaction mixture according to Tables 1 or 2 for duplex PCR and Table 3 for triplex and 4-plex PCR.
  3. Mix the reaction mixture thoroughly and dispense appropriate volumes into PCR tubes or the wells of a PCR plate or Rotor-Disc®.
  4. Add template DNA (start with 10 ng; use 30 pg–30 ng as the range) to the individual PCR tubes or wells. See Tables 1–3 for details.
  5. Program the real-time cycler according to Table 4.

**Note:** For correct instrument setup for multiplex analysis (e.g., setting up detection of multiple dyes from the same well), refer to the user manual for the real-time cycler and follow the manufacturer’s instructions. Be sure to activate the detector for each reporter dye used. Depending on your instrument, it may also be necessary to perform a calibration procedure for each of the reporter dyes before they are used for the first time.

6. Place the PCR tubes, Rotor-Disc, or plate in the real-time cycler, and start the cycling program.
7. Perform data analysis.

**Note:** Before performing data analysis, select the analysis settings for each probe (i.e., baseline settings and threshold values). Optimal analysis settings are a prerequisite for accurate quantification data.

**Note:** Only if using the Applied Biosystems® 7500, 7500 Fast, or ViiA™ 7 Real-Time PCR Systems, it is recommended to use the ‘manual C<sub>T</sub>’ function instead of the ‘auto C<sub>T</sub>’ function for data analysis. Use a value of 0.01 as a

starting point for the threshold setting. For all other cyclers, use the automatic  $C_T$  function as a starting point.

**Table 1. Reaction setup for duplex PCR**

Component	Volume ( $\mu$ l)
<b>Reaction mix</b>	
2x Type-it Probe PCR Master Mix	12.5
25x Type-it CNV Reference Probe Assay solution	1
25x primer–probe mix for GOI* (FAM labeled)	1
RNase-free water	Variable
<b>Template DNA (added at step 4)</b>	Start with 10 ng (use 30 pg–30 ng as the range)
<b>Total reaction volume</b>	<b>25<sup>†</sup></b>

\* **IMPORTANT:** For duplex PCR, a 25x primer–probe mix consists of 10  $\mu$ M forward primer, 10  $\mu$ M reverse primer, and 5  $\mu$ M probe in Buffer TE, resulting in a final concentration of 0.4  $\mu$ M forward and reverse primer and 0.2  $\mu$ M probe.

<sup>†</sup> If the real-time cycler requires a final reaction volume other than 25  $\mu$ l, adjust the amount of master mix and all other reaction components accordingly. If using 384-well plates on the ABI PRISM<sup>®</sup> 7900, use a reaction volume of 10  $\mu$ l.

**Table 2. Reaction setup for duplex PCR using TaqMan Copy Number Assays from Life Technologies<sup>‡</sup>**

Component	Volume ( $\mu$ l)
<b>Reaction mix</b>	
2x Type-it Probe PCR Master Mix	12.5
25x Type-it CNV Reference Probe Assay solution	1
20x TaqMan Copy Number Assay for GOI (FAM labeled) <sup>‡</sup>	1.25
RNase-free water	Variable
<b>Template DNA (added at step 4)</b>	Start with 10 ng (use 30 pg–30 ng as the range)
<b>Total reaction volume</b>	<b>25<sup>§</sup></b>

<sup>‡</sup> Delivered at a concentration of 20x.

<sup>§</sup> If the real-time cycler requires a final reaction volume other than 25  $\mu$ l, adjust the amount of master mix and all other reaction components accordingly. If using 384-well plates on the ABI PRISM 7900, use a reaction volume of 10  $\mu$ l.

**Table 3. Reaction setup for triplex and 4-plex PCR**

Component	Volume ( $\mu$ l)
<b>Reaction mix</b>	
2x Type-it Probe PCR Master Mix	12.5
25x Type-it CNV Reference Probe Assay solution	1
25x primer–probe mix for GOI 1*	1
25x primer–probe mix for GOI 2*	1
<b>Only for 4-plex PCR:</b> 25x primer–probe mix for GOI 3*	1
RNase-free water	Variable
<b>Template DNA (added at step 4)</b>	Start with 10 ng (use 30 pg–30 ng as the range)
<b>Total reaction volume</b>	<b>25<sup>†</sup></b>

\* **IMPORTANT:** For triplex and 4-plex RT-PCR, a 25x primer–probe mix consists of 10  $\mu$ M forward primer, 10  $\mu$ M reverse primer, and 5  $\mu$ M probe in Buffer TE, resulting in a final concentration of 0.4  $\mu$ M forward and reverse primer, and 0.2  $\mu$ M probe.

<sup>†</sup> If the real-time cycler requires a final reaction volume other than 25  $\mu$ l, adjust the amount of master mix and all other reaction components accordingly. If using 384-well plates on the ABI PRISM 7900, use a reaction volume of 10  $\mu$ l.

**Table 4. Cycling conditions**

Step	Time	Temperature
<b>PCR initial activation step</b>	5 min	95°C
<b>2-step cycling:</b>	30 s	95°C
Denaturation		
Annealing/extension	30 s <sup>‡</sup>	60°C
<b>Number of cycles</b>	40	

<sup>‡</sup> For real-time PCR systems that require a minimum annealing/extension time longer than 30 seconds, adjust the time accordingly, as per the requirements of the real-time PCR system in operation.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN®, HotStarTaq®, Rotor-Disc®, Rotor-Gene®, Type-it® (QIAGEN Group); ABI PRISM®, Applied Biosystems®, FAM™, HEX™, MAX™, ROX™, VIC®, ViiA™ (Life Technologies Corporation); MGB™ (Epoch Biosciences); TaqMan® (Roche Group). 1070459 02/2012 © 2012 QIAGEN, all rights reserved.

