

Stefanie Schröer, Dominic O'Neil, Claudia Dienemann, Thomas Deutschmann, Markus Sprenger-Haussels  
 QIAGEN GmbH, QIAGEN Strasse 1, 40724 Hilden, Germany


## Novel, Procedure for Extraction of Microbial gDNA from Environmental Samples

Studies of the environmental and human microbiome face the problem that sample materials, such as soil and stool, are extremely difficult sample types. This is due to their high abundance of inhibitory substances such as humic and fulvic acids in soil, and bile acids and proteoglycans in stool. These inhibitory substances are often co-purified and have the potential to inhibit enzymatic reactions. Extraction of microbial DNA is therefore challenging and often requires tedious and time-consuming workflows, especially with large sample sizes in high-throughput studies.

To address these challenges, we have developed a second generation of the QIAGEN®/MO BIO Inhibitor Removal Technology®, which efficiently removes PCR inhibitors from even the most difficult soil types and stool samples. High-throughput bead-beating homogenization allows rapid, thorough, and unbiased lysis by a combination of mechanical and chemical methods. Inhibitory substances are removed by a dedicated precipitation step, all in a fast and streamlined workflow. Subsequent DNA purification is automated on a 96-well silica membrane plate format on the QIAcube® HT instrument using vacuum technology. This technology enables high-throughput processing of environmental samples for downstream applications, such as PCR amplification and NGS analysis, in under 3 hours.


The automated protocol provides equivalent yields and purity to manual workflows, while enabling much higher throughput. 16S sequencing shows equivalent community composition to the manual methods, as shown by alpha and beta diversity measures. The new automated protocol allows for a complete picture of the diversity in the samples, thereby demonstrating that its suitability for high-throughput sample processing of both environmental and stool samples.

## Convenient, Automated, High-Throughput Workflow



## Reproducible, High-Quality DNA Yields from Soil and Stool


Soil and human stool samples were prepared using the new DNeasy 96 PowerSoil Pro QIAcube HT Kit. Microbial gDNA was either extracted from individually processed samples (A) or from aliquots of a lysate pool (B).




The DNeasy 96 PowerSoil Pro QIAcube HT yields high amounts of DNA from soil and stool samples with high reproducibility and low variance over the 96-well plate.

## Inhibitor Removal Technology Increases Purity of Isolated DNA

Soil samples were processed with the new DNeasy 96 PowerSoil Pro QIAcube HT Kit.



## Identification of Total Diversity and Community Representation in Soil and Stool Samples



## Conclusions

Relatively inexpensive, highly multiplexed 16S sequencing enables studies of large sample size, for which improved sample extraction techniques are needed in order to enable true high-throughput studies. A new semi-automated, high-throughput method is presented that:

- Allows simultaneous purification of microbial genomic DNA of up to 96 samples in a convenient semi-automated workflow using the QIAcube HT instrument
- Efficiently lyses bacteria in soil and stool samples
- Provides high yields of useful, pure DNA
- Removes inhibitors commonly present in soil and stool samples, allowing eluates to be used directly in downstream NGS applications.

These data show that the new high-throughput method on the QIAcube HT achieves successful DNA isolation from soil and stool samples in a semi-automated 96-well format.

The method presented here is intended for molecular biology applications. This method is not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, Sample to Insight®, QIAcube®, DNeasy®, Inhibitor Removal Technology®, PowerSoil® (QIAGEN Group). MiSeq® (Illumina, Inc.); NanoDrop™, Qubit®, (Thermo Fisher Scientific, Inc.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law. © 2019 QIAGEN, all rights reserved. PROM-13956-001