



July 2025

QIAseq[®] Targeted cfDNA Ultra Handbook

For ultrasensitive targeted next-generation sequencing (NGS) of cfDNA for Illumina[®] NGS systems

Table of Contents

Kit Contents	3
Shipping and Storage	6
Intended Use	7
Safety Information	8
Quality Control	8
Introduction	9
Principle and Procedure	10
UMI assignment	12
Target enrichment and final library construction	12
NGS adapter and index technologies	12
Principle of variant detection with UMIs	13
Equipment and Reagents to Be Supplied by User	15
Important Notes	16
Protocol: DNA End Prepare	23
Protocol: Adapter Ligation	26
Protocol: Target Enrichment	30
Protocol: Universal PCR	35
Recommendations: Library QC and Quantification	41
Recommended setting for checking QIAseq cfDNA Ultra library on QIAxcel Advanced and Connect systems	41
Protocol: Sequencing Setup on Illumina MiSeq, NextSeq 500/550, NextSeq 1000/2000, MiniSeq and NovaSeq	44
Troubleshooting Guide	48
Appendix A: Combining an Existing Panel with a Booster Panel	50
Appendix B: Data Analysis Using QIAGEN's CLC Genomic Workbench	51
Ordering Information	52
Document Revision History	54

Kit Contents

QIAseq Targeted cfDNA Ultra

Catalog no.	334051	334055
No. of samples	12	96
One pool of region-specific primers	60 µL	480 µL
ERA Buffer, 10x	60 µL	480 µL
ERA Enzyme	60 µL	480 µL
UPH Ligation Buffer, 2.5x	350 µL	2 x 1250 µL
DNA Ligase	75 µL	600 µL
HiFi Ultra Buffer, 5x	144 µL	2 x 768 µL
HiFi Ultra Polymerase	40 µL	315 µL
Nuclease-free Water	1.5 mL	10 mL
QIAseq Bead Binding Buffer (one bottle)	7 mL	54 mL
QIAseq Beads (one bottle)	10 mL	38.4 mL

QIAseq Targeted cfDNA Ultra Booster

Catalog no.	334085
No. of samples	96
One pool of region-specific primers	80 µL

Cat. no.	Product name	Total number of primers*	Panel size (bases)
YHS-001Z	Breast Cancer Ultra Panel	1045	21,960
YHS-002Z	Colorectal Cancer Ultra Panel	529	6,618
YHS-005Z	Lung Cancer Ultra Panel	948	26,448
YHS-101Z	Actionable Solid Tumor Ultra Panel	2284	44,012
YHS-003Z	Myeloid Neoplasms Ultra Panel	2302	32,716

* The number of primers in Custom and Booster panels is represented by the last digits of the catalog number. For example, a booster panel with catalog number SYHS-00100Z-48 has 48 primers.

QIAseq cfDNA Ultra Reagent Kit

QIAseq cfDNA Ultra Reagent Kit	
Catalog no.	334455
No. of samples	96
ERA Buffer, 10x	480 µL
ERA Enzyme	480 µL
UPH Ligation Buffer, 2.5x	2 x 1250 µL
DNA Ligase	600 µL
HiFi Ultra Buffer, 5x	2 x 768 µL
HiFi Ultra Polymerase	315 µL
Nuclease-free Water	10 mL
QIAseq Bead Binding Buffer (one bottle)	54 mL
QIAseq Beads (one bottle)	38.4 mL

QIAseq Targeted cfDNA Ultra 96-Unique Dual Indices

QIAseq Targeted cfDNA Ultra UDI Set*	Set A	Set B	Set C	Set D
Catalog no.	334165	334175	334185	334195
No. of samples	96	96	96	96
QUDI-96AA †	9 µL	N/A	N/A	N/A
QUDI-96BA †	N/A	9 µL	N/A	N/A
QUDI-96CA †	N/A	N/A	9 µL	N/A
QUDI-96DA †	N/A	N/A	N/A	9 µL
AdP-DP-Phased Adapter	144 µL	144 µL	144 µL	144 µL
SmP-IL5 TEPCRF Primer	250 µL	250 µL	250 µL	250 µL

* 10 bp dual indices.

† Index Primer Plate (DNA Ultra UDI Set A, B, C and D Final Plate); each plate contains 96 pairs of sample index primers plus universal primers, with each well corresponding to one pair of UDI sample index; each index is single-use.

QIAseq Targeted cfDNA Ultra 12-Unique Dual Indices

QIAseq Targeted cfDNA Ultra UDI (12)*

Catalog no.	334151
No. of samples	12
QUDI-12A†	9 µL
AdP-DP-Phased Adapter	18 µL
SmP-IL5 TEPCR-F Primer	31 µL

* 10 bp dual indices.

† Index Primer Plate (DNA Ultra UDI 12 Index Final Plate); each plate contains 12 pairs of sample index primers plus universal primers, with each well corresponding to one pair of UDI sample index; each index is single-use.

Shipping and Storage

The QIAseq Targeted cfDNA Ultra Kits (except QIAseq Beads and QIAseq Bead Binding Buffer) are shipped on dry ice and should be stored at -30°C to -15°C in a constant-temperature freezer upon arrival. The QIAseq Beads and QIAseq Bead Binding Buffer are shipped on cold packs and should be stored at $2-8^{\circ}\text{C}$ upon arrival.

The QIAseq Targeted cfDNA Ultra Index Kits are shipped on dry ice and should be stored at -30°C to -15°C upon arrival.

When stored correctly, the QIAseq Targeted cfDNA Ultra Kits are good until the expiration date printed on the kit label.

Intended Use

The QIAseq Targeted cfDNA Ultra Kits are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

All due care and attention should be exercised in the handling of the products. We recommend all users of QIAGEN® products to adhere to the NIH guidelines that have been developed for recombinant DNA experiments or to other applicable guidelines.

Safety Information

When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at www.qiagen.com/safety, where you can find, view and print the SDS for each QIAGEN kit and kit component.

Quality Control

In accordance with QIAGEN's ISO-certified Quality Management System, each lot of the QIAseq Targeted cfDNA Ultra is tested against predetermined specifications, to ensure consistent product quality.

Introduction

The QIAseq Targeted cfDNA Ultra enables streamlined Sample to Insight[®], targeted next-generation sequencing (NGS) of cell-free DNA (cfDNA). This highly optimized, automation-friendly solution facilitates ultrasensitive variant detection down to 0.1% by using integrated unique molecular indices (UMIs) and high-fidelity chemistry from cfDNA in biofluids within 8 hours and is coupled with an error-correction data analysis. The required amount of template for a single QIAseq Targeted cfDNA Ultra sequencing reaction ranges from 5 to 80 ng of cfDNA.

The NGS of DNA is a powerful tool for the detection of genetic variations, including somatic mutations, single-nucleotide polymorphisms, copy-number variation, and small insertions and deletions. Target-enrichment technology enhances DNA NGS by enabling users to sequence specific regions of interest — instead of the entire genome — which effectively increases sequencing depth and sample throughput while minimizing cost. Many commercially available target enrichment, library preparation and sequencing methods use regular DNA polymerase and amplification processes that introduce substantial bias and artifacts. This results in artifactual errors that greatly limit the detection of true low-frequency variants below 0.5% in heterogeneous samples, such as cfDNA. The QIAseq Targeted cfDNA Ultra overcomes these biases and artifacts with the use of a highly optimized, high-fidelity reaction chemistry and by incorporating UMIs into a single gene-specific, primer-based targeted enrichment process.

The QIAseq Targeted cfDNA Ultra has been optimized with a specially formulated enrichment chemistry to achieve highly efficient enrichment on both regular and GC-rich regions at high multiplex levels. In addition, the QIAseq Targeted cfDNA Ultra library can be sequenced with Illumina default sequencing primers and is compatible with most medium- and high-throughput Illumina sequencers.

Our data analysis tools have been developed to perform all steps necessary to generate a DNA sequence variant report from NGS data. Collectively, the QIAseq Targeted cfDNA Ultra is a Sample to Insight solution for precision variant detection of targeted genomic regions using NGS (Figure 1).




Figure 1. Overview of the Sample to Insight NGS workflow for cell-free DNA with the QIAseq Targeted cfDNA Ultra.

This complete Sample to Insight procedure begins with DNA extraction followed by library construction and target enrichment with the QIAseq Targeted cfDNA Ultra. After NGS, data analysis is performed using the QIAGEN CLC Genomics Workbench. Ultimately, detected variants can be interpreted with the QIAGEN Clinical Insight Interpret for clinical relevant information.

Principle and Procedure

The QIAseq Targeted cfDNA Ultra Kits are provided as single-tube primer mixes, with up to 20,000 primers per panel. The QIAseq Targeted cfDNA Ultra Kit is designed to enrich selected genes and regions using 5–80 ng cfDNA (Figure 2). Lower input amounts are possible; however, this will lead to fewer sequenced UMIs and reduced variant detection sensitivity.




Figure 2. QIAseq Targeted cfDNA Ultra workflow.

End Preparation

cfDNA samples are first end repaired and A tailed within a single, controlled multienzyme reaction. The prepared DNA fragments are then ligated with a sequencing platform-specific adapter containing UMIs.

UMI assignment

Prior to target enrichment and library amplification, each original DNA molecule is assigned a unique sequence or index, commonly referred to as a UMI. This assignment is accomplished by ligating fragmented DNA with an adapter containing a 12-base fully random sequence (i.e., the UMI). Statistically, this process provides 4^{12} possible indices per adapter, and each DNA molecule in the sample receives a unique UMI sequence.

Target enrichment and final library construction

Target enrichment is performed post-UMI assignment to ensure that DNA molecules containing UMIs are sufficiently enriched in the sequenced library. For enrichment, ligated DNA molecules are subject to several cycles of targeted PCR using one region-specific primer and one universal primer complementary to the adapter. A universal PCR is ultimately carried out to amplify the library and add platform-specific adapter sequences and additional sample indices.

NGS adapter and index technologies

The QIAseq Targeted cfDNA Ultra Panels use unique dual index (UDI) primers for sample indexing. The UDI primers significantly reduce the risk of index-bleeding issues associated with different Illumina sequencing instruments, as well as reducing the impact of low-level contamination during oligo synthesis, and kit manufacturing, as well as carry over on the Illumina sequencing instrument itself. Hence, each sample will be assigned two unique indices to mitigate errors from image analysis, sequencing error, demultiplexing, and oligo synthesis contamination, reducing reads mis-assignment to wrong samples.

Next-generation sequencing

The QIAseq Targeted cfDNA Ultra is compatible with most medium- and high-throughput sequencers including Illumina NGS systems (MiniSeq®, MiSeq®, NextSeq® 500/550, NextSeq 1000/2000, NovaSeq® 6000, and NovaSeq® X).

Principle of variant detection with UMIs

The principle of variant detection with UMIs is described in Figure 3. Indexed molecules may be amplified unevenly across the target regions due to intrinsic noise and sequence-dependent bias. However, target region coverage can be better achieved by counting the number of UMIs rather than counting the number of total reads for each region. Sequence reads having different UMIs represent different original molecules, while sequence reads having the same UMIs are the result of PCR duplication from one original molecule. Errors from PCR amplification and from the sequencing process may also be present in final reads that lead to false positive variants in sequencing results. These artifactual variants can be greatly reduced by calling variants across all reads within a unique UMI instead of picking up variants at the original read level.




Figure 3. Principle of variant detection with UMIs. Each original molecule is tagged by a UMI. True variants are mutations present in the majority of reads within a UMI, while false positives are mutations present in only one or a few reads within a UMI.

Data analysis

The data for the QIAseq Targeted cfDNA Ultra can be analyzed using the Biomedical Genomics Analysis plugin to the QIAGEN CLC Genomics Workbench. The plugin provides workflows and tools for all steps from the initial data processing and quality assurance through data analyses, annotation, and reporting. A detailed guide to UMI-directed variant detection in CLC Genomics Workbench can be found in *Biomedical Genomics Analysis Plugin User Manual*.

All detected variants can be further interpreted using QIAGEN’s Clinical Insight (QCI®) Interpret.

Equipment and Reagents to Be Supplied by User

When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles. For more information, consult the appropriate safety data sheets (SDSs), available from the product supplier.

In addition to the QIAseq Targeted cfDNA Ultra and Index Kit, the following are required:

- DynaMag™-96 Side Magnet (Thermo Fisher Scientific Inc., cat. no. 12331D)
- QIAxpert® Instrument (cat. no. 9002340), Qubit™ 4 Fluorometer (Thermo Fisher Scientific Inc., cat. no. Q33238), or Qubit Flex Fluorometer (Thermo Fisher Scientific Inc., cat. no. Q33327)
- Ethanol, 80% (made fresh each day before use)*
- Nuclease-free pipette tips and tubes
- 1.5 mL DNA LoBind® tubes (Eppendorf®, cat. no. 022431021)
- PCR tubes and caps (0.2 mL individual PCR tubes [VWR, cat. no. 20170-012], tube and caps 8 strips [VWR, cat. no. 93001 118]) or 96-well PCR plates and caps
- Ice
- Microcentrifuge
- Thermal cycler
- Multichannel pipettes
- Single-channel pipettes
- QIAxcel® Connect System (cat. no. 9003110), QIAxcel Advanced System (cat. no. 9001941), or Agilent® 2100 Bioanalyzer® (Agilent, cat. no. G2939AAR)
- QIAxcel DNA High Sensitivity Kit (cat. No. 929012), QIAxcel DNA High Resolution Kit (QIAGEN, cat. no. 929002), or Agilent High Sensitivity DNA Kit (Agilent, cat. no. 5067-4626)

*Do not use denatured alcohol, which contains other substances, such as methanol or methylethylketone.

Important Notes

For optimal results, all DNA samples should demonstrate consistent quality according to the following criteria.

DNA isolation and quality check

The most important prerequisite for DNA sequence analysis is consistent, high-quality DNA from every experimental sample. Therefore, sample handling and DNA isolation procedures are critical to the success of the experiment. Residual traces of proteins, salts or other contaminants may either degrade the DNA or decrease the efficiency of, if not block completely, the enzymatic activity necessary for optimal targeted enrichment. Sample purity can be checked with the QIAxpert System.

The QIAGEN kits listed in Table 1 are recommended for the preparation of cfDNA from serum or plasma samples. For best results, all DNA samples should be resuspended in DNase-free water or, alternatively, in DNase-free 10 mM Tris buffer, pH 8.0.

Important: Do not use DEPC-treated water.

Note: If genomic DNA must be isolated from biological samples for which kits are not available, please contact QIAGEN Technical Support for suggestions.

Table 1. Recommended kits for purification of cfDNA

Kit	Starting material	Cat. no.
QIAamp® Circulating Nucleic Acid Kit	Up to 5 mL of serum or plasma	55114
QIAamp MinElute ccfDNA Mini Kit	Up to 4 mL of serum or plasma	55204
QIAamp MinElute ccfDNA Midi Kit	Up to 10 mL of serum or plasma	55284
EZ1&2 ccfDNA Kit	Automated extraction from up to 8 mL serum or plasma	954854

For best results, all DNA samples should also demonstrate consistent quality according to the following criteria.

cfDNA quality

cfDNA quality can be checked using the QIAxcel System, Agilent Bioanalyzer or Agilent TapeStation. Good cfDNA samples should have a distinct cfDNA peak and a minimal amount of cellular genomic DNA that is of high-molecular weight.

cfDNA quantification

If the cfDNA is free of high-molecular weight cellular DNA, its concentration can be measured by fluorescence-based quantification, such as with the Qubit Fluorometer. If a significant amount of cellular DNA is present in cfDNA sample, it is recommended to measure the amount based on the cfDNA peak with either the QIAxcel System or the Bioanalyzer.

DNA input amount and sequencing depth

The number of UMIs captured from the original DNA sample is related to the DNA input amount and sequencing depth. Adequate sequencing of captured UMIs is necessary for UMI-based variant detection and requires relatively deep sequencing coverage. Table 2 provides guidance on variant detection with cfDNA amounts at different depths of coverage.

Table 2. Suggested cfDNA input amount and sequencing depth for variant detection*

Variant frequency (%)	Input (ng)	Read pairs per UMI	Mean read
0.5	30	4	24,000
0.2	30	4	24,000
0.1	60	4	48,000

* Variant detection is based on 80–90% sensitivity on the entire panel region of the QIAseq Targeted cfDNA Ultra.

Variant detection

The number of UMIs sequenced directly impacts the sensitivity of variant detection. Therefore, low-frequency mutation detection usually requires more DNA input and sequencing at deeper coverage (i.e., more reads per UMI) to generate a sufficient amount of UMIs.

Sequencing capacity and sample multiplex level

Sample multiplexing is one of the most important NGS tools for increasing throughput and reducing costs. It works by combining multiple samples to be processed together in a single sequencing run; as a consequence, sequencing reads need to be demultiplexed by reassigning each single read to its original source library. This is facilitated by the integration of sample index sequences into the individual library molecules.

The QIAseq Targeted cfDNA Ultra Panels are paired with a fully compatible UDI sample-indexing solution. Each QIAseq Targeted cfDNA Ultra UDI Kit includes one of the following:

- QIAseq Targeted cfDNA Ultra UDI Set A (96): DNA Ultra UDI Set A index primer plate, QUDI-96AA
- QIAseq Targeted cfDNA Ultra UDI Set B (96): DNA Ultra UDI Set B index primer plate, QUID-96BA
- QIAseq Targeted cfDNA Ultra UDI Set C (96): DNA Ultra UDI Set C index primer plate, QUDI-96CA
- QIAseq Targeted cfDNA Ultra UDI Set D (96): DNA Ultra UDI Set D index primer plate, QUDI-96DA
- QIAseq Targeted cfDNA Ultra UDI (12): DNA Ultra UDI 12 Index primer plate, QUDI-12A

The QIAseq Ultra UDI Kits use a fixed combination of two unique bar code motives per sample-index primer pair. Therefore, each single-index motive is only used once on any UDI index primer.

Use of UDI indexes effectively mitigates the risk of read misassignment due to index hopping, which is important for detecting variants at very low frequency. By filtering misassigned reads during the demultiplexing of individual samples, highly accurate output data is generated.

To multiplex more than 96 libraries in a single sequencing run, combine kits with different sets, QIAseq Targeted cfDNA Ultra UDI Set A, B, C or D. For example, combining the unique dual QIAseq Targeted cfDNA Ultra UDI Set A and B (96) kits will allow the generation of 192 libraries with different unique dual sample indexes for 192-plex sequencing.

Sample multiplexing level is determined by the size of the panel, required depth of coverage and sequencing platform read capacity. General guidelines are provided for the number of samples that can be multiplexed with different sequencing platforms, based on panel size and read depth (Table 3 and Table 4). Read depth can be fine-tuned after the first run.

Table 3. Number of multiplexed samples based on panel size with 25,000x mean coverage*

Instrument	Version	Capacity (paired-ends reads)	500 primers	1000 primers	2500 primers	5000 primers	10,000 primers
iSeq	v2 reagents	8 M	N/A	N/A	N/A	N/A	N/A
MiniSeq	Mid output	16 M	1	N/A	N/A	N/A	N/A
MiniSeq	High output	50 M	4	2	N/A	N/A	N/A
MiSeq	v2 reagents	30 M	2	1	N/A	N/A	N/A
MiSeq	v3 reagents	50 M	4	2	N/A	N/A	N/A
NextSeq 500	Mid output	260 M	20	10	4	2	1
NextSeq 500	High output	800 M	64	32	12	6	3
NextSeq 1000/2000	P1 flow cell	200 M	16	8	3	1	N/A
NextSeq 1000/2000	P2 flow cell	800 M	64	32	12	6	3
NextSeq 2000	P3 flow cell	2.4 B	192	96	38	19	9
NextSeq 2000	P4 flow cell	3.6B	288	144	57	28	13
NovaSeq 6000	SP (2 lanes per flow cell)	1.6 B	128	64	25	12	6
NovaSeq 6000	S1 (2 lanes per flow cell)	3.2 B	256	128	51	25	12
NovaSeq 6000	S2 (2 lanes per flow cell)	8.2 B	656	328	131	65	32
NovaSeq 6000	S4 (4 lanes per flow cell)	20 B	1600	800	320	160	80
NovaSeq X	1.5B flow cell	3.2 B	256	128	51	25	12
NovaSeq X	10B flow cell	20 B	1600	800	320	160	80
NovaSeq X	25B flow cell	52 B	4160	2080	832	416	208

* Based on 2 x 149 bp paired-end reads.

N/A: Not applicable no samples can be run.

Table 4. Number of multiplexed samples based on panel size with 50,000x mean coverage*

Instrument	Version	Capacity (paired-ends reads)	500 primers	1000 primers	2500 primers	5000 primers	10,000 primers
iSeq	v2 reagent	8 M	N/A	N/A	N/A	N/A	N/A
MiniSeq	Mid output	16 M	N/A	N/A	N/A	N/A	N/A
MiniSeq	High output	50 M	2	1	N/A	N/A	N/A
MiSeq	v2 reagents	30 M	1	N/A	N/A	N/A	N/A
MiSeq	v3 reagents	50 M	2	1	N/A	N/A	N/A
NextSeq 500	Mid output	260 M	10	5	2	1	N/A
NextSeq 500	High output	800 M	32	16	6	3	1
NextSeq 1000/2000	P1 flow cell	200 M	8	4	1	N/A	N/A
NextSeq 1000/2000	P2 flow cell	800 M	32	16	6	3	1
NextSeq 2000	P3 flow cell	2.4 B	96	48	19	9	4
NextSeq 2000	P4 flow cell	3.6B	144	72	28	13	6
NovaSeq 6000	SP (2 lanes per flow cell)	1.6 B	64	32	12	6	3
NovaSeq 6000	S1 (2 lanes per flow cell)	3.2 B	128	64	25	12	6
NovaSeq 6000	S2 (2 lanes per flow cell)	8.2 B	328	164	65	32	16
NovaSeq 6000	S4 (4 lanes per flow cell)	20 B	800	400	160	80	40
NovaSeq X	1.5B flow cell	3.2 B	128	64	25	12	6
NovaSeq X	10B flow cell	20 B	800	400	160	80	40
NovaSeq X	25B flow cell	52 B	2080	1040	416	208	104

* Based on 2 x 149 bp paired-end reads for.

N/A: Not applicable, no samples can be run.

NGS read-length recommendations

When using Illumina NGS systems, the QIAseq Targeted cfDNA Ultra UDI libraries require 149 bp paired-end reads and dual 10 bp indices.

Protocol: DNA End Prepare

Important points before starting

- This protocol covers all procedures required for the preparation of libraries for Illumina sequencers from cfDNA.
- Before setting up the reaction, it is critical to accurately determine the amount of the input cfDNA. When cfDNA is contaminated with cellular genomic DNA, measuring the concentration of the cfDNA peak with the QIAxcel, Bioanalyzer or similar methods is much more accurate than measuring the concentration of the whole sample.
- Reaction and cleanup procedures can be performed in either PCR tubes or 96-well plates.
- Set up reactions on ice.
- Do not vortex any reagents or reactions.

Procedure: end prepare

1. Thaw nucleic acid samples on ice. Gently mix, briefly centrifuge to collect residual liquid from the sides of the tubes, and then return to ice.
2. Prepare the reagents required for preparing the ends.
 - a. Thaw 10x ERA Buffer on ice or, if required, at room temperature (15–25°C). Immediately upon thawing, place the buffer on ice. Keep ERA Enzyme on ice.
 - b. Mix all reagents by flicking the tube, and centrifuge briefly.
3. On ice, prepare the end prepare mix according to Table 5. Briefly centrifuge, mix by pipetting up and down at least 12 times with pipetting volume close to 25 μ L and briefly centrifuge again.

Note: In general, increasing the amount of DNA input will improve variant detection sensitivity.

Important: Keep the reaction tubes or plates on ice during the entire reaction setup.

Note: If setting up more than one reaction, prepare a volume of master mix 10% greater than what is required for the total number of reactions.

Table 5. Reaction mix for end prepare

Component	Volume per reaction
cfDNA*	Variable
10x ERA Buffer	2.5 µL
ERA Enzyme	5 µL
Nuclease-free Water	Variable
Total	25 µL

* For cfDNA, add 5–80 ng. We recommend 30 ng or more for detection of variants below 0.5%.

4. Program the thermal cycler according to Table 6. Use the instrument’s heated lid.
5. Before adding the tubes or plate to a thermal cycler, start the program. When the thermal cycler reaches 4°C, pause the program.

Important: The thermal cycler must be prechilled and paused at 4°C.

Important: Do not use the heated lid during the 20°C incubation step. Alternatively, the lid temperature can be set at 65°C.

6. Transfer the tubes or plate prepared in step 3 to the prechilled thermal cycler and resume the cycling program.

Table 6. Incubation conditions for fragmentation and end prepare

Step	Incubation temperature (°C)	Incubation time
1	4	1 min
2	20	15 min
3	65	15 min
4	4	Hold

7. Upon completion, allow the thermal cycler to return to 4°C.
8. Place the samples on ice and immediately proceed to "Protocol: Adapter Ligation" on the next page.

Protocol: Adapter Ligation

Important points before starting

- The 25 μL product from "Protocol: DNA End Prepare" on page 23 is the starting material for this protocol.
- Set up reactions on ice.
- Do not vortex any reagents or reactions.
- The QIAseq Beads are used for ligation reaction cleanup. There is no need to bring beads to room temperature before use.
- Prepare fresh 80% ethanol each day before use.
- Ensure that the QIAseq Beads are thoroughly mixed at all times. This necessitates working quickly and resuspending the beads immediately before use. If a delay in the protocol occurs, simply vortex the beads.

Procedure: adapter ligation

1. Prepare the reagents required for the DNA ligation.
 - a. Thaw AdP-DP-Phased Adapter on ice. Thaw UPH Ligation Buffer, 2.5x, on ice or at room temperature. Immediately upon thawing, place the buffer on ice. Keep DNA Ligase on ice.
 - b. Mix by flicking the tube, and then centrifuge briefly.
2. Prepare the adapter ligation mix according to Table 7. Briefly centrifuge, mix by pipetting up and down at least 12 times with pipetting volume close to 50 μL and briefly centrifuge again.

Important: The AdP-DP-Phased Adapter does not contain a sample index; hence, one single adapter is used for all samples.

Note: If setting up more than one reaction, prepare a volume of master mix 10% greater than what is required for the total number of reactions.

Table 7. Reaction mix for adapter ligation

Component	Volume/reaction, <10 ng cfDNA	Volume/reaction, 10–50 ng cfDNA	Volume/reaction, >50 ng cfDNA
End prepare reaction (already in tube)	25 µL	25 µL	25 µL
UPH Ligation Buffer, 2.5x	20 µL	20 µL	20 µL
AdP-DP-Phased Adapter	0.5 µL	1 µL	1.5 µL
DNA Ligase	5 µL	5 µL	5 µL
Total	50.5 µL	51 µL	51.5 µL

3. Program the thermal cycler according to Table 8.

Important: Do not use the heated lid during the 20°C incubation step. Alternatively, the lid temperature can be set at 65°C.

- Before adding the tubes or plate to a thermal cycler, start the program. When the thermal cycler reaches 4°C, pause the program.
- Transfer the tubes or plate prepared in step 2 to the prechilled thermal cycler and resume the cycling program.

Table 8. Incubation conditions for DNA ligation

Step	Incubation temperature (°C)	Incubation time
1	4	1 min
2	20	15 min
3	65	15 min
4	4	Hold

6. Upon completion, allow the thermal cycler to return to 4°C.
7. Place the samples on ice and immediately proceed to "Procedure: ligation cleanup" (below). Alternatively, the samples can be stored at -20°C in a constant-temperature freezer for up to 3 days.

Procedure: ligation cleanup

8. Add 40 µL Nuclease-free Water to bring each sample to approximately 90 µL.
9. Add 108 µL QIAseq Beads. Mix well by vortexing or pipetting up and down at least 15 times with pipetting volume close to 190 µL. Incubate for 5 min at room temperature.
10. Place the tubes or plate on a magnetic rack for 10 min (tube) or 15 min (plate). Once the solution has cleared, with the beads still on the magnetic stand, carefully remove and discard the supernatant.

Important: Do not discard the beads as they contain the DNA of interest.

11. With the beads still on the magnetic stand, add 80 µL water to the beads, then add 80 µL QIAseq bead binding buffer.
12. Take the tubes or plate off the magnetic stand, mix well by vortexing or pipetting up and down at least 15 times with pipetting volume close to 160 µL.
13. Return the tubes or plate to the magnetic rack for 5 min. Once the solution has cleared, with the beads still on the magnetic stand, carefully remove and discard the supernatant.

Important: Do not discard the beads as they contain the DNA of interest.

14. With the beads still on the magnetic stand, add 200 µL 80% ethanol and leave beads in ethanol for 30 seconds. Carefully remove and discard the ethanol.
15. Repeat the ethanol wash.

Important: Completely remove all traces of the ethanol wash after this second wash. Remove the ethanol with a 200 μ L pipette first, and then use a 10 μ L pipette to remove any residual ethanol.

16. With the beads still on the magnetic stand, air dry at room temperature for at least 15 min.

Note: Visually inspect that the pellet is completely dry. Ethanol carryover to the target enrichment PCR step will significantly impact enrichment PCR efficiency. Over drying beads will not affect DNA elution.

17. Remove the beads from the magnetic stand, and elute the DNA from the beads by adding 12 μ L Nuclease-free Water. Mix well by pipetting.
18. Return the tubes or plate to the magnetic rack until the solution has cleared.
19. Transfer 10 μ L of the supernatant to clean tubes or plate.

Important: Make sure no beads remain in the elution. Bead carryover to the target enrichment PCR step will significantly impact enrichment PCR efficiency.

20. Proceed with "Protocol: Target Enrichment" (next page). Alternatively, the samples can be stored at -20°C in a constant-temperature freezer for up to 3 days.

Protocol: Target Enrichment

Important points before starting

- The 10 μ L cleaned ligation reaction from "Protocol: Adapter Ligation" on page 26 is the starting material for this protocol.
- Set up reactions on ice.
- Do not vortex any reagents or reactions.
- The QIAseq Beads are used for target enrichment PCR cleanup. There is no need to bring beads to room temperature before use.
- Prepare fresh 80% ethanol each day before use.
- Ensure that the QIAseq Beads are thoroughly mixed at all times. This necessitates working quickly and resuspending the beads immediately before use. If a delay in the protocol occurs, simply vortex the beads.

Procedure: target enrichment

1. Prepare the reagents required for the Target Enrichment PCR (TEPCR).
 - a. Thaw HiFi Ultra Buffer, 5x; QIAseq Targeted cfDNA Ultra Panel; and SmP-IL5 TEPCR-F Primer on ice or at room temperature. Immediately upon thawing, place the buffer on ice. Keep HiFi Ultra Polymerase on ice.
 - b. Mix all reagents by flicking the tube, and then centrifuge briefly.
2. Prepare the target enrichment mix according to Table 9. Briefly centrifuge, mix by pipetting up and down at least 12 times with pipetting volume close to 20 μ L, and briefly centrifuge again.

Note: If setting up more than one reaction, prepare a volume of master mix 10% greater than what is required for the total number of reactions.

Important: Do not add more HiFi Ultra Polymerase than it is indicated in Table 9. Adding more could result in low enrichment specificity.

Table 9. Reaction mix for Target Enrichment

Component	Volume per reaction			
	≤1100 primers	1101–5000 primers	5001–10,000 primers	>10,000 primers
Cleaned Adapter-ligated DNA (already in tube) from "Procedure: ligation cleanup"	10 µL	10 µL	10 µL	10 µL
HiFi Ultra Buffer, 5x	4 µL	4 µL	4 µL	4 µL
QIAseq Targeted cfDNA Ultra Panel	5 µL	5 µL	5 µL	5 µL
SmP-IL5 TEPCR-F Primer	2 µL	2 µL	2 µL	2 µL
HiFi Ultra Polymerase*	–	0.25 µL	0.5 µL	0.75 µL
HiFi Ultra Polymerase (10x diluted)* †	0.6 µL	–	–	–
Total	21.6 µL	21.25 µL	21.5 µL	21.75 µL

* It is recommended to make a master mix instead of pipetting small amount of HiFi Ultra Polymerase enzyme.

† Immediately before setting up the reaction, dilute HiFi Ultra Polymerase 10 folds with Nuclease-free Water. There is enough amount of HiFi Ultra Polymerase in the kit for the dilution.

3. Program a thermal cycler using the cycling condition in Table 10 (panel with ≤1100 primers/tube), Table 11 (panel with 1101–10,000 primers/tube) or Table 12 (panel with ≥10,001 primers/tube). Use the instrument's heated lid.

Table 10. Cycling conditions for target enrichment if number of primers ≤ 1100 /tube

Step	Time (<400 primers)	Time (400-1100 primers)	Temperature (°C)
Initial denaturation	2 min	2 min	98
8 cycles	20 s	20 s	98
	5 min	10 min	70
1 cycle	3 min	3 min	72
Hold	∞	∞	4

Table 11. Cycling conditions for target enrichment if number of primers 1101–10000/tube

Step	Time	Temperature (°C)
Initial denaturation	2 min	98
8 cycles	20 s	98
	12 min	69
1 cycle	3 min	72
Hold	∞	4

Table 12. Cycling conditions for target enrichment if number of primers ≥ 10001 /tube

Step	Time	Temperature (°C)
Initial denaturation	2 min	98
8 cycles	20 s	98
	10 min	68
1 cycle	3 min	72
Hold	∞	4

- Place the target enrichment reaction in the thermal cycler and start the run.
- After the reaction is complete, place the reactions on ice and proceed to "Procedure:

TEPCR cleanup" (next page). Alternatively, the samples can be stored at -20°C in a constant-temperature freezer for up to 3 days.

Procedure: TEP CR cleanup

1. Add 70 μL Nuclease-free Water to bring each sample to approximately 90 μL .
2. Add 108 μL QIAseq Beads. Mix well by vortexing or pipetting up and down at least 15 times with pipetting volume close to 190 μL . Incubate for 5 min at room temperature.
3. Place the tubes or plate on a magnetic rack for 5 min (tube) or 10 min (plate). After the solution has cleared, with the beads still on the magnetic stand, carefully remove and discard the supernatant.

Important: Do not discard the beads as they contain the DNA of interest.

4. With the beads still on the magnetic stand, add 200 μL 80% ethanol and leave beads in ethanol for 30 s. Carefully remove and discard the ethanol.
5. Repeat the ethanol wash.

Important: Completely remove all traces of the ethanol wash after this second wash. Remove the ethanol with a 200 μL pipette first, and then use a 10 μL pipette to remove any residual ethanol.

6. With the beads still on the magnetic stand, air dry at room temperature for at least 10 min.

Note: Visually inspect that the pellet is completely dry. Ethanol carryover to the next universal PCR step will significantly impact PCR efficiency. Over drying beads will not affect DNA elution.

7. Remove the beads from the magnetic stand, and elute the DNA from the beads by adding 52 μL Nuclease-free Water. Mix well by pipetting.
8. Return the tubes or plate to the magnetic rack until the solution has cleared.

9. Transfer 50 μL of the supernatant to clean tubes or plate. Transfer 5 μL to a clean tube or plate for universal PCR setup.

Important: Make sure no beads transfer to the elution. Bead carryover to the universal PCR step will significantly impact PCR efficiency.

Note: Only 5 μL of cleaned TEPCR volume will be used in the universal PCR. This volume is sufficient, as the DNA molecules were amplified in multiple cycles during TEPCR. The rest of the TEPCR volume can be stored at -20°C if needed.

10. Proceed with "Protocol: Universal PCR". Alternatively, the samples can be stored at -20°C in a constant-temperature freezer for up to 3 days.

Protocol: Universal PCR

Important points before starting

- The 5 μ L cleaned product from "Protocol: Target Enrichment" on page 30 is the starting material for this protocol.
- Set up reactions on ice.
- Do not vortex any reagents or reactions.
- The QIAseq Targeted cfDNA Ultra UDI (12) or QIAseq Targeted cfDNA Ultra UDI Set A (96), Set B (96), Set C (96) or Set D (96) is used for sample indexing.
- The Index Primer Plate contains predispensed index primer pairs and the universal PCR primer and is sealed with pierceable aluminum heat sealing film.
- Puncture the film using standard 200 μ L pipette tips to transfer the appropriate amount of index primer to the tubes or plate for the universal PCR.
- The QIAseq Beads are used for universal PCR cleanup. There is no need to bring beads to room temperature before use.
- Prepare fresh 80% ethanol each day before use.
- Ensure that the QIAseq Beads are thoroughly mixed at all times. This necessitates working quickly and resuspending the beads immediately before use. If a delay in the protocol occurs, simply vortex the beads.

Procedure: universal PCR

1. Prepare the reagents required for the universal PCR.
 - a. Thaw HiFi Ultra Buffer (5x), DNA Ultra UDI 12, and 96-Index Plate on ice or at room temperature. Immediately upon thawing, place the items on ice. Keep HiFi Ultra

Polymerase on ice.

- b. Mix by flicking the tube, and then centrifuge briefly.
2. Prepare the universal PCR in the tubes or plate containing the cleaned target-enriched from TEPCR cleanup reaction.

Note: If setting up more than one reaction, prepare a volume of master mix 10% greater than what is required for the total number of reactions.

3. For the QIAseq Targeted cfDNA Ultra UDI plates, pierce the foil seal associated with each well that will be used, and transfer 2 μL (each well contains a forward primer and a reverse primer, each with a unique index) to the cleaned target-enriched DNA from "Procedure: TEPCR cleanup" sample tubes or plate according to Table 13. Mix by pipetting up and down at least 12 times with pipetting volume close to 40 μL and briefly centrifuge again.

Important: Only one UDI pair should be used per universal PCR.

Important: The QIAseq Targeted cfDNA Ultra UDI index plates are stable for a maximum of 10 freeze-thaw cycles. If all 96 wells have not been used, cover the used wells with foil and return to the freezer. Do not reuse wells from the QIAseq Targeted cfDNA Ultra UDI index plates once the foil seals have been pierced. Reusing wells would risk significant cross-contamination.

Table 13. Reaction components for universal PCR if using QIAseq Targeted cfDNA Ultra UDI (12) or QIAseq Targeted cfDNA Ultra UDI Set A, B, C, or D (96)

Component	Volume/ reaction (µL)
Cleaned target-enriched DNA from "Procedure: TEPCR cleanup"	5
HiFi Ultra Buffer, 5x	8
Index primers from QIAseq Targeted cfDNA Ultra UDI index plate*	2
HiFi Ultra Polymerase	2
Nuclease-free water	23
Total	40

* Applies to QIAseq Targeted cfDNA Ultra UDI (12) or QIAseq Targeted cfDNA Ultra UDI Set A, B, C, and D (96).

DNA Ultra UDI 12 Index Plate (QUDI-12A) in QIAseq Targeted cfDNA Ultra UDI (12)

	1	2	3	4	5	6	7	8	9	10	11	12
A	LS-QUD1-1 LS-QUD1-9	LS-QUD1-10 LS-QUD1-18	LS-QUD1-11 LS-QUD1-19	LS-QUD1-12 LS-QUD1-20	LS-QUD1-13 LS-QUD1-21	LS-QUD1-14 LS-QUD1-22	LS-QUD1-15 LS-QUD1-23	LS-QUD1-16 LS-QUD1-24	LS-QUD1-17 LS-QUD1-25	LS-QUD1-18 LS-QUD1-26	LS-QUD1-19 LS-QUD1-27	LS-QUD1-20 LS-QUD1-28
B	LS-QUD1-21 LS-QUD1-29	LS-QUD1-22 LS-QUD1-30	LS-QUD1-23 LS-QUD1-31	LS-QUD1-24 LS-QUD1-32	LS-QUD1-25 LS-QUD1-33	LS-QUD1-26 LS-QUD1-34	LS-QUD1-27 LS-QUD1-35	LS-QUD1-28 LS-QUD1-36	LS-QUD1-29 LS-QUD1-37	LS-QUD1-30 LS-QUD1-38	LS-QUD1-31 LS-QUD1-39	LS-QUD1-32 LS-QUD1-40
C	LS-QUD1-33 LS-QUD1-41	LS-QUD1-34 LS-QUD1-42	LS-QUD1-35 LS-QUD1-43	LS-QUD1-36 LS-QUD1-44	LS-QUD1-37 LS-QUD1-45	LS-QUD1-38 LS-QUD1-46	LS-QUD1-39 LS-QUD1-47	LS-QUD1-40 LS-QUD1-48	LS-QUD1-41 LS-QUD1-49	LS-QUD1-42 LS-QUD1-50	LS-QUD1-43 LS-QUD1-51	LS-QUD1-44 LS-QUD1-52
D	LS-QUD1-45 LS-QUD1-53	LS-QUD1-46 LS-QUD1-54	LS-QUD1-47 LS-QUD1-55	LS-QUD1-48 LS-QUD1-56	LS-QUD1-49 LS-QUD1-57	LS-QUD1-50 LS-QUD1-58	LS-QUD1-51 LS-QUD1-59	LS-QUD1-52 LS-QUD1-60	LS-QUD1-53 LS-QUD1-61	LS-QUD1-54 LS-QUD1-62	LS-QUD1-55 LS-QUD1-63	LS-QUD1-56 LS-QUD1-64
E	LS-QUD1-57 LS-QUD1-65	LS-QUD1-58 LS-QUD1-66	LS-QUD1-59 LS-QUD1-67	LS-QUD1-60 LS-QUD1-68	LS-QUD1-61 LS-QUD1-69	LS-QUD1-62 LS-QUD1-70	LS-QUD1-63 LS-QUD1-71	LS-QUD1-64 LS-QUD1-72	LS-QUD1-65 LS-QUD1-73	LS-QUD1-66 LS-QUD1-74	LS-QUD1-67 LS-QUD1-75	LS-QUD1-68 LS-QUD1-76
F	LS-QUD1-69 LS-QUD1-77	LS-QUD1-70 LS-QUD1-78	LS-QUD1-71 LS-QUD1-79	LS-QUD1-72 LS-QUD1-80	LS-QUD1-73 LS-QUD1-81	LS-QUD1-74 LS-QUD1-82	LS-QUD1-75 LS-QUD1-83	LS-QUD1-76 LS-QUD1-84	LS-QUD1-77 LS-QUD1-85	LS-QUD1-78 LS-QUD1-86	LS-QUD1-79 LS-QUD1-87	LS-QUD1-80 LS-QUD1-88
G	LS-QUD1-81 LS-QUD1-89	LS-QUD1-82 LS-QUD1-90	LS-QUD1-83 LS-QUD1-91	LS-QUD1-84 LS-QUD1-92	LS-QUD1-85 LS-QUD1-93	LS-QUD1-86 LS-QUD1-94	LS-QUD1-87 LS-QUD1-95	LS-QUD1-88 LS-QUD1-96	LS-QUD1-89 LS-QUD1-97	LS-QUD1-90 LS-QUD1-98	LS-QUD1-91 LS-QUD1-99	LS-QUD1-92 LS-QUD1-100
H	LS-QUD1-93 LS-QUD1-101	LS-QUD1-94 LS-QUD1-102	LS-QUD1-95 LS-QUD1-103	LS-QUD1-96 LS-QUD1-104	LS-QUD1-97 LS-QUD1-105	LS-QUD1-98 LS-QUD1-106	LS-QUD1-99 LS-QUD1-107	LS-QUD1-100 LS-QUD1-108	LS-QUD1-101 LS-QUD1-109	LS-QUD1-102 LS-QUD1-110	LS-QUD1-103 LS-QUD1-111	LS-QUD1-104 LS-QUD1-112

DNA Ultra UDI Set A Plate (QUDI-96AA) in QIAseq Targeted cfDNA Ultra UDI Set A (96)

	1	2	3	4	5	6	7	8	9	10	11	12
A	LS-QUD1-101 LS-QUD1-109	LS-QUD1-110 LS-QUD1-118	LS-QUD1-111 LS-QUD1-119	LS-QUD1-112 LS-QUD1-120	LS-QUD1-113 LS-QUD1-121	LS-QUD1-114 LS-QUD1-122	LS-QUD1-115 LS-QUD1-123	LS-QUD1-116 LS-QUD1-124	LS-QUD1-117 LS-QUD1-125	LS-QUD1-118 LS-QUD1-126	LS-QUD1-119 LS-QUD1-127	LS-QUD1-120 LS-QUD1-128
B	LS-QUD1-121 LS-QUD1-129	LS-QUD1-122 LS-QUD1-130	LS-QUD1-123 LS-QUD1-131	LS-QUD1-124 LS-QUD1-132	LS-QUD1-125 LS-QUD1-133	LS-QUD1-126 LS-QUD1-134	LS-QUD1-127 LS-QUD1-135	LS-QUD1-128 LS-QUD1-136	LS-QUD1-129 LS-QUD1-137	LS-QUD1-130 LS-QUD1-138	LS-QUD1-131 LS-QUD1-139	LS-QUD1-132 LS-QUD1-140
C	LS-QUD1-133 LS-QUD1-141	LS-QUD1-134 LS-QUD1-142	LS-QUD1-135 LS-QUD1-143	LS-QUD1-136 LS-QUD1-144	LS-QUD1-137 LS-QUD1-145	LS-QUD1-138 LS-QUD1-146	LS-QUD1-139 LS-QUD1-147	LS-QUD1-140 LS-QUD1-148	LS-QUD1-141 LS-QUD1-149	LS-QUD1-142 LS-QUD1-150	LS-QUD1-143 LS-QUD1-151	LS-QUD1-144 LS-QUD1-152
D	LS-QUD1-145 LS-QUD1-153	LS-QUD1-146 LS-QUD1-154	LS-QUD1-147 LS-QUD1-155	LS-QUD1-148 LS-QUD1-156	LS-QUD1-149 LS-QUD1-157	LS-QUD1-150 LS-QUD1-158	LS-QUD1-151 LS-QUD1-159	LS-QUD1-152 LS-QUD1-160	LS-QUD1-153 LS-QUD1-161	LS-QUD1-154 LS-QUD1-162	LS-QUD1-155 LS-QUD1-163	LS-QUD1-156 LS-QUD1-164
E	LS-QUD1-157 LS-QUD1-165	LS-QUD1-158 LS-QUD1-166	LS-QUD1-159 LS-QUD1-167	LS-QUD1-160 LS-QUD1-168	LS-QUD1-161 LS-QUD1-169	LS-QUD1-162 LS-QUD1-170	LS-QUD1-163 LS-QUD1-171	LS-QUD1-164 LS-QUD1-172	LS-QUD1-165 LS-QUD1-173	LS-QUD1-166 LS-QUD1-174	LS-QUD1-167 LS-QUD1-175	LS-QUD1-168 LS-QUD1-176
F	LS-QUD1-169 LS-QUD1-177	LS-QUD1-170 LS-QUD1-178	LS-QUD1-171 LS-QUD1-179	LS-QUD1-172 LS-QUD1-180	LS-QUD1-173 LS-QUD1-181	LS-QUD1-174 LS-QUD1-182	LS-QUD1-175 LS-QUD1-183	LS-QUD1-176 LS-QUD1-184	LS-QUD1-177 LS-QUD1-185	LS-QUD1-178 LS-QUD1-186	LS-QUD1-179 LS-QUD1-187	LS-QUD1-180 LS-QUD1-188
G	LS-QUD1-181 LS-QUD1-189	LS-QUD1-182 LS-QUD1-190	LS-QUD1-183 LS-QUD1-191	LS-QUD1-184 LS-QUD1-192	LS-QUD1-185 LS-QUD1-193	LS-QUD1-186 LS-QUD1-194	LS-QUD1-187 LS-QUD1-195	LS-QUD1-188 LS-QUD1-196	LS-QUD1-189 LS-QUD1-197	LS-QUD1-190 LS-QUD1-198	LS-QUD1-191 LS-QUD1-199	LS-QUD1-192 LS-QUD1-200
H	LS-QUD1-193 LS-QUD1-201	LS-QUD1-194 LS-QUD1-202	LS-QUD1-195 LS-QUD1-203	LS-QUD1-196 LS-QUD1-204	LS-QUD1-197 LS-QUD1-205	LS-QUD1-198 LS-QUD1-206	LS-QUD1-199 LS-QUD1-207	LS-QUD1-200 LS-QUD1-208	LS-QUD1-201 LS-QUD1-209	LS-QUD1-202 LS-QUD1-210	LS-QUD1-203 LS-QUD1-211	LS-QUD1-204 LS-QUD1-212

DNA Ultra UDI Set B Plate (QUDI-96BA) in QIAseq Targeted cfDNA Ultra UDI Set B (96)

	1	2	3	4	5	6	7	8	9	10	11	12
A	LS-QUD1-201 LS-QUD1-209	LS-QUD1-210 LS-QUD1-218	LS-QUD1-211 LS-QUD1-219	LS-QUD1-212 LS-QUD1-220	LS-QUD1-213 LS-QUD1-221	LS-QUD1-214 LS-QUD1-222	LS-QUD1-215 LS-QUD1-223	LS-QUD1-216 LS-QUD1-224	LS-QUD1-217 LS-QUD1-225	LS-QUD1-218 LS-QUD1-226	LS-QUD1-219 LS-QUD1-227	LS-QUD1-220 LS-QUD1-228
B	LS-QUD1-221 LS-QUD1-229	LS-QUD1-222 LS-QUD1-230	LS-QUD1-223 LS-QUD1-231	LS-QUD1-224 LS-QUD1-232	LS-QUD1-225 LS-QUD1-233	LS-QUD1-226 LS-QUD1-234	LS-QUD1-227 LS-QUD1-235	LS-QUD1-228 LS-QUD1-236	LS-QUD1-229 LS-QUD1-237	LS-QUD1-230 LS-QUD1-238	LS-QUD1-231 LS-QUD1-239	LS-QUD1-232 LS-QUD1-240
C	LS-QUD1-233 LS-QUD1-241	LS-QUD1-234 LS-QUD1-242	LS-QUD1-235 LS-QUD1-243	LS-QUD1-236 LS-QUD1-244	LS-QUD1-237 LS-QUD1-245	LS-QUD1-238 LS-QUD1-246	LS-QUD1-239 LS-QUD1-247	LS-QUD1-240 LS-QUD1-248	LS-QUD1-241 LS-QUD1-249	LS-QUD1-242 LS-QUD1-250	LS-QUD1-243 LS-QUD1-251	LS-QUD1-244 LS-QUD1-252
D	LS-QUD1-245 LS-QUD1-253	LS-QUD1-246 LS-QUD1-254	LS-QUD1-247 LS-QUD1-255	LS-QUD1-248 LS-QUD1-256	LS-QUD1-249 LS-QUD1-257	LS-QUD1-250 LS-QUD1-258	LS-QUD1-251 LS-QUD1-259	LS-QUD1-252 LS-QUD1-260	LS-QUD1-253 LS-QUD1-261	LS-QUD1-254 LS-QUD1-262	LS-QUD1-255 LS-QUD1-263	LS-QUD1-256 LS-QUD1-264
E	LS-QUD1-257 LS-QUD1-265	LS-QUD1-258 LS-QUD1-266	LS-QUD1-259 LS-QUD1-267	LS-QUD1-260 LS-QUD1-268	LS-QUD1-261 LS-QUD1-269	LS-QUD1-262 LS-QUD1-270	LS-QUD1-263 LS-QUD1-271	LS-QUD1-264 LS-QUD1-272	LS-QUD1-265 LS-QUD1-273	LS-QUD1-266 LS-QUD1-274	LS-QUD1-267 LS-QUD1-275	LS-QUD1-268 LS-QUD1-276
F	LS-QUD1-269 LS-QUD1-277	LS-QUD1-270 LS-QUD1-278	LS-QUD1-271 LS-QUD1-279	LS-QUD1-272 LS-QUD1-280	LS-QUD1-273 LS-QUD1-281	LS-QUD1-274 LS-QUD1-282	LS-QUD1-275 LS-QUD1-283	LS-QUD1-276 LS-QUD1-284	LS-QUD1-277 LS-QUD1-285	LS-QUD1-278 LS-QUD1-286	LS-QUD1-279 LS-QUD1-287	LS-QUD1-280 LS-QUD1-288
G	LS-QUD1-281 LS-QUD1-289	LS-QUD1-282 LS-QUD1-290	LS-QUD1-283 LS-QUD1-291	LS-QUD1-284 LS-QUD1-292	LS-QUD1-285 LS-QUD1-293	LS-QUD1-286 LS-QUD1-294	LS-QUD1-287 LS-QUD1-295	LS-QUD1-288 LS-QUD1-296	LS-QUD1-289 LS-QUD1-297	LS-QUD1-290 LS-QUD1-298	LS-QUD1-291 LS-QUD1-299	LS-QUD1-292 LS-QUD1-300
H	LS-QUD1-293 LS-QUD1-301	LS-QUD1-294 LS-QUD1-302	LS-QUD1-295 LS-QUD1-303	LS-QUD1-296 LS-QUD1-304	LS-QUD1-297 LS-QUD1-305	LS-QUD1-298 LS-QUD1-306	LS-QUD1-299 LS-QUD1-307	LS-QUD1-300 LS-QUD1-308	LS-QUD1-301 LS-QUD1-309	LS-QUD1-302 LS-QUD1-310	LS-QUD1-303 LS-QUD1-311	LS-QUD1-304 LS-QUD1-312

Figure 4. Layout of DNA Ultra UDI Index Primer Plate in QIAseq Targeted cfDNA Ultra UDI (12) and QIAseq Targeted cfDNA Ultra UDI Sets A and B (96). Each well contains one pair of predispensed sample index primers plus universal primers for a single reaction in "Protocol: Universal PCR".

DNA Ultra UDI Set C Plate (GUDDI-96CA) in QIAseq Targeted cDNA Ultra UDI Set C (96)

	1	2	3	4	5	6	7	8	9	10	11	12
A	L5-QUD1193	L5-QUD1201	L5-QUD1209	L5-QUD1217	L5-QUD1225	L5-QUD1233	L5-QUD1241	L5-QUD1249	L5-QUD1257	L5-QUD1265	L5-QUD1273	L5-QUD1281
B	L7-QUD1194	L7-QUD1202	L7-QUD1210	L7-QUD1218	L7-QUD1226	L7-QUD1234	L7-QUD1242	L7-QUD1250	L7-QUD1258	L7-QUD1266	L7-QUD1274	L7-QUD1282
C	L5-QUD1196	L5-QUD1204	L5-QUD1212	L5-QUD1220	L5-QUD1228	L5-QUD1236	L5-QUD1244	L5-QUD1252	L5-QUD1260	L5-QUD1268	L5-QUD1276	L5-QUD1284
D	L7-QUD1198	L7-QUD1206	L7-QUD1214	L7-QUD1222	L7-QUD1230	L7-QUD1238	L7-QUD1246	L7-QUD1254	L7-QUD1262	L7-QUD1270	L7-QUD1278	L7-QUD1286
E	L5-QUD1197	L5-QUD1205	L5-QUD1213	L5-QUD1221	L5-QUD1229	L5-QUD1237	L5-QUD1245	L5-QUD1253	L5-QUD1261	L5-QUD1269	L5-QUD1277	L5-QUD1285
F	L7-QUD1199	L7-QUD1207	L7-QUD1215	L7-QUD1223	L7-QUD1231	L7-QUD1239	L7-QUD1247	L7-QUD1255	L7-QUD1263	L7-QUD1271	L7-QUD1279	L7-QUD1287
G	L5-QUD1200	L5-QUD1208	L5-QUD1216	L5-QUD1224	L5-QUD1232	L5-QUD1240	L5-QUD1248	L5-QUD1256	L5-QUD1264	L5-QUD1272	L5-QUD1280	L5-QUD1288
H	L7-QUD1202	L7-QUD1210	L7-QUD1218	L7-QUD1226	L7-QUD1234	L7-QUD1242	L7-QUD1250	L7-QUD1258	L7-QUD1266	L7-QUD1274	L7-QUD1282	L7-QUD1290

DNA Ultra UDI Set D Plate (GUDDI-96DA) in QIAseq Targeted cDNA Ultra UDI Set D (96)

	1	2	3	4	5	6	7	8	9	10	11	12
A	L5-QUD1289	L5-QUD1297	L5-QUD1305	L5-QUD1313	L5-QUD1321	L5-QUD1329	L5-QUD1337	L5-QUD1345	L5-QUD1353	L5-QUD1361	L5-QUD1369	L5-QUD1377
B	L7-QUD1290	L7-QUD1298	L7-QUD1306	L7-QUD1314	L7-QUD1322	L7-QUD1330	L7-QUD1338	L7-QUD1346	L7-QUD1354	L7-QUD1362	L7-QUD1370	L7-QUD1378
C	L5-QUD1291	L5-QUD1299	L5-QUD1307	L5-QUD1315	L5-QUD1323	L5-QUD1331	L5-QUD1339	L5-QUD1347	L5-QUD1355	L5-QUD1363	L5-QUD1371	L5-QUD1379
D	L7-QUD1292	L7-QUD1300	L7-QUD1308	L7-QUD1316	L7-QUD1324	L7-QUD1332	L7-QUD1340	L7-QUD1348	L7-QUD1356	L7-QUD1364	L7-QUD1372	L7-QUD1380
E	L5-QUD1293	L5-QUD1301	L5-QUD1309	L5-QUD1317	L5-QUD1325	L5-QUD1333	L5-QUD1341	L5-QUD1349	L5-QUD1357	L5-QUD1365	L5-QUD1373	L5-QUD1381
F	L7-QUD1294	L7-QUD1302	L7-QUD1310	L7-QUD1318	L7-QUD1326	L7-QUD1334	L7-QUD1342	L7-QUD1350	L7-QUD1358	L7-QUD1366	L7-QUD1374	L7-QUD1382
G	L5-QUD1295	L5-QUD1303	L5-QUD1311	L5-QUD1319	L5-QUD1327	L5-QUD1335	L5-QUD1343	L5-QUD1351	L5-QUD1359	L5-QUD1367	L5-QUD1375	L5-QUD1383
H	L7-QUD1296	L7-QUD1304	L7-QUD1312	L7-QUD1320	L7-QUD1328	L7-QUD1336	L7-QUD1344	L7-QUD1352	L7-QUD1360	L7-QUD1368	L7-QUD1376	L7-QUD1384

Figure 5. Layout of DNA Ultra UDI Index Primer Plate in QIAseq Targeted cDNA Ultra UDI Sets C and D (96). Each well contains one pair of predispersed sample index primers plus universal primers for a single reaction in "Protocol: Universal PCR".

4. Program a thermal cycler using the cycling conditions in Table 14 and Table 15 indicates the cycle number to use dependent on the number of primers in the pool. Use the instrument's heated lid.

Table 14. Incubation conditions for universal PCR

Step	Time	Temperature (°C)
Initial denaturation	2 min	98
Number of cycles (see Table 15)	20 s	98
	1 min	60
1 cycle	3 min	72
Hold	∞	4

Table 15. Amplification cycles for universal PCR

Primers per pool	Cycle number
6–24	26
25–96	24
97–288	22
289–1056	21
1057–3072	20
3073–5999	18
6000–12,000	17
≥12,001	16

5. After the reaction is complete, place the reactions on ice and proceed to "Procedure: cleanup of universal PCR" (below). Alternatively, the samples can be stored at -20°C in a constant-temperature freezer for up to 3 days.

Procedure: cleanup of universal PCR

6. Add 50 μL Nuclease-free Water to bring each sample to 90 μL .
7. Add 108 μL QIAseq Beads to the completed universal PCR; mix well by vortexing or pipetting up and down at least 15 times with pipetting volume close to 190 μL .
8. Incubate for 5 min at room temperature.
9. Place the tubes or plate on magnetic rack for 5 min (tube) or 10 min (plate) to separate beads from supernatant. Once the solution has cleared, with the tubes or plate still on the magnetic stand, carefully remove and discard the supernatant.
Important: Do not discard the beads as they contain the DNA of interest.
10. With the tubes or plate still on the magnetic stand, add 200 μL 80% ethanol and leave beads in ethanol for 30 seconds. Carefully remove and discard the ethanol.

11. Repeat the ethanol wash.

Important: Completely remove all traces of the ethanol wash after this second wash. Remove the ethanol with a 200 μL pipette first, and then use a 10 μL pipette to remove any residual ethanol.

12. With the tubes or plate still on the magnetic stand, air dry at room temperature for at least 10 min.

Note: Visually inspect that the pellet is completely dry. Over drying the beads will not affect DNA elution.

13. Remove the tubes or plate from the magnetic stand, and elute the DNA from the beads by adding 30 μL Nuclease-free Water. Mix well by pipetting or vortexing.

14. Return the tubes or plate to the magnetic rack until the solution has cleared.

15. Transfer 28 μL supernatant to clean tubes or plate.

16. Proceed to "Recommendations: Library QC and Quantification" (next page). Alternatively, the library can be stored at -30°C to -15°C in a constant-temperature freezer. Amplified libraries are stable for several months at -30°C to -15°C . Once quantification is performed, proceed to "Protocol: Sequencing Setup on Illumina MiSeq, NextSeq 500/550, NextSeq 1000/2000, MiniSeq and NovaSeq" on page 44

Recommendations: Library QC and Quantification

NGS library QC

After the library is constructed and purified, QC can be performed with QIAGEN's QIAxcel Systems or Agilent's Bioanalyzer or TapeStation to check for the correct size distribution of the library fragments and for the absence of primer dimers (approximately <200 bp) and concentration. Majority library fragments prepared for Illumina instruments demonstrate a size distribution between 200 and 400 bp (Figure 6 and Figure 7). Library overamplification is normal (Figure 7B), and this should not affect the sequencing results. Overamplified libraries are usually single-stranded libraries with correct size but appear as "larger fragments" due to secondary structures. Amounts of DNA under the appropriate peaks can be used to quantify the libraries. However, due to the superior sensitivity of qPCR, we recommend quantifying the libraries using the QIAseq Library Quant System, especially when there are overamplified libraries (See "Preferred library quantification method").

Recommended setting for checking QIAseq cfDNA Ultra library on QIAxcel Advanced and Connect systems

QIAxcel Connect

- QIAxcel DNA High Sensitivity Kit (1200) (cat. no. 929012)
- Use method Default High Sensitivity with 1–10 µL library

QIAxcel Advanced or QIAxcel Connect

- QIAxcel DNA High Resolution Kit (cat. no. 929002)
- QX Alignment Marker: 15 bp/5 kb (cat. no. 929524)

- QX DNA Size Marker: 100 bp – 2.5 kb (cat. no. 929559)
- Use the *Application Guide for Low-Concentration Libraries*. To access guides and system files for Library QC, contact QIAGEN Technical Services.





Figure 6. Sample QIAXcel image of QIAseq Targeted cfDNA Ultra libraries for Illumina instruments. The library assessed using QIAXcel illustrates the size of the majority of the library fragments are between 200 and 400 bp.

A: Library (without overamplification) prepared for Illumina instruments



B: Library (with overamplification) prepared for Illumina instruments




Figure 7. Sample Bioanalyzer images of QIAseq Targeted cfDNA Ultra libraries for Illumina instruments. The size of the majority of the library fragments are between 200 and 400 bp. **A:** Library without overamplification. **B:** Library with overamplification as indicated by the “larger fragment” peak.

Preferred library quantification method

The library yield measurements from the Bioanalyzer or TapeStation rely on fluorescence dyes that intercalate into DNA. These dyes cannot discriminate between molecules with or without adapter sequences, yet only complete QIAseq Targeted cfDNA Ultra libraries with full adapter sequences will be sequenced. Due to the superior sensitivity of qPCR, we recommend quantifying the libraries using, QIAGEN’s QIAseq Library Quant Assay Kit (cat. no. 333314), which contains laboratory-verified forward and reverse primers, together with a DNA standard. With this system, the correct dilution of the library can be determined for sequencing. Please refer to the relevant handbook (available at www.qiagen.com) for library quantification.

Protocol: Sequencing Setup on Illumina MiSeq, NextSeq 500/550, NextSeq 1000/2000, MiniSeq and NovaSeq

Important points before starting

- Recommendations for library dilution concentrations and library loading concentrations are based on QIAseq Library Quant System (see "Preferred library quantification method" on the previous page). If using Qubit for library quantification, load 50% less as starting point to avoid over clustering and adjust accordingly after the first run.
- Paired-end sequencing should be used for the QIAseq Targeted cfDNA Ultra on Illumina platform.
- To make sequencing preparation more convenient, download Illumina-compatible sample sheets for different sequencing instruments on www.qiagen.com, from the Resources tab of the QIAseq Targeted cfDNA Ultra page.
- Paired-end sequencing of 149 bp should be used for QIAseq Targeted cfDNA Ultra UDI libraries and dual 10 bp indices on Illumina platforms.
- For two-channel sequencing chemistry platforms such as MiniSeq, NextSeq, and NovaSeq, 10% PhiX can be included in the run to improve sequencing quality. For complete instructions on how to denature sequencing libraries, and set up a sequencing run, please refer to the system-specific Illumina documents.

Sequencing preparations for MiSeq with QIAseq Targeted cfDNA Ultra UDI Sets

1. When working with the QIAseq Targeted cfDNA Ultra UDI Sets, use Local Run Manager (LRM) v2 or later on the instrument to upload a sample sheet (see the Resources tab of the QIAseq Targeted cfDNA Ultra and download the appropriate template) and proceed with

sequencing: Read 1 is 149 bp, Read 2 is 149 bp, and each Index Read is 10 bp.

2. **Sample dilution and pooling:** Dilute libraries to 2 or 4 nM for MiSeq. Then, combine libraries with different sample indexes in equimolar amounts if similar sequencing depth is needed for each library.

Note: Recommendations for library dilution concentrations are based on QIAseq Library Quant System.

Note: If combining libraries with the same number of primers, pool equal volumes of the individual libraries at 4 nM together. If the libraries have different primer numbers, then combine the libraries at volume ratios according to their number of primers. For example, Library A has 5000 primers at 4 nM and Library B has 600 primers at 4 nM; combining 50 μ L of Library A with 6 μ L of Library B will result in similar coverage depth for both Libraries A and B in the same sequencing run.

3. **Library preparation and loading:** Prepare and load the library onto a MiSeq according to the MiSeq System Denature and Dilute Libraries Guide. The final library concentration is 10–12 pM on the MiSeq.

Note: Recommendations for library loading concentrations are based on QIAseq Library Quant System.

4. Upon completion of the sequencing run, proceed to "Appendix B: Data Analysis Using QIAGEN's CLC Genomic Workbench" (page "Appendix B: Data Analysis Using QIAGEN's CLC Genomic Workbench" on page 51).

Sequencing preparations for MiniSeq, NextSeq 500/550 and NextSeq 1000/2000 with QIAseq Targeted cfDNA Ultra UDI Sets

1. When working with the QIAseq Targeted cfDNA Ultra UDI Sets, use LRM v2 or later on the instrument to upload a sample sheet (see the Resources tab of the QIAseq Targeted cfDNA Ultra page and download the appropriate template) and proceed with

sequencing: Read 1 is 149 bp, Read 2 is 149 bp and each Index Read is 10 bp.

2. **Sample dilution and pooling:** Dilute libraries to 1 nM for MiniSeq, 0.5, 1, 2 or 4 nM for NextSeq 500/550 and 2 nM for NextSeq 1000/2000. Then, combine libraries with different sample indexes in equimolar amounts if similar sequencing depth is needed for each library.

Note: Recommendations for library dilution concentrations are based on QIAseq Library Quant System.

Note: If combining libraries with the same number of primers, pool equal volumes of the individual libraries at 4 nM together. If the libraries have different primer numbers, then combine the libraries at volume ratios according to their number of primers. For example, Library A has 5000 primers at 4 nM and Library B has 600 primers at 4 nM; combining 50 μ L Library A with 6 μ L Library B will result in similar coverage depth for both Libraries A and B in the same sequencing run.

3. **Library preparation and loading:** Prepare and load the library onto a MiniSeq, NextSeq 500/550 or NextSeq 1000/2000 according to the MiniSeq or NextSeq System Denature and Dilute Libraries Guide. The final library concentration is 1.0–1.2 pM on the MiniSeq or NextSeq 500/550 and 650 pM for NextSeq 1000/2000 onboard denature and dilute.

Note: Recommendations for library loading concentrations are based on QIAseq Library Quant System.

Note: For P3/P4 flow cell with XLEAP-SBS chemistry kit of NextSeq2000, load 488 pM that's about 25% less compared to standard SBS according to Illumina recommendation.

4. Upon completion of the sequencing run, proceed to "Appendix B: Data Analysis Using QIAGEN's CLC Genomic Workbench" on page 51".

Sequencing preparations for NovaSeq with QIAseq Targeted cfDNA Ultra UDI Index Sets

1. When working with the QIAseq Targeted cfDNA Ultra UDI Sets, upload a sample sheet (see the Resources tab of the QIAseq Targeted DNA Panel page and download the appropriate template) and proceed with sequencing: Read 1 is 149 bp, Read 2 is 149 bp and each Index Read is 10 bp.
2. Sample dilution and pooling: Dilute libraries to 4 nM for NovaSeq. Then, combine libraries with different sample indexes in equimolar amounts if similar sequencing depth is needed for each library.

Note: Recommendations for library dilution concentrations are based on QIAseq Library Quant System.

Note: If combining libraries with the same number of primers, pool equal volumes of the individual libraries at 4 nM together. If the libraries have different primer numbers, then combine the libraries at volume ratios according to their number of primers. For example, Library A has 5000 primers at 4 nM and Library B has 600 primers at 4 nM; combining 50 μ L Library A with 6 μ L Library B will result in similar coverage depth for both Libraries A and B in the same sequencing run.

3. Library preparation and loading: Prepare and load the library onto a NovaSeq according to the NovaSeq 6000 or NovaSeq x Sequencing System Guide. The final pooled library concentration recommendation is between 2–4 nM yielding a final loading concentration of between 140–200 pM on the NovaSeq.

Note: Recommendations for library loading concentrations are based on QIAseq Library Quant System.

4. Upon completion of the sequencing run, proceed to "Appendix B: Data Analysis Using QIAGEN's CLC Genomic Workbench" on page 51.

Troubleshooting Guide

This troubleshooting guide may be helpful in solving any problems that may arise. For more information, see also the Frequently Asked Questions page at our Technical Support Center: www.qiagen.com/FAQ/FAQList.aspx

The scientists in QIAGEN Technical Services are always happy to answer any questions you may have about either the information and/or protocols in this handbook or sample and assay technologies (for contact information, visit www.qiagen.com).

Comments and suggestions

Low library yield

- | | |
|---|---|
| a) Suboptimal reaction conditions due to low cDNA input | Make sure cDNA is free of cellular genomic DNA. If there is cellular genomic DNA, make sure input is based on cDNA peak, not all the DNA in the sample. |
| b) Inefficient targeted enrichment or universal PCR | Check to see if correct thermocycling condition is used during target enrichment or universal PCR. Ensure no ethanol or bead carryover to the PCR. |

Unexpected signal peaks

- | | |
|---|---|
| a) Short peaks less than 200 bp | These are primer–dimers from targeted enrichment PCR or universal PCR. The presence of primer dimers indicates either not enough DNA input or inefficient PCR. |
| b) Larger library fragments after universal PCR | After the universal PCR, library fragments are larger than the intended peak and can be a PCR artifact due to overamplification of the DNA library. Overamplification of the library will not affect the QIAseq Targeted cDNA Ultra sequencing performance. Decreasing the number of universal PCR cycle numbers can reduce over-amplification. |

Sequencing issues

- | | |
|--|---|
| a) Too low or too high cluster density | Accurate library quantification is the key for optimal cluster density on any sequencing instrument. PCR-based quantification method is recommended. Other methods may lead to the incorrect quantification of the library especially when there is over amplification. |
| b) Very low clusters passing filter | Make sure that the library is accurately quantified and that the correct amount is loaded onto the sequencing instrument. |

Comments and suggestions

Variant detection issues

Known variants not detected

Variant detection sensitivity is directly related to the input DNA and read depth. Check Table 2 to see if the required input DNA and read depth are met for the specific variant detection application.

Appendix A: Combining an Existing Panel with a Booster Panel

If additional primers need to be added into an existing panel, a Booster Panel with up to 100 primers can be ordered. To combine the existing panel with a Booster Panel, follow the volume ratio indicated in Table 16.

Table 16. Combining an existing panel (at 50 μ L) with a booster panel

No. of primers in existing panel	Volume of existing panel to combine (μ L)	Volume of booster panel to combine (μ L)
1–2000	50	5
2001–4000	50	3.75
4001–12,000	50	2.5
12,001–20,000	50	1.25

Appendix B: Data Analysis Using QIAGEN's CLC Genomic Workbench

After sequencing, data from the QIAseq Targeted cfDNA Ultra can be analyzed using the Biomedical Genomics Analysis plugin for CLC Genomics Workbench, which allows you to optimize analysis parameters to your specific panel. The parameters can then be locked for routine use. Contact your account manager for further details.

Ordering Information

Product	Contents	Cat. no.
QIAseq Targeted cfDNA Ultra (12)	All reagents (except indexes) for targeted DNA sequencing; fixed panel for 12 samples; less than 200 genes	334051
QIAseq Targeted cfDNA Ultra (96)	All reagents (except indexes) for targeted DNA sequencing; fixed panel for 96 samples; less than 200 genes	334055
QIAseq Targeted cfDNA Ultra Custom (96)	All reagents (except indexes) for targeted DNA sequencing; custom panel for 96 samples	Inquire
QIAseq Targeted cfDNA Ultra Booster (96)	Pool of primers used in combination with either cataloged or custom panels	334085
QIAseq cfDNA Ultra Reagent Kit (96)	Kit containing library preparation and target enrichment reagents (except indexes and primers) for targeted DNA sequencing; fixed reagents for 96 samples	334455
QIAseq Targeted cfDNA Ultra Unique Dual Indices		
QIAseq Targeted cfDNA Ultra UDI Set A (96)	Box containing unique molecularly indexed adapters and primers, enough for a total of 96 samples, for indexing up to 96 samples for targeted panel sequencing on Illumina platforms; Set A (of A, B, C, and D) required for multiplexing 384 samples in one run	334165
QIAseq Targeted cfDNA Ultra UDI Set B (96)	Box containing unique molecularly indexed adapters and primers, enough for a total of 96 samples, for indexing up to 96 samples for targeted panel sequencing on Illumina platforms; Set B (of A, B, C, and D) required for multiplexing 384 samples in one run	334175
QIAseq Targeted cfDNA Ultra UDI Set C (96)	Box containing unique molecularly indexed adapters and primers, enough for a total of 96 samples, for indexing up to 96 samples for targeted panel sequencing on Illumina platforms; Set C (of A, B, C, and D) required for multiplexing 384 samples in one run	334185

Product	Contents	Cat. no.
QIAseq Targeted cfDNA Ultra UDI Set D (96)	Box containing unique molecularly indexed adapters and primers, enough for a total of 96 samples, for indexing up to 96 samples for targeted panel sequencing on Illumina platforms; Set D (of A, B, C, and D) required for multiplexing 384 samples in one run	334195
QIAseq Targeted cfDNA Ultra UDI (12)	Box containing unique molecularly-indexed adapters and primers, enough for a total of 12 samples, for indexing up to 12 samples for targeted panel sequencing on Illumina platforms	334151
Related products		
QIAseq Library Quant Assay Kit	Reagents for quantification of libraries prepared for Illumina or Ion Torrent platforms; assay format	333314
QIAamp Circulating Nucleic Acid Kit (50)	For 50 preps: QIAamp Mini Columns, Tube Extenders (20 mL), QIAGEN Proteinase K, Carrier RNA, Buffers, VacConnectors, and Collection Tubes	55114
QIAamp MinElute ccfDNA Mini Kit (50)	For 50 DNA preps: QIAamp UCP MinElute Columns, QIAGEN Proteinase K, Magnetic Bead Suspension, Buffers, Bead Elution Tubes, Collection Tubes	55204
QIAamp MinElute ccfDNA Midi Kit (50)	For 50 preps: QIAamp UCP MinElute Columns, QIAGEN Proteinase K, Magnetic Bead Suspension, Buffers, Bead Elution Tubes, and Collection Tubes	55284
EZ1&2 ccfDNA Kit (48)	For 48 preps: 48 reagent cartridges (EZ1&2 ccfDNA), Magnetic Bead Suspension, Elution Buffer, Large-Volume Tubes, Disposable Tip Holders, Disposable Filter-Tips, Elution Tubes	984854

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Document Revision History

Date	Changes
12/2022	Initial revision
04/2023	With reduced HiFi polymerase during TEPCR for panel with 1100 primers or less; Added two panels, Actionable Solid Tumor and Myeloid Neoplasms; Included NextSeq 1000/2000 in sequencing set up; Corrected UDI index plate name in index kit component.
07/2025	Change in list of Illumina NGS systems: removed HiSeq with NovaSeq; Update in Important Notes: Addition of column "10,000 primers" for number of multiplexed samples; Update in Protocols; addition of cycling conditions in Protocol: Target Enrichment.

Limited License Agreement for QIAseq® Targeted cfDNA Ultra Kit

Use of this product signifies the agreement of any purchaser or user of the product to the following terms:

1. The product may be used solely in accordance with the protocols provided with the product and this Instructions for Use and for use with components contained in the panel only. QIAGEN grants no license under any of its intellectual property to use or incorporate the enclosed components of this panel with any components not included within this panel except as described in the protocols provided with the product, this Instructions for Use, and additional protocols available at www.qiagen.com. Some of these additional protocols have been provided by QIAGEN users for QIAGEN users. These protocols have not been thoroughly tested or optimized by QIAGEN. QIAGEN neither guarantees them nor warrants that they do not infringe the rights of third-parties.
2. Other than expressly stated licenses, QIAGEN makes no warranty that this panel and/or its use(s) do not infringe the rights of third-parties.
3. This panel and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. QIAGEN specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the panel agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. QIAGEN may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the panel and/or its components.

For updated license terms, see www.qiagen.com.

Trademarks: QIAGEN®, Sample to Insight®, QIAamp®, QIAseq®, QIAxcel®, QIAxpert®, QCI® (QIAGEN Group); Agilent®, Bioanalyzer®, TapeStation® (Agilent Technologies); DynaMag™, Qubit™ (Thermo Fisher Scientific Inc.); LoBind® (Eppendorf AG); Illumina®, HiSeq®, MiniSeq®, MiSeq®, NextSeq®, NovaSeq™ (Illumina, Inc.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

07/2025 HB-3094-003 © 2025 QIAGEN, all rights reserved.

