Januar 2025.

PartoSure[®] test

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Verzija 1

Za korišćenje u in vitro dijagnostici

0197

TTDT-1-20-IVDR

QIAGEN Sciences LLC,
19300 Germantown Road,
Germantown,
MD 20874, SAD

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3. P-0250 Poster - Mechanisms for preterm labor and fetal injury - Evaluation of PAMG -1 for the Prediction of Preterm Birth in Patients Symptomatic of Preterm Labor Lotfi, G.; Faraz, S.; Al Swalhee, N.; Nasir, R.; Somini, S.; Abdeldayem, R.; Koratkar, R.; Ammar, A. Dubai Health Authority, Duba Latfa Hospital, United Arab Emirates.
4. Dawes LK, Prentice LR, Huang Y, Groom KM. The Biomarkers for Preterm Birth Study-A prospective observational study comparing the impact of vaginal biomarkers on clinical practice when used in women with symptoms of preterm labor. *Acta Obstet Gynecol Scand.* 2020 Feb;99(2):249-258
5. E1294 - Placental Alpha Microglobulin 1 To Predict Spontaneous Preterm Birth In Symptomatic Woman Morales F. 1, Galarce V., Guerra M. (2019). 3 – Poster Presentations, *Journal of Perinatal Medicine*, 47(s1), eA327-eA550.
6. Ali Gokce, Erkan Kalafat, Yavuz Emre Sukur, Orhan Altinboga & Feride Soylemez (2020): Role of cervical length and placental alpha microglobulin-1 to predict preterm birth, *The Journal of Maternal-Fetal & Neonatal Medicine*.
7. Maryam Kashanian, Nooshin Eshraghi, Maryam Rahimi & Narges Sheikhsari (2020): Evaluation of placental alpha microglobulin-1(PAMG1) accuracy for prediction of preterm delivery in women with the symptoms of spontaneous preterm labor; a comparison with cervical length and number of contractions, *The Journal of Maternal-Fetal & Neonatal Medicine*.
8. Melchor JC, Khalil A, Wing D, Schleussner E, Surbek D. Prediction of preterm delivery in symptomatic women using PAMG-1, fetal fibronectin and pHGFBP-1 tests: systematic review and meta-analysis. *Ultrasound Obstet Gynecol.* 2018 Oct;52(4):442-451.
9. Pirjani R, Moini A, Almasi-Hashiani A, Farid Mojtahedi M, Vesali S, Hosseini L, Sepidarkish M. Placental alpha microglobulin-1 (PartoSure) test for the prediction of preterm birth: a systematic review and meta-analysis. *J Matern Fetal Neonatal Med.* 2019 Nov 17:1-13.
10. Nikolova T, Uotila J, Nikolova N, Bolotskikh VM, Borisova VY, Di Renzo GC. Prediction of spontaneous preterm delivery in women presenting with premature labor: a comparison of placenta alpha microglobulin-1, phosphorylated insulin-like growth factor binding protein-1, and cervical length. *Am J Obstet Gynecol.* 2018 Dec;219(6):610.e1-610.e9.

	11. Sergey V. Barinov, Gian Carlo Di Renzo, Antonina A. Belinina, Olga V. Koliado & Olga V. Remneva (2021): Clinical and biochemical markers of spontaneous preterm birth in singleton and multiple pregnancies, The Journal of Maternal-Fetal & Neonatal Medicine. 12. Chawanpaiboon, Saifon, Titapant, Vitaya and Pooliam, Julaporn. "Placental α-microglobulin-1 in cervicovaginal fluid and cervical length to predict preterm birth by Thai women with symptoms of labor" Asian Biomedicine, vol.15, no.3, 2021, pp.119-127. 13. Cnota W, Jagielska A, Janowska E, Banas E, Kierach R, Nycz-Reska M, Czuba B. Prediction of preterm birth using PAMG-1 test: a single centre experience - preliminary report. Ginekol Pol. 2022 Jan 24. 14. Pirjani R, Moini A, Almasi-Hashiani A, Farid Mojtahedi M, Vesali S, Hosseini L, Sepidarkish M. Placental alpha microglobulin-1 (PartoSure) test for the prediction of preterm birth: a systematic review and meta-analysis. J Matern Fetal Neonatal Med. 2021 Oct;34(20):3445-3457. 15. Konoplyannikov AG, Dikke GB, Karaganova EY. Combination of the placental alpha-1 microglobulin test and ultrasonic cervical length measurement to predict the time of preterm birth. J Matern Fetal Neonatal Med. 2022 Feb; 35(3):541-545. 16. Barinov SV, Di Renzo GC, Belinina AA, Koliado OV, Remneva OV. Clinical and biochemical markers of spontaneous preterm birth in singleton and multiple pregnancies. J Matern Fetal Neonatal Med. 2021 Feb 24:1-6. 17. Dochez V, Ducarme G, Gueudry P, Joueidi Y, Boivin M, Boussamet L, Pelerin H, Le Thuaut A, Lamoureux Z, Riche VP, Winer N, Thubert T, Marie E. Methods of detection and prevention of preterm labour and the PAMG-1 detection test: a review. J Perinat Med. 2020 Oct 2;49(2):119-126. 18. Kehl S, Weiss C, Pretscher J, Baier F, Faschingbauer F, Beckmann MW, Stumpfe FM. The use of PAMG-1 testing in patients with preterm labor, intact membranes and a short sonographic cervix reduces the rate of unnecessary antenatal glucocorticoid administration. J Perinat Med. 2021 Jul 19;49(9):1135-1140.
6.2 Sažetak podataka o učinku sa ekvivalentnog medicinskog sredstva, ako je primenljivo	Nije primenljivo.

6.3 Sažetak podataka o učinku iz sprovedenih studija uređaja pre dobijanja CE oznake	Sažetak kliničkih i analitičkih studija performansi je dat u nastavku: Klinički učinak medicinskog sredstva u evaluaciji ● Očekivane vrednosti PartoSure test je imunohromatografsko ispitivanje s lateralnim protokom namenjeno otkrivanju prisustva humanog alfa-1 mikroglobulina posteljice (PAMG-1) u amnionskoj tečnosti. PAMG-1 je izabran kao marker za preciznu procenu rizika od prevremenog porođaja zbog svojih jedinstvenih karakteristika, tj. visokog nivoa u amnionskoj tečnosti, niskog nivoa u krvi i izuzetno niskog pozadinskog nivoa (50-220 pikograma/ml) u cervikovaginalnom iscetku. Test koristi monoklonska antitela koja su u dovoljnoj meri osetljiva da otkriju 1 ng/ml PAMG-1 kada se eluira u bočicu rastvarača. Sa razblaživanjem od 3 do 4 puta jednom eluiran u bočicu rastvarača, granična koncentracija je 13–18 puta iznad maksimalne pozadinske koncentracije PAMG-1. Uzorak cervikovaginalnog iscetka za analizu se uzima pomoću štapića za uzorkovanje iz vagine i izdvaja se u rastvaraču. Prisustvo antigena PAMG-1 se zatim otkriva uranjanjem test trake sa lateralnim protokom u bočicu. Uzorak teče od upijajućeg jastučića na nitroceluloznu membranu, prolazi kroz reaktivnu zonu u kojoj se nalaze monoklonska anti-PAMG-1 antitela konjugovana sa zlatnim česticama. Kompleks antigena i antitela teče prema zoni za testiranje u kojoj se imobilise drugim anti-PAMG-1 antitelom. Taj događaj uzrokuje pojavu linije testa na traci. Nevezani kompleksi antigena i antitela nastavljaju da se slivaju niz test traku i imobilise ih drugo antitelo. To uzrokuje pojavu linije interne kontrole. PartoSure otkriva humani PAMG-1 u tragovima. Koncentracija PAMG-1 veća od 1 ng/ml ukazuje na povećani rizik od porođaja u narednih 7 dana. Karakteristike učinka: Preciznost i mogućnost ponavljanja su utvrđeni pomoću tri serije PartoSure testa na tri lokacije od strane tri različita korisnika na sve tri lokacije. Korišćeno je pet replikata sedam različitih koncentracija PAMG-1 iznad i ispod ograničenja detekcije PartoSure testa, uključujući apsolutnu nulu (0,0 ng/ml), nisku negativnu vrednost (0,2 ng/ml), negativnu vrednost različitu od nule (0,5 ng/ml), ograničenje detekcije (1,0 ng/ml), nisku pozitivnu vrednost (2,0 ng/ml).
--	--

Klinički učinak testa PartoSure je procenjen u više recenziranih publikacija koje opisuju upotrebu testa PartoSure u kliničkim studijama. Da bi bili uključeni u pregled, radovi su morali da navedu metode koje su uključivale izvođenje testa PartoSure u skladu sa svojim uputstvima za upotrebu (IFU) dostupnim u to vreme. Ukupno je u konačnu analizu uključeno 17 publikacija. Uz odgovarajući interval pouzdanosti od 95% izračunat je sledeći učinak:

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