

June 2015

GeneRead™ DNaseq Targeted Panels V2 Handbook

For targeted enrichment prior to next-
generation sequencing

(96-well plate compatible, all-bead purification)



Sample & Assay Technologies

QIAGEN Sample and Assay Technologies

QIAGEN is the leading provider of innovative sample and assay technologies, enabling the isolation and detection of contents of any biological sample. Our advanced, high-quality products and services ensure success from sample to result.

QIAGEN sets standards in:

- Purification of DNA, RNA and proteins
- Nucleic acid and protein assays
- microRNA research and RNAi
- Automation of sample and assay technologies

Our mission is to enable you to achieve outstanding success and breakthroughs. For more information, visit www.qiagen.com.

Contents

Kit Contents	5
Storage	7
Intended Use	8
Safety Information	8
Quality Control	8
Introduction	9
Principle and procedure	9
Equipment and Reagents to be Supplied by User	12
Important Notes	14
DNA preparation and quality control	14
DNA quantification and quality control	14
Protocol: PCR Setup	16
Protocol: Sample Pooling and Purification	22
Troubleshooting Guide	24
References	24
Appendix A: Library Construction Using GeneRead Library Prep Kits for Illumina	25
End repair of DNA	25
A-addition	26
Adapter ligation	26
Cleanup of adapter-ligated DNA with AMPure XP beads	27
PCR amplification of purified library	28
Cleanup of amplified library with AMPure XP beads	29
Appendix B: Library Construction Using GeneRead Library Prep Kits for Ion PGM Sequencer/Proton	31
End repair of DNA	31
Adapter ligation	32
Cleanup of adapter-ligated DNA with AMPure XP beads	33
PCR amplification of the purified library	34
Cleanup of amplified library with AMPure XP beads	35

Appendix C: FFPE DNA Quality and Quantity	37
Appendix D: Library Quantification and Quality Control	38
Appendix E: Data Analysis using QIAGEN's GeneRead DNaseq Sequence Variant Analysis Software	39
Appendix F: Combine Libraries for Multiplex Sequencing	40
Appendix G: Analyze the PCR amplicons and Library Using the QIAxcel Advanced with a QIAxcel DNA High Resolution Kit or Agilent 2100 Bioanalyzer	41
Appendix H: Downloading Individual Unaligned .BAM File with a Multiplex Sample on Ion PGM Sequencer	43
Ordering Information	45

Kit Contents

GeneRead DNaseq Targeted Panel V2	181900*
Catalog no.	181900*
Pools with enough primers for 12 or 96 samples, depending on pack size	1 [†] /4
Handbook	1

* Gene panel pools are labeled A1, A2, A3 and A4.

† E.g., the Tumor Actionable Mutations GeneRead DNaseq Targeted Panel V2.

GeneRead DNaseq Targeted HC Panel V2	181901[‡]
Catalog no.	181901[‡]
Pools with enough primers for 12 or 96 samples, depending on pack size	4
Handbook	1

‡ Gene panel pools are labeled A1, A2, A3 and A4.

§ E.g., the Human Comprehensive Cancer, Carrier Testing and Cancer Predisposition GeneRead DNaseq Targeted HC Panel V2.

GeneRead DNaseq Custom Panel V2	181902
Catalog no.	181902
Pools with primers for any gene or region in the human genome, up to 9600 primer pairs, for 480 samples	1/2/3/4 [¶]
Handbook	1

¶ Number of pools determined by covered region.

GeneRead DNaseq Mix-n-Match Panel V2	181905
Catalog no.	181905
Pools with bench-tested primers for a mix of any gene in catalog panels, up to 9600 primer pairs, for 96 samples	4
Handbook	1

Catalog no.	Product name	No. of total primer pairs	No. of pools
NGHS-001X	Human Breast Cancer Panel	2915	4
NGHS-002X	Human Colorectal Cancer Panel	1954	4
NGHS-003X	Human Myeloid Neoplasms Panel	2536	4
NGHS-004X	Human Liver Cancer Panel	2052	4
NGHS-005X	Human Lung Cancer Panel	3586	4
NGHS-006X	Human Ovarian Cancer Panel	2021	4
NGHS-007X	Human Prostate Cancer Panel	1837	4
NGHS-008X	Human Gastric Cancer Panel	2377	4
NGHS-009X	Human Cardiomyopathy Panel	2657	4
NGHS-011X	Human Carrier Testing Panel	6943	4
NGHS-013X	Human Cancer Predisposition Panel	6582	4
NGHS-501X	Human Comprehensive Cancer Panel	7951	4
NGHS-101X	Clinically Relevant Tumor Panel	602	4
NGHS-102X	Human BRCA1 and BRCA2 Panel	250	4
NGHS-201X	Tumor Actionable Mutations Panel	118	1

GeneRead DNaseq Panel PCR Kit V2	181940 (12)	181942 (96)
1	(24/24/16/12* samples)	(192/192/128/96* samples)
GeneRead DNaseq Panel 5x PCR Buffer	230 µl	1800 µl
GeneRead HotStarTaq® DNA Polymerase 6 U/µl	80 µl	600 µl
DNase-free water	1000 µl	1000 µl

* Number of samples that can be processed for 1, 2, 3 or 4 pool panels, respectively.

Required mastermix quantities depend on the number of PCR pools and number of samples.

# PCR pools	# samples per panel	Panel type	PCR Kit required	Quantity
1	12	Cataloged	181940	1
	96	Cataloged	181942	1
	480	Custom	181942	3
2	480	Custom	181942	3
3	480	Custom	181942	4
4	12	Cataloged	181940	1
	96	Cataloged and Mix-n-Match	181942	1
	480	Custom	181942	5

Storage

GeneRead DNaseq Panel Kits are shipped on dry ice and should be stored at -30°C to -15°C upon arrival. When stored properly at -30°C to -15°C , all reagents are stable for up to 6 months after delivery. GeneRead DNaseq Panel PCR Kits are shipped on cold packs. For long-term storage, keep tubes at -30°C to -15°C . If the entire volume will not be used at once, we recommend dividing into aliquots and storing at -30°C to -15°C . Avoid repeated freezing and thawing. If stored under these conditions, GeneRead DNaseq Panel PCR Kits are stable for 6 months after receipt.

Intended Use

GeneRead DNaseq Targeted Panels and GeneRead DNaseq Panel PCR Kits are intended for molecular biology applications. These products are not intended for the diagnosis, prevention or treatment of a disease.

All due care and attention should be exercised in the handling of the products. We recommend all users of QIAGEN products to adhere to the NIH guidelines that have been developed for recombinant DNA experiments, or to other applicable guidelines.

Safety Information

When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at www.qiagen.com/safety where you can find, view and print the SDS for each QIAGEN kit and kit component.

Quality Control

GeneRead DNaseq Panel Kits are tested, and each assay in the GeneRead DNaseq Targeted Panels is tested when designed, against predetermined specifications to ensure consistent product quality.

Introduction

DNA sequencing is a useful tool to detect genetic variations, including somatic mutations, SNPs and small insertions and deletions. Targeted enrichment technology enables next-generation sequencing (NGS) platform users to sequence specific regions of interest instead of the entire genome, effectively increasing sequencing depth and throughput with lower cost. GeneRead DNaseq Targeted Panels V2 use multiplex PCR-based targeted enrichment technology, in combination with a sophisticated primer design algorithm, to enable amplification and enrichment of any gene or targeted region in the human genome in order to detect genetic variation using NGS (Figure 1). Adjacent and potentially interacting primer pairs are separated into different pools for optimal performance. GeneRead DNaseq Targeted Panels V2 are designed to analyze a panel of genes and/or regions related to a disease state and can be used with any major NGS platform. The targeted enrichment process is essential for the efficient utilization of medium- and high-throughput sequencers such as Life Technologies[®], Ion Personal Genome Machine[®] (PGM[™]) Sequencer and Ion Proton[™], as well as Illumina[®]'s MiSeq[®] Personal Sequencer, HiSeq[®] 1000, HiSeq 1500, HiSeq 2000, HiSeq 2500 and GAllx.

GeneRead DNaseq Targeted Panels V2 have been optimized in combination with GeneRead DNaseq Panel PCR Kits V2 to provide superior sensitivity and linear multiplex amplification. The simplicity of the PCR method makes these panels accessible for routine use in every research laboratory.

Principle and procedure

GeneRead DNaseq Targeted Panels V2 are provided as sets of 1, 2, 3 or 4 pools, each containing primer mixes, with up to 9600 primer pairs per 4-pool set. The number of pools included is determined by the region covered by amplicons. Most panels cover full coding regions of genes, so one 4-pool set is provided. For special panels, like the Tumor Actionable Mutations Panel covering tumor mutation hotspots, one pool is provided. GeneRead DNaseq Targeted Panels V2 can enrich selected genes and/or regions using as little as 40 or 20 ng genomic DNA in 3 hours for a 4- or 1-pool panel, respectively (Figure 2). Briefly, genomic DNA samples are combined with primer mix and PCR reagent and PCR is performed in a standard thermocycler. After the reaction is complete, the reactions for each sample are pooled and the enriched DNA is purified. The purified DNA then is ready for NGS library construction and sequencing using the NGS platform of your choice. Since amplicons from each sample are pooled before library construction, each sample results in one

library only. The sequencing results can be analyzed using the GeneRead DNaseq Analysis Software at <http://ngsdataanalysis.sabiosciences.com>, which will automatically perform all steps necessary to generate a DNA sequence variant report from your NGS data (Figure 3). Alternatively, sequencing results can be analysed by the Biomedical Genomics Workbench platform. All detected variants can be interpreted by the Ingenuity Variant Analysis tool (<https://www.qiagen.com/us/products/catalog/sample-technologies/dna-sample-technologies/genomic-dna/generad-dnaseq-gene-panels-v2/#resources>).




Figure 1. Multiplex PCR-based targeted enrichment scheme. GeneRead DNaseq Targeted Panels V2 use multiplex PCR-based targeted enrichment technology in combination with a sophisticated primer design algorithm to maximize design coverage and minimize nonspecific amplification. The adjacent primer sets are distributed across an appropriate number of pools to minimize nonspecific amplification products.




Figure 2. The GeneRead DNaseq Targeted Panel V2 procedure.




Figure 3. Overview of the sample-to-insight NGS workflow with GeneRead DNaseq Targeted Panels V2. The complete sample-to-insight procedure begins with DNA extraction, followed by targeted enrichment with GeneRead DNaseq Targeted Panels V2, NGS library construction, sequencing and data analysis using either the QIAGEN NGS Data Analysis Web Portal or Biomedical Genomics Workbench. Follow-up experiments or sample verification against specific targets can be performed with qBiomarker Somatic Mutation PCR Assays. Detected variants can be interpreted with the Ingenuity® Variant Analysis tool.

Equipment and Reagents to be Supplied by User

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, consult the appropriate safety data sheets (SDSs), available from the product supplier.

In addition to the GeneRead DNaseq Targeted Panel V2 and GeneRead DNaseq Panel PCR Kit V2, the following supplies are required:

For genomic DNA isolation:

- See page 14 for specific recommendations.

For targeted enrichment:

- High-quality, nuclease-free water. **Do not use DEPC-treated water.**
- GeneRead DNA QuantiMIZE Array or Assay Kit if using FFPE samples (QIAGEN cat. nos. 180642/180654)
- Agencourt® AMPure® XP Kit (Beckman Coulter cat. no. A63880)
- Microcentrifuge
- 1.5 ml LoBind tubes
- 0.2 ml PCR tubes, 96-well reaction plates, or PCR strips and caps
- Thermal cycler
- Multichannel pipettor
- Single-channel pipettor
- DNase-free pipet tips and tubes
- 0.2 ml 96-well PCR plate
- QIAxcel Advanced with a QIAxcel DNA High Resolution Kit or Agilent® 2100 Bioanalyzer®
- Agilent High Sensitivity DNA Kit (Agilent cat. no. 5067-4626)
- DynaMag™-96 Side Magnet (Thermo Fisher cat no. 12331D)

For NGS library construction for Ion PGM and Proton Sequencers:

- GeneRead DNA Library L Core Kit (QIAGEN cat. no. 180462)
- GeneRead DNA L Amp Kit (QIAGEN cat. no. 180485)

- GeneRead Adapter L Set 1-plex (QIAGEN cat. no. 180922) or GeneRead Adapter L Set 12-plex (QIAGEN cat. no. 180994)
- Agencourt AMPure XP Kit (Beckman Coulter cat. no. A63880)
- GeneRead DNAseq Library Quant Kit for Ion PGM Sequencer (QIAGEN cat. no. 180601)
- 0.2 ml 96-well PCR plate
- 80% ethanol
- 96-well Thermal cycler
- A real-time thermal cycler compatible with 96-well/100-well/384-well plates
- DynaMag-96 Side Magnet (Thermo Fisher cat no. 12331D)

For NGS library construction for Illumina MiSeq/HiSeq:

- GeneRead DNA Library I Core Kit (12) (QIAGEN cat. no. 180432) or GeneRead DNA Library I Core Kit (48) (QIAGEN cat. no. 180434)
- GeneRead DNA I Amp Kit (QIAGEN cat. no. 180455)
- GeneRead Adapter I Set A 12-plex (QIAGEN cat. no. 180985) or GeneRead Adapter I Set B 12-plex (QIAGEN cat. no. 180986) for up to 12 samples, or both kits for up to 24 samples.
- GeneRead DNAseq Library Quant Kit for Illumina (QIAGEN cat. no. 180601)
- Agencourt AMPure XP Kit (Beckman Coulter cat. no. A63880)
- 0.2 ml 96-well PCR plate
- 80% ethanol
- Thermal cycler
- A real-time PCR machine compatible with 96-well/100-well/384-well plates
- DynaMag-96 Side Magnet (Thermo Fisher cat no. 12331D)

Important Notes

DNA preparation and quality control

High-quality DNA is essential for obtaining good sequencing results

The most important prerequisite for DNA sequence analysis is consistent, high-quality DNA from every experimental sample. Therefore, sample handling and DNA isolation procedures are critical to the success of the experiment. Residual traces of proteins, salts, or other contaminants may either degrade the DNA or decrease the efficiency of, if not block completely, the enzyme activities necessary for optimal targeted genome amplification and real-time PCR performance.

Recommended genomic DNA preparation method

The QIAGEN QIAamp® DNA Mini Kit (cat. no. 51304), QIAamp DNA FFPE Tissue Kit (cat. no. 56404) and GeneRead DNA FFPE Kit (cat. no. 180134) are highly recommended for the preparation of genomic DNA samples from fresh tissues and FFPE tissue samples. Ensure that samples have been treated for the removal of RNA, as RNA contamination will cause inaccuracies in DNA concentration measurements. **Do not** omit the recommended RNase treatment step to remove RNA. If genomic DNA samples need to be harvested from biological samples for which kits are not available, please contact Technical Support representatives for suggestions.

For best results, all DNA samples should be resuspended in DNase-free water or alternatively in DNase-free 10 mM Tris buffer pH 8.0. **Do not use DEPC-treated water.**

DNA quantification and quality control

For best results, all DNA samples should also demonstrate consistent quality according to the following criteria:

Concentration and purity determined by UV spectrophotometry

The concentration and purity of DNA should be determined by measuring the absorbance in a spectrophotometer. Prepare dilutions and measure absorbance

in 10 mM Tris·Cl, * pH 8.0. The spectral properties of nucleic acids are highly dependent on pH.

$A_{260}:A_{280}$ ratio should be greater than 1.8

Concentration determined by A_{260} should be $>2.5 \mu\text{g/ml}$ DNA.

DNA integrity

For best results, the genomic DNA should be greater than 2 kb in length with some fragments greater than 10 kb. This can be checked by running a fraction of each DNA sample on a 1% agarose gel.

FFPE DNA

If FFPE DNA will be used for GeneRead DNaseq Targeted Panels, the QIAGEN GeneRead DNA QuantiMIZE Array or Assay Kit is recommended for determining optimal DNA amount and PCR cycling conditions for each FFPE DNA sample (see Appendix C).

* When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs), available from the product supplier.

Protocol: PCR Setup

Procedure

1. Dilute DNA to 2.5 ng/ μ l with DNase-free water in a LoBind tube. For each sample, 20 ng (8 μ l, 2.5 ng/ μ l) DNA is required for the 1-pool or 2-pool panels, 30 ng (12 μ l, 2.5 ng/ μ l) for the 3-pool panels, or 40 ng (16 μ l, 2.5 ng/ μ l) for the 4-pool panels.

Note: Dilution of FFPE DNA samples should be determined by the QIAGEN GeneRead DNA QuantiMIZE Array or Assay Kit for optimal results.

2. Determine the number of reactions needed. For a 1-pool panel, one 40 μ l reaction for each sample is required. For 2-, 3- or 4-pool panels, 2, 3, or 4 x 20 μ l reactions for each sample are required. Prepare PCR strips or a PCR plate according to the number of reactions. Label with sample names and pool numbers.
3. Aliquot 8 μ l (1-pool panel) or 4 μ l (2-, 3- or 4-pool panels) of each DNA sample into each PCR reaction. Figure 4 shows an example of a 4-pool panel being set up for 3 samples.

Note: If the GeneRead DNA QuantiMIZE Array or Assay Kit is used, the same volume of FFPE DNA will be used after diluting DNA samples according to the GeneRead DNA QuantiMIZE Array or Assay Kit.




Figure 4. An example of how a 4-pool panel can be set up for 3 samples.

4. Prepare the PCR mix on ice according to Table 1 (1-pool panel) or Table 2 (2, 3 or 4-pool panel), and as shown in Figure 5. For each sample, 1, 2, 3 or 4 PCR mixes will be needed. Mix gently by pipetting up and down.

Table 1. Preparation of PCR mix for each primer mix pool (1-pool panel e.g., Tumor Actionable Mutations Panel)

Component	Per 1 sample (µl)	Per <i>n</i> samples (µl)
GeneRead DNaseq Panel PCR Buffer (5x)	8.8	8.8 x <i>n</i>
Primer mix pool (2x)	22	22 x <i>n</i>
GeneRead HotStarTaq DNA Polymerase (6 U/µl)	2.9	2.9 x <i>n</i>
DNase-free water	1.5	1.5 x <i>n</i>
Total volume	35.2	35.2 x <i>n</i>

Table 2. Preparation of PCR mix for each primer mix pool (2, 3 or 4-pool panels)

Component	Per 1 sample (µl)	Per <i>n</i> samples (µl)
GeneRead DNaseq Panel PCR Buffer (5x)	4.4	4.4 x <i>n</i>
Primer mix pool x* (2x)	11	11 x <i>n</i>
GeneRead HotStarTaq DNA Polymerase (6 U/µl)	1.5	1.5 x <i>n</i>
DNase-free water	0.7	0.7 x <i>n</i>
Total volume	17.6	17.6 x <i>n</i>

* The number of primer mix pools is determined by the panel size. Primer pools are named A1 through A4.




Figure 5. Preparation of PCR mastermix for each primer mix pool. For ease of setup, prepare enough PCR mastermix to support all the primer pools and samples (left). Then, prepare pool-specific mastermixes (right, MM A1, MM A2 ...).

5. Aliquot 32 μ l (1-pool panel) or 16 μ l (2, 3 or 4-pool panel) of each PCR mix, and add it to the well with DNA samples accordingly, as shown in Figure 6. Mix gently by pipetting up and down.




Figure 6. Example of a 4-pool panel. DNA samples have already been aliquoted into wells in step 3. **Do not add anymore DNA.**

6. Seal the wells with PCR tube caps. Place strips or plate in thermocycler and set up reaction parameters according to Tables 3, 4A and 4B.

Table 3. PCR program

Cycle	Temperature	Time
1	95°C	15 min
Number of cycles according to tables 4A or 4B	95°C	15 s
	60°C	4/8 min*
1	72°C	10 min
1	4°C	∞

* If number of primer pairs is <1200 in each pool, 4 min; if number of primer pairs is 1201–2500 per pool, 8 min. To determine number of primer pairs per pool, refer to tables 4A or 4B.

Table 4A. PCR cycles (if using a cataloged panel, refer to Table 4B)

Primer pairs per pool	No. of cycles for standard DNA*
1-11	25
12-23	24
24-47	23
48-95	22
96-191	21
192-287	20
288-399	19
400-1200	18
1201-2500	16

* Number of cycles for FFPE samples should be determined by QIAGEN GeneRead DNA QuantiMIZE Array or Assay Kit.

Table 4B. PCR cycles and cycling times for cataloged panels

Cat #	Panel name	# primer pairs	# pools	# primer pairs per pool	# standard PCR cycles needed	Cycling time (min)
NGHS-001X	Breast cancer	2915	4	729	18	4
NGHS-002X	Colorectal cancer	1954	4	489	18	4
NGHS-003X	Myeloid Neoplasms	2536	4	634	18	4
NGHS-004X	Liver cancer	2052	4	513	18	4
NGHS-005X	Lung cancer	3586	4	897	18	4
NGHS-006X	Ovarian cancer	2021	4	506	18	4
NGHS-007X	Prostate cancer	1837	4	460	18	4
NGHS-008X	Gastric cancer	2377	4	595	18	4
NGHS-009X	Cardiomyopathy	2657	4	665	18	4
NGHS-011X	Carrier testing	6943	4	1736	16	8
NGHS-013X	Cancer predisposition	6582	4	1646	16	8
NGHS-501X	Comprehensive cancer	7951	4	1988	16	8
NGHS-101X	Clinically relevant	602	4	151	21	4
NGHS-102X	BRCA1 and BRCA2	250	4	63	22	4
NGHS-201X	Actionable mutations	118	1	118	21	4

- After the reaction is complete, place on ice and proceed with sample pooling and purification using AMPure XP beads.

Note: If the samples are to be stored prior to purification, transfer them to a -20°C freezer. Samples are stable for 3 days.

Protocol: Sample Pooling and Purification

1. For a 2, 3 or 4-pool panel, combine all 2, 3 or 4 reactions from the same sample into one well of a PCR plate/strip. Mix thoroughly. The volume of each sample should be approximately 40 μ l for a 1- or 2-pool panel and 60 or 80 μ l for a 3- or 4-pool panel.
2. Transfer 40 μ l from each sample to a 1.5 ml LoBind tube or 96-well PCR plate for purification.
3. Add 36 μ l (0.9x volume) AMPure XP beads to 40 μ l PCR product. Mix well by pipetting.
4. Incubate for 5 min at room temperature.
5. Place the tube or 96-well PCR plate on magnetic rack to separate beads from supernatant. After the solution is clear (approximately 5 min), carefully transfer 70 μ l of the supernatant to a new tube or wells in 96-well plate without disturbing the beads. Discard the beads, which contain unwanted large DNA fragments.

Note: Do not discard the supernatant.

IMPORTANT: Transferring 70 μ l supernatant will leave behind about 6 μ l supernatant. This is to ensure that no beads are carried over into the supernatant. Any bead carryover will result in a significant amount of larger fragments present in the library, which will affect sequencing specificity.

6. Add 64 μ l (1.6x the original volume of PCR product, which was 40 μ l) AMPure XP beads to the supernatant, mix well by pipetting and incubate for 5 min at room temperature.
7. Place the tube or 96-well PCR plate on a magnetic rack and wait until solution is clear (approximately 5 min). Carefully remove and discard supernatant. Be careful not to disturb the beads, which contain the DNA target.

Note: Do not discard the beads.

8. Add 200 μ l fresh 80% ethanol to the tube or well while it is on the magnetic rack. Rotate the tube or move the plate side-to-side in the two positions of the magnet to wash the beads, then carefully remove and discard the supernatant.
9. Repeat previous step once.
10. Completely remove ethanol and dry beads for 15 min while the tube or plate is on the rack.

11. Elute DNA target beads in 28 μ l nuclease-free water. Mix well by pipetting. Place tube or plate on magnetic rack until solution is clear. Transfer 25 μ l supernatant to a clean LoBind tube or well in the PCR plate.
12. Determine the amount of your sample with a QIAxcel Advanced with a QIAxcel DNA High Resolution Kit or an Agilent 2100 Bioanalyzer, using the High Sensitivity DNA Kit. Normally 15–60 ng of PCR product will be obtained after purification.
13. Proceed to library construction according to the sequencing platform of your choice. Refer to Appendix A for the recommended library construction protocol for sequencing with Illumina MiSeq/HiSeq. Refer to Appendix B for the recommended library construction protocol for sequencing with Ion PGM Sequencer.

IMPORTANT: The protocols described in Appendix A and B are different from the protocols described in the *GeneRead Library Prep (I) Handbook*. For preparing libraries using panels, follow the steps described here only.

Note: If reactions are to be stored prior to library construction, transfer them to a -20°C freezer. Samples are stable for 3 days.

Troubleshooting Guide

For technical support, please call us at 1-888-503-3187 or 1-301-682-9200. For more information, see also the Frequently Asked Questions page at our Technical Support Center: www.qiagen.com/FAQ/FAQList.aspx. The scientists in QIAGEN Technical Services are always happy to answer any questions you may have about either the information and protocols in this handbook or sample and assay technologies (for contact information, see back cover or visit www.qiagen.com).

References

QIAGEN maintains a large, up-to-date online database of scientific publications utilizing QIAGEN products. Comprehensive search options allow you to find the articles you need, either by a simple keyword search or by specifying the application, research area, title, etc.

For a complete list of references, visit the QIAGEN Reference Database online at www.qiagen.com/RefDB/search.asp or contact QIAGEN Technical Services or your local distributor.

Appendix A: Library Construction Using GeneRead Library Prep Kits for Illumina

The following kits are required for NGS library construction for Illumina MiSeq/HiSeq:

- GeneRead DNA Library I Core Kit (QIAGEN cat. no. 180432 or 180434)
- GeneRead DNA I Amp Kit (QIAGEN cat. no. 180455)
- GeneRead Adapter I Set A 12-plex (QIAGEN cat. no. 180985) or GeneRead Adapter I Set B 12-plex (QIAGEN cat. no. 180986) for up to 12 samples, or both kits for up to 24 samples.
- Agencourt AMPure XP Kit (Beckman Coulter cat. no. A63880)
- GeneRead DNaseq Library Quant Kit for Illumina (QIAGEN cat. no. 180601)

End repair of DNA

Note: 10–200 ng PCR-enriched DNA should be used for library construction. Starting amounts <10 ng or >200 ng will decrease the efficiency of library construction.

- A1.** Prepare a reaction mix for end-repair according to Table 5, dispensing the reagents into a PCR tube or the well of a 96-well PCR plate on ice.

Table 5. Reaction mix for end-repair

Component	Volume (µl)
PCR-enriched DNA from previous step	20.5
End-Repair Buffer, 10x	2.5
End-Repair Enzyme Mix	2.0
Total	25

- A2.** Mix the components by pipetting up and down several times.
- A3.** Incubate in a thermal cycler for 30 min at 25°C, followed by 20 min at 75°C.
- A4.** Pulse-spin the microfuge tube and return to ice.

Note: If reactions are to be stored during library construction, transfer them to a -20°C freezer. Samples are stable for 3 days.

A-addition

A5. Prepare a reaction mix for A-addition according to Table 6, adding the components to the PCR tube or plate containing the end-repaired DNA from step A4.

Table 6. Reaction mix for A-addition

Component	Volume
End-repaired DNA (from step A4)	25 μl
A-addition Buffer, 10x	3 μl
Klenow Fragment (3'→5' exo-)	3 μl
Total	31 μl

A6. Mix the components by pipetting up and down several times.

A7. Incubate in a thermal cycler for 30 min at 37°C , followed by 10 min at 75°C .

Note: If reactions are to be stored during library construction, transfer them to a -20°C freezer. Samples are stable for 3 days.

Adapter ligation

A8. Prepare a reaction mix for adapter ligation according to Table 7, adding the components to the PCR tube or plate containing DNA that has undergone end-repair and A-addition (step A7).

Note: When using barcode adapters, open one adapter tube at a time and change gloves between pipetting the different barcode adapters to avoid cross-contamination.

IMPORTANT: Only one of the 24 adapters (Adapter Bc1–Bc127) should be used per ligation reaction.

Table 7. Reagents for adapter ligation of PCR product

Component	Volume/reaction (µl)
DNA from step A7	31
Ligation Buffer, 2x	45
Adapter†	1*
T4 DNA Ligase	4
DNase-free water	9
Total	90

* If other suppliers' adapters will be used, use 0.28 µM final concentration or add the correct amount of adapter according to supplier's directions. If more than 24 samples will be multiplexed, NEXTflex™ DNA Barcodes from Bioo Scientific can be used.

† This applies to GeneRead Adapter I Set A 12-plex and GeneRead Adapter I Set B 12-plex

A9. Mix the components by pipetting up and down several times.

A10. Program a thermocycler to incubate at 25°C for 10 min.

IMPORTANT: Do not use a thermocycler with a heated lid.

A11. After the reaction is complete, place the reactions on ice and proceed with purification using AMPure XP beads.

Note: If reactions are to be stored during library construction, transfer them to a -20°C freezer. Samples are stable for 3 days.

Cleanup of adapter-ligated DNA with AMPure XP beads

A12. Transfer 90 µl ligation reaction from A11 to a 1.5 ml LoBind tube if using tube, keep ligation reaction in PCR plate if using plate protocol.

A13. Add 108 µl (1.2 x volume) AMPure XP beads to 90 µl DNA solution. Mix well by pipetting up and down several times.

A14. Incubate for 5 min at room temperature.

A15. Place the tube or 96-well PCR plate on magnetic rack to separate beads from supernatant. After the solution is clear (approximately 5–10 min), carefully remove and discard supernatant. Be careful not to disturb the beads, which contain the DNA target.

Note: Do not discard the beads.

- A16.** Add 200 μ l freshly made 80% ethanol to the tube or plate while it is on the magnetic rack. Rotate the tube or move the plate side-to-side in the two positions of the magnet to wash the beads. Carefully remove and discard the supernatant.
- A17.** Repeat previous step once.
- A18.** Completely remove ethanol and dry beads for 10 min while the tube or plate is on the rack.
- A19.** Elute DNA target beads in 19 μ l nuclease-free water. Mix well by pipetting. Place tube or plate on the rack until solution is clear.
- A20.** Transfer 17 μ l supernatant to a clean PCR tube or well in PCR plate and proceed to PCR amplification.
- Note:** The median size of the library will be 280 bp.
- Note:** If reactions are to be stored after library construction, transfer them to a -20°C freezer. Samples are stable for 3 days.

PCR amplification of purified library

- A21.** Mix the components in Table 8 in a 0.2 ml PCR tube or 96-well PCR plate.

Table 8. Reaction components for PCR amplification

Component	Volume/reaction (μ l)
HiFi PCR Master Mix, 2x	25
Primer Mix (10 μ M each)*	1.5
Library DNA (from step A20)	17
RNase-free water	6.5
Total	50

* Use 0.3 μ M final concentration of the PCR primer mix. Alternatively, add the correct amount of primer according to supplier's directions.

- A22.** Set up the cycler using the cycling conditions in Table 9.

Table 9. Cycling conditions for amplification of the DNA library

Step	Temperature	Time
Initial denaturation	98°C	2 min
4 cycles	98°C	20 sec
	60°C	30 sec
	72°C	30 sec
1 cycle	72°C	1 min
Hold	4°C	∞

A23. After the reaction is complete, place the reactions on ice and proceed with cleanup with AMPure XP beads.

Note: If reactions are to be stored after library amplification, transfer them to a -20°C freezer. Samples are stable for 3 days.

Cleanup of amplified library with AMPure XP beads

A24. Transfer 50 µl PCR reaction from A23 to a 1.5 ml LoBind tube if using tube, keep PCR reaction in PCR plate if using plate protocol.

A25. Add 40 µl (0.8x volume) AMPure XP beads to 50 µl PCR solution. Mix well by pipetting up and down several times.

A26. Incubate for 5 min at room temperature.

A27. Place the tube or 96-well PCR plate on magnetic rack to separate beads from supernatant. After the solution is clear (approximately 5 min), carefully transfer 86 µl of the supernatant to a new tube or well in 96-well plate without disturbing the beads. Discard the beads, which contain unwanted large DNA fragments.

Note: Do not discard the supernatant.

A28. Add 20 µl (0.4 x the original volume of PCR product, which was 50 µl) AMPure XP beads to the supernatant, mix well by pipetting and incubate for 5 min at room temperature.

A29. Place the tube or 96-well PCR plate on a magnetic rack and wait until solution is clear (approximately 5 min). Carefully remove and discard

supernatant. Be careful not to disturb the beads, which contain the DNA target.

Note: Do not discard the beads.

- A30.** Add 200 μ l fresh 80% ethanol to the tube or well while it is on the magnetic rack. Rotate the tube or move the plate side-to-side in the two positions of the magnet to wash the beads, then carefully remove and discard the supernatant.
- A31.** Repeat previous step once.
- A32.** Completely remove ethanol and dry beads for 10 min while the tube or plate is on the rack.
- A33.** Elute DNA library beads in 30 μ l nuclease-free water or appropriate buffer. Mix well by pipetting. Place tube or plate on magnetic rack until solution is clear. Transfer 28 μ l supernatant to a clean LoBind 1.5 ml tube or PCR tube.
- A34.** The library can be stored in a -20°C freezer prior to quantification using the GeneRead DNaseq Library Quant Array. Amplified libraries are stable for 1 month at -20°C .

Appendix B: Library Construction Using GeneRead Library Prep Kits for Ion PGM Sequencer/Proton

The following kits are required for NGS library construction for the Ion PGM Sequencer or Ion Proton:

- GeneRead DNA Library L Core Kit (QIAGEN cat. no. 180462)
- GeneRead DNA L Amp Kit (QIAGEN cat. no. 180485)
- GeneRead Adapter L Set 1-plex (QIAGEN cat. no. 180922) or GeneRead Adapter L Set 12-plex (QIAGEN cat. no. 180994)
- Agencourt AMPure XP Kit (Beckman Coulter cat. no. A63880)
- GeneRead DNaseq Library Quant Kit for Ion PGM Sequencer (QIAGEN cat. no. 180601)

End repair of DNA

Note: 10–200 ng PCR-enriched DNA should be used for library construction. Starting amounts <10 ng or >200 ng will decrease the efficiency of library construction.

- B1.** Prepare a reaction mix for end-repair according to Table 10, dispensing the reagents into a PCR tube or the wells of a PCR plate on ice.

Table 10. DNA end-repair reaction components

Component	Volume (µl)
PCR-enriched DNA from previous step	20.5
End-Repair Buffer, 10x	2.5
End-Repair Enzyme Mix	2.0
Total	25

- B2.** Mix the components by pipetting up and down several times.

- B3.** Incubate in a thermal cycler for 20 min at 25°C, followed by 10 min at 70°C.

- B4.** Pulse-spin the microfuge tube and return to ice.

Note: If reactions are to be stored during library construction, transfer them to a –20°C freezer. Samples are stable for 3 days.

Adapter ligation

- B5.** Prepare a reaction mix for adapter ligation according to Table 11, adding the components to the PCR tube containing the end-repaired DNA from previous step. Mix thoroughly.

Note: If analyzing 1 sample only, use the GeneRead Adapter L Set 1-plex. If analyzing up to 12 samples in a single run, use the GeneRead Adapter L Set 12-plex. If analyzing more than 12 samples, refer to table 11 for instructions.

Note: When using barcode adapters, open one adapter tube at a time and change gloves between pipetting the different barcode adapters to avoid cross-contamination.

IMPORTANT: Only one of the 12 adapters (Adapter Bc1–Bc12) should be used per ligation reaction, in combination with the universal adapter BcGen.

Table 11. Reaction setup for adapter ligation

Component	Singleplex Volume/reaction (µl)	Multiplex Volume/reaction (µl)
End-repaired DNA (from step B4)	25	25
Ligation Buffer, 2x	40	40
Adapter (singleplex)	0.5*	—
Universal Adapter BcGen	—	0.5†
Barcode Adapter 1–12	—	0.5†
Ligation and Nick Repair Mix	4	4
dNTP Mix (10 mM)	1	1
DNase-free water	9.5	9
Total	80.0	80.0

* If other suppliers' adapter will be used, use 0.16 µM final concentration or add the correct amount of adapter according to supplier's directions. † If more than 12 samples will be multiplexed, Ion Xpress™ Barcode Adapters from Life Technologies can be used. For each reaction, use 0.5 µl Ion Xpress P1 Adapter, and 0.5 µl of one of the Ion Xpress Barcodes. Only one of the Ion Xpress Barcodes is used for each sample.

- B6.** Mix the contents by pipetting up and down several times.
- B7.** Program a thermocycler to incubate for 10 min at 25°C, followed by 5 min at 72°C.

IMPORTANT: Do not use a thermocycler with a heated lid.

- B8.** After the reaction is complete, place the reactions on ice and proceed with purification using AMPure XP beads.

Note: If reactions are to be stored during library construction, transfer them to a –20°C freezer. Samples are stable for 3 days.

Cleanup of adapter-ligated DNA with AMPure XP beads

- B9.** Transfer 80 µl ligation reaction from B8 to a 1.5 ml LoBind tube if using tube, keep ligation reaction in PCR plate if using plate protocol.

- B10.** Add 112 μ l (1.4x volume) AMPure XP beads to 80 μ l DNA solution. Mix well by pipetting up and down several times.
- B11.** Incubate for 5 min at room temperature.
- B12.** Place the tube or 96-well PCR plate on magnetic rack to separate beads from supernatant. After the solution is clear (approximately 5–10 min), carefully remove and discard supernatant. Be careful not to disturb the beads, which contain the DNA target.
- Note:** Do not discard the beads.
- B13.** Add 200 μ l freshly made 80% ethanol to the tube or plate while it is on the magnetic rack. Rotate the tube or move the plate side-to-side in the two positions of the magnet to wash the beads. Carefully remove and discard the supernatant.
- B14.** Repeat previous step once.
- B15.** Completely remove ethanol and dry beads for 10 min while the tube or plate is on the rack.
- B16.** Elute DNA target beads in 19 μ l nuclease-free water. Mix well by pipetting. Place tube or plate on the rack until solution is clear.
- B17.** Transfer 17 μ l supernatant to a clean PCR tube or well in PCR plate and proceed to PCR amplification.

Note: The median size of the library will be 220 bp.

Note: If reactions are to be stored after library construction, transfer them to a -20°C freezer. Samples are stable for 3 days.

PCR amplification of the purified library

Note: This step is required to ensure that only properly-made libraries proceed to the next-generation sequencing step.

- B18.** Mix the components in Table 12 in a 0.2 ml PCR tube or plate.

Table 12. Reaction components for PCR amplification

Component	Volume (µl)
HiFi PCR Master Mix, 2x	25
Primer Mix (10 µM each)	1.5
Library DNA (from step B17)	17
RNase-free water	6.5
Total	50

B19. Set up the cycler using the cycling conditions in Table 13.

Table 13. Cycling conditions for amplification of adapter-ligated DNA

Step	Temperature	Time
Initial denaturation	98°C	2 min
5 cycles	98°C	20 sec
	60°C	30 sec
	72°C	30 sec
1 cycle	72°C	1 min
Hold	4°C	∞

B20. After the reaction is complete, place the reactions on ice and proceed with cleanup using the AMPure XP beads.

Note: If reactions are to be stored during library amplification, transfer them to a -20°C freezer. Samples are stable for 3 days.

Cleanup of amplified library with AMPure XP beads

B21. Transfer 50 µl PCR reaction from B20 to a 1.5 ml LoBind tube if using tube, keep reaction in PCR plate if using plate protocol.

B22. Add 70 µl (1.4 x volume) AMPure XP beads to 50 µl PCR reaction. Mix well by pipetting up and down several times.

- B23.** Incubate for 5 min at room temperature.
- B24.** Place the tube or 96-well PCR plate on magnetic rack to separate beads from supernatant. After the solution is clear (approximately 5 min), carefully remove and discard supernatant. Be careful not to disturb the beads, which contain the DNA target.
- Note:** Do not discard the beads.
- B25.** Add 200 μ l freshly made 80% ethanol to the tube or plate while it is on the magnetic rack. Rotate the tube or move the plate side-to-side in the two positions of the magnet to wash the beads. Carefully remove and discard the supernatant.
- B26.** Repeat previous step once.
- B27.** Completely remove ethanol and dry beads for 10 min while the tube or plate is on the rack.
- B28.** Elute DNA library beads in 30 μ l nuclease-free water or appropriate buffer. Mix well by pipetting. Place tube or plate on magnetic rack until solution is clear. Transfer 28 μ l supernatant to a clean LoBind 1.5 ml tube or PCR tube.
- B29.** The library can be stored in a -20°C freezer prior to quantification using the GeneRead DNaseq Library Quant Array. Amplified libraries are stable for 1 month at -20°C .

Appendix C: FFPE DNA Quality and Quantity

Genomic DNA present in FFPE archives is usually damaged and fragmented to an uncertain extent. Commonly used DNA quantification methods including spectrometers or fluorometers do not differentiate between amplifiable and non-amplifiable DNA. Therefore, they cannot reliably measure the amplifiable amounts of DNA that are able to participate in the multiplex PCR-based targeted enrichment step in the NGS workflow.

The QIAGEN GeneRead DNA QuantiMIZE System is a qPCR-based approach that determines the quantity and quality of DNA that is amenable to PCR-based targeted enrichment prior to NGS. The system provides a sensitive and accurate approach to qualify and quantify DNA isolated from biological samples, mainly for FFPE samples. Please refer to the corresponding user manual for determining FFPE DNA quantity and quality with GeneRead DNA QuantiMIZE System.

Appendix D: Library Quantification and Quality Control

Quality control for the targeted enrichment and library construction process can be performed using QIAGEN's GeneRead DNAseq Library Quant Array (cat. no. 180601). With this array, the correct dilution of the library can also be determined for sequencing. Please refer to the corresponding user manual for library quantification and quality control.

Appendix E: Data Analysis using QIAGEN's GeneRead DNaseq Sequence Variant Analysis Software

After sequencing, results can be analyzed using QIAGEN's Cloud-Based GeneRead DNaseq Sequence Variant Analysis Software. Our data analysis software will perform mapping to the reference genome, read trimming (removing primer sequences) and variant identification. Please refer to the corresponding document for data analysis. Alternatively, sequencing results can be analysed by the Biomedical Genomics Workbench platform. All detected variants can be interpreted by the Ingenuity Variant Analysis tool (<https://www.qiagen.com/us/products/catalog/sample-technologies/dna-sample-technologies/genomic-dna/generead-dnaseq-gene-panels-v2/#resources>).

Appendix F: Combine Libraries for Multiplex Sequencing

Libraries can be combined into one sequencing run, as long as each library uses a different barcode.

Ion PGM Sequencer libraries

Barcoded libraries can be constructed using the GeneRead Adapter L Set 12-plex. If more than 12 samples will be multiplexed, the Ion Xpress Barcode Adapters from Life Technologies (cat. nos. 4471250, 4474009, 4474518, 4474519, 4474520, 4474521, or 4474517) can be used to replace the “GeneRead Adapter L Set 12-plex” in Table 11, page 33.

After the library is constructed, follow Appendix D to determine the library dilution factor (which dilutes libraries to 4–8 pM) and dilute each individual library according to this factor. Combine libraries in equimolar amounts, and mix well. At least 25 μ l of the mixture is required. Proceed to template preparation using the mixture.

Illumina libraries

Barcoded libraries can be constructed using GeneRead Adapter I Set A 12-plex or GeneRead Adapter I Set B 12-plex as described in Appendix A, Table 7. If more than 24 samples will be multiplexed, replace the GeneRead Adapters with NEXTflex™ DNA Barcodes from Bio Scientific, as described in Appendix A, Table 7.

After the library is constructed, follow Appendix D to determine the library concentration. Dilute individual libraries to 4 nM, then combine libraries in equimolar amounts, and mix well. At least 5 μ l of the mixture is needed. Proceed to denature libraries using fresh NaOH and generate clusters using this mixture.

Appendix G: Analyze the PCR amplicons and Library Using the QIAxcel Advanced with a QIAxcel DNA High Resolution Kit or Agilent 2100 Bioanalyzer

The QIAxcel Advanced with a QIAxcel DNA High Resolution Kit or Agilent 2100 Bioanalyzer can be used at multiple steps of the targeted enrichment workflow as a quality control checkpoint. If the Agilent 2100 Bioanalyzer is used, the High Sensitivity DNA kit must be used.

After the multiplex PCR run, sample pooling, and purification, the PCR product can be analyzed using the QIAxcel Advanced or an Agilent 2100 Bioanalyzer. A sample image from both instruments is shown in Figure 7. The amplicons should be in the correct size range (usually around 160 bp), and the measured amount of DNA under the appropriate peak should be greater than 10 ng total. The OD reading (NanoDrop®) method is not recommended for determination of DNA concentration because the product may be below the detection limitation of the instrument.




Figure 7. QIAxcel (top) and Bioanalyzer (bottom) trace of pooled and purified multiplex PCR product. A peak of around 160 bp is observed, which represents the amplicons. The primer

dimers will be removed during purification steps after adapter ligation and library amplification and, therefore, will have no effect on the final library (see Figure 8).

After the library is constructed, amplified and purified, the QIAxcel Advanced or Bioanalyzer can be used to check the fragment size and concentration. For Ion PGM Sequencer libraries, a peak around 220 bp is expected. There should be no significant peak at around 100 bp, which represents adapter dimers. For Illumina libraries, a peak around 280 bp is expected (Figure 8) and no significant peak should be observed around 120 bp, which represents adapter dimers. Amounts of DNA under the appropriate peaks can be used to quantify libraries. However, due to the superior sensitivity of qPCR, we recommend quantifying libraries with the GeneRead DNAseq Library quantification kits (cat. no. 180601).




Figure 8. Sample Agilent Bioanalyzer image of a MiSeq Sequencer library. A peak of around 280 bp is observed.

Appendix H: Downloading Individual Unaligned .BAM File with a Multiplex Sample on Ion PGM Sequencer

H1. When the run is finished, navigate to the report page on the Torrent Browser. Locate the “Output Files” section near the end of the report.




Figure 9. Report page and output files.

H2. Click the “UBAM” button in the row corresponding to individual barcoded samples and column labeled “Files” in the table. These are the unaligned reads in BAM format, with the barcode separated for each sample. Save the .bam file to your local disk. The file is usually several hundred megabytes to several gigabytes, depending on the size of the sequencing chip being used.

Output Files




Figure 10. Unaligned reads in BAM format.

- H3.** Locate the file that was just downloaded to the local disk and upload the individual .BAM file to QIAGEN's NGS Sequence Variant Analysis Web Portal for analysis.

Ordering Information

Product	Contents	Cat. no.
GeneRead DNaseq Targeted Panels V2	Sets of 1 or 4 pools containing wet-bench verified primer sets for targeted enrichment of a focused panel of <100 genes	181900
GeneRead DNaseq Targeted HC Panel V2	Sets of 4 pools containing wet-bench verified primer sets for targeted enrichment of a focused panel of >100 genes	181901
GeneRead DNaseq Custom Panel V2	Pools containing primer sets for targeted enrichment of a customized panel of genes or genomic regions	181902
GeneRead DNaseq Mix-n-Match Panel V2	Pools containing wet-bench verified primer sets for targeted enrichment of a custom panel of genes	181905
GeneRead DNaseq Panel PCR Kit V2	PCR chemistry for use with the GeneRead DNaseq Panel V2 System	Varies
Related products		
GeneRead DNaseq Library Quant Array	Reagents for NGS sample library quantification following targeted enrichment with the GeneRead DNaseq Targeted Panel System	180601
GeneRead qPCR SYBR® Green Mastermix	Mastermix for use with the GeneRead Library Quant Arrays and Kit	Varies
GeneRead DNA QuantiMIZE Array Kit	qPCR arrays for optimizing amount of input DNA and PCR cycling conditions for targeted enrichment of FFPE DNA	180642

Product	Contents	Cat. no.
GeneRead DNA QuantiMIZE Assay Kit	qPCR assays for optimizing amount of input DNA and PCR cycling conditions for targeted enrichment of FFPE DNA	180654
QIAamp DNA Mini Kit (50)	For 50 DNA preps: 50 QIAamp Mini Spin Columns, QIAGEN Proteinase K, Collection Tubes (2 ml), reagents and buffers	51304
QIAamp DNA FFPE Tissue Kit (50)	For 50 DNA preps: 50 QIAamp MinElute® Columns, Proteinase K, Collection Tubes (2 ml), buffers	56404
GeneRead DNA FFPE Kit (50)	QIAamp MinElute columns, Proteinase K, UNG, Collection Tubes (2 ml), Buffers, Deparaffinization Solution, RNaseA	180134
Ion PGM Sequencer/Proton Library Prep		
GeneRead DNA Library L Core Kit (12)	For 12 reactions: Buffers and reagents for end-repair, ligation, and nick repair, for use with Ion PGM Sequencer/Proton Instruments from Life Technologies	180462
GeneRead DNA L Amp Kit (100)	For 100 reactions: Buffers and reagents for library amplification, for use with Ion PGM Sequencer/Proton Instruments from Life Technologies	180485
GeneRead Adapter L Set 12-plex (72)	For 72 reactions: 12 barcoded adapters for ligation to DNA library, for use with instruments from Life Technologies	180994
GeneRead Adapter L Set 1-plex (12)	For 12 reactions: Adapters for DNA ligation, for use with Life Technologies instruments	180922

Product	Contents	Cat. no.
Illumina Library Prep		
GeneRead DNA Library I Core Kit (48)	For 48 reactions: Buffers and reagents for end-repair, A-Addition, and ligation, for use with Illumina instruments	180434
GeneRead DNA I Amp Kit (100)	For 100 reactions: Buffers and reagents for library amplification, for use with Illumina instruments	180455
GeneRead Adapter I Set A 12-plex (144)	For 144 reactions: 12 barcoded adapters for ligation to DNA library, for use with Illumina instruments	180985
GeneRead Adapter I Set B 12-plex (144)	For 144 reactions: 12 barcoded adapters for ligation to DNA library, for use with Illumina instruments	180986

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, QIAamp®, QIAquick®, QIAxcel®, GeneRead™, HotStarTaq®, Ingenuity®, MinElute® (QIAGEN Group); Agilent®, Biaoanalyzer® (Agilent Technologies, Inc.); AMPure®, Agencourt® (Beckman Coulter, Inc.); HiSeq®, MiSeq®, Illumina® (Illumina, Inc.); Ion Xpress™, Ion Proton™, Personal Genome Machine®, PGM™, SYBR®, Life Technologies® (Life Technologies Corporation); NanoDrop® (NanoDrop Technologies, Inc.); NEBNext®, (New England BioLabs, Inc.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Limited License Agreement for GeneRead DNaseq Targeted Panels V2

Use of this product signifies the agreement of any purchaser or user of the product to the following terms:

1. The product may be used solely in accordance with the protocols provided with the product and this handbook and for use with components contained in the kit only. QIAGEN grants no license under any of its intellectual property to use or incorporate the enclosed components of this kit with any components not included within this kit except as described in the protocols provided with the product, this handbook, and additional protocols available at www.qiagen.com. Some of these additional protocols have been provided by QIAGEN users for QIAGEN users. These protocols have not been thoroughly tested or optimized by QIAGEN. QIAGEN neither guarantees them nor warrants that they do not infringe the rights of third-parties.
2. Other than expressly stated licenses, QIAGEN makes no warranty that this kit and/or its use(s) do not infringe the rights of third-parties.
3. This kit and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. QIAGEN specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the kit agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. QIAGEN may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the kit and/or its components.

For updated license terms, see www.qiagen.com.

© 2015 QIAGEN, all rights reserved.

www.qiagen.com

Australia ■ techservice-au@qiagen.com

Austria ■ techservice-at@qiagen.com

Belgium ■ techservice-bnl@qiagen.com

Brazil ■ suportetecnico.brasil@qiagen.com

Canada ■ techservice-ca@qiagen.com

China ■ techservice-cn@qiagen.com

Denmark ■ techservice-nordic@qiagen.com

Finland ■ techservice-nordic@qiagen.com

France ■ techservice-fr@qiagen.com

Germany ■ techservice-de@qiagen.com

Hong Kong ■ techservice-hk@qiagen.com

India ■ techservice-india@qiagen.com

Ireland ■ techservice-uk@qiagen.com

Italy ■ techservice-it@qiagen.com

Japan ■ techservice-jp@qiagen.com

Korea (South) ■ techservice-kr@qiagen.com

Luxembourg ■ techservice-bnl@qiagen.com

Mexico ■ techservice-mx@qiagen.com

The Netherlands ■ techservice-bnl@qiagen.com

Norway ■ techservice-nordic@qiagen.com

Singapore ■ techservice-sg@qiagen.com

Sweden ■ techservice-nordic@qiagen.com

Switzerland ■ techservice-ch@qiagen.com

UK ■ techservice-uk@qiagen.com

USA ■ techservice-us@qiagen.com

