April 2016

STAR Q Punch AS

User Manual

Version 1

9002651 (EC)

Hamilton Bonaduz AG
Via Crusch 8
7402 Bonaduz
SWITZERLAND

R1

1100450

Contents

1	Introduction	6
1.1	About this user manual	6
1.2	General information	7
1.2.1	Technical assistance	7
1.2.2	Policy statement	7
1.2.3	Version management	7
1.3	Intended use of STAR Q Punch AS instrument	7
1.4	STAR Q Punch AS card formats	8
1.5	Requirements for STAR Q Punch AS users	8
2	Safety Information	9
2.1	Proper use	10
2.2	Electrical safety	12
2.3	Biological and chemical safety	13
2.4	Mechanical hazards	14
2.5	Waste disposal	14
2.6	Maintenance safety	14
2.7	Symbols on STAR Q Punch AS	16
3	General Description	19
3.1	Platform	19
3.2	Card and plate gripper	22
3.3	Punch module	23
3.3.1	XY-table and plate slide	24
3.3.2	Punch head and actuator	24
3.3.3	Vacuum system	24
3.3.4	Ionizer	24
3.3.5	Camera and light	25
3.4	Card magazines and magazine carriers	25
3.5	Card recovery magazines	25

3.6	Plate adapter frames	26
3.7	Autoload function	26
3.8	Carriers.....	27
3.9	CO-RE gripper	28
3.10	Pipetting	29
3.10.1	Air displacement pipetting	29
3.10.2	Tip recognition	30
3.10.3	Liquid level detection.....	31
3.10.4	Monitored air displacement	31
3.11	Disposable tips.....	32
3.11.1	Tip packaging.....	32
3.12	Computer requirements	33
3.13	STAR Q Punch AS Software	34
3.13.1	easyPunchimaging software.....	34
3.13.2	Security and data protection	34
3.14	Electronics and interfaces	35
3.15	Power/voltage	36
3.16	Training.....	37
3.17	STAR Q Punch AS functionalities	37
3.17.1	Card transport.....	37
3.17.2	Card analysis.....	37
3.17.3	Punching of cards.....	37
4	Sample Processing Using STAR Q Punch AS.....	38
4.1	STR reaction PCR setup from forensic reference samples on FTA cards	39
4.1.1	Setup of the run specifications.....	39
4.1.2	Loading a run	48
4.1.3	Output files	57
4.1.4	Re-entry procedure.....	59
5	Troubleshooting	63
5.1	Error handling.....	63
5.1.1	Barcode reading error	63

5.1.2	Pipetting error	65
5.1.3	Error handling for easyPunchimaging software.....	67
5.2	Log files.....	68
6	Maintenance	69
6.1	Maintenance intervals	70
6.1.1	Preventative maintenance	70
6.2	Materials required for maintenance procedures.....	70
6.3	Maintenance procedures	71
6.3.1	Daily maintenance.....	72
6.3.2	Punch module maintenance	74
6.3.3	Card gripper and dust collector maintenance.....	75
6.3.4	Weekly maintenance	77
6.4	Printing a maintenance report	81
6.5	Instrument decontamination.....	81
6.5.1	Surface decontamination using liquids	82
6.5.2	Decontamination using gases.....	83
6.5.3	Decontamination using UV light.....	84
7	Glossary	85
	Appendix A — Ordering Information.....	88
	Appendix B — Technical Specifications	89
	STAR Q Punch AS instrument specifications	89
	Operating data.....	90
	QIAGEN standard computer specifications	91
	Labware for STAR Q Punch AS	92
	Pipetting specifications for disposable tips (individual 1000 µl channels)	93
	CO-RE gripper 1000 µl specifications.....	93
	Autoload: barcode and reader specifications	95
	Barcode symbologies	95
	Reading accuracy	95
	Barcode specifications.....	96
	Tube barcodes.....	97

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1 Introduction

Thank you for choosing the STAR Q Punch AS instrument. We are confident it will become an integral part of your laboratory.

Before using the STAR Q Punch AS instrument, it is essential that you read this user manual carefully and pay particular attention to the safety information. The instructions and safety information in the user manual must be followed to ensure safe operation of the instrument and to maintain the instrument in a safe condition.

1.1 About this user manual

This user manual provides information about the STAR Q Punch AS instrument in the following sections:

1. Introduction
2. Safety Information
3. General Description
4. Sample Processing Using STAR Q Punch AS
5. Troubleshooting
6. Maintenance
7. Glossary

The appendices include the following:

- Ordering Information
- Technical Specifications
- Instrument Dimensions
- Chemical Compatibility
- Verification
- Legal information
- Regulatory information
- Safety information in French (FR)
- Safety information in German (DE)

1.2 General information

1.2.1 Technical assistance

At QIAGEN we pride ourselves on the quality and availability of our technical support. Our Technical Services department is staffed by experienced scientists with extensive practical and theoretical expertise in sample and assay technologies and the use of QIAGEN products. If you have any questions or experience any difficulties regarding STAR Q Punch AS or QIAGEN products in general, do not hesitate to contact us.

QIAGEN customers are a major source of information regarding advanced or specialized uses of our products. This information is helpful to other scientists as well as to the researchers at QIAGEN. We therefore encourage you to contact us if you have any suggestions about product performance or new applications and techniques.

For technical assistance and more information, please see our Technical Support Center at www.qiagen.com/goto/TechSupportCenter or call one of the QIAGEN Technical Service Departments or local distributors (see back cover or visit www.qiagen.com).

1.2.2 Policy statement

It is the policy of QIAGEN to improve products as new techniques and components become available. QIAGEN reserves the right to change specifications at any time.

In an effort to produce useful and appropriate documentation, we appreciate your comments on this user manual. Please contact QIAGEN Technical Services.

1.2.3 Version management

This document is the *STAR Q Punch AS User Manual*, version 1.

1.3 Intended use of STAR Q Punch AS instrument

The STAR Q Punch AS instrument is designed to perform automated sample punching and short tandem repeat (STR) PCR setup in molecular biology applications.

The STAR Q Punch AS instrument is intended to be used only in combination with QIAGEN kits indicated for use with STAR Q Punch AS instruments for applications described in the respective QIAGEN kit handbooks.

If the STAR Q Punch AS instrument is used with other than QIAGEN kits, it is the user's responsibility to validate the performance of such product combination for any particular application.

The STAR Q Punch AS instrument is intended for use by professional users, such as technicians and physicians trained in molecular biological techniques and the operation of the STAR Q Punch AS instrument.

1.4 STAR Q Punch AS card formats

The STAR Q Punch AS (EC) instrument with Autoload including 4 x 1000 µl channels, card and plate gripper and punch module can be used with FTA™ cards and rigid frames. (FTA is an acronym for fast technology for analysis of nucleic acids.) The following formats are supported:

Manufacturer	Description	Dimensions
GE Healthcare	EasiCollect™ card	51 mm x 51 mm
COPAN®	NUCLEIC-CARD™ collection card	51 mm x 51 mm

1.5 Requirements for STAR Q Punch AS users

This table covers the general level of competence and training necessary for transportation, installation, use, maintenance and servicing of the STAR Q Punch AS instrument.

Task	Personnel	Training and experience
Delivery	No special requirements	No special requirements
Installation	QIAGEN service personnel or service technicians of an authorized agent	Trained and authorized by QIAGEN
Routine use	Laboratory technicians or equivalent	Appropriately trained and experienced personnel familiar with use of computers and automation in general
Preventive maintenance	QIAGEN service personnel or service technicians of an authorized agent	Trained and authorized by QIAGEN
Servicing	QIAGEN service personnel or service technicians of an authorized agent	Trained and authorized by QIAGEN

2 Safety Information

Before using the STAR Q Punch AS instrument, it is essential that you read this user manual carefully and pay particular attention to the safety information. The instructions and safety information in the user manual must be followed to ensure safe operation of STAR Q Punch AS and to maintain the instrument in a safe condition.

Note: Translations of the Safety Information in French and German are available in Appendix J – Safety Information (French, FR) and Appendix K – Safety Information (German, DE).

The following types of safety information appear throughout this user manual.

WARNING

The term WARNING is used to inform you about situations that could result in **personal injury** to you or others.

Details about these circumstances are provided to avoid personal injury to you or other persons.

CAUTION

The term CAUTION is used to inform you about situations that could result in **damage to the instrument** or other equipment.

Details about these circumstances are provided to avoid damage to the instrument or other equipment.

The advice given in this manual is intended to supplement, not supersede, the normal safety requirements prevailing in the user's country.

2.1 Proper use

CAUTION

Loss of data

Do not open the front cover during a run. An aborted run, that is one stopped by opening the front cover, cannot be recovered. To open the window during a run, click **Pause** in the run screen, wait until the instrument stops and then open the window.

CAUTION

Loss of data and damage to the instrument

Use only labware defined in this manual with the STAR Q Punch AS instrument. Failure to do so may result in damage to the instrument and incorrect results.

WARNING/ CAUTION

Risk of personal injury and material damage

Improper use of STAR Q Punch AS may cause personal injuries or damage to the instrument. STAR Q Punch AS must only be operated by appropriately trained and experienced personnel.

Servicing of STAR Q Punch AS must only be performed by QIAGEN Field Service Specialists or service technicians of an authorized agent.

WARNING

Magnetic field

The magnetic plate slide used with the STAR Q Punch AS instrument emits a strong magnetic field that can be harmful to pacemaker wearers.

Pacemaker wearers should maintain a minimum distance of 50 cm (20 in.) from the plate slide at all times.

CAUTION

Damage to the instrument

Direct sunlight may bleach parts of the instrument and cause damage to parts. STAR Q Punch AS must be located out of direct sunlight, away from heat sources and away from sources of vibration and electrical interference.

CAUTION**Damage to the instrument**

Avoid spilling water or chemicals onto STAR Q Punch AS. Damage caused by water or chemical spillage will void your warranty.

**WARNING/
CAUTION****Risk of personal injury and material damage**

Do not attempt to move STAR Q Punch AS during operation.

Never lift a fully installed instrument for transportation from one location to another. It must be reinstalled in the new location by an authorized service engineer.

**WARNING/
CAUTION****Risk of personal injury and material damage**

STAR Q Punch AS weighs approximately 135 kg (297 lbs). Necessary safety precautions should be taken when transporting the instrument.

**WARNING/
CAUTION****Explosive atmosphere**

STAR Q Punch AS is not designed for use in an explosive atmosphere.

WARNING**Risk of explosion**

STAR Q Punch AS is intended for use with reagents supplied with QIAGEN kits. Use of other reagents and substances may lead to fire or explosion.

In case of emergency, switch off STAR Q Punch AS at the power switch and unplug the power supply from the power outlet.

The instrument should be positioned in the laboratory in a way permitting personnel to access the front and sides of the instrument in order to operate, maintain, open and close the protective covers. Consider the dimensions of the instrument (see "Appendix B — Technical Specifications," page 89) and calculate sufficient room for a person to move and work comfortably.

The STAR Q Punch AS instrument is equipped with dark panels to reduce the influence of ambient light. Direct exposure of the instrument's front towards light sources (artificial or natural) must be avoided.

The optical path must not be obstructed by the operator, specifically during the loading phase, to ensure detection reliability.

2.2 Electrical safety

Note: Disconnect the line power outlet before servicing.

WARNING

Electrical hazard

Any interruption of the protective conductor (earth/ground lead) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument dangerous.

Intentional interruption is prohibited.

Lethal voltages inside the instrument

When the instrument is connected to line power, terminals may be live.

Opening covers or removing parts is likely to expose live parts.

Avoid spilling liquid onto or into the instrument. In case of spilling liquid over the instrument, immediately disconnect the instrument from the mains power.

To ensure satisfactory and safe operation of STAR Q Punch AS, follow the guidelines below:

- The line power cord must be connected to a line power outlet that has a protective conductor (earth/ground).
- Keep mains plug easily accessible in case the equipment needs to be disconnected quickly from mains power.
- Use only power supplies and cords supplied with the system.
- If the instrument becomes electrically unsafe, prevent other personnel from operating it and contact QIAGEN Technical Services. The instrument may be electrically unsafe when:
 - The line power cord appears to be damaged.
 - It has been stored for a prolonged period of time in conditions which are outside of the storage conditions outlined in "Appendix B — Technical Specifications," page 89.
 - It has been subjected to severe transport stresses.
 - Liquid has entered the instrument.

2.3 Biological and chemical safety

When handling biological material, use safe laboratory procedures as outlined in publications such as Biosafety in Microbiological and Biomedical Laboratories, HHS (<http://www.cdc.gov/biosafety/publications/bmbl5/index.htm>).

WARNING

Biological materials

Handle biological material with the greatest of care and in accordance with the required safety regulations. Always wear safety glasses, 2 pairs of gloves and a lab coat.

The responsible body (e.g., laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe and that the instrument operators are not exposed to hazardous levels of infectious agents as defined in the applicable Safety Data Sheets (SDSs) or OSHA,* ACGIH† or COSHH‡ documents.

For more information, visit www.qiagen.com/safety.

Venting for fumes and disposal of wastes must be in accordance with all national, state and local health and safety regulations and laws.

WARNING

Hazardous chemicals

Always wear safety glasses, gloves and a lab coat.

The responsible body (e.g., laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe and that the operators are not exposed to hazardous levels of toxic substances (chemical or biological) as defined in the applicable Safety Data Sheets (SDSs) or OSHA,* ACGIH† or COSHH‡ documents.

For more information, visit www.qiagen.com/safety.

Venting for fumes and disposal of wastes must be in accordance with all national, state and local health and safety regulations and laws.

* OSHA: Occupational Safety and Health Administration (United States of America).

† ACGIH: American Conference of Government Industrial Hygienists (United States of America).

‡ COSHH: Control of Substances Hazardous to Health (United Kingdom).

2.4 Mechanical hazards

WARNING

Moving parts

To avoid contact with moving parts during operation of STAR Q Punch AS, the instrument must be operated with the cover closed.

Do not remove the cover panels since there are no user-serviceable parts inside. If there is a problem with STAR Q Punch AS, contact QIAGEN Technical Services immediately.

WARNING

Risk of personal injury

Do not touch the LED light source during run time and for 1 hour after finishing a run as it might be hot.

2.5 Waste disposal

CAUTION

Disposal of plasticware

Used plasticware may contain hazardous chemicals, or contagious/biohazardous materials. Such wastes must be collected and disposed of properly according to local safety regulations.

2.6 Maintenance safety

Perform the maintenance as described in Section 6. QIAGEN charges for repairs that are required due to incorrect maintenance.

WARNING/ CAUTION

Risk of personal injury and material damage

Only perform maintenance that is specifically described in this user manual.

WARNING/ CAUTION

Risk of electric shock

Do not open any panels on the STAR Q Punch AS instrument.

Only perform maintenance that is specifically described in this user manual.

CAUTION**Damage to the instrument**

Do not use solvents or reagents containing acids, alkalis or abrasives to clean the STAR Q Punch AS instrument. Do not use disinfecting materials which contain hypochlorite or other bleaching fluids. Use non-corrosive, neutral liquids.

CAUTION**Damage to the instrument**

Autoclaving cannot be used for instrument components or accessories (sample carriers, magazines, plate slide or punch heads).

CAUTION**Damage to the instrument**

Ethylene oxide fumigation may increase service and maintenance requirements (O-rings exchange, greasing of spindles, etc.) and may make shorter maintenance intervals necessary.

CAUTION**Damage to the instrument**

Hydrogen peroxide fumigation leads to bleaching or discoloration of many instrument materials and may increase service and maintenance requirements (O-rings exchange, greasing of spindles, etc.) and may make shorter maintenance intervals necessary.

CAUTION**Damage to the instrument**

Do not use formaldehyde fumigation or chlorine oxides (chemical compounds of chlorine and oxygen such as bleach). They are not suitable for the STAR Q Punch AS instrument because of chemical reaction and corrosion.

CAUTION**Damage to the instrument**

UV irradiation causes many synthetic materials to become brittle. This may increase service and maintenance requirements and may make shorter maintenance intervals necessary.


2.7 Symbols on STAR Q Punch AS

The following symbols may be found on the instrument or in this user manual.

Symbol	Location	Description
	Type plate on the instrument	Manufacturer
	Type plate on the instrument	Waste Electrical and Electronic Equipment (WEEE)
	Type plate on the instrument	FCC mark of the United States Federal Communications Commission
	Type plate on the instrument	RoHS mark for China (the restriction of the use of certain hazardous substances in electrical and electronic equipment)
	Type plate on the instrument	RCM mark for Australia/New Zealand
	Type plate on the instrument	Serial number
	Type plate on the instrument	Global trade item number
	On the instrument	General warning sign
	On the instrument	Warning, dangerous voltage

Symbol	Location	Description
	On the LED light source	Warning, hot surface
	On the barcode reader	Warning, laser
	On the instrument	Warning, biological hazard
	In this user manual	Warning, strong magnetic field
	In this user manual	Warning, no access for persons with pacemaker

Location and explanation of warning and attention labels:


3 General Description

The STAR Q Punch AS instrument is designed for punching and processing of sample collection cards. The hardware and software are capable of handling cards, identifying samples on the cards and punching the filter paper pieces into a 96-well PCR plate. The STAR Q Punch AS instrument uses the general STAR Q resources such as the channels, CO-RE gripper and Autoload.

3.1 Platform

The STAR Q Punch AS instrument work surface, the deck, is used for placing loadable carriers. The carriers hold reagent containers, microplates, filter tips and card magazines.

- | | |
|--|--|
| 1 Card and plate gripper (mounted on the left side of the arm; refer to Section 3.2) | 6 Carriers: Card magazines left, microplates right |
| 2 Pipetting channels (mounted on the right side of the arm) | 7 Autoload drive and barcode scanner |
| 3 Front cover frame and window | 8 Waste station |
| 4 Punch module (refer to Section 3.3) | 9 Pipetting axes and coordinates (x, y, z) |
| 5 Deck | |

CAUTION

Loss of data

Do not open the front cover during a run. An aborted run, that is one stopped by opening the front cover, cannot be recovered. To open the window during a run, click **Pause** in the run screen, wait until the instrument stops and then open the window.

The instrument is fully covered by an acrylic glass hood. The front cover consists of a hinged transparent window made of acrylic glass. The window is equipped with a magnetic switch that is monitored during a run. Opening the cover will abort the run.

The instrument deck is divided into 30 equal tracks (T) for loading carriers in predetermined positions. This eliminates the need for precise measurement of positions. The deck has partitions of 22.5 mm, equivalent to 1-T. The labware carriers are adapted to the partitions: 6-T carriers for microplates or Compression-induced O-Ring Expansion (CO-RE) tips and other labware. An additional partition provides space for the tip waste. The deck partitions accommodate a maximum of 30 1-T carriers for sample tubes (when applicable), or a maximum of five 6-T carriers for microplates and CO-RE tips. A total of 25 SBS (Standard format of the Society for Biomolecular Screening) positions can be placed onto the STAR Q Punch AS deck.

The instrument's internal coordinate (x, y, z) system is shown in the picture above, located at its origin. Please note that the ZERO position is 100 mm below the metal deck sheet.


Note: The punch module is serviced from the left of the system. Please make sure that there is about 1 m of space to the left of the instrument for servicing.

The STAR Q Punch AS instrument has the following features:

- Card and plate gripper specific for the cards used (see "STAR Q Punch AS card formats," page 8)
- Punch module holding the punch head
- STAR Q Punch AS Autoload
- Four 1000 µl pipetting channels on a modular pipetting head MPH/iSWAP® arm

Magazines specific to card type, the corresponding card recovery magazines and carriers, multiflex MFX modules and additional carriers are provided.

The card and plate gripper is mounted on the STAR Q Punch AS arm and requires the MPH position on this arm. The punch module requires tracks 1–6 on the left of the instrument plus space in the negative T area. The instrument carriers occupy positions 1–4 (depicted with blue squares).

- | | |
|--------------------------|--------------------------------------|
| 1 MPH arm (not shown) | 5 Carrier positions 1–4 |
| 2 Card and plate gripper | 6 Carrier |
| 3 Punch module | 7 Autoload drive and barcode scanner |
| 4 Deck | 8 Waste station |

The STAR Q Punch AS magazines are designed to hold 20 cardboard framed cards per magazine. Card-specific magazines are loaded on a corresponding magazine carrier or MFX module on a multiflex carrier. Magazine carriers are designed to allow for barcode reading on the magazines from the right side. The barcodes must follow the specifications described in “Autoload function,” page 26, and in “Appendix B — Technical Specifications,” page 89.

Card recovery magazines are used to recover cards that generated errors during automated processing and require manual inspection.

3.2 Card and plate gripper

The card and plate gripper is mounted in the MPH position of the arm. The card and plate gripper moves in the XYZ-position to pick up cards or plates and position them on the punch module.

Note: Do not touch the pointed ends of the card and plate gripper fingers as they are a potential source of injury.

Note: Bending of gripper fingers can impair performance of the instrument.


Star Q Punch AS card and plate gripper.

3.3 Punch module

The punch module is mounted on the left side of the system. It occupies tracks 1–6 plus space in the negative T area. The module consists of 5 different parts:

- A magnetic XY-table with a plate slide
- An exchangeable punch head with punch actuator
- A vacuum system
- An ionizer
- A camera and light

- 1 XY-table and plate slide
- 2 Punch head and actuator
- 3 Vacuum

- 4 Ionizer
- 5 Camera and LED Light source

3.3.1 XY-table and plate slide

The XY-table contains linear motors for moving the plate slide and a light foil that is used as a backlight for imaging of the card and the plate.

WARNING

Magnetic field

The magnetic plate slide used with the STAR Q Punch AS instrument emits a strong magnetic field that can be harmful to pacemaker wearers.

Pacemaker wearers should maintain a minimum distance of 50 cm (20 in.) from the plate slide at all times.

Note: Handle the plate slide with care! Bending of the plate slide can impair performance of the instrument.

The plate slide is used to move plates on the XY-table using magnetic connections. It contains two positions for SBS format plates up to 43 mm high plus two containers for cleaning punches.

3.3.2 Punch head and actuator

The standard punch head used has a punch size of 1.2 mm. The punch head can be released easily. All punch heads contain an ejection pin that allows the punch to be ejected properly.

The punch head and the ejection pin are driven by the punch actuator.

3.3.3 Vacuum system

The vacuum system is connected to the punch head and removes any dust generated during the punching procedure.

3.3.4 Ionizer

The ionizer efficiently removes static forces from the plates to reduce the punch's attraction to static charges on the plates.

3.3.5 Camera and light

An industrial grade camera and two LEDs as a light source are mounted on top of the punch module. The camera is used for taking pictures of cards and plates for later analysis.

WARNING

Risk of personal injury

Do not touch the LED light source during run time and for 1 hour after finishing a run as it might be hot.

3.4 Card magazines and magazine carriers

Card magazines are used for loading cards onto the card and plate gripper. The card gripper can pick up cards from the magazines and replace them after punching. Magazines and carriers are provided for the EasiCollect and NUCLEIC-CARD formats (51mm x 51mm) (see "STAR Q Punch AS card formats," page 8).

Note: Performance of the cards depends on the quality of the cards. Bent or kinked cards can reduce the robustness of the process.

One magazine holds 20 cards. Predefined carriers and Multiflex modules are provided for each magazine. One carrier holds 5 magazines, one MFX module holds one. The number and position is predefined on the deck layout.

Note: Use only magazines developed for the STAR Q Punch AS instrument. Use of any other magazine can impair performance and will lead to expiry of warranty of the instrument.

Note: Before use, remove any covers from the magazine!

3.5 Card recovery magazines

If processing of cards leads to errors, cards are automatically recovered in recovery magazines for manual processing. These magazines are universal for all card types. Cards cannot be picked up from card recovery magazines by the instrument. Specific card recovery magazines are supplied with dedicated MFX modules.

3.6 Plate adapter frames

Specific adapter frames are used to process semi-skirted PCR plates (e.g., plates compatible with Life Technologies thermal cyclers). For plate types to be used, see "Labware for STAR Q Punch AS," page 92.



Plate adapter frame and semi-skirted PCR plate.

3.7 Autoload function

The Autoload drive is a device enabling automatic loading of carriers onto the STAR Q Punch AS instrument deck. The features of the Autoload are:

- Moves in X-direction,
- Shunts carriers onto and off the deck
- Reads the barcodes of microplates and tip trays in carriers

The Autoload option including the barcode reader is shown reading horizontally (1) for plates and vertically (2) for carriers.

Carrier identification by barcodes and reading of barcodes on plates is possible in combination with the Autoload drive. Barcodes of tip trays and racks are read automatically when loading tips to exclude the risk of mixing up tip types. This method is suitable for all tip types.

The following barcode types can be read by the Autoload:

- ISBT Standard
- Code 128 (Subset B and C)
- Code 39
- Codabar
- Code 2 of 5 interleaved
- UPC A

Barcodes must be black bars on white background. We recommend using the barcode type Code 128 (Subset B and C).

Note: Barcodes must have a minimum readability (i.e., good contrast, size, correct orientation and distance between bars) to be fully functional. Make sure that the barcode is in the correct orientation for tubes (when applicable) and plates. For details of barcodes, see the specifications given in "Autoload: barcode and reader specifications," page 95.

3.8 Carriers


Labware such as card magazines and plates are placed on a carrier loaded onto the STAR Q Punch AS instrument deck. Standard carriers are loaded to the deck by the Autoload drive.

The Multiflex carrier consists of a multiple-use carrier base offering space for up to 5 modules. The modular design of the Multiflex carrier allows optimization of space.

3.9 CO-RE gripper

The CO-RE gripper is a plate handling tool picked up by two pipetting channels. On the STAR Q Punch AS instrument, the CO-RE gripper is used during load check. During the run, the gripper on the card and plate gripping module is used for plate transport.

The traverse height of the channels with the gripping jaws is the same as with tips: 145 mm above the deck.

CO-RE gripper on the pipetting head.

The holder of the two gripping paddles is on the waste block.

Gripping paddles in park position in the waste block.

3.10 Pipetting

The STAR Q Punch AS instrument has four 1000 µl pipetting channels working in parallel for simultaneous transfer of liquids. Each pipetting channel moves independently on the Y-axis, as well as the Z-axis. Each channel uses its own high-precision motors and electronics to reach the pipetting position. The minimum distance between two 1000 µl pipetting channels on the arm is 9 mm.

1000 µl pipetting channel.


Note: The system reports an error if pipetting positions outside the access range of the instrument are used.

The four 1000 µl channels support pipetting with disposable 50 µl CO-RE tips.

The 1000 µl pipetting channels have a set “traverse height” of 245 mm above the origin, or 145 mm between the top of the disposable tip and the deck of the instrument. This means that when a channel moves from one location on the deck to another, it automatically does so at that particular height. This is a safety precaution to prevent the channels from colliding with items that may be on the deck.

3.10.1 Air displacement pipetting

The pipetting operation is based on the air displacement pipetting principle, comparable to the functioning of handheld pipets. Air displacement means that the liquid is aspirated into and dispensed from a disposable tip by the movement of a plunger. Only air is present between the plunger and the liquid surface and no system liquid of any kind is involved in the instrument.

- | | |
|------------------|-----------------|
| 1 Plunger | 4 Air |
| 2 Barrel | 5 Liquid |
| 3 Tip | |

Pipetting with the STAR Q Punch AS instrument is identical to pipetting with a handheld plunger pipette. The pipetting head's barrel and plunger are not intended to be cleaned.

3.10.2 Tip recognition

The tips used in a pipetting procedure are matched with the pipetting channel to prevent damage to the instrument. The STAR Q Punch AS instrument uses two methods for recognition of the tip type:

- Tip trays and racks have color-coded barcode labels for automated recognition
- The color and text on the barcode label distinguishes the tip volume by eye

Tip volume	Barcode label	Color code
50 µl		Light red

3.10.3 Liquid level detection

The liquid level in a container can be detected before aspiration of the liquid by the liquid level detection (LLD) feature of the STAR Q Punch AS instrument. LLD is based on either capacitive (cLLD) or pressure (pLLD) signal detection. Usually cLLD is used for conductive liquids. The sensitivity of the cLLD depends on the vessel size, volume and the conductivity (or polarity) of the liquid to be detected.

For non-conductive liquids, or where there is an insufficient coupling between container bottom and carrier, pLLD is used; pLLD only works with new and empty tips for the aspiration of liquids and is available on the individual pipetting channels only.


In cases of LLD under demanding circumstances, e.g., foaming liquids, cLLD and pLLD can be used at the same time.

3.10.4 Monitored air displacement

STAR Q Punch AS is equipped with an aspiration monitoring feature. During the aspiration process, the pressure within the pipetting channel is measured in real time. By analyzing the shape of the $p(t)$ curve, the system can distinguish the following situations:

- A correct aspiration had taken place
- Air is aspirated into the tip (because, for example, the container has not been filled properly)
- A clot has blocked the tip

The diagram shows the functioning of monitored air displacement (MAD) based on pressure.

Aspiration monitoring based on pressure.

The MAD feature is available on the individual pipetting channels only.

3.11 Disposable tips

CO-RE tips are the only tips designed to work with the STAR Q Punch AS instrument. Channels lock on to the tip rather than the standard method of forcing a tip on a mandrel. This design assures superior tip alignment, seals the tip to the pipette channel, eliminates tip distortion and mitigates aerosol contamination.

The use of non-CO-RE tips is not supported on the STAR Q Punch AS instrument. Non-CO-RE tips do not have the annular space to receive the O-ring. Use of non-CO-RE tips may result in unintended deformation of the O-ring and the tip material itself. It may lead to misalignment of the tip, improper sealing, reduced O-ring life, improper drop-off and potential aborted runs from randomly dropped tips. Therefore, QIAGEN does not support the use of non-CO-RE tips on the STAR Q Punch AS instrument.

Important: Only CO-RE disposable tips should be used with the pipetting channels of the STAR Q Punch AS instrument. Other tips may cause contaminated or damaged pipetting channels.

Note: Pipetting specifications are only guaranteed when using CO-RE tips.

Disposable CO-RE tips are produced under clean room conditions (class 8), based on ISO 14644 standards. "Biological purity tested" tips are free of DNA, DNase/RNase, PCR inhibitors and endotoxin (non-pyrogenic). In addition to these criteria, "Biological purityPLUS" tips are sterile, according to ISO 11135, and free of ATP.

- Filter tips prevent aerosol contamination.
- Conductive (black) tips are specifically designed for cLLD

3.11.1 Tip packaging

All CO-RE disposable tip types are available in trays with a sealed paper lid (blister pack). The tip trays are barcode labeled for automatic identification during the loading process. QIAGEN instruments check for loading of the correct tip type via the barcode reader.

The low-volume CO-RE tip, 1–50 µl, is available as a conductive (black) tip for use with cLLD, with filter for use with 1000 µl channels. The 1–50 µl tip is compatible with the MAD feature.

CO-RE tips, 1–50 µl, come in trays of 96 tips. One blister pack contains 5 trays.

The trays are compatible with the tip carriers.

Multiflex carrier with two modules holding card recovery magazines and three tip rack modules and tip trays.

3.12 Computer requirements

The STAR Q Punch AS instrument is controlled by dedicated STAR Q Punch AS Software which manages all functions of the daily work routine.

QIAGEN offers standard computers for use with the software. For specifications of the QIAGEN standard computers, refer to "QIAGEN standard computer specifications," page 91.

Note: Computer requirements as well as operating systems described in this manual are subject to change without notice.

The STAR Q Punch AS instrument configuration and the processor board implemented have the following interfaces for linking the instrument to the control computer:

- Serial interface RS-232C with dual processor board
- Unified Serial Bus (USB) interface with dual processor board and LAN dual processor board
- Ethernet with LAN dual processor board

The communication interface used on the computer is set by the configuration editor. For further information about the recommended computer model refer to "QIAGEN standard computer specifications," page 91 and "Electronics and interfaces," page 35.

To avoid data loss, use of an uninterruptible power supply (UPS) for the computer is recommended.

3.13 STAR Q Punch AS Software

The STAR Q Punch AS method is based on the STAR Q Punch AS Software. All pipetting parameters, aliquoting steps, labware configurations, consumables management, etc. are configured and handled by the STAR Q Punch AS Software.

3.13.1 easyPunchimaging software

The easyPunchimaging software is used by the STAR Q Punch AS imaging method for image processing. In routine operation, the easyPunchimaging software appears briefly during analysis.

If additional support is required with the application, the easyPunchimaging software can export images and settings so that they can be sent to QIAGEN Technical Services for support.

3.13.2 Security and data protection

Use the necessary precautions against software viruses. Use only manufacturer's original installation DVD/CD-ROM sets for installation of the operating system.

Running other software in parallel to the STAR Q Punch AS Software may negatively affect the operation of the STAR Q Punch AS instrument.

Any manipulation to the software data files or other information determining or affecting the functions of the STAR Q Punch AS instrument can result in erroneous test results or instrument failure.


Only the STAR Q Punch AS Software may be used to control the instrument.

For reasons of data security and integrity, the use of an uninterrupted power supply (UPS) is recommended, since a loss of power may cause data to be lost or corrupted.

To avoid computer breakdowns, configure a hard disk of sufficient space in the computer. Make sure that there is always enough storage capacity on the hard drive. Delete the log files from time to time. Generated data within the log files directory, e.g., traces, pipetting files and the image files generated during the run and stored in the imaging database should be backed up on the laboratory's host device and deleted from the control computer's hard disk at weekly intervals. See "Log files," page 68, for more information.

3.14 Electronics and interfaces

All the electrical connections are placed on the left side of the instrument, as shown below:

- | | |
|----------------|------------------------|
| 1 TCC1 | 7 Extensions I/O |
| 2 TCC2 | 8 DIV 2 |
| 3 External CAN | 9 DIV 3 |
| 4 Power 1 | 10 DIV 1 |
| 5 Power 2 | 11 USB communication |
| 6 Main power | 12 RS232 communication |

The main power connection is located near the front of the instrument (on the right side in the picture above). The communication connections to the computer are positioned beneath this main power connection. The STAR Q Punch AS instrument can communicate either via USB (the preferred option) or via RS232.

Important: Never use both connections, USB and RS232, together!

Note: Do not attempt to install a 3rd party device via the electrical connections of the STAR Q Punch AS instrument on your own! Always consult a local QIAGEN representative for installation.

The connectors labeled “Power 1” and “Power 2” provide different power supply voltages. There is also a CAN bus for communication.

“TCC1” and “TCC2” are not applicable in this configuration.

The connectors “DIV 1”, “DIV 2”, “DIV 3” and “Extensions I/O” deliver several digital input/output signals as well as pulse-width-modulated (PWM) outputs, CAN bus and the TTL levels.

3.15 Power/voltage

Make sure that the instrument is connected to a 100/115/230 VAC (50 or 60 Hz) socket. The instrument can be used with any voltage within the range 100–230 VAC. The main plug is on the left-hand side of the instrument at the front. We recommend using an uninterruptible power supply (UPS) for the instrument. Plug the main cables for the computer and the instrument into the same electrical outlet. Connect them only to a grounded outlet.

Ensure that the instrument is correctly grounded when connected to the electrical outlet.

The fuses for the instrument are placed in the main power socket (see picture below). During installation/IQ, the appropriate fuse is selected by the QIAGEN Field Service Specialist.

Note: When replacing a fuse make sure to use the appropriate fuse (see “Appendix B — Technical Specifications,” page 89) and place it in the main power switch before switching on the instrument.

3.16 Training

Training in the operation of the STAR Q Punch AS instrument and general use of the STAR Q Punch AS Software will be provided by QIAGEN service personnel or service technicians of an authorized agent at initial installation and setup.

3.17 STAR Q Punch AS functionalities

3.17.1 Card transport

The card and plate gripper is used to pick up cards from the specific magazines and move them around on the deck, e.g., to the imaging position or to the punching position.

3.17.2 Card analysis

The camera is used for taking pictures for analysis by the software. This analysis includes barcode reading and sample identification, as well as determination of the correct punching position.

3.17.3 Punching of cards

The system punches cards that are placed into the punch head. The punch size is predefined with a diameter of 1.2 mm per punch.

To remove potential contamination from the punch head, cleaning punches can be performed either on cards containing sample or on separate sample-free cards. (These cards are kept in a specific magazine and are reused until no punching position is left). If the cleaning punch is to be done on a card containing sample, the software determines sample-free areas on the card suitable for cleaning punches.

4 Sample Processing Using STAR Q Punch AS

The purpose of the QIAGEN STAR Q Punch AS Graphical User Interface (GUI) described in this user manual is to perform the following sample processing procedure:

- STR reaction PCR setup from forensic reference samples on FTA cards

Minimum: 1 reaction

Maximum: 180 reactions

All listed workflows are in association with QIAGEN reagents. Refer to “Appendix A — Ordering Information,” page 88 to order kits.

Important notes

- Do not exchange positions of sample and reagents or switch microplates after they have been identified by the instrument. This could result in incorrect test data or an instrument crash.
- Microplates must be placed on the carrier such that well A1 is located in the top-left position.
- When handling liquid, ensure that there is no foam on the surface of the liquid. Note that foam may cause pipetting problems.
- Do not overfill containers.
- Do not mix tip size and type (e.g., with or without filter, or different volumes) in the same tip rack. Take care if using tips which cannot be distinguished by the tip recognition feature (refer to “Tip recognition,” page 30).
- Do not fill partially consumed tip racks with tips from other racks. Tip racks should be loaded into the tip rack carriers as they are provided in the original package. The tip racks are individually labeled with a barcode for identification.
- Do not try to open the front cover of the instrument during a run because the system will abort and this will cause a loss of data.
- When the system is paused, do not wait too long before resuming the run. Loss of liquid from a full tip may result in invalid data.
- Discard used tips. Do not reuse them.
- Do not empty the tip waste during a run.
- Do not leave tips on the pipetting channels for a long period of time, e.g. overnight. This may cause damage to the CO-RE O-rings. A daily maintenance procedure will remove the tips.

4.1 STR reaction PCR setup from forensic reference samples on FTA cards


Note: Refer to information in the appropriate kit handbook before starting a run.

4.1.1 Setup of the run specifications

1. Launch the QIAGEN user interface from the computer desktop.

2. Enter the run specifications by choosing **PCR Setup**.

3. In the **Labware** panel, select the FTA card type from the drop-down list.
For a list of available FTA card types, see "STAR Q Punch AS card formats," page 8.

Labware	
Card type	GEHC EasiCollect
Plate type	GEHC EasiCollect
Use card barcode	<input type="checkbox"/>
Use magazine barcode	<input type="checkbox"/>
Use plate barcode	<input type="checkbox"/>

4. In the **Labware** panel, select the labware type for the PCR output plate.

Labware	
Card type	GEHC EasiCollect
Plate type	ABI MicroAmp 96well PCR
Use card barcode	ABI MicroAmp 96well PCR
Use magazine barcode	Bio-Rad 96well PCR
Use plate barcode	<input type="checkbox"/>


Note: The system is configured for the following PCR output plates:

- “Bio-Rad® 96well PCR” plate (Bio-Rad Hard-Shell® 96-well PCR plate, cat. no. HSP9901)
- “ABI MicroAmp® 96well PCR” plate (Applied Biosystems® MicroAmp 96-well plate, cat.no. N8016154)

Note: The half-skirted Applied Biosystems MicroAmp 96-well plate must be used in combination with an adapter (see “Plate adapter frames,” page 26).

Plate adapter frame.

Left: “Bio-Rad 96well PCR” plate; right: “ABI MicroAmp 96well PCR” mounted on the adapter.

5. In the **Labware** panel, check the boxes to use barcodes on the FTA sample cards, the FTA card magazines or the PCR output plates.

Labware	
Card type	GEHC EasiCollect
Plate type	ABI MicroAmp 96well PCR
Use card barcode	<input checked="" type="checkbox"/>
Use magazine barcode	<input checked="" type="checkbox"/>
Use plate barcode	<input checked="" type="checkbox"/>

6. In the **Punch and Plate Setup** panel, check the box to add the positive and negative controls after the samples. Dedicated positions for the positive and negative controls can be entered manually.

Punch and Plate Setup	
Add controls after samples	<input checked="" type="checkbox"/>
Positive control positions	F11:F11
Negative control positions	G11:G11
Ladder positions	H2:H4:H6:H8:H10:H12
Create duplicate of plates	<input checked="" type="checkbox"/>
Number of punches per card and plate	1
Number of cleaning punches	3
Use separate cards for cleaning punches	<input checked="" type="checkbox"/>

Note: When more than one PCR output plate is processed, the positions for the positive and negative controls must be entered for all PCR output plates as shown above. The positions for the positive and negative controls on the different PCR output plates may differ.

7. In the **Punch and Plate Setup** panel, manually enter the positions for the ladder; these positions are left empty during the PCR setup.

Punch and Plate Setup	
Add controls after samples	<input checked="" type="checkbox"/>
Positive control positions	F11:F11
Negative control positions	G11:G11
Ladder positions	H2:H4:H6:H8:H10:H12
Create duplicate of plates	<input checked="" type="checkbox"/>
Number of punches per card and plate	1
Number of cleaning punches	3
Use separate cards for cleaning punches	<input checked="" type="checkbox"/>

8. **Optional:** In the **Punch and Plate Setup** panel, check the box to set up the PCR in duplicate on identical independent PCR output plates.

Punch and Plate Setup	
Add controls after samples	<input checked="" type="checkbox"/>
Positive control positions	F11:F11
Negative control positions	G11:G11
Ladder positions	H2:H4:H6:H8:H10:H12
Create duplicate of plates	<input checked="" type="checkbox"/>
Number of punches per card and plate	1
Number of cleaning punches	3
Use separate cards for cleaning punches	<input checked="" type="checkbox"/>

9. In the **Punch and Plate Setup** panel, enter the number of punches per FTA sample card and the number of cleaning punches.

Check the box to use separate sample-free FTA cards for cleaning punches. If the box is unchecked, cleaning punches are performed on the FTA sample card outside the sample area.

Punch and Plate Setup	
Add controls after samples	<input checked="" type="checkbox"/>
Positive control positions	F11:F11
Negative control positions	G11:G11
Ladder positions	H2:H4:H6:H8:H10:H12
Create duplicate of plates	<input checked="" type="checkbox"/>
Number of punches per card and plate	1
Number of cleaning punches	3
Use separate cards for cleaning punches	<input checked="" type="checkbox"/>

Note: If the option **Use separate cards for cleaning punches** is chosen, the sample-free FTA cards used for the cleaning punches must be placed in a separate FTA card magazine on the deck.


10. In the **Sample information** panel, check the box to use a worklist.

Sample information	
Use worklist	<input checked="" type="checkbox"/>
Worklist mode	<input checked="" type="radio"/> Worklist with card barcodes <input type="radio"/> Worklist with magazine barcodes
Worklist file path	C:\Input\Worklists\PcrSetupftaCards_20CardsWithCardBarcode.csv 
Number of cards	20
Number of magazines	1
Sample type	Saliva

11. If the **Use worklist** box is checked, click on the blue folder icon to browse to the **Worklist file path**.

Sample information	
Use worklist	<input checked="" type="checkbox"/>
Worklist mode	<input checked="" type="radio"/> Worklist with card barcodes <input type="radio"/> Worklist with magazine barcodes
Worklist file path	C:\Input\Worklists\PcrSetupftaCards_20CardsWithCardBarcode.csv 
Number of cards	20
Number of magazines	1
Sample type	Saliva

12. When the browser opens, navigate to the **Worklists** folder, click on a worklist and then click **Open** to select it.

Worklist mode has two categories of worklist:

- **Worklist with card barcodes**
- **Worklist with magazine barcodes**

Worklist with card barcodes

The worklist with card barcodes must contain FTA sample card barcodes. The FTA card magazine barcodes or the positions in the FTA card magazines are not taken into account. The FTA sample cards must be loaded into the FTA card magazines beginning at the FTA card magazine 1, position 1 and without gaps. The order of the FTA sample cards must not be as given in the worklist. The system will determine if an FTA sample card barcode is contained in the worklist regardless of the position in the worklist.

Example of a **Worklist with card barcodes**:

SampleID	SampleCardBarcode	MagazineBarcode	PositionInTheMagazine	Comments
Sample_1	Sample01			Comment_1
Sample_2	Sample02			Comment_2
Sample_3	Sample03			Comment_3
Sample_4	Sample04			Comment_4
Sample_5	Sample05			Comment_5
Sample_6	Sample06			Comment_6
Sample_7	Sample07			Comment_7
Sample_8	Sample08			Comment_8
Sample_9	Sample09			Comment_9
Sample_10	Sample10			Comment_10
Sample_11	Sample11			Comment_11
Sample_12	Sample12			Comment_12
Sample_13	Sample13			Comment_13
Sample_14	Sample14			Comment_14
Sample_15	Sample15			Comment_15
Sample_16	Sample16			Comment_16
Sample_17	Sample17			Comment_17
Sample_18	Sample18			Comment_18
Sample_19	Sample19			Comment_19
Sample_20	Sample20			Comment_20


Worklist with magazine barcodes

The worklist with magazine barcodes must contain FTA card magazine barcodes and the position of the FTA sample cards within the FTA card magazine. The FTA sample card barcodes are not taken into account. The FTA sample cards must be loaded in the FTA card magazine and the position given in the worklist. The **Worklist with magazine barcodes** category allows for “cherry picking” dedicated FTA sample cards in a set of up to nine FTA card magazines.

Example of a **Worklist with magazine barcodes**:

SampleID	SampleCardBarcode	MagazineBarcode	PositionInTheMagazine	Comments
Sample_1		MBC01	1	Comment_1
Sample_2		MBC01	2	Comment_2
Sample_3		MBC01	3	Comment_3
Sample_4		MBC01	4	Comment_4
Sample_5		MBC01	5	Comment_5
Sample_6		MBC01	6	Comment_6
Sample_7		MBC01	7	Comment_7
Sample_8		MBC01	8	Comment_8
Sample_9		MBC02	9	Comment_9
Sample_10		MBC02	10	Comment_10
Sample_11		MBC02	11	Comment_11
Sample_12		MBC02	12	Comment_12
Sample_13		MBC02	13	Comment_13
Sample_14		MBC02	14	Comment_14
Sample_15		MBC03	5	Comment_15
Sample_16		MBC04	1	Comment_16
Sample_17		MBC04	7	Comment_17
Sample_18		MBC05	10	Comment_18
Sample_19		MBC05	14	Comment_19
Sample_20		MBC06	20	Comment_20

After a worklist is loaded, the sample list is shown in the right-hand **Sample list** panel. All samples are automatically selected. Uncheck the **Selected** box if a particular sample should not be processed.

13. If **Use worklist** is not selected in the **Sample information** panel and no worklist is used, manually enter the number of FTA sample cards to be used and the number of magazines.

Sample information

Use worklist ☐

Worklist mode

- ☒ Worklist with card barcodes
- ☐ Worklist with magazine barcodes

Worklist file path

Number of cards

Number of magazines

Sample type

Note: The FTA sample cards must be loaded into the FTA card magazines beginning at FTA card magazine 1, position 1 and without gaps.

14. In the **Sample information** panel, select the sample type from the drop-down list.

Sample information

Use worklist ☐

Worklist mode

- ☒ Worklist with card barcodes
- ☐ Worklist with magazine barcodes

Worklist file path

Number of cards

Number of magazines

Sample type

Kit information

Blood

Saliva

15. In the **Kit information** panel, select the kit to be used from the drop-down list.

Kit information	
Kit type	Investigator 24plex GO! Kit
Kit barcode	Investigator 24plex GO! Kit Investigator ESSplex SE GO! Kit Investigator IDplex GO! Kit
Kit part number	
Kit lot number	
Kit expiry date (yyymmdd)	
Volume positive control	2
Volume master mix	22

16. In the **Kit information** panel, scan the kit barcode into the **Kit barcode** field. The grayed-out fields will automatically populate with the relevant kit details.

Kit information	
Kit type	Investigator 24plex GO! Kit
Kit barcode	03824261812150011512345
Kit part number	0382426
Kit lot number	0011512345
Kit expiry date (yyymmdd)	181215
Volume positive control	2
Volume master mix	22

Note: The instrument will identify an out-of-date kit by highlighting the **Kit expiry date** field in red and warning the user. It is possible to continue with an out-of-date kit.

Example of a kit barcode:

Investigator[®]
24plex GO! Kit (1000)

Cat. No. 382428

Store at -30 to -15°C

Protect from light!

For molecular biology applications in forensic,
human identity and paternity testing

DOM: 2014-12
(YY-MM)

0 3 8 2 4 2 8 1 8 1 2 3 1 1 5 1 2 3 4 5 6 7

Product of Germany

17. In the **Kit information** panel, select the volume for the positive control and the master mix.

Note: Refer to the appropriate kit handbook for information on assay and control volumes.

Kit information	
Kit type	Investigator 24plex GOI Kit
Kit barcode	03824261812150011512345
Kit part number	0382426
Kit lot number	0011512345
Kit expiry date (yyymmdd)	181215
Volume positive control	2
Volume master mix	22

18. Enter a user name in the **User** panel

User	
User name	QIAGEN


19. When setup is completed, start the run by clicking **Start Run**.

The screenshot shows the QIAGEN STAR Q Punch AS software interface. The left sidebar contains four main panels: 'Punch and Plate Setup', 'Sample information', 'Kit information', and 'User'. The 'Punch and Plate Setup' panel includes options for adding controls, selecting control positions, and setting the number of punches. The 'Sample information' panel includes options for selecting a worksheet mode and setting the number of cards and magazines. The 'Kit information' panel includes fields for kit type, barcode, part number, lot number, expiry date, and volumes for positive control and master mix. The 'User' panel includes a field for the user name. The 'Start Run' button is located at the bottom right of the interface.

4.1.2 Loading a run

The **Preparation** window provides, besides general information, a list of the reagents needed as well as the tube or container type in which the reagents need to be loaded. Where a reagent must be prepared manually prior to the start of the run, the composition of the reagent is displayed.

The **Preparation** window serves as a general overview; detailed information regarding the loading of the reagents is given at a later stage in the loading process.

1. Click **Ok** to begin loading.


CAUTION

Loss of data and damage to the instrument

Use only labware defined in this manual with the STAR Q Punch AS instrument. Failure to do so may result in damage to the instrument and incorrect results.

2. Load the punch module in the location indicated in the **Load punch module** window.

Note: The use of different PCR output plates requires the use of different waste troughs.


- Use the 24 mm waste trough in combination with the Applied Biosystems MicroAmp 96-well plate
- Use the 14.5 mm waste trough in combination with the Bio-Rad 96-well PCR plate

Left: 14.5 mm waste trough; right: 24 mm waste trough

3. Press **Ok** to load.

- Load the FTA card magazines containing the FTA cards as indicated in the **Loading of card magazine carrier** window.



Note: If the option to **Use separate cards for cleaning punches** was chosen, the sample-free FTA cards must be placed in FTA card magazine number 5 on the FTA card magazine carrier on track 13–18.

Note: Insert the carriers into the tracks between the front and rear slide blocks of the Autoload tray until they touch the stop hooks on the far side of the tray.

Slide blocks for carriers.

Stop hooks for carriers.

Note: For the barcode to be correctly read, the FTA card magazine barcode must be placed outside the frame printed on the FTA card magazine.

Example of a non-barcoded FTA card magazine:


Example of a barcoded FTA card magazine:


Note: Performance of the cards depends on the quality of the cards. Bent or kinked cards can reduce the robustness of the process.

Example of FTA card magazine containing bent and straight FTA cards:

Note: Correct functioning of the imaging system of the STAR Q Punch AS system is reliant on the quality of the FTA sample card. Make sure that sufficient sample is available in the correct position. Insufficient sample material may make it impossible to punch the sample. A sample area too large may make it impossible for cleaning punches to be made on the same FTA sample card.


Sample area of a FTA sample card:

5. Press **Ok** to continue.

The STAR Q Punch AS instrument will load the FTA card magazine carriers.


6. Load the FTA card recovery magazines and the 50 µl tips in the positions indicated in the **Loading of card recovery and tip carrier** window.

Note: A tip rack must be present in each location. If no tips are required in a position, load an empty tip rack in that position with barcodes facing to the right.


7. Press **Ok** to continue.

The **Edit Tip Count** window allows the user to enter the locations of 50 µl tips being loaded. The system tracks the number of 50 µl tips remaining from the previous run and displays the remaining 50 µl tips in brown.

At the bottom of the window, the **Remove All** button allows the user to remove all tips from their location. Use the zoom and scaling buttons for easier viewing of the locations of the tips.

To add new tips, keep the left mouse button pressed and apply a rubber band over a tip box or a subset of tips. Tips can also be added individually by selecting them with a mouse click.


- Click **OK** when this is done.

The instrument will load the tip carriers.

The system will perform a load check.


If insufficient tips are detected, a warning message will appear and the tip loading procedure will be repeated.

9. Load the PCR output plates, the tubes containing the master mix and the positive control as indicated in the **Loading of plate and tube module carrier** window.


Note: Depending on the number of samples being loaded, one or more tubes containing master mix will be required.

10. Press **Ok** to load. The system will perform a load check.

11. Loading is complete. Press **Ok** to start the run.

After the run the **Maintenance** window displays what to do with each of the components on the STAR Q Punch AS instrument deck.

12. Press **Ok** to finish the run and close the **Run Control** environment.

4.1.3 Output files

The STAR Q Punch AS system generates two types of output files:

- *.csv file: This file is used as a general output file summarizing the run results.
- *.txt file: This file is used as an input file for the subsequent capillary electrophoresis (CE).

Example of a *.csv output file:

	A	B	C	D	E	F
1	<u>EasyPunchCardID</u>	<u>EasyPunchPlateID</u>	<u>SampleID</u>	<u>SampleCardBarcode</u>	<u>PunchType</u>	<u>SourceMagazineBarcode</u>
2	1499	1369	8432_2016-03-14-094942_ID001		SP	
3	1500		8432_2016-03-14-094942_ID002		SP	
4	1501	1370	8432_2016-03-14-094942_ID003		SP	
5	1502		8432_2016-03-14-094942_ID004		SP	
6	1503	1371	8432_2016-03-14-094942_ID005		SP	
7	1504		8432_2016-03-14-094942_ID006		SP	
8			LAD			
9			LAD			
10			LAD			
11			LAD			
12			LAD			
13			LAD			
14			PTC			
15			NTC			

Part A (columns A to F).

	A	G	H	I	J	K
1	<u>EasyPunchCardID</u>	<u>SourceMagazineLabware</u>	<u>SourcePositionInMagazine</u>	<u>DestMagazineBarcode</u>	<u>DestMagazineLabware</u>	<u>DestPosition</u>
2	1499	NucleiCardMagazine_01	1		NucleiCardMagazine_01	1
3	1500	NucleiCardMagazine_01	2	Recovery1	NucleiCardRecoveryMagazine_01	1
4	1501	NucleiCardMagazine_01	3		NucleiCardMagazine_01	3
5	1502	NucleiCardMagazine_01	4	Recovery1	NucleiCardRecoveryMagazine_01	2
6	1503	NucleiCardMagazine_01	5		NucleiCardMagazine_01	5
7	1504	NucleiCardMagazine_01	6	Recovery1	NucleiCardRecoveryMagazine_01	3
8						
9						
10						
11						
12						
13						
14						
15						

Part B (columns G to K).

Example of a .txt output file (CE input file):

1	3500 Plate Layout File Version 1.0					
2						
3	Plate Name	Application Type	Capillary Length (cm)	Polymer	Number of Wells	Owner Name
4	HID 36	POP4	96			
5						
6	Well	Sample Name	Assay	Results Group	File Name Convention	Sample Type User Defined
7	A01	8432_2016-03-14-094942_ID001			Sample	
8	B01	8432_2016-03-14-094942_ID002			Sample	
9	C01	8432_2016-03-14-094942_ID003			Sample	
10	D01	8432_2016-03-14-094942_ID004			Sample	
11	E01	8432_2016-03-14-094942_ID005			Sample	
12	F01	8432_2016-03-14-094942_ID006			Sample	
13	G01	PTC	Positive Control			
14	H01	NTC	Negative Control			

Part A.



1					
2					
3	<u>Barcode</u> Number Comments				
4					
5					
6	ed Field 1	User Defined Field 2	User Defined Field 3	User Defined Field 4	User Def
7					
8					
9					
10					
11					
12					
13					
14					

Part B.

4.1.4 Re-entry procedure

If there are hardware errors during punching, it is possible to re-start the run at the point where the error occurred.

Examples of error messages due to a hardware error:


Follow these steps if a hardware error occurs.

1. Switch OFF the instrument at the power connection on the bottom-left front of the instrument.

2. Close the user messages on the GUI by clicking **Cancel** and/or **Ok** (see examples of error messages above).
3. If it does not occur automatically, close the **HAMILTON Run Control** environment. Use the Windows Task Manager if necessary.


The **HAMILTON Run Control** window:


Note: Make sure the instrument is switched OFF.

4. Open the instrument cover.
5. If necessary, carefully remove any FTA card from the card gripper manually. If appropriate, replace the FTA card in its original position in the FTA card magazine. Make sure the FTA card is not bent.
6. Don't change anything on the STAR Q Punch AS instrument deck. Close the instrument cover!
Note: The instrument cover must be fully closed before the instrument is switched on again.
7. Switch ON the instrument at the power connection on the bottom-left front of the instrument.


8. Return to the GUI. In the **Punch Error** window, choose whether the run should be continued or aborted.

9. In the **Punch Re-Entry** window, select whether the run should be continued with the FTA card being processed when the error occurred or with the FTA card following this card.

The carrier for the PCR output plates moves out. The PCR output plate must be returned manually to its original loading position as indicated in the **Load plates** window.

10. Click **Ok** to continue the run.

The carrier for the PCR output plates will be loaded again and the run is continued.

5 Troubleshooting

5.1 Error handling

If an error occurs, the process can be continued using the error handling procedure. A detailed description is available in the STAR Q Punch AS Software **Help** function. Click on **Error Settings** within the single step dialogs of the STAR Q Punch AS line-specific commands and select **Help**.

5.1.1 Barcode reading error

If, for example, a barcode of a carrier cannot be read, the **Load Carrier – Error** dialog opens:

Position	Error	Assigned recovery
2	Barcode Error	
4	Barcode Error	
5	Barcode Error	
6	Barcode Error	
7	Barcode Error	
8	Barcode Error	
9	Barcode Error	
10	Barcode Error	

In the **Error** panel, all positions producing an error are listed. The following table describes the columns in the table:

Column	Description
(First column)	A red dot identifies a position with an error
Position	The position number of the barcoded item
Error	Short description of the error
Assigned recovery	Selected action

Different errors can lead to the same short error description. A detailed error description is shown for the selected position in the **Description:** field of the **Error** panel.

The following table describes the buttons in the **Recovery** panel and the error recovery program selected:

Button	Error recovery program
Continue	Ignore the error message (here, failure of the barcode reading) Selecting Continue is not indicated in the case of a barcode reading error; a manual entry has to be made so that barcode data exist for further processing.
Repeat	Try again to read the barcode Often a repetition of the reading will solve the problem because the reading speed is reduced.
Barcode...	Enter the barcode manually Click Barcode... to open a dialog box where a barcode can be entered (no entry is also allowed).

1. Assign a recovery option to the selected positions.

The selected action is displayed in the **Assigned recovery** column and the **Execute** button becomes active.

2. Click **Execute**.


The instrument proceeds with the selected recovery option.

Click **Abort** to exit the error recovery procedure without further recovery options.

5.1.2 Pipetting error

If an error occurs with the pipetting channels, a dialog opens to show the error state and recovery options for every channel. Different channels can have different errors.

In case of an LLD error, such as no liquid in the container while aspirating, a dialog similar to the following example opens:

In the **Error** panel, all pipetting channels which produced an error are listed in a table. The following table describes the columns in the table:

Column	Description
(First column)	A red dot identifies a channel with an error
Channel	The number of the channel in the pipetting head
Error	Short description of the error
Assigned recovery	Selected action

Different errors can lead to the same short error description. A detailed error description is shown for the selected channel in the **Description:** field of the **Error** panel.

The following table describes the buttons in the **Recovery** panel and the error recovery program selected:

Button	Error recovery program
Repeat	Executes again the command that caused the error
Available	Aspirates the available volume from the source and fills up the missing volume with air
Continue	Continues as if no error was recognized
Waste	The tip with the error is ejected to the waste and the channel is excluded
Bottom	Activates the channel to move to the bottom of the container; the available volume is aspirated without LLD
Air	Air is pipetted instead of liquid and the method will continue
Exclude	Disables any further action on the selected channel
Move Up	Not an error recovery procedure, but useful e.g., to manually remove a clot. This causes the following actions (this action can be repeated): <ul style="list-style-type: none"> • Moves the barcode reader of the Autoload (if present) to the far right • Moves the selected channel up by 10 mm
Open cover/ Close cover	Not an error recovery procedure; triggers the following actions: <ul style="list-style-type: none"> • Open Cover: enables opening the front cover during error recovery • Close Cover: enables closing the front cover before executing error recovery

1. Assign a recovery option to the selected channels.

The selected action is displayed in the **Assigned recovery** column and the **Execute** button becomes active.

2. Click **Execute**.

The instrument proceeds with the selected recovery option.

Click **Abort** to exit the error recovery procedure without further recovery options.


An error recovery option must be assigned for every error. Selecting a channel followed by any possible recovery procedure assigns the selected recovery procedure to all error-affected channels with the same error.

If any recovery procedure is assigned to a channel (even one which is not desired to be assigned), the red dot in the first column (see picture above) changes color to green.

Some recovery buttons are disabled to prevent further faulty steps. For example, **Continue** cannot be selected for an error-affected aspiration step to prevent any later dispense action with insufficient volume. When a recovery procedure is assigned to the last faulty channel, the **Execute** button becomes active and the system can proceed.

5.1.3 Error handling for easyPunchimaging software

If either the dongle is missing from the computer or the license file is missing, the easyPunchimaging software will not start up and the following message is displayed:

Make sure the USB dongle is connected and the corresponding license file is installed in the correct location.

If the dongle is connected and the license file is in the correct location, contact QIAGEN Technical Services.

Note: For all other instrument error messages, please contact QIAGEN Technical Services

5.2 Log files

In the event of an error, please contact QIAGEN Technical Services.

QIAGEN Technical Services may need your assistance with troubleshooting. The STAR Q Punch AS instrument logs all hardware and software events. These log or trace files might be required for appropriate action to be taken. In addition, screenshots of the error message and/or photos of the problem on the instrument will assist QIAGEN to support you.

Trace files are located in the following system path:

C:\Program Files(x86)\HAMILTON\LogFiles

Two trace files are required for troubleshooting:

- System trace file

[NameMethod]_[24 characters]_Trace.trc

This trace file is generated every time a method (run) is started. Send the system trace file of the erroneous run to QIAGEN Technical Services to assist with troubleshooting.


- Communication trace file

HxUsbCommYYYYMMDD.trc

(If a simulation run was performed, the file name is ComTrace_SimulatorYYYYMMDD.trc)

A communication trace file is generated every day. There is only one communication file per day. This file is also required for troubleshooting.

An example showing the communication trace file and a system trace file:

6 Maintenance

**WARNING/
CAUTION**

Risk of personal injury and material damage

Only perform maintenance that is specifically described in this user manual.

**WARNING/
CAUTION**

Risk of electric shock

Do not open any panels on the STAR Q Punch AS instrument.

Only perform maintenance that is specifically described in this user manual.

The STAR Q Punch AS instrument must be maintained on a regular basis. Regular maintenance consists primarily of surface cleaning and requires opening up the front cover of the instrument. Unscrewing of parts, removal of the deck, side covers or other fixed parts voids the warranty.

The operator is responsible for changing consumable parts including disposable tips and waste bags. QIAGEN Technical Services or service technicians of an authorized agent change spare parts (PCBs, cables, channels, etc.).

Periodic maintenance routines are run to ensure safe and reliable operation of the STAR Q Punch AS instrument and its accessories. A QIAGEN Field Service Specialist or service technician of an authorized agent performs servicing and preventive maintenance on the instrument at least twice a year.

In addition, QIAGEN Technical Services may be called to repair a damaged component of the instrument or to resolve a functional problem which the user cannot resolve such as adjusting and calibrating the channels.

If an error is encountered during a maintenance procedure, attempt to rectify the problem and restart the maintenance procedure. If this fails, contact QIAGEN Technical Services for assistance.

6.1 Maintenance intervals

We recommend the following maintenance intervals:

- Daily before the instrument is shut down
- Weekly at the end of the week before the instrument is shut down
- Six-monthly preventive service maintenance carried out by a QIAGEN Field Service Specialist or service technician of an authorized agent

The STAR Q Punch AS instrument is preconfigured to generate warning messages to perform daily and weekly maintenance procedures.

Note: If the operator does not run either daily or weekly maintenance before shutting down the instrument, these routines must be implemented when the next run is started.

Note: If any parts of the instrument, carriers or racks have become contaminated, the weekly maintenance procedure must be performed (see “Weekly maintenance,” page 77).

6.1.1 Preventative maintenance

Preventive maintenance including verification should be carried out at regular intervals by QIAGEN Technical Services or service technicians of an authorized agent. A service agreement ensures regular maintenance and verification for a specified period of time. We recommend that maintenance and verification take place twice a year (see “Appendix E — Verification,” page 105).

6.2 Materials required for maintenance procedures

CAUTION

Damage to the instrument

Do not use solvents, or reagents containing acids, alkalis or abrasives to clean the STAR Q Punch AS instrument. Do not use disinfecting materials which contain hypochlorite or other bleaching fluids. Use non-corrosive, neutral liquids.

When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles. Use cleaning, disinfecting and decontaminating fluids in accordance with manufacturer's instructions.

The following materials are required for maintenance procedures:

- Paper towels
- Lint-free cloths or Q-tips
- Set of 8 maintenance needles
- Ethanol (70%)
- Deionized water
- Compressed air

6.3 Maintenance procedures

The STAR Q Punch AS Software guides the operator through the regularly scheduled maintenance procedures. The instrument will display a message at start up for optional, warning or mandatory maintenance.


**WARNING/
CAUTION**

Risk of personal injury and material damage


Do not clean the instrument in the vicinity of unshielded flames or devices which could create sparks. Do not use hot air blowers to dry the instrument. The liquids used for cleaning may be flammable.

Example of a maintenance warning message:

1. To initiate maintenance procedures, double-click the **Maintenance** icon on the desktop.

The **Maintenance and Verification Run** window opens. Information is listed on the process and status of all maintenance and verification procedures for the connected/selected instrument.

2. Select the desired maintenance routine by checking the box in the **Execute** column.
3. Press the **Run Process**  button.

The STAR Q Punch AS Software displays instructions detailing all procedures required to perform the selected maintenance routine.

A maintenance routine is completed when the procedure has been fully implemented and the results are within the specifications.

Aborting a maintenance procedure will lead to a failed status, and maintenance will need to be started again.

6.3.1 Daily maintenance

CAUTION


Disposal of plasticware

Used plasticware may contain hazardous chemicals, or contagious/biohazardous materials. Such wastes must be collected and disposed of properly according to local safety regulations.

The following tasks are performed during daily maintenance:

- Check if the deck is clean
- Empty the tip waste
- Check the tightness of the four 1000 µl channels
- Verify the functioning of the cLLD (1000 µl pipetting channel)

After initialization of the instrument, the operator is reminded to execute daily maintenance:


1. Click **Yes** to start the daily maintenance procedure. Click **No** to cancel the procedure.

The hinged acrylic glass window that shields the front of the instrument can be opened for user intervention. The pipetting arm moves to the left side giving the operator access to the deck to check if cleaning is needed.

2. If the deck is clean continue with the daily maintenance.

If the deck needs to be cleaned, interrupt the daily maintenance to carry out weekly maintenance (see "Weekly maintenance," page 77).

3. Empty the tip waste. Dispose of it with the laboratory's contaminated waste.
4. Continue the procedure with the tightness check of the pipetting channels.


The pipetting arm travels to the right-hand side to pick up the teaching needles. Two checks are done with the pipetting channels, the over-pressure and the under-pressure check.

5. Check that the needles are picked-up again for cLLD.

All channels are checked for the proper functioning of the cLLD.

The following message is displayed at completion of daily maintenance:

The daily maintenance process status is saved on the instrument and a report file is created. Refer to "Printing a maintenance report," page 81.

6.3.2 Punch module maintenance

In addition to the daily general instrument maintenance described above, the STAR Q Punch AS also requires maintenance of specific components such as the punch module. This maintenance is automatically started within the daily maintenance procedure and the software guides the user through the process.

CAUTION

Damage to the instrument

Autoclaving cannot be used for instrument components or accessories (sample carriers, magazines, plate slide or punch heads).

CAUTION

Damage to the instrument

Do not use solvents, or reagents containing acids, alkalis or abrasives to clean the STAR Q Punch AS instrument. Do not use disinfecting materials which contain hypochlorite or other bleaching fluids. Use non-corrosive, neutral liquids.

First, the arm moves to the right so access to the punch module is free. The user is then asked to execute the following steps:

1. Empty the contents of punch waste troughs into an appropriate laboratory waste receptacle.
2. Rinse the troughs with water. Make sure that the troughs are dry before replacing them.

3. Remove the plate slide and clean the plate slide and the glass plate with ethanol.


Make sure to wipe the glass plate dry after cleaning.

Note: Leftover liquid on the XY-table leads to dust accumulation impairing the movement of the plate slide.

Note: Do not heat the plate slide over 80°C as this will damage the magnets.

After confirmation, the system moves the punch head into position for it to be removed and cleaned.

4. Press the blue button located to the left of the punch head to release it.

5. Use compressed air to remove all dust from the punch head.

Note: Dust released during cleaning is a potential hazard. Take all necessary safety precautions.

Note: Do not use liquids to sterilize the punch head. This might result in clogging of the dye or the vacuum system thus impairing proper functionality.

6. Before reassembly of the cleaned punch head or fitting of a different punch head, enter the serial number and part number of the punch head when prompted.

This step makes sure that the correct punch head is used for your application. The punch cycles of the head are saved automatically.

The ionizer and the motors for the punch module and the card and plate gripper are switched on.

Note: Replace mechanically deformed carriers or the plate slide to prevent damage to the STAR Q Punch AS instrument.

6.3.3 Card gripper and dust collector maintenance

In addition to the weekly general instrument maintenance described below, the STAR Q Punch AS also requires maintenance of the card gripper and dust collector. This maintenance is automatically started within the weekly maintenance procedure and the software guides the user

through the process. The instrument performs the regular daily maintenance first before it requests the user to perform the following maintenance procedure for the card gripper and dust collector.

CAUTION

Damage to the instrument

Autoclaving cannot be used for instrument components or accessories (sample carriers, magazines, plate slide or punch heads).

CAUTION

Damage to the instrument

Do not use solvents, or reagents containing acids, alkalis or abrasives to clean the STAR Q Punch AS instrument. Do not use disinfecting materials which contain hypochlorite or other bleaching fluids. Use non-corrosive, neutral liquids.

First, the arm moves to the right so access to the punch module is free. The user is then asked to remove all carriers.

1. Remove the dust collector from the vacuum system and empty it.

Note: Dust released during cleaning is a potential hazard. Take all necessary safety precautions.

2. Clean the card gripper fingers with a lint-free cloth soaked in disinfectant.

The instrument tests the x-y-z position of the gripper and compares it to the saved value.

6.3.4 Weekly maintenance

CAUTION

Disposal of plasticware

Used plasticware may contain hazardous chemicals, or contagious/biohazardous materials. Such wastes must be collected and disposed of properly according to local safety regulations.

WARNING

Electrical hazard


When the instrument is connected to line power, terminals may be live. Opening covers or removing parts is likely to expose live parts. Avoid spilling liquid onto or into the instrument. In case of spilling liquid over the instrument, immediately disconnect the instrument from the mains power.

Note: If any parts of the instrument, carriers or racks have become contaminated, the weekly maintenance procedure must be performed.

The following tasks are performed during weekly maintenance:


- Clean the deck and carriers
- Check the condition of the carriers
- Empty and clean the tip waste
- Check the tightness of the four 1000 µl channels
- Verify the function of the cLLD of the four 1000 µl channels
- Clean the four 1000 µl channels: stop disk, O-ring and tip eject sleeve
- Clean the covers and Autoload protecting ribbon

After initialization of the instrument, the operator is reminded to execute weekly maintenance:

1. Unload the deck. The Autoload function carries out this step automatically.
2. Remove all carriers and clean them with disinfectant spray. Leave them to dry.
If carriers are heavily soiled, soak them in a solution of cleaning liquid (refer to the product data sheet for further information).
3. Examine each carrier for scratches on the barcode and any signs of damage. If damage is apparent, replace with new carriers.
The Autoload moves to the right-hand side of the instrument.
4. Open the front cover and wipe the deck with a cloth saturated in disinfectant. The slide blocks in particular must be checked for cleanliness.
5. Close the front cover.
The Autoload moves to the left-hand side of the instrument.
Note: The tip waste, the tip eject plate and the plastic bag are always to be regarded as contaminated.
6. Remove the tip eject plate of the tip waste station, spray disinfectant directly onto the surface and wipe.
7. Remove the frame that holds the plastic bag in place, and discard the plastic in the laboratory's contaminated waste. Pull a new plastic bag over the frame and re-attach it.
8. Put the clean tip eject plate back in place.
9. Remove the maintenance needle and spray disinfectant directly onto the surface and wipe. Replace the clean and dry maintenance needle.
10. Make sure that the reading of barcodes is reliable by checking the laser scanner window of the barcode reader and cleaning it with a lint-free cloth or Q-tips lightly soaked in ethanol (70%).
Note: The laser scanner window must be completely dry and free from dust and fibers before the instrument can be reused.

11. Continue the procedure with the tightness check of the pipetting channels.

The pipetting arm travels to the right-hand side to pick up the teaching needles. Two checks are done with the pipetting channels, the over-pressure and the under-pressure check.

12. Check that the needles are picked-up again for cLLD.

All channels are checked for the proper functioning of the cLLD.

13. Clean the tip eject sleeve on the outer part of the pipetting channels with a lint-free cloth soaked in disinfectant.

Note: Take care that no liquid gets inside the tip channel.

Whenever it is necessary to move channels on the x-arm, move them gently by pushing close to their Y-slide. Never force them as this may lead to damage. If possible, switch on the instrument as this will result in a smoother motion when channels have to be moved on the x-arm.

14. Clean the stop disk and the O-rings of the pipetting head on the outer part of the pipetting channels with a lint-free cloth soaked in disinfectant.


Note: Take care that no liquid gets inside the tip channel.

Whenever it is necessary to move channels on the x-arm, move them gently by pushing close to their Y-slide. Never force them as this may lead to damage. If possible, switch on the instrument as this will result in a smoother motion when channels have to be moved on the x-arm.

15. Spray the front and side cover with disinfectant and wipe dry.
16. Clean the Autoload protecting ribbon with a cloth soaked in disinfectant and wipe without exerting pressure.
17. Clean the x-guide shaft behind the upper front cover with a dry cloth at least once a month.
18. Make sure that carriers are completely clean and dry before re-using.

The following message is displayed at completion of weekly maintenance:

The weekly maintenance process status is saved on the instrument and a report file is created. Refer to "Printing a maintenance report," page 81.

6.4 Printing a maintenance report

The maintenance process status can be printed. To print such a report:

1. From the **File** menu, select **Open Report**.
All maintenance and verification processes found in the default **Report Path** are listed.
2. If necessary, change the report path using the browse button <...>.
3. Select a report and press the **Open** button. The **Report Viewer** displays the selected report file.
4. From the **File** menu, select **Print** to print the report file.

6.5 Instrument decontamination

CAUTION

Damage to the instrument

Do not use solvents, or reagents containing acids, alkalis or abrasives to clean the STAR Q Punch AS instrument. Do not use disinfecting materials which contain hypochlorite or other bleaching fluids. Use non-corrosive, neutral liquids.

CAUTION

Damage to the instrument

Autoclaving cannot be used for instrument components or accessories (sample carriers, magazines, plate slide or punch heads).

In general, good laboratory practice must be observed for decontamination.

1. Spray the front and side cover of the instrument with disinfectant.
2. Open the front cover and spray the deck with disinfectant.
3. Remove the tip eject plate of the tip waste station and clean it.
4. Spray disinfectant directly onto the surface of the tip waste station.
5. Remove the frame that holds the plastic bag in place and discard the plastic bag in the laboratory's contaminated waste.
6. Replace the tip eject plate.
7. Clean the tip eject sleeve on the outer part of the pipetting channels with a lint-free cloth soaked in disinfectant.

8. Clean all carriers with disinfectant liquid and leave them to dry.

If the carriers are heavily soiled, soak them in a disinfectant solution (see the product data sheet for further information).

Note: Autoclaving, using superheated steam under pressure, may only be used for decontamination followed by final disposal of consumables and racks; they will be destroyed by the process.

6.5.1 Surface decontamination using liquids

Use the disinfecting procedure described above for decontamination. Other procedures have not been tested by QIAGEN.

In some cases other decontamination procedures may be desirable, e.g., for reliably destroying infectious materials or DNA/RNA. Many of these decontamination procedures are very aggressive and can cause damage to the STAR Q Punch AS instrument.

If using other decontamination procedures not listed here, be aware that they may increase service and maintenance requirements and may make shorter maintenance intervals necessary.

If using other liquids or sprays for surface decontamination, follow the manufacturer's instructions. Pay particular attention to potential corrosiveness, e.g., acidic or alkaline solutions and oxidizing agents. Use of such agents may increase service and maintenance requirements (O-rings exchange, greasing of spindles, etc.) and may make shorter maintenance intervals necessary. .

When the instrument deck and carriers are cleaned using enzyme solutions such as DNase and RNase, make sure to thoroughly remove any remainders by wiping deck and carriers with deionized water to avoid leaving aggressive substances on the surface.

6.5.2 Decontamination using gases

CAUTION

Damage to the instrument

Ethylene oxide fumigation may increase service and maintenance requirements (O-rings exchange, greasing of spindles, etc.) and may make shorter maintenance intervals necessary.

CAUTION

Damage to the instrument

Hydrogen peroxide fumigation leads to bleaching or discoloration of many instrument materials and may increase service and maintenance requirements (O-rings exchange, greasing of spindles, etc.) and may make shorter maintenance intervals necessary.

CAUTION

Damage to the instrument

Do not use formaldehyde fumigation or chlorine oxides (chemical compounds of chlorine and oxygen such as bleach). They are not suitable for the STAR Q Punch AS instrument because of chemical reaction and corrosion.

Should decontamination by fumigation be necessary, QIAGEN recommends using ethylene oxide. QIAGEN does not carry out such fumigation procedures; use a 3rd party contractor for such service.

Fumigation using hydrogen peroxide (H_2O_2) is possible but not recommended. QIAGEN does not carry out such fumigation procedures; use a 3rd party contractor for such service.

6.5.3 Decontamination using UV light

CAUTION

Damage to the instrument

UV irradiation causes many synthetic materials to become brittle. This may increase service and maintenance requirements and may make shorter maintenance intervals necessary.

A UV light is mounted on the STAR Q Punch AS instrument. We recommend UV irradiation for 15 minutes per day. The arm will move left and right while the UV lamp is on. Make sure that the path is not blocked. If using other UV irradiation options, use wavelength, intensity and duration according to manufacturer's instructions.

The instrument cover is made of Makrolon®, not of standard acrylic glass, due to its better resistance to UV irradiation.

7 Glossary

Term	Description
Adjustments	Detailed positional setting for the hardware.
Air displacement pipetting	Air displacement means that the liquid is aspirated into and dispensed from a disposable tip by the movement of a plunger.
Aliquots	Aliquots are identical small volumes of liquid.
Autoload	Hardware assembly that enables automatic loading of the STAR Q Punch AS instrument. It consists of a loading head movable in the Y-direction, which draws the carriers into the instrument and reads the barcodes on the carriers.
Barcode reader	Device for reading sample/plate barcodes; part of the Autoload function.
Carrier	Unit for loading plates and tips on the STAR Q Punch AS instrument deck. A carrier is handled by the Autoload function.
cLLD	Capacitive liquid level detection.
Container	A container defines a vessel or a single well of a plate.
Deck	The work surface of the STAR Q Punch AS instrument presents the largest possible area. The placing of the carriers on the deck is defined by the tracks (see T/track, below) providing the tracks are in the operating range of the pipetting area.
Deck layout	A collection of labware placed on a deck.
Dispense	To distribute quantities of liquid from a pipetting device.
Firmware	Lower level program code that is carried out on the processors of the STAR Q Punch AS instrument.
Front cover	Protective covering for the STAR Q Punch AS instrument featuring a hinged front window made of transparent acrylic glass. The work surface of the STAR Q Punch AS instrument is shielded from user intervention and other outside influences such as dust. It also protects the user from the movements of the STAR Q Punch AS instrument.
FTA	Special surface treatment technology for forensics cards.
Good laboratory practices	Also written as GLP, are set of appropriate laboratory behaviors which should be observed.

Term	Description
GUI	Graphical user interface: allows users to interact with electronic devices through graphical icons and visual indicators as opposed to text-based interfaces.
Labware	Movable items placed on the deck, such as carriers, containers or racks.
Liquid	Includes reagents, controls, standards and wash fluids.
LLD (liquid level detection)	Positive detection of liquid achieved either by pressure or capacitive signal detection and transfer.
Loading/unloading	The process by which plates, tubes and tip carriers are brought on and off the deck either manually or automatically by the Autoload function.
MAD	Monitored air displacement: aspiration monitoring feature. During the aspiration process, the pressure within the pipetting channel is measured in real time.
Magazine	A magazine is designed to hold different numbers of cards (20 cardboard framed cards per magazine).
MXF	Multiflex; a module that can be flexibly arranged on specific carriers for the STAR Q Punch AS instrument
Microplate	A plate with 96 wells (8 x 12) 9 mm wide, standard SBS format.
Pause	Interruption of processing. The current processing steps are completed.
Pipetting	Transfer of liquids, usually a defined volume, from one container to another.
Pipetting arm	Assembly equipped with the pipetting tool and/or plate handler, as well as the common X-drive.
Pipetting channel	Hardware including the function of picking up a tip, aspirating, dispensing, tip eject, liquid level detection and the Y/Z-movements.
Pipetting module	Firmware (lower level program code) which controls a pipetting channel. Included in the category are the Y and Z pipetting movements, and LLD.
pLLD	Pressure-based liquid level detection.
Processing step	Defines what must be carried out on the STAR Q Punch AS instrument, as well as the location where it must be carried out and possible interaction with other system components or labware. The action is defined in accordance with the loading and the tasks.

Term	Description
Rack	Grouping of containers, such as magazines, microplates, etc.
Rack identification	Barcode for rack identification.
Run	<p>Execution of the processing steps defined in the relevant method with the aim of processing one or more liquids and containers. The run is a series of timed commands to carry out processing on the instrument according to the processing plan.</p> <p>The run can be interrupted to load more elements. Then processing continues according to a newly calculated processing plan; the run starts again. Loading is not a part of the run.</p>
Sample	Refers to a liquid in an unambiguously identified container which is to be processed.
SBS format	Standard format for microplates defined by the Society for Biomolecular Screening.
STR	Short tandem repeat.
T/track	The STAR Q Punch AS instrument has equal partitions of 22.5 mm, equivalent to 1-T on the deck. Labware carriers are adapted to those partitions.
TCC	Temperature controlled carrier (not applicable in this configuration).
Tip	Disposable tip for pipetting.
Tip rack	Frame that holds the tips.
Tip waste	Container for used and ejected tips.
Trace	Record of the status during processing.
User	User of the software. Access rights for different types of users can be defined, such as operators, laboratory managers, etc.
Verification kit	Balance, liquid and disposable tips to verify the function (volume check) of the pipetting heads.
Well	The individual container of a microplate.
Work area	The area on the STAR Q Punch AS instrument accessed during processing. Elements to be pipetted or handled can be placed in this area.

Appendix A — Ordering Information

Product	Contents	Cat. no.
Instrument		
STAR Q Punch AS (EC)	Instrument with Autoload including 4 x 1000 µl channels, card and plate gripper and punch module. This product configuration is used for EasiCollect and COPAN cards.	9002651
Investigator kits		
Investigator® 24plex GO! (1000)	Primer Mix, Fast Reaction Mix 2.0 including Taq DNA polymerase, Control DNA, allelic ladder 24plex, DNA size standard 550 (BTO)	382428
Investigator IDplex GO! Kit (1000)	Primer mix, Fast Reaction Mix including HotStarTaq® Plus DNA Polymerase, Control DNA, allelic ladder IDplex GO!, DNA size standard 550 (BTO)	381638
Investigator ESSplex SE GO ! Kit (1000)	Primer mix, Fast Reaction Mix including HotStarTaq Plus DNA Polymerase, Control DNA, allelic ladder ESSplex SE GO!, DNA size standard 550 (BTO)	381568
Consumables		
CO-RE Filter-Tips, 50 µl (5760)	50 µl CO-RE Filter Tip, sterile, black, for STAR Q Instruments	990065
Biohazard Waste Bags (25)	25 biohazard waste bags	990123

Appendix B — Technical Specifications

STAR Q Punch AS instrument specifications

Dimensions (w x d x h)	Width: 112.4 cm (44.25 inches) [96/384 Probe Head: 145 cm (57 inches)] Depth: 79.5 cm (31.3 inches) [Autoload: 100.6 cm (39.6 inches)] Height: 90.3 cm (35.6 inches)
Work area dimensions	Width (x): 67.5 cm (26.6 inches) Depth (y): 49.7 cm (19.6 inches) Height (z): 14.5 cm (5.7 inches) The maximum height for labware used on the deck is 14.0 cm (5.5 inches)
Weight	135 kg (297 lbs.) [96/384 Probe Head: 150 kg (331 lbs)]
Deck capacity	30 tracks (T) allow combinations of: Maximum of 30 tube carriers (1 T) holding 24 or 32 tubes per carrier Maximum of 5 carriers (6 T) holding 5 tip racks or 5 plate positions per carrier
Modal precision	x-y-z positional accuracy of 0.1 mm
Tip size	Low and intermediate volume: 50 µl
Power consumption	Standby power consumption: 100 VA

Operating data

Maximum power consumption	600 VA or 1000 VA (depending on configuration)
Voltage	100 VAC/115 VAC/230 VAC
Frequency	50/60 Hz \pm 3 Hz
Delayed action fuse 600 VA 1000 VA	115V: 6.3A (T6.3AL250) 230V: 3.15A (T3.15A250) 115V: 10A (T10AL250) 230V: 5A (T5AL250)
Installation category	II
Pollution degree	2
Temperature range	15–35°C
Relative humidity	30 % – 85 % (non-condensing, indoors)
Noise level	<65 dBA (according to EN27779) <46 dBA in standby mode
Altitude	Maximum 2000 m above sea level
Heat: the power consumed will be transferred to head	Example: 600 or 1000 Watts of Heat = 600 or 1000 Joules/second
Recommended computer	Intel® Core™ 2 Duo, \geq 4 GB RAM, 500 GB Hard Drive, 16x DVD +/-RW, DirectX 250 MB graphic card, Microsoft® Windows® 7 Professional. (not included in shipment)
Communication	USB or RS232 with dual processor board Ethernet or USB with LAN dual processor board
Location	Indoor use only

Storage and transportation

Temperature range	–25°C to 70°C
Relative humidity	10–90 % (non-condensing, indoors)

QIAGEN standard computer specifications

Configuration (overview)	<p>OptiPlex® XE DT: N Series Base Desktop</p> <p>Minimum requirements:</p> <p>Intel Core 2 Duo E8400 (3GHz, 1333 MHz, 6MB) 4096 MB (2x 2048 MB) 1333MHz DDR3 Dual Channel RAM 16x DVD+/-RW 2x 250 GB, RAID1, HDD</p> <ul style="list-style-type: none">• 2x RS232, 6x USB, 2x Ethernet, 1x Digital Video Port <p>Flat Screen 21.5 inches (1920 x 1080) Webcam Mouse and Dell® Soundbar Country specific keyboard Microsoft Windows 7 Ultimate 64 Bit Microsoft Office 2013 for Home and Business Recovery DVD 2 year Dell warranty (next business day support)</p>
-----------------------------	--

Labware for STAR Q Punch AS

Manufacturer	Catalog no.	Description
GE Healthcare	WB120206	EasiCollect card with indicator for saliva and blood
Copan	4474001 4473974 (with indicator)	NUCLEIC-CARD collection card for blood and saliva
Bio-Rad	HSP9901	Hard-Shell Low-Profile Thin-Wall 96-well PCR Plate, fully-skirted
ABI	N8016154	Applied Biosystems MicroAmp Optical 96-well Reaction Plate
ABgene®	AB1100	PCR Plate, 96-Well, Semi-Skirted, Raised Deck
Sarstedt®	72.694	2.0 ml skirted tube
Sarstedt	72.730	0.5 ml skirted tube

Pipetting specifications for disposable tips (individual 1000 µl channels)

Disposable tip size	Volume	Trueness R (%)	Precision CV (%)
50 µl	0.5 µl	10.0 %	6.0 %
50 µl	1 µl	5.0 %	4.0 %
50 µl	5 µl	2.5 %	1.5 %
50 µl	50 µl	2.0 %	0.75 %

Operating data for tips

Temperature range	15–25°C
-------------------	---------

Storage data for tips

Temperature	Maximum 55°C
-------------	--------------

Liquid level detection

Individual channels	Capacitive Liquid Level Detection (cLLD) and pressure (pLLD) on aspiration, cLLD on dispense, minimum volume 10 µl, depending on container type
---------------------	---

CO-RE gripper 1000 µl specifications

Labware format

Microtiter footprint	plate height ≤43 mm
----------------------	---------------------

Absolute positioning

Accuracy	x, y, z = 0.5 mm
Reproducibility	x, y, z = 0.25

Movement range

x	Track 1 – n (depending on instrument type)
y	Depending on number of channels and used front channel
z	Lowest position = 15 mm over metal deck sheet

Gripper opening

Arm	Modular pipetting arm
Minimum opening	9 mm
Maximum opening	Dependent upon travel range on arm

Gripping force	5 N – 16 N (default 9 N)
----------------	--------------------------

Transport mass	300 g filled deep well plate
----------------	------------------------------

Autoload: barcode and reader specifications

Carriers, containers, racks and tip racks can be identified by a barcode. A reader mounted on the Autoload slide scans the barcode. The system must allow specification of ranges (barcode mask) for plausibility checking of barcode information.

Barcode symbologies

The following barcode symbologies can be detected by the system:

- ISBT standard
- Code 128 (subset B and C)
- Code 39
- Codabar
- Code 2 of 5 Interleaved
- UPC A

QIAGEN recommends using barcode type Code128 (subset B and C)

Reading accuracy

The rate of inaccurate readings of sample plates and container barcodes is less than 1 ppm.

The above-mentioned specification is valid under the following conditions:

- Barcode symbology module: ISBT standard
- Code density: 0.0065 inch (0.1651 mm)
- Print Quality: see "Barcode specifications," page 96
- Recognized errors are defined as an accurate reading

Barcode specifications

Type	Modular pipetting arm
Length of string	9 mm
Maximum opening	Dependent upon travel range on arm
Code density	Minimum module width (x dimension) including a print tolerance: ≥ 0.0065 inch (0.1651 mm)
Tolerance	Maximum module width (x dimension) including a print tolerance: ≤ 0.02 inch (0.508 mm) Best reading performance with x dimension: ≥ 0.01 inch (0.254 mm)
Check character:	
ISBT standard	One character
Code 128	One character
Code 39	None
Codabar	None
Code 2 of 5 interleaved	None
UPC A	One character
Quiet zone	≥ 10 times the x dimension, but at least 3 mm
Print quality	The barcode print must be of a high quality. A printed barcode with an ANSI/CEN/ISO grade A or B is required. Offset, typographic, intaglio and flexographic printing are suitable. Mechanical dot matrix and thermo matrix printing are not suitable. The surface may be treated, sealed or plastic-coated.

Tube barcodes

Specifications

	Dimension	Minimum	Maximum
A	Label length	–	80 mm
B	Code length	–	74 mm
C	Quiet zone	3 mm	–
D	Label width	12 mm	–
E	Code width	12 mm	–
F	Distance from code to label edge	–	1 mm

Positioning tube labels

The label must be glued within a range of between 20 mm to 100 mm from the bottom of the tube. The label must fit tightly at an angle of approximately 90° to the tube. The label must fit tightly over its whole length.




Plate barcodes

Specifications


	Dimension	Minimum	Maximum
A	Label length	–	66 mm
B	Code length	–	30 mm
C	Quiet zone	3 mm	–
D	Label width	10 mm	–
E	Code width	7 mm	–
	Distance from code to label edge (if necessary)	–	1 mm

Positioning plate labels

The plate barcode must fit on side A of the plate.

The barcode label must be centered and parallel to the edge of the plate.

The barcode label must not protrude above or below the edge of the plate.

Specifications for STAR Q Punch AS specific sample barcodes

1D barcodes

1D barcode types	ISBT Standard Code 128 (Subset A,B and C) Code 39 Codabar Code 2 of 5 interleaved UPC A/E JAN/EAN 8/13
Check character: ISBT Standard (128) Code 128 Code 39 Codabar Code 2 of 5 interleaved UPC A/E JAN/EAN 8	One character One character None, optional check character to be checked within programming None, optional check character to be checked within programming None, optional check character to be checked within programming One character
Length of string	Maximum 20 characters excluding start stop and check characters, depending on the code length (see label dimensions)
Code density	Minimum module width (X dimension) including a print tolerance: ≥ 0.0065 inch (0.1651 mm)
Tolerance	Maximum module width (X dimension) including a print tolerance: ≤ 0.0 inch (0.508 mm)
Quiet zone	≥ 10 times the x dimension, but at least 3 mm
Print quality	The barcode print must be of a high quality. A printed barcode with an ISO/IEC 15416: 2002 grade 4 (A) or 3 (B) is required. Offset, typographic, intaglio and flexographic printing are suitable. Mechanical dot matrix and thermo matrix printing are not suitable. The surface may be treated, sealed or plastic-coated.

2D barcodes

Length of string	Maximum 2000 characters depending on density and code length (see label dimensions)
Code density	Minimum module width (X dimension) including a print tolerance: ≥ 0.010 inch (0.254 mm)
Tolerance	Maximum module width (X dimension) including a print tolerance: ≤ 0.04 inch (1.0 mm)
Check character: ISBT standard DataMatrix (EC200) PDF417	ISO/IEC 16022: 2006 ISO/IEC 15438: 2006
Error correction: DataMatrix (ECC200) QR-Code	ISO/IEC 16022: 2006 ISO/IEC 24720: 2008 ISO/IEC 18004: 2009
Quiet zone	≥ 4 cells, but at least 3 mm, each side vertically and horizontally
Print quality	The barcode print must be of a high quality. A printed barcode with an ISO/IEC 15416: 2002 grade 4 (A) or 3 (B) is required. Offset, typographic, intaglio and flexographic printing are suitable. Mechanical dot matrix and thermo matrix printing are not suitable. The surface may be treated, sealed or plastic-coated.

Card specifications

Manufacturer	Description	Dimensions
GE Healthcare	EasiCollect card	51 mm x 51 mm
Copan	NUCLEIC-CARD collection card	51 mm x 51 mm

STAR Q Punch AS with Autoload function

Appendix D — Chemical Compatibility

The table for chemical compatibility is based on information from different manufacturers. The results refer to laboratory tests with raw materials. The outcomes with these materials are liable to effects which cannot be observed under laboratory conditions (e.g., temperature, pressure, tension, chemical influences of substances, design features, etc.). The results listed should be considered only as a guideline. In case of doubt we recommend significant tests. The chemical resistance is not sufficient for an evaluation of a particular material for a product. Particular regulations, e.g., explosion prevention in the case of flammable liquids, have to be taken into account.

Key:

1.4034	X46Cr13 steel	FFKM	Kalrez®
1.4301	X2CrNi19-11 steel	FKM	Viton®
1.4305	X8CrNiS18-9 steel	PE	Polyethylene
1.4404	X2CrNiMo17-12-2 steel	PEEK	Polyetheretherketone
1.4435	X2CrNiMo18-14-3 steel	PP	Polypropylene
EPDM	Ethylene-propylene elastomer	PTFE	Polytetrafluorethylene
EPT	Ethylene-propylene terpolymer	ZrO ₂	Zirconium oxide

CO-RE head consists of 1.4305, EPDM, PEEK, ZrO₂ and PTFE

- 1 = No effect, little or no noticeable change
- 2 = Slight corrosion or discoloration
- 3 = Moderate corrosion or other change in physical properties or dimensions;
not recommended for continuous contact
- 4 = Severe corrosion or physical change; prolonged contact not recommended
- 0 = No data
- L = Danger of pitting corrosion (a localized form of corrosion that leads to the creation of small holes in the metal)

Chemical compatibility of materials and level of resistance to corrosion

Chemical	Material													
	1.4034	1.4301	1.4305	1.4404	1.4435	PE	PP	PTFE	PEEK	FKM	FFKM	EPT	ZrO ₂	CO-RE head
Acetic acid, 20%	2	1	1	1	1	1	1	1	1	2	1	1	0	1
Acetic acid, glacial	2	1	1	1	1	1	1	1	1	4	1	1	0	1
Acetone	1	1	1	1	1	2	1	1	1	4	1	1	0	1
Acetonitrile	1	1	1	1	1	1	3	1	0	2	0	3	0	3
Ammonium hydroxide, 5%	1	1	1	1	1	1	1	1	1	1	1	1	0	1
Chloroform	1	1	1	1	1	3	3	1	1	1	1	4	0	4
Deionized water	1	1	1	1	1	1	1	1	1	1	1	1	0	1
Dimethyl formamide	1	1	1	1	1	1	1	1	1	3	1	1	0	1
Dimethyl sulfoxide	1	1	1	1	1	1	1	1	0	0	0	1	0	1
Ethanol	1	1	1	1	1	1	1	1	1	1	1	1	0	1
Ethyl acetate	1	1	1	1	1	2	1	1	1	4	1	1	0	1
Hexane	1	1	1	1	1	3	2	1	1	1	1	4	0	4
Hydrochloric acid, 5%	4L	2L	3L	2L	2L	1	1	1	1	1	1	1	1	1
Hydrochloric acid, 20%	4L	3L	3L	2L	2L	1	1	1	1	1	1	1	1	1
Hydrogen peroxide, 10%	1	1	1	1	1	2	2	1	1	2	2	2	1	2
Isopropyl alcohol	1	1	1	1	1	1	1	1	1	1	1	1	0	1
Methanol	1	1	1	1	1	1	1	1	1	2	1	1	0	1
Methylene chloride	1	1	1	1	1	4	3	1	2	2	1	4	0	4
Nitric acid, 5–10%	1	1	1	1	1	1	1	1	1	1	1	3	1	3
Nitric acid, 70%	1	1	1	1	1	3	4	1	1	2	1	3	1	3
Phosphate buffer	1	1	1	1	1	1	1	1	0	1	1	1	0	1
Phosphoric acid, 85%	3	2	3	2	2	1	1	1	0	1	1	1	1	3
Potassium hydroxide conc.	3	1	2	1	1	1	1	1	1	3	1	1	1	2
Sodium acetate	1	1	1	1	1	1	1	1	0	4	1	1	0	1
Sodium borate	1	1	1	1	1	1	1	1	0	1	1	1	0	1
Sulfuric acid, 1–75%	4	2	3	2	2	1	1	1	2	1	1	1	1	3
Urine	1	1	1	1	1	1	1	1	1	1	1	1	0	1
Triethylamine	1	1	1	1	1	0	4	1	0	4	0	4	0	4
Toluene	1	1	1	1	1	3	3	1	1	1	1	4	0	4
Sodium hydroxide 5%	1	1	1	1	1	1	1	1	1	1	1	1	1	1
Formic acid 5%	3	1	2	1	1	1	1	1	1	2	1	1	0	2
Sodium hypochlorite 10%	3L	2L	2L	1L	1L	1	1	1	0	1	1	1	0	1

Appendix E — Verification

To perform the verification of the STAR Q Punch AS instrument, Field Verification 2 is used. Field Verification 2 can only be performed by trained personnel.

Verification specifications are different from the specifications given in the technical specifications of this user manual.

The volume verification specifications are different from the pipetting specifications for disposable tips given in the technical specifications of this user manual. The field verification contains validated procedures and equipment defined by QIAGEN to demonstrate and to verify the correct function of the instrument according to specifications given by QIAGEN suitable for the field. The field verification is therefore a reference defined by QIAGEN to compare the instrument's performance according to given procedures valid for a broad operating range.

Based on that, specifications as applicable in the technical specifications will be achieved by maintaining defined environmental conditions in the laboratory, by keeping the operating range as small as possible, by optimizing the methods such as adapting the liquid classes, knowing the sample liquids and the characteristics of used labware, etc. See also "Appendix B — Technical Specifications," page 89.

For the four 1000 µl channels, a dye-pipetting procedure followed by gravimetric and photometric analysis is used to verify the trueness and precision.

Devices and accessories such as the barcode reader and cover safety can also be verified with Field Verification 2.

Note: This is an additional feature of the STAR Q Punch AS instrument not included in the STAR Q Punch AS Software.

Appendix F — Waste Electrical and Electronic Equipment (WEEE)

This section provides information about disposal of waste electrical and electronic equipment by users.

The crossed-out wheeled bin symbol (see below) indicates that this product must not be disposed of with other waste; it must be taken to an approved treatment facility or to a designated collection point for recycling, according to local laws and regulations.

The separate collection and recycling of waste electronic equipment at the time of disposal helps to conserve natural resources and ensures that the product is recycled in a manner that protects human health and the environment.

Recycling can be provided by QIAGEN upon request at additional cost. In the European Union, in accordance with the specific WEEE recycling requirements and where a replacement product is being supplied by QIAGEN, free recycling of its WEEE-marked electronic equipment is provided.

To recycle electronic equipment, contact your local QIAGEN sales office for the required return form. Once the form is submitted, you will be contacted by QIAGEN either to request follow-up information for scheduling collection of the electronic waste or to provide you with an individual quote.

Appendix G – Warranty Statement

Thank you for your purchase of QIAGEN instrumentation. Your instrument has been carefully tested to ensure optimum operating efficiency and reproducibility of results. QIAGEN warrants that all new instrumentation manufactured by QIAGEN will correspond to the product specifications and be free from defects in workmanship and materials for a period of twelve (12) months from the original date of shipment. Repair or replacement of defective parts will be provided to the purchaser during this time period provided the QIAGEN instrumentation is operated under conditions of normal and proper use, but not for damage caused by the customer. If any part or subassembly proves to be defective, it will be repaired or replaced at QIAGEN's sole option, subsequent to inspection at the factory, or in the field by an authorized factory representative, provided that such defect manifested under normal and proper use.

Limitation of warranties and remedies

THE FOREGOING WARRANTY IS QIAGEN'S SOLE AND EXCLUSIVE WARRANTY, AND REPAIR OR REPLACEMENT OF DEFECTIVE PARTS IS THE SOLE AND EXCLUSIVE REMEDY. THERE ARE NO OTHER WARRANTIES OR GUARANTEES, EXPRESS OR IMPLIED. THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE EXPRESSLY EXCLUDED, TO THE FULLEST EXTENT PERMITTED BY LAW. (NOTE: SOME STATES DO NOT PERMIT DISCLAIMERS OF IMPLIED WARRANTIES SO THIS LIMITATION MAY NOT APPLY TO YOU). WITH THE EXCEPTION OF THE ABOVE-REFERENCED REPAIR OR REPLACEMENT REMEDY, QIAGEN SHALL HAVE NO OBLIGATION OR LIABILITY OF ANY NATURE WHATSOEVER WITH RESPECT TO THE QIAGEN INSTRUMENTATION, WHETHER ARISING IN CONTRACT, TORT, STRICT LIABILITY, OR OTHERWISE, INCLUDING BUT NOT LIMITED TO, LIABILITY FOR INDIRECT, CONSEQUENTIAL, INCIDENTAL AND/OR SPECIAL, PUNITIVE, MULTIPLE AND/OR EXEMPLARY DAMAGES AND/OR OTHER LOSSES (INCLUDING LOSS OF USE, LOST REVENUES, LOST PROFITS AND DAMAGE TO REPUTATION), EVEN IF SUCH DAMAGES WERE FORESEEN OR FORSEEABLE, OR WERE BROUGHT TO QIAGEN'S ATTENTION. IN NO EVENT SHALL QIAGEN'S LIABILITY TO YOU EXCEED THE PURCHASE PRICE OF THE PRODUCT.

Liability clause

QIAGEN shall be released from all obligations under its warranty in the event repairs or modifications are made by persons other than its own personnel, except in cases where the Company has given its written consent to perform such repairs or modifications. All materials replaced under this warranty will be warranted only for the duration of the original warranty

period, and in no case beyond the original expiration date of original warranty unless authorized in writing by an officer of the Company. Read-out devices, interfacing devices and associated software will be warranted only for the period offered by the original manufacturer of these products. Representations and warranties made by any person, including representatives of QIAGEN, which are inconsistent or in conflict with the conditions in this warranty shall not be binding upon the Company unless produced in writing and approved by an officer of QIAGEN.

Appendix H — Declaration of Conformity

Name and address of the legal manufacturer:

Hamilton Bonaduz AG
Via Crusch 8
7402 Bonaduz
Switzerland

An up-to-date Declaration of Conformity can be requested from QIAGEN Technical Support.

Appendix I — FCC Declaration

The United States Federal Communications Commission (FCC) in 47 CFR 15.105 declares that the users of this equipment must be informed of the following facts and circumstances:

Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

The following text appears on the right inside front corner of the “Output” drawer of the instrument:

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

QIAGEN is not responsible for any radio or television interference caused by unauthorized modifications of this instrument or the substitution or attachment of connection cables and equipment other than those specified by QIAGEN. The correction of interference caused by such unauthorized modification, substitution or attachment will be the responsibility of the user.

Appendix J – Safety Information (French, FR)

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Index

- 2D barcodes, 101
- About this user manual, 6
- Achtung, 121
- Actuator, 24
- Adapter frame, 26, 40
- Attention, 111
- Autoload, 26
 - specifications, 95
- Avertissement, 111
- Barcode
 - 2D barcodes, 101
 - autoload, 26
 - card magazine, 51
 - error, 63
 - kit, 46
 - specifications, 95
 - tip recognition, 30
 - worklist, 43
- Camera, 25
- Capillary electrophoresis, 57
- Card
 - Copan, 8
 - EasiCollect, 8
 - FTA, 8
 - NUCLEIC-CARD, 8
- Card gripper, 22, 75
- Cards
 - analysis, 37
 - magazines, 25
 - punching, 37
 - recovery magazines, 25
 - specifications, 101
 - transport, 37
- Carrier, 27
 - magazine, 21
- Carrier slide blocks, 50
- Carrier stop hooks, 50
- Cautions, 9
- Computer, 33
 - specifications, 91
- Copan, 8
- Data security, 34
- Deck, 19
 - layout, 21
- Decontamination, 81
- Dust collector, 75
- EasiCollect, 8
- Error
 - Help function, 63
 - imaging software, 67
 - pipetting, 65
- Error recovery
 - barcode, 64
 - pipetting, 66
 - re-entry, 59
- File
 - output, 57
- FTA, 8
- Glossary, 85
- Gripper, 28
 - specifications, 93
- Imaging software, 34, 67
- Indicator, 52
- Intended use, 7
- Interfaces, 35
- Ionizer, 24
- Labware, 92
- Light source, 25
- Liquid level detection, 31
- Magazines, 21, 25
- Magnetic field safety, 10, 24
- Maintenance
 - card gripper, 75
 - daily, 72
 - dust collector, 75
 - intervals, 70
 - materials, 70
 - preventative, 70
 - punch module, 74
 - report, 81
 - warning, 71
 - warning messages, 70
 - weekly, 77
- NUCLEIC-CARD, 8
- Ordering information, 88
- Output files, 57
- Pipetting, 29
 - liquid level detection, 31
 - MAD, 31
- Pipetting arm, 29
- Plate gripper, 22
- Plate slide, 24
- Platform, 19
- Policy statement, 7

- Punch head, 24
- Punch module, 23, 74
- Recovery magazines, 25
- Re-entry procedure, 59
- Safety information
 - biological, 13
 - chemicals, 13
 - electrical, 12
 - labels, 18
 - magnetic field, 10, 24
 - maintenance safety, 14
 - mechanical hazards, 14
 - proper use, 10
 - symbols, 16
 - waste disposal, 14
- Safety Information (French, FR), 111
- Safety Information (German, DE), 121
- Sample area, 52
- Slide blocks, 50
- Software, 34
 - imaging, 34
- Specifications
 - 2D barcodes, 101
 - Autoload, 95
 - barcode, 95
 - cards, 101
 - chemical compatibility, 103
 - computer, 91
 - gripper, 93
 - instrument, 89
 - instrument dimensions, 102
 - labware, 92
 - operating, 90
 - sample barcodes, 100
 - storage and transport, 91
 - tips, 93
- Stop hooks, 50
- STR reaction
 - loading, 48
 - PCR setup, 39
- Symbols, 16
- Teaching needles, 73, 79
- Technical assistance, 7
- Tip recognition, 30
- Tips, 29, 32
 - liquid level detection, 32
 - packaging, 32
 - specifications, 93
- Trace files, 68
- Tracks, 20
- Training, 37
- User
 - requirements, 8
- Vacuum system, 24
- Verification, 105
- Version management, 7
- Voltage, 36
- Warnings, 9
- Warnung, 121
- Waste Electrical and Electronic Equipment (WEEE), 106
- Waste trough, 49
- Worklist, 42
 - barcode, 43
- XY-table, 24

This page intentionally left blank

This page intentionally left blank

Trademarks: QIAGEN®, Sample to Insight®, HotStarTaq®, Investigator® (QIAGEN Group); ABgene® (Thermo Fisher Group); Applied Biosystems®, MicroAmp® (Life Technologies Corporation); Bio-Rad®, Hard-Shell® (Bio-Rad Laboratories, Inc.); Copan®, NUCLEIC-CARD™ (Copan Italia S.P.A.); Core™, Intel® (Intel Corporation); Dell®, OptiPlex® (Dell, Inc.); EasiCollect™, FTA™ (GE Healthcare); Microsoft®, Windows® (Microsoft Corporation); iSWAP® (Hamilton Company); Kalrez®, Viton® (DuPont Performance Elastomers L.L.C.); Makrolon® (Bayer Aktiengesellschaft); Sarstedt® (Sarstedt AG and Co.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

HB-2123-001 1100450 04-2016 © 2016 QIAGEN, all rights reserved.

Ordering www.qiagen.com/shop | Technical Support support.qiagen.com | Website www.qiagen.com