



April 2023

QIAseq[®] Targeted cfDNA Ultra Handbook

For ultrasensitive targeted next-generation sequencing (NGS) of cfDNA for Illumina[®] NGS systems

Contents

Kit Contents	4
Shipping and Storage	6
Intended Use	6
Safety Information	7
Quality Control	7
Introduction	8
Principle and procedure	9
Equipment and Reagents to Be Supplied by User	14
Important Notes	15
Protocol: DNA End Prepare	21
Protocol: Adapter Ligation	24
Protocol: Target Enrichment	28
Protocol: Universal PCR	32
Recommendations: Library QC and Quantification	38
NGS library QC	38
Recommended setting for checking QIAseq cfDNA Ultra library on QIAxcel Advanced and Connect systems	39
Preferred library quantification method	40
Protocol: Sequencing Setup on Illumina MiSeq, NextSeq 500/550, NextSeq 1000/2000, MiniSeq and NovaSeq	41
Troubleshooting Guide	45
Appendix A: Combining an Existing Panel with a Booster Panel	47

Appendix B: Data Analysis Using QIAGEN's CLC Genomics Workbench	48
Ordering Information	49
Document Revision History	52

Kit Contents

QIAseq Targeted cfDNA Ultra Catalog no. No. of samples	334051 12	334055 96	Custom 334075 96
One pool of region-specific primers	60 µL	480 µL	480 µL
ERA Buffer, 10x	60 µL	480 µL	480 µL
ERA Enzyme	60 µL	480 µL	480 µL
UPH Ligation Buffer, 2.5x	350 µL	2x1250 µL	2x1250 µL
DNA Ligase	75 µL	600 µL	600 µL
HiFi Ultra Buffer, 5x	144 µL	2x768 µL	2x768 µL
HiFi Ultra Polymerase	40 µL	315 µL	315 µL
Nuclease-free Water	1.5 mL	10 mL	10 mL
QIAseq Bead Binding Buffer (one bottle)	7 mL	54 mL	54 mL
QIAseq Beads (one bottle)	10 mL	38.4 mL	38.4 mL

QIAseq Targeted cfDNA Ultra Booster Catalog no. No. of samples	334085 96
One pool of region-specific primers	80 µL

Cat. no.	Product name	Total number of primers*	Panel size (bases)
YHS-001Z	Breast Cancer Ultra Panel	1045	21,960
YHS-002Z	Colorectal Cancer Ultra Panel	529	6618
YHS-005Z	Lung Cancer Ultra Panel	948	26,448
YHS-101Z	Actionable Solid Tumor Ultra Panel	2284	44,012
YHS-003Z	Myeloid Neoplasms Ultra Panel	2302	32,716

* The number of primers in Custom and Booster panels is represented by the last digits of the catalog number. For example, a custom panel with catalog number CYHS-00100Z-1256 has 1256 primers.

QIAseq Targeted cfDNA Ultra 96-Unique Dual Indices

QIAseq Targeted cfDNA Ultra UDI Set* Catalog no. No. of samples	Set A 334165 96	Set B 334175 96	Set C 334185 96	Set D 334195 96
QUDI-96AA †	9 µL	N/A	N/A	N/A
QUDI-96BA †	N/A	9 µL	N/A	N/A
QUDI-96CA †	N/A	N/A	9 µL	N/A
QUDI-96DA †	N/A	N/A	N/A	9 µL
AdP-DP-Phased Adapter	144 µL	144 µL	144 µL	144 µL
SmP-IL5 TEPCR-F Primer	250 µL	250 µL	250 µL	250 µL

* 10 bp dual indices.

† Index Primer Plate (DNA Ultra UDI Set A, B, C and D Final Plate); each plate contains 96 pairs of sample index primers plus universal primers, with each well corresponding to one pair of UDI sample index; each index is single use.

QIAseq Targeted cfDNA Ultra 12 Unique Dual Indices

QIAseq Targeted cfDNA Ultra UDI (12)* Catalog no. No. of samples	334151 12
QUDI-12A †	9 µL
AdP-DP-Phased Adapter	18 µL
SmP-IL5 TEPCR-F Primer	31 µL

* 10 bp dual indices.

† Index Primer Plate (DNA Ultra UDI 12 Index Final Plate); each plate contains 12 pairs of sample index primers plus universal primers, with each well corresponding to one pair of UDI sample index; each index is single use.

Shipping and Storage

The QIAseq Targeted cfDNA Ultra Kits (except QIAseq Beads and QIAseq Bead Binding Buffer) are shipped on dry ice and should be stored at -30°C to -15°C in a constant-temperature freezer upon arrival. The QIAseq Beads and QIAseq Bead Binding Buffer are shipped on cold packs and should be stored at $2-8^{\circ}\text{C}$ upon arrival.

The QIAseq Targeted cfDNA Ultra Index Kits are shipped on dry ice and should be stored at -30°C to -15°C upon arrival.

When stored correctly, the QIAseq Targeted cfDNA Ultra Kits are good until the expiration date printed on the kit label.

Intended Use

The QIAseq Targeted cfDNA Ultra Kits are intended for molecular biology applications. These products are not intended for the diagnosis, prevention or treatment of a disease.

All due care and attention should be exercised in the handling of the products. We recommend all users of QIAGEN® products to adhere to the NIH guidelines that have been developed for recombinant DNA experiments or to other applicable guidelines.

Safety Information

When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at www.qiagen.com/safety, where you can find, view and print the SDS for each QIAGEN kit and kit component.

Quality Control

In accordance with QIAGEN's ISO-certified Quality Management System, each lot of the QIAseq Targeted cfDNA Ultra is tested against predetermined specifications, to ensure consistent product quality.

Introduction

The QIAseq Targeted cfDNA Ultra enables streamlined Sample to Insight[®], targeted next-generation sequencing (NGS) of cell-free DNA (cfDNA). This highly optimized, automation-friendly solution facilitates ultrasensitive variant detection down to 0.1% by using integrated unique molecular indices (UMIs) and high fidelity chemistry from biofluids within 8 hours and is coupled with an error-correction data analysis. The required amount of template for a single QIAseq Targeted cfDNA Ultra sequencing reaction ranges from 5 to 80 ng of cfDNA.

The NGS of DNA is a powerful tool for the detection of genetic variations, including somatic mutations, single-nucleotide polymorphisms, copy-number variation, and small insertions and deletions. Target-enrichment technology enhances DNA NGS by enabling users to sequence specific regions of interest — instead of the entire genome — which effectively increases sequencing depth and sample throughput while minimizing cost. Many commercially available target enrichment, library preparation and sequencing methods use regular DNA polymerase and amplification processes that introduce substantial bias and artifacts. This results in artifactual errors that greatly limit the detection of true low-frequency variants below 0.5% in heterogeneous samples, such as cfDNA. The QIAseq Targeted cfDNA Ultra overcomes these biases and artifacts with the use of a highly optimized, high-fidelity reaction chemistry and by incorporating UMIs into a single gene-specific, primer-based targeted enrichment process.

The QIAseq Targeted cfDNA Ultra has been optimized with a specially formulated enrichment chemistry to achieve highly efficient enrichment on both regular and GC-rich regions at high multiplex levels. In addition, the QIAseq Targeted cfDNA Ultra library can be sequenced with Illumina default sequencing primers and is compatible with most medium- and high-throughput Illumina sequencers.

Our data analysis tools have been developed to perform all steps necessary to generate a DNA sequence variant report from NGS data. Collectively, the QIAseq Targeted cfDNA Ultra

is a Sample to Insight solution for precision variant detection of targeted genomic regions using NGS (Figure 1).




Figure 1. Overview of the Sample to Insight NGS workflow for cell-free DNA with the QIAseq Targeted cfDNA Ultra.

This complete Sample to Insight procedure begins with DNA extraction followed by library construction and target enrichment with the QIAseq Targeted cfDNA Ultra. After NGS, data analysis is performed using the QIAGEN CLC Genomics Workbench. Ultimately, detected variants can be interpreted with the QIAGEN Clinical Insight Interpret for clinical relevant information.

Principle and procedure

The QIAseq Targeted cfDNA Ultra Kits are provided as single-tube primer mixes, with up to 20,000 primers per panel. The QIAseq Targeted cfDNA Ultra Kit is designed to enrich selected genes and regions using 5–80 ng cfDNA (Figure 2). Lower input amounts are possible; however, this will lead to fewer sequenced UMIs and reduced variant detection sensitivity.




Figure 2. QIAseq Targeted cfDNA Ultra workflow.

End preparation

cfDNA samples are first end repaired and A tailed within a single, controlled multienzyme reaction. The prepared DNA fragments are then ligated with a sequencing platform-specific adapter containing UMIs.

UMI assignment

Prior to target enrichment and library amplification, each original DNA molecule is assigned a unique sequence or index, commonly referred to as a UMI. This assignment is accomplished by ligating end-prepared DNA with an adapter containing a 12-base, fully random sequence (i.e., the UMI). Statistically, this process provides 4^{12} possible indices per adapter, and each DNA molecule in the sample receives a unique UMI sequence.

Target enrichment and final library construction

Target enrichment is performed after UMI assignment to ensure that DNA molecules containing UMIs are sufficiently enriched in the sequenced library. For enrichment, ligated DNA molecules are subject to several cycles of targeted PCR using one region-specific primer and one universal primer complementary to the adapter. A universal PCR is ultimately carried out to amplify the library and add platform-specific adapter sequences and sample indices.

NGS adapter and index technologies

The QIAseq Targeted cfDNA Ultra Panels use unique dual index (UDI) primers for sample indexing. The UDI primers significantly reduce the risk of index-bleeding issues associated with different Illumina sequencing instruments, as well as reducing the impact of low-level contamination during oligo synthesis and kit manufacturing, as well as carryover on the Illumina sequencing instrument itself. Hence, each sample is assigned two unique indices to overcome the error introduced by image analysis, sequencing error, demultiplexing, and oligo synthesis contamination, which reduces the mis-assignment of reads to wrong samples.

Next-generation sequencing

The QIAseq Targeted cfDNA Ultra is compatible with most medium- and high-throughput sequencers including Illumina NGS systems (MiniSeq®, MiSeq®, NextSeq® 500/550, NextSeq 1000/2000, HiSeq® 2500, HiSeq 3000/4000, and NovaSeq® 6000).

Principle of variant detection with UMIs

The principle of variant detection with UMIs is described in Figure 3. Due to intrinsic noise and sequence-dependent bias, indexed molecules may be amplified unevenly across the target regions. However, target region coverage can be better achieved by counting the number of UMIs rather than counting the number of total reads for each region. Sequence reads having different UMIs represent different original molecules, while sequence reads having the same UMIs are the result of PCR duplication from one original molecule. Errors from PCR amplification and from the sequencing process may also be present in final reads that lead to false positive variants in sequencing results. These artifactual variants can be greatly reduced by calling variants across all reads within a unique UMI instead of picking up variants at the original read level.




Figure 3. Principle of variant detection with UMIs. Each original molecule is tagged by a UMI. True variants are those mutations present in the majority of reads within a UMI, while false positives are mutations present in only one or a few reads within a UMI.

Data analysis

The data for the QIAseq Targeted cfDNA Ultra can be analyzed using the Biomedical Genomics Analysis plugin to the QIAGEN CLC Genomics Workbench. The plugin provides workflows and tools for all steps from the initial data processing and quality assurance through data analyses, annotation and reporting. A detailed guide to UMI-directed variant detection in CLC Genomics Workbench can be found in *Biomedical Genomics Analysis Plugin User Manual*.

All detected variants can be further interpreted using QIAGEN's Clinical Insight (QCI®) Interpret.

Equipment and Reagents to Be Supplied by User

When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles. For more information, consult the appropriate safety data sheets (SDSs), available from the product supplier.

In addition to the QIAseq Targeted cfDNA Ultra and Index Kit, the following are required:

- Ethanol, 80% (made fresh each day before use) *
- Nuclease-free pipette tips and tubes
- 1.5 mL LoBind® tubes (Eppendorf®, cat. no. 022431021)
- PCR tubes (0.2 mL individual PCR tubes [VWR, cat. no. 20170-012], tube strips [VWR, cat. no. 93001 118]) or 96-well PCR plates and caps
- Ice
- Microcentrifuge
- Thermal cycler
- Multichannel pipettes
- Single-channel pipettes
- QIAxcel® Connect System (cat. no. 9003110) or QIAxcel Advanced System (cat. no. 9001941) or Agilent® 2100 Bioanalyzer® (Agilent, cat. no. G2939BA) or Agilent TapeStation® (Agilent, cat. no. G2991AA)
- QIAxcel DNA High Sensitivity Kit (cat. No. 929012), QIAxcel DNA High Resolution Kit (QIAGEN, cat. no. 929002) or Agilent High Sensitivity DNA Kit (Agilent, cat. no. 5067-4626)
- DynaMag™-96 Side Magnet (Thermo Fisher Scientific Inc., cat. no. 12331D)
- QIAxpert® Instrument (cat. no. 9002340) or Thermo Fisher Scientific Qubit Fluorometer

* Do not use denatured alcohol, which contains other substances, such as methanol or methyl ethyl ketone.

Important Notes

For optimal results, all DNA samples should demonstrate consistent quality according to the following criteria.

DNA isolation and quality check

The most important prerequisite for DNA sequence analysis is consistent, high-quality DNA from every experimental sample. Therefore, sample handling and DNA isolation procedures are critical to the success of the experiment. Residual traces of proteins, salts or other contaminants may either degrade the DNA or decrease the efficiency of, if not block completely, the enzymatic activity necessary for optimal targeted enrichment. Sample purity can be checked with the QIAxpert System.

The QIAGEN kits listed in Table 1 are recommended for the preparation of cfDNA from serum or plasma samples. For best results, all DNA samples should be resuspended in DNase-free water or, alternatively, in DNase-free 10 mM Tris buffer, pH 8.0.

Important: Do not use DEPC-treated water.

Note: If genomic DNA must be isolated from biological samples for which kits are not available, please contact QIAGEN Technical Support for suggestions.

Table 1. Recommended kits for purification of cfDNA

Kit	Starting material	Cat. no.
QIAamp® Circulating Nucleic Acid Kit	Up to 5 mL of serum or plasma	55114
QIAamp MinElute ccfDNA Mini Kit	Up to 4 mL of serum or plasma	55204
QIAamp MinElute ccfDNA Midi Kit	Up to 10 mL of serum or plasma	55284
EZ1&2 ccfDNA Kit	Automated extraction from up to 8 mL serum or plasma	954854

For best results, all DNA samples should also demonstrate consistent quality according to the following criteria.

cfDNA quality

cfDNA quality can be checked using the QIAxcel System, Agilent Bioanalyzer or Agilent TapeStation. Good cfDNA samples should have a distinct cfDNA peak and a minimal amount of cellular genomic DNA that is of high-molecular weight.

cfDNA quantification

If the cfDNA is free of high-molecular weight cellular DNA, its concentration can be measured by fluorescence-based quantification, such as with the Qubit Fluorometer. If a significant amount of cellular DNA is present in cfDNA sample, it is recommended to measure the amount based on the cfDNA peak with either the QIAxcel System or the Bioanalyzer.

DNA input amount and sequencing depth

The number of UMIs captured from the original DNA sample is related to the DNA input amount and sequencing depth. Adequate sequencing of captured UMIs is necessary for UMI-based variant detection and requires relatively deep sequencing coverage. Table 2 provides guidance on variant detection with cfDNA amounts at different depths of coverage.

Table 2. Suggested cfDNA input amount and sequencing depth for variant detection*

Variant frequency (%)	Input (ng)	Read pairs per UMI	Mean read
0.5	30	4	24,000
0.2	30	4	24,000
0.1	60	4	48,000

* Variant detection is based on 80–90% sensitivity on the entire panel region of the QIAseq Targeted cfDNA Ultra.

Variant detection

The number of UMIs sequenced directly impacts the sensitivity of variant detection. Therefore, low-frequency mutation detection usually requires more DNA input and sequencing at deeper coverage (i.e., more reads per UMI) to generate a sufficient amount of UMIs.

Sequencing capacity and sample multiplex level

Sample multiplexing is one of the most important NGS tools for increasing throughput and reducing costs. It works by combining multiple samples to be processed together in a single sequencing run; as a consequence, sequencing reads need to be demultiplexed by reassigning each single read to its original source library. This is facilitated by the integration of sample index sequences into the individual library molecules.

The QIAseq Targeted cfDNA Ultra Panels are paired with a fully compatible UDI sample-indexing solution. Each QIAseq Targeted cfDNA Ultra UDI Kit includes one of the following:

- QIAseq Targeted cfDNA Ultra UDI Set A (96): DNA Ultra UDI Set A index primer plate, QUDI-96AA
- QIAseq Targeted cfDNA Ultra UDI Set B (96): DNA Ultra UDI Set B index primer plate, QUID-96BA
- QIAseq Targeted cfDNA Ultra UDI Set C (96): DNA Ultra UDI Set C index primer plate, QUDI-96CA

- QIAseq Targeted cfDNA Ultra UDI Set D (96): DNA Ultra UDI Set D index primer plate, QUDI-96DA
- QIAseq Targeted cfDNA Ultra UDI (12): DNA Ultra UDI 12 Index primer plate, QUDI-12A

The QIAseq Ultra UDI Kits use a fixed combination of two unique bar code motives per sample-index primer pair. Therefore, each single-index motive is only used once on any UDI index primer.

Use of UDI indexes effectively mitigates the risk of read misassignment due to index hopping, which is important for detecting variants at very low frequency. By filtering misassigned reads during the demultiplexing of individual samples, highly accurate output data is generated.

To multiplex more than 96 libraries in a single sequencing run, combine kits with different sets, QIAseq Targeted cfDNA Ultra UDI Set A, B, C or D. For example, combining the unique dual QIAseq Targeted cfDNA Ultra UDI Set A and B (96) kits will allow the generation of 192 libraries with different unique dual sample indexes for 192-plex sequencing.

Sample multiplexing level is determined by the size of the panel, required depth of coverage and sequencing platform read capacity. General guidelines are provided for the number of samples that can be multiplexed with different sequencing platforms, based on panel size and read depth (Table 3 and Table 4). Read depth can be fine-tuned after the first run.

Table 3. Number of multiplexed samples based on panel size with 25,000x mean coverage*

Instrument	Version	Capacity (paired-ends reads)	500 primers	1000 primers	2500 primers	5000 primers
iSeq	v2 reagents	8 M	N/A	N/A	N/A	N/A
MiniSeq	Mid output	16 M	1	N/A	N/A	N/A
MiniSeq	High output	50 M	4	2	N/A	N/A
MiSeq	v2 reagents	30 M	2	1	N/A	N/A
MiSeq	v3 reagents	50 M	4	2	N/A	N/A
NextSeq 500	Mid output	260 M	20	10	4	2
NextSeq 500	High output	800 M	64	32	12	6
NextSeq 1000/2000	P1 flow cell	200 M	16	8	3	1
NextSeq 1000/2000	P2 flow cell	800 M	64	32	12	6
NextSeq 2000	P3 flow cell	2.4 B	192	96	38	19
HiSeq 2500 rapid run	Dual flow cell v2	1200 M	96	48	19	9
HiSeq 3000	8 lanes per flow cell	5 B	400	200	80	40
HiSeq 4000	8 lanes per flow cell	10 B	800	400	160	80
NovaSeq 6000	SP (2 lanes per flow cell)	1.6 B	128	64	25	12
NovaSeq 6000	S1 (2 lanes per flow cell)	3.2 B	256	128	51	25
NovaSeq 6000	S2 (2 lanes per flow cell)	8.2 B	656	328	131	65
NovaSeq 6000	S4 (4 lanes per flow cell)	20 B	1600	800	320	160

* Based on 2 x 149 bp paired-end reads.

N/A: Not applicable no samples can be run.

Table 4. Number of multiplexed samples based on panel size with 50,000x mean coverage*

Instrument	Version	Capacity (paired-ends reads)	500 primers	1000 primers	2500 primers	5000 primers
iSeq	v2 reagent	8 M	N/A	N/A	N/A	N/A
MiniSeq	Mid output	16 M	N/A	N/A	N/A	N/A
MiniSeq	High output	50 M	2	1	N/A	N/A
MiSeq	v2 reagents	30 M	1	N/A	N/A	N/A
MiSeq	v3 reagents	50 M	2	1	N/A	N/A
NextSeq 500	Mid output	260 M	10	5	2	1
NextSeq 500	High output	800 M	32	16	6	3
NextSeq 1000/2000	P1 flow cell	200 M	8	4	1	N/A
NextSeq 1000/2000	P2 flow cell	800 M	32	16	6	3
NextSeq 2000	P3 flow cell	2.4 B	96	48	19	9
HiSeq 2500 rapid run	Dual flow cell v2	1200 M	48	24	9	4
HiSeq 3000	8 Lanes per flow cell	5 B	200	100	40	20
HiSeq 4000	8 Lanes per flow cell	10 B	400	200	80	40
NovaSeq 6000	SP (2 lanes per flow cell)	1.6 B	64	32	12	6
NovaSeq 6000	S1 (2 lanes per flow cell)	3.2 B	128	64	25	12
NovaSeq 6000	S2 (2 lanes per flow cell)	8.2 B	328	164	65	32
NovaSeq 6000	S4 (4 lanes per flow cell)	20 B	800	400	160	80

* Based on 2 x 149 bp paired-end reads for.
N/A: Not applicable, no samples can be run.

NGS read-length recommendations

When using Illumina NGS systems, the QIAseq Targeted cfDNA Ultra UDI libraries require 149 bp paired-end reads and dual 10 bp indices

Protocol: DNA End Prepare

Important points before starting

- This protocol covers all procedures required for the preparation of libraries for Illumina sequencers from cfDNA.
- Before setting up the reaction, it is critical to accurately determine the amount of the input cfDNA. When cfDNA is contaminated with cellular genomic DNA, measuring the concentration of the cfDNA peak with the QIAxcel, Bioanalyzer or similar methods is much more accurate than measuring the concentration of the whole sample.
- Reaction and cleanup procedures can be performed in either PCR tubes or 96-well plates.
- Set up reactions on ice.
- Do not vortex any reagents or reactions.

Procedure: end prepare

1. Thaw nucleic acid samples on ice. Gently mix, briefly centrifuge to collect residual liquid from the sides of the tubes, and then return to ice.
2. Prepare the reagents required for preparing the ends.
 - 2a. Thaw 10x ERA Buffer on ice or, if required, at room temperature (15–25°C). Immediately upon thawing, place the buffer on ice. Keep ERA Enzyme on ice.
 - 2b. Mix all reagents by flicking the tube, and centrifuge briefly.

3. On ice, prepare the end prepare mix according to Table 5. Briefly centrifuge, mix by pipetting up and down 10–12 times and briefly centrifuge again.

Note: In general, increasing the amount of DNA input will improve variant detection sensitivity.

Important: Keep the reaction tubes or plate on ice during the entire reaction setup.

Note: If setting up more than one reaction, prepare a volume of master mix 10% greater than what is required for the total number of reactions.

Table 5. Reaction mix for end prepare

Component	Volume per reaction
cfDNA*	Variable
10x ERA Buffer	2.5 µL
ERA Enzyme	5 µL
Nuclease-free Water	Variable
Total	25 µL

* For cfDNA, add 5–80 ng. We recommend 30 ng or more for detection of variants below 0.5%.

4. Program the thermal cycler according to Table 6. Use the instrument's heated lid.
5. Before adding the tubes or plate to a thermal cycler, start the program. When the thermal cycler reaches 4°C, pause the program.

Important: The thermal cycler must be prechilled and paused at 4°C.

Important: Do not use the heated lid during the 20°C incubation step. Alternatively, the lid temperature can be set at 65°C.

6. Transfer the tubes or plate prepared in step 3 to the prechilled thermal cycler and resume the cycling program.

Table 6. Incubation conditions for fragmentation and end prepare

Step	Incubation temperature (°C)	Incubation time
1	4	1 min
2	20	15 min
3	65	15 min
4	4	Hold

7. Upon completion, allow the thermal cycler to return to 4°C.
8. Place the samples on ice and immediately proceed to “Protocol: Adapter Ligation” (page 24).

Protocol: Adapter Ligation

Important points before starting

- The 25 μ L product from “Protocol: DNA End Prepare” (page 21) is the starting material for this protocol.
- Set up reactions on ice.
- Do not vortex any reagents or reactions.
- The QIAseq Beads are used for ligation reaction cleanup. There is no need to bring beads to room temperature (15–25°C) before use.
- **Important:** Prepare fresh 80% ethanol each day before use.
- Ensure that the QIAseq Beads are thoroughly mixed at all times. This necessitates working quickly and resuspending the beads immediately before use. If a delay in the protocol occurs, simply vortex the beads.

Procedure: adapter ligation

1. Prepare the reagents required for the DNA ligation.
 - 1a. Thaw AdP-DP-Phased Adapter on ice. Thaw UPH Ligation Buffer, 2.5x, on ice or at room temperature (15–25°C). Immediately upon thawing, place the buffer on ice. Keep DNA Ligase on ice.
 - 1b. Mix by flicking the tube, and then centrifuge briefly.
2. Prepare the adapter ligation mix according to Table 7. Briefly centrifuge, mix by pipetting up and down 10–12 times and briefly centrifuge again.

Important: The AdP-DP-Phased Adapter does not contain a sample index; hence, one single adapter is used for all samples.

Note: If setting up more than one reaction, prepare a volume of master mix 10% greater than what is required for the total number of reactions.

Table 7. Reaction mix for adapter ligation

Component	Volume/reaction, <10 ng cfDNA	Volume/reaction, 10–50 ng cfDNA	Volume/reaction, >50 ng cfDNA
End prepare reaction (already in tube)	25 µL	25 µL	25 µL
UPH Ligation Buffer, 2.5x	20 µL	20 µL	20 µL
AdP-DP-Phased Adapter	0.5 µL	1 µL	1.5 µL
DNA Ligase	5 µL	5 µL	5 µL
Total	50.5 µL	51 µL	51.5 µL

3. Program the thermal cycler according to Table 8.

Important: Do not use the heated lid during the 20°C incubation step. Alternatively, the lid temperature can be set at 65°C.

- Before adding the tubes or plate to a thermal cycler, start the program. When the thermal cycler reaches 4°C, pause the program.
- Transfer the tubes or plate prepared in step 2 to the prechilled thermal cycler and resume the cycling program.

Table 8. Incubation conditions for DNA ligation

Step	Incubation temperature (°C)	Incubation time
1	4	1 min
2	20	15 min
3	65	15 min
4	4	Hold

- Upon completion, allow the thermal cycler to return to 4°C.
- Place the samples on ice and immediately proceed to “Procedure: ligation cleanup” (page 26). Alternatively, the samples can be stored at –20°C in a constant-temperature freezer for up to 3 days.

Procedure: ligation cleanup

8. Add 40 μL nuclease-free water to bring each sample to approximately 90 μL .
9. Add 108 μL QIAseq Beads. Mix well by vortexing or pipetting up and down several times. Incubate for 5 min at room temperature (15–25°C).
10. Place the tubes or plate on a magnetic rack for 10 min. Once the solution has cleared, with the beads still on the magnetic stand, carefully remove and discard the supernatant.
Important: Do not discard the beads as they contain the DNA of interest.
11. With the beads still on the magnetic stand, add 80 μL **water** to the beads, then add 80 μL QIAseq bead binding buffer.
12. Take the tubes or plate off the magnetic stand, mix well by vortexing or pipetting up and down several times.
13. Return the tubes or plate to the magnetic rack for 5 min. Once the solution has cleared, with the beads still on the magnetic stand, carefully remove and discard the supernatant.
Important: Do not discard the beads as they contain the DNA of interest.
14. With the beads still on the magnetic stand, add 200 μL 80% ethanol and leave beads in ethanol for 30 seconds. Carefully remove and discard the ethanol.
15. Repeat the ethanol wash.
Important: Completely remove all traces of the ethanol wash after this second wash. Remove the ethanol with a 200 μL pipet first, and then use a 10 μL pipet to remove any residual ethanol.
16. With the beads still on the magnetic stand, air dry at room temperature for 15 min.
Note: Visually inspect that the pellet is completely dry. Ethanol carryover to the target enrichment PCR step will significantly impact enrichment PCR efficiency. Over drying beads will not affect DNA elution.
17. Remove the beads from the magnetic stand, and elute the DNA from the beads by adding 12 μL nuclease-free water. Mix well by pipetting.
18. Return the tubes or plate to the magnetic rack until the solution has cleared.

19. Transfer 10 μL of the supernatant to clean tubes or plate.

Important: Make sure no beads remain in the elution. Bead carryover to the target enrichment PCR step will significantly impact enrichment PCR efficiency.

20. Proceed with "Protocol: Target Enrichment" (page 28). Alternatively, the samples can be stored at -20°C in a constant-temperature freezer for up to 3 days.

Protocol: Target Enrichment

Important points before starting

- The 10 μ L cleaned ligation reaction from “Protocol: Adapter Ligation” (page 24) is the starting material for this protocol.
- Set up reactions on ice.
- Do not vortex any reagents or reactions.
- The QIAseq Beads are used for target enrichment PCR cleanup. There is no need to bring beads to room temperature (15–25°C) before use.
- **Important:** Prepare fresh 80% ethanol each day before use.
- Ensure that the QIAseq Beads are thoroughly mixed at all times. This necessitates working quickly and resuspending the beads immediately before use. If a delay in the protocol occurs, simply vortex the beads.

Procedure: target enrichment

1. Prepare the reagents required for the Target Enrichment PCR (TEPCR).
 - 1a. Thaw HiFi Ultra Buffer, 5x; QIAseq Targeted cfDNA Ultra Panel; and SmP-IL5 TEPCR-F Primer on ice or at room temperature (15–25°C). Immediately upon thawing, place the buffer on ice. Keep HiFi Ultra Polymerase on ice.
 - 1b. Mix all reagents by flicking the tube, and then centrifuge briefly.
2. Prepare the target enrichment mix according to Table 9. Briefly centrifuge, mix by pipetting up and down 7–8 times, and briefly centrifuge again.

Note: If setting up more than one reaction, prepare a volume of master mix 10% greater than what is required for the total number of reactions.

Important: Do not add more HiFi Ultra Polymerase than it is indicated in the Table 9. Adding more could result in low enrichment specificity.

Table 9. Reaction mix for Target Enrichment

Component	Volume per reaction 1–1100 primers	Volume per reaction 1101–5000 primers	Volume per reaction 50001–10000 primers	Volume per reaction >10000 primers
Cleaned Adapter-ligated DNA (already in tube) from “Procedure: Ligation Cleanup” (page 26)	10 µL	10 µL	10 µL	10 µL
HiFi Ultra Buffer, 5x	4 µL	4 µL	4 µL	4 µL
QIAseq Targeted cfDNA Ultra Panel	5 µL	5 µL	5 µL	5 µL
SmP-IL5 TEPCR-F Primer	2 µL	2 µL	2 µL	2 µL
HiFi Ultra Polymerase*	–	0.25 µL	0.5 µL	0.75 µL
HiFi Ultra Polymerase (10x diluted)* †	0.6 µL	–	–	–
Total	21.6 µL	21.25 µL	21.5 µL	21.75 µL

* It is recommended to make a master mix instead of pipetting small amount of HiFi Ultra Polymerase enzyme.

† Immediately before setting up the reaction, dilute HiFi Ultra Polymerase 10 folds with nuclease-free H₂O.

3. Program a thermal cycler using the cycling condition in Table 10 (panel with <1600 primers/tube) or Table 11 (panel with ≥1600 primers/tube).

Table 10. Cycling conditions for target enrichment if number of primers <1600/tube

Step	Time	Temperature (°C)
Initial denaturation	2 min	98
8 cycles	20 s 10 min	98 68
1 cycle	3 min	72
Hold	∞	4

Table 11. Cycling conditions for target enrichment if number of primers ≥ 1600 /tube

Step	Time	Temperature (°C)
Initial denaturation	2 min	98
8 cycles	20 s	98
	10 min	65
1 cycle	3 min	72
Hold	∞	4

- Place the target enrichment reaction in the thermal cycler and start the run.
- After the reaction is complete, place the reactions on ice and proceed to “Procedure: TEPCR cleanup” (below). Alternatively, the samples can be stored at -20°C in a constant-temperature freezer for up to 3 days.

Procedure: TEPCR cleanup

- Add 70 μL nuclease-free water to bring each sample to approximately 90 μL .
- Add 108 μL QIAseq Beads. Mix well by vortexing or pipetting up and down several times. Incubate for 5 min at room temperature ($15\text{--}25^{\circ}\text{C}$).
- Place the tubes or plate on a magnetic rack for 5 min. After the solution has cleared, with the beads still on the magnetic stand, carefully remove and discard the supernatant.

Important: Do not discard the beads as they contain the DNA of interest.

- With the beads still on the magnetic stand, add 200 μL 80% ethanol and leave beads in ethanol for 30 s. Carefully remove and discard the ethanol.
- Repeat the ethanol wash.

Important: Completely remove all traces of the ethanol wash after this second wash. Remove the ethanol with a 200 μL pipet first, and then use a 10 μL pipet to remove any residual ethanol.

- With the beads still on the magnetic stand, air dry at room temperature ($15\text{--}25^{\circ}\text{C}$) for 10 min.

Note: Visually inspect that the pellet is completely dry. Ethanol carryover to the next universal PCR step will significantly impact PCR efficiency. Over drying beads will not affect DNA elution.

7. Remove the beads from the magnetic stand, and elute the DNA from the beads by adding 52 μL nuclease-free water. Mix well by pipetting.
8. Return the tubes or plate to the magnetic rack until the solution has cleared.
9. Transfer 50 μL of the supernatant to clean tubes or plate. Transfer 5 μL to a clean tube or plate for universal PCR setup.

Important: Make sure no beads transfer to the elution. Bead carryover to the universal PCR step will significantly impact PCR efficiency.

Note: Only 5 μL of cleaned TEPCR volume will be used in the universal PCR. This volume is sufficient, as the DNA molecules were amplified in multiple cycles during TEPCR. The rest of the TEPCR volume can be stored at -20°C if needed.

10. Proceed with "Protocol: Universal PCR" (page 32). Alternatively, the samples can be stored at -20°C in a constant-temperature freezer for up to 3 days.

Protocol: Universal PCR

Important points before starting

- The 5 μ L cleaned product from “Procedure: target enrichment” (page 28) is the starting material for this protocol.
- Set up reactions on ice.
- Do not vortex any reagents or reactions.
- **Important:** The QIAseq Targeted cfDNA Ultra UDI (12) or QIAseq Targeted cfDNA Ultra UDI Set A (96), Set B (96), Set C (96) or Set D (96) is used for sample indexing.
- The Index Primer Plate contains predispensed index primer pairs and the universal PCR primer and is sealed with pierceable aluminum heat sealing film.
- Puncture the film using standard 200 μ L pipet tips to transfer the appropriate amount of index primer to the tubes or plate for the universal PCR.
- The QIAseq Beads are used for universal PCR cleanup. There is no need to bring beads to room temperature (15–25°C) before use.
- **Important:** Prepare fresh 80% ethanol each day before use.
- Ensure that the QIAseq Beads are thoroughly mixed at all times. This necessitates working quickly and resuspending the beads immediately before use. If a delay in the protocol occurs, simply vortex the beads.

Procedure: universal PCR

1. Prepare the reagents required for the universal PCR.
 - 1a. Thaw HiFi Ultra Buffer (5x), DNA Ultra UDI 12, and 96-Index Plate on ice or at room temperature (15–25°C). Immediately upon thawing, place the items on ice. Keep HiFi Ultra Polymerase on ice.
 - 1b. Mix by flicking the tube, and then centrifuge briefly.

2. Prepare the universal PCR in the tubes or plate containing the cleaned target-enriched from TEPCR cleanup reaction.
Note: If setting up more than one reaction, prepare a volume of master mix 10% greater than what is required for the total number of reactions.
3. For the QIAseq Targeted cfDNA Ultra UDI plates, pierce the foil seal associated with each well that will be used, and transfer 2 μ L (each well contains a forward primer and a reverse primer, each with a unique index) to the cleaned target-enriched DNA from “Procedure: TEPCR cleanup” (page 30) sample tubes or plate according to Table 12. Mix by pipetting up and down 7–8 times and briefly centrifuge again.

Important: Only one UDI pair should be used per universal PCR.

Important: The QIAseq Targeted cfDNA Ultra UDI index plates are stable for a maximum of 10 freeze–thaw cycles. If all 96 wells have not been used, cover the used wells with foil and return to the freezer. Do not reuse wells from the QIAseq Targeted cfDNA Ultra UDI index plates once the foil seals have been pierced. Reusing wells would risk significant cross-contamination.

Table 12. Reaction components for universal PCR if using QIAseq Targeted cfDNA Ultra UDI (12) or QIAseq Targeted cfDNA Ultra UDI Set A, B, C or D (96)

Component	Volume/reaction (μ L)
Cleaned target-enriched DNA from “Procedure: TEPCR Cleanup” (page 30)	5
HiFi Ultra Buffer, 5x	8
Index primers from QIAseq Targeted cfDNA Ultra UDI index plate*	2
HiFi Ultra Polymerase	2
Nuclease-free water	23
Total	40

* Applies to QIAseq Targeted cfDNA Ultra UDI (12) or QIAseq Targeted cfDNA Ultra UDI Set A, B, C and D (96).

DNA Ultra UDI 12 Index Plate (QUDI-12A) in QIAseq Targeted cDNA Ultra UDI (12)

	1	2	3	4	5	6	7	8	9	10	11	12
A	IL5-QUDI-1 IL7-QUDI-1	IL5-QUDI-9 IL7-QUDI-9										
	IL5-QUDI-2 IL7-QUDI-2	IL5-QUDI-10 IL7-QUDI-10										
B	IL5-QUDI-3 IL7-QUDI-3	IL5-QUDI-11 IL7-QUDI-11										
	IL5-QUDI-4 IL7-QUDI-4	IL5-QUDI-12 IL7-QUDI-12										
C	IL5-QUDI-5 IL7-QUDI-5	IL5-QUDI-19 IL7-QUDI-19										
	IL5-QUDI-6 IL7-QUDI-6	IL5-QUDI-20 IL7-QUDI-20										
D	IL5-QUDI-7 IL7-QUDI-7	IL5-QUDI-21 IL7-QUDI-21										
	IL5-QUDI-8 IL7-QUDI-8	IL5-QUDI-22 IL7-QUDI-22										

DNA Ultra UDI Set A Plate (QUDI-96AA) in QIAseq Targeted cDNA Ultra UDI Set A (96)

	1	2	3	4	5	6	7	8	9	10	11	12
A	IL5-QUDI-1 IL7-QUDI-1	IL5-QUDI-9 IL7-QUDI-9	IL5-QUDI-17 IL7-QUDI-17	IL5-QUDI-25 IL7-QUDI-25	IL5-QUDI-33 IL7-QUDI-33	IL5-QUDI-41 IL7-QUDI-41	IL5-QUDI-49 IL7-QUDI-49	IL5-QUDI-57 IL7-QUDI-57	IL5-QUDI-65 IL7-QUDI-65	IL5-QUDI-73 IL7-QUDI-73	IL5-QUDI-81 IL7-QUDI-81	IL5-QUDI-89 IL7-QUDI-89
	IL5-QUDI-2 IL7-QUDI-2	IL5-QUDI-10 IL7-QUDI-10	IL5-QUDI-18 IL7-QUDI-18	IL5-QUDI-26 IL7-QUDI-26	IL5-QUDI-34 IL7-QUDI-34	IL5-QUDI-42 IL7-QUDI-42	IL5-QUDI-50 IL7-QUDI-50	IL5-QUDI-58 IL7-QUDI-58	IL5-QUDI-66 IL7-QUDI-66	IL5-QUDI-74 IL7-QUDI-74	IL5-QUDI-82 IL7-QUDI-82	IL5-QUDI-90 IL7-QUDI-90
B	IL5-QUDI-3 IL7-QUDI-3	IL5-QUDI-11 IL7-QUDI-11	IL5-QUDI-19 IL7-QUDI-19	IL5-QUDI-27 IL7-QUDI-27	IL5-QUDI-35 IL7-QUDI-35	IL5-QUDI-43 IL7-QUDI-43	IL5-QUDI-51 IL7-QUDI-51	IL5-QUDI-59 IL7-QUDI-59	IL5-QUDI-67 IL7-QUDI-67	IL5-QUDI-75 IL7-QUDI-75	IL5-QUDI-83 IL7-QUDI-83	IL5-QUDI-91 IL7-QUDI-91
	IL5-QUDI-4 IL7-QUDI-4	IL5-QUDI-12 IL7-QUDI-12	IL5-QUDI-20 IL7-QUDI-20	IL5-QUDI-28 IL7-QUDI-28	IL5-QUDI-36 IL7-QUDI-36	IL5-QUDI-44 IL7-QUDI-44	IL5-QUDI-52 IL7-QUDI-52	IL5-QUDI-60 IL7-QUDI-60	IL5-QUDI-68 IL7-QUDI-68	IL5-QUDI-76 IL7-QUDI-76	IL5-QUDI-84 IL7-QUDI-84	IL5-QUDI-92 IL7-QUDI-92
C	IL5-QUDI-5 IL7-QUDI-5	IL5-QUDI-13 IL7-QUDI-13	IL5-QUDI-21 IL7-QUDI-21	IL5-QUDI-29 IL7-QUDI-29	IL5-QUDI-37 IL7-QUDI-37	IL5-QUDI-45 IL7-QUDI-45	IL5-QUDI-53 IL7-QUDI-53	IL5-QUDI-61 IL7-QUDI-61	IL5-QUDI-69 IL7-QUDI-69	IL5-QUDI-77 IL7-QUDI-77	IL5-QUDI-85 IL7-QUDI-85	IL5-QUDI-93 IL7-QUDI-93
	IL5-QUDI-6 IL7-QUDI-6	IL5-QUDI-14 IL7-QUDI-14	IL5-QUDI-22 IL7-QUDI-22	IL5-QUDI-30 IL7-QUDI-30	IL5-QUDI-38 IL7-QUDI-38	IL5-QUDI-46 IL7-QUDI-46	IL5-QUDI-54 IL7-QUDI-54	IL5-QUDI-62 IL7-QUDI-62	IL5-QUDI-70 IL7-QUDI-70	IL5-QUDI-78 IL7-QUDI-78	IL5-QUDI-86 IL7-QUDI-86	IL5-QUDI-94 IL7-QUDI-94
D	IL5-QUDI-7 IL7-QUDI-7	IL5-QUDI-15 IL7-QUDI-15	IL5-QUDI-23 IL7-QUDI-23	IL5-QUDI-31 IL7-QUDI-31	IL5-QUDI-39 IL7-QUDI-39	IL5-QUDI-47 IL7-QUDI-47	IL5-QUDI-55 IL7-QUDI-55	IL5-QUDI-63 IL7-QUDI-63	IL5-QUDI-71 IL7-QUDI-71	IL5-QUDI-79 IL7-QUDI-79	IL5-QUDI-87 IL7-QUDI-87	IL5-QUDI-95 IL7-QUDI-95
	IL5-QUDI-8 IL7-QUDI-8	IL5-QUDI-16 IL7-QUDI-16	IL5-QUDI-24 IL7-QUDI-24	IL5-QUDI-32 IL7-QUDI-32	IL5-QUDI-40 IL7-QUDI-40	IL5-QUDI-48 IL7-QUDI-48	IL5-QUDI-56 IL7-QUDI-56	IL5-QUDI-64 IL7-QUDI-64	IL5-QUDI-72 IL7-QUDI-72	IL5-QUDI-80 IL7-QUDI-80	IL5-QUDI-88 IL7-QUDI-88	IL5-QUDI-96 IL7-QUDI-96

DNA Ultra UDI Set B Plate (QUDI-96BA) in QIAseq Targeted cDNA Ultra UDI Set B (96)

	1	2	3	4	5	6	7	8	9	10	11	12
A	IL5-QUDI-97 IL7-QUDI-97	IL5-QUDI-105 IL7-QUDI-105	IL5-QUDI-113 IL7-QUDI-113	IL5-QUDI-121 IL7-QUDI-121	IL5-QUDI-129 IL7-QUDI-129	IL5-QUDI-137 IL7-QUDI-137	IL5-QUDI-145 IL7-QUDI-145	IL5-QUDI-153 IL7-QUDI-153	IL5-QUDI-161 IL7-QUDI-161	IL5-QUDI-169 IL7-QUDI-169	IL5-QUDI-177 IL7-QUDI-177	IL5-QUDI-185 IL7-QUDI-185
	IL5-QUDI-98 IL7-QUDI-98	IL5-QUDI-106 IL7-QUDI-106	IL5-QUDI-114 IL7-QUDI-114	IL5-QUDI-122 IL7-QUDI-122	IL5-QUDI-130 IL7-QUDI-130	IL5-QUDI-138 IL7-QUDI-138	IL5-QUDI-146 IL7-QUDI-146	IL5-QUDI-154 IL7-QUDI-154	IL5-QUDI-162 IL7-QUDI-162	IL5-QUDI-170 IL7-QUDI-170	IL5-QUDI-178 IL7-QUDI-178	IL5-QUDI-186 IL7-QUDI-186
B	IL5-QUDI-99 IL7-QUDI-99	IL5-QUDI-107 IL7-QUDI-107	IL5-QUDI-115 IL7-QUDI-115	IL5-QUDI-123 IL7-QUDI-123	IL5-QUDI-131 IL7-QUDI-131	IL5-QUDI-139 IL7-QUDI-139	IL5-QUDI-147 IL7-QUDI-147	IL5-QUDI-155 IL7-QUDI-155	IL5-QUDI-163 IL7-QUDI-163	IL5-QUDI-171 IL7-QUDI-171	IL5-QUDI-179 IL7-QUDI-179	IL5-QUDI-187 IL7-QUDI-187
	IL5-QUDI-100 IL7-QUDI-100	IL5-QUDI-108 IL7-QUDI-108	IL5-QUDI-116 IL7-QUDI-116	IL5-QUDI-124 IL7-QUDI-124	IL5-QUDI-132 IL7-QUDI-132	IL5-QUDI-140 IL7-QUDI-140	IL5-QUDI-148 IL7-QUDI-148	IL5-QUDI-156 IL7-QUDI-156	IL5-QUDI-164 IL7-QUDI-164	IL5-QUDI-172 IL7-QUDI-172	IL5-QUDI-180 IL7-QUDI-180	IL5-QUDI-188 IL7-QUDI-188
C	IL5-QUDI-101 IL7-QUDI-101	IL5-QUDI-109 IL7-QUDI-109	IL5-QUDI-117 IL7-QUDI-117	IL5-QUDI-125 IL7-QUDI-125	IL5-QUDI-133 IL7-QUDI-133	IL5-QUDI-141 IL7-QUDI-141	IL5-QUDI-149 IL7-QUDI-149	IL5-QUDI-157 IL7-QUDI-157	IL5-QUDI-165 IL7-QUDI-165	IL5-QUDI-173 IL7-QUDI-173	IL5-QUDI-181 IL7-QUDI-181	IL5-QUDI-189 IL7-QUDI-189
	IL5-QUDI-102 IL7-QUDI-102	IL5-QUDI-110 IL7-QUDI-110	IL5-QUDI-118 IL7-QUDI-118	IL5-QUDI-126 IL7-QUDI-126	IL5-QUDI-134 IL7-QUDI-134	IL5-QUDI-142 IL7-QUDI-142	IL5-QUDI-150 IL7-QUDI-150	IL5-QUDI-158 IL7-QUDI-158	IL5-QUDI-166 IL7-QUDI-166	IL5-QUDI-174 IL7-QUDI-174	IL5-QUDI-182 IL7-QUDI-182	IL5-QUDI-190 IL7-QUDI-190
D	IL5-QUDI-103 IL7-QUDI-103	IL5-QUDI-111 IL7-QUDI-111	IL5-QUDI-119 IL7-QUDI-119	IL5-QUDI-127 IL7-QUDI-127	IL5-QUDI-135 IL7-QUDI-135	IL5-QUDI-143 IL7-QUDI-143	IL5-QUDI-151 IL7-QUDI-151	IL5-QUDI-159 IL7-QUDI-159	IL5-QUDI-167 IL7-QUDI-167	IL5-QUDI-175 IL7-QUDI-175	IL5-QUDI-183 IL7-QUDI-183	IL5-QUDI-191 IL7-QUDI-191
	IL5-QUDI-104 IL7-QUDI-104	IL5-QUDI-112 IL7-QUDI-112	IL5-QUDI-120 IL7-QUDI-120	IL5-QUDI-128 IL7-QUDI-128	IL5-QUDI-136 IL7-QUDI-136	IL5-QUDI-144 IL7-QUDI-144	IL5-QUDI-152 IL7-QUDI-152	IL5-QUDI-160 IL7-QUDI-160	IL5-QUDI-168 IL7-QUDI-168	IL5-QUDI-176 IL7-QUDI-176	IL5-QUDI-184 IL7-QUDI-184	IL5-QUDI-192 IL7-QUDI-192

Figure 4A. Layout of DNA Ultra UDI Index Primer Plate in QIAseq Targeted cDNA Ultra UDI (12) and QIAseq Targeted cDNA Ultra UDI Sets A and B (96). Each well contains one pair of predisposed sample index primers plus universal primers for a single reaction in "Protocol: Universal PCR" (page 32).

DNA Ultra UDI Set C Plate (QUDI-96CA) in QIAseq Targeted cfDNA Ultra UDI Set C (96)

	1	2	3	4	5	6	7	8	9	10	11	12
A	IL5-QUDI-193	IL5-QUDI-201	IL5-QUDI-209	IL5-QUDI-217	IL5-QUDI-225	IL5-QUDI-233	IL5-QUDI-241	IL5-QUDI-249	IL5-QUDI-257	IL5-QUDI-265	IL5-QUDI-273	IL5-QUDI-281
	IL7-QUDI-193	IL7-QUDI-201	IL7-QUDI-209	IL7-QUDI-217	IL7-QUDI-225	IL7-QUDI-233	IL7-QUDI-241	IL7-QUDI-249	IL7-QUDI-257	IL7-QUDI-265	IL7-QUDI-273	IL7-QUDI-281
B	IL5-QUDI-194	IL5-QUDI-202	IL5-QUDI-210	IL5-QUDI-218	IL5-QUDI-226	IL5-QUDI-234	IL5-QUDI-242	IL5-QUDI-250	IL5-QUDI-258	IL5-QUDI-266	IL5-QUDI-274	IL5-QUDI-282
	IL7-QUDI-194	IL7-QUDI-202	IL7-QUDI-210	IL7-QUDI-218	IL7-QUDI-226	IL7-QUDI-234	IL7-QUDI-242	IL7-QUDI-250	IL7-QUDI-258	IL7-QUDI-266	IL7-QUDI-274	IL7-QUDI-282
C	IL5-QUDI-195	IL5-QUDI-203	IL5-QUDI-211	IL5-QUDI-219	IL5-QUDI-227	IL5-QUDI-235	IL5-QUDI-243	IL5-QUDI-251	IL5-QUDI-259	IL5-QUDI-267	IL5-QUDI-275	IL5-QUDI-283
	IL7-QUDI-195	IL7-QUDI-203	IL7-QUDI-211	IL7-QUDI-219	IL7-QUDI-227	IL7-QUDI-235	IL7-QUDI-243	IL7-QUDI-251	IL7-QUDI-259	IL7-QUDI-267	IL7-QUDI-275	IL7-QUDI-283
D	IL5-QUDI-196	IL5-QUDI-204	IL5-QUDI-212	IL5-QUDI-220	IL5-QUDI-228	IL5-QUDI-236	IL5-QUDI-244	IL5-QUDI-252	IL5-QUDI-260	IL5-QUDI-268	IL5-QUDI-276	IL5-QUDI-284
	IL7-QUDI-196	IL7-QUDI-204	IL7-QUDI-212	IL7-QUDI-220	IL7-QUDI-228	IL7-QUDI-236	IL7-QUDI-244	IL7-QUDI-252	IL7-QUDI-260	IL7-QUDI-268	IL7-QUDI-276	IL7-QUDI-284
E	IL5-QUDI-197	IL5-QUDI-205	IL5-QUDI-213	IL5-QUDI-221	IL5-QUDI-229	IL5-QUDI-237	IL5-QUDI-245	IL5-QUDI-253	IL5-QUDI-261	IL5-QUDI-269	IL5-QUDI-277	IL5-QUDI-285
	IL7-QUDI-197	IL7-QUDI-205	IL7-QUDI-213	IL7-QUDI-221	IL7-QUDI-229	IL7-QUDI-237	IL7-QUDI-245	IL7-QUDI-253	IL7-QUDI-261	IL7-QUDI-269	IL7-QUDI-277	IL7-QUDI-285
F	IL5-QUDI-198	IL5-QUDI-206	IL5-QUDI-214	IL5-QUDI-222	IL5-QUDI-230	IL5-QUDI-238	IL5-QUDI-246	IL5-QUDI-254	IL5-QUDI-262	IL5-QUDI-270	IL5-QUDI-278	IL5-QUDI-286
	IL7-QUDI-198	IL7-QUDI-206	IL7-QUDI-214	IL7-QUDI-222	IL7-QUDI-230	IL7-QUDI-238	IL7-QUDI-246	IL7-QUDI-254	IL7-QUDI-262	IL7-QUDI-270	IL7-QUDI-278	IL7-QUDI-286
G	IL5-QUDI-199	IL5-QUDI-207	IL5-QUDI-215	IL5-QUDI-223	IL5-QUDI-231	IL5-QUDI-239	IL5-QUDI-247	IL5-QUDI-255	IL5-QUDI-263	IL5-QUDI-271	IL5-QUDI-279	IL5-QUDI-287
	IL7-QUDI-199	IL7-QUDI-207	IL7-QUDI-215	IL7-QUDI-223	IL7-QUDI-231	IL7-QUDI-239	IL7-QUDI-247	IL7-QUDI-255	IL7-QUDI-263	IL7-QUDI-271	IL7-QUDI-279	IL7-QUDI-287
H	IL5-QUDI-200	IL5-QUDI-208	IL5-QUDI-216	IL5-QUDI-224	IL5-QUDI-232	IL5-QUDI-240	IL5-QUDI-248	IL5-QUDI-256	IL5-QUDI-264	IL5-QUDI-272	IL5-QUDI-280	IL5-QUDI-288
	IL7-QUDI-200	IL7-QUDI-208	IL7-QUDI-216	IL7-QUDI-224	IL7-QUDI-232	IL7-QUDI-240	IL7-QUDI-248	IL7-QUDI-256	IL7-QUDI-264	IL7-QUDI-272	IL7-QUDI-280	IL7-QUDI-288

DNA Ultra UDI Set D Plate (QUDI-96DA) in QIAseq Targeted cfDNA Ultra UDI Set D (96)

	1	2	3	4	5	6	7	8	9	10	11	12
A	IL5-QUDI-289	IL5-QUDI-297	IL5-QUDI-305	IL5-QUDI-313	IL5-QUDI-321	IL5-QUDI-329	IL5-QUDI-337	IL5-QUDI-345	IL5-QUDI-353	IL5-QUDI-361	IL5-QUDI-369	IL5-QUDI-377
	IL7-QUDI-289	IL7-QUDI-297	IL7-QUDI-305	IL7-QUDI-313	IL7-QUDI-321	IL7-QUDI-329	IL7-QUDI-337	IL7-QUDI-345	IL7-QUDI-353	IL7-QUDI-361	IL7-QUDI-369	IL7-QUDI-377
B	IL5-QUDI-290	IL5-QUDI-298	IL5-QUDI-306	IL5-QUDI-314	IL5-QUDI-322	IL5-QUDI-330	IL5-QUDI-338	IL5-QUDI-346	IL5-QUDI-354	IL5-QUDI-362	IL5-QUDI-370	IL5-QUDI-378
	IL7-QUDI-290	IL7-QUDI-298	IL7-QUDI-306	IL7-QUDI-314	IL7-QUDI-322	IL7-QUDI-330	IL7-QUDI-338	IL7-QUDI-346	IL7-QUDI-354	IL7-QUDI-362	IL7-QUDI-370	IL7-QUDI-378
C	IL5-QUDI-291	IL5-QUDI-299	IL5-QUDI-307	IL5-QUDI-315	IL5-QUDI-323	IL5-QUDI-331	IL5-QUDI-339	IL5-QUDI-347	IL5-QUDI-355	IL5-QUDI-363	IL5-QUDI-371	IL5-QUDI-379
	IL7-QUDI-291	IL7-QUDI-299	IL7-QUDI-307	IL7-QUDI-315	IL7-QUDI-323	IL7-QUDI-331	IL7-QUDI-339	IL7-QUDI-347	IL7-QUDI-355	IL7-QUDI-363	IL7-QUDI-371	IL7-QUDI-379
D	IL5-QUDI-292	IL5-QUDI-300	IL5-QUDI-308	IL5-QUDI-316	IL5-QUDI-324	IL5-QUDI-332	IL5-QUDI-340	IL5-QUDI-348	IL5-QUDI-356	IL5-QUDI-364	IL5-QUDI-372	IL5-QUDI-380
	IL7-QUDI-292	IL7-QUDI-300	IL7-QUDI-308	IL7-QUDI-316	IL7-QUDI-324	IL7-QUDI-332	IL7-QUDI-340	IL7-QUDI-348	IL7-QUDI-356	IL7-QUDI-364	IL7-QUDI-372	IL7-QUDI-380
E	IL5-QUDI-293	IL5-QUDI-301	IL5-QUDI-309	IL5-QUDI-317	IL5-QUDI-325	IL5-QUDI-333	IL5-QUDI-341	IL5-QUDI-349	IL5-QUDI-357	IL5-QUDI-365	IL5-QUDI-373	IL5-QUDI-381
	IL7-QUDI-293	IL7-QUDI-301	IL7-QUDI-309	IL7-QUDI-317	IL7-QUDI-325	IL7-QUDI-333	IL7-QUDI-341	IL7-QUDI-349	IL7-QUDI-357	IL7-QUDI-365	IL7-QUDI-373	IL7-QUDI-381
F	IL5-QUDI-294	IL5-QUDI-302	IL5-QUDI-310	IL5-QUDI-318	IL5-QUDI-326	IL5-QUDI-334	IL5-QUDI-342	IL5-QUDI-350	IL5-QUDI-358	IL5-QUDI-366	IL5-QUDI-374	IL5-QUDI-382
	IL7-QUDI-294	IL7-QUDI-302	IL7-QUDI-310	IL7-QUDI-318	IL7-QUDI-326	IL7-QUDI-334	IL7-QUDI-342	IL7-QUDI-350	IL7-QUDI-358	IL7-QUDI-366	IL7-QUDI-374	IL7-QUDI-382
G	IL5-QUDI-295	IL5-QUDI-303	IL5-QUDI-311	IL5-QUDI-319	IL5-QUDI-327	IL5-QUDI-335	IL5-QUDI-343	IL5-QUDI-351	IL5-QUDI-359	IL5-QUDI-367	IL5-QUDI-375	IL5-QUDI-383
	IL7-QUDI-295	IL7-QUDI-303	IL7-QUDI-311	IL7-QUDI-319	IL7-QUDI-327	IL7-QUDI-335	IL7-QUDI-343	IL7-QUDI-351	IL7-QUDI-359	IL7-QUDI-367	IL7-QUDI-375	IL7-QUDI-383
H	IL5-QUDI-296	IL5-QUDI-304	IL5-QUDI-312	IL5-QUDI-320	IL5-QUDI-328	IL5-QUDI-336	IL5-QUDI-344	IL5-QUDI-352	IL5-QUDI-360	IL5-QUDI-368	IL5-QUDI-376	IL5-QUDI-384
	IL7-QUDI-296	IL7-QUDI-304	IL7-QUDI-312	IL7-QUDI-320	IL7-QUDI-328	IL7-QUDI-336	IL7-QUDI-344	IL7-QUDI-352	IL7-QUDI-360	IL7-QUDI-368	IL7-QUDI-376	IL7-QUDI-384

Figure 4B. Layout of DNA Ultra UDI Index Primer Plate in QIAseq Targeted cfDNA Ultra UDI Sets C and D (96). Each well contains one pair of predispersed sample index primers plus universal primers for a single reaction in “Protocol: Universal PCR (page 32).

- Program a thermal cycler using the cycling conditions in Table 13 and Table 14 indicates the cycle number to use dependent on the number of primers in the pool.

Table 13. Incubation conditions for universal PCR

Step	Time	Temperature (°C)
Initial denaturation	2 min	98
Number of cycles (see Table 14)	20 s	98
	1 min	60
1 cycle	3 min	72
Hold	∞	4

Table 14. Amplification cycles for universal PCR

Primers per pool	Cycle number
6–24	26
25–96	24
97–288	22
289–1056	21
1057–3072	20
3073–5999	18
6000–12,000	17
≥12,001	16

5. After the reaction is complete, place the reactions on ice and proceed to “Procedure: cleanup of universal PCR” (below). Alternatively, the samples can be stored at -20°C in a constant-temperature freezer for up to 3 days.

Procedure: cleanup of universal PCR

6. Add 50 μL nuclease-free water to bring each sample to 90 μL .
7. Add 108 μL QIAseq Beads to the completed universal PCR; mix well by vortexing or pipetting up and down several times.
8. Incubate for 5 min at room temperature ($15\text{--}25^{\circ}\text{C}$).
9. Place the tubes or plate on magnetic rack for 5 min to separate beads from supernatant. Once the solution has cleared, with the tubes or plate still on the magnetic stand, carefully remove and discard the supernatant.
Important: Do not discard the beads as they contain the DNA of interest.
10. With the tubes or plate still on the magnetic stand, add 200 μL 80% ethanol and leave beads in ethanol for 30 seconds. Carefully remove and discard the ethanol.

11. Repeat the ethanol wash.

Important: Completely remove all traces of the ethanol wash after this second wash. Remove the ethanol with a 200 μL pipet first, and then use a 10 μL pipet to remove any residual ethanol.

12. With the tubes or plate still on the magnetic stand, air dry at room temperature (15–25°C) for 10 min.

Note: Visually inspect that the pellet is completely dry. Over drying the beads will not affect DNA elution.

13. Remove the tubes or plate from the magnetic stand, and elute the DNA from the beads by adding 30 μL nuclease-free water. Mix well by pipetting or vortexing.

14. Return the tubes or plate to the magnetic rack until the solution has cleared.

15. Transfer 28 μL supernatant to clean tubes or plate.

16. Proceed to “Recommendations: Library QC and Quantification” (page 38). Alternatively, the library can be stored at -30°C to -15°C in a constant-temperature freezer. Amplified libraries are stable for several months at -30°C to -15°C . Once quantification is performed, proceed to “Protocol: Sequencing Setup on Illumina MiSeq, NextSeq 500/550, MiniSeq and NovaSeq” (page 41).

Recommendations: Library QC and Quantification

NGS library QC

After the library is constructed and purified, QC can be performed with QIAGEN's QIAxcel Systems or Agilent's Bioanalyzer or TapeStation to check for the correct size distribution of the library fragments and for the absence of primer dimers (approximately <200 bp) and concentration. Majority library fragments prepared for Illumina instruments demonstrate a size distribution between 200 and 400 bp (Figure 5 and Figure 6). Library overamplification is normal (Figure 6B), and this should not affect the sequencing results. Overamplified libraries are usually single-stranded libraries with correct size but appear as "larger fragments" due to secondary structures. Amounts of DNA under the appropriate peaks can be used to quantify the libraries. However, due to the superior sensitivity of qPCR, we recommend quantifying the libraries using the QIAseq Library Quant System, especially when there are overamplified libraries (See "Preferred library quantification method", page 40).

Recommended setting for checking QIAseq cfDNA Ultra library on QIAxcel Advanced and Connect systems

QIAxcel Connect

- QIAxcel DNA High Sensitivity Kit (1200) (cat. no. 929012)
- Use method Default High Sensitivity with 1–10 μ L library

QIAxcel Advanced or QIAxcel Connect

- QIAxcel DNA High Resolution Kit (cat. no. 929002)
- **QX Alignment Marker:** 15 bp/5 kb (cat. no. 929524)
- **QX DNA Size Marker:** 100 bp – 2.5 kb (cat. no. 929559)
- Use the *Application Guide for Low-Concentration Libraries*. To access guides and system files for Library QC, contact QIAGEN Technical Services.





Figure 5. Sample QIAxcel image of QIAseq Targeted cfDNA Ultra libraries for Illumina instruments. The library assessed using QIAxcel illustrates the size of the majority of the library fragments are between 200 and 400 bp.

A: Library (without overamplification) prepared for Illumina instruments



B: Library (with overamplification) prepared for Illumina instruments




Figure 6. Sample Bioanalyzer images of QIAseq Targeted cfDNA Ultra libraries for Illumina instruments. The size of the majority of the library fragments are between 200 and 400 bp. A: Library without overamplification. B: Library with overamplification as indicated by the “larger fragment” peak.

Preferred library quantification method

The library yield measurements from the Bioanalyzer or TapeStation rely on fluorescence dyes that intercalate into DNA. These dyes cannot discriminate between molecules with or without adapter sequences, yet only complete QIAseq Targeted cfDNA Ultra libraries with full adapter sequences will be sequenced. Due to the superior sensitivity of qPCR, we recommend quantifying the libraries using, QIAGEN’s QIAseq Library Quant Assay Kit (cat. no. 333314), which contains laboratory-verified forward and reverse primers, together with a DNA standard. With this system, the correct dilution of the library can be determined for sequencing. Please refer to the relevant handbook (available at www.qiagen.com) for library quantification.

Protocol: Sequencing Setup on Illumina MiSeq, NextSeq 500/550, NextSeq 1000/2000, MiniSeq and NovaSeq

Important points before starting

- **Important:** Recommendations for library dilution concentrations and library loading concentrations are based on QIAseq Library Quant System (see “Preferred library quantification method”, page 40). If using Qubit for library quantification, load 50% less as starting point to avoid over clustering and adjust accordingly after the first run.
- **Important:** Paired-end sequencing should be used for the QIAseq Targeted cfDNA Ultra on Illumina platform.
- **Important:** To make sequencing preparation more convenient, download Illumina-compatible sample sheets for different sequencing instruments on www.qiagen.com, from the Resources tab of the QIAseq Targeted cfDNA Ultra page.
- **Important:** Paired-end sequencing of 149 bp should be used for QIAseq Targeted cfDNA Ultra UDI libraries and dual 10 bp indices on Illumina platforms.
- **Important:** For two-channel sequencing chemistry platforms such as MiniSeq, NextSeq and NovaSeq, 10% PhiX can be included in the run to improve sequencing quality. For complete instructions on how to denature sequencing libraries, and set up a sequencing run, please refer to the system-specific Illumina documents.

Sequencing preparations for MiSeq with QIAseq Targeted cfDNA Ultra UDI Sets

1. When working with the QIAseq Targeted cfDNA Ultra custom QIAseq 96-Unique Dual Index Sets, use Local Run Manager (LRM) v2 or LRM v3 on the instrument to upload a sample sheet (see the Resources tab of the QIAseq Targeted cfDNA Ultra and download the appropriate template) and proceed with sequencing: Read 1 is 149 bp, Read 2 is 149 bp, and each Index Read is 10 bp.

2. Sample dilution and pooling: Dilute libraries to 2 or 4 nM for MiSeq. Then, combine libraries with different sample indexes in equimolar amounts if similar sequencing depth is needed for each library.

Note: Recommendations for library dilution concentrations are based on QIAseq Library Quant System.

Note: If combining libraries with the same number of primers, pool equal volumes of the individual libraries at 4 nM together. If the libraries have different primer numbers, then combine the libraries at volume ratios according to their number of primers. For example, Library A has 5000 primers at 4 nM and Library B has 600 primers at 4 nM; combining 50 μ L of Library A with 6 μ L of Library B will result in similar coverage depth for both Libraries A and B in the same sequencing run.

3. Library preparation and loading: Prepare and load the library onto a MiSeq according to the *MiSeq System Denature and Dilute Libraries Guide*. The final library concentration is 10–12 pM on the MiSeq.

Note: Recommendations for library loading concentrations are based on QIAseq Library Quant System.

4. Upon completion of the sequencing run, proceed to Appendix B (page 48).

Sequencing preparations for MiniSeq, NextSeq 500/550 and NextSeq1000/2000 with QIAseq Targeted cfDNA Ultra UDI Sets

1. When working with the QIAseq Targeted cfDNA Ultra custom QIAseq 96-Unique Dual Index Sets, use LRM v2 on the instrument to upload a sample sheet (see the Resources tab of the QIAseq Targeted cfDNA Ultra page and download the appropriate template) and proceed with sequencing: Read 1 is 149 bp, Read 2 is 149 bp and each Index Read is 10 bp.
2. Sample dilution and pooling: Dilute libraries to 1 nM for MiniSeq, 0.5, 1, 2 or 4 nM for NextSeq 500/550 and 2 nM for NextSeq 1000/2000. Then, combine libraries with different sample indexes in equimolar amounts if similar sequencing depth is needed for each library.

Note: Recommendations for library dilution concentrations are based on QIAseq Library Quant System.

Note: If combining libraries with the same number of primers, pool equal volumes of the individual libraries at 4 nM together. If the libraries have different primer numbers, then combine the libraries at volume ratios according to their number of primers. For example, Library A has 5000 primers at 4 nM and Library B has 600 primers at 4 nM; combining 50 μ L Library A with 6 μ L Library B will result in similar coverage depth for both Libraries A and B in the same sequencing run.

3. Library preparation and loading: Prepare and load the library onto a MiniSeq, NextSeq 500/550 or NextSeq 1000/2000 according to the *MiniSeq or NextSeq System Denature and Dilute Libraries Guide*. The final library concentration is 1.0–1.2 pM on the MiniSeq or NextSeq 500/550 and 650 pM for NextSeq 1000/2000 onboard denature and dilute.

Note: Recommendations for library loading concentrations are based on QIAseq Library Quant System.

4. Upon completion of the sequencing run, proceed to Appendix B (page 48).

Sequencing preparations for NovaSeq with QIAseq Targeted cfDNA Ultra UDI Index Sets

1. When working with the QIAseq Targeted cfDNA Ultra custom QIAseq 96-Unique Dual Index Sets, upload a sample sheet (see the Resources tab of the QIAseq Targeted DNA Panel page and download the appropriate template) and proceed with sequencing: Read 1 is 149 bp, Read 2 is 149 bp and each Index Read is 10 bp.
2. Sample dilution and pooling: Dilute libraries to 10 nM for NovaSeq. Then, combine libraries with different sample indexes in equimolar amounts if similar sequencing depth is needed for each library.

Note: Recommendations for library dilution concentrations are based on QIAseq Library Quant System.

Note: If combining libraries with the same number of primers, pool equal volumes of the individual libraries at 10 nM together. If the libraries have different primer numbers, then combine the libraries at volume ratios according to their number of primers. For example, Library A has 5000 primers at 10 nM and Library B has 600 primers at 10 nM; combining 50 μ L Library A with 6 μ L Library B will result in similar coverage depth for both Libraries A and B in the same sequencing run.

3. Library preparation and loading: Prepare and load the library onto a NovaSeq according to the *NovaSeq 6000 Sequencing System Guide* (part #1000000019358). The final pooled library concentration recommendation is between 1.0 and 1.5 nM yielding a final loading concentration of between 200 and 300 pM on the NovaSeq.

Note: Recommendations for library loading concentrations are based on QIAseq Library Quant System.

4. Upon completion of the sequencing run, proceed to Appendix B (page 48).

Troubleshooting Guide

This troubleshooting guide may be helpful in solving any problems that may arise. For more information, see also the Frequently Asked Questions page at our Technical Support Center: www.qiagen.com/FAQ/FAQList.aspx. The scientists in QIAGEN Technical Services are always happy to answer any questions you may have about either the information and/or protocols in this handbook or sample and assay technologies (for contact information, visit www.qiagen.com).

Comments and suggestions

Low library yield

- a) Suboptimal reaction conditions due to low cfDNA input Make sure cfDNA is free of cellular genomic DNA. If there is cellular genomic DNA, make sure input is based on cfDNA peak, not all the DNA in the sample.
- b) Inefficient targeted enrichment or universal PCR Check to see if correct thermocycling condition is used during target enrichment or universal PCR. Ensure no ethanol or bead carryover to the PCR.

Unexpected signal peaks

- a) Short peaks less than 200 bp These are primer–dimers from targeted enrichment PCR or universal PCR. The presence of primer dimers indicates either not enough DNA input or inefficient PCR.
- b) Larger library fragments after universal PCR After the universal PCR, library fragments are larger than the intended peak and can be a PCR artifact due to overamplification of the DNA library. Overamplification of the library will not affect the QIAseq Targeted cfDNA Ultra sequencing performance. Decreasing the number of universal PCR cycle numbers can reduce over-amplification.

Sequencing issues

- a) Too low or too high cluster density Accurate library quantification is the key for optimal cluster density on any sequencing instrument. PCR-based quantification method is recommended. Other methods may lead to the incorrect quantification of the library especially when there is over amplification.
- b) Very low clusters passing filter Make sure that the library is accurately quantified and that the correct amount is loaded onto the sequencing instrument.

Comments and suggestions

Variant detection issues

Known variants not detected

Variant detection sensitivity is directly related to the input DNA and read depth. Check Table 2 to see if the required input DNA and read depth are met for the specific variant detection application.

Appendix A: Combining an Existing Panel with a Booster Panel

If additional primers need to be added into an existing panel, a Booster Panel with up to 100 primers can be ordered. To combine the existing panel with a Booster Panel, follow the volume ratio indicated in Table 15.

Table 15. Combining an Existing Panel (at 50 μ L) with a Booster Panel

No. of primers in existing panel	Volume of existing panel to combine (μ L)	Volume of booster panel to combine (μ L)
1–2000	50	5
2001–4000	50	3.75
4001–12,000	50	2.5
12,001–20,000	50	1.25

Appendix B: Data Analysis Using QIAGEN's CLC Genomics Workbench

After sequencing, data from the QIAseq Targeted cfDNA Ultra can be analyzed using the Biomedical Genomics Analysis plugin for CLC Genomics Workbench, which allows you to optimize analysis parameters to your specific panel. The parameters can then be locked for routine use. Contact your account manager for further details.

Ordering Information

Product	Contents	Cat. no.
QIAseq Targeted cfDNA Ultra (12)	All reagents (except indexes) for targeted DNA sequencing; fixed panel for 12 samples; less than 200 genes	334051
QIAseq Targeted cfDNA Ultra (96)	All reagents (except indexes) for targeted DNA sequencing; fixed panel for 96 samples; less than 200 genes	334055
QIAseq Targeted cfDNA Ultra Custom (96)	All reagents (except indexes) for targeted DNA sequencing; custom panel for 96 samples	334075
QIAseq Targeted cfDNA Ultra Booster (96)	Pool of primers used in combination with either cataloged or custom panels	334085

Product	Contents	Cat. no.
QIAseq Targeted cfDNA Ultra Unique Dual Indices		
QIAseq Targeted cfDNA Ultra UDI Set A (96)	Box containing unique molecularly indexed adapters and primers, enough for a total of 96 samples, for indexing up to 96 samples for targeted panel sequencing on Illumina platforms; Set A (of A, B, C, and D) required for multiplexing 384 samples in one run	334165
QIAseq Targeted cfDNA Ultra UDI Set B (96)	Box containing unique molecularly indexed adapters and primers, enough for a total of 96 samples, for indexing up to 96 samples for targeted panel sequencing on Illumina platforms; Set B (of A, B, C, and D) required for multiplexing 384 samples in one run	334175
QIAseq Targeted cfDNA Ultra UDI Set C (96)	Box containing unique molecularly indexed adapters and primers, enough for a total of 96 samples, for indexing up to 96 samples for targeted panel sequencing on Illumina platforms; Set C (of A, B, C, and D) required for multiplexing 384 samples in one run	334185
QIAseq Targeted cfDNA Ultra UDI Set D (96)	Box containing unique molecularly indexed adapters and primers, enough for a total of 96 samples, for indexing up to 96 samples for targeted panel sequencing on Illumina platforms; Set D (of A, B, C, and D) required for multiplexing 384 samples in one run	334195
QIAseq Targeted cfDNA Ultra UDI (12)	Box containing unique molecularly-indexed adapters and primers, enough for a total of 12 samples, for indexing up to 12 samples for targeted panel sequencing on Illumina platforms	334151

Product	Contents	Cat. no.
Related products		
QIAseq Library Quant Assay Kit	Reagents for quantification of libraries prepared for Illumina or Ion Torrent platforms; assay format	333314
QIAamp Circulating Nucleic Acid Kit (50)	For 50 preps: QIAamp Mini Columns, Tube Extenders (20 mL), QIAGEN Proteinase K, Carrier RNA, Buffers, VacConnectors, and Collection Tubes	55114
QIAamp MinElute ccfDNA Mini Kit (50)	For 50 DNA preps: QIAamp UCP MinElute Columns, QIAGEN Proteinase K, Magnetic Bead Suspension, Buffers, Bead Elution Tubes, Collection Tubes	55204
QIAamp MinElute ccfDNA Midi Kit (50)	For 50 preps: QIAamp UCP MinElute Columns, QIAGEN Proteinase K, Magnetic Bead Suspension, Buffers, Bead Elution Tubes, and Collection Tubes	55284
EZ1&2 ccfDNA Kit (48)	For 48 preps: 48 reagent cartridges (EZ1&2 ccfDNA), Magnetic Bead Suspension, Elution Buffer, Large-Volume Tubes, Disposable Tip Holders, Disposable Filter-Tips, Elution Tubes	984854

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Document Revision History

Date	Changes
12/2022	Initial revision
04/2023	With reduced HiFi polymerase during TEPCR for panel with 1100 primers or less; Added two panels, Actionable Solid Tumor and Myeloid Neoplasms; Included NextSeq 1000/2000 in sequencing set up; Corrected UDI index plate name in index kit component.

This page intentionally left blank

This page intentionally left blank

Limited License Agreement for QIAseq Targeted DNA Pro

Use of this product signifies the agreement of any purchaser or user of the product to the following terms:

1. The product may be used solely in accordance with the protocols provided with the product and this handbook and for use with components contained in the kit only. QIAGEN grants no license under any of its intellectual property to use or incorporate the enclosed components of this kit with any components not included within this kit except as described in the protocols provided with the product, this handbook, and additional protocols available at www.qiagen.com. Some of these additional protocols have been provided by QIAGEN users for QIAGEN users. These protocols have not been thoroughly tested or optimized by QIAGEN. QIAGEN neither guarantees them nor warrants that they do not infringe the rights of third-parties.
2. Other than expressly stated licenses, QIAGEN makes no warranty that this kit and/or its use(s) do not infringe the rights of third-parties.
3. This kit and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. QIAGEN specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the kit agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. QIAGEN may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the kit and/or its components.

For updated license terms, see www.qiagen.com.

Trademarks: QIAGEN[®], Sample to Insight[®], QIAamp[®], QIAseq[®], QIAxcel[®], QIAxperi[®], QCI[®] (QIAGEN Group); Agilent[®], Bioanalyzer[®], TapeStation[®] (Agilent Technologies); DynaMag[™] (Thermo Fisher Scientific Inc.); LoBind[®] (Eppendorf AG); Illumina[®], HiSeq[®], MiniSeq[®], MiSeq[®], NextSeq[®], NovaSeq[™] (Illumina, Inc.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

04/2023 HB-3094-002 © 2023 QIAGEN, all rights reserved.

