

# Automation of the QIAGEN QIAseq Targeted DNA Pro Panels on the Hamilton NGS STAR MOA Generates High-Quality Libraries for Targeted Next-Generation Sequencing

**Authors:** Guido Brandt<sup>1</sup>, Andreas Finke<sup>1</sup>, Thomas Zoll<sup>1</sup>, Ronja Tauschek<sup>1</sup>, Cornelia Mechlen<sup>2</sup>, Berthold Lechtenberg<sup>2</sup>, Nadine Thaler<sup>1</sup>, Birgit Ottenwalder<sup>1</sup>, Peter Hahn<sup>2</sup>, Alexander Vial<sup>2</sup>, Christian Vogler<sup>1</sup>, Dominik Laubscher<sup>1</sup>, Kristina Ueffing<sup>1</sup>

<sup>1</sup> Hamilton Bonaduz AG, Via Crusch 8, 7402 Bonaduz, Switzerland

<sup>2</sup> QIAGEN GmbH, Qiagen Strasse 1, 40724 Hilden, Germany

## Introduction

Library preparation is essential for Next-Generation Sequencing (NGS) and one of the costliest steps in the workflow. It is time-consuming and prone to sample loss or reduced DNA quality due to handling errors. Furthermore, efficient genetic testing using targeted NGS applications demand precise detection of genetic variations, such as somatic mutations, single nucleotide polymorphisms, copy number variation, and small insertions/deletions, in defined genomic regions. The QIAseq Targeted DNA Pro Panels provide a powerful and flexible solution for targeted DNA sequencing. With their optimized workflow, the QIAseq Targeted DNA Pro Panels ensure smooth automation of library preparation with high-sample throughput on the Hamilton NGS STAR MOA (Fig. 1).

- Increased throughput – process up to 96 samples in parallel
- Enhanced reproducibility – minimized manual variability and human error
- Time efficiency – reduced hands-on time and more walk-away time
- Streamlined workflow – reduced sample loss and fewer process errors
- Automated barcode verification ensures error-free setup and traceability

## System Description


The NGS STAR MOA (Fig. 2) is designed for high-quality NGS library preparation in high-throughput workflows. The On-Deck Thermal Cycler (ODTC), one ANSI/SLAS cooling position (CPAC), and the Heater Shaker modules



Fig. 1: Hamilton NGS STAR MOA.

(HHSs) ensure optimal sample handling. A magnetic stand and tip carriers, together with carriers for samples and reagents, complete the ideal deck of the NGS STAR for NGS library preparation.

The NGS STAR MOA enables the fully-automated processing of up to 96 samples, depending on the kit used. This reduces the amount of manual work to a minimum. The correct placement of samples, reagents, plates, and tips is guaranteed by using automated barcode verification. In addition, the user can define a worklist with the combination of indexes and samples. Automated error handling and the easy-to-use software framework ensure a smooth setup of the workflow, which can also be started and stopped at specific steps within the process.



- 1 On-Deck Thermal Cycler (ODTC)
- 2 Magnetic Stand
- 3 Heater Shaker Modules (HHSs)
- 4 High Performance Cooling Modules (CPAC)

Fig. 2: Deck Layout of the NGS STAR MOA.

### Method Description

The QIAseq Targeted DNA Pro Panels - MOA v1.0 (further referred to as QIAseq Targeted DNA Pro Panels) method automates the QIAseq Targeted DNA Pro protocol (document 12/2022 HB-2979-002) on the NGS STAR MOA. The method covers the library preparation and the enrichment processes. It facilitates targeted

NGS library preparation for sequencing on Illumina sequencers and allows for the processing of a maximum of 96 samples per run, with an input amount of 10 - 80 ng gDNA or cfDNA or 100 - 250 ng FFPE DNA (Fig. 3).




Fig. 3: Graphical Overview of the QIAseq Targeted DNA Pro Panel Method Workflow.

### Enhancement Option

Add the Hamilton FLUOREYE (Fig. 4) for even more walk-away time.



Fig. 4: The Hamilton FLUOREYE

### Biological Qualification Results

To assess the biological performance of the QIAseq Targeted DNA Pro Panels method on the NGS STAR MOA, library preparation of 96 human genomic DNA samples was executed using a custom panel with 4110 primers. As input DNA, five distinct genotypes with two different input amounts each (10 ng and 40 ng) were used (see Method Requirements at the end of the Application Note). The run was conducted using 8 PCR cycles for Target Enrichment and 25 PCR cycles for Universal PCR. DNA concentration and total yield of the libraries obtained from the 96-sample biological verification run were determined using the Thermo Fisher Scientific Qubit 4 Fluorometer with the Quant-iT 1x dsDNA HS Assay Kit. Size distribution of library DNA was measured with the Agilent TapeStation 4150 using the High-Sensitivity D5000 ScreenTape with High-Sensitivity D5000 Reagents (Table 1).

Out of the generated libraries from the 96-sample biological verification run, three libraries per DNA type and input amount were selected for sequencing. They were sequenced using the NovaSeq6000 S1 Reagent Kit v1.5 on an Illumina NovaSeq 6000 at QIAGEN GmbH


(Hilden, Germany). Sequencing data was analyzed by a QIAGEN specialist using CLC Genomics Workbench 24.0 and smCounter2. On average, each sample generated over 59 million ( $\pm 7$  million) reads, with 95% of clusters passing the filter.

**Tab. 1:** Summary of Library Size and Yields and Sequencing Results.


| DNA Type                                | HD734          |                | HD752          |                | HD827          |                | HD829          |                | NA12878        |                |
|---|----------------|----------------|----------------|----------------|----------------|----------------|----------------|----------------|----------------|----------------|
| Input DNA Amount [ng]                   | 10             | 40             | 10             | 40             | 10             | 40             | 10             | 40             | 10             | 40             |
| Insert Size [bp]                        | 200 - 1000     |                |                |                |                |                |                |                |                |                |
| Sample Number                           | 6              | 6              | 6              | 6              | 6              | 6              | 6              | 6              | 24             | 24             |
| Library Yield [ng $\pm$ SD]             | 106 $\pm$ 18   | 174 $\pm$ 52   | 82 $\pm$ 41    | 189 $\pm$ 46   | 107 $\pm$ 49   | 226 $\pm$ 50   | 76 $\pm$ 7     | 139 $\pm$ 78   | 87 $\pm$ 27    | 250 $\pm$ 148  |
| Library Size [bp $\pm$ SD]              | 472 $\pm$ 16   | 466 $\pm$ 20   | 467 $\pm$ 30   | 469 $\pm$ 18   | 477 $\pm$ 39   | 464 $\pm$ 14   | 465 $\pm$ 27   | 449 $\pm$ 30   | 479 $\pm$ 23   | 479 $\pm$ 37   |
| Total Number of Reads [Mio $\pm$ SD]    | 57 $\pm$ 4     | 59 $\pm$ 8     | 60 $\pm$ 4     | 62 $\pm$ 10    | 51 $\pm$ 6     | 66 $\pm$ 9     | 49 $\pm$ 9     | 63 $\pm$ 7     | 51 $\pm$ 16    | 71 $\pm$ 9     |
| Mean Read Fragments per UMI [ $\pm$ SD] | 8.6 $\pm$ 0.5  | 5.8 $\pm$ 2.3  | 11.1 $\pm$ 2.4 | 6.3 $\pm$ 1.3  | 9.1 $\pm$ 0.6  | 4.8 $\pm$ 0.9  | 13.3 $\pm$ 1.4 | 6.2 $\pm$ 1.1  | 12.3 $\pm$ 4.3 | 4.5 $\pm$ 0.6  |
| Specificity [% $\pm$ SD]                | 79.8 $\pm$ 0.6 | 81.3 $\pm$ 3.7 | 77.7 $\pm$ 0.7 | 81.0 $\pm$ 1.9 | 79.9 $\pm$ 1.7 | 82.6 $\pm$ 1.5 | 77.0 $\pm$ 2.2 | 81.3 $\pm$ 1.0 | 77.5 $\pm$ 3.2 | 84.0 $\pm$ 0.4 |

The accuracy of Variant Allelic Frequency (VAF) estimation for HD827, HD829, and HD734 gDNA libraries depicts a strong match between expected and detected frequencies (Fig. 5 and Fig. 6). Libraries made from

10 ng and 40 ng input DNA performed consistently well. This accuracy highlights the method's ability to detect important variants, even with low DNA input, proving the workflow's reliability in various experiments.



**Fig. 5:** Accuracy of Variant Allelic Frequency (VAF) estimation for (A) 10 ng and (B) 40 ng HD827 gDNA NGS libraries.



**Fig. 6:** Accuracy of Variant Allelic Frequency (VAF) estimation of selected mutations of the for (A) HD829 and (B) HD734 gDNA NGS libraries.

In conclusion, automating the QIAseq Targeted DNA Pro Panels method on the Hamilton NGS STAR MOA offers an efficient and reliable solution for high-throughput NGS library preparation. It reduces hands-on time

and errors, ensuring consistent library quality for various sequencing applications. This system's robust performance across different samples and library sizes optimizes NGS workflows for all lab scales.

All Hamilton Consumables can be ordered from our Web Shop.



HAMILTON'S CONSUMABLES  
**WEBSHOP**  
Web: [eshop.hamilton.com](http://eshop.hamilton.com)

| Requirements   |                          |             |
|--|--------------------------|-------------|
| System Requirements  | Provider                 | Part Number |
| NGS STAR MOA / NGS STAR MOA UV incl. Deck Components       | Hamilton Bonaduz AG      | 806800      |
| ODTC 96  | Hamilton Bonaduz AG      | 10147734    |
| iSWAP  | Hamilton Bonaduz AG      | 190220      |
| Labware Requirements                                       |                          |             |
| 50 µL CO-RE Filter Tips                                    | Hamilton Bonaduz AG      | 235948      |
| 300 µL CO-RE Filter Tips                                   | Hamilton Bonaduz AG      | 235903      |
| 1000 µL CO-RE Filter Tips                                  | Hamilton Bonaduz AG      | 235905      |
| PCR ComfortLid   | Hamilton Bonaduz AG      | 814300      |
| PCR FramePlate 96-well                                     | Hamilton Bonaduz AG      | 814302      |
| 20 mL Reagent Container, natural color, no lid             | Hamilton Bonaduz AG      | 10161052    |
| 60 mL PP Reagent Trough with Lid                           | Hamilton Bonaduz AG      | 56694-01    |
| 0.5 mL Screw Cap Micro Tubes                               | Sarstedt                 | 72.730.006  |
| 2 mL Screw Cap Micro Tubes                                 | Sarstedt                 | 72.694.006  |
| Abgene 96-Well 0.8 mL Polypropylene DeepWell Storage Plate | Thermo Fisher Scientific | AB0859      |
| Method Requirements  |                          |             |
| CEPH1463   | Coriell Institute        | NA12878     |
| Myeloid DNA Reference Standard                             | Horizon Discovery        | HD829       |
| OncoSpan gDNA  | Horizon Discovery        | HD827       |
| Tru-Q 7 (1.3% Tier) Reference Standard                     | Horizon Discovery        | HD734       |
| Tru-Q 0 (100% Wildtype) Reference Standard                 | Horizon Discovery        | HD752       |
| QIAseq Targeted DNA Pro (96)                               | QIAGEN                   | 333655      |
| QIAseq Targeted DNA Pro HC (96)                            | QIAGEN                   | 333665      |
| QIAseq Targeted DNA Pro UDI Set A                          | QIAGEN                   | 333455      |
| QIAseq Targeted DNA Pro UDI Set B                          | QIAGEN                   | 333465      |
| QIAseq Targeted DNA Pro UDI Set C                          | QIAGEN                   | 333475      |
| QIAseq Targeted DNA Pro UDI Set D                          | QIAGEN                   | 333485      |
| High Sensitivity D5000 Reagents                            | Agilent                  | 5067-5593   |
| High Sensitivity D5000 ScreenTape                          | Agilent                  | 5067-5592   |
| NovaSeq6000 S1 Reagent Kit v1.5                            | Illumina                 | 20028317    |
| QuantIT 1X dsDNA HS Assay Kit                              | Thermo Fisher Scientific | Q33232      |

Throughout this document, protected product names may be used without being specifically marked as such. Research use only. Not for use in diagnostics procedures. All rights reserved. All other trademarks are the sole property of their respective owners. © 2025 Hamilton Company. All rights reserved. All trademarks are owned and/or registered by Hamilton Company in the U.S. and/or other countries. Lit. No. AN-2502-01-02/2025



To find a representative in your area, please visit:

[hamiltoncompany.com/contact](http://hamiltoncompany.com/contact)  
[infoservice@hamiltonrobotics.com](mailto:infoservice@hamiltonrobotics.com)

**United States**  
+1-775-858-3000  
**United Kingdom, Ireland**  
+44 121 272 92 80  
**Brazil**  
+55 11 95914 5000  
**China**  
+86 21 6164 6567

**Denmark, Norway, Sweden, Finland**  
+46 8410 27 373  
**France**  
+33 184 008 420  
**Germany, Austria, Switzerland**  
+49 89 248 804 804

**Benelux**  
+31 40 209 178 0  
**Italy**  
+39 039 930 06 06  
**Japan**  
+81 3 6435 6850  
**Spain, Portugal**  
+34 930 186 262

