


# QIAcube<sup>®</sup> Connect MDx User Manual

For use with software version 2.x



**IVD** For in vitro diagnostic use



**REF** 9003070



QIAGEN, GmbH, QIAGEN Strasse 1, 40724 Hilden, GERMANY

**MAT**<sub>R3</sub>

# Table of Contents

1. Introduction .....	4
1.1. About this user manual .....	4
1.2. General information .....	5
1.3. Intended use of the QIAcube Connect MDx .....	6
2. Safety Information .....	7
2.1. Proper use .....	7
2.2. Electrical safety .....	9
2.3. Environment .....	10
2.4. Biological safety .....	10
2.5. Chemical safety .....	11
2.6. Waste disposal .....	12
2.7. Mechanical hazards .....	12
2.8. Heat hazard .....	13
2.9. Maintenance safety .....	13
2.10. Radiation safety .....	14
2.11. Symbols on the QIAcube Connect MDx .....	15
3. General Description .....	17
3.1. QIAcube Connect MDx principle .....	17
3.2. External features of the QIAcube Connect MDx .....	19
3.3. Internal features of the QIAcube Connect MDx .....	24
3.4. Disposables .....	29
4. Installation Procedures .....	30
4.1. Installation environment .....	30
4.2. Unpacking the QIAcube Connect MDx .....	32
4.3. Installing the QIAcube Connect MDx .....	33
4.4. Repackaging and shipping of the QIAcube Connect MDx .....	42
4.5. Configuration of the QIAcube Connect MDx .....	45
5. Operating Procedures .....	57
5.1. Use of the QIAcube Connect MDx software .....	59
5.2. Switching the QIAcube Connect MDx on and off .....	63
5.3. Logging in and logging out .....	63
5.4. Setting up a protocol run .....	64
5.5. Starting a protocol run .....	82
5.6. Stopping a protocol run .....	86
5.7. Saving run reports to the USB flash drive .....	87
5.8. Independent heater/shaker operation .....	90
5.9. Independent centrifuge operation .....	91
5.10. Managing protocols .....	94
5.11. User Management .....	104
6. Cleaning and Maintenance .....	110
6.1. Cleaning agents .....	111
6.2. Decontaminating the QIAcube Connect MDx surface .....	111
6.3. Regular maintenance .....	113
6.4. Daily maintenance .....	114
6.5. Monthly maintenance .....	115
6.6. Periodic maintenance .....	116
6.7. Optional maintenance .....	124
6.8. Decontaminating the QIAcube Connect MDx .....	126
6.9. QIAcube Connect MDx repair .....	127

7. Troubleshooting .....	128
7.1. Creating a support package .....	129
7.2. Operation .....	129
8. Glossary .....	141
9. Technical Specifications .....	142
9.1. Operating conditions .....	142
9.2. Transport conditions .....	142
9.3. Storage conditions .....	142
9.4. Mechanical data and hardware features .....	143
Appendix A – Legal .....	144
Declaration of Conformity .....	144
Waste Electrical and Electronic Equipment (WEEE) .....	144
EMC Declaration .....	146
California Proposition 65 .....	147
Liability Clause .....	147
Appendix B – QIAcube Connect MDx Accessories .....	148
Ordering information .....	148
Document Revision History .....	150

# 1. Introduction

Thank you for choosing the QIAcube Connect MDx. We are confident it will become an integral part of your laboratory. Before using the QIAcube Connect MDx, it is essential that you read this user manual carefully and pay attention to the safety information. The instructions and safety information in the user manual must be followed to ensure safe operation of the instrument and to maintain the instrument in a safe condition.

## 1.1. About this user manual

This user manual provides information about the QIAcube Connect MDx in the following sections:

- Introduction
- Safety Information
- General Description
- Installation Procedures
- Operating Procedures
- Cleaning and Maintenance
- Troubleshooting
- Glossary
- Document Revision History

The appendices contain the following information:

- Technical Specifications
- Appendix A – Legal
- Appendix B – QIAcube Connect MDx Accessories



## 1.2. General information

### 1.2.1. Technical assistance

At QIAGEN®, we pride ourselves on the quality and availability of our technical support. Our Technical Services Departments are staffed by experienced scientists with extensive practical and theoretical expertise in molecular biology and the use of QIAGEN products. If you have any questions or experience any difficulties regarding the QIAcube Connect MDx or QIAGEN products in general, do not hesitate to contact us.

QIAGEN customers are a major source of information regarding advanced or specialized uses of our products. This information is helpful to other scientists as well as to the researchers at QIAGEN. We therefore encourage you to contact us if you have any suggestions about product performance or new applications and techniques.

For technical assistance, contact QIAGEN Technical Services.

Website: [support.qiagen.com](https://support.qiagen.com)

When contacting QIAGEN Technical Services about errors, please have the following information ready:

- QIAcube Connect MDx serial number, type, and version
- Error code (if applicable)
- Timepoint when the error occurred for the first time
- Frequency of error occurrence (i.e., intermittent or persistent error)
- Copy of log files (support package generated according to Section 7.1)

For up-to-date information about QIAcube Connect MDx instrument, visit [www.qiagen.com/QIAcube-Connect-MDx](https://www.qiagen.com/QIAcube-Connect-MDx).

For technical assistance and more information, please see our Technical Support Center at [www.qiagen.com/support/technical-support](https://www.qiagen.com/support/technical-support) or call one of the QIAGEN Technical Service Departments or local distributors (visit [www.qiagen.com](https://www.qiagen.com)).

### 1.2.2. Policy statement

It is the policy of QIAGEN to improve products as new techniques and components become available. QIAGEN reserves the right to change specifications at any time.

To produce useful and appropriate documentation, we appreciate your comments on this user manual. Please contact QIAGEN Technical Services.

## 1.3. Intended use of the QIAcube Connect MDx

The QIAcube Connect MDx instrument is designed to perform fully automated isolation and purification of nucleic acids in molecular diagnostic and/or molecular biology applications. The system is intended for use by professional users, such as technicians and physicians trained in molecular biological techniques and the operation of the instrument.

The QIAcube Connect MDx instrument (in the IVD mode) is intended to be used only in combination with QIAGEN and PreAnalytiX kits indicated for use with the QIAcube Connect MDx instrument for the applications that will be described in the kit handbooks.

### 1.3.1. Limitations of use

Only use the instrument in combination with the accessories specified in Appendix B – QIAcube Connect MDx Accessories. Other limitations of the applications are specified in the respective kit handbooks.

### 1.3.2. Requirements for QIAcube Connect MDx users

The table below covers the general level of competence and training necessary for transportation, installation, use, maintenance, and servicing of the QIAcube Connect MDx.

Task	Personnel	Level of competence and expertise
Delivery	No special requirements	No special requirements
Installation, routine use, and maintenance	Laboratory technicians or equivalent	Appropriately trained and experienced personnel familiar with use of computers and automation in general
Servicing and required annual maintenance	QIAGEN Field Service specialists or service technicians of an authorized agent	Trained and authorized by QIAGEN

## 2. Safety Information

Before using the QIAcube Connect MDx, it is essential that you read this user manual carefully and pay attention to the safety information. The instructions and safety information in the user manual must be followed to ensure safe operation of the instrument and to maintain the instrument in a safe condition.

Possible hazards that could harm the user or result in damage to the instrument are clearly stated at the appropriate places throughout this consolidated operating guide.

If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

The following types of safety information appear in this manual.

### **WARNING**



The term **WARNING** is used to inform you about situations that could result in personal injury to you or others.

Details about these circumstances are given in a box like this one.

### **CAUTION**



The term **CAUTION** is used to inform you about situations that could result in **damage to an instrument** or other equipment.

Details about these circumstances are given in a box like this one.

The advice given in this user manual is intended to supplement, not supersede, the normal safety requirements prevailing in the user's country.

Please be aware that you may be required to consult your local regulations for reporting serious incidents that have occurred in relation to the device to the manufacturer and the regulatory authority in which the user and/or the patient is established.

### 2.1. Proper use

#### **WARNING/ CAUTION**



#### **Risk of personal injury and material damage**

Improper use of the QIAcube Connect MDx may cause personal injuries or damage to the instrument. The QIAcube Connect MDx must only be operated by qualified personnel who have been appropriately trained. Servicing of the QIAcube Connect MDx must only be performed by a QIAGEN field service specialist.

Perform the maintenance as described in Section 6, Cleaning and Maintenance. QIAGEN charges for repairs that are required due to incorrect maintenance.

### **WARNING**



#### **Risk of personal injury and material damage**

The QIAcube Connect MDx is too heavy to be lifted by one person. To avoid personal injury or damage to the instrument, do not lift the instrument alone.

Contact QIAGEN Technical Services to relocate the instrument.

### **WARNING**



#### **Risk of personal injury and material damage**

Do not attempt to move the QIAcube Connect MDx during operation.

**CAUTION**    **Damage to the instrument**



Avoid spilling water or chemicals onto the QIAcube Connect MDx. Instrument damage caused by water or chemical spillage will void your warranty.

In case of emergency, power OFF the QIAcube Connect MDx at the power switch located in front of the instrument and unplug the power cord from the power outlet.

**CAUTION**    **Damage to the instrument**



Only use QIAGEN spin columns and the QIAcube Connect MDx-specific consumables with the QIAcube Connect MDx. Instrument damage caused by use of other types of spin columns or chemistries will void your warranty. during operation.

**WARNING**    **Risk of personal injury and material damage**



Do not use damaged rotor adapters. The rotor adapters can only be used once. High *g* forces exerted in the centrifuge can cause damage to re-used rotor adapters.

**CAUTION**    **Damage to the instrument**



Empty the tip disposal container prior to use to prevent a tip jam in the waste drawer. Failure to empty the waste container may block the robotic arm that could cause run failure or instrument damage.

**WARNING**    **Risk of personal injury and material damage**



To prevent plastic crash, load the tubes properly. After a plastic crash, sharp plastic particles could be inside the centrifuge. Be careful when handling items inside the centrifuge.

**CAUTION**    **Damage to the instrument**



Only use the correct volume of liquids.  
Exceeding the recommended volume of liquids may damage the centrifuge rotor or instrument.

**WARNING**    **Risk of fire or explosion**



When using ethanol or ethanol-based liquids on the QIAcube Connect MDx, handle such liquids carefully and in accordance with the required safety regulations. If liquid has been spilled, wipe it off and leave the QIAcube Connect MDx hood open to allow flammable vapors to disperse.

**WARNING**    **Risk of explosion**



The QIAcube Connect MDx is intended for use with reagents and substances supplied with QIAGEN kits as outlined in respective information for use. Use of other reagents and substances may lead to fire or explosion.

If hazardous material is spilled on or inside the QIAcube Connect MDx, the user is responsible for carrying out appropriate decontamination.

**Note:** Do not place items on top of the QIAcube Connect MDx hoods.

**CAUTION**    **Damage to the instrument**



Do not lean against the touchscreen when it is pulled out.

## 2.2. Electrical safety

**Note:** Disconnect the line power cord from the power outlet before servicing.

### **WARNING** Electrical hazard



Any interruption of the protective conductor (earth/ground lead) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument dangerous.

Intentional interruption is prohibited.

#### **Lethal voltages inside the instrument**

When the instrument is connected to line power, terminals may be live and opening covers or removing parts is likely to expose live parts.

### **WARNING** Damage to electronics



Before powering ON the instrument, make sure that the correct supply voltage is used.

Use of incorrect supply voltage may damage the electronics.

To check the recommended supply voltage, refer to the specifications indicated in the type plate of the instrument.

### **WARNING** Risk of electric shock



Do not open any panels on the QIAcube Connect MDx.

#### **Risk of personal injury and material damage**

Only perform maintenance that is specifically described in this user manual. Any other maintenance or repair may only be carried out by an authorized Field Service Specialist.

To ensure satisfactory and safe operation of the QIAcube Connect MDx, follow these guidelines:

- The line power cord must be connected to a line power outlet that has a protective conductor (earth/ground).
- Place instrument in a location so that the power cord is accessible and can be connected/disconnected.
- Use only the power cord delivered by QIAGEN.
- Do not adjust or replace internal parts of the instrument.
- Do not operate the instrument with any covers or parts removed.
- If liquid has spilled inside the instrument, power OFF the instrument, disconnect it from the power outlet, and contact QIAGEN Technical Services.

If the instrument becomes electrically unsafe, prevent other personnel from operating it, and contact QIAGEN Technical Services.

The instrument may be electrically unsafe when:

- The instrument or the line power cord appears to be damaged.
- It has been stored under unfavorable conditions for a prolonged period.
- It has been subjected to severe transport stresses.
- Liquids come in contact directly with electrical components of the QIAcube Connect MDx.
- The power cord has been exchanged with a power cord that was not intended to be used with the QIAcube Connect MDx.

## 2.3. Environment

Parameters such as temperature range and humidity range are described in Technical Specifications.

### 2.3.1. Operating conditions

#### **WARNING** Explosive atmosphere



The QIAcube Connect MDx is not designed for use in an explosive atmosphere.

#### **CAUTION** Damage to the instrument



Direct sunlight may bleach parts of the instrument and cause damage to plastic parts. The QIAcube Connect MDx must be located out of direct sunlight.

#### **CAUTION** Damage to the instrument



Do not use the QIAcube Connect MDx in the vicinity of sources of strong electromagnetic radiation (e.g., unshielded, deliberately operated high-frequency sources or mobile radio devices), because these can interfere with the proper operation.

## 2.4. Biological safety

Specimens and reagents containing materials from humans should be treated as potentially infectious. Use safe laboratory procedures as outlined in publications such as Biosafety in Microbiological and Biomedical Laboratories, HHS ([www.cdc.gov/biosafety](http://www.cdc.gov/biosafety)).

### 2.4.1. Samples

Samples may contain infectious agents. You should be aware of the health hazard presented by such agents and should use, store, and dispose of such samples according to the required safety regulations.

### **WARNING** Samples containing infectious agents



Samples used with the QIAcube Connect MDx may contain infectious agents. Handle such samples with the greatest of care and in accordance with the required safety regulations.

Always wear safety glasses, gloves, and a lab coat.

The responsible body (for example, a laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe, and that the instrument operators are suitably trained and not exposed to hazardous levels of infectious agents as defined in the applicable Material Safety Data Sheets (MSDSs) or the OSHA<sup>1</sup>, \* ACGIH,<sup>†</sup> or COSHH<sup>‡</sup> documents.

Venting for fumes and disposal of waste must be in accordance with all national, state, and local health and safety regulations and laws.

\* OSHA — Occupational Safety and Health Organization (United States of America)

† ACGIH — American Conference of Government Industrial Hygienists (United States of America)

‡ COSHH — Control of Substances Hazardous to Health (United Kingdom)

## 2.5. Chemical safety

### **WARNING** Hazardous chemicals



Some chemicals used with the QIAcube Connect MDx may be hazardous or may become hazardous after completion of purification.

Always wear safety glasses, gloves, and a lab coat.

The responsible body (for example, a laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe, and that the instrument operators are suitably trained and not exposed to hazardous levels of infectious agents as defined in the applicable Material Safety Data Sheets (MSDSs) or the OSHA<sup>1</sup>, \* ACGIH,<sup>†</sup> or COSHH<sup>‡</sup> documents.

Venting for fumes and disposal of waste must be in accordance with all national, state, and local health and safety regulations and laws.

\* OSHA — Occupational Safety and Health Organization (United States of America)

† ACGIH — American Conference of Government Industrial Hygienists (United States of America)

‡ COSHH — Control of Substances Hazardous to Health (United Kingdom)

### 2.5.1. Toxic fumes

If working with volatile solvents or toxic substances, you must provide an efficient laboratory ventilation system to remove vapors that may be produced.

### **WARNING** Toxic fumes



Do not use bleach to clean or disinfect the QIAcube Connect MDx or the labware, as bleach in contact with salts from the buffers can produce toxic fumes.

### **WARNING** Toxic fumes



Do not use bleach to disinfect used labware. Bleach in contact with salts from the buffers used can produce toxic fumes.

## 2.6. Waste disposal

Used labware, such as sample tubes, QIAGEN spin columns, filter-tips, buffer bottle and enzyme tubes, or rotor adapters, may contain hazardous chemicals or infectious agents from the purification process. These hazardous wastes must be collected and disposed of properly according to local safety regulations.

For more information about how to dispose of the QIAcube Connect MDx, see Appendix A – Legal, Waste Electrical and Electronic Equipment (WEEE)).

### **WARNING** Hazardous chemicals and infectious agents



Waste may contain toxic or infectious material and must be disposed of properly. Refer to your local safety regulations for proper disposal procedures.

## 2.7. Mechanical hazards

The hood of the QIAcube Connect MDx must remain closed during operation of the instrument. Only open the hood when instructed to do so by the instruction for use.

Do not lean on the worktable when the robotic arm of the instrument is moving to reach loading position with its lid open. Wait until the robotic arm has stopped moving before loading or unloading the instrument workdeck.

### **WARNING** Moving parts



Avoid contact with moving parts during operation of the QIAcube Connect MDx. Do not place your hands under the robotic arm when it is lowering. Do not attempt to move any tip racks or tubes while the instrument is operating.

### **WARNING** Moving parts



To avoid contact with moving parts during the operation of the QIAcube Connect MDx, the instrument must be operated with the hood closed.

If the hood sensor or lock is not functioning properly, contact QIAGEN Technical Services.

### 2.7.1. Centrifuge

Make sure that the rotor and buckets are installed correctly. All buckets must be mounted before starting a protocol run, regardless of the number of samples to be processed. If the rotor or buckets show signs of mechanical damage or corrosion, or the rotor positioning PIN is loose or damaged, do not use the QIAcube Connect MDx; contact QIAGEN Technical Services.



**CAUTION**    **Damage to the instrument**



The QIAcube Connect MDx must not be used if the centrifuge lid is broken, or if the lid lock is damaged. Make sure that no loose material is inside the centrifuge during operation.

Make sure that the rotor is installed correctly and that all buckets are properly mounted, regardless of the number of samples to be processed. Load the rotor only as instructed by the software.

Only use rotors, buckets, and consumables designed for use with the QIAcube Connect MDx. Damage caused by use of other consumables will void your warranty.

We recommend replacing the centrifuge rotor and buckets after 20,000 cycles, which is equivalent to 9 years of usage with two runs per day for 220 days each year. For more information contact QIAGEN Technical Services.

In case of breakdown caused by power failure, the centrifuge lid can be opened manually to remove the samples (see Section 7.2.2).

**WARNING**    **Moving parts**



In case of breakdown caused by power failure, remove the power cord and wait 10 minutes before attempting to manually open the centrifuge lid.

**CAUTION**    **Damage to the instrument**



After a power failure, do not move the z-module (robotic arm) manually in front of the instrument. Damage may occur if the QIAcube Connect MDx hood is closed and collides with the z-module.

**WARNING**    **Risk of personal injury and material damage**



Raise the centrifuge lid carefully. The lid is heavy and may cause injury if it falls.

**WARNING**    **Risk of overheating**



To ensure proper ventilation, maintain a minimum clearance of 10 cm at the sides and rear of the QIAcube Connect MDx.

Slits and openings that ensure the ventilation of the instrument must not be covered.

## 2.8. Heat hazard

The QIAcube Connect MDx worktable contains a heated shaker.

**WARNING**    **Hot surface**



The shaker can reach temperatures of up to 70°C (158°F). Avoid touching it when it is hot, in particular, shortly after a run has been carried out.

## 2.9. Maintenance safety

**WARNING/  
CAUTION**    **Risk of personal injury and material damage**



Only perform maintenance that is specifically described in this user manual.

**WARNING** Risk of explosion



When cleaning the QIAcube Connect MDx with alcohol-based disinfectant, leave the QIAcube Connect MDx hood open to allow flammable vapors to disperse.

Only clean the QIAcube Connect MDx when worktable components have cooled down.

**WARNING** Risk of fire



Do not allow cleaning fluid or decontamination agents to come into contact with the electrical parts of the QIAcube Connect MDx.

**WARNING** Risk of personal injury and material damage



To prevent the rotor nuts from loosening during operation of the centrifuge, securely tightened the nuts using the rotor key supplied with the QIAcube Connect MDx.

**WARNING** UV irradiation hazard



A mechanical lock ensures that the hood must be closed for operation of the UV LED.

If the hood sensor or lock is not functioning correctly, contact QIAGEN Technical Services.

**WARNING** Risk of personal injury and material damage



Make sure that lids from spin columns and 1.5 mL microcentrifuge tubes are in the correct position and pushed all the way down to the bottom of the slots on the sides of the rotor adapter. Incorrectly positioned lids can break off during centrifugation.

**WARNING** Risk of personal injury and material damage



Be sure the lid is completely removed from the spin column. Spin columns with partially removed lids may not be removed properly from the rotor, causing the protocol run to crash.

**CAUTION** Damage to the instrument



Do not use bleach, solvents, or reagents containing acids, alkalis, or abrasives to clean the QIAcube Connect MDx.

**CAUTION** Damage to the instrument



Do not use spray bottles containing alcohol or disinfectant to clean surfaces of the QIAcube Connect MDx. Spray bottles should be used only to clean items that have been removed from the worktable and if permitted by local laboratory operating practices.

## 2.10. Radiation safety

**WARNING** Risk of personal injury



Do not expose your skin to UV-C light from the UV LED lamp.









**WARNING** Risk of personal injury



Hazard Level 2 laser light: Do not stare into the light beam when using handheld barcode scanner.

## 2.11. Symbols on the QIAcube Connect MDx

Symbol	Location	Description
	Next to the shaker	Heat hazard – the temperature of the shaker can reach up to 70°C (158°F).
	Near the centrifuge; near the robotic arm	Mechanical hazard — avoid contact with moving parts.
	On the instrument, near bottle rack	Fire hazard – usage of ethanol in bottle rack.
	In front of worktable	Biological hazard – some samples used with this instrument may contain infectious agents and must be handled with gloves.
	Inside the waste drawer	Biohazard — the waste drawer may be contaminated with biohazardous material and must be handled with gloves.
	Type plate on the back of the instrument	CE mark for European Conformity
	Type plate at the back of the instrument	CSA listing mark for Canada and the USA
	Type plate on the back of the instrument	FCC mark of the United States Federal Communications Commission
	Type plate on the back of the instrument	RCM mark for Australia and New Zealand
	Type plate on the back of the instrument	RoHS mark for China (the restriction of the use of certain hazardous substances in electrical and electronic equipment)
 	Type plate on the back of the instrument	Waste Electrical and Electronic Equipment (WEEE) mark for Europe
	Type plate on the back of the instrument	Legal manufacturer

Symbol	Location	Description
	At the back of the instrument	Consult instructions for use
	Type plate on the back of the instrument	See warnings and precautions
	Type plate on the back of the instrument	In vitro diagnostic medical device
	Type plate on the back of the instrument	Unique device identifier
	Type plate on the back of the instrument	Unique Device Identifier (UDI) as a 2D bar code in Data Matrix format
	Type plate on the back of the instrument	Global Trade Item Number
	Type plate on the back of the instrument	Serial number
	Type plate on the back of the instrument	Catalog number

## 3. General Description

The QIAcube Connect MDx performs fully automated isolation and purification of nucleic acids in molecular diagnostic and molecular biology applications.

It can process up to 12 samples per single run. The QIAcube Connect MDx is designed to automate selected QIAGEN DSP and non-DSP kits and the PAXgene® Blood RNA Kit. The QIAcube Connect MDx controls integrated components, including a centrifuge, a heated shaker, a pipetting system, a UV LED, and a robotic gripper.

The QIAcube Connect MDx gives the option to start a protocol either in the IVD mode of the software (for validated IVD applications only) or in the Research mode of the software (for Molecular Biology Applications (MBA) only). The usage of IVD protocols is only possible and strictly restricted to the IVD mode of the software. This user manual focuses on operating the QIAcube Connect MDx in the IVD- software mode. For detailed instructions on how to operate the QIAcube Connect MDx using the Research mode of the software (with MBA protocols or any customized protocols), please refer to the *QIAcube Connect User Manual* (available on the QIAcube Connect product webpage under the **Resources** tab: [www.qiagen.com/HB-2594](http://www.qiagen.com/HB-2594)).

The QIAcube Connect MDx is preinstalled with various protocols for processing QIAGEN spin columns for purification of RNA, genomic DNA, and viral nucleic acids. In the Research mode of the software, further protocols are available, for example, plasmid DNA and protein purification as well as DNA and RNA cleanup. The user first selects the software mode for the type of application to be performed using the touchscreen then chooses an application or scans a kit barcode and loads labware, samples, and reagents onto the QIAcube Connect MDx worktable. The user then closes the instrument hood and starts the protocol, which provides all necessary commands for sample lysis and purification using QIAGEN spin columns. A fully automated load check helps to ensure correct loading of the worktable.


Offering an augmented user interface, users stay connected to their instrument through the built-in screen and also remotely with a computer or a mobile device (e.g., a tablet) and the QIASphere® App, enabling quick response times and ability to monitor runs while being away from the instrument.

### 3.1. QIAcube Connect MDx principle

Sample preparation using QIAcube Connect MDx follows the same steps as the manual procedure (i.e., lyse, bind, wash, and elute as outlined below). Depending on the chosen application, the procedure may be alternated or steps may not be required. Moving forward to laboratory automation, no change in purification chemistry is required as you simply continue to use trusted QIAGEN spin-column kits.


1. Samples are lysed in the orbital shaker, which can be heated if required by the protocol.
2. Each lysate is transferred to a spin column in a rotor adapter. If the lysate needs to be homogenized or cleared, it is first transferred to the middle position of the rotor adapter.
3. Nucleic acids or proteins bind to the silica membrane or purification resin of the QIAGEN spin column and are washed to remove contaminants.
4. The spin column is transferred to a microcentrifuge tube for elution of purified nucleic acids or protein.

For an example workflow using a QIAGEN spin-column kit, see the flowchart on the next page.



Exemplary procedure flowchart.

### 3.2. External features of the QIAcube Connect MDx



Front view of the QIAcube Connect MDx.



Pulled-out touchscreen.



Rear view of the QIAcube Connect MDx.



Rear view of the QIAcube Connect MDx.

- |   |              |   |  |
|---|--------------|---|--|
| 1 | Touchscreen  | 5 | 2 USB ports on the left side of the touchscreen; 2 USB ports behind the touchscreen (Wi-Fi module plugged into 1 USB port) |
| 2 | Hood         | 6 | RJ-45 Ethernet port  |
| 3 | Waste drawer | 7 | Power cord socket  |
| 4 | Power switch | 8 | Cooling air outlet   |
|   |              | 9 | External barcode scanner (not depicted)  |

## Touchscreen

The QIAcube Connect MDx is controlled using a swivel-mounted touchscreen (1). The touchscreen allows the user to operate the instrument and to guide the user through run setup and worktable loading. During sample processing, the touchscreen shows the protocol status and remaining time.



Pulled-out touchscreen.

## Hood

The QIAcube Connect MDx hood (2) protects users from the moving robotic arm and from potentially infectious material placed on the worktable. The hood can be manually opened to gain access to the worktable. During operation of the QIAcube Connect MDx, the hood must remain closed and should only be opened when instructed to do so by the software. Therefore, a hood lock has been implemented to avoid inappropriate opening.

### **WARNING** Moving parts



To avoid contact with moving parts during the operation of the QIAcube Connect MDx, the instrument must be operated with the hood closed.

If the hood sensor or lock is not functioning properly, contact QIAGEN Technical Services.

## Power switch

The power switch (4) is located at the front right of the QIAcube Connect MDx and is used to power ON and OFF the instrument.

## RJ-45 Ethernet port

The RJ-45 Ethernet port (6) located at the back of the instrument beside the power cord socket, is only used to connect the QIAcube Connect MDx to local area network via cable.



## USB ports

The QIAcube Connect MDx has four USB ports (5). Two are located on the left of the touchscreen and two are located behind the touchscreen.

The USB ports located on the left of the touchscreen allow connection of the QIAcube Connect MDx to a USB flash drive. Data files, such as support package, protocols, or report files, can be transferred via the USB port from the QIAcube Connect MDx to the USB flash drive. Protocols can be uploaded via the USB flash drive. The USB ports can also be used to plug in the supplied external barcode scanner.

USB ports located below the touchscreen, allow insertion of a Wi-Fi adapter to enable Wi-Fi connection to a local network.

**Important:** Use only the USB flash drive provided by QIAGEN. Do not connect other USB flash drive devices to USB ports. Insert only one USB stick for data transfer. Otherwise no USB stick can be recognized.

**Important:** Do not remove the USB flash drive while downloading or transferring data or software to or from the instrument.

**Important:** Always turn off the QIAcube Connect MDx before plugging or unplugging the Wi-Fi USB device. Plug-and-play of the Wi-Fi USB device while the instrument is turned on is not supported. The Wi-Fi USB device may have been delivered to you (availability may differ from country to country based on regulations and approvals). If you have not received a Wi-Fi USB device from QIAGEN, ensure that the Wi-Fi adapter supports IEEE 802.11-2016 standards, including WIFI 4 (802.11n), WIFI 5 (802.11a/c), WPA/WPA2 (802.11i), and WPA3 (SAE). It is recommended to use a Wi-Fi adapter with the RTL8723BU chipset. The adapter should be in accordance with your local laws and regulations. For more information contact QIAGEN Technical Services.

## Waste drawer

Used disposable filter-tips are discarded through two slots in the worktable and collected in the waste drawer (3). Used columns (e.g., QIAshredder columns) are also discarded into this drawer.

### CAUTION Damage to the instrument



Empty the tip disposal container prior to use to prevent a tip jam in the waste drawer. Failure to empty the waste container may block the robotic arm that could cause run failure or instrument damage.

### CAUTION Hazardous chemicals and infectious agents



Waste may contain toxic or infectious material and must be disposed of properly. Refer to your local safety regulations for proper disposal procedures.

## **WARNING** Hazardous chemicals



Some chemicals used with the QIAcube Connect MDx may be hazardous or may become hazardous after completion of purification.

Always wear safety glasses, gloves, and a lab coat.

The responsible body (for example, a laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe, and that the instrument operators are suitably trained and not exposed to hazardous levels of infectious agents as defined in the applicable Material Safety Data Sheets (MSDSs) or the OSHA<sup>1</sup>, \* ACGIH,<sup>†</sup> or COSHH<sup>‡</sup> documents.

Venting for fumes and disposal of waste must be in accordance with all national, state, and local health and safety regulations and laws.

\* OSHA — Occupational Safety and Health Organization (United States of America)

† ACGIH — American Conference of Government Industrial Hygienists (United States of America)

‡ COSHH — Control of Substances Hazardous to Health (United Kingdom)

## **WARNING** Samples containing infectious agents



Samples used with the QIAcube Connect MDx may contain infectious agents. Handle such samples with the greatest of care and in accordance with the required safety regulations.

Always wear safety glasses, gloves, and a lab coat.

The responsible body (for example, a laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe, and that the instrument operators are suitably trained and not exposed to hazardous levels of infectious agents as defined in the applicable Material Safety Data Sheets (MSDSs) or the OSHA<sup>1</sup>, \* ACGIH,<sup>†</sup> or COSHH<sup>‡</sup> documents.

Venting for fumes and disposal of waste must be in accordance with all national, state, and local health and safety regulations and laws.

\* OSHA — Occupational Safety and Health Organization (United States of America)

† ACGIH — American Conference of Government Industrial Hygienists (United States of America)

‡ COSHH — Control of Substances Hazardous to Health (United Kingdom)

## Power cord socket

The power cord socket (7) is located at the rear right of the QIAcube Connect MDx and allows connection of the QIAcube Connect MDx to a power outlet via the supplied power cord.

### **WARNING** Electrical hazard



Any interruption of the protective conductor (earth/ground lead) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument dangerous.

Intentional interruption is prohibited.

#### **Lethal voltages inside the instrument**

When the instrument is connected to line power, terminals may be live and opening covers or removing parts is likely to expose live parts.

### **WARNING** Damage to electronics



Before powering ON the instrument, make sure that the correct supply voltage is used.

Use of incorrect supply voltage may damage the electronics.

To check the recommended supply voltage, refer to the specifications indicated in the type plate of the instrument.

### **WARNING** Risk of electric shock



Do not open any panels on the QIAcube Connect MDx.

#### **Risk of personal injury and material damage**

Only perform maintenance that is specifically described in this user manual. Any other maintenance or repair may only be carried out by an authorized Field Service Specialist.

## Cooling air outlet

Cooling air outlets are located at the rear left side of the QIAcube Connect MDx and allow cooling of the internal components of the QIAcube Connect MDx.

### **WARNING** Risk of overheating



To ensure proper ventilation, maintain a minimum clearance of 10 cm at the sides and rear of the QIAcube Connect MDx.

Slits and openings that ensure the ventilation of the instrument must not be covered.

## External barcode scanner


The QIAcube Connect MDx is equipped with a 2D handheld barcode scanner to allow kit barcode and sample barcode scanning.

### **WARNING** Risk of personal injury



Hazard Level 2 laser light: Do not stare into the light beam when using handheld barcode scanner.

### 3.3. Internal features of the QIAcube Connect MDx



Internal view of the QIAcube Connect MDx.

- |   |                          |   |   |
|---|--------------------------|---|---|
| 1 | Centrifuge lid           | 6 | Microcentrifuge tube slots  |
| 2 | Centrifuge               | 7 | 3 Slots for tip racks   |
| 3 | Shaker                   | 8 | Disposal slots for tips and columns   |
| 4 | Reagent bottle rack      | 9 | Robotic arm (includes gripper, pipetting system, optical sensor, ultrasonic sensor, and UV LED) |
| 5 | Tip sensor and hood lock |   |   |

#### Centrifuge

The centrifuge is equipped with 12 swing-out buckets, each of which can hold a disposable rotor adapter. Up to 12 samples can be processed per run. For ease of use and high process safety, a gray line marks the side of the bucket that must face toward the center of the rotor. All centrifuge buckets must be mounted before starting a run, regardless of the number of samples to be processed. Be sure to follow the loading instructions provided by the user interface to ensure correct loading of the centrifuge.

The centrifuge can also be operated individually via the touchscreen (refer to Section 5.9 Independent centrifuge operation).

**Note:** Be sure to follow the centrifuge loading instructions provided by the user interface.

#### **WARNING** Risk of personal injury and material damage



To prevent plastic crash, load the tubes properly. After a plastic crash, sharp plastic particles could be inside the centrifuge. Be careful when handling items inside the centrifuge.


#### **WARNING** Risk of personal injury and material damage



Raise the centrifuge lid carefully. The lid is heavy and may cause injury if it falls.

## Shaker

The heated orbital shaker enables fully automated lysis of up to 12 samples. Two kinds of shaker adapters are available for 2 mL microcentrifuge tubes (imprinted “2”) and 2 mL screw-cap tubes (imprinted “S2”). Sample tubes are placed into a rack that fits onto the shaker adapter. The lid of each microcentrifuge tube or shaker rack plug of each screw-cap tube is held in a slot at the edge of the shaker rack. This ensures that microcentrifuge tubes cannot be displaced during sample processing and allows shaker loading to be checked. The shaker can also be operated individually via the touchscreen (see Section 5.8 Independent heater/shaker operation).



Sample tube lids are held in slots at the edge of the shaker rack

Shaker rack with 2 mL microcentrifuge tubes.

**Note:** For shaker loading follow instructions provided by the software.

### **WARNING** Risk of personal injury and material damage



To prevent plastic crash, load the tubes properly. After a plastic crash, sharp plastic particles could be inside the centrifuge. Be careful when handling items inside the centrifuge.

### **WARNING** Hot surface



The shaker can reach temperatures of up to 70°C (158°F). Avoid touching it when it is hot, in particular, shortly after a run has been carried out.

## Reagent bottle rack

The reagent bottle rack holds up to six 30 mL QIAcube Connect MDx-specific reagent bottles and, for ease of use and high process safety, fits onto the QIAcube Connect MDx worktable only in the correct orientation. Liquid is aspirated from the bottles by the pipetting system. A labeling strip should be attached to the reagent bottle rack. For increased convenience and ease of use, the labeling strip fits only in the correct orientation on the reagent bottle rack. The use of the labeling strip ensures that the rack is correctly positioned on the worktable for liquid-level detection.

**Note:** Reagent bottles designed for use with the QIAcube Connect MDx and supplied by QIAGEN must be used. Otherwise, errors may occur during liquid detection.



Reagent bottle rack with white labeling strips at the sides.

**WARNING** Risk of fire or explosion



When using ethanol or ethanol-based liquids on the QIAcube Connect MDx, handle such liquids carefully and in accordance with the required safety regulations. If liquid has been spilled, wipe it off and leave the QIAcube Connect MDx hood open to allow flammable vapors to disperse.

**WARNING** Risk of explosion



The QIAcube Connect MDx is intended for use with reagents and substances supplied with QIAGEN kits or other than outlined in respective instruction for use. Use of other reagents and substances may lead to fire or explosion.

**WARNING** Risk of personal injury and material damage



To prevent plastic crash, load the tubes properly. After a plastic crash, sharp plastic particles could be inside the centrifuge. Be careful when handling items inside the centrifuge.

**WARNING** Hot surface



The shaker can reach temperatures of up to 70°C (158°F). Avoid touching it when it is hot, in particular, shortly after a run has been carried out.

### Tip sensor

During sample preparation, the tip sensor checks that the tip adapter has picked up a tip and checks whether it is a 200 µL or a 1000 µL type filter-tip.

### Microcentrifuge tube slots

In addition to the 12 tubes that the shaker can accommodate, up to 3 additional microcentrifuge tubes can be used in the microcentrifuge accessory position. These slots are used by applications in which, for example, proteinase K or another enzyme is required for the purification protocol.

**Note:** There is no liquid level detection on those slots. Ensure loading of the exact volume specified on the user interface.

### Tip rack slots

Three tip racks can be placed on the QIAcube Connect MDx worktable. Tips can be purchased in prefilled tip racks holding 200 µL filter-tips or 1000 µL filter tips, regular or wide-bore.

**Note:** Only filter-tips designed for use with the QIAcube Connect MDx and supplied by QIAGEN must be used. Do not refill the racks manually.

### Disposal slots for tips and columns

Disposable filter-tips are alternately discarded through each of the round-shaped tip disposal slots into the waste drawer. This prevents discarded tips from piling up in the waste drawer.

Used columns (e.g., QIAshredder columns) are discarded through the square-shaped disposal slot into the waste drawer.

### Robotic arm

The robotic arm provides accurate and precise positioning of the robotic gripper and pipetting system on the QIAcube Connect MDx worktable and includes an optical and ultrasonic sensor, as well as a UV LED.

#### **WARNING** Moving parts



To avoid contact with moving parts during the operation of the QIAcube Connect MDx, the instrument must be operated with the hood closed.

If the hood sensor or lock is not functioning properly, contact QIAGEN Technical Services.

#### **WARNING** Moving parts



Avoid contact with moving parts during operation of the QIAcube Connect MDx. Do not place your hands under the robotic arm when it is lowering. Do not attempt to move any tip racks or tubes while the instrument is operating.

### Robotic gripper

The robotic gripper transfers spin columns. During transfer of a spin column, a stabilizing rod holds the rotor adapter in place, ensuring it remains properly seated in the centrifuge bucket. The robotic gripper is behind the panel covering the robotic arm.



Robotic gripper

Stabilizing rod

**Robotic gripper automates spin-column processing.**

### Pipetting system

The QIAcube Connect MDx is equipped with a single-channel pipetting system that moves in the X, Y, and Z directions. The pipetting channel, fitted with a tip adapter, is connected to a precision syringe pump, which enables accurate transfer of liquids. The tip adapter allows aspiration and dispensing of liquid through an attached disposable tip. Disposable filter-tips (200, 1000, and 1000 µL wide-bore) are used for sample processing to minimize the risk of cross-contamination.

### Optical sensor

During the load check, the optical sensor checks that the number of rotor adapters corresponds to the number of samples in the shaker and that both shaker and rotor are correctly loaded. The optical sensor also checks the type of tips loaded on the worktable and whether there are enough tips for the protocol run.

### Ultrasonic sensor

During the load check, the ultrasonic sensor checks that the buffer bottles in the reagent bottle rack contain enough buffer for the protocol run.

**Note:** The ultrasonic sensor has a black beam collimator. If, for any reason, this beam collimator falls off or is missing, the instrument will display an error message to inform the user that the beam collimator is missing and that runs cannot be started. To replace the beam collimator, it needs to be manually adjusted to its original position (see the image below). If you are still having difficulties and the error message remains, contact your local technical service department for further support.



Black beam collimator (see red circle) of the ultrasonic sensor.

### Speakers

The system is equipped with speakers with different audio signals to notify different instrument states such as:

- Run completed
- Error
- Aborted run



## Interior LED

The QIAcube Connect MDx is equipped with an interior LED that illuminates the worktable for ease of use. The LED can also indicate the instrument status (e.g., error) by blinking.

## UV LED

The QIAcube Connect MDx is equipped with a UV LED light for supporting decontamination. During the maintenance decontamination procedure, the UV LED is moved over the worktable. The hood and the waste drawer must be closed before starting the procedure and must not be opened during the procedure.

### **WARNING** Risk of personal injury




Do not expose your skin to UV-C light from the UV LED lamp.

## 3.4. Disposables

### Rotor adapter

A disposable rotor adapter holds a QIAGEN spin column and a microcentrifuge tube in a centrifuge bucket during sample processing. If required by the protocol, an additional column (e.g., QIAshredder column) can be placed in the middle position of the rotor adapter. For ease of use and high process safety, the rotor adapters are designed so that they fit into a centrifuge bucket only in the correct orientation. Spin-column and microcentrifuge tube lids are held securely in slots at the edge of the rotor adapter.



Assembly of a rotor adapter.

The wash position of the rotor adapter is open at the bottom, enabling wash buffers to flow through and collect at the bottom of the rotor adapter during centrifugation. The other two positions in the rotor adapter are closed. Be sure to follow the loading instructions provided by the user interface.

### **WARNING** Risk of personal injury and material damage



Do not use damaged rotor adapters. The rotor adapters can only be used once. High *g* forces exerted in the centrifuge can cause damage to used rotor adapters.

### **WARNING** Risk of personal injury and material damage



To prevent plastic crash, load the tubes properly. After a plastic crash, sharp plastic particles could be inside the centrifuge. Be careful when handling items inside the centrifuge.

## 4. Installation Procedures

This section provides instructions on installation environment requirements as well as unpacking, installing, configuring, and packing the QIAcube Connect MDx.

### 4.1. Installation environment

#### 4.1.1. Site requirements

The QIAcube Connect MDx must be located out of direct sunlight, away from heat sources and away from sources of vibration and electrical interference. Refer to Technical Specifications for the operating conditions (temperature and humidity). The site of installation should be free of excessive drafts, excessive moisture, and excessive dust, and should not be subject to large temperature fluctuations.

Use a level workbench that is large enough and strong enough to accommodate the QIAcube Connect MDx. Refer to Technical Specifications for the weight and dimensions of the QIAcube Connect MDx.

Ensure that the workbench is dry, clean, and vibration-proof, and has additional space for accessories.

Do not use this device in close proximity to sources of strong electromagnetic radiation (e.g., unshielded intentional RF sources), as these can interfere with proper operation.

The QIAcube Connect MDx must be placed within approximately 1.5 m of a properly grounded (earthed) AC power outlet. The power line to the instrument should be voltage regulated and surge protected. Ensure that the QIAcube Connect MDx is positioned so that it is easy to access the power connector at the back of the instrument and the power switch on the front at all times and that it is easy to power the instrument OFF and disconnect it.

**Note:** It is recommended to plug the instrument directly into its own power outlet and not to share the power outlet with another lab equipment. Do not place the QIAcube Connect MDx on a vibrating surface or near vibrating objects.

#### **WARNING** Explosive atmosphere



The QIAcube Connect MDx is not designed for use in an explosive atmosphere.

#### **WARNING** Risk of overheating



To ensure proper ventilation, maintain a minimum clearance of 10 cm at the sides and rear of the QIAcube Connect MDx.

Slits and openings that ensure the ventilation of the instrument must not be covered.

#### **WARNING** Risk of personal injury and material damage



The QIAcube Connect MDx is too heavy to be lifted by one person. To avoid personal injury or damage to the instrument, do not lift the instrument alone.

#### **CAUTION** Damage to the instrument



Direct sunlight may bleach parts of the instrument and cause damage to plastic parts. The QIAcube Connect MDx must be located out of direct sunlight.

### 4.1.2. Power requirements

The QIAcube Connect MDx operates at: 100–240 V AC, 50/60 Hz, 650 VA.

**Note:** The apparent power can exceed 650 VA for up to 2 seconds during the centrifuge acceleration and can reach an approximate value of 1200 VA. The QIAcube Connect MDx can be connected to an uninterruptible power supply (UPS).

Minimum specifications for UPS:

<b>Power capacity</b>	1200 VA
<b>AC voltage</b>	220–240 VAC 100–120 VAC
<b>Frequency</b>	50/60 Hz
<b>Waveform</b>	Pure sine wave

Ensure that the voltage rating of the QIAcube Connect MDx is compatible with the AC voltage available at the installation site. Mains supply voltage fluctuations are not to exceed 10% of nominal supply voltages.

#### **WARNING** Damage to electronics



Before powering ON the instrument, make sure that the correct supply voltage is used.

Use of incorrect supply voltage may damage the electronics.

To check the recommended supply voltage, refer to the specifications indicated in the type plate of the instrument.

#### **WARNING** Electrical hazard



Any interruption of the protective conductor (earth/ground lead) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument dangerous.

Intentional interruption is prohibited.

##### **Lethal voltages inside the instrument**

When the instrument is connected to line power, terminals may be live and opening covers or removing parts is likely to expose live parts.

### 4.1.3. Grounding requirements

To protect operating personnel, the National Electrical Manufacturers' Association (NEMA) recommends that the QIAcube Connect MDx be correctly grounded (earthed). The instrument is equipped with a 3-conductor AC power cord that, when connected to an appropriate AC power outlet, grounds (earths) the instrument. To preserve this protection feature, do not operate the instrument from an AC power outlet that has no ground (earth) connection.

#### **WARNING** Electrical hazard



Any interruption of the protective conductor (earth/ground lead) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument dangerous.

Intentional interruption is prohibited.

##### **Lethal voltages inside the instrument**

When the instrument is connected to line power, terminals may be live and opening covers or removing parts is likely to expose live parts.

## 4.2. Unpacking the QIAcube Connect MDx

### **WARNING** Risk of personal injury and material damage



The QIAcube Connect MDx is too heavy to be lifted by one person. To avoid personal injury or damage to the instrument, do not lift the instrument alone.

1. Before unpacking the QIAcube Connect MDx, move the package to the site of installation and check that the arrows on the package point upward. In addition, check whether the package is damaged. In case of damage, contact QIAGEN Technical Services.
2. Open the top of the transportation box to remove the *QIAcube Connect MDx Quick-Start Guide*, packing list document, certificate of manufacturing, and power cords before lifting the box.



#### **Documents and power cord location.**

3. Remove the black foam protector lid and lift the box.
4. When lifting the QIAcube Connect MDx, slide your fingers under both sides of the workstation and keep your back straight.  
**Important:** Do not hold the touchscreen display while unpacking or lifting the QIAcube Connect MDx, this might damage the instrument.
5. Carefully remove the instrument from the protective film shipping bag, including the silica gel pack.
6. Check that the QIAcube Connect MDx is not damaged and that there are no loose parts. If anything is damaged, contact QIAGEN Technical Services. Make sure that the QIAcube Connect MDx has equilibrated to ambient temperature before operating it.
7. Retain the package in case you need to transport the QIAcube Connect MDx in the future. Refer to Section 4.4, Repackaging and shipping of the QIAcube Connect MDx, for more details. Using the original package minimizes the possibility of damage during transportation of the QIAcube Connect MDx.

## 4.3. Installing the QIAcube Connect MDx

This section describes important actions that must be performed before operating the QIAcube Connect MDx. These actions include:

- Removal of the QIAcube Connect MDx accessories and shipping material (also see Section 4.2 Unpacking the QIAcube Connect MDx)
- Installation of the AC power cord
- Installation of external barcode scanner
- Installation of the centrifuge rotor and buckets
- If an installation qualification (IQ/OQ) is required in your laboratory setup, this service can be ordered together with the instrument. For details, please contact QIAGEN Technical Services.

### 4.3.1. Removal of the QIAcube Connect MDx accessories and shipping material

1. Remove the USB flash drive, rotor key, rotor nut, Allen key, S2 shaker adapter and shaker rack plugs from the waste drawer.
2. Read the packing list to check that you have received all items. If anything is missing, contact QIAGEN Technical Services.
3. To remove the foam protectors above the centrifuge, gently pull the top foam protector towards you (see the picture below). After removal of the top foam protector, gently pull the middle foam protector towards you and find the barcode reader enclosed within the middle foam (see the picture below). Gently remove the bottom foam protector above the centrifuge.
4. To remove the foam protector around the robotic arm, gently pull the foam protector towards you (see the picture below). Carefully push back the robotic arm to reveal and remove the small foam pad underneath (see the picture below). After removing the protector for the robotic arm, ensure to close the QIAcube Connect MDx hood.



Foam protector above centrifuge.



Barcode scanner enclosed within middle foam protector.



Foam protector for robotic arms.



Foam protector under robotic arms.

5. Carefully peel off the protective film from the QIAcube Connect MDx hood.

#### 4.3.2. Installation of AC power cord

1. Take the power cord previously removed from the foam packing material on top of the QIAcube Connect MDx.

**Note:** Only use the power cord that is supplied with the QIAcube Connect MDx.

2. Ensure that the power switch is set to OFF: outer position is OFF and inner position is ON.
3. Check that the voltage rating on the label at the back of the QIAcube Connect MDx matches the voltage available at the installation site.
4. Plug the power cord into the instrument power cord socket.
5. Plug the power cord into a grounded power outlet.

#### **WARNING** Damage to electronics



Before powering ON the instrument, make sure that the correct supply voltage is used.

Use of incorrect supply voltage may damage the electronics.

To check the recommended supply voltage, refer to the specifications indicated in the type plate of the instrument.

#### **WARNING** Electrical hazard



Any interruption of the protective conductor (earth/ground lead) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument dangerous.

Intentional interruption is prohibited.

#### **Lethal voltages inside the instrument**

When the instrument is connected to line power, terminals may be live and opening covers or removing parts is likely to expose live parts.

### 4.3.3. Installation of external barcode scanner

1. Remove the barcode scanner from the box.
2. Plug the USB connector of the scanner into one of the USB ports located on the left of the QIAcube Connect MDx touchscreen.

### 4.3.4. Installing the centrifuge rotor and buckets

The centrifuge rotor and buckets are preinstalled in the QIAcube Connect MDx. When setting up the QIAcube Connect MDx for the first time, switch on the instrument (see Section 5.2) and remove the transportation foam inserts from the centrifuge after it opens. In case the centrifuge rotor and buckets have been manually removed (e.g., during maintenance), follow the instructions below to re-install them.

1. The rotor can be mounted in only one orientation. The pin on the rotor shaft fits into a notch on the underside of the rotor directly underneath rotor position 1.
2. Line up position 1 of the rotor with the pin on the rotor shaft and carefully lower the rotor onto the shaft.
3. Install the rotor nut on top of the rotor and tighten securely using the rotor key supplied with the QIAcube Connect MDx. Make sure that the rotor is securely seated. If the rotor nut is not tightened properly, it can become loose during operation of the centrifuge and can cause serious damage to the instrument. Such damage is not covered by the warranty.
4. Insert the rotor buckets. The side of the rotor bucket that must face toward the rotor shaft is marked with a gray line. Hold the bucket at an angle with the gray line facing the center of the rotor and hang the bucket on the rotor. Check that all buckets are properly suspended and can swing freely.

**Important:** All centrifuge buckets must be mounted before starting a run.

Before starting next protocol run, follow the instructions in Section 6.6.3 Operating the centrifuge after cleaning.

#### **WARNING** Risk of personal injury and material damage



To prevent the rotor nuts from loosening during operation of the centrifuge, securely tightened the nuts using the rotor key supplied with the QIAcube Connect MDx.

#### **WARNING** Risk of personal injury and material damage



Raise the centrifuge lid carefully. The lid is heavy and may cause injury if it falls.

#### **CAUTION** Damage to the instrument



The QIAcube Connect MDx must not be used if the centrifuge lid is broken, or if the lid lock is damaged. Make sure that no loose material is inside the centrifuge during operation.

Make sure that the rotor is installed correctly and that all buckets are properly mounted, regardless of the number of samples to be processed. Load the rotor only as instructed by the software.

Only use rotors, buckets, and consumables designed for use with the QIAcube Connect MDx. Damage caused by use of other consumables will void your warranty.

We recommend replacing the centrifuge rotor and buckets after 20,000 cycles, which is equivalent to 9 years of usage with two runs per day for 220 days each year. For more information contact QIAGEN Technical Services.




### 4.3.5. Installing the shaker adapter

A shaker adapter must be installed before the shaker can be used. Two types of shaker adapter are available:

- Adapter for 2 mL microcentrifuge safe-lock tubes (marked with “2”)
- Adapter for 2 mL screw-cap tubes (marked with “S2”).

The shaker adapter to be used is outlined by the loading instructions on the user interface (number in the overview picture and text below the picture: “use shaker type...”).



#### Indication of shaker type on GUI.

The QIAcube Connect MDx is delivered with the shaker adapter for 2 mL microcentrifuge safe-lock tubes already installed. If you need to install the shaker adapter for 2 mL screw-cap tubes, follow these steps:

1. Remove the shaker rack.
2. Remove the shaker adapter for 2 mL microcentrifuge safe-lock tubes by unscrewing the retaining screws. Use the Allen key supplied with the QIAcube Connect MDx.
3. Place the shaker adapter for 2 mL screw-cap tubes onto the shaker.
4. Tighten the 2 retaining screws using the Allen key.

**Note:** Make sure to use the correct adapter as displayed on the touchscreen during run setup. This helps to ensure optimal instrument performance. Using an incorrect shaker adapter can adversely affect pipetting performance and protocol results.



#### 4.3.6. Software upgrade

**Note:** Only administrators can update the software.

**Note:** System Configuration is only available to users with Administrator role.

**Note:** To ensure that your QIAcube Connect MDx has the latest software versions installed, please visit the QIAcube Connect MDx webpage at [www.qiagen.com/QIAcube-Connect-MDx](http://www.qiagen.com/QIAcube-Connect-MDx). Currently installed software version can be found from the **Configuration** menu, under the **System configuration** tab.


**Important:** Use only the USB stick provided by QIAGEN. Do not connect other USB stick devices to USB ports.

**Important:** Use only QIAcube Connect MDx–related files downloaded from [www.qiagen.com](http://www.qiagen.com) or provided by QIAGEN Technical Services.

**Note:** Checksum confirmation is required to secure software integrity after web download was successfully completed and before subsequent handling of the software. For detailed information on confirmation of software integrity during download and file transfer, please check the “QIAGEN software integrity verification process” description document, which is provided next to the software package on the QIAGEN webpage. The checksum provided on the download page is the checksum of the **.zip** package. Ensure to perform the checksum comparison before unzipping.

If an updated software version is available for download, it can be accessed at [www.qiagen.com/QIAcube-Connect-MDx](http://www.qiagen.com/QIAcube-Connect-MDx); see the **Resources** tab. The download creates a ZIP file. If you require translation of the user interface, please also check for the corresponding language packages and download them as well.

The software can only be updated by users assigned with the Administrator role. It is recommended to download all run reports before updating the software and to create a support package, because run reports and support packages will be lost during software update (see Section 5.7 Saving run reports to the USB flash drive and Section 7.1 Creating a support package). In addition, it is strongly recommended to create a backup of all protocol files according to Section 5.10.4. The software update package contains the newest standard protocol package. If you use customized protocols or your process relies on a certain protocol version, those have to be restored from the backup package after the software update (see Section 5.10.1 Installing new protocols via USB stick).

1. On the menu bar, press the **Configuration** icon (.
2. Create a protocol backup according to Section 5.10.4.
3. Press the **System** tab.

4. The currently installed software version is shown at the right.

System configuration screen.

5. On a computer running Microsoft® Windows®, download and transfer the software ZIP file to the main folder of the USB flash drive provided with the QIAcube Connect MDx and extract the ZIP file there.

**Note:** After extraction, make sure the following files are in the main folder of the USB flash drive:

- qcc1.sig
- qcc2.sig
- qcc3.sig
- qcc4.sig
- qiacube1.bin
- qiacube2.bin
- qiacube-connect-<version>.tar.gz
- qiacube-connect-<version>.tar.gz.md5sum.
- **Optional:** A folder called "Language\_Upload" containing the language files corresponding to the software version


The update will not work if one of the files is missing or has been renamed. Make sure that only the files for one software version are in the main folder of the USB flash drive.

6. If a translated user interface is required, download the respective language pack from [www.qiagen.com/QIAcube-Connect-MDx](http://www.qiagen.com/QIAcube-Connect-MDx) (see the **Resources** tab) as well and extract it to the same USB flash drive.


7. Connect the USB flash drive to the instrument using one of the USB ports at the left of the touchscreen.

**Important:** Make sure that all required run reports, support packages and protocols have been backed up before proceeding to the next step. See Section 5.7 Saving run reports to the USB flash drive, Section 7.1 Creating a support package, and Section 5.10.4 Saving protocols.

8. Press **Update Software** to start the software update. Follow the instructions on the screen.
9. An administrator login is required.




10. The software update tool is displayed. If language packages are detected by the system on the same USB device, the drop-down list under "select language" will be active and allow to choose the desired display language. All detected language packs will be installed simultaneously to the software and the selected one will be the display language after reboot. Up to date language packages are available for download on [www.qiagen.com/QIAcube-Connect-MDx](http://www.qiagen.com/QIAcube-Connect-MDx); see the **Resources** tab. Ensure that the language files are stored in a folder called "Language\_Upload" on the USB device.



**Software update tool with active language selection.**


11. If no translation of the user interface is required and no language package is detected on the USB flash drive, the option will be grayed out.
12. Press **Start update** to start the software update. Press **Cancel** if you do not want to update the software. In this case, the instrument will initialize without updating the software.
13. The software update tool screen also offers to reset the system to factory settings. This process is recommended before an instrument is discarded. Ensure to backup all required data before resetting.

14. Wait until the update has completed.
15. When the update has finished the user is prompted to remove the USB stick. Doing so will trigger automatic reboot of the system.




**Software update tool when update is ready.**

16. After a while without user action, the screen will actively instruct to remove the USB flash drive from the USB port.



17. After removing the USB, the instrument will reboot. The following message appears:




The instrument will initialize with the updated software.


18. After the software update, the system automatically checks the firmware versions and updates or downgrades the firmware if it does not match the expected firmware versions.
19. This happens automatically after the reboot and does not require any user interaction. It is outlined by the following screen:



20. If the firmware update has finished successfully, the login screen will be displayed. Wait until you see the login screen.



21. If the software update has finished successfully, it is recommended to delete the installation files from the USB device as they may disturb other file transfer actions. Use a computer running Microsoft Windows to delete the previously downloaded software installation files from the USB flash drive.
22. After successful reboot, the language selected above or the default language “English” is displayed. If you decide later to work with additional languages, which have not been stored on the USB flash drive during software upgrade, the upload of the respective language pack according to Section 4.5.1 (step 7) is required.
23. The software installation contains all standard protocol files. If you want to reinstall the protocols from the protocol backup follow Section 5.10.1.
24. If you do not install all required protocols the following message will occur after scanning a kit barcode.



#### 4.4. Repackaging and shipping of the QIAcube Connect MDx

When repackaging the QIAcube Connect MDx for shipping, the original packaging materials must be used. If the original packaging materials are not available, contact QIAGEN Technical Services. Make sure that the instrument has been properly prepared (see Section 6 Cleaning and Maintenance) prior to packing and that it poses no biological or chemical hazard.

**WARNING** Risk of personal injury and material damage




The QIAcube Connect MDx is too heavy to be lifted by one person. To avoid personal injury or damage to the instrument, do not lift the instrument alone.

Before transporting the QIAcube Connect MDx, the instrument must first be decontaminated. Refer to Section 6.8, Decontaminating the QIAcube Connect MDx, for more details. Then prepare the instrument as follows:

1. Prepare the packing material. Materials required are the cardboard carton, the pallet with foam blocks, the foam lid, and the foam protector for the robotic arm.  
**Note:** The centrifuge lid must be open to allow access to the inside of the centrifuge. If the lid is not open, perform steps 2–5 below. If the lid is already open, proceed to step 6.
2. Close the instrument hood.
3. On the menu bar, press the **Tools** button.
4. In the **Tools** menu, press the **Run Modules** tab.

5. In the **Run Modules** menu, press the **Centrifuge** tab and then press **Open Centrifuge Lid** (grayed out button in the picture below).



6. Undo the rotor nut on top of the rotor using the rotor key, and carefully lift the rotor off the rotor shaft.
7. Place the rotor into the black foam lid provided.
8. Close the hood.
9. On the menu bar, press the **Tools** button.
10. In the **Tools** menu, press the **Run Modules** tab.
11. In the **Run Modules** menu, press the **Centrifuge** tab and then press **Close Centrifuge Lid**.
12. When the centrifuge lid is closed, power OFF the QIAcube Connect MDx and open the hood.
13. Insert the foam protector into the front of the instrument.
14. Press the foam down between the centrifuge and the robotic arm.



Foam protector inserted between the centrifuge and the robotic arm.

15. Push the foam until the rear end touches the back wall of the instrument. Ensure that the arm is held firmly in place and cannot move.
16. Make sure that the QIAcube Connect MDx hood can be closed properly. The hood should lightly brush against the foam.
17. Place the accessories into the waste drawer. The following accessories should be packed in air cushion bags:
  - Rotor key
  - Allen key
  - Rotor nut
  - USB flash drive
  - Wi-Fi USB device – if this was supplied with your QIAcube Connect MDx
  - Shaker rack plugs
  - Shaker adapter
18. Place the handheld scanner into dedicated handheld scanner box.
19. Place the QIAcube Connect MDx onto the pallet and put the black foam lid over the top of the instrument. Place the box onto the instrument.

**Important:** When lifting the QIAcube Connect MDx, slide your fingers under both sides of the instrument and keep your back straight.

**Important:** Do not hold the touchscreen display while lifting the QIAcube Connect MDx, as this might damage the instrument.

**WARNING** Risk of personal injury and material damage



The QIAcube Connect MDx is too heavy to be lifted by one person. To avoid personal injury or damage to the instrument, do not lift the instrument alone.

20. Place the accessories into the black foam lid. The following accessories should be packed in air cushion bags:
  - Rotor with swing-out buckets
  - Power cord
21. Seal the outside edges of the carton with tape to protect against moisture.

**Note:** Using the original package minimizes potential damage during transportation of the QIAcube Connect MDx.




## 4.5. Configuration of the QIAcube Connect MDx

When using the QIAcube Connect MDx for the first time, it is recommended to define the required settings. Other settings can be made later when needed.


For details on using the touchscreen and software, refer to Section 5.1 Use of the QIAcube Connect MDx software.

To configure the QIAcube Connect MDx, follow the steps below.


1. Close the instrument hood.
2. Press the power switch to the inner position to power ON the instrument. The start-up screen appears and the beeper sounds (if enabled in the sound settings). The instrument automatically performs the initialization tests. If the centrifuge lid is closed, it will open.
3. Initially, only one user account is available: the preinstalled default user. Press **OK** on the touchscreen to confirm the message.



4. Initially, if no maintenance has been recorded yet, the maintenance status is initialized using a default file. Press **OK** to confirm the message. The maintenance counter starts after the first run has been performed.




5. The Login screen appears after the initialization.



The login screen features the QIAGEN logo at the top left and a 'Help' icon at the top right. The main heading is 'Please log in to start'. Below this, there are two input fields: 'User ID' with 'Admin' entered, and 'Password' which is empty. Under the password field are two radio buttons: 'Research Mode' (selected) and 'IVD Mode'. A 'Log in' button is positioned at the bottom center.

Login screen.

6. Initially, only one default user is available. In this case, enter "Admin" in both the user ID and password fields using the on-screen keyboard. Touch in the entry field to open the on-screen keyboard.
7. Select which software mode (IVD or Research) to launch and press **Log in**.
- For details on software modes, please also refer to Section 5.1 Use of the QIAcube Connect MDx software.
8. After the first login each user is asked to accept the end user license agreement. Press **Accept**.



The screen displays the 'End user license agreement' with a yellow header. The text of the agreement is as follows:


**QIAGEN End User License Agreement**

IMPORTANT: PLEASE READ THIS END USER LICENSE AGREEMENT CAREFULLY. ACCESSING OR USING QIAGEN SOFTWARE OR ANY COMPONENT OF LICENSED MATERIALS (DEFINED BELOW) OR CLICKING THE "ACCEPT" BUTTON BELOW CONSTITUTES ACCEPTANCE OF THIS AGREEMENT. THE TERMS AND CONDITIONS OF THIS USER AGREEMENT GOVERN YOUR RIGHTS TO THE QIAGEN SOFTWARE, LICENSED MATERIALS AND SERVICES TO BE SUPPLIED BY QIAGEN HEREUNDER.

YOU REPRESENT THAT (1) YOU HAVE READ, UNDERSTAND, AND AGREE TO BE BOUND BY THIS USER AGREEMENT, (2) YOU ARE OF LEGAL AGE TO FORM A BINDING CONTRACT WITH QIAGEN AND (3) YOU HAVE THE AUTHORITY TO ENTER INTO THIS USER AGREEMENT PERSONALLY OR ON BEHALF OF THE COMPANY NAMED AS THE CUSTOMER ON THE ORDERING DOCUMENT (DEFINED BELOW), AND TO BIND THAT COMPANY TO THIS USER AGREEMENT. THE TERM "**CUSTOMER**" REFERS TO THE INDIVIDUAL OR LEGAL ENTITY, AS APPLICABLE, IDENTIFIED AS THE CUSTOMER ON THE ORDERING DOCUMENT. **IF YOU DO NOT AGREE TO BE BOUND BY THIS USER AGREEMENT, YOU MAY NOT ACCESS OR USE QIAGEN SOFTWARE OR ANY COMPONENT OF LICENSED MATERIALS.**

At the bottom, there are two buttons: 'Decline' and 'Accept'.

- After the first login, the system will ask you to change the password for the User Admin. The new password must contain 8–40 characters.

A screenshot of a 'Change Password' dialog box. The dialog has a title bar with the QIAGEN logo on the left and a 'Help' icon on the right. The main area contains three text input fields labeled 'Old password', 'New password', and 'Confirm new password'. At the bottom right are 'OK' and 'Cancel' buttons.

Change password screen.

- For details on how to enter text or numbers, refer to Section 5.1 Use of the QIAcube Connect MDx software.
- Users with Administrator role are allowed to change or reset passwords for every other user including other administrators. We recommend creating at least one additional administrator as a backup for the preinstalled administrator user Admin.
- The Setup screen appears.



Setup screen.

- If you need to return to the Setup screen from another screen, press the **Setup** icon (  ).

### 4.5.1. System configurations

This section describes how to set the following system configurations:

- Name for the QIAcube Connect MDx
- Current date and time
- System language

These settings can only be made by users with Administrator role. When using the QIAcube Connect MDx for the first time, it is recommended to set the current date and time.

**Important:** Changes to the system configurations performed in the Research mode of the software will be automatically applied as well in the IVD mode of the software.


1. On the menu bar, press the **Configuration** icon (⚙️).
2. Press the **System** tab. This tab is only available for users assigned the Administrator role.
3. Optionally, enter a name for the QIAcube Connect MDx. The device name serves as the network/host name when connecting the instrument to the network.

The name can have up to 24 characters: letters A–Z, a–z, digits 0–9, and a hyphen (-).

The name must start with a letter and cannot end with a hyphen (-).

4. In the Date and Time fields, select the current date and enter the current time for the instrument. These are used to track the start and end time of a run and are also part of the run report. Date and time are not synchronized using the network.  
To change the date, touch the **Calendar** icon (📅) and select the date.

5. Use the left and right arrow icons to change the month Then touch the current day and press **OK** to confirm.



Calendar window with date picker.

6. In the Language Setting field, a selection of languages according to countries local requirements is available upon delivery. Select the language of your choice in the drop-down menu to run the software in a translated version. A restart of the instrument is necessary to activate the new language setting.
7. Language packages can also be uploaded via the **Load from USB** button. This process may be required after a software update or service intervention of if new language packages become available. A language package only delivers the translation of the graphical user interface. To see a completely translated user interface in the IVD mode of the software, the translated DSP/IVD protocols have to be uploaded as well. Therefore, the translation of the user interface is a two-

step process. First, the translated graphical user interface is uploaded via the **Load from USB** button next to the Language Setting field. The following steps describe the process. For more details or support, contact QIAGEN Technical Services.


- a. On a computer running Microsoft Windows, download the language packages from the Resources tab at [www.qiagen.com/QIAcube-Connect-MDx](http://www.qiagen.com/QIAcube-Connect-MDx).
- b. Unzip the downloaded package. This will result in the subfolder **Language\_Upload**. Transfer this folder to the main directory of the USB flash drive.
- c. Use the USB flash drive that was shipped with the QIAcube Connect MDx to transfer the language package to the instrument.

**Note:** Do not rename or modify the language files or folder name. Otherwise, they cannot be used.


- d. Connect the USB flash drive to the QIAcube Connect MDx using one of the USB ports at the left of the touchscreen.
- e. If not already done, select the **Configuration** icon (⚙️).
- f. The language package can be uploaded using the **Load from USB** button next to the Language Setting field.

**Note:** Language packages are only compatible to a certain software version. Please ensure to upload the language pack version (visible in file name) corresponding to your installed software version. Compatible files will be specified by a message box that appears after pressing **Load from USB**.

- g. Confirm by pressing **Yes**.



- h. After successful upload the following message appears. Finish the process by pressing **OK**.



**Note:** The Research mode of the software is translated after upload of a language package. However, no translated protocols are available in the Research mode of the software. This means that some screens of the protocol setup and at the end of a protocol run appear in mixed language (English and configured language).

### 4.5.2. Settings configuration

This section describes the optional settings that can be defined by each user:

- Audio volume
- Display brightness
- Internal light color intensity (red, green, blue)
- Internal light (enabled/disabled)

The settings apply to the current user.

1. On the menu bar, press the **Configuration** icon (⚙️).
2. Press the **Settings** tab. This tab is available for all users.



#### Settings screen.

3. To adjust the audio volume, display brightness, or RGB settings for the internal light, touch the desired position on the virtual slider on the screen. For audio volume, a sound is played with the set volume.
4. Check the box next to Use internal light to turn on the LED light inside the instrument. Uncheck the box to turn it off.
5. Adjust the color of your internal light to your needs.

**Note:** Be sure to use only well visible colors as the internal LED has a warning function in case of error.

6. To return to the Setup screen, press the **Setup** icon (⏪).

### 4.5.3. Network configuration and QIASphere Base connection

The QIAcube Connect MDx can be connected to a network, allowing real-time status monitoring on a computer or a mobile device (e.g., a tablet) using the QIASphere App, see [www.qiagen.com/qiasphere](http://www.qiagen.com/qiasphere). This configuration requires connecting the QIAcube Connect MDx as well as the QIASphere Base to your network. Follow the instructions below to configure a wired or wireless network connection first before you connect your QIAcube Connect MDx instrument to your network, either by local area network (LAN)/Ethernet cable, or by wireless connection. For more details on the QIASphere Base network configuration and how to connect both devices with each other, refer to the *QIASphere User Manual*, which is available at [www.qiagen.com](http://www.qiagen.com).

#### Using QIASphere communication on the QIAcube Connect MDx

On the instrument side the user enables QIASphere communication by checking the box “enable QIASphere communication” and entering the QIASphere Base details. A default QIASphere certificate is available on the instrument, but it can be exchanged by QIAGEN Technical Services.

QIASphere Base settings

Enable QIASphere communication: ☒

QIASphere Base URL:

QIASphere Base password:

Nov 09, 2023, 10:58 User: Admin Admin


#### QIASphere settings.

**Important:** The QIASphere Base communicates with the QIASphere App, and it also communicates with the QIASphere cloud. The QIASphere App helps to monitor the instrument status, for example:

- Instrument is running.
- Instrument is available.
- Instrument requires maintenance.

If your QIASphere Base is connected to the QIASphere Cloud, run reports generated by the QIAcube Connect MDx (including sample IDs) are transferred to the QIASphere Cloud. This transfer is encrypted. In case this transfer of information would still not be in accordance with your local regulations or with your laboratory regulations, the connection between the QIASphere Base and the QIASphere cloud needs to be actively deactivated. Disable the connection by unchecking the box “enable QIASphere communication”. To disconnect the connection between the QIASphere Base and the QIASphere Cloud, refer to the instructions provided in the *QIASphere User Manual*.

Messages coming from QIASphere are displayed under the **QIASphere** button in the software header.



#### QIASphere Notification Center.

The notifications can also be selected or deleted from this screen. Notifications being too long to be displayed in the preview will be displayed with "...". The full message becomes visible upon pressing the **Information** icon (i).

Only users assigned with the Administrator role can change the network configuration. It is recommended to consult your network administrator when configuring the network. For communication with QIASphere Base, the outbound TCP port 443 (https) is used; ping is supported. If QIASphere connection is activated, the QIACube Connect MDx sends the following information to the QIASphere Base and to the network:

- Export files
  - Run reports files
  - Hardware counter file
  - Audit trail / Statistics events file
  - Log files
- System status
- Device configuration (MDx)
- Maintenance status
- Protocol list



During a run, the QIAcube sends the following additional information to the QIAsphere Base and to the network:

- Application
- Kit name
- Material information
- Protocol name
- Sample count
- Start Time
- Estimation End Time
- Run status (running, finished successfully)

### Configuring a wired network connection of QIAcube Connect MDx

1. Connect the QIAcube Connect MDx to an LAN by using an Ethernet cable and the RJ45 ethernet port on the rear side of the QIAcube Connect MDx.
2. On the menu bar, press the **Configuration** icon (⚙️).
3. Press the **LAN** tab.
4. To automatically configure the network via DHCP, check the DHCP enabled box. Leave all fields empty when using this setting. The assigned IP address will be displayed below the box.

The screenshot displays the QIAcube Connect MDx configuration interface. At the top, a menu bar includes icons for QIAcube, Setup, Tools, Configuration (selected), QIAsphere, Network, Lock, and Logout. Below this, a row of tabs shows System, Users, Protocols, Settings, LAN (selected), Wi-Fi, and QIAsphere. The main content area is titled "Device network settings" and contains a form with the following fields: IP Address, Subnet Mask, DNS Server, Gateway, and Hardware Address. The DHCP enabled checkbox is checked, and the IP address 10.100.226.119 is displayed below it. The Hardware Address is 50:2D:F4:23:1E:76. At the bottom right of the form are buttons for "Test connection" and "Apply". A footer at the bottom left shows the date and time "May 23, 2024, 13:08", and the bottom right shows the user "User: Admin Admin".

Device network settings screen.

5. To manually configure the network, uncheck the DHCP enabled box. Enter the IP address, Subnet mask, and Gateway in the respective fields, using the IPv4 format shown in the picture below. Entering the DNS server is optional. These settings will not be validated by the QIAcube Connect MDx.

The screenshot shows the QIAcube Connect MDx configuration interface. At the top, there is a menu bar with icons for Setup, Tools, Configuration (selected), QIAsphere, Network, Lock, and Logout. Below the menu bar, there are tabs for System, Users, Protocols, Settings (selected), LAN, Wi-Fi, and QIAsphere. The main content area is titled "Device network settings" and contains a form with the following fields:

IP Address 10.10.10.127	Subnet Mask 255.255.255.0	<input type="checkbox"/> DHCP enabled
DNS Server 8.8.8.8	Gateway 10.10.10.0	Hardware Address 22:46:BF:43:A2:FE

At the bottom right of the form, there are two buttons: "Test connection" and "Apply".

At the bottom of the screen, there is a status bar with the following information:

- Jun 21, 2024, 09:19
- Mode: Research
- User: Admin Admin

Device network settings screen with manual network configuration entries.

6. Press **Apply** to confirm and save the settings made.

### Configuring a wireless connection of the QIAcube Connect MDx

The QIAcube Connect MDx supports only WPA-PSK and WPA2-PSK. Additionally, the SSID of the Wi-Fi network must be visible. Connection to a hidden SSID is not supported. The Wi-Fi USB device may have been delivered to you (availability may differ from country to country based on regulations and approvals). If you have not received a Wi-Fi USB device from QIAGEN, ensure that the Wi-Fi adapter supports IEEE 802.11-2016 standards, including WIFI 4 (802.11n), WIFI 5 (802.11a/c), WPA/WPA2 (802.11i), and WPA3 (SAE). It is recommended to use a Wi-Fi adapter with the RTL8723BU chipset. The adapter should be in accordance with your local laws and regulations. For more information contact QIAGEN Technical Services.

Before Wi-Fi can be configured, the Wi-Fi USB device must be plugged into one of the USB ports behind the touchscreen.

**Important:** Always turn off the QIAcube Connect MDx before you plug in or unplug the Wi-Fi USB device. Plug-and-play of the Wi-Fi USB device while the instrument is turned on is not supported.

1. On the menu bar, press the **Configuration** icon (⚙️).
2. Press the **Wi-Fi** tab.
3. Press **Scan** to scan for available networks. The networks are listed based on their signal strength.
4. Select one of the available networks from the list. The details for the selected network are shown at the right.
5. Enter the wireless network password and press **Connect** to connect to the network. The connected network is marked in the list.

**Note:** If a network is previously configured and a connection is established successfully at least once, the instrument will

connect to that network automatically.

6. To disconnect from Wi-Fi, press **Disconnect**.

Proceed with connecting QIAcube Connect MDx to the QIASphere Base as described below.

### Connecting QIAcube Connect MDx to the QIASphere Base

Follow first the instructions of the *QIASphere User Manual* ([www.qiagen.com/qiasphere](http://www.qiagen.com/qiasphere)) to connect the QIASphere Base to the same local network as the QIAcube Connect MDx is connected to. During this procedure, the QIASphere Base receives an IP address which is required in the following configuration:

1. On the QIAcube Connect touchscreen, press **Configuration**, then press the **QIASphere Base** tab.
2. Ensure that the Enable QIASphere communication box is checked.
3. Enter QIASphere Base IP address in the QIASphere Base URL field.
4. Enter the QIASphere Base password and press **Apply**.
5. Your QIAcube Connect MDx is now connected to the QIASphere Base, and you can proceed with setting up QIASphere according to the *QIASphere User Manual*.

The screenshot displays the QIAcube Connect MDx Configuration interface. At the top, a navigation bar includes icons for QIAGEN, Setup, Tools, Configuration (highlighted), QIASphere, Network, Lock, and Logout. Below this is a secondary menu with tabs for System, Users, Protocols, Settings, LAN, Wi-Fi, and QIASphere (selected). The main content area is titled 'QIASphere Base settings' and contains three fields: 'Enable QIASphere communication:' with a checked checkbox, 'QIASphere Base URL:' with the text 'https://10.136.76.12', and 'QIASphere Base password:' with a masked password field. At the bottom right of the settings area are two buttons: 'Test connection' and 'Apply'.

May 23, 2024, 13:17








User: Admin Admin

#### QIASphere screen for connection with QIAcube Connect MDx.

**Note:** To disable QIASphere Base connection, uncheck Enable QIASphere communication box.

The header of the user interface shows you at all times the network and QIASphere status.

#### Network status icons:

-  No connection (gray LAN icon with “disabled” symbol)
-  LAN connection (blue LAN icon with “check mark” symbol)
-  Wi-Fi connection (blue Wi-Fi icon with “check mark” symbol)
- The blue symbol  is displayed if there are no new notifications in the QIASphere notifications tab or all notifications in the **QIASphere notifications** tab were read and the connection is fine.
- The blue symbol with the red dot  is displayed if new QIASphere notifications are available.
- The yellow “warning” symbol  is displayed if there is a connection problem. Please tap the icon to see the corresponding warning notification.
- The gray icon  with a “disabled” symbol is displayed if the Enabled QIASphere communication checkbox is unmarked.

## 5. Operating Procedures

This section describes how to operate the QIAcube Connect MDx.

Before proceeding, we recommend that you familiarize yourself with the features of the instrument by referring to Sections 3.2 and 3.3.

**Important:** The QIAcube Connect MDx is designed for use with QIAGEN spin columns only. Geometry of spin columns manufactured by other suppliers may not be compatible with the QIAcube Connect MDx.

### **CAUTION** Damage to the instrument



Only use QIAGEN spin columns and the QIAcube Connect MDx-specific consumables with the QIAcube Connect MDx. Damage caused by use of other types of spin columns or chemistries will void your warranty. during operation.

The hood of the QIAcube Connect MDx must remain closed during operation of the instrument. Only open the hood when instructed to do so by the software.

### **WARNING** Moving parts



To avoid contact with moving parts during the operation of the QIAcube Connect MDx, the instrument must be operated with the hood closed.

If the hood sensor or lock is not functioning properly, contact QIAGEN Technical Services.

### **WARNING** Risk of personal injury and material damage



Do not attempt to move the QIAcube Connect MDx during operation.

### **WARNING/CAUTION** Risk of personal injury and material damage



Improper use of the QIAcube Connect MDx may cause personal injuries or damage to the instrument. The QIAcube Connect MDx must only be operated by qualified personnel who have been appropriately trained. Servicing of the QIAcube Connect MDx must only be performed by a QIAGEN field service specialist.

### **WARNING** Risk of personal injury and material damage



Do not use damaged rotor adapters. The rotor adapters can only be used once. High *g* forces exerted in the centrifuge can cause damage to re-used rotor adapters.

### **CAUTION** Damage to the instrument



Empty the tip disposal container prior to use to prevent a tip jam in the waste drawer. Failure to empty the waste container may block the robotic arm that could cause run failure or instrument damage.

### **WARNING** Risk of personal injury and material damage



To prevent plastic crash, load the tubes properly. After a plastic crash, sharp plastic particles could be inside the centrifuge. Be careful when handling items inside the centrifuge.

**CAUTION****Damage to the instrument**

Only use the correct volume of liquids.

Exceeding the recommended volume of liquids may damage the centrifuge rotor or instrument.

**WARNING****Risk of fire or explosion**

When using ethanol or ethanol-based liquids on the QIAcube Connect MDx, handle such liquids carefully and in accordance with the required safety regulations. If liquid has been spilled, wipe it off and leave the QIAcube Connect MDx hood open to allow flammable vapors to disperse.

**WARNING****Risk of explosion**

The QIAcube Connect MDx is intended for use with reagents and substances supplied with QIAGEN kits as outlined in respective information for use. Use of other reagents and substances may lead to fire or explosion.

**CAUTION****Damage to the instrument**

Do not lean against the touchscreen when it is pulled out.

**WARNING****Samples containing infectious agents**

Samples used with the QIAcube Connect MDx may contain infectious agents. Handle such samples with the greatest of care and in accordance with the required safety regulations.

Always wear safety glasses, gloves, and a lab coat.

The responsible body (for example, a laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe, and that the instrument operators are suitably trained and not exposed to hazardous levels of infectious agents as defined in the applicable Material Safety Data Sheets (MSDSs) or the OSHA<sup>1</sup>, \* ACGIH,<sup>†</sup> or COSHH<sup>‡</sup> documents.

Venting for fumes and disposal of waste must be in accordance with all national, state, and local health and safety regulations and laws.

\* OSHA — Occupational Safety and Health Organization (United States of America)

† ACGIH — American Conference of Government Industrial Hygienists (United States of America)

‡ COSHH — Control of Substances Hazardous to Health (United Kingdom)

**WARNING****Hazardous chemicals**

Some chemicals used with the QIAcube Connect MDx may be hazardous or may become hazardous after completion of purification.

Always wear safety glasses, gloves, and a lab coat.

The responsible body (for example, a laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe, and that the instrument operators are suitably trained and not exposed to hazardous levels of infectious agents as defined in the applicable Material Safety Data Sheets (MSDSs) or the OSHA<sup>1</sup>, \* ACGIH,<sup>†</sup> or COSHH<sup>‡</sup> documents.

Venting for fumes and disposal of waste must be in accordance with all national, state, and local health and safety regulations and laws.

\* OSHA — Occupational Safety and Health Organization (United States of America)

† ACGIH — American Conference of Government Industrial Hygienists (United States of America)

‡ COSHH — Control of Substances Hazardous to Health (United Kingdom)

**WARNING** Moving parts



Avoid contact with moving parts during operation of the QIAcube Connect MDx. Do not place your hands under the robotic arm when it is lowering. Do not attempt to move any tip racks or tubes while the instrument is operating.

**WARNING** Hot surface



The shaker can reach temperatures of up to 70°C (158°F). Avoid touching it when it is hot, in particular, shortly after a run has been carried out.

**WARNING/  
CAUTION** Risk of personal injury and material damage



Only perform maintenance that is specifically described in this user manual.

## 5.1. Use of the QIAcube Connect MDx software

The QIAcube Connect MDx gives the option to start a protocol either in the IVD mode of the software (validated IVD applications only) or in the Research mode of the software (MBA and customized protocols only). The usage of IVD protocols is only possible and strictly restricted to the IVD mode of the software. This user manual focuses on operating the QIAcube Connect MDx in the IVD-software mode. For detailed instructions on how to operate the QIAcube Connect MDx using the Research mode of the software, please refer to the *QIAcube Connect User Manual* (available on the QIAcube Connect product webpage under the **Resources** tab).

To change the software mode, the user must first log out from the current software mode before logging in to another mode. At the bottom of the touchscreen, the system shows which software mode is in use.

---

Feb 22, 2023, 13:17

Mode: IVD 

User: Admin Admin

IVD mode footer.

---

Feb 22, 2023, 13:18

Mode: Research 

User: Admin Admin

Research mode footer.



















The QIAcube Connect MDx is operated through a touchscreen, which guides you step-by-step through the correct loading of the work deck and selection of the protocol.

**Note:** The instrument's touchscreen does not support swiping and multigestures.

General functions of the QIAcube Connect MDx touchscreen are described on the next page.


**Note:** A red dot on the pressed button indicates a longer reaction time.

## Touchscreen buttons and icons overview for general use


Button/Icon	Function
	Enables the user to scroll up through a list.
	Enables the user to scroll down through a list.
	The software automatically proceeds to the next screen.
	Returns to the previous screen.
	Returns to previous screen without saving any changes.
	Enables the user to change certain settings (e.g., to edit a user account).
	Enables the user to delete certain settings (e.g., to delete a user).
<b>Text fields</b>	Enables to edit text or value. A pop-up keyboard enables these changes.
<b>Row in tables</b>	Can be pressed to select the respective row. Either item will be selected, or row will be highlighted.
	Press to show additional information to respective item
	Press to show important information that must be followed during run setup for respective item.
	Press to show additional information to respective item
	Press to show important information that must be followed during run setup for respective item.
	Navigation back to the Setup screen
	Tools/Maintenance functions
	Configuration
	QIASphere connection is enabled and works properly.
	QIASphere connection is enabled, but there is some network or configuration issue. Click on the icon to see detailed information. It is highly recommended to solve the QIASphere connection issue or otherwise to disable of the QIASphere connection to avoid software instability.
	QIASphere connection is disabled.
	Log out from the instrument



To enter text or numbers, touch the respective field. The corresponding on-screen keyboard will be displayed. See the examples below:




Keyboard for entering a sample name.



Keyboard for editing a protocol parameter.











For protocol parameters, the value range is displayed. In the example shown in the screenshot above, values from 50 to 100 µL can be entered, but only in steps of 10 µL.

If the entered value is not correct, the field border will change to red and the allowed input range is shown in red. It is not possible to proceed to the next screen in this case. Press the field again and correct the value according to the range displayed next to the field.



Buttons and icons in the on-screen keyboard are described below.

Touchscreen buttons and icons in the on-screen keyboard

Button/Icon	Function
	Remove character to the left.
	Clear all from the field.
	To type the next letter in upper-case. After the letter has typed, the keyboard will show lower-case letters again.
	Switch to upper-case letters. Allows to type a number of upper-case letters. Press the symbol again to return to lower-case letters.
	Show special characters.
	Return to letters.
 	Navigate through the text in the input field.
	Confirm and close.
	Discard and close.

## 5.2. Switching the QIAcube Connect MDx on and off

### Switching on the QIAcube Connect MDx

1. Close the instrument hood.
2. Switch on the instrument by pressing the power switch inside (the button will stay at the inner position). You will hear a sound (if sound settings are enabled), and the start-up screen appears. The instrument automatically performs initialization tests. If the centrifuge lid is closed, it will open. If initialization is finished without error, it is ensured that the instrument is properly installed and operating as intended.

### Switching off the QIAcube Connect MDx

Switch off the instrument by pressing the power switch inside (the button will move back to the outer position).


**Note:** After powering OFF the QIAcube Connect MDx, you must wait for a few seconds before you power ON the instrument again. The system might fail to start if you do not allow the QIAcube Connect MDx to rest for a few seconds before powering ON.

## 5.3. Logging in and logging out

### Logging in

1. Close the instrument hood.
2. Switch on the instrument.

After the initialization is complete, the Login screen appears.




Login screen.

3. Enter the user ID and password using the on-screen keyboard.
4. Select software mode (IVD or Research).

**Note:** The selected software mode is shown at the bottom of the screen as long as you are logged in.

5. Press **Log in**.
6. The Setup screen appears.

In case of a failed login, an exclamation icon (❗) and information screen will appear. Touch the respective field to enter the user ID and password again, ensuring that you enter the information correctly. The user ID is case sensitive.



Information screen from failed login, for example, due to entering incorrect password.

## Logging out

1. Press **Logout** at the top right of the screen.
2. To log out, confirm the message with **OK**. To stay logged in, press **Cancel**.

**Note:** The system will log out automatically if you are inactive for a certain time. The administrator can set the number of minutes before automatic logout (see Section 5.11.2).

3. The login screen appears.

**Note:** In case of automatic log out, only the same operator who previously used the instrument or an administrator can log in again. In the case, that another user logs in, the configuration settings from the previous user will be applied.

## 5.4. Setting up a protocol run

All released QIAGEN standard protocols are installed on the QIAcube Connect MDx upon delivery. Available QIAGEN IVD protocols can also be downloaded from the Resources tab of the QIAcube Connect MDx product page: [www.qiagen.com/QIAcube-Connect-MDx](http://www.qiagen.com/QIAcube-Connect-MDx)

For instructions on installing downloaded protocols, see Section 5.10.1 Installing new protocols via USB stick.

**Important:** Before starting any protocol, thoroughly read the relevant QIAGEN kit handbook.

Protocol setup starts from the Setup (➤➤➤) screen.




Entering the setup screen.

The touchscreen software will guide you through the protocol run setup and worktable loading steps. The display screens vary depending on the protocol in use and may look different from the screens shown in this section.

**Note:** If you need to pause the setup, you can click the **lock** icon (🔒) to lock the screen. To unlock the screen, you need to enter your credentials. Only the same operator who previously used the instrument or an administrator can unlock the screen.


1. To start the run setup in the IVD mode, scan the 2D barcode of the QIAGEN kit. Some kits also have a 1D barcode on their label. Please ensure to use the 2D barcode (QR-code) to enter kit information. Press **Scan kit** and then use the handheld scanner. It is also possible to just scan the barcode without pressing **Scan kit**.



Setup screen.

The following information from the kit barcode will be added to the run report created at the end of the run (if no kit barcode is scanned in the Research mode of the software, no kit information will be in the run report):

- Kit name
- Material number
- Lot number
- Expiry date
- The kit is expired (using an expired kit, will flag all samples of the run as Invalid and the following warning will appear).



**Note:** Only use QIAGEN IVD kits within their expiration date. The run is not anymore valid if you use a kit with an expired shelf life, and accordingly, the run results cannot be used for diagnostic use.


Moreover, it is not recommended to start a run at the end of the day (overnight run) and collect the eluates on the next day. As the QIAcube Connect MDx has no technical capability to cool the eluates after the completion of a run, the quality of the eluates might be compromised after a longer storage time at room temperature.

If scanning of the kit barcode fails, you can also type the kit barcode via the user interface. The code is structured as follows:

#### Kit barcode structure

Position	Length	Value	Description
1 to 2	2	01	Identifier "GTIN"
3 to 16	14		GTIN, not used by the system, see label
17 to 18	2	17	Identifier "Expiry date"
19 to 24	6		Expiry date (YYMMDD), see label. If not used: 000000
25 to 26	2	10	Identifier "Lot"
27 to "]"	4 to 10		Lot number, variable length, see label
	1	]	Marker for the end of the lot number
	3	240	Identifier "Product code"
after "240"	0 to 15		Material number (REF), either contains a catalog number or a material number, see label.

The sample barcode label below means that for this label, the barcode string would be 010405322800290117181231101151234567]24061704:



#### Example of kit barcode label.

- After scanning of the 2D barcode, the software will automatically proceed to the next screen. If you scanned a kit barcode, the software may skip the Kit, Material, and/or Protocol Selection screens. The software will skip the selection screen if the required information is provided by kit barcode scanning.

Prepare the samples to be processed using the protocol indicated in the protocol selection screen. Refer to the kit handbook for more information, if required. The needed sample pretreatment can be found in the respective kit handbooks.

To enter information in the next screens, follow the instructions provided in the sections below. Depending on your selections, the number and order of the screens displayed on your instrument may vary.


Each section below contains an exemplary screenshot image. Follow the instructions in that section with the corresponding screen displayed on your instrument.

In general, press Next to proceed to the next screen or press **Back** to return to the previous screen. **Next** will only be active if all required information has been entered in the current screen.

**Note:** Many screens have **Up** (▲) and **Down** (▼) arrow icons to scroll. Make sure to scroll all texts down to the end and follow the instructions completely.

#### 5.4.1. Kit selection (only in Research mode of software)

This step only available in the Research mode of the software and replaces the scanning the 2D barcode of DSP kits in the IVD mode of the software.



##### Select kit screen.

1. Use the **Up** and **Down** arrow icons (▲ and ▼) to scroll through the list of kits.  
Select the kit to be used for your run by touching the corresponding row. Only one kit can be selected per run. Information about the selected kit is shown in the right pane.
2. Press **Next** to proceed to defining the sample material.

## 5.4.2. Material selection

QIAcube Connect MDx

Setup Tools Configuration QIAsphere Network Lock Logout

### Select material

QIAamp DNA Mini Step 2 of 9

Bacteria (Gram+) or yeast	<input type="radio"/>
Bacterial pellet	<input type="radio"/>
Blood or body fluid	<input checked="" type="radio"/>
PAXgene Saliva	<input type="radio"/>
Tissue	<input type="radio"/>

Cancel Back Next

May 23, 2024, 13:22 User: Admin Admin

### Select material screen.

1. Select the sample material by touching the corresponding row. Only one sample material type can be selected per run.
2. Press **Next** to proceed to define the protocol.

## 5.4.3. Protocol selection

QIAcube Connect MDx

Setup Tools Configuration QIAsphere Network Lock Logout

### Select protocol

QIAamp DNA Mini Standard Step 3 of 9

Manual lysis	<input type="radio"/>
Standard	<input checked="" type="radio"/>

Cancel Back Next

Version: 1

This protocol is for purification of total DNA from up to 200 µl of whole blood, plasma, serum, buffy coat, or body fluids, or up to 5 x 10e6 lymphocytes in 200 µl PBS using the QIAamp DNA Mini Kit.

**Input material:** 200 µl blood or body fluid

May 23, 2024, 13:22 User: Admin Admin

### Select protocol screen.

1. Select the protocol by touching the corresponding row. Only one protocol can be selected per run.  
**Important:** Ensure to read all essential and critical information in the right pane (scroll down if necessary) before proceeding to the next step.
2. Press **Next** to proceed to define the run parameters.



#### 5.4.4. Parameter definition

Depending on the selected protocol, certain parameters must be defined. Some protocols do not allow the modification of parameters. These parameters are fixed as they are validated for the procedure. For protocols with editable parameters, default settings are defined but can be changed. Follow the instructions in the information pane on the right regarding changing the values and which increments can be used.

Parameter name	Value
1st elution vol (Default: 100 µl)	100 µl
2nd elution vol (Default: 100 µl)	100 µl

QIAamp DNA Mini  
Standard Step 4 of 9

1st elution vol 50-100 µl in increments of 10 µl  
2nd elution vol 50-100 µl in increments of 10 µl

Cancel

Back Next

May 23, 2024, 13:23 User: Admin Admin

##### Define parameters screen.

1. If required, press the Value field to change a parameter value using the on-screen keyboard. Refer to Section 5.1 for details about the on-screen keyboard.
2. Press **Next** to proceed to define the sample number. The software automatically proceeds to the next screen. Follow the instructions in the corresponding section below.

### 5.4.5. Sample number definition

QIAcube Connect MDx

Setup Tools Configuration QIAsphere Network Lock Logout

### Select the number of samples

1 2 3 4 5 6

7 8 9 10 11 12

Cancel

QIAamp DNA Mini Standard Step 5 of 9

The numbers of samples available for selection allow the centrifuge to be correctly balanced.

Back Next

May 23, 2024, 13:23 User: Admin Admin

#### Define sample numbers screen.

1. To select the number of samples for the run, press the corresponding number on the screen. Sample numbers (1 and 11) that would lead to a misbalance during centrifugation cannot be chosen.
2. Press **Next** to proceed to loading the reagents.

### 5.4.6. Loading buffer bottles

QIAcube Connect MDx

Setup Tools Configuration QIAsphere Network Lock Logout

### Load buffer bottles

Pos.	Name	Volume	Info
2	Buffer AL	≥ 5 ml	i
3	100% ethanol	≥ 5 ml	i
4	Buffer AW1	≥ 8 ml	i
5	Buffer AW2	≥ 8 ml	i
6	Buffer AE	≥ 5 ml	i

Extra kits Add Delete Show

Cancel

QIAamp DNA Mini Standard Step 6 of 9

**Bottle Type**  
30 ml, Cat. No 990393

1 2


3 4

5 6

Back Next

May 23, 2024, 13:24 User: Admin Admin

#### Load buffer bottles screen.


The Load buffer bottles screen guides you through setting up the required buffers for the run. Make sure to read all essential and critical information before proceeding to the next step. If required, click  **Add** to add extra kits via 2D barcode scanning.

**Note:** To avoid run issues and to guarantee proper sitting of buffer bottle rack, the buffer bottle rack must be equipped with rack labeling strips. Please ensure that the buffer bottle rack is sitting properly by pressing it down.

Depending on the selected protocol, loading of buffer bottles might not be required. In this case, the software will indicate that this step can be skipped.


1. Prepare the reagents as shown on the screen. Refer to the respective kit handbook for further information to ensure that the correct buffers are used on each position (see the blue-shaded circle on the touchscreen). When pouring, ensure that the buffer does not foam or contain large air bubbles.

**Note:** Use a volume of reagents as close as possible to the volumes required for the selected protocol and the respective number of samples processed (as indicated on the reagent table on the touchscreen). Do not use less than 5 mL of buffer.



Load buffer bottle screen circles the position on the buffer bottle rack for selected buffer in the list.

Before proceeding to the next step, ensure that you read all essential and critical information indicated in the **Information** icon (i) marked in red. Press the icon to open the information.



Example of message box displayed by hitting the information icon (i).

2. Make sure that the buffer bottles contain the minimum volumes described in the volume column. Each bottle can hold a maximum volume of 30 mL, which can also be used for subsequent runs. However, it is recommended not to use much more than the minimal volume. The physical mark of the bottle (30 mL) must not be exceeded. Later, when the run is started, the instrument will determine the filling volume.
3. Ensure to label the buffer bottles properly and according to safety requirements. Buffer bottles can be stored according to the storage conditions described in the kit handbooks. However, prolonged times of buffer bottles standing open on the instrument must be avoided. For subsequent runs, fresh buffer should be refilled. We recommend reusing the buffer bottles only until one kit is used up. As soon as a new QIAGEN kit is opened, new buffer bottles should be used.
4. Place each opened buffer bottle in the correct position of the reagent bottle rack as shown on the screen. The buffer bottle rack positions are numbered for easy identification.
5. Once all the buffer bottles are placed on the reagent bottle rack, place the rack on the worktable. Ensure that the rack is oriented correctly with number 1 on the top. The rack will only fit on the worktable if its orientation is correctly set.

**Important:** Be sure to place the reagent bottle rack correctly into the designated worktable slot. Bottle racks that are in a leaning position could cause errors during liquid detection.

**Important:** Ensure that buffer bottles are open. A closed buffer bottle will be detected by the instrument and prevent the run from being started.

6. Press **Next** to proceed to loading the tips and enzymes. The software automatically proceeds to the next screen. Follow the instructions in the corresponding section below.

**WARNING** Risk of fire or explosion




When using ethanol or ethanol-based liquids on the QIAcube Connect MDx, handle such liquids carefully and in accordance with the required safety regulations. If liquid has been spilled, wipe it off and leave the QIAcube Connect MDx hood open to allow flammable vapors to disperse.

### 5.4.7. Loading tips and enzymes

**Important:** When the Loading tips and enzymes screen appears, the robotic arm will automatically move slowly – even when the instrument hood is open – so you can access all of the loading positions. Always stand clear of the instrument while the robotic arm is moving. Wait until the robotic arm has completed its movements before you start to load or unload the tip racks or enzymes. After you are finished loading and you proceed from this screen, the robotic arm will automatically move back to its original position (above tip rack position 3).

If more than one rack of the same tip type is loaded, the instrument will use the tip rack placed in position 1 first, then continue to position 2, and then position 3. To use a partially filled rack first, load it in position 1.

Depending on the selected protocol, loading of tips and enzymes might not be required. In this case, the software will indicate that this step can be skipped.



#### Loading tip racks and enzymes screen.

If for any reason the robotic arm prevents you from reaching a loading position, do not move the robotic arm manually. Instead, proceed as follows:

- Press **Move left** or **Move right**. The robotic arm will start to move. The hood can remain open during this movement.
- Ensure that you stand clear of the instrument while the robotic arm is moving. Wait until the robotic arm has completed its movements.

Follow the instructions below to load enzyme, reagents and tips:

1. Prepare enzymes and/or reagent(s) listed on the screen. Refer to the corresponding kit handbook for more information. Before proceeding to the next step, make sure that you read and follow all essential and critical information indicated in the **Information** icon () marked in red.
2. Make sure that you are using the correct tube type. Press the **Information** icon () in the respective row to show the details.

3. Supported enzyme tube types are 1.5 mL micro-centrifuge tube (Sarstedt®, cat. no. 72.706), 2 mL screw-cap tube without skirted base (QIAGEN, cat. no. 990382), and 2 mL processing tube (QIAGEN, supplied with PAXgene Blood RNA Kit).
4. Make sure that you provide the correct volume as indicated on screen. The volume shown on the screen is the exact volume to be loaded. Do not overfill.
5. Place the opened tube in the worktable position as indicated in the table on the screen. It is important that you load the tube in the correct worktable position.
6. Place the lid of the microcentrifuge tube securely in the lid slot next to the tube.
7. Make sure that the required number of tips for each tip type is loaded as indicated on the screen. It is possible to use up the tip racks, if the minimum required number of tips for each type is loaded. However, it is recommended to load more than the minimum amount of tips.

**Note:** The loading position shown on the screen is the recommended position for tip racks. The position can also be changed. Later, when starting the run, the instrument will check if the correct tip racks are placed on the worktable and if there are enough tips for the protocol run.

There are 3 different kinds of tip racks that can be used on the QIAcube Connect MDx, depending on the protocol selected. A blue rack for 200 µL filter-tips, a light gray rack for 1000 µL filter-tips, and a dark gray rack for 1000 µL wide-bore filter-tips. The instrument uses the notches on the filter-tip rack to identify the type. To prevent mix up that may lead to a run issue, do not manually refill tip racks. Only use tips designed for use with the QIAcube Connect MDx.

**Important:** Do not use damaged filter tips. Do not load damaged tip racks onto the worktable.

**Note:** When using partially filled tip racks, consider the order in which the racks are loaded. The tip rack placed in position 1 will be used first.

8. Ensure that you empty the waste drawer containing used disposable labware before every run to avoid waste accumulation.
9. Press **Next** to proceed to loading the centrifuge or shaker rack, depending on the selected protocol.

**Important:** After loading, the robotic arm will automatically move back to its original position (above tip rack position 3). Ensure to stand clear of the instrument while the robotic arm is moving. Wait until the robotic arm has completed its movements.

#### 5.4.8. Loading the centrifuge


The Loading centrifuge screen guides you through setting up the required rotor adapters and centrifuge for the run. Make sure to read all essential and critical information before proceeding to the next step.

Depending on the selected protocol, loading of the centrifuge might not be required. In this case, the software will indicate that this step can be skipped.

Rotor adapters can be placed into a rotor adapter holder, enabling convenient and easy preparation and loading of columns. Place the columns, tubes, or samples into the appropriate positions in each rotor adapter as instructed by the software. It is recommended to label the elution tubes appropriately (e.g., with the corresponding sample ID).

For some protocols (e.g., PAXgene Blood RNA Kit), the software might instruct you to cut off the lid of a special spin column – the QIAshredder (pink) – for the middle position of the rotor adapter. Do this before loading the spin column (QIAshredder).

Be sure that the lid is completely removed from the spin column. Spin columns with partially removed lids may not be gripped properly by the robotic gripper and may cause the protocol run to crash.



**Comparing correct and incorrect removed column lids.**

Ensure that the tubes and spin columns have been pushed firmly into the appropriate rotor adapter position.

Place the lids into the correct lid position of the rotor adapter as indicated on your screen in the Lid position table column and the rotor adapter illustration. Make sure that the lids are pushed all the way down to the bottom of the slots on the sides of the rotor adapter. Incorrectly positioned lids can break off during centrifugation and cause the protocol run to crash.

**A**



1.5 mL microcentrifuge tube lid is in correct position

**B**



**Correctly loaded rotor adapter.** **A** Rotor adapter is correctly loaded, and the 1.5 mL microcentrifuge tube lid is in the correct position. **B** Correctly loaded rotor adapter seen from the side.

**C**



1.5 mL microcentrifuge tube lid not pushed down into slot

**D**



**Incorrectly loaded rotor adapter.** **C** Rotor adapter is incorrectly loaded with a 1.5 mL microcentrifuge tube. The lid of the tube is not pushed all the way down to the bottom of the slot of the rotor adapter and could break off during centrifugation (compare with part A of the above figure); **D** Incorrectly loaded rotor adapter seen from the side (compare with part B of the above figure).

E



1.5 mL microcentrifuge tube lid is in the wrong slot of the rotor adapter

**Rotor adapter is incorrectly loaded with a 1.5 mL microcentrifuge tube.** The lid of the tube is positioned in the wrong slot of the rotor adapter. During column transfer, the lid of the spin column could crash onto the lid of the 1.5 mL microcentrifuge tube, causing the protocol run to crash.

**WARNING Risk of personal injury and material damage**



To prevent plastic crash, load the tubes properly. After a plastic crash, sharp plastic particles could be inside the centrifuge. Be careful when handling items inside the centrifuge.

**WARNING Risk of personal injury and material damage**



Do not use damaged rotor adapters. The rotor adapters can only be used once. High  $g$  forces exerted in the centrifuge can cause damage to re-used rotor adapters. To prevent plastic crash, load the tubes properly. After a plastic crash, sharp plastic particles could be inside the centrifuge. Be careful when handling items inside the centrifuge.

**WARNING Risk of personal injury and material damage**



Make sure that lids from spin columns and 1.5 mL microcentrifuge tubes are in the correct position and pushed all the way down to the bottom of the slots on the sides of the rotor adapter. Incorrectly positioned lids can break off during centrifugation. Do not use damaged rotor adapters. The rotor adapters can only be used once. High  $g$  forces exerted in the centrifuge can cause damage to used rotor adapters.

**WARNING Risk of personal injury and material damage**



Be sure the lid is completely removed from the spin column. Spin columns with partially removed lids may not be removed properly from the rotor, causing the protocol run to crash. Make sure that lids from spin columns and 1.5 mL microcentrifuge tubes are in the correct position and pushed all the way down to the bottom of the slots on the sides of the rotor adapter. Incorrectly positioned lids can break off during centrifugation.

Depending on the selected protocol, samples can be loaded onto the shaker or directly into the centrifuge. Follow the instructions below the screenshot corresponding to your screen. Your screen may look different depending on the selected protocol.



## Loading the centrifuge in case samples are loaded onto shaker

This section describes a workflow that includes the shaker unit (e.g., for lysis). Tubes with samples need to be loaded onto the shaker unit (refer to Section 5.4.9 Loading the shaker), and the centrifuge must be prepared according the following description.

Pos.	Labware	Lid position
1	QIAamp spin column	L1
2	-	-
3	1.5 ml Elution tube	L3

6 x Rotor adapter

QIAamp DNA Mini Standard Step 8 of 9

Label Elution Tubes with Sample ID or Rotor Position

Cancel Back Next

May 23, 2024, 13:25 User: Admin Admin

**Load the centrifuge rotor adapter screen when samples are loaded onto the shaker.** Position 2 of the rotor adapter is empty.

The number and tube positions of rotor adapters required for the protocol run will be shown in the table and illustration on your screen. The table shows how to load and position each rotor adapter. The Pos. column indicates the position in the rotor adapter and the Lid position column indicates where to place the lid of a particular tube. These positions are outlined in the rotor adapter illustration below the table.

For each rotor adapter:

1. Load each tube/spin column in its correct position as indicated in the table on the screen. Touch the table row to highlight the particular tube position in the illustration below the table.
2. Ensure that the tubes and spin columns are pushed firmly into the appropriate rotor adapter position.
3. Make sure that the lids are pushed all the way down to the bottom of the slots on the sides of the rotor adapter. Ensure that you place the lids in the correct lid positions.
4. Label the elution tubes appropriately (e.g., with the same sample ID as the sample on the respective input position on the shaker (see Section 5.4.9 Loading the shaker) or the rotor position number). Be sure to use a sticky label, which is safely attached.
5. If required and described in the table, cut off the lid, and break off the bottom of the spin column.
6. Repeat steps 1–5 until all rotor adapters have been prepared.
7. Place the loaded rotor adapters into the centrifuge buckets as shown on the right side of the screen. For ease of use and high process safety, the rotor adapters only fit into the centrifuge buckets in one orientation.
8. Press **Next** to proceed with loading the samples into the shaker. Follow the instructions in Section 5.4.9 Loading the shaker. Depending on the selected protocol, the sequence of the next screens can be different.
9. If samples are loaded into the shaker, the next section can be skipped.

## Loading samples into the centrifuge

This section describes a workflow that does not include the shaker unit (e.g., for lysis). Samples are directly loaded into the centrifuge.

The procedures for loading samples into the centrifuge are shown on both sides of the screen. Such protocols are only available in the Research mode of the software.

Prepare the required number of rotor adapters as shown on the screen.

**Load the centrifuge rotor adapter**

Pos.	Labware	Lid position
1	MB RNA spin Column	L1
2	450 µl sample	-
3	1.5 ml Elution tube	L3

6 x Rotor adapter

Diagram of rotor adapter showing positions 1, 2, and 3 with lid positions L1, L2, and L3.

**RNeasy PowerMicrobiome**  
IRT with DNase Step 9a of 9

Rotor Adapter Position for sample loading: 2

Lyse samples mechanically according to RNeasy PowerMicrobiome Kit instruction manual. After centrifugation in step 4, Transfer at least 450 µl of sample into Rotor Adapter position 2 (Rotor Adapter Middle Position). If the sample volume is less than 450 µl, add Solution PM1 up to the final volume. Important: It is critical that you use Solution PM1 and not water to adjust the volume. Solution PM1 is available as an accessory product (cat. no. 26000-50-1).

Buttons: Cancel, Back, Next

May 23, 2024, 13:28 User: Admin Admin

**Load the centrifuge rotor adapter screen when samples are loaded into the centrifuge.** Samples are loaded in position 2 of the rotor adapter.

The number and tube positions of the rotor adapters required for the run are shown in the table and in the illustration. The table shows how to load each rotor adapter. The Pos. column indicates the position in the rotor adapter, and the Lid position column indicates where to place the lid of a particular tube.

For each rotor adapter:


1. For samples: Prepare and load the samples as outlined on the screen. Make sure to load the correct sample amount. Ensure that you read all essential and critical information indicated in the blue **Information** box on the right side of the screen.
2. Load each tube/spin column in its correct position as indicated in the table on the screen. Touch the table row to highlight the particular tube position in the illustration below the table.
3. Ensure that the tubes and spin columns are pushed firmly into the appropriate rotor adapter position.
4. Make sure that the lids are pushed all the way down to the bottom of the slots on the sides of the rotor adapter. Ensure that you place the lids in the correct lid positions.

**Note:** If required and described in the table, cut off the lid, and break off the bottom of the spin column.

5. Repeat steps 1–4 until all rotor adapters have been prepared.
6. Press **Next** to proceed loading of rotor adapters into centrifuge. Load the rotor adapters into centrifuge. Place the prepared rotor adapters into the centrifuge buckets as shown on the right side of screen. For ease of use and safety, the

rotor adapters fit in the centrifuge buckets in only one orientation. To prevent sample mix up, ensure to load the particular Sample ID into the defined centrifuge position.

- 7. If required, change the default value in the Sample ID fields using the on-screen keyboard. You can enter the value manually or scan the sample barcode using the external barcode scanner. The sample ID is initially created using the YYYYMMDD-HHMM-no. format. Be sure that the same ID can be found on the respective elution tube on a sticky label, which is safely attached.



[View sample details screen.](#)

**Note:** Sample IDs are part of the run reports and can be part of log files and audit trail. These IDs are not encrypted.

**Important:** Please note that the sample ID field must not contain any personal data.

- 8. Press **Next** to proceed to start the run.

### 5.4.9. Loading the shaker

The Load shaker screen guides you through the loading of the shaker.

Depending on the selected protocol, loading of the shaker might not be required. In this case, the software will indicate that this step can be skipped. According to the requirements of the selected protocol, samples and/or other tubes must be loaded on the shaker.

**Load shaker**

Pos	Sample ID	Tube Type	Value
01	20240610-0224-01	2 ml screw-cap	100 µl
02	20240610-0224-02	2 ml screw-cap	100 µl
07	20240610-0224-07	2 ml screw-cap	100 µl
08	20240610-0224-08	2 ml screw-cap	100 µl

Cancel

**DNeasy Blood and Tissue**  
Standard Step 9 of 9

01 07  
02 08

S2

⚠ Use shaker type S2

i Sample information

◀ Back Next ▶

Jun 10, 2024, 14:25 Mode: Research User: Admin Admin

#### Loading of shaker; lid positions filled with shaker rack plugs.

In this step, the software displays the shaker positions, tubes, and volume to be loaded in the table and in the schematic on the right. Make sure to load the correct shaker rack type as described on the right side of the screen. The shaker adapter can only be loaded in the correct orientation. Make sure to read all essential and critical information stated under Sample information before proceeding to the next step.

1. Ensure that you are using the correct shaker rack type indicated by **S2** in the shaker illustration
2. If required, change the default Sample IDs in the respective fields using the on-screen keyboard. You can enter the value manually or scan the sample barcode using a barcode scanner. The sample ID is initially created using the YYYYMMDD-HHMM-no. format.

**Note:** Sample IDs are part of the run reports and can be part of log files and audit trail. They are not encrypted.

**Important:** Please note that the sample ID field must not contain any personal data.

**Note:** Sample IDs will be lost if the user moves back to kit selection screen in order to change the selection

3. Prepare the correct tubes. Ensure that you read all essential and critical information indicated in the **Information** icon (i) marked in red. Information about the labware to use can also be found in the respective kit handbook. If you use sticky labels on the tubes, ensure to use a thin label that allows to insert the tubes completely into the shaker position.
4. Load the tubes into the positions on the shaker rack according to the allocation of the sample IDs as displayed by the touchscreen. The shaker rack positions are numbered for easy identification. Touch the table row to highlight the position in the schematic at the right.

- Depending on the tube type, a shaker rack plug or the tube lid must be placed into the slot next to the tube, as shown on the screen or indicated in the **Information** icon (i). Make sure that the lid/shaker rack plug is securely seated in the slot. Do not place a lid or shaker rack plug next to an empty shaker rack position.

**Note:** Depending on the protocol used, positions 1 and 7 may be used differently from the other positions. Ensure that you follow the instructions in the table and in the schematic to load the shaker correctly. In this exemplary screen below, no lid or shaker rack plug is required for these positions.

**Load shaker**

Pos	Sample ID	Tube Type	Value
01	Empty tube required for run	2 ml safe-lock... Cut off lid	Empty tube
02	sample 1	2 ml safe-lock...	350 µl
03	20240523-0133-03	2 ml safe-lock...	350 µl
04	20240523-0133-04	2 ml safe-lock...	350 µl
08	20240523-0133-08	2 ml safe-lock...	350 µl
09	20240523-0133-09	2 ml safe-lock...	350 µl
10	20240523-0133-10	2 ml	350 µl

**AIIPrep DNA RNA Micro**  
Standard part A Step 9 of 9

Use shaker type 2


Sample information ...

Back Next

May 23, 2024, 13:34 User: Admin Admin

**Example of protocol using shaker positions 1 and 7 differently. For this example, do not load lids or shaker rack plugs for these positions.**


The pictures below indicate how the shaker loading illustrated on the user interface (left-hand side) translates into the real shaker layout (right hand side).




Loading the shaker rack with sample tubes that have attached lids.



Lids of sample tubes must be securely placed into the slots at the edge of the shaker rack.



Loading the shaker rack with sample tubes that have screw-cap lids.



Shaker rack plugs must be placed into the slots at the edge of the shaker rack.

6. Press **Next** to proceed to starting the run or loading the centrifuge, depending on your selected protocol.

**WARNING** Risk of personal injury and material damage




To prevent plastic crash, load the tubes properly. After a plastic crash, sharp plastic particles could be inside the centrifuge. Be careful when handling items inside the centrifuge.

**Important:** Do not use 1.5 mL microcentrifuge tubes on the shaker. These microcentrifuge tubes can cause filter-tips to stick during sample transfer. Using this tube on the shaker can damage the pipetting system and can cause a centrifuge crash.

## 5.5. Starting a protocol run

A confirmation message will appear when the final step on the last setup screen has been completed.

1. In this window the user can add a comment to the run. The comment will be part of the run report.




2. Close the hood to continue.

**Note:** Ensure that the waste drawer is empty and closed before starting the run.

3. Press **Start** to begin the run. If required, press **Back** to return to the previous setup screen.


**Important:** Do not try to open the instrument hood during a run.

**Important:** Once you set up a run on the instrument and have pressed the Start button, we strongly recommend to stay for a few minutes near the instrument until the load check has been completed. This will allow you to add missing reagents or consumables in case the instrument detects missing items. The status of the load check and the recommendation to wait until it is finished is also shown on the user interface.




**Information to stay near the instrument during the load check.**

If the load check fails, pop-up dialogs will show you the problem. If the user presses **OK** in the dialog, the system will return the last page of the loading workflow. This allows the user to check the loading and start the run again without losing the entered data.




If the load check is successful, the run starts immediately.

**Note:** For very first run of a certain application and sample number, the estimated run duration is not available. If the same application has been used before (with the same amount of samples), an approximate run time estimation is available.



**Run status screen during a protocol run.**


During the run, you can press the **Run Details** button to display the run steps and details. To return to the run view, press **Close**.



**Run details screen.**



4. When the protocol run is completed, the output position and content will be shown on the right side of the screen. For some protocols, further treatment of samples is described on the right side of the screen. Remove the eluates/samples from the instrument directly after the run is finished and make sure that you follow the proper procedures for storing and handling the eluates/samples.



**Run completed screen.**

5. Press **Done** to create the report file. The run report is a PDF file and contains the following information:


- the protocol information (name and the version of the completed application process file)
- the serial number of the instrument
- the software version
- the sample IDs and their positions
- the time, date, and user at the run start
- the time and date of the run end
- the user ID that confirmed the end of run
- the kit material number, lot number, and expiration date
- the error and warning descriptions
- the run validity (valid or invalid)
- the run status (completed or aborted)
- the maintenance status (due, last executed)
- the software mode (IVD or Research)
- the run ID
- the elution volume
- the final eluate position

**Important:** It is recommended to perform regular maintenance as described in Section 6.3 Regular maintenance before starting the next run.

**Note:** For data protection reasons only the user ID (not the user name) and sample IDs are displayed on the run report. Please ensure that these IDs do not contain any clear names if this is not according to your regulations.

## 5.6. Stopping a protocol run

If required, a run can be stopped by pressing the **Stop Run** button on the run status screen (see Section 5.5 Starting a protocol run). To confirm the run stopped, click **Yes** in the Stop Run dialog box.



### Stop run screen.

If a run is stopped, perform the daily maintenance as described in Section 6.4, Daily maintenance, and ensure that no plastic parts are present in the centrifuge before starting the next run. Additionally, it is recommended to reboot the system before starting the next run.


**Note:** If a protocol run is stopped, the run cannot be restarted. Under the **Run Details**, you will find the step at which the protocol stopped.


**Note:** The run will also stop immediately if the hood is opened during a run. Do not open the hood during a run.

## 5.7. Saving run reports to the USB flash drive

Run reports are saved on the instrument after each run is confirmed by clicking the **Done** button.

To transfer run reports to the USB flash drive, proceed as follows:

1. Press the **Tools** icon () on the menu bar.
2. Press the **Data Exchange** tab. The number of available run reports is shown on the screen.



**Data exchange screen.**

3. This screen offers the opportunity to upload a custom logo to be shown on the run reports. You can upload logos of your company or university in several picture formats by clicking the **Upload** button.
4. If not yet connected, connect the USB flash drive, which is delivered with the instrument, to one of the USB ports at the left of the touchscreen.

**Important:** Only use USB flash drive that is delivered with the instrument and ensure that the USB flash drive has enough remaining space prior to initiate the procedure of saving the run reports.

5. To save all available run reports to the USB flash drive, select **All reports**. To save only the recent reports, select **Recent reports**. To enter the number of reports to be saved, touch in the Recent reports field and enter the desired number.
6. If you would like to delete reports from the instrument after download, check the box "Delete reports after download".

**Important:** The deleted reports cannot be restored from the instrument. Ensure to store the files from the USB stick to a safe place.

Press **Download run report(s)** to save the reports to the USB flash drive. A confirmation message that run reports are successfully saved to the USB flash drive will appear. The USB flash drive can be removed from the instrument.

**Important:** Do not remove the USB flash drive while the files are downloading. Wait until the download is completed.

**Note:** Run reports downloaded to the USB flash drive contain unencrypted sample IDs and the User ID. Ensure that these IDs do not contain clear names of users and patients to comply to your local data protection regulations.

7. On the data exchange screen, you can also create and download a support package, which contains additional information like the logfiles. You can also check the available disk space and free disk space if required.

**Important:** If you use Delete Files, ensure to store the files from a previously created backup package at a safe place.

### 5.7.1. Low space on flash memory

The number of run reports on the instrument is limited. If the remaining space on the internal storage reaches 10% of the total capacity, you are prompted to perform a system backup. Follow this recommendation. When the system reaches a completely full disk, no more runs can be started.

System alert - Disk space is almost full

The available disk space has reached a critical limit (9999 KB). Message ID: 2700  
Step 1: Backup files by pressing 'Download files'. Message ID: 2704

Step 1: Download data to USB stick

Download files

Step 2: Check downloaded files on USB stick

Check done

Step 3: Delete files on disk to free up disk space

Delete files

USB

Connected

Cancel

OK

#### Warning on critical remaining disk space.

You will be guided through the 3 steps of downloading a backup, checking the backup for completeness, and freeing up disk space by deleting system files. In the end, the following screen will be displayed: persistent files deleted. Press **OK** to leave the message box and return to normal operation.

System is OK - Normal operation

Persistent files deleted. Available disk space is now 51000 KB.  
Press 'OK' to leave the message box and return to normal operation.  
Message ID: 2053

Step 1: Download data to USB stick

Download files

Step 2: Check downloaded files on USB stick

Check done

Step 3: Delete files on disk to free up disk space

Delete files

USB

Connected

OK

On a regular basis, you can also check the free disk space by clicking on the **Configuration** icon located on the top of the screen. At the right-hand side, system information is displayed showing the free disk space.

The screenshot shows the QIAGEN configuration interface. At the top, there is a navigation bar with icons for Setup, Tools, Configuration, QIAsphere, Network, Lock, and Logout. Below this is a sub-menu with tabs for System, Users, Protocols, Settings, LAN, Wi-Fi, and QIAsphere. The 'System' tab is selected, displaying the following information:

- Device name:** IVDNew
- Date:** Nov 16, 2023
- Time:** 11 : 42
- Language setting:** American English (with a dropdown arrow and a 'Load from USB' button)
- Serial number:** 13
- Free disk space:** 48096 KB
- Software version:** 2.0.0.0 CI 66dceb2
- USB:** Connected

At the bottom of the system information panel, there are two buttons: 'Reset Centrifuge Counter' and 'Update Software'. A footer bar at the very bottom shows the date 'Nov 16, 2023, 11:42' and the user 'User: Admin Admin'.

#### Free disk space information.

If you do not follow the recommendations above, the system memory will be full after a few more runs. If the flash memory disk is full, you cannot start a run. The system will show an alert. You will be guided through the 3 steps of downloading a backup, checking the backup for completeness, and freeing up disk space by deleting system files. In the end the following screen will be displayed: persistent files deleted. Press **OK** to leave the message box and return to normal operation.

**Note:** the backup package contains sensitive user data. Please ensure data safety according to your local regulations.

The screenshot shows a system alert dialog box with a red header bar containing an exclamation mark icon and the text 'System alert - Disk space is full'. The main content area contains the following text:

There is not enough disk space available to continue instrument operation. Free disk space now!

Step 1: Backup disk files by pressing 'Download files'.

Below this text is a table with three steps and their corresponding actions:

Step 1: Download data to USB stick	Download files
Step 2: Check downloaded files on USB stick	Check done
Step 3: Delete files on disk to free up disk space	Delete files

At the bottom right of the dialog box, there is a 'USB Connected' status indicator and an 'OK' button.

#### Full disk space system alert.

## 5.8. Independent heater/shaker operation

The heater/shaker can be operated individually if the QIAcube Connect MDx is not running a protocol. The heating and shaking functions are not interlinked and can be used independently or in combination.

### **WARNING** Risk of personal injury and material damage




Do not attempt to move the QIAcube Connect MDx during operation.

### **WARNING** Hot surface



The shaker can reach temperatures of up to 70°C (158°F). Avoid touching it when it is hot, in particular, shortly after a run has been carried out.

1. Press the **Tools** icon (  ) on the menu bar.
2. Press the **Run Modules** tab.
3. Press the **Heater Shaker** tab.




**Heater Shaker operation screen.**

4. Press the appropriate field to select the Frequency, Temperature, and Duration using the on-screen keyboard.
5. Load the shaker rack tubes containing samples.
6. Close the hood to start run. Press **Start**.

**Note:** The remaining time and current temperature and status of the heater are displayed on the screen. Wait until the operation is completed. The ongoing operation is indicated by a moving circle.

7. To stop the run progress, press **Stop**.



Heater Shaker operation screen.

## 5.9. Independent centrifuge operation

The centrifuge can be operated individually if the QIAcube Connect MDx is not running a protocol.

Do not attempt to move the QIAcube Connect MDx during operation.

### CAUTION Damage to the instrument



The QIAcube Connect MDx must not be used if the centrifuge lid is broken, or if the lid lock is damaged. Make sure that no loose material is inside the centrifuge during operation.

Make sure that the rotor is installed correctly and that all buckets are properly mounted, regardless of the number of samples to be processed. Load the rotor only as instructed by the software.

Only use rotors, buckets, and consumables designed for use with the QIAcube Connect MDx. Damage caused by use of other consumables will void your warranty.

We recommend replacing the centrifuge rotor and buckets after 20,000 cycles, which is equivalent to 9 years of usage with two runs per day for 220 days each year. For more information contact QIAGEN Technical Services.

### WARNING Moving parts



In case of breakdown caused by power failure, remove the power cord and wait 10 minutes before attempting to manually open the centrifuge lid.

### CAUTION Damage to the instrument



After a power failure, do not move the z-module (robotic arm) manually in front of the instrument. Damage may occur if the QIAcube Connect MDx hood is closed and collides with the z-module.

### WARNING Risk of overheating




To ensure proper ventilation, maintain a minimum clearance of 10 cm at the sides and rear of the QIAcube Connect MDx.


Slits and openings that ensure the ventilation of the instrument must not be covered.

**WARNING Risk of personal injury and material damage**



Raise the centrifuge lid carefully. The lid is heavy and may cause injury if it falls.

1. Press the **Tools** icon (  ) on the menu bar.
2. Press the **Run Modules** tab.
3. Press the **Centrifuge** tab.












**Centrifuge operation screen.**

4. Press the appropriate field to select the **Speed** and **Duration** using the on-screen keyboard.
5. If the centrifuge lid is not open, press **Open Centrifuge Lid**.
6. If required, load opened 1.5 mL microcentrifuge elution tubes and/or QIAGEN spin columns into the rotor adapters and place the lids into the appropriate slots in the rotor adapter.
7. Ensure that the tubes and spin columns are pushed firmly into the appropriate rotor adapter position.
8. Make sure that the lids are pushed all the way down to the bottom of the slots on the sides of the rotor adapter. If required, cut off the lid.
9. Place the rotor adapters into the centrifuge.
- Important:** If fewer than 12 samples are to be processed, ensure that you load the correct centrifuge positions as described in the Loading scheme table below. One or 11 samples cannot be loaded.
10. Close the hood and press **Start** to start centrifugation.

**Note:** The **Close centrifuge lid** button is not needed to start a centrifuge run, as the lid will close automatically. It is only required in case you need to prepare the QIAcube Connect MDx for shipment.



Centrifuge loading scheme			
No. of samples	Centrifuge loading scheme	No. of samples	Centrifuge loading scheme
2		7	
3		8	
4		9	
5		10	
6		12	Load all positions

## 5.10. Managing protocols

Commonly used QIAGEN standard protocols are installed on the QIAcube Connect MDx upon delivery. The range of QIAGEN standard protocols is continuously expanding, and these protocols are downloadable at no cost. See the **Resources** tab at [www.qiagen.com/QIAcube-Connect-MDx](http://www.qiagen.com/QIAcube-Connect-MDx) for the protocols running in the IVD mode of the software. For protocols running in the Research mode of the software please refer to [www.qiagen.com/QIAcube-Connect](http://www.qiagen.com/QIAcube-Connect). QIAGEN's Application Lab Specialists can also customize these protocols or develop new protocols depending on your needs. Customized protocols can only be used in the Research mode of the software and are not validated and must not be used for diagnostic purposes. Protocols that are no longer required can be removed from the QIAcube Connect MDx. Protocols can only be managed by users assigned the Administrator role.

### 5.10.1. Installing new protocols via USB stick

This process is used to install new protocols and translated protocols if required by the language settings – see Section 4.5.1 System configurations – or to reinstall a protocol backup.

1. On a computer running Microsoft Windows, download the new protocols from the **Resources** tab at [www.qiagen.com/QIAcube-Connect-MDx](http://www.qiagen.com/QIAcube-Connect-MDx). A previously created protocol backup package is located in the folder **Protocol\_Download** on the used USB stick.

Use the USB flash drive that was shipped with the QIAcube Connect MDx to transfer the protocol files to the instrument.


**Note:** Checksum confirmation is required to secure software integrity after web download was successfully completed and before subsequent handling of the software. For detailed information on confirmation of software integrity during download and file transfer, please check the "QIAGEN software integrity verification process" description document, which is provided on the same download page as the protocol package.

2. Unzip the download folder. This results in separated **.zip** files for each protocol or in a folder called **Protocol\_Upload** containing the individual **.zip** files.
3. Create a new folder on the USB flash drive with the name **Protocol\_Upload** and copy the downloaded protocol zip file(s) or the protocol zip file(s) from the folder **Protocol\_Download** to this directory. The folder **Protocol\_Download** is created by the instrument during a protocol backup according to Section 5.10.4 Saving protocols. Do not rename, unzip, or modify the individual protocol files. Otherwise, they cannot be used. If the download package already contained the folder **Protocol\_Upload**, simply copy it to the main directory of the USB flash drive. Always make sure to use the correct directory (folder **Protocol\_Upload** in the main directory of the USB flash drive); otherwise, QIAcube Connect MDx will not find the protocols.

**Note:** Do not rename or modify the protocol files. Otherwise, they cannot be used.


4. Connect the USB flash drive to the QIAcube Connect MDx using one of the USB ports at the left of the touchscreen.
5. Select the **Configuration** icon ()

- Press the **Protocols** tab.



#### Protocols configuration screen.

- Press **Copy from USB**.
- A message is displayed indicating how many protocols are found on the USB flash drive and how many will be installed or overwritten.




#### Screen for doubled protocols.

**Note:** Already installed protocols with the same name but an older version number will be overwritten. The older version can only be re-installed if the newer version is deleted before (see Section 5.10.3).


- Press **OK** to start uploading.

All compatible protocol zip file(s) in the **Protocol\_Upload** folder will be installed. Protocols may be incompatible if they are not available on the device type (QIAcube Connect vs. QIAcube Connect MDx) or not released for the particular instrument serial number.

10. During transfer the progress is visualized by a moving bar.




11. Wait until the transfer is complete. A message is displayed when the transfer is complete.
12. If duplicate or incompatible protocols have been detected the message box outlines the overwritten protocols. Even though the message box title is "Upload failure", only the protocols outlined in the message box have failed. All other protocols from the package are ready for use.




New protocols will be available immediately after upload.

### 5.10.2. Installing protocols via QIAsphere


This process is used to install new protocols and translated protocols via QIAsphere. If a new protocol package is available, QIAsphere will send a notification to your device. Protocol packages can be created and sent to the instrument by an administrator user using the QIAsphere App. For details, refer to the *QIAsphere User Manual*. The notification appears under the QIAsphere button and is highlighted by a red dot.



1. Press the **Info** button (i) .
2. The available protocol package is described in the message box. If you wish to download the specified protocol package press **Download**.




3. The following dialog appears. Confirm the protocol download with **OK**.




4. After successful download a red dot on the QIAsphere icon will indicate a new notification.

5. Check the QIAsphere notifications, and press the **Info** button (i) again.




6. If you are logged in with an administrator role, the following dialog box specifying the package version appears. Confirm the protocol update with **Yes**.


**Note:** If you confirm with **Yes**, all existing protocols will be overwritten. Only the protocols available in the new package will be installed on the instrument.




7. After checking the content of the upload package, another message box is displayed. This message box outlines whether protocols will be removed or overwritten and which protocol is affected.



8. The progress of the installation will be shown by a moving bar in the QIAsphere Notification Center.



9. After successful installation a message box with a summary of installed, removed, or overwritten protocols appear.






### 5.10.3. Deleting protocols


**Important:** Before deleting, you should back up the protocols on the USB flash drive provided with the instrument. See Section 5.10.4 Saving protocols.

1. Select the **Configuration** icon (⚙️).
2. Press the **Protocols** tab.



**Protocols configuration screen.**


3. To delete all protocols installed on the instrument, press **Delete All**. To delete selected protocols press **Delete selected**.
4. Click + to extend the list of protocol folder.



5. Click + to extend the list of protocols.


6. Select single or multiple protocol(s), and click **Delete** to remove them at once.

**Note:** Deletion of parent folder is not possible.



#### Deletion of selected protocols.

7. A message box with a summary of the deletion appears. Confirm by pressing **Yes** if you wish to proceed.




### 5.10.4. Saving protocols


Protocols can be downloaded from the instrument to USB flash drive to transfer them to another instrument or to save them for backup purposes, for example, before a software update. Use the USB flash drive provided by QIAGEN.

1. Connect the USB flash drive that was shipped together with the instrument to the QIAcube Connect MDx using one of the USB ports at the left of the touchscreen.
2. Select the **Configuration** icon (⚙️).
3. Press the **Protocols** tab.


4. From the Download Protocols section, press **Copy to USB**.



5. A message box appears indicating how many protocols have been found on the device.
6. Confirm by pressing **Yes**.



7. The successful download will be confirmed by a message box. Wait for this confirmation as the process can take a few minutes. Confirm with **OK**.



## 5.11. User Management

The QIAcube Connect MDx is provided with a User Management feature. This feature enables you to set up multiple users with two different roles: administrator and operator. For each operator, the software mode (IVD or Research) to be used can be set up. Access to both software modes can be selected for an operator, or a restricted access to only one software mode. When you use the QIAcube Connect MDx for the first time, a default user named Admin is already pre-installed and configured with both roles assigned. The user management feature is only available for users assigned the Administrator role.

### 5.11.1. Setting up a new user

1. Press the **Configuration** icon (⚙️) on the menu bar.
2. Press the **Users** tab.

The configured users are shown in the table. Each row contains the data for one user.

User ID	First Name	Last Name	Role(s)	Edit
Admin	Admin	Admin	Administrator, Operator	
RG	Ralf	Secret	Operator	

☒ Show only activated user profiles New ...

Jun 04, 2024, 10:02 User: Admin Admin

**List of configured users in the user management.**

**Note:** It is recommended to create at least two users with the administrator role.

3. Press **New** to add a new user.
4. Enter the respective data for the new user. Keep the "Activate User" box checked.

The User ID, First name, and Last name fields are mandatory. These fields may contain up to 30 letters and numerical characters. The user ID must be unique for each user profile. It must contain at least one letter and cannot contain blank spaces. The user ID is used for logging in and is printed on run reports. The first and last name are displayed on the touchscreen for the currently logged in user.

The Password field is mandatory and must contain 8–40 letters or numerical characters. Enter the same password into the Confirm password field.

Select the user role: **Administrator** and/or **Operator**. The operator may only use the instrument, while the administrator is only allowed to configure the system. One user can have both roles assigned at the same time. This is the recommended setting for an Administrator who wants to start application runs as well. The default user Admin has both user roles assigned. Select the software mode (Research and/or IVD) the user shall have access to.

The E-mail address field is optional. The system does not confirm whether the e-mail address entered is valid.

**Important:** A newly created user with administrator rights can only configure the system and cannot start a run. If this is required both roles have to be chosen.

5. Press **OK** to save the new user.

5.11.2. Changing data for an existing user

- 1. As an admin user, press the **Configuration** icon (⚙️) on the menu bar.
- 2. Press the **Users** tab.

The configured users are shown in the table. Each row contains the data for one user.

QIAGEN

Setup

Tools

Configuration

QIASphere

Network

Lock

Logout

System

Users

Protocols

Settings



LAN

Wi-Fi

QIASphere

User List

User Config

User Id	First Name	Last Name	Role(s)	Edit
Admin	Admin	Admin	Administrator, Operator	
RG	Ralf	Secret	Operator	

☒ Show only activated user profiles

New ...

Jun 04, 2024, 19:02

User: Admin Admin

List of configured users in the user management.

- 3. In the user profile row, press the **Edit** icon ().
- 4. A screen will appear showing the current information of the user. Edit the information as necessary.

»» ⚙️ 🔧 ⚙️ Simulator

Last Login: 2023-02-20 Next change: 58 days

Anonymous user ID

First name

Last name

User\_1

Admin

Admin

E-mail

n/a

Enter password

••••••••

Confirm password

••••••••

☒ Administrator

☒ Operator

☒ Research Mode

☒ IVD Mode

☒ Activate User

☐ Change Password

Cancel

OK

Jun 02, 2024, 12:57

Guest (ID)

User: Admin Admin

Edit User screen.

5. Users with the Administrator role are allowed to change or reset the passwords of all other users, including other admins. We recommend creating at least one additional administrator as a backup for the preinstalled administrator user Admin. Passwords are never displayed in this process, so the administrator cannot view passwords.

If you touch the password field, the existing password will be cleared, and a new password must be entered and confirmed.

6. To confirm the changes, press **OK**. To close the dialog and discard the changes, press **Cancel**.
7. The administrator can also change the user configuration on the Users tab. The administrator can set a number of login attempts, the number of days between password changes (Note: Setting to 0 leads to daily password change), and the number of minutes before automatic logout.

**Note:** The input range to define the number of login attempts is 2 to 10.


**Note:** If the number of minutes before automatic logout is set to 0, automatic logoff is disabled.


The screenshot displays the QIAGEN Connect MDx configuration interface. At the top, there is a navigation bar with icons for Setup, Tools, Configuration (selected), QIAsphere, Network, Lock, and Logout. Below this is a secondary navigation bar with tabs for System, Users (selected), Protocols, Settings, LAN, Wi-Fi, and QIAsphere. Under the Users tab, there are two sub-tabs: User List and User Config (selected). The main content area is titled "Settings for all users" and contains three input fields with labels and descriptions:


Input Field	Label	Description
10	[-]	Number of login attempts before user is locked
360	Days	Number of days between password changes
0	Minutes	Number of minutes before user logout (0 = no forced logout)

At the bottom of the screen, there is a status bar showing the date and time "Jun 28, 2023, 09:36", the mode "Mode: Research", and the user "User: Admin Admin".


### 5.11.3. Temporarily deactivating and reactivating a user

1. To temporarily deactivate a user, press the **Edit** icon () in the user profile row. Uncheck the Activate User box. It is not possible to deactivate the currently logged-in administrator.
2. To reactivate a user profile, deactivated users can be displayed in the users list by unchecking the box "Show only activated user profiles".



3. Press the **Edit** icon () in the user profile row. If required, change the password of the user. Check the Activate user box.

**Note:** If a user tries to login with the wrong password, the user profile will be automatically deactivated after the set number of failed login attempts.




The screenshot shows the "Edit User" dialog box. It contains fields for Anonymous user ID (User\_1), First name (Admin), and Last name (Admin). There is also an E-mail field with the value n/a. Below these fields are two password fields: Enter password and Confirm password, both containing masked characters. To the right of the password fields are four checkboxes: Administrator (checked), Operator (checked), Research Mode (checked), and IVD Mode (checked). At the bottom left, there are two checkboxes: Activate User (checked) and Change Password (unchecked). At the bottom right, there are two buttons: Cancel and OK. The top right corner of the dialog box shows "Last Login: 2024-06-10" and "Next change: 22 days".

### 5.11.4. Resetting a user password


If a user tries to log in with a wrong password more times than the set number of failed login attempts, the user will be inactivated. In this case the user can be re-activated by another administrator according to Section 5.11.3.

If no alternative user with administrator role is available, open the Help Center (**Help** button at upper right-hand side) and follow the screen instructions. The Help Center only works if the user name of an administrator (default: Admin) is known.



The screenshot shows the top of a web application with the QIAGEN logo on the left and a 'Help' button (a question mark icon) on the right. Below the header is a light blue box with the text 'Please log in to start'. Inside this box are two input fields: 'User ID' with 'Admin' entered, and 'Password'. Below these are two radio buttons: 'Research Mode' (selected) and 'IVD Mode'. At the bottom of the box is a dark blue 'Log in' button.

Login screen with Help button.



The screenshot shows the 'Help' screen. It has the QIAGEN logo at the top left. The main text reads: 'Another administrator can reset your password. If you are the only administrator, contact QIAGEN Technical Support and provide the generated challenge code and serial number. You will receive a one-time password.' Below this text is a button labeled 'Generate challenge code'. Underneath is the text 'Challenge code:' followed by a hyphen. Below that is a button labeled 'Show details'. Then is the text 'One-time password' followed by an empty input field. At the bottom are two buttons: 'Back' and 'Next'.

Note your instrument serial number and the challenge code generated in the Help Center and contact QIAGEN Technical Services. Ensure you have an email address available which is already known to QIAGEN Technical Services. You will receive a one-time password.




### 5.11.5. Changing password

Users with Administrator role are allowed to change the password for every user by editing the user profile. Refer to Section 5.11.2, Changing data for an existing user, for more details. Passwords are never displayed in this process, so the administrator cannot view password.

Users with Operator role can change their own password according to the following instructions:

1. Press the **Configuration** icon (⚙️) on the menu bar.
2. For users with the role Operator, the **Password** tab is automatically active.



The screenshot shows the QIAsphere Configuration screen. The top menu bar includes icons for Setup, Tools, Configuration (active), Simulator, QIAsphere, Network, Lock, and Logout. Below the menu bar, there are tabs for System, Password (active), Protocols, Settings, LAN, Wi-Fi, and QIAsphere. The Password tab contains three input fields: 'Old password', 'New password', and 'Confirm new password'. At the bottom of the form are 'OK' and 'Cancel' buttons. The status bar at the bottom shows the date and time 'Sep 18, 2023, 12:28' and the user 'User: John Doe'.

#### Change password screen.

3. Enter the old password into the Old password field. Touch the field to open the on-screen keyboard.
4. Enter a new password into the New password field, and re-enter the new password in the Confirm new password field.

**Note:** The new password must be different from the last three used passwords.

Press **OK** to save the new password. Press **Cancel** to discard any changes and to keep the old password.

To return to the setup screen, press the **Setup** icon (⏪).

## 6. Cleaning and Maintenance

### **WARNING/ CAUTION**



#### **Risk of personal injury and material damage**

Only perform maintenance that is specifically described in this user manual.

The following maintenance procedures must be carried out to ensure reliable operation of the QIAcube Connect MDx:

- Regular maintenance: after each protocol run
- Daily maintenance: after the last protocol run of the day and after switching from Research software mode to IVD mode
- Monthly maintenance: every month
- Periodic maintenance: when necessary; at least every 6 months
- Annual (preventive) maintenance performed by QIAGEN-authorized Service Specialists (for more information contact QIAGEN Technical Services)

Optionally, these procedures may be performed to check and ensure the reliability of operation of the QIAcube Connect MDx:

- UV Run: reduces contaminants (e.g., nucleic acids and *E. coli*)
- Tightness Test: ensures the tightness of the tip adapter (e.g., after O-Ring change)

The software provides step-by-step guidance under **Tools/Maintenance** for the maintenance procedures listed above, except for the regular maintenance.

Following these procedures ensures that the QIAcube Connect MDx is free from dust and liquid spills.

Select the cleaning agent according to the objective of the cleaning procedure, the sample material used and the downstream assay.

### **WARNING**

#### **Risk of fire or explosion**



When using ethanol or ethanol-based liquids on the QIAcube Connect MDx, handle such liquids carefully and in accordance with the required safety regulations. If liquid has been spilled, wipe it off and leave the QIAcube Connect MDx hood open to allow flammable vapors to disperse.

Before using any cleaning or decontamination methods except those recommended by the manufacturer, users should check with the manufacturer that the proposed method will not damage the equipment.

## 6.1. Cleaning agents

The following disinfectants and detergents are recommended for cleaning the QIAcube Connect MDx.

**Note:** If you want to use disinfectants different from those recommended, ensure that their compositions are like those described below.

General cleaning of the QIAcube Connect MDx:

- Mild detergents (e.g., Mikrozyd® AF sensitive)
- 70% ethanol (for cleaning the worktable only; not for cleaning the QIAcube Connect MDx hood)

## 6.2. Decontaminating the QIAcube Connect MDx surface

Ethanol-based disinfectants can be used for disinfection of surfaces, such as the worktable or inside the centrifuge: for example, 25 g ethanol and 35 g 1-propanol per 100 g liquid or Mikrozyd Liquid (Schülke & Mayr GmbH, e.g. cat. no. 109203 or 109160).

Disinfectants based on glyoxal and quaternary ammonium salt can be used for submerging worktable items, the centrifuge rotor, and the waste drawer: for example, 10 g glyoxal, 12 g lauryldimethylbenzylammonium chloride, 12 g myristyldimethylbenzylammonium chloride, and 5–15% nonionic detergent per 100 g liquid; Lysetol® AF (Gigasept® Instru AF in Europe, cat. no. 107410, or DECON-QUAT® 100, Veltek Associates, Inc., cat. no. DQ100-06-167-01, in the USA).

### General instructions

- Do not use spray bottles to spray cleaning or disinfectant liquids onto surfaces of the QIAcube Connect MDx workstation. Spray bottles should be used only for items that have been removed from the workstation.
- If solvents or saline, acidic, or alkaline solutions are spilled on the QIAcube Connect MDx or if QIAGEN buffers splash the instrument hood, wipe the spilled liquid away immediately.
- Follow manufacturer's safety instruction for handling cleaning agents.
- Follow manufacturer's instruction for soaking time and concentration of the cleaning agents. Immersing for longer than the recommended soak time can harm the instrument.
- Do not use alcohol or alcohol-based disinfectants to clean the QIAcube Connect MDx hood. Exposing the QIAcube Connect MDx hood to alcohol or alcohol-based disinfectants will cause surface cracking. Clean the QIAcube Connect MDx hood with distilled water only or a mild detergent.
- Do not submerge buffer bottles in 70% alcohol. The blue ring is not ethanol resistant.
- Take care that no liquid runs down the touchscreen. Liquid may be drawn through the dust protection sealing by capillary forces and cause malfunction of the display. To clean the touchscreen, moisten a soft lint-free cloth with water, ethanol, or a mild detergent and carefully wipe the display. Wipe dry with a paper towel.

## Removal of RNase contamination

RNaseZap® RNase Decontamination Solution (Ambion, Inc., cat. no AM9780) can be used for cleaning surfaces and submerging worktable items, centrifuge rotor, and waste drawer. RNaseZap can also be used to perform decontamination by spraying the respective worktable items outside the instrument. Use the RNase removal agent according to the manufacturer's instructions. Be aware that spraying of cleaning agents may not be allowed according to your local regulations. We recommend using lint-free towels moistened with the cleaning agent.

## Removal of nucleic acid contamination

DNA-ExitusPlus™ (AppliChem, cat. no. A7089,0100) can be used for cleaning surfaces and submerge worktable items, centrifuge rotor, and waste drawer. DNA-ExitusPlus can also be used to perform decontamination by spraying the respective worktable items outside the instrument (use the nucleic acid decontamination agent according to the manufacturer's instructions). Although the supplier of DNA-ExitusPlus only recommends to clean items when there are unwanted dried residual traces of the reagent, we recommend wiping items with a wet lint-free cloth and sterile water at any rate. This is especially important for the rotor and swing-out buckets so that the buckets do not get stuck during centrifugation and positioning.

### **CAUTION** Damage to the instrument



Do not use bleach, solvents, or reagents containing acids, alkalis, or abrasives to clean the QIAcube Connect MDx.

### **CAUTION** Damage to the instrument



Do not use spray bottles containing alcohol or disinfectant to clean surfaces of the QIAcube Connect MDx. Spray bottles should be used only to clean items that have been removed from the worktable and if permitted by local laboratory operating practices.

### **WARNING** Risk of fire



Do not allow cleaning fluid or decontamination agents to come into contact with the electrical parts of the QIAcube Connect MDx.

### **WARNING** Risk of electric shock



Do not open any panels on the QIAcube Connect MDx.

### **Risk of personal injury and material damage**

Only perform maintenance that is specifically described in this user manual. Any other maintenance or repair may only be carried out by an authorized Field Service Specialist.

### **WARNING** Hazardous chemicals and infectious agents



Waste may contain toxic or infectious material and must be disposed of properly. Refer to your local safety regulations for proper disposal procedures.

### **WARNING/CAUTION** Risk of personal injury and material damage



Improper use of the QIAcube Connect MDx may cause personal injuries or damage to the instrument. The QIAcube Connect MDx must only be operated by qualified personnel who have been appropriately trained. Servicing of the QIAcube Connect MDx must only be performed by a QIAGEN field service specialist.

**WARNING** Risk of explosion



When cleaning the QIAcube Connect MDx with alcohol-based disinfectant, leave the QIAcube Connect MDx hood open to allow flammable vapors to disperse.

Only clean the QIAcube Connect MDx when worktable components have cooled down.

**WARNING** Risk of fire or explosion



When using ethanol or ethanol-based liquids on the QIAcube Connect MDx, handle such liquids carefully and in accordance with the required safety regulations. If liquid has been spilled, wipe it off and leave the QIAcube Connect MDx hood open to allow flammable vapors to disperse.


**WARNING** Toxic fumes



Do not use bleach to clean or disinfect the QIAcube Connect MDx or the labware, as bleach in contact with salts from the buffers can produce toxic fumes.

## 6.3. Regular maintenance

After running a protocol, perform the regular maintenance procedure described below.



**Run completed screen.**

1. Open the waste drawer and empty tips and columns (if necessary) into a suitable laboratory waste container.
2. Remove used disposable labware and unwanted samples and reagents from the worktable. Discard them according to your local safety regulations.

**Note:** If the robotic arm prevents you from reaching a position, do not move the robotic arm manually. Instead, proceed as follows:

Press **Move left** or **Move right** on the Run Completed screen, as needed. The robotic arm will start to move. The hood can remain open during this movement.

**Important:** Ensure that you stand clear of the instrument while the robotic arm is moving. Wait until the robotic arm has completed its movements.


3. Replace the lids of the reagent bottles and close tightly. Store the bottles according to the instructions in the relevant kit handbook.

You can now run another protocol or switch off the QIAcube Connect MDx.

## 6.4. Daily maintenance

After running the last protocol of the day and after switching from the Research software mode to the IVD mode, perform the daily maintenance procedure. The software guides you through each step to be performed:

1. To start the daily maintenance, press the **Tools** icon (🔧) on the menu bar.
2. Then press the **Maintenance** tab, and select the **Daily** subtab. The screen shows the “Last Executed” and the “Next Due” daily maintenance dates. When daily maintenance was executed, “Next Due” date stays empty until the next protocol was executed. The status of the maintenance tasks (due, last executed) will also appear in the run report.



**Daily maintenance screen.**

3. Press **Start**. Follow the instructions on the screen. Details are provided in the next steps below.  
The robotic arm will automatically move slowly to the left – even if the instrument hood is open – to provide access to the loading positions. Always stand clear of the instrument while the robotic arm is moving. Wait until the robotic arm has completed its movements before you start to unload.
4. Remove used disposable labware, adapters, and unwanted samples and reagents from the worktable. If required, discard them according to your local safety regulations.
5. Close the buffer bottles tightly and store according to the instructions in the relevant kit handbook. We recommend reusing the buffer bottles only until the kit is used up. As soon as a new QIAGEN kit is opened, new buffer bottles should be used.
6. Press **Done** to confirm that the steps have been completed.
7. Empty the waste drawer and check that the inlay is clean. If necessary, clean the inlay of the waste drawer with alcohol-based disinfection wipes, or by soaking using one of the cleaning agents listed above, and then rinse with distilled water.

- Wipe and clean the worktable with alcohol-based disinfection wipes. Incubate as appropriate, wipe thoroughly with distilled water, and wipe dry with lint-free paper towels.

**Note:** Do not use alcohol or alcohol-based disinfectants to clean the hood.


- Press **Done** only when the steps listed above have been successfully completed. The date of the last performed daily maintenance is updated automatically.

The robotic arm will automatically move back to its original position (above tip rack position 3).

## 6.5. Monthly maintenance

Perform the daily maintenance procedure (refer to Section 6.4 "Daily maintenance" on the previous page) before you perform the monthly maintenance procedure. Select the appropriate cleaning agent according to the sample material and downstream assay (refer to Section 6.1 Cleaning agents).

- To start the monthly maintenance, press the **Tools** icon (🔧) on the menu bar.
- Then press the **Maintenance** tab and select the **Monthly** subtab. The screen shows the "Last Executed" and the "Next Due" monthly maintenance dates.



### Monthly maintenance screen.

- Close the hood.
- Press **Start**. Follow the instructions on the screen. Details are provided in the next steps below.  
The robotic arm will move to the cleaning position.
- Thoroughly clean the worktable with alcohol-based disinfection wipes. Incubate as appropriate, rinse thoroughly with distilled water, and wipe dry with paper towels.

**Important:** Do not use alcohol or alcohol-based disinfectants to decontaminate the QIAcube Connect MDx hood.

- Clean the touchscreen with alcohol-based disinfection wipes and wipe dry afterwards.

**Important:** Take care that no liquid runs down the touchscreen. Liquid may be drawn through the dust protection sealing by capillary forces and cause malfunction of the display. To clean the touchscreen, moisten a soft lint-free cloth with 70%

ethanol or a mild disinfectant and carefully wipe the display. Depending on the disinfectant, wipe the screen with distilled water. Wipe dry with a paper towel.

7. Clean the outer hood with a soft lint-free cloth moistened with water or mild detergent.
8. Clean the shaker adapter (gray), shaker tray (metal adapter), buffer bottle rack (and waste drawer in liner if not done during daily maintenance) with alcohol-based disinfection wipes.
9. Incubate the shaker adapter (gray), shaker tray (metal adapter), buffer bottle rack, and waste drawer in liner (if not done during daily maintenance) by soaking as appropriate. Rinse thoroughly with distilled water and wipe dry with lint-free paper towels. If the shaker rack plugs are used, treat them in the same way.
10. Press **Done** only when the steps listed above have been successfully completed. The date of the last performed monthly maintenance is updated automatically.

**Important:** Inspect the waste drawer during maintenance. Contact QIAGEN Technical Services if any broken parts are observed.

11. It is recommended to transfer the run reports from the instrument to the USB flash drive and remove the run reports from the instrument to free up disk space. For details refer to Section 5.7 Saving run reports to the USB flash drive.

## 6.6. Periodic maintenance

The periodic maintenance consists of cleaning the robotic arm modules and centrifuge. It is recommended to be conducted, when necessary, but at least every 6 months.

Select the appropriate cleaning agent according to the sample material and downstream assay (refer to Section 6.1 Cleaning agents).




### 6.6.1. Cleaning the robotic arm modules

Cleaning of the robotic arm modules must be performed periodically or could be performed if required. For example, the robotic arm modules must be cleaned if liquids were spilled due to crash.

**Note:** Perform the monthly maintenance procedure before you perform the robotic arm cleaning procedure.

1. To start cleaning the robotic arm modules, press the **Tools** icon (🔧) on the menu bar. Press the **Maintenance** tab and select the **Robotic arm** subtab. The screen shows the “Last Executed” maintenance date of the robotic arm modules.



**Robotic arm maintenance screen.**

2. Press **Start** to begin the cleaning of robotic arm modules. Follow the instructions on the screen. Details are provided in the next steps below.
3. Make sure that used labware, adapters, and reagents are removed from the worktable. Close the hood.
4. Press **Next** to move to cleaning position.
5. Remove the waste drawer and open hood.
6. Moisten a soft lint-free cloth with water and carefully clean the optical sensor, tip adapter, gripper unit, rotor adapter stabilization rod, and the spin column lid holder. Wipe these items dry as indicated on the touchscreen of the instrument.
7. Close the hood and press **Done** to finish cleaning of robotic arm. The date of the last performed cleaning of robotic arm is updated automatically.

## 6.6.2. Cleaning the centrifuge


Cleaning of the centrifuge must be performed periodically (at least every 6 month) or could be performed if required. For example, the centrifuge must be cleaned in case of plastic crash or spillage of liquids due to crash.


### **WARNING** Risk of personal injury and material damage



To prevent plastic crash, load the tubes properly. After a plastic crash, sharp plastic particles could be inside the centrifuge. Be careful when handling items inside the centrifuge.

**Note:** Perform the monthly maintenance procedure before you perform the cleaning of centrifuge procedure.


1. To start cleaning the centrifuge, press the **Tools** icon () and press the **Centrifuge** subtab under the **Maintenance** tab. The screen shows the Last Executed centrifuge maintenance date and a cycle count.




### **Centrifuge maintenance screen.**

2. Press **Start** to begin the centrifuge cleaning procedure. Follow the instructions on the screen. Details are provided in the next steps below.

- The centrifuge lid must be open to allow access to the inside of the centrifuge. The lid should be opened only after the centrifuge has come to a complete stop. If the lid does not open automatically, close the hood and press the **Open Centrifuge Lid** button.



- Switch off the instrument and perform cleaning as outlined in the following sections (below):
  - Cleaning the rotor and buckets
  - Cleaning the centrifuge chamber
  - Maintenance of the rotor nut
  - Installing the centrifuge rotor and buckets
- When cleaning is completed, turn on the instrument and log in. Press the **Tools** icon (  ) and then the **Maintenance** tab. Select the **Centrifuge** subtab.
- Press **Start** again then press **Done** to confirm cleaning. The date of the last performed cleaning of centrifuge is updated automatically.

### Cleaning the rotor and buckets

- Make sure the QIAcube Connect MDx is switched off.
- Remove all disposable rotor adapters, including tubes and spin columns, from the buckets.
- Remove the buckets from the rotor. Undo the rotor nut on top of the rotor using the rotor key, and carefully lift the rotor off the rotor shaft.



**Rotor key.**

- Submerge the rotor, buckets, and rotor nut in cleaning agent. Incubate as appropriate.

5. Rinse thoroughly with distilled water. Use a brush (e.g., a toothbrush or tube brush) to clean any parts that are difficult to access, such as the bucket mount and the rotor head. Wipe surfaces dry with a soft lint-free cloth. If available, dry the buckets and rotor with pressurized air.



Brushing a bucket.



Brushing the rotor.


**Important:** Make sure the paper towels and brush used are lint-free.

**Important:** Make sure that all residual salt is removed.

**Important:** Make sure to remove all traces of cleaning agent from the centrifuge buckets. Residual agent can cause the buckets to jam.

6. Carefully check the rotor for damage. If the rotor is damaged or shows signs of wear or corrosion, do not use the rotor. Contact QIAGEN Technical Services.
7. Apply a few drops of mineral oil (Anti-Corrosion Oil (rotor), cat. no. 9018543) on a soft, lint-free cloth, and wipe the bucket mount and rotor claw. A thin, invisible oil film should cover the bucket mount and rotor claw, but no droplets or smear should be apparent.

**Important:** Before applying oil to the rotor buckets on the rotor, make sure that the rotor and all buckets are completely dry.



Rotor head.



Bucket mounts.

### Cleaning the centrifuge chamber

**Note:** Ensure that the instrument is switched off during cleaning.

1. Moisten a soft lint-free cloth with cleaning agent and clean the inside of the centrifuge and the centrifuge gasket. Incubate as appropriate.
2. Clean the inside of the centrifuge and the gasket with distilled water and wipe dry with lint-free paper towels. If available, use a vacuum cleaner.

**Important:** Make sure the gaskets remain in the proper positions.


3. Clean the centrifuge lid with a soft lint-free cloth moistened with cleaning agent. Incubate as appropriate, clean with water and wipe dry with lint-free paper towels.
4. Check the centrifuge gasket for damage. If the gasket is damaged or shows signs of wear, contact QIAGEN Technical Services.

### Maintenance of the rotor nut

**Note:** Ensure that the instrument is switched off during cleaning.

**Note:** Always perform the cleaning procedure after disassembly of the rotor and at least twice a year.


After cleaning the rotor thread, apply a few drops of mineral oil (Anti-Corrosion Oil (rotor), cat. no. 9018543) on a lint-free cloth, and wipe the thread. A thin, invisible oil film should cover the rotor thread but no droplets or smear should appear.



Rotor thread

Pin

Rotor thread.



Inner thread

Inner thread of the rotor nut.

After cleaning the inner thread of the rotor nut, wipe the thread using Anti-Corrosion Oil as described above.

**Note:** Contact QIAGEN Technical support if the pin on the rotor thread fell out. Do not reinsert the pin! Do not run the centrifuge!

### Installing the centrifuge rotor and buckets

**Note:** Ensure that the instrument is switched off during cleaning.

1. Mount the rotor.
2. The rotor can be mounted in only one orientation. The pin on the rotor shaft fits into a notch on the underside of the rotor directly underneath rotor position 1. Line up position 1 of the rotor with the pin on the rotor shaft and carefully lower the rotor onto the shaft.
3. Install the rotor nut on top of the rotor and tighten securely using the rotor key supplied with the QIAcube Connect MDx. Make sure that the rotor is securely seated.



Rotor key.



Rotor nut.

If the rotor nut is not tightened properly, it can become loose during operation of the centrifuge and can cause serious damage to the instrument. Such damage is not covered by the warranty.

### **WARNING** Risk of personal injury and material damage



To prevent the rotor nuts from loosening during operation of the centrifuge, securely tightened the nuts using the rotor key supplied with the QIAcube Connect MDx.

4. Insert the rotor buckets. The side of the rotor bucket that must face toward the rotor shaft is marked with a gray line. Hold the bucket at an angle with the gray line facing the center of the rotor and hang the bucket on the rotor. Check that all buckets are properly suspended and can swing freely.

**Important:** All centrifuge buckets must be mounted before starting any centrifuge run.

Before starting next protocol run, follow the instructions in Section 6.6.3 Operating the centrifuge after cleaning.


### 6.6.3. Operating the centrifuge after cleaning

The centrifuge must be operated independently before starting further runs to check if residual plastic parts are still in the centrifuge.

**Note:** Rotor adapters and other consumables are not required.

**Important:** Ensure that the rotor and all centrifuge buckets are mounted properly before starting any centrifuge run.

1. Switch the instrument on and log in.
2. To start a centrifuge run, press the **Tools** icon (🔧) on the menu bar and then the **Run Modules** tab. The **Centrifuge** tab is open by default.



3. In the Set speed and Set duration fields, set the speed to 10,000 x g and the duration to 1 min (1:0 m:s), respectively.
4. Press **Start** to begin the centrifuge run.
5. Carefully listen to the sound during centrifugation. See below for more details regarding the sound.

### Unusual sound during centrifugation

If any grinding, rattling or crunching sounds are heard during the centrifugation, there could be still loose plastic particles inside the centrifuge. Repeat the cleaning procedure as described in Section 6.6.2 Cleaning the centrifuge.

**Note:** It may be necessary to repeat the procedure several times to remove all plastic particles.

## No unusual sound during centrifugations

If no unusual sound from loose plastic particles can be heard during centrifugation, the next protocol run can be started.

**Note:** The **Open centrifuge lid** and **Close centrifuge lid** buttons are not needed to start a centrifuge run, as the lid will close automatically. Instead, they are needed in case you need to prepare the QIAcube Connect MDx for shipment or during troubleshooting.

## 6.7. Optional maintenance

### 6.7.1. UV run


The UV run is recommended to be performed daily to support decontamination of the device. It helps to reduce possible contaminants (e.g., nucleic acids and *E. coli*) of the QIAcube Connect MDx worktable. The efficiency of inactivation depends, for example, on layer thickness and sample type. QIAGEN cannot guarantee complete eradication of specific contaminants.

During UV decontamination, the robotic arm will move slowly across the worktable. The default cycle number is 1 (approx. 12 minutes) for maintenance. In case you visually see splashes on the worktable after a run, clean them first according to the instruction above (see Section 6.4), then, increase the cycle number based on used sample material or contaminants (e.g., nucleic acids or *E. coli*).

**Note:** Before starting the UV irradiation procedure, ensure that daily maintenance (see Section 6.4) is performed and thereby all samples, eluates, reagents and disposable labware are removed from the worktable and the worktable has been wiped.

During each cycle, an average summed dose rate of 28 to 46 mW\*s/cm<sup>2</sup> can be achieved by the UV LED light.

1. To start the UV irradiation, press the **Tools** icon on the menu bar. Press the **Maintenance** tab and select the **UV Run** subtab. The screen shows the Last Executed UV run date and the duration.



UV run screen.



2. In the Cycle field, change the number of cycles based on used sample material or contaminants (e.g., nucleic acids or *E. coli*). The default cycle number is 1 (approx. 12 minutes).

3. Ensure that all disposable labware has been removed from the worktable.

**Important:** Ensure that the waste drawer is closed. Do not open it during the UV run. Ensure that the rotor and the rotor buckets are installed in the centrifuge.

4. Close the hood and press **Start** to begin the UV run.

5. Press **Done** once the UV run is completed. The date of the last performed UV run is updated automatically.

#### **WARNING** Risk of personal injury




Do not expose your skin to UV-C light from the UV LED lamp.

### 6.7.2. Tightness test

To ensure that the tightness of the tip adapter is sufficient for accurate pipetting, the tightness test of the tip adapter can be performed. This test must be performed after replacing a tip adapter O-Ring to verify if replacement is successful.

**Note:** Perform the daily maintenance and cleaning of the robotic arm procedure before you perform the tightness test. See Sections 6.4 Daily maintenance and 6.6.1 Cleaning the robotic arm modules.

1. To start the tightness test, press the **Tools** icon on the menu bar. Press the **Maintenance** tab and select the **Tightness** subtab. The screen shows the Last Executed tightness test date.



#### **Tightness test screen.**

2. Press **Start** to begin with the tightness test procedure. Follow the instructions on the screen. Details are provided in the next steps below.

3. Open the hood and load a 1000 µL tip rack with at least one 1000 µL tip into tip rack position 1.

4. Place an empty 2 mL safe-lock microcentrifuge tube (cat. no. 990381) in position 1 of the shaker (shaker type 2).

5. Place a buffer bottle filled with ≥10 mL 96–100% ethanol in position 1.

6. Close the hood, and press Next to start tightness test.
7. After the load check, the robotic arm will pick up a tip, aspirate ethanol, and move to the tube. The tip will remain in place above the tube for 2 minutes. The tip will be discarded into the waste afterwards.
8. Wait until the test has been completed and then press Next.
9. Open the QIAcube Connect MDx hood and remove the buffer bottle and tips to store them accordingly.
10. Remove the tube and visually check if liquid is present:  
 If no liquid is present, press **Yes** to record that the test passed.  
 If liquid is present, press **No** to record that the test failed.

QIAGEN

Setup Tools Configuration QIAcube Network Lock Logout

Data exchange Run modules Maintenance

Daily UV Run Monthly Robotic Arm Centrifuge Tightness O-Ring

### Tightness Test - Result

been passed.  
- If liquid is present, press 'No' to record that test has been failed.

In this case repeat the test. If test fails again it is recommended to replace O-Ring first (using tab 'O-Ring') or contact QIAGEN Technical Services.

Yes No

Press [Done] to finalize tightness test procedure.

Step 3 | 3 Done

Sep 12, 2023, 09:49 Mode: Research User: Admin Admin

11. In case the test failed, repeat the test. If test fails again, it is recommended to replace the O-Ring first (see Section 7.2.5 O-Ring exchange) or contact QIAGEN Technical Services.
12. Press **Done** to finalize tightness test procedure. The date of the last performed tightness test is updated automatically.

## 6.8. Decontaminating the QIAcube Connect MDx

If the QIAcube Connect MDx is contaminated with infectious material, it should be decontaminated. If hazardous material is spilled on or inside the QIAcube Connect MDx, the user has responsibility for carrying out appropriate decontamination.

The QIAcube Connect MDx should also be decontaminated before shipping (e.g., back to QIAGEN). In this case, a decontamination certificate must be completed to confirm that the decontamination procedure has been carried out.

To decontaminate the QIAcube Connect MDx, follow the daily, monthly, and periodic maintenance procedure in Sections 6.4–6.6, using the recommended disinfection agents. In addition, perform a UV run with at least 5 cycles as outlined in Section 6.7.1.

## 6.9. QIAcube Connect MDx repair

Contact your local QIAGEN Field Service Specialist or your local distributor for more information about flexible Service Support Agreements from QIAGEN. Preventive maintenance contracts are available as well to ensure performance of the inspection at least once a year.

### **WARNING/ CAUTION**



#### **Risk of personal injury and material damage**

Improper use of the QIAcube Connect MDx may cause personal injuries or damage to the instrument. The QIAcube Connect MDx must only be operated by qualified personnel who have been appropriately trained. Servicing of the QIAcube Connect MDx must only be performed by a QIAGEN field service specialist.

## 7. Troubleshooting

This section provides information about what to do if an error occurs when using the QIAcube Connect MDx.

If further assistance is required, contact QIAGEN Technical Services using the contact information below:

**Website:** [support.qiagen.com](https://support.qiagen.com)

When contacting QIAGEN Technical Services about an error with the QIAcube Connect MDx, note the steps leading to the error and any information appearing in any dialog boxes. This information will help the QIAGEN Technical Services solve the problem.

When contacting QIAGEN Technical Services about errors, please have the following information ready:


- Protocol name and version (found in the report file)
- Software version (see Section 4.5.1)
- Serial number of the instrument can be found at the right of the **System** tab in the configuration screen.
- Sample input material
- Detailed description of the error situation, especially on the worktable status after an interrupted run.
- Download a support package from the instrument.


This information will help you and your QIAGEN Technical Service Specialist to deal most efficiently with your issue.

**Note:** Information about the latest software and protocol versions can be found at [www.qiagen.com](https://www.qiagen.com). In some cases, updates may be available for addressing specific problems.

## 7.1. Creating a support package

The support package is a zip file that can be sent to QIAGEN Technical Services for diagnosis and trouble shooting.

1. On the menu bar, press the **Tools** icon () .
2. Press the **Data Exchange** tab.
3. Connect the USB flash drive to one of the two USB ports next to the touchscreen.



**Data exchange screen.**

4. Press **Create and download** in the **Support package** section. The support package will be created and saved on the USB flash drive . It is a .zip with the filename “QIAcube-SN-YYYYMMDDhmm.zip” (where SN is the serial number of your device). The support package will contain all relevant data from the past 6 weeks including protocols, run reports, the audit trail, and log files.
5. To read the audit trail, the **.csv** file has to be imported to a suitable application (e.g., Microsoft Excel) using UTF-8 formatting.

## 7.2. Operation

### Comments and suggestions

Centrifuge	
Bucket does not swing back into place	Clean the centrifuge and rotors as described in Section 6.6.2.
Imbalance detected	Make sure the rotor is symmetrically loaded according to instructions on run setup screens. Remove the rotor and check the centrifuge chamber for loose plasticware. Power OFF the QIAcube Connect MDx, wait for a few minutes, and power it ON again. If the error persists, contact QIAGEN Technical Services.
Imbalance detected; loud noise heard during centrifugation	Ensure that loose parts are cleared from the worktable prior to centrifuge operation to avoid loose parts jamming or damaging the centrifuge.
Rotor location pin fell out	Contact QIAGEN Technical Services, do not use the centrifuge.
Shaker	
Incorrect repositioning of shaker	The shaker should re-position itself toward the right side once shaking is completed. Remove any obstructions that prevent the shaker from returning to the correct position.

## Comments and suggestions

### Robotic arm

Robotic arm does not return to set position	Ensure that the instrument is placed on a stable, flat and level surface as described in Section 4.1.1. In other cases, contact QIAGEN Technical Services.
Instrument loading	In case of incorrect instrument loading, thoroughly read the error message. It will guide you to the missing/wrong item.

### Pipettor

Pipet tips not picked up by automatic pipettor	Make sure that the tip rack is not damaged and is correctly positioned on the worktable.
Pipet tips not disposed correctly	Empty the waste drawer and ensure it is not broken. Check that tip disposal slot is not damaged or obstructed. Perform regular maintenance, as described in Section 6.3.
Droplets observed on worktable	<p>The pipettor is dripping liquid. Check that the reagent bottles contain the correct buffers and are correctly placed in the reagent bottle rack. Be sure to use the right plasticware. Check the volumes in the sample tubes and tubes of accessory buffer(s), if applicable. Do not exceed the recommended amount of starting material to avoid blocking disposable filter-tips. Do not use manually filled tip racks.</p> <p>Check tightness of the pipettor as described in Section 6.7.2. If leakage could be detected change O-Ring as described in Section 7.2.5. If issue persists, contact QIAGEN Technical Services.</p>

### Mechanical

Frame of instrument is distorted (e.g., uneven, unstable or not level)	Ensure that the instrument is placed on a stable, flat and level surface as described in Section 4.1.1.
Hood sensor error: instrument will not function	Ensure the hood is properly closed. The instrument will not function if the hood is open.
Broken instrument hood	Ensure that only the cleaning products as described in Section 6.6.1 are used on the hood.
Waste drawer jams but can still be inserted	Empty the waste drawer. Perform daily maintenance, as described in Section Daily maintenance.
Incorrectly inserted waste drawer	Handle the waste drawer with both hands when inserting or removing the drawer.
Pipet tips not disposed of correctly	Make sure the top of the tip disposal slot (refer to Section 3.3) is not broken.
Scratches appear on the instrument	Always use the cleaning products as described in Section Cleaning the robotic arm modules. Do not use bleach or ethanol, as they can damage the surface of the instrument.

### Electronic

Display does not turn on	Do not touch the display with excessive force, do not use corrosive chemicals to clean the display surface. Contact QIAGEN Technical Services for repair.
Error when copying files to USB	Power OFF the QIAcube Connect MDx, wait for a few minutes, and power it ON again. Save the file(s) to the USB flash drive again. Check the USB flash drive on a PC to ensure it is functional. If possible, format the USB flash drive before use on the instrument. If the error persists, contact QIAGEN Technical Services.
USB device not detected	<p>Be sure to use only use the USB stick delivered with the instrument.</p> <p>Power OFF the QIAcube Connect MDx, wait for a few minutes, and power it ON again. Insert the USB flash drive into the USB port. Check the USB flash drive on a PC to ensure it is functional. Check if only one USB flash drive is inserted. Otherwise, the instrument will not detect a USB flash drive. If the error persists, contact QIAGEN Technical Services.</p>
Login screen not visible when launching instrument	If the touchscreen does not display the login screen, but instead a software update message is shown, power OFF the QIAcube Connect MDx, wait for a few minutes. Ensure that the USB flash drive is not inserted in the USB port. Power ON the QIAcube Connect MDx again. The login screen should be visible. If the error persists, contact QIAGEN Technical Services.
System only knows default user	The user files were disturbed. Log in with the default user (Admin) and the default password to create a new user file. If the error persists, contact QIAGEN Technical Services.
Error displayed when inserting the USB flash drive into a Windows PC	Ignore the message. In most cases, no scan is needed; use the USB flash drive as usual. Reformat the USB on the Windows PC after securing all data, which has been stored on it

### 7.2.1. Protocol interruption

If an error occurs during a protocol run, it is possible to continue sample preparation manually.

**Important:** It is not recommended for QIAGEN DSP/IVD protocols to finish the run manually; the run will be invalid and the sample result from the manual continuation of the protocol must not be used for diagnostic purposes. It is under your responsibility to continue sample processing manually, as this invalidates the whole procedure.

The error code, the description and the step at which the protocol stopped are displayed in the touchscreen.

To continue sample processing:

1. Note the step at which the protocol stopped. This is displayed in the touchscreen under **Run Details**.
2. Remove the samples and reagents from the QIAcube Connect MDx.
3. Refer to the appropriate protocol in the relevant kit handbook, find the last performed protocol step (e.g., Wash step) and continue sample processing manually.

### 7.2.2. Centrifuge

#### Opening the centrifuge lid in the event of a breakdown

In case of power failure, the centrifuge lid can be manually opened so that the samples can be removed. To open the centrifuge lid, follow the instructions below.

#### **WARNING** Moving parts



In case of breakdown caused by power failure, remove the power cord and wait 10 minutes before attempting to manually open the centrifuge lid.


#### **WARNING** Risk of personal injury and material damage



Raise the centrifuge lid carefully. The lid is heavy and may cause injury if it falls.

1. Power OFF the QIAcube Connect MDx.
2. Unplug the power cord from the power outlet. Wait 10 minutes for the rotor to stop.
3. Open the instrument hood.
4. Carefully move the robotic arm to the right side of the worktable, furthest away from the centrifuge lid.

5. Remove the screw protection on top of the centrifuge lid. Using the rotor key, turn the screw counterclockwise.



**Turning the screw in the centrifuge lid.**

6. Remove the waste drawer. The centrifuge release cord will be visible on the left side of the waste drawer compartment.




**Removed waste drawer.**



**Centrifuge release cord.**

7. Pull the cord firmly to release the lid from the lock.
8. Manually raise the centrifuge lid.
9. Hold the raised lid and remove the samples and rotor adapters from the rotor.



**Removing rotor adapters.**

Contact QIAGEN Technical Services for instructions on how to reset the lid.



### Liquid spills in the centrifuge

The rotor adapter is designed for use with QIAGEN automated protocols. Do not fill the rotor adapters with liquid.

Liquid spills may occur if QIAGEN spin columns become blocked due to sample overloading. Do not exceed the maximum amount of starting material.

Incorrect installation of the centrifuge buckets may also cause rotor adapters to leak. Check that the buckets are installed properly and can swing freely.

If there is a liquid spill in the centrifuge, clean according to the instructions in Section 6.

### 7.2.3. Reagent volume detection and ultrasonic pipe

To help prevent errors during detection of reagent volumes, make sure that both rack labeling strips are attached to the reagent bottle rack. These strips ensure that the reagent bottle rack is positioned correctly on the worktable for liquid detection during the load check.

The instrument does not start a load check if the ultrasonic pipe (black cap) of the ultrasonic sensor is missing. Check if cap has been installed before starting a load check.



Black beam collimator (see red circle) of the ultrasonic sensor.

### 7.2.4. Touchscreen

Every time the user presses a button on the touchscreen, a small red sign is displayed at the place where the touchscreen recognizes the contact. If the point of touch and the recognized contact are at different positions, a recalibration of the touchscreen can be performed. The calibration function can be reached during the instrument's start-up procedure.

It is recommended to use a touch pen or an unused tip for optimal calibration results. In case you use a tip, discard the tip after calibration.

To re-calibrate the touchscreen:


1. Turn OFF the QIAcube Connect MDx.
2. Wait a few minutes and then turn ON the instrument again.
3. In the second screen press the QIAGEN logo.

**Note:** If you do not press the logo, the instrument will continue initialization.




Startup screen.

4. Press **Calibrate Touch Screen**.




Initial screen for the touchscreen calibration.

5. An information message is displayed. The message will close automatically after 10 seconds if you are unable to close it by pressing **OK**.




6. After 10 seconds or pressing **OK**, the following screen appears.

7. Touch the plus sign in the upper left of the screen.



8. One after another, plus signs (+) will be shown at different positions on the screen. For each of these, press the center of the + sign. After one position is touched, the next position will be shown. The graphic below shows the positions and the order in which the + signs will appear.



[Expected spots for touchscreen calibration.](#)

9. After all five positions have been touched, the following screen appears.



10. Press **Quit** to continue initialization using the new calibration settings.

11. To cancel the calibration process, turn OFF the QIAcube Connect MDx.

### 7.2.5. O-Ring exchange

O-Ring replacement must be performed if the tightness test (see Section 6.7.2 Tightness test) failed or if the following issues are observed:

- Uneven volume transfers
- Dripping on the worktable

In any case, it is recommended to consult QIAGEN Technical Services. The replacement procedure requires the O-Ring change tool and an O-Ring. See Appendix B – QIAcube Connect MDx Accessories for ordering details.




**O-Ring tool with prepared new O-Ring.**

The O-Ring replacement is semi-automatically and includes movement of the robotic arm.

**Note:** Perform the daily maintenance and cleaning of the robotic arm procedure before you replace the O-Ring.

1. To start the O-Ring replacement, press the **Tools** icon on the menu bar. Press the **Maintenance** tab and select the **O-Ring** subtab. The screen shows the last O-Ring replacement date.



#### O-Ring maintenance screen.

2. Close the hood and press **Start** to begin the O-Ring replacement procedure. Follow the instructions on the screen. Details are provided in the next steps below.
3. To prepare the O-Ring tool, perform the following steps:
  - a. Slide the new O-Ring over the small end of the peg.



#### How to attach the new O-Ring.

- b. Push the gray lever until you reach the black stop, and insert the small end of the peg into the hole.

- c. Press the peg down using the back end of the tweezers until the O-Ring sits (in the middle) on the larger end of the peg.



Press the peg down using the back end of the tweezer.

- d. Open the gray lever and insert the peg with small end first into hole as shown.



Insert the peg to the hole.


4. Press **Next** on the screen and begin loading of the O-Ring tool into the QIAcube Connect MDx.
5. Load the O-Ring tool by opening the gray lever into tip rack position 1 (nearest to user).




Load the O-ring tool.

6. Close the hood and press **Next** to begin cutting of the old O-Ring.
7. To cut and remove the old O-Ring, perform the following steps:

- a. To cut the O-Ring, open the hood and rotate the gray lever counter-clockwise until you reach the black stop.



Status after opening the hood.



Rotate the gray lever counter-clockwise.

- b. Open the gray lever and remove the O-Ring (by using the tweezers) from pipetting channel.


**Note:** If required, repeat cutting process until O-Ring is cut completely and can be removed.



Open the gray lever and remove the O-Ring using the tweezers.

8. Close the hood and press **Next** to pick up the prepared new O-Ring.

9. Open the hood and visually check if new O-Ring sits firmly on tip adapter.



Check if the new O-Ring sits firmly.

**Note:** If the O-Ring was not successfully picked up, complete the O-Ring replacement procedure and restart.

10. Close the hood press **Next**.
11. Open the hood and remove the O-Ring change tool.
12. Wipe and clean the O-Ring change tool with alcohol-based disinfection wipes. Incubate as appropriate, rinse thoroughly with distilled water and wipe dry with lint-free paper towels.
13. Press **Done** to complete the O-Ring replacement. The date of the last performed O-Ring replacement is updated automatically.

**WARNING/  
CAUTION**



**Risk of personal injury and material damage**

Improper use of the QIAcube Connect MDx may cause personal injuries or damage to the instrument. The QIAcube Connect MDx must only be operated by qualified personnel who have been appropriately trained. Servicing of the QIAcube Connect MDx must only be performed by a QIAGEN field service specialist.



## 8. Glossary

Term	Description
Centrifuge	A component of the QIAcube Connect MDx that accommodates a rotor with 12 swing-out buckets. Each bucket holds a disposable rotor adapter.
Disposal slots	Slots in the QIAcube Connect MDx worktable through which used tips and columns (e.g., QIAshredder columns) are discarded into the waste drawer.
Hood	The main door at the front of the QIAcube Connect MDx. When open, it provides complete access to the worktable.
Error code	A 3- or 4-digit number that indicates an error of the QIAcube Connect MDx.
Filter-tip	An item of labware that is picked up by the tip adapter during operation of the QIAcube Connect MDx. Liquid is aspirated into and dispensed from a filter-tip.
Initialization	An operation performed automatically when the QIAcube Connect MDx is switched on and if required before each protocol run to check the operation of the QIAcube Connect MDx.
Microcentrifuge tube slots	Three slots located in the labware tray that accommodate accessory buffers, in 1.5 mL or 2 mL microcentrifuge tubes.
Pipetting system/ Pipettor unit	The component of the QIAcube Connect MDx that aspirates and dispenses liquid. The pipetting system moves up and down above the worktable and contains a syringe pump that is connected to a tip adapter.
Power switch	A button located at the front of the QIAcube Connect MDx in the bottom-right corner. It allows the user to switch the QIAcube Connect MDx on and off; inner position is ON and outer position is OFF.
Protocol	A set of instructions for the QIAcube Connect MDx that allows the instrument to automate a nucleic acid or protein purification procedure. Protocols are run using the touchscreen.
Reagent bottle rack	A rack that can accommodate six 30 mL bottles on the QIAcube Connect MDx worktable.
Robotic gripper	A component of the QIAcube Connect MDx robotic arm that moves spin columns during sample processing.
Rotor adapter	A disposable plastic adapter that fits into a centrifuge bucket and holds a QIAGEN spin column and microcentrifuge tube during sample processing.
Tip adapter	A metal probe installed on the pipettor head. During operation of the QIAcube Connect MDx, the tip adapter picks up filter-tips from the worktable.
Tip rack	A plastic rack that accommodates filter-tips on the worktable.
Touchscreen	The user interface that allows the user to operate the QIAcube Connect MDx.
Waste drawer	A drawer that collects used filter-tips and disposable columns.
Worktable	The surface of the QIAcube Connect MDx where samples, reagents, and filter-tips are loaded.

## 9. Technical Specifications

QIAGEN reserves the right to change specifications at any time.

### 9.1. Operating conditions

Description	Requirement
Power	100–240 V AC, 50/60 Hz, 650 VA Mains supply voltage fluctuations are not to exceed 10% of nominal supply voltages. <b>Note:</b> The apparent power can exceed 650 VA for up to 2 seconds during the centrifuge acceleration and can reach an approximate value of 1200 VA.
Fuse	2x T8A L 250V
Overvoltage category	II
Air temperature	18–28°C (64.4–82.4°F)
Relative humidity	15–75% (noncondensing)
Altitude	Up to 2000 m (6500 ft.)
Place of operation	For indoor use only
Pollution level	2
Environmental class	IEC 60721-3-3

### 9.2. Transport conditions

Description	Requirement
Air temperature	–25°C to 60°C (–13°F to 140°F) in manufacturer's package
Relative humidity	Max. 75% (noncondensing)
Environmental class	2K2 & 2M2 (IEC 60721-3-2)

### 9.3. Storage conditions

Description	Requirement
Air temperature	5°C to 40°C (41°F to 104°F) in enclosed location
Relative humidity	Max. 75% (noncondensing)
Environmental class	1K2 (IEC 60721-3-1)

## 9.4. Mechanical data and hardware features

Dimensions (hoods closed)	<b>Width:</b> 65 cm (25.6 in.) <b>Height:</b> 58 cm (22.8 in.) <b>Depth:</b> 62 cm (24.4 in.)
Dimensions (hood open)	<b>Width:</b> 65 cm (25.6 in.) <b>Height:</b> 86 cm (34 in.) <b>Depth:</b> 62 cm (24.4 in.)
Mass	QIAcube Connect MDx: 73 kg (160.9 lb.) Accessories: 3 kg (6.6 lb.)
Centrifuge	10,640 rpm maximum 12,000 x g maximum Swing-out rotor, maximum 45degrees 12 rotor positions
Shaker	<b>Speed:</b> 100–2000 RPM <b>Amplitude:</b> 2 mm Heating range of ambient temperature to 70°C (158°F) Ramp-up time of <5 minutes from ambient temperature to 55°C (±3°C) Difference in the temperature detected by the internal sensor and the temperature of the sample liquid is approximately –2°C
Pipetting system	Syringe size 1 mL Pipetting range 5–900 µL
Capacity	Up to 12 samples per run
Touchscreen	10.4" TFT Touchscreen, active area 211.2 x 158.4 mm, resolution 800*600 SVGA
USB flash drive	USB2.0 Compatible with devices and Operating Systems that support the USB standard <b>Operating temperature range:</b> 0°C to 70°C <b>Storage temperature range:</b> –40°C to 85°C (–40°F to 185°F) <b>Formatting:</b> FAT32
UV LED	<b>Wavelength:</b> 278 nm <b>Optical power:</b> 200–300 mW
Barcode scanner	<b>Scan Pattern:</b> Area Image (838 x 640-pixel array) <b>Motion Tolerance:</b> Up to 610 cm/s (240 in/s) for 13 mil UPC at optimal focus <b>Symbol Contrast:</b> 20% minimum reflectance difference <b>Decode Capability:</b> Reads standard 1D, PDF, 2D, Postal, and OCR symbologies
Software	QIAGEN protocols are preinstalled on the QIAcube Connect MDx or can be downloaded at <a href="http://www.qiagen.com">www.qiagen.com</a>

## Appendix A – Legal

### Declaration of Conformity

Name and address of the legal manufacturer:

QIAGEN GmbH  
QIAGEN Strasse 1  
40724 Hilden  
Germany

An up-to-date declaration of conformity can be requested from QIAGEN Technical Services.

### Waste Electrical and Electronic Equipment (WEEE)

This section provides information about disposal of waste electrical and electronic equipment by users.

The crossed-out wheeled bin symbol (see below) indicates that this product must not be disposed of with other waste; it must be taken to an approved treatment facility or to a designated collection point for recycling, according to local laws and regulations.

The separate collection and recycling of waste electronic equipment at the time of disposal helps to conserve natural resources and ensures that the product is recycled in a manner that protects human health and the environment.



Recycling can be provided by QIAGEN upon request at additional cost. In the European Union, in accordance with the specific WEEE recycling requirements and where a replacement product is being supplied by QIAGEN, free recycling of its WEEE-marked electronic equipment is provided.

To recycle electronic equipment, contact your local QIAGEN sales office for the required return form. Once the form is submitted, you will be contacted by QIAGEN either to request follow-up information for scheduling collection of the electronic waste or to provide you with an individual quote.

## EMC Declaration

The IVD medical equipment complies with the emission and immunity requirements described in IEC 61326-2-6.

The "United States Federal Communications Commission" (USFCC) (in 47 CFR 15. 105) declared that the users of this product must be informed of the following facts and circumstances.

"This device complies with part 15 of the FCC:

Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation."

"This Class B digital apparatus complies with Canadian ICES-0003."

The following statement applies to the products covered in this manual, unless otherwise specified herein. The statement for other products will appear in the accompanying documentation.

**Note:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules and meets all requirements of the Canadian Interference-Causing Equipment Standard ICES-003 for digital apparatus. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that the interference will not occur in an installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

QIAGEN GmbH Germany is not responsible for any radio television interference caused by unauthorized modifications of this equipment or the substitution or attachment of connection cables and equipment other than those specified by QIAGEN GmbH, Germany. The correction of interference caused by such unauthorized modification, substitution or attachment will be the responsibility of the user.

## California Proposition 65

Using this product can expose you to chemicals including lead acetate, which is known to the state of California to cause cancer and DEHP, which is known to the State of California to cause birth defects and/or other reproductive harm. For more information, go to [www.p65warnings.ca.gov](http://www.p65warnings.ca.gov).

## Liability Clause

QIAGEN shall be released from all obligations under its warranty in the event repairs or modifications are made by persons other than its own personnel, except in cases where the Company has given its written consent to perform such repairs or modifications.

All materials replaced under this warranty will be warranted only for the duration of the original warranty period, and in no case beyond the original expiration date of original warranty unless authorized in writing by an officer of the Company. Read-out devices, interfacing devices, and associated software will be warranted only for the period offered by the original manufacturer of these products. Representations and warranties made by any person, including representatives of QIAGEN, which are inconsistent or in conflict with the conditions in this warranty shall not be binding upon the Company unless produced in writing and approved by an officer of QIAGEN.

The QIAcube Connect MDx is equipped with an Ethernet port and a Wi-Fi USB device (optional). The Purchaser of the QIAcube Connect MDx is solely responsible for preventing any and all computer viruses, worms, trojans, malware, hacks, or any other type of cybersecurity breaches. QIAGEN assumes no liability for computer viruses, worms, trojans, malware, hacks, or any other type of cybersecurity breaches.

## Appendix B – QIAcube Connect MDx Accessories

For more information and an up-to-date list of available protocols, visit [www.qiagen.com/QIAcube-Connect-MDx](http://www.qiagen.com/QIAcube-Connect-MDx) and go to the **Resources** tab.

### Ordering information

Product	Contents	Cat. no.
QIAcube Connect MDx	Instrument and 1-year warranty on parts and labor	9003070
QIAcube Connect MDx System FUL-2	Instrument and service agreement package: includes installation, application training, Full Agreement for two years with a two business day response time and two inspection service visits.	9003071
QIAcube Connect MDx System FUL-3	Instrument and service agreement package: includes installation, application training, Full Agreement for three years with a two business day response time and three inspection service visits.	9003072
QIAcube Connect MDx System PRV-1	Instrument and service agreement package: includes installation, application training and one preventive maintenance visit. One-year warranty on labor, travel and parts also included.	9003073
QIAcube Connect MDx Device PRV-1	Instrument and service agreement package: includes one inspection service visit. One-year warranty on labor, travel and parts are also included. Does not include installation and training.	9003074
QIAcube Connect MDx System PRM-1	Instrument and service agreement package: includes installation, application training, Premium Agreement for 1 year with a next business day response time and 1 inspection service visit.	9003075
QIAcube Connect, Premium Agreement	On-site repair service with next business day response time. Includes 1 Inspection Service and travel, labor, and parts	9245209
QIAcube Connect, Full Agreement	On-site repair service with two business day response time. Includes 1 Inspection Service and travel, labor, and parts	9245208
QIAcube Connect, Core Agreement	On-site repair service and one on-site Inspection Service, including travel, labor and spare parts for a period of 1 year. Response time of five business days. Includes 10% discount on additional repair service during the agreement period.	9245260
QIAcube Connect, Installation & Training	On-site installation and setup of instrument hardware and system software. Training demonstrating routine maintenance, basic troubleshooting, and more for up to 4 laboratory staff members.	9245211
Starter Pack, QIAcube	200 µL filter-tips (1024); 1000 µL filter-tips (1024); 30 mL reagent bottles (12); rotor adapters (240); 1.5 mL elution tubes (240); rotor adapter holder (1)	990395
QIAcube Connect, IQ/OQ Service	On-site installation qualification and operational quality service.	9245232
<b>Other consumables</b>		
Filter-Tips, 1000 µL (1024)	Disposable Filter-Tips, racked; (8 x 128)	990352



Product	Contents	Cat. no.
Filter-Tips, 1000 µL, wide-bore (1024)	Disposable Filter-Tips, wide-bore, racked; (8 x 128); not required for all protocols	990452
Filter-Tips, 200 µL (1024)	Disposable Filter-Tips, racked; (8 x 128); not required for all protocols	990332
Rotor, centrifuge	Rotor for the QIAcube centrifuge	9017848
Swing-out Buckets	Swing-out buckets for the rotor of the QIAcube centrifuge	9017849
Rotor Adapters (10 x 24)	For 240 preps: 240 Disposable Rotor Adapters and 240 microcentrifuge tubes (1.5 mL); for use with QIAcube instruments	990394
Rotor Adapter Holder	Holder for 12 disposable rotor adapters; for use with QIAcube instruments	990392
Reagent Bottle Rack	Rack for accommodating 6 x 30 mL reagent bottles on the QIAcube instruments worktable	9026197
Reagent Bottles, 30 mL (6)	Reagent Bottles (30 mL) with lids; pack of 6; for use with the QIAcube instruments reagent bottle rack	990393
Shaker Rack Plugs (12)	For use with 2 mL screw cap tubes	9017854
Sample Tubes RB (2 mL)	1000 safe-lock microcentrifuge tubes (2 mL) for use with the QIAcube instruments	990381
Sample Tubes CB (2 mL)	1000 conical screw-cap tubes without skirted base (2 mL) for use with the QIAcube instruments	990382
1.5 mL elution tubes	Set of 50; for use with QIAcube instruments	1050875
USB flash drive	USB flash drive; for use with QIAcube instruments	9026881
O-Ring Change Tool	O-Ring change tool for use with QIAcube instruments	9026181
O-Ring Set	Set of 10 O-Rings for use with QIAcube instruments	9018472

For up-to-date licensing and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

# Document Revision History

Revision	Changes
R3, March 2025	Update of all screenshots to reflect new software version 2.0. General review of document to keep it up to date.
R2, April 2022	Added statement on reporting incidents. Updated the URL for QIAcube Connect MDx webpage. Updated the Ordering Information section.
R1, December 2020	Initial release.

Trademarks: QIAGEN®, Sample to Insight®, QIAcube®, QIASphere®, Qproteome® (QIAGEN Group); DNA-ExitusPlus™ (AppliChem); RNaseZap® (Ambion, Inc.); Sarstedt® (Sarstedt AG and Co.); Microsoft®, Windows® (Microsoft Corporation); PAXgene® (PreAnalytiX GmbH); Gigasept®, Lysetol®, Mikrozyd® (Schülke & Mayr GmbH); DECON-QUAT® (Veltek Associates, Inc.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

R3 03/2025 HB-2794-003 © 2025 QIAGEN, all rights reserved.

