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产品	规格	货号
Rotor-Gene Q 2plex	Real-time PCR cyclers with 2 channels (green, yellow), laptop computer, software, accessories, 1-year warranty on parts and labor, optional installation and training	Inquire
Rotor-Gene Q 2plex HRM	Real-time PCR cycler and HRM instrument with 2 channels (green, yellow) plus HRM channel, laptop computer, software, accessories, 1-year warranty on parts and labor, optional installation and training	Inquire
Rotor-Gene Q 5plex	Real-time PCR cycler with 5 channels (green, yellow, orange, red, crimson), laptop computer, software, accessories, 1-year warranty on parts and labor, optional installation and training	Inquire
Rotor-Gene Q 5plex HRM	Real-time PCR cycler and HRM instrument with 5 channels (green, yellow, orange, red, crimson) plus HRM channel, laptop computer, software, accessories, 1-year warranty on parts and labor, optional installation and training	Inquire
Rotor-Gene Q 6plex	Real-time PCR cycler with six channels (blue, green, yellow, orange, red, crimson), laptop computer, software, accessories, 1-year warranty on parts and labor, optional installation and training	Inquire
Warranty PLUS 1 Basic, Rotor-Gene Q	2-year warranty, all labor, travel, and parts	9241779
Rotor-Disc 100 Starter Kit	Kit includes: 2 Rotor-Disc 100 packs, Rotor-Disc Heat Sealer, Rotor-Disc Heat Sealing Film, Rotor-Disc 100 Rotor and Locking Ring, Rotor-Disc 100 Loading Block, Rotor-Disc Pipetting Aid	Inquire
Rotor-Disc 100 (30)	Pack of 30 individually wrapped discs for 3000 reactions	981311
Rotor-Disc 72 Starter Kit	Kit includes: 3 Rotor-Disc 72 packs, Rotor-Disc Heat Sealer, Rotor-Disc Heat Sealing Film, Rotor-Disc 72 Rotor and Locking Ring, Rotor-Disc 72 Loading Block, Rotor-Disc Pipetting Aid	Inquire
Rotor-Gene Probe PCR Kit (400)*	For 400 x 25 µl reactions: 2x Master Mix, RNase-Free Water	204374
Rotor-Gene Probe RT-PCR Kit (400)	For 400 x 25 µl reactions: 2x Master Mix, RT Mix, RNase-Free Water	204574
QuantiFast Probe Assay (400)*	For 400 x 25 µl reactions: dual-labeled, probe-based, predesigned 20x lyophilized assays; includes master mix and reagents for real-time one-step or two-step RT-PCR	Varies
Rotor-Gene Multiplex PCR Kit (400)*	For 400 x 25 µl reactions: 2x Master Mix, RNase-Free Water	204774
Rotor-Gene Multiplex RT-PCR Kit (400)*	For 400 x 25 µl reactions: 3 x 1.7 ml 2x Rotor-Gene Multiplex RT-PCR Master Mix, 100 µl Rotor-Gene RT Mix, 2 x 2 ml RNase-Free Water	204974
QuantiTect Primer Assay (200)*	For 200 x 50 µl reactions: 10x QuantiTect Primer Assay (lyophilized)	Varies

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产品	规格	货号
QuantiTect Reverse Transcription Kit (50)*	For 50 x 20 µl reactions: Buffers, Quantiscript® Reverse Transcriptase, RT Primer Mix, RNase-Free Water	205311
miScript SYBR Green PCR Kit (200)*†	For 200 x 50 µl reactions: 2x Master Mix, Universal Primer, RNase-Free Water	218073
QuantiTect Virus +ROX Vial Kit (200)*	For 200 x 50 µl reactions: 5x Master Mix, RT Mix, RNase-Free Water, Nucleic Acid Dilution Buffer	211033
QuantiFast Pathogen PCR +IC Kit (400)*	For 400 x 25 µl reactions: Master Mix, lyophilized Internal Control Assay, lyophilized Internal Control DNA, ROX Dye Solution, High-ROX Dye Solution, RNase-Free Water, Nucleic Acid Dilution Buffer, Buffer TE	211354
QuantiFast Pathogen RT-PCR +IC Kit (400)*	For 400 x 25 µl reactions: Master Mix, RT Mix, lyophilized Internal Control Assay, lyophilized Internal Control RNA, ROX Dye Solution, High-ROX Dye Solution, RNase-Free Water, Nucleic Acid Dilution Buffer, Buffer TE	211454
Type-it Fast SNP Probe PCR Kit (800)*	For 800 x 25 µl reactions: 2x Master Mix, RNase-Free Water, Q-Solution®	206045
EpiTect MethyLight PCR Kit (200)*	For 200 x 50 µl reactions: 2x Master Mix, RNase-Free Water	59436
EpiTect HRM PCR Kit (100)	For 100 x 25 µl reactions: 2x Master Mix, RNase-Free Water	59445
Type-it HRM PCR Kit (100)	For 100 x 25 µl reactions: 2x Master Mix, RNase-Free Water	206542
Pathway-Powered RT ² Profiler PCR Arrays	PCR Array includes SYBR Green-optimized primer assays for a panel of relevant pathway- or disease-focused genes	330231
Custom RT ² Profiler PCR Arrays	Custom PCR Arrays containing any set of human, mouse, rat, dog, rhesus macaque, or fruit fly genes	330131
RT ² SYBR Green ROX FAST Mastermix	For 200 x 20 µl reactions: 2x SYBR Green ROX FAST qPCR Mastermix, HotStart DNA Taq Polymerase, dNTP mix, dyes	330520

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Trademarks: QIAGEN®, QIAgility®, QIAAsymphony®, cadon®, EpiTect®, GeneGlobe®, mericon®, HRM®, Q-Bond®, Q-Solution®, Quantiscript®, QuantiTect®, QuantiFast®, Rotor-Disc®, Rotor-Gene®, ScreenClust™, Type-it® (QIAGEN Group); Alexa Fluor®, Marina Blue®, SYBR®, SYTO® Texas Red® (Life Technologies); CAL Fluor® (Biosearch Technologies, Inc.); Cy® (GE Healthcare); EvaGreen® (Biotium, Inc.); FAM™, JOE™, ROX™, TET™, VIC® (Applied Biosystems or its subsidiaries); LC Green® (Idaho Technology, Inc.); LightCycler®, TaqMan® (Roche Group); Quasar® (Biosearch Technologies, Inc.); Yakima Yellow® (Nanogen, Inc.); MAX™ (Integrated DNA Technologies, Inc.).

For HRM versions only: Hi-Res Melting™ Mark under license from Idaho Technology, Inc for applicable countries.

For applicable countries:

The purchase of this product includes a limited, non transferable license to one or more of U.S. Patent Nos. 6,787,338; 7,238,321; 7,081,226; 6,174,670; 6,245,514; 6,569,627; 6,303,305; 6,503,720; 5,871,908; 6,691,041; 7,387,887; U.S. Patent Application Nos. 2003/0224434, and 2006-0019253, and all continuations and divisionals, and corresponding claims in patents and patent applications outside the United States, owned by the University of Utah Research Foundation, Idaho Technology, Inc., Evotec Biosystems GmbH, and/or Roche Diagnostics GmbH, for internal research use or for non-in vitro diagnostics applications. No right is conveyed, expressly, by implication or estoppel, for any reagent or kit, or under any other patent or patent claims owned by the University of Utah Research Foundation, Idaho Technology, Inc., and/or Roche Diagnostics GmbH, or by any other Party. For information on purchasing licenses for in-vitro diagnostics applications or reagents contact Roche Molecular Systems, 4300 Hacienda Drive, Pleasanton, CA 94588, USA.

For applicable countries:

This real-time thermal cycler is licensed under U.S. Patent No. 6,814,934 and corresponding claims in any foreign counterpart patent thereof owned by Applied Biosystems LLC, in all fields, including research and development, all applied fields, and human and animal in-vitro diagnostics. No rights are conveyed expressly, by implication or estoppel to any patents on real-time methods, including but not limited to 5' nuclease assays, or to any patent claiming a reagent or kit. For further information on purchasing additional rights, contact the Director of Licensing at Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California, 94404, USA.

Purchase of this product (Rotor-Gene SYBR Green Kits and miScript SYBR Green PCR Kit) is accompanied by a limited, non-transferable immunity from suit to use it with detection by a dsDNA-binding dye as described in U.S. Patents Nos. 5,994,056 and 6,171,785 and corresponding patent claims outside the United States for the purchaser's own internal research. No real-time apparatus or system patent rights or any other patent rights, and no right to use this product for any other purpose are conveyed expressly, by implication or by estoppel.

For all kits propagating melting curve analysis: The purchase of this product includes a limited, non-transferable license under U.S. Patent No. 5,871,908 and all continuations and divisionals, and corresponding claims in patents and patent applications outside the United States, owned by Roche Diagnostics GmbH, for internal research use or for non-in vitro diagnostics applications with authorized reagents with regard to Melting Curve Analysis. No right is conveyed, expressly, by implication or estoppel, under any other patent or patent claims owned by Roche Diagnostics GmbH, or by any other Party.

This product (Rotor-Gene Probe Kits, Rotor-Gene Multiplex PCR Kit, QuantiTect Virus +ROX Vial Kit, Type-it Fast SNP Probe PCR Kit, and EpiTect MethylLight PCR +ROX Vial Kit) is an Authorized 5' Nuclease Core Kit without Licensed Probe. Its purchase price includes a limited, non-transferable immunity from suit under certain patents owned by Roche Molecular Systems, Inc. or F. Hoffmann-La Roche Ltd, for using only this amount of the product in the practice of the 5' nuclease process solely for the purchaser's own internal research when used in conjunction with Licensed Probe. No right under any other patent claims (such as apparatus or system claims) and no right to use this product for any other purpose is hereby granted expressly, by implication or by estoppel. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

The use of this product (miScript SYBR Green PCR Kit) is covered by at least one claim of U.S. Patent No. 7,687,247 owned by Life Technologies Corporation. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product, (b) its components, or (c) materials made by the employment of this product or its components to a third party or otherwise use this product or its components or materials made by the employment of this product or its components for Commercial Purposes. Commercial Purposes means any activity for which a party receives or is due to receive consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. The buyer cannot use this product or its components or materials made using this product or its components for therapeutic, diagnostic or prophylactic purposes. Further information on purchasing licenses under the above patents may be obtained by contacting the Licensing Department: Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, CA 92008. Email: outlicensing@lifetech.com

Use of this product (RT² Profiler PCR Arrays) is covered by one or more of the following US patents and corresponding patent claims outside the US: 5,994,056 and 6,171,785. The purchase of this product includes a limited, nontransferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

This product (RT² Profiler PCR Arrays) is provided under an agreement between Molecular Probes, Inc. and SABiosciences and the manufacture, use, sale, or import of this product is subject to one or more of U.S. Patent Nos. 5,436,134; 5,658,751 and corresponding international equivalents, owned by Invitrogen Corp. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer, where such research does not include testing, analysis or screening services for any third party in return for compensation on a per test basis. The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product of its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402 Tel: (541) 465-8300, Fax: (541) 335-0504.

The purchase of this product (RT² Profiler PCR Arrays, Rotor-Gene SYBR Green Kits, miScript SYBR Green PCR Kit) includes a limited, non-transferable license under specific claims of U.S. Patent Nos. 6,174,670, 6,569,627 and 5,871,908, owned by the University of Utah Research Foundation or Evotec Biosystems GmbH and licensed to Idaho Technology, Inc. and Roche Diagnostics GmbH, to use only the enclosed amount of product according to the specified protocols. No right is conveyed, expressly, by implication, or by estoppel, to use any instrument or system under any claim of U.S. Patent Nos. 6,174,670, 6,569,627 and 5,871,908, other than for the amount of product contained herein.

The purchase of this product (QuantiTect Virus +ROX Vial Kit, RT² SYBR Green ROX FAST Mastermix) includes a limited, non-transferable right to use the purchased amount of the product to perform Applied Biosystems' patented Passive Reference Method for the purchaser's own internal research. No right under any other patent claim and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

RESEARCH USE ONLY: A license under U.S. Patents 5,786,146, 6,017,704, 6,265,171 and 6,200,756 or their foreign counterparts is required to practice MSP, Methylation-Specific PCR. The purchase price of this product (EpiTect MethylLight PCR Kit) includes limited, nontransferable rights to use only this amount of the product to practice MSP solely for the internal research activities of the purchaser. To obtain rights to practice MSP in any in vitro diagnostic application or to perform or offer commercial services of any kind using MSP, contact OncoMethylome Sciences.

Certain specific embodiments of the process of multiplex PCR may be covered by patents of third parties in certain countries and may require a license.

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© 2011 QIAGEN, all rights reserved.