



March 2025

Quick-Start Protocol

QIAamp[®] DNA Host-Free Microbiome Kit

The QIAamp DNA Host-Free Microbiome Kit is shipped at room temperature (15–25°C). Upon receipt, store the Host Depletion Solution and the QIAamp UCP Mini Spin Columns at 2–8°C, Benzonase at –20°C, and all other kit components dry at room temperature.

Further information

- Safety Data Sheets: www.qiagen.com/safety
- Technical assistance: support.qiagen.com

Notes before starting

- Dissolve lyophilized Host Depletion Solution in 1.8 mL PBS. Mix gently by inverting. Do not spin down. For long-term storage, store single-use aliquots at –15°C to –25°C. Thawed aliquots can be stored at 2–8°C for 6 weeks. Note that precipitate may form during storage but do not compromise the activity of the product.
- Tissue samples: mince or grind sample before starting the protocol (e.g., using TissueRuptor II, cat. no. 990890).
- Swab samples: swirl the swab in 1 mL transport media or PBS for at least 20 s and dry off by pressing against the wall of the tube multiple times. Spin down (10 min, 10,000 × g) and remove the supernatant.

- If a precipitate has formed in Buffer ATL or Buffer APL2, dissolve by incubating at 56°C.
- Ensure that Buffers AW1 and AW2 have been prepared according to the instructions on the bottle.
- Add 100 µL Reagent DX to 15 mL Buffer ATL. If smaller amounts are needed, transfer 1.5 mL of Buffer ATL into a sterile 2 mL vial and add 10 µL Reagent DX. Mix well. After preparation, the mixture is stable for 6 months at room temperature (15–25°C).
- Spin the PowerBead Pro Tube briefly to ensure that the beads have settled at the bottom.

Equipment and reagents to be supplied by user

- Ethanol (96–100 %)
- Phosphate-buffered saline (PBS)
- Microcentrifuge (with rotor for 2 mL tubes)
- Shaker–incubator
- Equipment for sample disruption and homogenization, one of the following:
 - Vortex Genie 2 and Vortex Adapter for 24 (1.5–2 mL) tubes (cat. no. 13000-V1-24)
 - TissueLyser III (cat. no. 9003240) with adapter sets for use with the PowerBead Pro Tubes (TissueLyser Adapter Set 2 x 24, cat. no. 69982), or 2 mL Tube Holder (cat. no. 11993), in conjunction with Plate Adapter Set, cat. no. 11990)

Procedure

1. Add 220 µL of Buffer RDD, 3 µL Benzonase and 30 µL Host Depletion Solution to up to 100 mg of tissue or the pre-processed swab sample in a 2 mL microcentrifuge tube (in that order). Mix well and incubate at 37°C for 45 min at 600 rpm in a heating block or water bath.

2. Add 20 μL Proteinase K, mix well, and incubate at 56°C for 20 min at 600 rpm in a heating block or water bath.
3. Add 200 μL Buffer ATL (containing Reagent DX). Mix well to avoid loss of sample material, and transfer into a PowerBead Pro Tube.
4. Mechanically disrupt samples using one of the following methods:
 - a. Secure the PowerBead Pro Tube horizontally on a Vortex Adapter for 24 (1.5–2 mL) tubes (cat. no. 13000-V1-24). Orient the tube caps to point toward the center of the vortex adapter. Vortex at maximum speed for 10 min.
 - b. Use TissueLyser III. Place the PowerBead Pro Tube into the TissueLyser Adapter Set 2 x 24 (cat. no. 69982) or 2 mL Tube Holder (cat. no. 11993) and Plate Adapter Set (cat. no. 1190). Fasten the adapter into the instrument and shake for 5 min at 30 Hz. Reorient the adapter so that the side that was closest to the machine body is now furthest from it. Shake again for 5 min at 30 Hz.
5. Centrifuge the PowerBead Pro Tube at $10,000 \times g$ for 1 min. Transfer the supernatant to a fresh microcentrifuge tube.
6. Add 200 μL APL2. Mix by pulse vortexing for 30 s.
7. Incubate at 70°C for 10 min and briefly spin the tube to remove condensation.
8. Carefully apply up to 700 μL of the mixture from step 7 to the QIAamp UCP Mini Column. Close the cap and centrifuge at $6000 \times g$ for 1 min.
9. Discard the flow-through. Put the column back into the collection tube to repeat step 8 with any remaining mixture from step 7.
10. Transfer the QIAamp UCP Mini Column to a fresh collection tube. Carefully open the cap and add 500 μL Buffer AW1 without wetting the rim. Close the cap and centrifuge at $6000 \times g$ for 1 min. Discard the filtrate.

11. Add 500 μ L Buffer AW2 to the QIAamp UCP Mini Column without wetting the rim. Centrifuge at full speed (20,000 \times g) for 3 min. Discard the filtrate.
12. Centrifuge at full speed (20,000 \times g) for 1 min.
13. Place the QIAamp UCP Mini Column into a fresh 1.5 mL tube and apply 50 μ L RNase-Free Water directly onto the center of the membrane. Close the lid and incubate at room temperature for 5 min.
14. Centrifuge at 6000 \times g for 1 min to elute the DNA.

Document Revision History

Date	Description
03/2025	Initial release

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN®, Sample to Insight®, QIAamp®. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

03/2025 HB-3652-001 1136340 © 2024 QIAGEN, all rights reserved.

Ordering www.qiagen.com/shop | Technical Support support.qiagen.com | Website www.qiagen.com